PREFACE

Thank you for choosing FitTrack. Please read our user manual carefully before using **Atria 2.0**.

This manual will introduce you to the method, functions and operation of **Atria 2.0**. This will help familiarize yourself with the product.

Disclaimer. All copyrights in this manual are owned by FitTrack Inc. No part of this manual may be used, reproduced or trans-

mitted in any form or by any means for any other purpose without the express written permission of the Directors. This manual is based on available information and its contents are subject to change without notice. We have tried our best to make the contents accurate and reliable at the time of writing this manual, but we are not responsible for any loss or damage caused by omissions, inaccuracies or misalignment in this manual.

PRODUCT DESCRIPTION

Atria 2.0 is a fitness tracker that integrates heart rate, blood oxygen monitoring, exercise data collection, sleep monitoring and body temperature. Users can monitor their body condition and cultivate a healthy lifestyle. **Atria 2.0** can be combined with the **FitTrack Health** app to personalize and synchronize data. In addition, the watch has functions such as push notifications.

INSTRUCTIONS

BASIC FUNCTIONS

Power on: When the watch is off, press and hold the power key for 8 seconds to turn on, and the main page of the watch will display after turning on.

Touch screen operation: Swipe to the left to the menu page, then tap the icons to enter the corresponding sub-interface.

Bright screen: When the watch screen has timed-out...

Option 1: Press the power key on the right side of the watch to relight the screen.

Option 2: In the FitTrack Health App, turn on “Turn Wrist Light” feature. When you turn your wrist towards you, the screen will relight.

Shutdown: Long press the power key on the side, then click the “Shutdown” icon to shut down.

BLUETOOTH CONNECTION/DISCONNECTION

Follow the steps below to connect to the App. If the following steps are not followed, **Atria 2.0** will not function properly.

(1) Connect Atria 2.0 to FitTrack app.

- Step 1: Search for "FitTrack Health" in the App Marketplace (Android) or App Store (iOS)
- Step 2: Download the latest version of FitTrack Health
- Step 3: Ensure Bluetooth on your phone is turned on

- Step 4:
 - + For new users: Sign Up using your email address and go through onboarding process.
 - + For returning users: Log in using your e-mail address address. If you have an existing FitTrack watch or tracker connected, please remove it from the device list in the app and from your phone's bluetooth list. Note: only one watch can be connected to one account at a time.
- Step 5: From the home page, click + on the upper right corner

to enter into device pairing screen

+ Click Search for Devices.

+ Select your **Atria 2.0** device and Connect device to bind watch.

(2) Disconnect Atria 2.0 from FitTrack Health app:

- Step 1: Open **FitTrack Health app**

- Step 2: From the home page, click + on the upper right corner to enter into device pairing screen.

- Step 3: Swipe left on the device that you wish to remove or

- click Remove All Device to remove all existing devices.
- *Step 4 (iOS Users only) If you are using an iOS system, please remove the device from your phone's bluetooth settings after removing the device from the app.*

FEATURES INTRODUCTION AND OPERATING INSTRUCTIONS

SPORT MODE

How to operate Sports Mode:

- From the home screen, swipe left once to enter into the

Sports Mode list.

- Select the Sport by clicking on the sports icon.

Note: Some modes are able to customize targets and reminders.

- Click Start when ready to begin.

- **Atria 2.0** will then display the movement time, calories, steps etc...

While in sports mode: Swipe up and down to view more details.

Pause: Sports mode will automatically be paused when swiping right once. The watch will display **Pause Time**. Total time will

continue to calculate, however the workout time will be paused.

Exit: Swipe the screen to the right on the data display interface and click the stop icon to exit the mode.

*Sync **Atria 2.0** to FitTrack Health to review more details*

STEP/DISTANCE/CALORIE

(1) Functions: The watch will display the user's total number of steps on the day, walking distance, calories burned on the

same day, and the data will be cleared at 00:00 (12AM) o'clock every day local time.

BLOOD OXYGEN MONITORING

(1) Function introduction: *Atria 2.0* will measure the user's blood oxygen under this interface.

HEART RATE MONITORING

Function introduction: *Atria 2.0* will measure your heart rate under the heart rate measurement interface.

BODY TEMPERATURE MONITORING

Function introduction: *Atria 2.0* will measure your body temperature under the Body Temp. interface.

SCIENTIFIC SLEEP

Function introduction: *Atria 2.0* will monitor the sleep cycle for the night, sync to the FitTrack Health app to view more details from the night before. Sleep measurement time is ONLY between 20:00 – 9:00 next day. If you fall asleep anytime during this timeframe, *Atria 2.0* will monitor your sleep.

STATEMENT

It is strictly forbidden to copy, transmit, distribute and store any part of this manual in any form without our prior written permission. The correctness and completeness of the content have been requested during the preparation of this manual. The single statement does not guarantee that there are no errors or omissions in this manual. The company will not be responsible for any errors or omissions that may occur in this

manual. As there are continuous new improvements, we reserve the right to modify and improve any of the products described in this manual without prior notice. Under no circumstances shall we be liable for any loss of data and income or for any special, incidental, incidental or consequential damages, regardless of the cause of the loss.

Need Help? For more product and troubleshooting information, visit <https://getfittrack.com/pages/guide>

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FCC Warning Statement

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- - Reorient or relocate the receiving antenna.
- - Increase the separation between the equipment and receiver.
- - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- - Consult the dealer or an experienced radio/TV technician for help.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

The device has been evaluated to meet general RF exposure requirement.