EARMOR[®]

**XDZ-H20T In-ear Bluetooth Audiphone
User manual**

XDZ-H20T

In-ear Bluetooth Audiphone

**This product complies
with user compatibility standards**

INTRODUCTION

EARMOR[®] XDZ-H20T Bluetooth Audiphone use in-ear hearing aids design, which can obviously Amplifies ambient sound, allowing users to improve the perception of the outside world. With the Bluetooth connection function, you can connect devices in a noisecancelling state, freely listen to music, answer calls or reject calls, enhancing the diversity of communication and entertainment. The headset can last 8 hours when fully charged, and has a voice prompt function for instant confirmation of low battery or any operating situation. Compared with all kinds of traditional sound pickup and noise reduction headphones, the product is smaller in size, lighter in weight, and integrates protective gear compatibility, portability and comfort.

1.PRODUCT


Item: In-ear Bluetooth Audiphone

2.ACCESSORIES INCLUDED IN THE PACKAGE

XDZ-H20T In-ear Bluetooth Audiphone Package includes earbuds, charging case, battery and main unit

NAME OF MAIN PARTS (IN-EAR BLUETOOTH AUDIPHONE)


1. Power/switch mode button
2. Pickup microphone
3. Charging copper pin
4. Earplugs
5. Ear clip
6. Indicator light

NAME OF THE RECHARGING BOX PARTS

1. Anti-lost rope fixing post
2. Power button
3. Charging copper pin
4. Battery indicator
5. Hook hole
6. Audiphone Indicator (The red light is the charging state; the 6.red and blue lights flash alternately for the power-on state)

*The fixed columns on both sides of the charging box can be rotated 360°, and the anti-lost rope can be wrapped and stored here when not in use.

3. WORKING PRINCIPLE

The sound sensor converts the sound information into an electrical signal, which is processed by the analog chip and transmitted to the receiver through the sound guide tube for output. The signal processing adopts the principle of analog signal amplification.

4.SCOPE OF APPLICATION

After testing and inspection, it can be used to compensate hearing loss for patients with air conduction hearing loss.

Model	XDZ-H20T
Maximum OSPL90	$\leq 115\text{dB}+3\text{dB}$
high frequency average OSPL90	$\leq 105\text{dB}\pm 4\text{dB}$
full sound gain	$15\pm 5\text{dB}$
Equivalent Input Noise Level	$\leq 29+3\text{dB}$
total harmonic distortion	$\leq 3\%$
Frequency response range (not narrower than)	200Hz-5000Hz
Rated power supply current consumption	$\leq 14\text{mA}$
Rated voltage	d.c.3.7V

Model specification parameter table

Note: The above parameter index inspection equipment is AWA6128AH Hearing aid tester
Standard: IEC60118-7 2005; by using 2CC coupled cavity HA-1

Model	Recharging box
Working current	$\leq 140\text{mA}$
d.c.3.7V	d.c.3.7V
Charging box input voltage	d.c.5V
Hearing aid charging current	$\leq 140\text{mA}$
Audiphone charging voltage	d.c.5V

Recharging box model parameter table

5. PRODUCT PERFORMANCE

1. Appearance and Structure

- 1.1 All parts should be fully assembled and fixed reliably.
- 1.2 Text, symbols or marks should be clear, correct and firm.
- 1.3 The surface should be free of burrs; burrs, dents, scratches and other defects.
- 1.4 The connection of each component should be reliable and effective.
- 1.5 The operation of each adjustment button, switch, and other mechanisms should be flexible and effective.

2. Electroacoustic performance

Meet the requirements of the model specification parameter table

3. Electroacoustic performance after environmental test

Meet the requirements of the model specification parameter table

4. Electromagnetic Compatibility

- a) Should meet GB/T 25102.13-2010 "Electroacoustic Hearing Aids Part 13: Electromagnetic Compatibility (EMC)" requirements of capacitive regulations
- b) Should meet YY 0505-2012 "Medical Electrical Equipment Part 1-2: General Requirements for Safety Collateral Standard: Electromagnetic Compatibility Required Summation Test".

5. Environmental test requirements

The environmental test of hearing aids shall comply with the provisions in 4.4 and Tables 2.1 and 2.2 of GB/T 14199-2010. The transport test shall comply with The provisions of Chapter 4 in GB/T 14710-2009.

6. safety requirements

It should comply with the provisions of GB9706.1-2007 "Medical Electrical Equipment Part 1: General Requirements for Safety".

6. CONTRAINDICATIONS

- 1) Congenital atresia or deformity of the external auditory canal, a large number of foreign bodies in the external auditory canal;
- 2) Deafness status is uncertain, and the degree of hearing loss often fluctuates;
- 3) The main complaints are those with obvious symptoms such as headache, dizziness, earache, tinnitus, etc.;
- 4) Sudden deafness in the past 3 months; progressive rapid deafness; obvious unilateral deafness or other acute ear diseases
- 5) Middle ear fluid and often pus;
- 6) Central deafness; non-organic deafness;
- 7) Patients with acute otitis externa and tympanitis;
- 8) Patients with acute and chronic suppurative otitis media (in the period of purulent infection);
- 9) Those who are allergic to the material of this product

7. PRECAUTIONS, WARNINGS AND TIPS

1. Before fitting a audiphone, it should undergo a professional inspection and hearing test, and use it under the guidance of a doctor or a professional audiphone dispenser.

2. Pay attention to the moisture-proof, fog-proof and high-temperature proof of the product.

3. Do not let the audiphone immerse in any liquid. Remove the audiphone before swimming or washing your hair and bathing; after swimming or washing, wait

Put on audiphone after the ear canal is dry.

4. After removing the audiphone at night, it should be placed in a box with desiccant, protected from light, high temperature, high temperature environment, and avoid touching.

5. The outer surface of the audiphone should be kept clean, and the earwax in the earplug should be cleaned frequently.

6. The function and volume control switch should be used correctly, try to turn on and off lightly, and avoid excessive force and excessive force.

7. Audiphone should be regularly sent to the fitting center for maintenance and inspection to ensure that the audiphone are in good condition.

8. If you feel uncomfortable in your ears during use, please stop using it immediately.

9. Please do not put the audiphone where the baby can reach. If swallowed, please seek medical attention immediately.

10. Applicable objects: It is suitable for people with air conduction hearing loss of about 30-50dB.

11. Please use the audiphone with the audiphone charging box.

12. When not using audiphone, please charge the product regularly to avoid battery damage; the recommended charging interval for audiphone is ≤ 28 days

The recommended charging interval for the audiphone charging case is ≤ 90 days.

8. INSTALLATION AND USE INSTRUCTIONS AND DIAGRAMS

a. Installation instructions for ear buckles, ear plugs, and anti-lost cord


Use only ear clips and ear tips designed for the XDZ-H20T. Do not insert the earphones directly into your ears without the XDZ-H20T's ear clips and earbuds installed.

1. Install the ear clip of the appropriate size on the earphone according to the corresponding direction (Note: When installing, make sure that the fixing hole in the ear clip The ear buckle fixing post of the machine fits perfectly)

2. To install the Anti-Lost Cord, attach one end of the Anti-Lost Cord to the earbud stem

3. Install the earbuds on the earbud stems

4. Make sure the earbuds are fully seated on the stems






2. Instructions for wearing the headset

Turn off the power before wearing the headset. Earbuds should be clean and in good condition prior to installation.

1. Choose the right ear tip size for your ear canal
2. Hold the XDZ-H20T with the earbuds attached
3. Insert the earbud into the ear canal while gently pulling up on the ear with the other hand
4. Rotate the headset towards the inside of the ear canal to ensure that the ear clip is fully seated in the ear socket
5. With the XDZ-H20T turned off, check to make sure the earbuds are the right size, and gently pull the earbuds to make sure they won't easily come out of your ears.

If it comes out easily, please repeat the above steps to wear.

3. Headphone Operation Instructions

1. Power on from the warehouse: the headset will automatically turn on after taking it out

Manual power on: When the earphone is turned off to restart, you need to press and hold the button on the left and right ears for 3 seconds to power on each, and when you hear the prompt tone

Release the button when "power on", the headset is turned on and in outdoor mode.

2. Volume switch: click any button left or right to adjust the mode, the order is: large -- small -- mute -- large -- small.

3. Device pairing: Take out the headset and automatically turn it on, then enter the pairing mode (white and blue lights flash alternately), the pairing name

"XDZ-H20T", after the device selects the headset and the pairing is successful, the LED stops flashing, and the prompt sound "Connected" appears.

4. Listening to music: After the device is successfully connected, double-click any button left or right to play/pause the music.

5. Transparency mode: During music playback, you can click any button on the left or right to adjust the mode.

6. Call mode: click any button on the left or right to answer the call, double-click any button to hang up the call, and long press any button when an incoming call

key for 2 seconds to reject the call.

7. Shutdown when entering the warehouse: The headset will automatically shut down after being placed in the charging compartment

Manual shutdown: Press and hold the unilateral button for 3 seconds to synchronously shutdown, and a prompt sound "power off" will appear.

Automatic shutdown: no operation/use, it will automatically shutdown after two hours, and the prompt sound "power off" will appear.

FAULT SELF-CHECK

Items	Possible Causes	Solutions
No sound from audiphone	dead battery	recharge
	The sound outlet or microphone is blocked	Clean the sound outlet or microphone port
The sound is weak and blurry	low battery	recharge
	Volume setting is too low	Turn up the volume
	loose earbuds	Re-tighten audiphone

4. Instructions for charging operation


Place the audiphone into the charging box, use the Type-C data cable to connect the charging box to the power supply, and automatically charge when the circuit is connected.

The battery indicator also flashes.

If the battery indicator does not flicker, please check whether the charging box, data cable, and power supply are properly connected.

When charging is complete, the battery indicator light is always on.

NOTE: Make sure the charging contacts of the headset point down to make contact with the charging point in the socket.

Items	possible reason	Solutions
After the audiphone is connected, indicator light is off	Poor contact between the audiphone and the charging case	Reinsert the audiphone into the charging case
	Audiphone reversed	Check if audiphone are reversed
Light does not shine	Power is not plugged in	Connect the USB cable to charge

WARN

Do not disassemble, puncture, crush, heat or throw batteries into fire

Do not use after severe swelling or immersion in water

9. MAINTENANCE AND MAINTENANCE

- Earplug cleaning and disinfection method: Take off the earplugs from the audiphone, rinse with water, or wipe with alcohol, and wait for the earplugs to dry
After that, you can install it on the audiphone.
- Cleaning method of the product: It can be wiped with a soft lint-free cloth or a sponge soaked with detergent. Do not pour liquids on the instrument during cleaning
On the instrument, make sure that no liquid can enter the inside of the instrument. Be especially careful when wiping to make sure no water can get in.

10. STORAGE/TRANSPORTATION CONDITIONS AND METHODS

- the normal use of environmental conditions
 - Ambient temperature 0°C-40°C;
 - Relative humidity: ≤ 80%
 - Atmospheric pressure: 86kPa-106kPa;
- The product should be stored in a cool and dry place, and the environment should meet the following conditions:
 - Ambient temperature: -10°C -55°C
 - Relative temperature: ≤ 90%
 - Atmospheric pressure: 50kPa-150kPa;
- Handle with care to avoid collision or strong impact on the product.

11. DISPOSAL OF WASTE AND RESIDUES

Please do not discard discarded earbuds, mainframes, and charging boxes at will. The disposal of wastes and residues should comply with the corresponding national laws and regulations.

12. PARTS LIST

- A pair of XDZ-H20T in-ear Bluetooth audiphone (with a pair of medium ear clips and a pair of medium ear tips)
- Portable lithium battery charging case *1
- Silicone anti-lost rope *1
- Charging cable *1
- A pair of small ear clips & earplugs
- A pair of large ear clips & earplugs
- Manual *1
- conformity label

ATTACHMENT: ELECTROMAGNETIC COMPATIBILITY INFORMATION

	Note, refer to accompanying documents		manufacturer
	Type B applied part equipment		up
	CE marked		protected from rain
	Production batch number		stacking layers
	Separate disposal of waste electrical and electronic equipment sign (please observe local laws and regulations)		

NOTICE:

1. In-the-ear audiphone shall comply with the requirements of YY0505-2012 and GB/T25102.13-2010 for neighbor compatibility; Users should install and use according to the electromagnetic compatibility information provided in the random document;
2. Portable and mobile RF communication equipment may affect the performance of in-ear audiphone, avoid strong electromagnetic interference when using, such as close to mobile phones, microwave oven, etc.;
3. The guidelines and manufacturer's declaration are detailed in the appendix.

WARNING:

1. The in-ear audiphone should not be used in close stacking with other devices. If it must be used in close stacking, it should be verified by
It can run normally under the configuration used;
2. With the exception of cables sold by manufacturers of in-the-ear audiphone as spare parts for internal components, the use of unspecified accessories and cables may
Causes an increase in emissions from in-ear audiphone or a decrease in noise immunity.

Attachment: Electromagnetic Compatibility Information

In-the-ear audiphone are intended to be used in the electromagnetic environment specified below.

launch test	compliance	Electromagnetic Environment - Guidelines
GB4824 RF emissions	Group 1	In-the-ear audiphone use RF energy only for their internal function, therefore, its RF emissions are low and may not cause any interference to nearby equipment In-the-ear hearing aids are suitable for use in the home and in all installations directly connected to the residential public low-voltage supply network for domestic use.
GB4824 RF emissions	Class B	
Gb17625.1 Harmonic emission	Not applicable	
GB17625.2 voltage fluctuation/ flicker launch	Not applicable	

Guidelines and Manufacturer's Declaration - Electromagnetic Immunity

In-the-ear audiphone are intended to be used in the electromagnetic environment specified below, and the purchaser or user of the in-the-ear audiphone should ensure that it is used in this electromagnetic environment:

Immunity test	IEC60601 test level	match level	Electromagnetic Environment - Guidelines
Electrostatic discharge (ESD) GB/T17626.2	±6kV contact discharge ±1kV air discharge	±6kV contact discharge ±6kV air discharge	Floors should be wood, concrete or tile, and if floors are covered with synthetic materials, the relative humidity should be at least 30%
electrical fast transient burst GB/T 17626.4	±2kV to power line ±1kV for input/output lines	Not applicable	Not applicable
surge GB/T17626.5	±1kV differential mode voltage ±2kV	Not applicable	Not applicable
Voltage dips, short interruptions, and voltage variations on power input lines GB/17626.11	<5% Ut for 0.5 weeks (>95% dip on Ut) 40% Ut for 5 weeks (60% dip on Ut) 70% Ut for 25 weeks (on Ut, 30 % dip) < 5% Ut for 5s (> 95% dip on Ut)	Not applicable	Not applicable
Power frequency magnetic field (50/60HZ) GB/T 17626.8	3A/m	3A/m	The power frequency magnetic field shall have the power frequency magnetic field level of a typical location in a typical commercial or hospital environment

Note: Ut refers to the AC mains voltage before applying the test voltage

Guidelines and Manufacturer's Declaration - Electromagnetic Immunity

In-the-ear audiphone are intended to be used in the electromagnetic environment specified below, and the purchaser or user of the in-the-ear audiphone should ensure that it is used in this electromagnetic environment:

Immunity test	IEC60601 test level	match level	Electromagnetic Environment - Guidelines
RF conduction GB/T17626.6	3 Vrms 150 kHz to 80MHz	Not applicable 3V/m	Portable and mobile RF communications equipment should not be used closer to any part of the in-the-ear audiphone, including cables, than the recommended isolation distance. The distance should be determined by the public frequency corresponding to the frequency of the transmitter formula calculation. Recommended isolation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80MHz to 800MHz $d = 2.3 \sqrt{P}$ 800MHz to 2.5GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters are determined by surveying the electromagnetic site and should be lower than the compliance level in each frequency range. Interference may occur near equipment marked with the following symbols:
RF radiation GB/T17626.3	3V/m 80 MHz to 2.5 GHz	3V/m	

Note 1: At 80MHz and 800MHz, the formula for the higher frequency band applies.

Note 2: These guidelines may not be suitable for all situations, electromagnetic propagation is affected by absorption and reflection from buildings, objects and people.

a Field strengths of fixed transmitters, such as: base stations for wireless (cellular/cordless) telephones and land mobile radios, amateur radio, AM (amplitude modulation) and FM (frequency modulation) radio broadcasts, and television broadcasts, are not theoretically capable of Accurate prediction. To assess the electromagnetic environment of fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the in-ear audiphone is located is higher than the RF compliance level for the above application, the in-the-ear audiphone should be observed to verify proper operation. If abnormal performance is observed, supplementary measures may be necessary, such as reorienting or repositioning the intraocular audiphone. b In the whole frequency range of 150KHz-80MHz, the field strength should be lower than 3V/m.

Recommended separation distances between portable and mobile RF communication equipment and in-the-ear audiphone					
In-the-ear audiphone are intended for use in electromagnetic environments where radiated RF disturbances are controlled. Buyers and users of in-the-ear audiphone can prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and in-the-ear audiphone as recommended below, based on the maximum output power of the communication device.					
The rated maximum output power of the transmitter/w.	Corresponding to the isolation distance of different frequencies of the transmitter /m				
	150kHz-80MHz $d=1.2\sqrt{P}$	80MHz-800MHz $d=1.2\sqrt{P}$	800MHz-2.5GHz $d=2.3\sqrt{P}$		
0.01	/	0.12	0.23		
0.1	/	0.38	0.73		
1	/	1.2	2.3		
10	/	3.8	7.3		
100	/	12	23		
determined by the formula in the Transmitter Frequency column, where P is the transmitter's maximum power rating in watts (w) as provided by the transmitter manufacturer Note: At 80MHz and 800MHz, the formula for the higher frequency range is used. NOTE: These guidelines may not apply in all situations, electromagnetic propagation is affected by absorption and reflection from buildings, objects and people.					
The field strength of the RF test signal when the audiphone complies with the audiphone immunity specified in GB/T25102.13.					
	Neighbor Compatibility at IRIL \leq 55Db Frequency Range GHz Field Strength in V/m				
Frequency range GHz	0.08-0.8	0.8-0.96	0.96-1.4	1.4-2.0	2.0-3.0
Microphone mode	/	3	/	2	/

This manual is only suitable for rechargeable battery models

Electromagnetic compatibility: should meet the requirements specified in GB/T25102.13-2010.

[Product name] Bluetooth in-ear audiphone

[Product Model] XDZ-H20T

[Working power supply]

[Registration certificate number]

[Product technical requirement number]

[Production license number]

[Registrant / Manufacturer / After-sales / Brand Operation]

Guangzhou Actions Technology Co., Ltd.

Address : Room 601, Block A, No. 94, Liwan Road, Liwan District, Guangzhou City, Guangdong Province, China

Tel : +86 20 81179170 Fax : +86 20 81179171

Email : info@opsmen.com/cs@opsmen.com

MAINTENANCE SERVICE

The EARMOR XDZ-H20T Electronic Hearing Protection Earplug Kit can be repaired or replaced free of charge within one year from the date of purchase

Repairs and shipping charges for component failure are excluded.

This repair service does not apply to the following situations:

1. Vandalism, disassembly, modification of this product
2. Product damage caused by wrong operation

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

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