Assert-IQ™

Model DM5000 Model DM5300 Model DM5500 Insertable Cardiac Monitor

User's Manual



CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

WARNING: This product can expose you to chemicals including ethylene oxide, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

[™] Indicates a trademark of the Abbott group of companies.

‡ Indicates a third-party trademark, which is property of its respective owner.

Bluetooth and Bluetooth logo are registered trademarks of Bluetooth SIG, Inc.

Pat. http://www.abbott.com/patents

© 2022 Abbott. All Rights Reserved.

Description

This manual describes the following Abbott Medical device:

Table 1. Assert-IQ™ insertable cardiac monitor

Name	Model Number	Description	MRI Status
Assert-IQ ICM 3	DM5000	Insertable cardiac	MR Conditional
Assert-IQ ICM3 +	DM5300	monitor	
Assert-IQ ICM3 EL+	DM5500		

Table 2. Assert-IQ accessories

Name	Model Number	Description
Insertion tool	DM5310 (for device models DM5000 and DM5300) DM5510 (for device model DM5500)	Tool to insert the Assert-IQ ICM
Incision tool	DM5320	Tool to cut the skin during insertion procedure

The Abbott Medical Assert-IQ insertable cardiac monitor (ICM) is designed to detect arrhythmias and wirelessly transmit data to the Merlin.net™ Patient Care Network (PCN).

The ICM constitutes the inserted portion of the system. The Merlin™ Patient Care System (PCS) with software version 25.3 (or greater), magnet, myMerlin™ mobile application (app), and Merlin.net PCN constitute the external portion of the system.

The Merlin PCS and magnet are used to interrogate and program the device in the clinic. Remote transmissions are performed using the app. The app also allows patients to record and send EGMs of symptomatic events to the clinic. All remotely transmitted data is made available on Merlin.net where clinicians can log in, review data, make a diagnosis and edit the monitoring parameters.

Indications for Use

The Assert-IQ™ ICM is indicated for the monitoring and diagnostic evaluation of patients who experience unexplained symptoms that may be cardiac-related such as: dizziness, palpitations, chest pain, syncope, and shortness of breath, as well as patients who are at risk for cardiac arrhythmias such as bradycardia, tachycardia, and sinus pauses. The Assert-IQ ICM is also indicated for patients who have been previously diagnosed with atrial fibrillation (AF) or who are susceptible to developing AF. The Assert-IQ ICM is intended to be inserted subcutaneously in the left pectoral region, also described as the left anterior chest wall.

Intended Use

The Assert-IQ™ ICM is intended to help physicians and clinicians monitor, diagnose and document the heart rhythm in patients who are susceptible to cardiac arrhythmias and unexplained symptoms by detecting arrhythmias and transmitting data for review.

Table 3. Accessories and their intended uses

Accessory	Model Number	Intended purpose
Incision tool	DM5320	The Assert-IQ™ ICM incision tool is intended to puncture the skin, as a first step during the insertion procedure of the Assert-IQ ICM.
Insertion tool	DM5310 (for device models DM5000 and DM5300) DM5510 (for device model DM5500)	The Assert-IQ™ ICM insertion tool is intended to create a subcutaneous pocket and deliver the Assert-IQ ICM through the skin incision into the subcutaneous pocket.

MRI Safety Information

Testing has demonstrated that the Assert-IQ ICM is conditionally safe for use in the MRI environment when used according to instructions in the MRI-Ready Monitor Systems manual.

CAUTION: Do not bring any external control devices, such as the Merlin™ Patient Care System (PCS) Model 3650 or a mobile device, into the scanner magnet room (Zone IV). These devices are considered MR Unsafe.

Contraindications

There are no known contraindications for the insertion of the Assert-IQ™ ICM. However, the patient's particular medical condition may dictate whether or not a subcutaneous, chronically inserted device can be tolerated.

Warnings and Precautions

Sterilization

- The device and the incision and insertion tools have been sterilized with ethylene oxide prior to shipment. They are intended for single use only and should not be resterilized.
- If the sterile package has been compromised, contact Abbott Medical.
- Do not insert the device if the dot on the ethylene oxide label is purple. Purple indicates that the
 package has not been sterilized. Return the device to Abbott Medical.

Package Inspection

- Check the "use-before" date on the package label. Do not insert the device if its "use-before" date has expired.
- Do not use the device if the packaging is wet, punctured, opened or damaged because the integrity of the sterile packaging may be compromised. Return the device to Abbott Medical.

Storage and Handling

- Store the device between 59°F and 86°F (15°C and 30°C).
- During transportation and handling, the device can be exposed to temperature excursions between -4°F and 140°F (-20°C and 60°C).
- Exposure to temperatures outside this range may result in damage to the device or device malfunction.

Temperature Equilibrium

 After cold storage, allow the device to reach room temperature before programming or inserting the device because cold temperature may affect initial device function.

Insertion

 Insert the device in the subcutaneous space, just under the skin. Insert the device no deeper than 0.80 in to ensure reliable data transmission.

Device Replacement

Replace the device within one month of receiving a low battery alert, if necessary or desired. Replace
the device immediately upon receiving a low battery alert if frequent EGMs are being stored and
remotely transmitted.

Explant and Disposal

- Interrogate the device to acquire stored EGMs and episode data and turn monitoring off before
 explanting, cleaning or shipping the device to prevent unwanted EGM and episode storage.
- Explant the device with standard surgical tools upon receiving an End-of-Service (EOS) alert.
- Locate the device by palpating it or by using fluoroscopy and mark its location (directly above the
 implanted device or the header side of the device). Following aseptic techniques, make an incision
 either on the implant scar or at the header end of the ICM wide enough to accommodate the ICM.
- Return all explanted devices to Abbott Medical.

• Never incinerate the device because of the potential for explosion. Explant the device before cremation.

Environmental and Medical Therapy Hazards

Instruct patients to avoid devices that generate a strong electric or magnetic interference (EMI). EMI
could cause device malfunction or damage, resulting in inappropriate episode storage or inhibition of
episode storage. Moving away from the source or turning it off will usually allow the device to return to
its normal mode of operation.

Hospital and Medical Environments

- Electrosurgical cautery may cause device malfunction or damage. If electrocautery is necessary, keep
 the current path and groundplate as far away from the device as possible.
- External defibrillation may damage the device. Minimize current flowing through the device by following these precautions when using external defibrillation on a patient with a device:
 - Position defibrillation paddles as far from the device as possible (minimum of 2.4 in)
 - Use the lowest clinically appropriate energy output
 - Confirm the device function following any external defibrillation
- Do not direct high radiation sources such as cobalt 60 or gamma radiation at the device. If a patient
 requires radiation therapy in the vicinity of the device, place lead shielding over the device to prevent
 radiation damage and confirm its function after treatment.

- Lithotripsy may permanently damage the device. Avoid it unless the therapy site is not near the device.
- Avoid diathermy, even if the device is programmed off, as it may damage tissue around the device or may permanently damage the device.
- The device should not be exposed to therapeutic levels of ultrasound energy, as the device can inadvertently concentrate the ultrasound field and cause harm that might not be immediately detectable. Diagnostic ultrasound treatment is not known to affect the function of the device.
- Transcutaneous Electrical Nerve Stimulation (TENS) may interfere with device function. To reduce interference, place the TENS electrodes close to one another and as far from the device as possible. Monitor cardiac activity during TENS use.
- Radiofrequency (RF) ablation in a patient with a device may cause device malfunction or damage.
 Minimize RF ablation risks by:
 - Disabling monitoring
 - Avoiding direct contact between the ablation catheter and the inserted device
 - Positioning the groundplate so that the current pathway does not pass near the inserted device, i.e., place the groundplate under the patient's buttocks or legs

Home and Occupational Environments

 High-voltage power transmission lines may generate enough EMI to interfere with device operation if approached too closely.

- Communication equipment such as microwave transmitters or high-power amateur transmitters may generate enough EMI to interfere with device operation if approached too closely.
- Home appliances in good working order and properly grounded do not usually produce enough EMI to interfere with device operation. There are reports of device disturbances caused by electric hand tools or electric razors used directly over the device insertion site.
- Wireless communication devices such as computers that operate on a wireless network, cellular phones, smart phones, tablets, and even cordless telephones may generate enough EMI to interfere with device operation.
- A variety of industrial equipment produce EMI of sufficient field strength and modulation characteristics
 to interfere with proper operation of the device. These include, but are not limited to: arc welders;
 induction furnaces; very large or defective electric motors; and internal combustion engines with poorly
 shielded ignition systems.

Electronic Article Surveillance (EAS)

Advise patients that the Electronic Article Surveillance/Anti-theft (EAS) systems such as those at the point of sale and entrances/exits of stores, libraries, banks, etc., emit signals that may interact with the device. It is very unlikely that these systems will interact with their device significantly. However, to minimize the possibility of interaction, advise patients to simply walk through these areas at a normal pace and avoid lingering near or leaning on these systems.

Metal Detectors

• Advise patients that metal detector security systems such as those found in airports and government buildings emit signals that may interact with their device. It is very unlikely that these systems will interact with their device significantly. To minimize the possibility of interaction, advise patients to simply walk through these areas at a normal pace and avoid lingering. Even so, the device contains metal that may set off the airport security system alarm. If the alarm does sound, the patient should present security personnel with their patient identification card. If security personnel perform a search with a handheld wand, the patient should ask that they perform the search quickly, stressing that they should avoid holding the wand over the device for a prolonged period.

Mobile Devices

- The device has been tested for compatibility with handheld wireless transmitters in accordance with the requirements of ISO 14117:2019. This testing covered the operating frequencies (385 MHz 3 GHz) and pulsed modulation techniques of all of the digital cellular phone technologies in worldwide use today. Based on the results of this testing, the device should not be affected by the normal operation of cellular phones when used more than 6 in from the device.
- To minimize the possibility of interaction, advise patients not to carry a cellular phone in a breast pocket or on a belt within 6 in of the device, and to use a cellular phone on the side of their body opposite from the device.

Potential Adverse Events

Possible adverse events (in alphabetical order) associated with the device, include the following:

- Allergic reaction
- Bleeding
- Chronic nerve damage
- Erosion
- Excessive fibrotic tissue growth
- Extrusion
- Formation of hematomas or cysts
- Infection
- Keloid formation
- Migration

Clinician Use Information

Physician Training

Physicians should be familiar with sterile device insertion procedures and with follow-up evaluation and management of patients with an insertable cardiac monitor (or should refer the patient to such a physician). See the Monitoring Devices Help manual for instructions on programming the insertable cardiac monitor.

Package Contents

The device is supplied in a sterile tray for introduction into the operating field. The tray contains:

- One insertable cardiac monitor
- One incision tool
- One insertion tool

The outer box contains:

Literature

Pre-Insertion Device Setup

Prior to opening the sterile package, apply a magnet to the Assert-IQ™ device for at least three seconds and then remove. Interrogate the device with the Merlin™ PCS and follow the on-screen prompts. The user interface will provide a prompt to insert the device once the required information is entered. See the Monitoring Devices Help manual for instructions on programming the insertable cardiac monitor.

Opening the Sterile Package

To open the package:

- 1. Peel back the outer tray cover, starting with the corner labeled with an arrow.
- 2. Observing sterile technique, lift up the end of the inner tray that rests in the recess in the outer tray or flip over the outer tray so that the inner tray falls onto the table.
- 3. Peel off the inner tray cover, starting with the corner labeled with an arrow.
- 4. Use the recessed areas to facilitate removing the tools from the tray.

Choosing the Insertion Location

The Assert-IQ™ ICM is inserted under the skin in the left pectoral region. Common insertion locations are listed in the table below. Mapping may be beneficial as part of the pre-insertion process, especially for the Anterolateral, inframammary position between the 5th and 6th ribs.

Table 4. Insertion locations

4th intercostal space, 45° relative to the sternum, along axis of the heart

4th intercostal space, parallel to the sternum

Anterolateral, inframammary between the 5th and 6th ribs

NOTE:

- An implant site parallel to the midline, closer to the sternum and away from the lower half of the pectoral region and breast area may help minimize device movement.
- · Consider patient comfort when selecting an insertion location.

General Site Mapping

Either clinical ECG equipment or the Merlin™ PCS may be used to perform surface mapping. See the Monitoring Devices Help manual for a description of how to perform mapping using the Merlin PCS.

When performing pre-insertion surface mapping, use conductive ECG patches spaced approximately 4 cm apart. This distance approximates the Assert-IQ™ ICM inter-electrode spacing.

Evaluate potential insertion locations to optimize the following signal characteristics:

- High amplitude R-wave that demonstrates minimal variation in different patient positions, such as sitting versus lying down.
- High R-wave to T-wave ratio


Inserting the Device

To insert the device:

- 1. Prepare the insertion site using conventional antiseptic and local anesthetic procedures.
- 2. Pull back the skin and make an angled cut with the incision tool.


NOTE: Create the incision above the final insertable cardiac monitor position to reduce the force of gravity acting on the incision site during healing.

Figure 1. Pull back the skin and make incision



3. Hold the insertion tool as shown in the figure below, at approximately a 45 degree angle to the incision. Insert the blunt dissection tip of the insertion tool just past the skin, then adjust the angle to guide the insertion tool almost parallel with the patient's chest. This creates a subcutaneous device pocket parallel to the skin. Advance the insertion tool as far as it can go, until the flared edge contacts the incision site. Keep the flared edge in contact with the incision site for the subsequent steps.


Figure 2. Insert the introducer



NOTE: The insertion tool is designed to form an ideally sized pocket with some distance between the device and the incision. Avoid forming a pocket larger than the introducer.


4. Hold the insertion tool firmly at the incision site by gripping the ribbed finger recesses as shown in the following figure. Withdraw the plunger until the plunger stops and the preloaded device drops completely into the insertion channel.

Figure 3. Withdraw the plunger




5. Continue to hold the insertion tool firmly at the incision site. Advance the plunger to insert the device. The plunger has a feature to stop the introducer when the device is at the proper depth. Do not force the plunger past the stop feature.

Figure 4. Advance the plunger



6. Apply pressure to the incision site so that the device does not move, and then remove the insertion tool.

Figure 5. Remove the insertion tool



7. Measure R-wave amplitude and observe the signal quality on the Merlin™ PCS.

Consider body position and arm movement as part of this assessment. If signals are small or unstable, remove the device, reload the device in the insertion tool, and reposition the device. See the Monitoring Devices Help manual for information about adjusting sensitivity parameters to match implant conditions. If R-wave amplitudes are below 0.2 mV, consider repositioning the device.

NOTE: The device can be repositioned through the same incision, if desired. When repositioning the device, significantly change the insertion angle to ensure that the new device pocket does not merge with the previous pocket.


8. Close the incision using a preferred method (such as Dermabond, sutures, or staples). Device migration and loss may occur through an insufficiently closed incision. If topical skin adhesive is used, avoid getting adhesive on the sensing electrodes of the device.

Reloading the Device

To reload the device:

Place the device into the insertion tool so that the lasermark (2) on the device reads toward the top of the insertion tool (1). Ensure the device header faces the proximal end of the insertion tool, with the large electrode facing up. Angle the device and position it under the first tab (3). Then press the device under the second tab (4) until it snaps into place. Refer to the figure below.

Figure 6. Device placement




- 1. Insertion tool
- 2. Lasermark
- First tab
- 4. Second tab

Explanting the Device

Explant the device using standard surgical tools.

NOTE: Grasp the device with the tool held parallel to the device. Do not grasp the device at an angle, as this could damage the device.

Figure 7. Explant technique



Patient Education

Abbott Medical provides a booklet for patients to explain the device and its operation, in addition to a patient manual to explain the myMerlin™ mobile application. You can use this to supplement your discussions with the patient, and caregiver or other interested persons. To obtain other available patient education materials, contact Abbott Medical.

Patient Identification Card

A patient identification (ID) card should be provided to all patients with an Assert-IQ™ ICM. The ID card indicates that the patient has an inserted cardiac monitor.

Radiopaque Identification

Each device has an x-ray absorptive marker for non-invasive identification. The two-letter model code is visible on a radiograph.

Table 5. X-ray ID code for Assert-IQ™ device

Device Model	X-ray ID Model Code
DM5000, DM5300	DS
DM5500	DL

Additional Information

The Volt™ PFA Generator tab contains the serial number and device version information of the Volt™ PFA Generator.

For additional information on this device, see the Monitoring Devices Help manual or the Merlin™ PCS on-screen help.

To enable communication with the Merlin.net™ Patient Care Network, the insertable cardiac monitor must be configured for use with the myMerlin™ mobile application on a compatible mobile device. See the myMerlin mobile application user's manual for instructions.

Physical Specifications

Device Specifications

Table 6. Device specifications DM5300

Specification ¹	Data
Dimensions (h x l x t) (mm)	46 x 9.4 x 3.1

¹ The dimensions, weight, and displacement volume are nominal values based on engineering model measurements.

Table 6. Device specifications DM5300

Specification	Data
Weight(g)	3.0
Displacement volume (cm ³)	1.2
Surface area of can electrode (mm²)	105.9
Surface area of header electrode (mm ²)	37
Shortest distance between electrodes (mm)	38.10
Can and electrode material	Titanium
Header material	Polyurethane and epoxy
Coating	Parylene
Table 7. Device specifications DM5000 and DM5500	
Specification	Data
Dimensions (h x l x t) (mm)	49 x 9.4 x 4.4
Weight(g)	3.0

Table 7. Device specifications DM5000 and DM5500

Specification	Data
Displacement volume (cm ³)	1.9
Surface area of can electrode (mm ²)	105.9
Surface area of header electrode (mm²)	37
Shortest distance between electrodes (mm)	39.78
Can and electrode material	Titanium
Header material	Polyurethane and epoxy
Coating	Parylene

Battery Specifications

Table 8. Battery specifications

Parameter	Data
Manufacturer	Eagle Picher

Table 8. Battery specifications

Parameter	Data
Model	ICM Battery
Chemistry	CFx
Number of cells	One cell
Battery voltage (beginning of service)	3.40 V
Elective replacement voltage (ERI)	2.81 V
End-of-service voltage (EOS)	2.67 V
Longevity	DM5000 and DM5300 3 years, DM5500 5 years, under the following usage scenarios: • Average of 1 auto-detected episode per day
	 Average of 1 patient-activated symptom episode per month
	 Up to 6 month shelf storage time
	NOTE: At a maximum shelf storage time of 18 months, longevity is reduced by approximately 5 months.

FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radiofrequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Operation is subject to the following two conditions:

- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired
 operation.

Modifications not expressly approved by the manufacturer could void the user's authority to operate the equipment under FCC rules.

Statement of Compliance with License-Exempt RSS Standard (Canada)

This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

Le présent appareil est conforme aux CNR d'Industrie Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes : (1) l'appareil ne doit pas produire de brouillage; (2) l'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement.

CAN ICES-3 (B)/NMB-3(B)

Identification Information for Product Registration

This device has a label that contains, among other information, a product identifier in the following format:

Table 9. Registration identification information

Identifier Type	Registration Identifier
FCC registration number	2A76T-ICM2ABT
Industry Canada (IC) registration number	7067A-ICM2ABT

Wireless Technology Information

The following table summarizes the technical details of the Bluetooth $^{\circ}$ low energy (Bluetooth LE) technology as it is implemented in the device.

Table 10. Bluetooth low energy information

Parameter	Data
Antenna type	Embedded antenna in header
Antenna dimensions	0.375 in x 0.05 in
Modulation	GFSK
Magnetic field strength (at 2 m distance)	92.30 uA/m
Electric field strength (at 2 m distance)	34.79 mV/m
Output power (EIRP*)	0.2 mW maximum
Range	2 m typical
Center frequency	2.44 GHz

Table 10. Bluetooth low energy information

Parameter	Data
Frequency (range)	2.4000 to 2.4835 GHz
Bandwidth	2 MHz per channel
Bandwidth(-15dB)	2.398 to 2.4855 GHz
Data flow	Bi-directional and semi-duplex
Duty cycle	Variable, but low (<12%)
Bandwidth	2 MHz per channel
Protocol	Bluetooth LE
*EIRP = Equivalent isotrop	cally radiated power

The radio receiver in the device is using the same frequency and bandwidth as the transmitter.



Cables and transducers:

Cables and transducers are not used during normal use of the device nor while programming the device.

Quality of Service for Wireless Technology

Bluetooth® low energy (Bluetooth LE) wireless technology enables communication between the monitor and the clinician programmer, smart phone, or tablet. The requirements for the quality of service (QoS) vary depending on the use environment (operating room, recovery room, and home environment).

After the clinician programmer, smart phone, or tablet is paired with a monitor, the Bluetooth symbol is visible on the clinician programmer, smart phone, or tablet. When the Bluetooth LE connection is not active, the symbol appears dimmed.

Other requirements include a semi-duplex transmission with a required acknowledge, a transmission latency in each direction (2x), and a receive-to-transmit mode (RX-to-TX) time. Data is resent if not sent successfully.

Wireless Security Measures

The insertable cardiac monitor is designed to include several features to enhance the security of wireless communications. Design features include the following:

- The monitor encrypts its wireless communication.
- The monitor is designed to limit its communications to authenticated paired smart phones or tablets.
- The monitor is designed to pair with a single smart phone or tablet at a time.
- The monitor uses a proprietary pairing protocol in addition to the pairing procedure specified in Bluetooth® low energy protocols.
- The monitor authenticates the pairing requests using a standard cloud-based authentication.
- The monitor uses an authorization protocol, which limits a paired smart phone or tablet's access to data appropriate for its functionality.
- The monitor creates a unique key for the paired unit and verifies it at the onset of every communication. If the unique key is not verified, the monitor denies access.

Technical Support

Abbott Medical maintains 24-hour phone lines for technical questions and support:

- 1 818 362 6822
- 1 800 722 3774 (toll-free within North America)
- + 46 8 474 4147 (Sweden)
- +61 2 9936 1200 (Australia)
- medical.abbott/manuals

For additional assistance, call your local Abbott Medical representative.

Symbols

The symbols below and harmonized symbols may be found on the product or product label. For harmonized symbols, refer to the Universal Symbols Glossary at medical, abbott/manuals.

Symbol

Description



Follow instructions for use on this website

medical.abbott/manuals

Symbol	Description		
9 3	Manufacturing facility		
	Product literature		
+	Accessories		
	Incision tool		
	Insertion tool		
☐ Abbott Jot Dx™ ੈ D	Insertable Cardiac Monitor		

Symbol	Description
	Warning; sharp element
⊕ OFF	Shipped settings off
FCC ID: 2A76T-ICM2ABT	Federal Communication Commission Number (FCC ID: #)
IC: 7067A-ICM2ABT	Industry Canada Radio Communications License (IC: #)

Symbol

Description



The device contains a battery and the label is affixed to this device in accordance with European Council Directives 2002/96/EC and 2006/66/EC.

These directives call for separate collection and disposal of electrical and electronic equipment and batteries. Sorting such waste and removing it from other forms of waste lessens the contribution of potentially toxic substances into municipal disposal systems and into the larger ecosystem.

Return the device to Abbott Medical at the end of its operating life.



Australian Communications and Media Authority (ACMA) and New Zealand Radio Spectrum Management (RSM) Regulatory Compliance Mark (RCM)



This equipment is certified for type certification pursuant of Article 38-24 of the Japan Radio Law

Symbol	Description
	Korea Certification mark for electrical devices
MEME	Malaysian Communication and Multimedia Commission (MCMC) certification mark for products meeting applicable MCMC Technical Codes
	Prescription only
Made in USA	Made in USA
-	Unique Device Identifier
REF	Catalog number

Symbol	Description
31	Date
W	Healthcare center or physician
† ?	Patient identification
<u>=</u> &	Patient identification card label
C	Physician telephone
0	Location of implant

Abbott Medical 15900 Valley View Court Sylmar, CA 91342 USA +1 818 362 6822

www.abbott.com

2022-10 ARTEN600265177 02



