
# Booklet-Page 1

## Welcome to Your HRV Coherence Pro

The HRV Coherence Pro helps you activate a healthy, high performance state called Heart Coherence, where the heart, mind and emotions are in sync and operating at optimal efficiency.

The sensor measures and analyzes your heart rate variability (HRV), and translates it into a real time coherence score. This sensor helps you learn what heart coherence feels like and how to increase it.

The Inner Balance App - Your personal coherence trainer to experience more ease, calm, clarity and focus. A few minutes a day can improve mood, decisions, sleep, meditation, intuitive guidance and more.

The Global Coherence App - See and track how much heart coherence you are adding to the global energetic field. Also use to measure team, group or family coherence.

## Getting Started

1. Download the Inner Balance app and watch the onboarding video. App Store for iOS, Google Play for Android or Amazon appstore for Amazon Fire.

2. Plug the Ear Sensor Cable into the round Ear Sensor Cable Port.

3. Clip the Ear Sensor to your earlobe.

4. Connect to the Inner Balance app: - Bluetooth (wireless): Press the Power Button on the Sensor Pod. Press the start button on your app. The app will start scanning. Choose the number that matches the Bluetooth Pairing Number on the back of your Sensor Pod.

-USB-C (wired): Plug the USB-C Charging/Connector Cable into the rectangular port located below the power button and the other end of the cable into your Android device.

5. Once connected, the white Add Heart Light on the Sensor Pod changes from blinking to solid and the Ring Light around the Power Button begins to pulse blue.

6. Secure the Sensor Pod to your shirt via the metal Clothing Clip.


7. When you complete your session, explore the guides, resources, journaling, tutorials, tracking your progress and much more.

Enjoy your journey!

Your HRV Coherence Pro sensor works with all HeartMath apps on both iOS and Android devices, as well as the emWave2, emWave Pro and Pro Plus desktop software for Mac and Windows.

To learn more: <https://www.heartmath.com/apps/>

## Sensor Diagram

## LIGHTS: Add Heart Light

OFF	When powered off
BLINKS WHITE	"Ready to Connect" to app (Bluetooth)
SOLID WHITE	Connected

### A few notes:

- If the sensor lights do not come on, charge with the included USB-C Charging/Connector cable (2 hrs = full charge).
- HeartMath apps have a battery power meter to see your device's charging status.
- Bluetooth is disabled when sensor is plugged into power.

## LIGHTS: RING Light (Power)

OFF	When powered off
PULSES BLUE	During a session
PULSES RED	Low battery
BLINKS RED	While charging
SOLID GREEN	Fully Charged

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