READ THIS MANUAL COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT

Keep this manual in a safe location for future reference

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CONTENT

1. IMPORTANT SAFETY INSTRUCTIONS.....	2	8. Other functions.....	9
2. Introduction	4	9. Care and Maintenance	9
3. Package Contents.....	5	10. Troubleshooting Guide.....	9
4. Parts and Controls.....	5	11. FCC Compliance Statement	10
5. Initial Use	5	12. Technical Specifications	10
6. Taking Measurements	6	13. Warranty	11
7. Evaluating Results.....	7		

SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE

If the instruction manual is damaged or if you no longer have the instruction manual in your possession, please contact Beurer customer service. Please see warranty for service contact.

Dear Customer,

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for applications in the areas of heating, weight management, blood pressure, body temperature, pulse, gentle therapy, massage, and beauty.

Please read these instructions for use carefully and keep them for later use. Make them accessible to other users and observe the information they contain.

With kind regards,
Your Beurer team

1. IMPORTANT SAFETY INSTRUCTIONS

Signs and symbols

The following signs appear in the Safety Section (page 3) and in this manual on page 3.

READ AND UNDERSTAND THIS ENTIRE MANUAL INCLUDING THE SAFETY SECTION AND ALL INSTRUCTIONS AND WARNINGS COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT. FOLLOW ALL SAFETY INSTRUCTIONS AND WARNINGS TO AVOID HAZARDOUS SITUATIONS AND TO MAKE CORRECT USE OF THIS PRODUCT.

This is the safety alert symbol. It alerts you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury or death.

WARNING

WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.


CAUTION

CAUTION indicates a hazardous situation which, if not avoided, may result in minor or moderate injury.

NOTICE

NOTICE addresses practices not related to personal injury, such as product and/or property damage.

	Caution	IP21	Protected against access to hazardous parts with a finger and against vertically falling water drops
	Follow instructions for use	Storage 	Permissible storage temperature and humidity
	Type BF applied part	Operating 	Permissible operating temperature and humidity
	Do not use the device with a pacemaker		

⚠ WARNING – To reduce the risk of fire, electric shock, or serious personal injury:

- Not to diagnose or treat any medical condition.
- The scale must not be used by people with medical implants (e.g. pacemakers), as this may affect their functionality.
- Do not use during pregnancy.
- To prevent falls, do not step onto the scale with wet feet or if the surface of the scale is damp.
- Do not stand on the outer edge of the scale, otherwise it may tip.
- Keep children away from packaging materials (risk of suffocation).

⚠ CAUTION – To reduce the risk of personal injury or product/property damage:

- This product is intended solely for household use and not for medical or commercial purposes.
- Place the scale on a hard, flat, level surface; a hard surface is required for accurate measurements.
- Protect the device from impacts, moisture, dust, chemicals, large temperature fluctuations, and nearby sources of heat (ovens, heaters).
- Do not use on children under 10 years of age.

NOTICE

- The device should be cleaned from time to time. Do not use any abrasive cleaning products and never submerge the device in water.
- Measurement variations are normal, as this scale is not calibrated for use in a professional medical environment.
- Make sure that no liquids come into contact with the scale. Never submerge the scale in water. Never rinse it in running water.
- Do not place any objects on the scale when not in use.
- Protect the device from knocks, damp, dust, chemicals, marked temperature fluctuations and nearby sources of heat (ovens, heaters).
- Do not press the buttons violently or with pointed objects.
- Do not expose the scale to high temperatures or strong electromagnetic fields.
- The scale's capacity is 400 lb (181.4 kg, 28.5 st). The results for weight and bone mass are displayed in 0.2 lb (100g) increments.
- Measurements of body fat, body water and muscle content are displayed in 0.1% increments.
- Calorie requirement is indicated in steps of 1 kcal.
- The scale is initially set to weigh and measure in "lb" and "in".
- Repairs must only be performed by Beurer Customer Service.

Battery Handling Safety Precautions

- Use only the size and type of batteries specified.
- Be sure to follow the correct polarity when installing the batteries. Reversed batteries may cause damage to the device.
- Do not mix different types of batteries together (e.g. Alkaline and Carbon-zinc or rechargeable batteries) or old batteries with fresh ones. Always replace batteries as a simultaneous set.

- If the batteries in the device are depleted or the device will not be used for a long period of time, remove the batteries to prevent damage or injury from possible battery leakage.
- Do not try to recharge batteries not intended to be recharged; they can overheat and rupture (follow battery manufacturer's directions).
- Do not dispose of batteries in fire, batteries may explode or leak.
- Clean the battery contacts and also those of the device prior to battery installation.
- Remove discharged batteries from the product and dispose/recycle in compliance with all applicable laws.
- Keep batteries away from children and pets. Batteries may be harmful if swallowed. Should a child or pet swallow a battery, seek medical assistance immediately.

2. Introduction

This digital body analysis scale is intended for measuring body weight and providing an analysis of your personal fitness data. It is intended for household use only. In addition to measuring body weight, the scale has the following functions that can be used by up to eight persons:

- Body weight in lb/kg/st
- Body fat percentage BF in %
- Body water percentage \approx in %
- Muscle percentage ✚ in %
- Bone mass ☞ in kg
- Basal metabolic rate and active metabolic rate

The scale also includes the following features:

- Switches between pounds (lb), kilograms (kg) and stones (st)
- Automatic shut-off
- Low battery indicator
- Automatic user recognition
- Stores last 30 measurements for each of eight users
- *Bluetooth*[®] communication via our free 'beurer HealthCoach' app

System requirements for the 'beurer HealthCoach' app

Bluetooth[®] 4.0, iOS from version 8.0, Android[™] from version 4.4 with *Bluetooth*[®] low energy technology.

List of compatible devices:

This scale works on the principle of Bioelectrical Impedance Analysis (BIA). This involves a quick calculation of body content using a harmless electric current that is not even noticeable.

When measurement of the electrical resistance (impedance) is considered alongside other data (age, height, gender, and activity level), body fat percentage and other variables in the body can be estimated. Muscle tissue and water conduct electricity well and therefore have a lower resistance. In contrast, bones and fat tissue have low conductivity as fat cells and bones only minimally conduct current due to very high levels of resistance.

Calculated values only approximate the actual values for your body. Only a doctor can precisely establish body fat, body water, muscle percentage and bone structure using medical methods.

Tips

- If possible, always weigh yourself at the same time every day (ideally in the morning), after having used the restroom, on an empty stomach and with as little clothing as possible, in order to obtain comparable, consistent results.

- Body fat percentage can only be calculated when you are barefoot; the soles of the feet may be lightly dampened where appropriate. Results may be unsatisfactory if the soles of your feet are completely dry or feature areas of rough skin; this reduces conductivity.
- Stand up straight and still while weighing.
- Wait a few hours after significant physical exertion.
- Wait about 15 minutes after waking so that your body's water content can normalize.
- Remember that only long-term trends are important. Brief deviations in weight within a few days are normally caused by changes in body water content.

Limitations

When calculating body fat and additional values, deviating and unexpected results may occur in the case of:


- Children under 10 years of age
- Athletes and bodybuilders
- Pregnant women
- Persons with a fever, persons undergoing dialysis, persons demonstrating edema symptoms, and persons suffering from osteoporosis
- Persons taking cardiovascular medicine (affecting the heart and vascular system)
- Persons taking vasodilating or vascular constricting medication
- Persons with significant anatomical deviations in their legs compared with their overall body height (leg length significantly shorter or longer).

3. Package Contents

- 1 x BF70 Body Analysis Scale
- 3 x 1.5 V AAA-size Alkaline Batteries
- 1 x Instruction manual
- 1 x Quick start guide

4. Parts and Controls

1. Battery Compartment (at rear)
2. Reset Button (at rear)
3. *Bluetooth*® Pairing Icon
4. User Initials (if recognized)
5. Body weight
6. Other Measurement Data (BMI, Body Fat Percentage, etc.)

5. Initial Use

Inserting the Batteries

Remove the Battery Compartment Cover and insert three AAA-size alkaline batteries (included) according to the indicated polarity. If the scale does not function, remove the batteries and reinsert them, ensuring proper polarity. Replace Battery Compartment Cover.

Positioning the Scale

Place the scale on a hard, flat surface with no carpet or rugs; a hard surface is required for accurate measurements.

Downloading and Using the App

Download and install the free 'beurer HealthCoach' app from Google Play (for Android devices) or the Apple App Store (for iOS devices). Start the app and enter your personal user data (height, birthdate, gender, etc.). Be sure to select "BF70" as your device before performing any measurements. If you do not use the app to enter your user data, the scale cannot recognize you and store your measurements. Enter the requested user data.

User data	Parameters
Initials/abbreviation of name	max. 3 characters or numbers
Body height	3 ft 3.5 in to 7 ft 2.5 in (100 to 220 cm)
Age	10 to 100 Years
Sex	man (♂), woman (♀)
Degree of activity	1 to 5

To select the proper Activity Level, consult the table below.

Activity Level	Physical Activity
1	None
2	Low: Small amount of light physical exertion (e.g. going for walks, easy gardening, simple gymnastics)
3	Medium: Physical exertion at least 2 to 4 times a week for 30 minutes each
4	High: Physical exertion at least 4 to 6 times a week for 30 minutes each
5	Very High: Intense physical exertion, training or hard physical work at least 1 hour every day

The scale has eight user memory locations that allow, for example, you and your family members to save personal settings. You can also activate users via the 'beurer HealthCoach' app on other mobile devices or by changing the user profile in the app (see the app settings).

Keep your smartphone close to the scale in order to maintain an active *Bluetooth*[®] connection.

To set automatic user recognition, the first measurement must be assigned to your personal user data; see below for instructions.

6. Taking Measurements

NOTE: Do not let your feet, legs, calves or thighs touch each other. Otherwise, the measurement cannot be taken correctly. The measurement will also be incorrect if taken while you are wearing socks.

1. Step onto the scale **with bare feet**; stand still with your weight distributed evenly on both sets of sensors.
2. The scale first displays the weight. If you are barefoot, BIA diagnostics will follow.
3. If a user has been assigned, the user ID displayed in the bottom left corner (e.g., "TOM").

The following values appear:

1. Weight in kg with BMI
2. Body fat in % with interpretation
3. Body water in % 
4. Muscle percentage in % 
5. Bone mass in kg 
6. Basal metabolic rate in kcal (BMR)
7. Active metabolic rate in kcal (AMR)

If the user is not recognized, only the weight is shown and "---" appears in the User ID area of the LCD display.

Weight Measurement Only

Step on the scale wearing shoes. Stand still with your weight distributed equally on both feet. The scale will display your weight after a few seconds. If a stored user is recognized, their BMI, BMR, and AMR are also displayed along with their initials.

NOTE: Other values cannot be calculated if user is wearing shoes.

The scale will shut off by itself after a few moments of inactivity.

7. Evaluating Results

Body Fat Percentage

The following body fat percentages are a guideline only; for more information, ask your doctor.

Men

Age	Very good	Good	Average	Poor
10-14	<11 %	11-16 %	16.1-21 %	>21.1 %
15-19	<12 %	12-17 %	17.1-22 %	>22.1 %
20-29	<13 %	13-18 %	18.1-23 %	>23.1 %
30-39	<14 %	14-19 %	19.1-24 %	>24.1 %
40-49	<15 %	15-20 %	20.1-25 %	>25.1 %
50-59	<16 %	16-21 %	21.1-26 %	>26.1 %
60-69	<17 %	17-22 %	22.1-27 %	>27.1 %
70-100	<18 %	18-23 %	23.1-28 %	>28.1 %

Women

Age	Very good	Good	Average	Poor
10-14	<16 %	16-21 %	21.1-26 %	>26.1 %
15-19	<17 %	17-22 %	22.1-27 %	>27.1 %
20-29	<18 %	18-23 %	23.1-28 %	>28.1 %
30-39	<19 %	19-24 %	24.1-29 %	>29.1 %
40-49	<20 %	20-25 %	25.1-30 %	>30.1 %
50-59	<21 %	21-26 %	26.1-31 %	>31.1 %
60-69	<22 %	22-27 %	27.1-32 %	>32.1 %
70-100	<23 %	23-28 %	28.1-33 %	>33.1 %

Athletes will usually have lower values. Depending on the sport, training intensity, and your physical constitution, results may still be below the stated ranges. Please note that there may be health risks in the case of extremely low values.

Body Water Content

Body water content in percent is normally within the following ranges:

Men

Poor	Good	Very good
<50%	50-65%	>65%

Women

Poor	Good	Very good
<45%	45-60%	>60%

Body fat contains relatively little water; therefore, body water content may be below the standard values in persons with a high body fat percentage. On the other hand, some athletes may exceed the standard values due to having low levels of fat and a high muscle percentage.

The body water calculation performed using this scale is not suitable for drawing medical conclusions, such as in the case of age-related water retention. Consult your doctor where necessary. A high body water content is generally desirable.

Muscle Percentage

Muscle percentage is normally within the following ranges:

Men

Age	Low	Normal	High
10-14	<44%	44-57%	>57%

Women

Age	Low	Normal	High
10-14	<36%	36-43%	>43%

15-19	<43%	43-56%	>56%
20-29	<42%	42-54%	>54%
30-39	<41%	41-52%	>52%
40-49	<40%	40-50%	>50%
50-59	<39%	39-48%	>48%
60-69	<38%	38-47%	>47%
70-100	<37%	37-46%	>46%

15-19	<35%	35-41%	>41%
20-29	<34%	34-39%	>39%
30-39	<33%	33-38%	>38%
40-49	<31%	31-36%	>36%
50-59	<29%	29-34%	>34%
60-69	<28%	28-33%	>33%
70-100	<27%	27-32%	>32%

Bone Mass

Like the rest of our body, bones are subject to natural growth, shrinking, and aging processes. Bone mass increases rapidly during childhood and reaches its peak between the ages of 30 and 40. As we age, our bone mass begins to decrease. You can combat this reduction to an extent with the help of a healthy diet (particularly calcium and vitamin D) and regular physical exercise. The stability of the skeletal system may be increased through muscle buildup. Note that this scale does not measure calcium content of bones; instead, it calculates the weight of all the components that make up the bones (organic matter, inorganic matter and water). It is very difficult to affect your bone mass, although it does fluctuate within the scope of influencing factors (weight, height, age, gender). There are no recognized guidelines or recommendations.

IMPORTANT:

- Do not confuse bone mass with bone density.
- Bone density can only be determined by a medical examination, so you cannot draw conclusions on changes to bones and bone hardness (e.g. osteoporosis) using this scale.

BMR

The basal metabolic rate (BMR) is the amount of energy required by the body at complete rest in order to maintain its basic functions. This value primarily depends on weight, height and age. It is displayed in kcal/day and calculated using the scientifically recognized Harris-Benedict equation. This is the amount of energy that is required by your body under all circumstances and must be resupplied to the body in the form of food. If your energy intake is below this level long term, it can affect your health.

AMR

The active metabolic rate (AMR) is the amount of energy that an active body consumes each day. An individual's energy requirement increases when the level of physical activity increases; the scale calculates this by means of the activity level (1- 5) that has been entered.

To maintain the current weight, the energy that the body uses must be replaced in the form of food and drink. If less energy is taken in than is used over a long period of time, the body takes the difference from its fat stores, and weight loss results. However, if the amount of energy taken in exceeds the calculated active metabolic rate (AMR) for a longer period, the body cannot burn off the excess energy. The excess is stored in the body as fat, leading to weight gain.

Interpreting the Results

NOTE: Only long-term trends are important. Brief deviations in weight within a few days are normally caused by loss of fluids.

Interpretation of these results is based on changes in total body weight, percentage of body fat, body water and muscle content, as well as the length of time over which these changes occur. Rapid changes within the scope of a few days are considered separately from medium-term changes (in the scope of weeks) and long-term changes (months). As a basic rule, short-term changes in weight are almost entirely changes in water content, whereas medium-term and long-term changes may indicate changes in fat percentage and muscle percentage.

- If you experience short-term weight loss but your body fat percentage increases or stays the same, you have only lost water, such as after a training session, visit to the sauna, or a diet aimed only at fast weight loss.

- On the other hand, if your weight increases in the medium term and your body fat percentage drops or stays the same, you may have increased your muscle mass.
- If your weight and body fat percentage both fall at the same time, then your diet is working – you are losing fat mass.
- Ideally, you should support your diet with physical activity, fitness or strength training. This enables you to increase your muscle percentage in the medium term.
- Body fat, body water, and muscle percentages should not be totaled, since muscle tissue also contains components made of body water.

8. Other functions

User assignment

It is possible to assign up to eight users for weight-only measurement (with shoes) or analysis (bare feet). For a new measurement, the scale assigns a measurement to the user whose most recently saved measurement is within 4 lb or 2 kg of the active measurement and within 2% of the body fat.

Unknown measurements

If a measurement cannot be assigned to a known user, it is stored as an unknown measurement. Up to 20 unknown measurements can be saved. Use the app to directly assign unknown measurements to a user.

Saving the assigned measurements on the scale

If the app is open and your mobile device is paired with the scale via *Bluetooth*[®], newly assigned measurements are sent to the app and not saved on the scale. If the app is not open, the newly assigned measurements are saved on the scale. A total of 30 measurements per user can be saved. The saved measurements are transferred automatically to the app if you open the app within *Bluetooth*[®] range. Data synchronization takes place within 10 seconds. It is also possible to transfer data when the scale is switched off.

The *Bluetooth*[®] word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer North America LP is under license. Other trademarks and trade names are those of their respective owners.

Deleting scale data

To delete all measurements and old user data on the scale, switch on the scale and press and hold the reset button for three seconds; “dEL” will appear in the LCD Display.

9. Care and Maintenance

Replacing Batteries

Your scale is equipped with a low battery indicator. If the batteries are weak, “Lo” will appear in the LCD Display and the scale will automatically shut off. In this case, replace the batteries with three new AAA-size alkaline batteries.

Cleaning

Clean the scale periodically using a damp cloth and, if necessary, a small amount of mild detergent.

10. Troubleshooting Guide

If the scale encounters an error during measurement, the following is displayed:

Display	Cause	Remedy
“- - -”	Unknown measurement as it is outside the user assignment limit or a unique assignment is not possible.	Assign unknown measurement in app or repeat user assignment.

"- - -" runs across the LCD.	The fat percentage is outside the measurable range (less than 3 % or greater than 65 %).	Repeat the measurement barefoot or, if necessary, moisten the soles of your feet slightly.
Err	The maximum load-bearing capacity of 400 lb (181.4 kg) was exceeded.	The load must not exceed 400 lb (181.4 kg).
Err or incorrect weight is displayed.	No flat, stable surface.	Place the scale on a flat, stable surface.
Err or incorrect weight is displayed.	Not standing still.	Stand as still as possible.
Incorrect weight is displayed.	Scale zero setting is incorrect.	Wait until the scale switches itself off. Activate scale, wait for "0.0 lb" to appear and repeat measurement.
No Bluetooth® connection (Ⓜ symbol is missing).	Device outside the range.	The minimum range in a free field is approx. 82 ft (25 m). Walls and ceilings reduce the range. Other radio waves may interfere with the transmission. Therefore, do not position the scale near devices such as WLAN routers, microwaves, induction hobs etc.
FULL	User memory space is full. No more measurements can be saved.	Open the app. The data is transferred automatically. This may take up to one minute.
Lo	The batteries in the scale are empty.	Replace the batteries of scale.

11. FCC Compliance Statement

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

12. Technical Specifications

Battery type	3 x 1.5 V AAA-size Alkaline Batteries
Measuring range	4 lb to 400 lb (2 kg to 181.4 kg)
Weight unit	0.2 lb increments (100 g)
Body fat, water, muscle percentage	unit 0.1 % increments
Electrical current of B.I.A.	max. 5 µA / 1 V
Fat percentage range	5 % to 55 %

Water percentage range	25 % to 75 %
Muscle percentage range	25 % to 80 %
Temperature operating range	41 °F to 95 °F (5 °C to 35 °C), 20-95%RH
Storage temperature range	-4 °F to 140 °F (-20 °C to 60 °C), 20-95%RH
Subject to technical modifications.	

13. Warranty

Limited Lifetime Warranty For Original Purchaser

Your Beurer Body Analysis Scale, Model BF70, is warranted to be free from defects in materials and workmanship for the life of the product under normal conditions of intended use and service. This warranty extends only to the original retail purchaser and does not extend to retailers or subsequent owners.

We will, at our option, repair or replace the Beurer Body Analysis Scale, Model BF70, without additional charge, for any part or parts covered by these written warranties. No refunds will be given. Repair or replacement is our only responsibility and your only remedy under this written warranty. If replacement parts for defective materials are not available, Beurer reserves the right to make product substitutions in lieu of repair or replacement. For warranty service contact our customer service department at **1-800-536-0366** or at **info@beurer.com** to provide a description of the problem. If the problem is deemed to be within the scope of the limited lifetime warranty, you will be asked to mail the product at your costs in its original package with proof of purchase, your name, address and phone number. If the problem is not deemed to be within the scope of the limited lifetime warranty, we will provide a quotation for repair respectively replacement and return shipping fee.

This warranty does not cover damage caused by misuse or abuse; accident; the attachment of unauthorized accessory; alteration to the product; improper installation; misapplication; lack of reasonable care with respect to the product; unauthorized repairs or modifications; improper use of electrical/power supply; normal wear; loss of power; dropped product; malfunction or damage of an operating part as a result of failure to comply with instructions for use or to provide manufacturer's recommended maintenance; transit damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of Beurer. This warranty is void if the product is ever used in a commercial or business environment. The maximum liability of Beurer under this warranty is limited to the purchase price actually paid by the customer for the product covered by the warranty, as confirmed by proof of purchase, regardless of the amount of any other direct or indirect damage suffered by the customer.

This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adaptation to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. ANY IMPLIED WARRANTIES, OBLIGATIONS, OR LIABILITIES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED IN DURATION TO THE DURATION OF THIS APPLICABLE WRITTEN WARRANTY. Some states do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you.

IN NO EVENT SHALL BEURER BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS, IMPLIED OR ANY OTHER THEORY OF LIABILITY, WHATSOEVER. Some states do not allow the exclusion or limitation of special, incidental, or consequential damages, so the above limitation may not apply to you.

Beurer does not authorize anyone, including, but not limited to, retailers, the subsequent consumer purchaser of the product from a retailer or remote purchaser, to obligate Beurer in any way beyond the terms set forth herein.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resealed products, including but not limited to sale of such products on Internet auction sites and/or products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate in connection with any products or parts thereof which are repaired, replaced, altered, or modified, without the prior explicitly written consent of Beurer.

This warranty gives you specific legal rights, and you may also have other rights which may vary from state to state.

Manufactured for:
Beurer North America LP,
900 N Federal Hwy Ste 300
Hallandale Beach, FL 33009 USA
Made in China

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