

MPI, Carriage Integration, USM, IS4000 MPI#: 844505-01 REV. L 1 of 20

SURGICAL

INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.

ITEM	PART#	DESCRIPTION	QTY
A	723618	TOOL, ASSY, SPRING BRAKE RELEASE, IS4000	NA


Note: Brake Release Tool (\mathbf{A}) may be used to release the Brake and position the Arm throughout this MPI.

1. Verify **Screws** holding **Carriage Top Plate** are tight. Record result in Results Sheet.




CAUTION: Once Carriage Cover Aide is removed the fragile ACCP LEDs will be

is removed the fragile ACCP LEDs will be exposed. Handle Carriage with care to avoid damage.



2. Carefully remove Carriage Cover Aide from Carriage.




Intuitive Surgical Confidential – This document contains confidential and proprietary information, and except with written permission of Intuitive Surgical, such information shall not be published or disclosed to others or used for any purpose other than on behalf of Intuitive Surgical. This document shall not be copied in whole or in part, including electronically. MPI Template: 800103



MPI, Carriage Integration, USM, IS4000 MPI #: 844505-01 REV. L 2 of 20

INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.

- 3. Carefully remove **Kapton tape** and **Protection Clip** securing the fiber cables to the FFCs.
- 4. Verify **FFC's** are in order 1, 2, 3, 4, and 5.
- 5. Fold **FFC** 90° flush with bottom edge of **Bracket**.




Intuitive Surgical Confidential – This document contains confidential and proprietary information, and except with written permission of Intuitive Surgical, such information shall not be published or disclosed to others or used for any purpose other than on behalf of Intuitive Surgical. This document shall not be copied in whole or in part, including electronically. MPI Template: 800103

Align Carriage to Insertion Assembly as shown.

CAUTION: Be careful working with Fiber Cables. Avoid extreme bends or pinching of cables.







MPI, Carriage Integration, USM, IS4000 REV. L


MPI #: 844505-01


4 of 20

SURGICAL

INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.


ITEM	PART#	DESCRIPTION	QTY
A	902947	SCREW, SHCS, M2X5, STST	2
В	902460	SCREW, SHCS, M1.6X12, STST	1
С	721916	EFT, TORQUE WRENCH, 1.0 IN-LBF, 1/4" FEM HEX ADAPTER	NA
D	722555	EFT, TORQUE WRENCH, 1.6 IN-LBF, 1/4" FEM HEX ADAPTER	NA







Install Screw (B) with Loctite 222. Torque Screws (A) to 1 in-lbs using (C). Torque Screws (B) to 1.6 in-lbs using (D).





MPI, Carriage Integration, USM, IS4000

MPI #: 844505-01

REV. L 5 of 20


SURGICAL

INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.

ITEM	PART#	DESCRIPTION	QTY
A	318276	CLAMP,FFC,IS4000 CARRIAGE	1
В	902463	SCREW, PHFH, 0-80 X .125, STST	2
С	721916	EFT, TORQUE WRENCH, 1.0 IN-LBF, 1/4" FEM HEX ADAPTER	NA

Route FFC's against Carriage. Secure by installing Clamp (A) using Qty: 2 Screws (B) Loctite 222. Torque to 1.0 in-lbs using (C).









MPI, Carriage Integration, USM, IS4000

MPI #: 844505-01 R


REV. L 6 of 20

SURGICAL


INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.

ITEM	PART#	DESCRIPTION	QTY
A	963004	TAPE,POLYIMIDE (KAPTON),ROLL,.001 X.50W	NA

10. Route **Fiber Cables** on **Carriage** and secure with **Kapton Tape** (A).







MPI, Carriage Integration, USM, IS4000

MPI #: 844505-01


REV. L

SURGICAL

INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.

11. Apply a service loop to route and connect the **Fiber Cables** to **ACCL**, **PCA**.

Note: "Arrow up" Cable is 'Receive' and connects to the 'RX1' connector.





MPI, Carriage Integration, USM, IS4000 MPI #: 844505-01 REV. L

SURGICAL


INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.


ITEM	PART#	DESCRIPTION	QTY
A	963004	TAPE,POLYIMIDE (KAPTON),ROLL,.001 X.50W	NA

12. On both sides of the Carriage Assembly, tuck the Motor Wires in between the ACCP and ACCL, PCAs and secure with a strip of Kapton Tape (A) to avoid being damaged. See photos below:

Note: The service loop of the Fiber Cables are placed under the ACCP, PCA before the Motor Wires.

Tuck wires to the left as shown.








CAUTION: Be careful working with Fiber Cables. Avoid extreme bends or pinching of cables.


MPI #: 844505-01 REV. L

INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.

13. Bend #1 and #2 FCCs 90° and align to groove closer to the spar. Repeat for #3, #4 and #5 FFCs, but align to groove further from spar.

FFC #s 1 and 2







FFC #s 3, 4, 5



MPI, Carriage Integration, USM, IS4000 MPI #: 844505-01 REV. L 10 of 20


SURGICAL


INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.

ITEM	PART#	DESCRIPTION	QTY
A	318283	CLAMP,FFC,CENTER, IS4000 CARRIAGE	1
В	902974	SCREW,SHCS, M1.6 X3,SST	1
С	722555	EFT, TORQUE WRENCH, 1.6 IN-LBF, 1/4" FEM HEX ADAPTER	NA

14. Install Clamp (A) to Carriage with Screw (B) to secure FFCs using Loctite 222. Torque to 1.6 in-lbs using (C).









MPI, Carriage Integration, USM, IS4000 MPI#: 844505-01 REV. L 11 of 20

SURGICAL

INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.

ITEM	PART#	DESCRIPTION	QTY
A	963004	TAPE,POLYIMIDE (KAPTON),ROLL,.001 X.50W	NA

FFC Connections		
FFC	Connector	
1	J9	
2	J11	
3	J10	
4	J8	
5	J12	


15. Connect **FFCs** as follows:


16. Bend and tape **FFCs** as shown. Secure with **Kapton Tape** (A).




ESD Precautions Required











MPI, Carriage Integration, USM, IS4000 MPI #: 844505-01 REV. L

INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.

17. Connect Fan Assembly to J18: Fan1 connection of ACCS Assembly.

Verify that the **Fan** and **Motor Stator** wires are fully seated in the connectors



ESD Precautions Required




New highlights, picture and text clarifying verification step

on top of the **PCA**. Record result in Result Sheet.



Fan1






MPI, Carriage Integration, USM, IS4000 MPI#: 844505-01 REV. L 13 of 20

INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.


ITEM	PART#	DESCRIPTION	QTY
A	318286	LOWER SUPPORT, CARRIAGE,IS4000	1

19. Install **Fan Assembly** onto **Lower Support** (**A**), as shown below. **Fan** slips into shoulder and secured with **PEM Locating Pins** on the **Lower Support**. Install **Lower Support** into Carriage as shown.



Angle up


Insert first.



Fan Shoulder



Carefully hinge downward



on ACCP PCA.

CAUTION: Be careful when installing
Lower Support. Make sure you
do not bump and damage LEDs



ESD Precautions
Required



MPI, Carriage Integration, USM, IS4000


Front View

MPI #: 844505-01 REV. L 14 of

INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.

ITEM	PART#	DESCRIPTION	QTY
A	902974	SCREW,SHCS, M1.6 X3,SST	4
В	722555	EFT, TORQUE WRENCH, 1.6 IN-LBF, 1/4" FEM HEX ADAPTER	NA

20. Secure Lower Support using Qty: 4 Screws (A) with Loctite 222. Torque to 1.6 in-lbs using (B).








MPI, Carriage Integration, USM, IS4000 MPI#: 844505-01 REV. L 15 of 20

INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.

ITEM	PART#	DESCRIPTION	QTY
A	971125	ADHESIVE,UV CURE GEL,LOCTITE 3979	AR
В	319736	LIGHTPIPE,RIGHT, IS4000,CARRIAGE	1

Glue Flanges

21. Apply a thin bead of **Loctite 3979** (A) to flanged surface on the **Right Light Pipe** (B).





MPI, Carriage Integration, USM, IS4000 MPI#: 844505-01 REV. L 16 of 20

SURGICAL


INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.


ITEM	PART#	DESCRIPTION	QTY
A	319736	LIGHTPIPE,RIGHT, IS4000,CARRIAGE	1
В	319733	COVER, IS4000 CARRIAGE	1
С	723349	AIDE,UVCURE, LOCTITE 7700-HD	NA
D	710349	EFT, TIMER, 1 SEC-99 MIN 59SEC, 30 SEC ALARM	NA
Е	319735	LIGHTPIPE,LEFT, IS4000,CARRIAGE	1

- 22. Place **Right Light Pipe** (**A**) into the **Cover** (**B**) **Housing** from the inside and press into the mating feature on the **Cover** (as shown).
- 23. Using **UV Cure Aide** (**C**) and cure adhesive for ~30 seconds. (adhesive turns pink when cured). **Note: Make sure to wear safety glasses during this step.** Use **Timer** (**D**) during curing process.
- 24. Repeat Steps 21 23 for second Light Pipe (E).



Safety Glasses Required








MPI, Carriage Integration, USM, IS4000 MPI#: 844505-01 REV. L 17 of 20

INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.

ITEM	PART#	DESCRIPTION	QTY
A	319739	FOAM PAD,FAN SUPPORT,IS4000 CARRIAGE	1

- 25. Remove paper to expose adhesive backing on Foam Pad (A).
- 26. Apply the **Foam Pad** while carefully aligning the edges (see arrows) on the **Cover**.
- 27. After placing **Foam Pad** into the well, press around the perimeter to ensure adhesion.





MPI, Carriage Integration, USM, IS4000 MPI #: 844505-01 REV. L 18 of 20

Torque

1.6 in-lbs

INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.

SURGICAL

SCRIPTION	OTY	28. Secure FFC with Kapton Tape (A).

ITEM	PART#	DESCRIPTION	QTY
A	963004	TAPE,POLYIMIDE (KAPTON),ROLL,.001 X.50W	NA
В	902974	SCREW,SHCS, M1.6 X3,SST	4
С	722555	EFT, TORQUE WRENCH, 1.6 IN-LBF, 1/4" FEM HEX ADAPTER	NA


29. Attach Cover to Carriage using Qty: 4 Screws (B) and Loctite 425. Torque to 1.6 in-lbs using (C).

CAUTION: Avoid damaging the cables while installing the Cover.

Apply Kapton Tape over this FFC and pay close attention to avoid damaging when installing the Cover.



Required






30. Connect **Motor Module Cables** according to table.



PITCH

ESD Precautions
Required

Yaw Motor	MTR1
Yaw Hall Sensor	HALL1
Yaw Brake	BRK1
Pitch Motor	MTR2
Pitch Hall Sensor	HALL2
Pitch Brake	BRK2





MPI, Carriage Integration, USM, IS4000 MPI #: 844505-01 REV. L 20 of 20

INTUITIVE SURGICAL CONFIDENTIAL AND PROPRIETARY. Users of any printed hard copy of this document must verify that it represents the most current revision prior to use.


ITEM	PART#	DESCRIPTION	QTY
A	319140	LIGHT PIPE, LINK 1, USM,IS4000	1
В	318349	STRUCTURAL COVER,LINK 1, USM,IS4000	1
С	901195	SCREW,SHCS, M1.6X4,STST	4
D	901212	SCREW,SHCS, M47 X 20,STST	2
Е	902036	SCREW,SHCS, M4 X 12 MM,STST	2

31. Install Light Pipe (A) into Link 1 Cover (B) with Screws (C) and Loctite 425.





32. Install Link 1 Cover onto Link 1 with Screws (D) and (E).



Intuitive Surgical Confidential - This document contains confidential and proprietary information, and except with written permission of Intuitive Surgical, such information shall not be published or disclosed to others or used for any purpose other than on behalf of Intuitive Surgical. This document shall not be copied in whole or in part, including electronically. MPI Template: 800103