## Smart Watch with Wrist-based Heart Rate / SpO2 / Blood Pressure/ Stress / HRV Monitor

1963YH


Smart Watch

### Get Started

What included:

- 1) Smart Watch
- 2) USB Charger
- 3) User Manual
- 4) Gift Box

1. USB charger

2. Home Button

3. PPG Sensor

4. Temperature Sensor

### 1. Download the APP

Go to Apple Store or Android Play Store, please search "JClife" to download our "JClife" application. Or scan the following QR code and jump directly to the APP download interface and enter into APP Stores.

QJClife

Ensure Bluetooth is enabled on your smart phone. Important: Please ensure that your Android OS is version 4.3 or above and iOS version is 7.0 or above and the device has Bluetooth 4.2.

### 2. Activate the Device

Before use, please connect the USB charger to charging the device for activation. The progress bar will change from empty to full to indicate full activation as below. After activation, the device will start charging the battery. Please use the device after it's fully charged.


### 3. Battery Charge

The battery indicator is displayed on the home screen of the device while charging. If the battery level is less than 10%, the device will display the empty battery symbol, prompting to charge the device. The green progress bar will jump during charging, and it will stop until the battery is fully charged.

### 4. Activity Mode

Touch the screen to check detailed daily activity data.

### 5. Heart Rate Monitor

Slide the watch screen to find the Heart Rate icon, press it to start heart rate monitoring. The time interval of automatic HR monitoring can be set via APP.

Please do wear tightly for accurate heart rate monitoring, otherwise the device will gently vibrate to remind you to wear tightly.

### 6. SpO2 Monitoring

Slide the watch screen to find the SpO2 icon, press it to start SpO2 monitoring.


### 7. Breath Training

Slide the watch screen to find the Breath icon, press it to start breath training, the training level & duration can be set via APP.



### 8. Blood Pressure, HRV & Stress Level Monitor

Slide the watch screen to find the Vitals icon, press it to measure blood pressure, HRV & stress level, please keep still while measuring, and after finished, you can slide to check the report.

### 9. Multi-sport Mode

In Multi-sport mode, you can choose Run, Cycling, Badminton, Football, Tennis, Yoga, Dance, Basketball, Hiking, Gym to record your exercise data, such as Exercise Time, Steps, Distance, Speed, Heart Rate, Calories, etc.


### 10. More Functions

Slide the watch screen to find more functions - Find Phone, Remote Camera Control, Music Control, Timer, Stopwatch, SOS Help and Language Exchange.

### 11. Setting

Slide the watch screen to find Setting menu, you can set the time unit, distance unit, wrist on/off, alarm on/off, weather display on/off, silent mode on/off, info.

### 12. Sleep Mode

The device can automatically monitor whether you have gone to sleep and record your sleep quality at the same time. You can also set a sleep alarm on the APP, the device will vibrate on time and wake you up gently.

### 13. Phone GPS Supported

Please click the GPS icon on the App to start the exercise. The APP will start the GPS of the mobile phone and record the speed, distance, calories and activity time. After the exercise starts, the device will automatically turn on the heart rate monitoring. After exiting the exercise, the heart rate will automatically stop monitoring. You can view the exercise data on the APP.


### 14. Phone Notifications

(Only effective when Bluetooth is connected between the device and the mobile APP) After turning on the smart reminder function on the APP, the device will gently remind you not to miss important information such as incoming calls, text messages, WeChat, Facebook, QQ and other social software messages on the mobile phones.

### 15. Personal Reminder

Personal reminders can be set on the APP to develop good lifestyle habits, such as regular alarm clocks, drinking water, taking pills, social distance alert, and other reminders.

### 16. Sedentary Reminder


Set sedentary reminder on the APP, the device will gently vibrate to remind you to exercise more within the set time period.

### 17. Wrist Sense


Alternatively, after setting on the APP, simply raise your wrist to activate the display automatically.

### 18. APP Basic Operation

1. Check the APP home page about this device and the features of this APP.
2. Set personal information such as photo, age, gender, height, weight, left / right hand habits to wear, and sports goals, etc. Connect the JClife APP, select the device you are currently wearing, click next, you can see the personal information interface.

After connected successfully, set your personal information and save it.

Set your smart alarm & notifications to get your personal reminders and then save.

### 19. Personal Goal Setting

Customize your workout plan on the APP, the device will gently vibrate to remind you to exercise more.


1. Tap Goal Setting
2. Choose your target number of steps during your daily target time period
3. Choose the sleep hours you may require everyday
4. Tap confirm

### 20. Automatic Health Monitoring

The device supports real-time heart rate monitoring. You can set the interval for automatic measurement on the APP. The default interval time is 60 minutes, and you can set the interval every 5, 15, 20 minutes and so on.

Note: The shorter the interval for automatic monitoring, the faster the power consumption.

### 21. Check History Data

Clicking the heart rate or SpO2, HRV, Stress icon will display the history data of the corresponding heart rate, SpO2, HRV and Stress. You can also click the history data curve with your finger to view the heart rate or SpO2, HRV, Stress value per minute. At the same time, you can also view the weekly / monthly / yearly data on the APP to understand your average data over a period of time.

### 22. Specifications

Size: Main body, 42\*35\*10.9mm  
Weight: 50g  
Battery: 230mAh rechargeable lithium polymer battery  
Display: 1.3" TFT  
Data Memory: 30 days  
Walking Distance: Max 9999.99km  
Calories Burned: Max 9999.9kcal  
Activity Time: 99h59min

### 23. Precautions for Use

▲ Warning! Children swallowing risk! Children are not allowed to play with the device. Keep children away from the device.

▲ Warning! Risk of fire / explosion or chemical burns! This product contains a lithium-ion battery. Do not remove, disassemble, throw into fire or other battery.

△ Caution! Possible property damage! Do not drop the device and protect it from impact. Do not expose the device to extreme temperatures or extreme temperature fluctuations. Protect the device from direct sunlight and dust. When cleaning, do not use any strong chemical, abrasive or foam cleaning materials. Protect the display from hard objects.

△ Caution! Do not dispose of the device in domestic waste! If you have any questions, please contact your municipal waste disposal authority. Do not dispose of used batteries in domestic waste. Instead, take them to a special waste disposal site or a retailer's battery collection point. Packaging is recyclable or can be recycled into raw materials. Please dispose of unnecessary packing materials properly.

### 24. FCC Statement

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
  - Increase the separation between the equipment and receiver.
  - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
  - Consult the dealer or an experienced radio/TV technician for help
- This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

### IC Statement

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

- (1) This device may not cause interference; and
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

The term "IC" before the certification/registration number only signifies that the Industry Canada technical specifications were met. This product meets the applicable Industry Canada technical specifications.

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