

user manual

Model name:HB-04

FCC ID: 2AD7CHB-04

T-Touch Operation Manual

Specification Data

Measurement Range: 33°C~42°C

Accuracy: $\pm 0.1^{\circ}\text{C}$

resolution 0.01°C

Measurement time : around 6 Mins

Monitoring time: up to 7 days

wireless Syncing: Bluetooth 4.0

Valid transmission distance: 10m

Battery: Polymer lithium battery

Charge Time: 3 hours.

Charging State: red light means charging, blue light means full

Dimension: 67×48×17mm

Operation Instruction:

1. Please connect the temperature probe with the device before turn ON the device. LED light will turn on when the device started. Open your smart phone and connect with the device.
2. Please put the probe under armpit and observe the temperature on mobile after at least 6 minutes.
3. Please turn OFF the device after monitoring to save energy.
4. Please using the USB cable connects power for charging the battery.

Notice:

1. Please make sure temperature cable connected with the device before monitoring.
2. Please don't do strenuous exercises before monitoring.
3. Please clamp the temperature probe under your armpit
4. Please don't use the thermometer in high magnetic environment
5. Please keep this product away from chemical solvent or diluents.
6. Children, younger than 3 years, use the device should under protecting of their guardians.
7. The device is for system IOS7 and Android version 4.3 or above, requires smart phone with BT4.0

FCC Caution.

§ 15.19 Labelling requirements.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

§ 15.21 Changes or modification warning

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

§ 15.105 Information to the user.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio

communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

*** RF warning for Portable device:**

The device has been evaluated to meet general RF exposure requirement.

The device can be used in portable exposure condition without restriction.