**USER MANUAL
&
INSTALLATION INSTRUCTIONS**

**WIRELESS SENSOR SIRES OF
ERGOLAB**

**Designation of type
ErgoLAB**

User Manual

1、 Overview.....	2
2、 Technical Index.....	3
3、 Indicator Light.....	3
4、 Use Steps.....	4
5、 Parameter Setting.....	6
5.1、 Channel Choosing.....	6
5.2、 State Monitoring.....	7
5.3、 Start Recording.....	9
6、 Data Recording.....	10
6.1、 Data Collecting.....	10
6.2、 Interface Presenting.....	11
6.3、 Data Analyzing.....	12
7、 Tips.....	13
Installation Instructions	13

1、 Overview

ErgoLAB wireless sensor adopts RF technology. RF is short for Radio Frequency which is RF current, it is short for electromagnetic wave of high frequency and AC change. The frequency is very high, which can be used for sending radio frequency, generally refers to dozens to hundreds frequency band, that is VHF-UHF frequency band. And Higher frequency is called micro wave. Broadly speaking, micro wave on radio-frequency spectrum refers to radio wave of which the frequency is 300MHz-300GHz, specifically refers to 1-30 frequency band. RF has broad and irreplaceable function in the field of wireless communication. Electric

information source(simulation or numeral) is modulated(amplitude or frequency modulated) with high frequency current to form RF signal and sent to the sky by antenna; RF signal is received in long distance and counter-modulated to restore to electric information source, this process is called wireless transmission.

ErgoLAB wireless sensor adopts 2.4GHz RF technology, the volume is small and easy to wear. It takes samples at 16 bits high-resolution and high speed, the sampling frequency reaches up to 4096Hz.Indoor data tracking is broad and prompt and can suppress noise.

2、 Technical Index

3

www.kingfar.net

Supplying battery	3.7V,200mA
Service time of battery	4-5h
Charging time	1h
Sampling rate	Max: 4096Hz
Working frequency	2402 MHz
Transmission power	0dbm
Operating environment	Temperature:-40°C~45°C Humidity: <60%
size	36mm*22mm*11mm
weight	20g

3、Indicator Light

Green light: after starting the sensor, green light keeps on, when collecting data, green indicator light will flash at the frequency of one second.

Red light: the sensor will self-test, red light will be put out 1s after it lights, when the capacity of battery is lower than 30%, red light will flash at the frequency of 1 second.


Figure 1

4、 Use Steps

It needs to combine with certain software when obtaining the data with this wireless sensor:

5

www.kingfar.net

- 1) Set the switch of wireless sensor to “on”, then green light lights, which indicates the wireless sensor is transmitting data;
- 2) Open ErgoLAB, choose wireless sensor collection module;

Figure 2

- 3) In software, choose “START” to receive the data of wireless sensor.
- 4) In software, choose “STOP” to stop receiving the data of wireless sensor, the data will save automatically.

Referring to the detailed process in the following part 5 and part 6 :

5、 Parameter Setting

We need to set the parameter of sensor before receiving the data of ErgoLAB wireless sensor.

5.1 Channel Selection

Wireless sensor has its own range of radio-frequency, in the software, we need to choose the channel of sensor, then refresh, search this sensor before receiving the data of sensor.

Complete with “channel” of software.


Figure 3

7

5.2 State Monitor

After searching this sensor, it presents relevant information of this sensor, sampling frequency of sensor which needs to detect, battery state and steady state of signal;

Sampling frequency of sensor and current battery state:

Sampling Rate	Battery
64	86%

Figure 4

We can continue with experiment or charge wireless sensor according to battery state.

It is suggested that we should consider

8

charging when capacity of battery is lower than 30%.

Monitoring steady state of sensor's signal:

Figure 5

If following cases were found:

Figure 6

In such case, data continues to be collected, but the collection of signal may interrupt, so it needs to monitor battery capacity of wireless

sensor and pay attention if connection of sensor is unobstructed.

5.3 Start Recording

When state monitoring is completed, it can normally search wireless sensor, and receive the data that sensor sends, we can click “start” to start collecting and recording data of sensor.

Figure 7

When the color of “start” turns from blue to red, it means it is in the process of recording data.

6、 Data Recording

6.1 Collecting Data


Figure 8

This is the interface of collecting data, up side shows collected data signal, down side shows data is being collected. If collection needs to be stopped, just click the red button:


Figure 9

When red button occurs “stop” tip, click it to stop recording data. Data of wireless sensor will not be received. Received data will be automatically saved in the software.

6.2 Presentation Interface

After data record is saved, recorded data of wireless sensor in this period of time can be replayed in software and processed.


Figure 10

12

6.3 Data Analysis

Collected data of wireless sensor can be compiled and analyzed in software.


Figure 11


Figure 12

7、 Tips

If red light of sensor flashes frequently, it indicates that the sensor is collecting data; if red indicator lights, it indicates the battery is low and needs to be charged promptly.

INSTALLATION INSTRUCTIONS

ErgoLAB wireless sensor

Overview

ErgoLAB wireless sensor adopts RF technology, freely recording relevant data information. Appearance size of this sensor is only 43mm*26mm, its volume is small and easy to wear.

Technical index

ErgoLAB wireless sensor adopts 2.4GHz RF technology, working on frequency band of 2.4 GHz ISM,

ErgoLAB wireless sensor adopts 3.7V Li battery to supply power, it takes samples at 16 bits high-resolution and high speed, the sampling frequency reaches up to 4096Hz.

Attention

Warranty period of this product is 2 years. If quality problem is proved, Kingfar group will be responsible for maintaining or replacing. If not, Kingfar will not have this obligation.

User should send written-notice to Kingfar in 30 days once finding quality problem. And

attach the detailed description of quality problems.

Attention before installing

Please make sure the sensor can work normally and correctly collect data before installing on test person. When installing on the skin, lines of sensor should be arranged carefully, lead wires can not be crossed each other, or it will affect the test result.

Warning

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired

operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

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