# Penguin ENT Scope

## 

Photograph Mode: Press briefly to take a snapshot. LED blinks once.

Recording Mode: Hold down the button for 1.5 seconds to start the video recording. LED blinking continuously. Press again to stop recording.

# Brightness

Press briefly to brighten/dim illumination. Hold down the button for 2 seconds to switch MSI modes.







Mode2: Vessel Vessel dark blue in blue LED Red



Mode3; Water Water dark blue in blue LED Blue●

## O Power & Light indication

**Power On:** Press and hold the button for 2 seconds, the LED turns on in last mode color.

Wi-Fi Connecting: The LED blinks slowly.

Wi-Fi Connected: The LED stays on

continuously.

USB Connection: The LED stays on

continuously.

Low Battery Warning: The LED blinks fast. Power Off: Press and hold the button for 2 seconds to turn off the device, the LED turns off.

Charging: The LED blinks fast.

Charged: The LED stays continuously.

Light indicator will remain the same mode color from the last status before power off.

#### Battery

Built-in rechargeable battery.

#### **USB Cable**

2-meter USB Type-C/Type-C cable for charging and UVC video streaming to PC/Pad/Mobile.

#### Runtime

The battery supports 10+ hours of continuous operation per full charge.







# Limited Warranty and Repairment

### 1. Within 1-Year of purchase:

The buyer may return the product to OPCOM for free repair or replacement if it has defects in material, workmanship, or malfunctions within the one-year warranty.

Please note that the warranty does not cover defects or malfunctions from misuse or mishandling.

## 2. Warranty Beyond 1 Year:

If the product experiences issues beyond the initial one-year warranty period, the buyer can return it to OPCOM for repair or replacement, subject to a service charge.

Before initiating a return, please contact our website's customer service to obtain a Return Merchandise Authorization (RMA) number and send the product to the original shipping location. Please note that the cost of returning the product is the buyer's responsibility.

#### FCC STATEMENT

- This device complies with Part 15 of the FCC Rules.
   Operation is subject to the following two conditions:
   (1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesired operation.
- Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna. Increase the separation between the equipment and receiver.

Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

Consult the dealer or an experienced radio/TV technician for help.