

## WARNINGS

- Do not use this device if you have a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device. Doing so could cause electric shock, burns, electrical interference or death.
- Do not use during pregnancy or if you are breastfeeding. If you are in the care of a physician, consult your physician before using this device.
- Do not use this device if you have epilepsy or a history of seizures.
- Do not use this device if you have a Temporomandibular Joint Disorder, Bell's Palsy, impaired cranial nerve function, or facial pain.
- Do not place Thync Strips on body in locations other than those directed.
- Do not use this device while driving, operating machinery, or during any activity in which electrical stimulation can put you at risk for injury.
- Do not use the device while in a shower, bath, pool, or other body of water.
- Do not place Thync Strips over open wounds, sores or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions. If you experience an adverse reaction, discontinue use.
- Do not use device if the housing has been damaged.
- Do not use device in the presence of strong electromagnetic fields.
- Do not use this device on children under the age of 18.
- Do not place this device across your chest. The introduction of an electrical current to the chest may cause rhythm disturbances to your heart, which could be lethal.
- Do not place device over the carotid sinus nerves, the front of the neck, or around the mouth.

# Thync

Shift Your State of Mind

# Thync

For more information, visit [thync.com](https://thync.com)

[thync.com](https://thync.com)

## ADVERSE REACTIONS

If you experience adverse reactions, stop using the system and consult with your physician. Users with sensitive skin may experience skin irritation in the area where the Thync Strip is applied. You may experience a headache and other painful sensations during or following the application of electrical stimulation.

## PRECAUTIONS

- Use caution if Thync Strip is placed over areas of skin that lack normal sensation.
- Clean and dry the area of skin to which the Thync Strip will be applied before applying it.
- Use of accessories not approved by the manufacturer may cause harm or injury.
- Do not disassemble the Module.
- Thync Strips should only be applied to normal, intact, healthy skin over locations as directed.
- Keep dry. Clean by wiping with dry cloth. Do not immerse.
- Operation in close proximity to short wave or microwave therapy equipment may produce instability in the Module output.
- Keep this system out of reach of children.

## LIMITED WARRANTY

Thync, Inc. warrants this Thync hardware product against material defects in materials and workmanship for a period of one year from the date of purchase by the original purchaser. This express written warranty is exclusive and in lieu of any other warranty or remedy, express or implied. For elaboration on this warranty and related claims and remedies, please go to our website at [thync.com](http://thync.com)

## ENVIRONMENTAL AND TECHNICAL

Store Between	-40 and 70 degrees C 10 – 90% Relative Humidity Atmospheric Pressure 50 – 106 kPa
Module Weight	17.9 g
Power Supply	x V Battery type model
Channels	1 channel output through female snaps
Output Current	0-23 mA peak into 500-ohm load
Average Output Current	< 10 mA into 500-ohm load
Intensity Control	Adjustable in 100 steps
Average Current Density	< 2 mA/cm^2
Waveform	Asymmetrical Biphasic Pulse
Operating Temperature	0 to 40 degrees C

## IC RSS-Gen

This device complies with Industry Canada's licence-exempt RSSs. Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

Le présent appareil est conforme aux CNR d'Industrie Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes : (1) l'appareil ne doit pas produire de brouillage, et (2) l'utilisateur de l'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement.

## FCC STATEMENT

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

## FCC NOTICE

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:  
(1) This device may not cause harmful interference, and  
(2) This device must accept any interference received, including interference that may cause undesired operation.

## FCC CAUTION

Changes or modifications not expressly approved by Thync, Inc., could void the user's authority to operate the equipment.