











User Manual

FCC ID: 2AGFWB80

Model Name: B80

1. Usage and requirements

- A : Long press bracelet's touch screen to power on.
- B : Touch screen to switch bracelet's interfaces.

 <p>Time interface: Show the current time.</p>	 <p>Pedometer interface:Show the current steps.</p>
 <p>Calories interface: Show the current calories.</p>	 <p>Distance page:Show the current sport distance.</p>
 <p>Blood pressure measurement interface: Stay in this interface for 1 second and green light on the back of bracelet will light, about 10 seconds later, bracelet's shocking remind measurement to end.</p>	 <p>Finding interface:When bracelet is connected with mobile phone,long press touch key for 3 seconds, mobile phone will ring.</p>
 <p>Power off page:Long press for 3 seconds,then it will show "Yes", and then click it for 3 seconds more,it for 3 seconds more,it will show "Byebye" to power off,long press to power on.</p>	 <p>Heart rate:Stay in this interface for 1 second,it will measure heart rate and light on the back of bracelet will light, measurement will stop about 10 second later.</p>
 <p>stopwatch interface:Press screen 3 seconds and wake up stopwatch function.Turn on device with short press,and pause with short press again, exit with long press.</p>	 <p>Brightness : Long press to brightness level setting mode, short press to set up different level and long press again to exit.</p>

Precautions :

- ※ Can wear this fitness tracker to wash hands, face,short time swimming in shallow water, can't touch hot water and for deep diving.
- ※ Don't throw and collide which would easy to crash r damage the device.
- ※ Don't disassemble the device by yourself which would cause it doesn't work.
- ※ Don't throw it in fire, easy to explode!
- ※ Please use it in -5°C to 45°C.
- ※ Use soft cotton flannel or glasses cloth to clean the screen, don't use soap shower gel or other detergents to clean.
- ※ This product can't use as disease diagnosis, treatment and prevention. Just for reference only.
- ※ Don't let kids or pets to touch this device in case of damage.

1.1 Time Page

Switch time interface, As following pictures:



Long press to select different watch-face in watch interface, short press touch screen to switch and long press to confirm exit. (After your cellphone connected with the APP, time will be sync with your smartphone's time.)


1.2 APP download/Connect the bracelet

Compatible for mobile phones : Android system 4.4 version or iOS system 9.0 version or above ;
Compatible for blue-tooth version : Bluetooth 4.0 version above.
 Use mobile phone scan QR code on users' manual or search "Hero band II " APP by APP store or Google Play to download and install it .




2.Bracelet connection

Before using bracelet ,users need register APP on mobile phone firstly .Chinese users support using mobile phone numbers to register ,foreign and Hong Kong, Macao and Taiwan regions support using e-mail to log in. Open APP, enter into "Mine" and click "connect a bracelet" to search a bracelet named"Hero band II " , keep mobile phone 's blue-tooth on when they are connected .



2.1 Notification settings

When finished connection, users enter into 【Mine】 by APP, click 【Device Management】 , and Android mobile phones to enter into this function firstly, APP will give a tips for allowing it could read mobile phone's notification permission, click confirm, then on 【Notification access】 users need turn "Hero band II " on; When iOS mobile phone connected with bracelets, it will give a tip for gairing, please allow it, otherwise bracelet will not receive notification. **As below:**



2.2 Lift wrsit to light up screen

Enter into 【Device Management】 turn on "Lift the wrist brighten screen" or "Turn the wrist brighten screen":

2.3 Message Reminder

Bracelets support message notification showing screen of Skype/Whatsapp/Facebook/We-chat/QQ/SMS and so on.Long press bracelet screen for 3 second to back into main interface when bracelets received messages.

2.4 Reminder settings

Turn sedentary / Take medicine / Drink / Meetings remind on and Set their working time as user's need.

2.5 Firmware upgrade

When mobile phone was connected with bracelet, click "Firmware upgrade": If it upgrade failed, users need charge the bracelet and try it again.

3.1 APP Sport / Health / Care

Sport / Health / Care three parts:




3.2 Sport

On this function,there two data of sport/sleep.

3.2.1 Sport (Steps / Distance / Calories)

As following pictures:



3.3.1 HR-Heart rate monitor

Normal heart rate value is 60-100 times/min,there will be a little difference between people with different age,gender and other factors.

3.3.2 BD-Blood pressure monitor

Normal blood pressure value:Systolic blood pressure 90-139mmHg, diastolic blood pressure 60-89mmHg.

3.4 Care/Remote view

Search friends and add them to check or share the sport data with them.

3.4.1 HR-Heart rate monitor

Normal heart rate value is 60-100 times/min,there will be a little difference between people with different age,gender and other factors.

3.4.2 BD-Blood pressure monitor

Normal blood pressure value:Systolic blood pressure 90-139mmHg, diastolic blood pressure 60-89mmHg.

3.4.3 Care/Remote view

Search friends and add them to check or share the sport data with them.

3.4.4 HR-Heart rate monitor

Normal heart rate value is 60-100 times/min,there will be a little difference between people with different age,gender and other factors.

3.4.5 BD-Blood pressure monitor

Normal blood pressure value:Systolic blood pressure 90-139mmHg, diastolic blood pressure 60-89mmHg.

3.4.6 Care/Remote view

Search friends and add them to check or share the sport data with them.

FCC Caution.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.