



User Manual
Bluetooth Thermometer
T101



1、 Upload battery



2、Download APP

- 1、Go to APP Store
Search iLece



3、Build connection

Follow operation insturction below



4、Prepare protect film



5、Prepare cover film



6、Wearing(iLece Wearing Method)



1 Please take one of sleeve off, and stick the iLece fever alert under armpits.



2. Put down the dress without wearing the sleeve.



iLece Fever Alert

Note :1 Adult should help child to wear and operate the device
2 Please holding your arm tight at least for 5 mins.
3. Put the mobile phone and the iLece fever alert on the same side closely.

Privacy Policy &Disclaimer

We use the information you provide to complete our production and service. We will not share your personal information we collected with unrelated service providers without your permission.

iLece Fever Alert is a kind of home-consumer electronic equipment,assist in the understanding of body temperature, environment temperature etc. iLece Fever alert is full complied with national standards. We do not recommend use in clinic diagnosis. The iLece Fever alertt is not a substitute for medical device. iLece IOT Tech company holds no liability for damages because of hardware,APP,mobile phone failure,bug,or improper use.

ACCESSORIES

1. Two CR2016 battery
2. Medical sticker (1 box)
3. User manual

TECHNICAL SPECIFICATIONS:

Model :	T 100	Battery :	CR2016
size :	34*34*4.2mm	Battery work time :	150 days
Weight :	5g (without battery)	Temperature :	-15—50℃
Measurement range :	0℃—50℃	Humidity :	10%—85%
Accuracy :	± 0.05 (35.00℃-38.50℃) ± 0.1 (25.00℃-34.99℃) (38.51℃-42.00℃) ± 0.5 (00.00℃-24.99℃) (42.01℃-50.00℃)		
Temperature unit :	℃ 或 °F		

FCC Statement

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the

instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation