



TBSI W16 System Manual v.1.0

Triangle BioSystems International

2224 Page Rd. Suite 108

Durham, NC 27703

Phone: (919) 361-2663 • Fax: (919) 544-3061

www.trianglebiosystems.com

Contains/Contient FCC: 2AJQ2W11; IC: 22071-W11

Notice: FCC 15.105(a) and 15.21: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense. The user shall not make modifications or changes not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesired operation.

Avis: FCC (Federal Communications Commission) 15.105(a) et 15.21: ce matériel a été testé et est soumis aux restrictions d'utilisation d'un appareil numérique de classe A conformément à la partie 15 de la réglementation de la FCC. Ces restrictions ont été conçues pour assurer une protection satisfaisante concernant les interférences nuisibles produites par l'appareil dans le cadre d'une utilisation commerciale. Dans son fonctionnement, cet équipement génère une fréquence radio qui, dans le cas d'une installation et/ou d'une utilisation non-conforme avec les instructions de la notice explicative, peuvent causer des interférences néfastes avec d'autres fréquences radio. Toute utilisation de cet équipement dans une zone résidentielle risque fortement de causer des interférences nocives qui, le cas échéant, devront être rectifiées à la charge de l'utilisateur. Il est fortement déconseillé de procéder à des modifications ou changements qui ne sont pas explicitement approuvés par une autorité compétente en charge de sa conformité. Dans le cas d'une modification non-autorisée, l'utilisateur pourra voir son autorisation d'utilisation révoquée. Cet appareil est conforme à la partie 15 des règles de la FCC. Son utilisation est soumise aux deux conditions suivantes: (1) Cet appareil ne doit pas causer d'interférences nuisibles, et (2) Cet appareil doit accepter toute interférence reçue, y compris les interférences pouvant entraîner un fonctionnement indésirable.

Contents

16 CHANNEL WIRELESS NEURAL HEADSTAGE SYSTEM	4
Introduction	4
Features	4
System Configuration.....	3
System Block Diagram.....	3
Wireless System Parts List	4
Headstage Transmitter Specifications	5
Headstage Mechanical Overview	6
Wireless Headstage Pinout	6
RF Receiver Specifications	7
RF Receiver Mechanical Overview.....	8
RF Receiver Signal Interface Options.....	8
Receiver Signals - Analog Output.....	9
Receiver Signals– Digital Interface.....	9
System Setup and Testing	10
Internal and External Batteries	13
Mouse and Rat Saddle Harness	14
Headstage Accessory Options.....	14
LED Headstage Options.....	15
Electrodes	15
Gain and Phase Response	15
Troubleshooting.....	16
Application Notes.....	16
Contact Us.....	17

16 CHANNEL WIRELESS NEURAL HEADSTAGE SYSTEM



Figure 1: System Setup

Introduction

The Triangle BioSystems International (TBSI) W16 wireless recording headstage system allows researchers to continuously and simultaneously monitor up to 15 single-ended biopotentials as single units, EEG and ECOG. Differential signals such as an EMG recording can easily be accommodated in software. No longer do experiments have to be constrained by tethering a test subject. The complete system is comprised of a wireless headstage transmitter with integrated battery, RF signal receiver/baseband demodulator, power supply and all required cables. With an effective range of 1 meter, this system provides a wireless connection between the implanted electrodes and the data recording system.

TBSI's custom ASIC technology and proprietary radio design provide up to 7kHz bandwidth in a wireless headstage that is both small and light weight (4 g). This design also incorporates preamplifier circuitry to create an extremely compact headstage.

Features

- Wireless operations across 1 meter
- 15 single ended recording channels
- Headstage transmitter weight: 4.0 grams
- Default Total System voltage gain of 800, other gain offerings available
- Rechargeable battery with 3.5 – 4.2 hours of battery life
- Bandpass filtering per channel at 0.8 Hz to 7 kHz typical
- Covered or dipped headstage options



Figure 2: Headstage



Figure 3: Receiver with Dual Antenna

System Configuration

The receiver antennas must be clipped to a surface within 1 meter from the animal at all times to achieve optimal signal quality. They should be positioned above the animal's cage. Be careful not to obstruct the line-of-sight path between the animal and the receiver antenna with any material except for glass or plastic.

See the section entitled [System Setup and Testing](#) for more details on proper system assembly and positioning.



Figure 4: Standard TBSI Wireless Recording Setup

System Block Diagram

The wireless neural headstage system consists of a wireless transmitter headstage, an RF receiver and baseband demodulation subsystems as shown below:




Figure 5: System Block Diagram

Wireless System Parts List



Wireless System Parts List		
Item #	Part Number	Description
1	Various	Headstage
2	Various	Wireless Radio Receiver (Analog or Digital Output)
3	110-0001-00	W-Series Receiver Antenna
4	300-0002-10	Receiver Antenna Extension Cable and Mounting Clamp, 5 ft
5	100-0000-10	Receiver 6 VDC power supply
6	100-00001-00	Headstage Battery Charger
7	200-0021-10	W16 Signal Input Test Cable, All Channels Combined, Single GND
8	300-0000-00	Magnetic Wand, Headstage On/Off, 12 In
9	202-0001-10	Male DB 37 Connector and housing
10	200-0011-00	4 foot earth GND cable for receiver banana jack (Not Pictured)

Headstage Transmitter Specifications

PARAMETER	MIN	TYP	MAX	UNITS	NOTES
Power					
Current Consumption		13		mA	
Battery Capacity		60		MAH	Default version, other versions available
Analog Input Specs					
Input Voltage Range			4	mVp-p	
Gain Selection	790	800	810		Default system gain. Other gains available and is set by manufacturer.
Bandwidth	0.8		8500	Hz	Bandwidth at -3 dB output signal level
Input Impedance		6.5		MΩ	At 1kHz
Total Noise		3.8		μVrms	25 μVpp, Input referred, 1 Hz – 8.5 kHz band
Analog Sampling Rate		50.78125		kHz	Per channel
Mechanical Specs					
Length		22.2		mm	Edge to Edge (including connectors)
Width		16.5		mm	Edge to Edge
Height		14.2		mm	Edge to Edge
Weight		4.0	4.2	g	With connector and dipped package
Input connector					18 Pin Omnetics, .025"
RF Transmission Specs					
Center frequency		3.05		GHz	+/- 25 MHz bandwidth
Transmit range		1		m	Typ., straight-line, unobstructed

Headstage Mechanical Overview





Figure 6: Mechanical Overview

Wireless Headstage Pinout



Figure 7: .025" Omnetics Female
TBSI P/N: A9094
Mating Male Electrode Connector
TBSI P/N: A8784

Pin #	Pin Names and Connection
1	Channel 15
2	Channel 13
3	Channel 11
4	Channel 9
5	Channel 7
6	Channel 5
7	Channel 3
8	Channel 1
9	ACgnd
10	Channel 14
11	Channel 12
12	Channel 10
13	Channel 8
14	Channel 6
15	Channel 4
16	Channel 2
17	NC
18	ACgnd



**Figure 8: .050" Omnetics Female
TBSI P/N: A11862
Mating Male Electrode Connector
TBSI P/N: A11449**

Pin #	Pin Names and Connection
1	Channel 15
2	Channel 13
3	Channel 11
4	Channel 9
5	Channel 7
6	Channel 5
7	Channel 3
8	Channel 1
9	NC
10	ACgnd
11	NC
12	Channel 14
13	Channel 12
14	Channel 10
15	Channel 8
16	Channel 6
17	Channel 4
18	Channel 2
19	NC
20	ACgnd

RF Receiver Specifications

- 1 meter typ. range between headstage transmitter and receiver antennae
- Front-end gain: 60 dB
- Intermediate gain: 10-20 dB
- Analog channel bandwidth: 8.5 kHz
- DC offset: < 100 μ Vdc
- Phase delay typical: 50 μ sec at 10 kHz
- Signal lock indicator: LED on front panel

RF Receiver Mechanical Overview




Figure 9: Receiver Mechanical Overview

RF Receiver Signal Interface Options

The W16 RF Receiver Base Station is available in different system configurations to accommodate a range of interface requirements. The RF signal from the wireless headstage is demodulated and processed by the RF Receiver Base Station. Each headstage neural channel is directly fed as an analog output signal to one DB37 connector. The analog neural channels are interfaced to the user's neural data acquisition system, where the analog signals are converted to digital signals, fed into a PC system, and analyzed by the user's neural software.

Receiver Signals - Analog Output




Figure 10: Receiver Box Diagram

Pin #	Description
20	Channel 1 Output
21	Channel 2 Output
22	Channel 3 Output
23	Channel 4 Output
24	Channel 5 Output
25	Channel 6 Output
26	Channel 7 Output
27	Channel 8 Output
28	Channel 9 Output
29	Channel 10 Output
30	Channel 11 Output
31	Channel 12 Output
32	Channel 13 Output
33	Channel 14 Output
34	Channel 15 Output
1,36	GND
2-19, 35, 37	No Connection

Receiver Signals- Digital Interface

TBSI offers an internal DAQ board option and [NeuroWare™](#) software for the 5ch, 16ch, 32ch and 64ch receivers which can be used for signal data acquisition and analysis. Digital output and event input information is sent to the PC via a USB connection. The recorded data file can also be read by MatLab® script or NeuroExplorer®.




Figure 11: NeuroWare Software Image

System Setup and Testing

The steps below should be followed when setting up your W16 system to optimize its performance.

1) Charging the Headstage Lithium Ion Battery

The wireless headstage includes an integrated rechargeable battery which should be recharged when not in use. Recharging is accomplished by the following procedure:

- Turn off the headstage transmitter using the magnetic wand (blue light will turn off as indication).
- Connect the battery charging adapter to a wall socket.
- The white charging plug is keyed and mates with the white socket on the side of the headstage. Connect the charging cable to the headstage, making sure that the connection is sound.
- The indicator light on the charging unit is red while charging and green when the connected battery is halfway recharged. A dead battery will become fully charged within approximately one and a half hours after the light has turned green.



Figure 12: Headstage and Charger

2) Correctly position the receiver near the animal's cage:

It is critical for transmission performance to correctly position the receiver next to or inside the animal's cage. First, screw the antennas onto both antenna clamps gold input SMA connectors located on the side of each clamp. Attach the SMA connector of the antenna cable to the mating SMA connector on the receiver box. For best results, position the antennae at 45 degree angles with respect to the cage area and aim them towards center of the cage.




Figure 13: Antennae Positions Image

3) Connect power supply to receiver

Plug the DC connector of the AC-DC power supply adapter into the back of receiver and plug the adapter into an AC outlet. The receiver system is powered by an AC-DC adapter that plugs into a US power outlet, 100-240 VAC, 47-63 Hz and provides 6 VDC, 2.5A. International adapters are available to convert the standard US plug to your AC outlet. **WARNING: DO NOT USE ANY OTHER AC-DC POWER SUPPLY ADAPTER OTHER THAN THE TBSI ADAPTER.** Using any other power adapter will void the warranty.



Figure 14: Connecting the Power Supply to the Receiver

4) Connect signal cable to receiver analog output DB37 connector

If your system was configured for analog system output, connect a suitable analog output cable to the DB37 receiver connector. Five signal wires and a ground wire can be added to the DB37 mating connector to check for signal output. The lengths of the wires are not critical. The analog output channel positions are described in the section entitled [Receiver Signal Interface – Analog Output](#).

If your receiver has an integrated DAQ, plug the USB cord into the receiver and PC and follow the instructions provided in the [Neuroware Manual](#) or your data acquisition software for signal acquisition and viewing.




Figure 15: Receiver Plug

5) Turn on Receiver

Flip the on/off button of the receiver to control power to the receiver.

6) Signal Lock LED and measure analog output signals

Once the receiver and headstage (see #7 below) have both been switched on, the **“Signal Lock”** LED on the receiver front panel should light up. After confirming the signal lock integrity, you can view the analog signals (output from the DB37 connector on the back of the receiver) with an oscilloscope. Please note the default system gain is 800, therefore you can expect the analog output values to be about 3.2 Vp-p of the signal from the function generator.

7) Test headstage signal transmission with the W16 Signal Input Test Cable











Figure 16: Headstage and Signal Transmission Test Cable

- Line up the two guide pins on the headstage with the guide pin holes on the test cable. Make sure that the labels on the connectors face the same direction.
- W16 Signal Input Test Cable (200-0021-10) supplied with receiver
- Signal leads – connect to positive side of function generator
- GND Leads – connected to signal Ground side of function generator: - Input a 1 kHz sine wave signal at maximum 4 mVp-p amplitude. The function generator signal output voltage cannot exceed 4 mVp-p, otherwise the headstage input voltage range will saturate.
- Once the test cable has been attached, use the magnetic wand to turn on the headstage. Briefly hold the black tip of the wand next to the top of the headstage until you see the blue “power on” LED. The headstage will transmit data continuously while this LED is on.

Internal and External Batteries

The TBSI W16 recording headstage is typically configured with a 60 mAH internal battery. The W16 headstage can also be built with an internal 75 mAH internal battery for increased recording capacity. Alternatively, the headstage can be configured for use with an external battery. Several different external batteries are available and are listed in the table below. For rodent applications we offer both a mouse and a rat saddle harness with Velcro pad (See the section entitled [Mouse and Rat Saddle Harness](#) below).

Removing the internal battery will reduce the weight and size of the headstage by an amount comparable to the battery specifications as listed in the following table.

TBSI W-Series External Batteries				
				
Capacity	10 mAH	40 mAH	60 mAH	75 mAH
Weight	0.6 g	1.5 g	1.6 g	4.4 g
Size (L xW x H)	12 x 11 x 3 mm	15.3 x 15 x 6 mm	16.2 x 11 x 6 mm	25.2 x 20 x 4.7 mm
Hours¹	0.8	3.1	4.6	5.8
Product Number	109-0001-10	109-0002-10	109-0003-10	109-0004-10
				
Capacity	180 mAH	200 mAH	220 mAH	500 mAH
Weight	5.3 g	4.7 g	5.3 g	11.6 g
Size (L xW x H)	28 x 20 x 4.5 mm	34 x 11.6 x 6.2 mm	30 x 25 x 3.8 mm	40 x 29.8 x 4.6 mm
Hours¹	13.9	15.4	16.9	38.5
Product Number	109-0005-10	109-0012-10	109-0013-10	109-0025-10

1. Approx. hours of continuous operation on a full-charge battery.



Figure 17: 100-0001-04
External Battery Charger 2 prong
male connector



Figure 18: External Battery and
Headstage Connector



Figure 19: International Plug Adapters

Mouse and Rat Saddle Harness

Saddle harnesses can be used for external battery options:



Figure 20: Mouse Saddle Harness
004-0001-10



Figure 21: Rat Saddle Harness



Figure 22: Rat Saddle Harness with Mounted 500 mAh External Battery

Headstage Accessory Options

Three biosensor options are available as add-ons for any wireless headstage. These accessories replace channels typically used for neural recording on the headstage. Installation of any of these features increases the headstage length by 0.4 inches, its height by 0.2 inches, and its weight by 0.4 grams. The features can be installed individually or together in any combination. See the [TBSI Biosensors Brochure](#) for more information.

Accelerometer

This add-on will monitor x, y, and z acceleration vectors and output the information via three analog signal channels. The animal head orientation or velocity can later be calculated with software such as Matlab.

Ultrasonic Microphone

This single ultrasonic microphone enables recording of audio frequencies between 1 kHz and 25 kHz (W5 and W32 Headstages) or 1-40 kHz (W16 and W64 Headstages) and uses one recording channel.

Temperature Sensor

A Thermistor temperature sensor option is available. Many sensor sizes are available including very small ones that can be implanted within or near the brain. The sensor must be attached to a mating connector and externalized to the head of the animal much like traditional recording electrodes. Sensors are customized to your particular research needs.




Figure 23: Headstage Sensors and Dimensions

LED Headstage Options

Red, blue or green permanently affixed LEDs or socket mounted LEDs are available for the W16 headstage. The LEDs are placed facing upward on top of the headstage and are suitable for use with most video tracking software applications. The LEDs will turn on when the headstage is turned on. Installation of LEDs on the headstage may reduce battery life by as much as 15%.



Figure 24: Headstage

Electrodes

Many electrode manufacturers offer electrode arrays compatible with TBSI wireless headstages. [Cambridge Neurotech](#) offers an extensive line of silicone probes suitable for chronic signal recordings from our wireless headstages. Compatible electrode arrays are also available from Microprobes, Inc., NeuroNexus and others electrode providers.

Gain and Phase Response

Please refer to the section of this document entitled "[Headstage Transmitter Specifications](#)" for numerical data on bandwidth and input referred noise.




Figure 26: Gain and Phase Graphs

Troubleshooting

Problem: No neural signals are visible on any of the analog outputs at the DB37 connector.

Suggestion: Verify the AC power connection is in place and the green "Power" LED is illuminated on the front of the RF receiver box. Also, verify that the yellow "Signal Lock" light is illuminated, which confirms that the receiver is receiving the signal transmitted from the headstage.

Problem: Visible neural signal is missing information.

Suggestion: Keep the animal within the 1 meter range of the receiver. If you exceed this range, the radio signal from the headstage will not be strong enough to maintain reliable signal monitoring of the animal. Also, be sure to keep the area under the RF receiver unit's antennas free from metallic objects, which will reduce signal range and introduce noise.

Problem: Not all channels are visible on the neural signal.

Suggestion: Make sure the headstage connection to the animal is secure.

Application Notes

- 1) The wireless headstage is a low-power device, it is critical that the RF receiver be carefully located for the system to operate. DO NOT place the receiver on the outside of a metal wire cage. Try to minimize the distance between the receiver and the wireless headstage. Please refer to the section in this document entitled [System Setup and Testing](#) for orientation suggestions.
- 2) If you do not intend to record from all of the available signal channels you must AC-ground any unused channels at the electrode connector or EIB interface. Failure to do so can create artifact noise from the floating channels on the other signal channels.




Figure 27: Headstage Tie-Down Mechanism

Contact Us

Triangle BioSystems International
2224 Page Rd. Suite 108
Durham, NC 27703
Phone: (919) 361-2663 • Fax: (919) 544-3061
www.trianglebiosystems.com

Need more help? You can reach TBSI customer support by phone at the number listed above (M-F Eastern US Time) or email us at support@trianglebiosystems.com. If you need replacement parts or accessories for your system or to learn the latest about available TBSI products visit our website, call or email us at sales@trianglebiosystems.com.

CAUTION

This product is not registered with the FDA and is not for clinical use on human or veterinary patients. It is intended for research use only.



Contains/Contient FCC: 2AJQ2W11; IC: 22071-W11

Notice: FCC 15.105(a) and 15.21: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense. The user shall not make modifications or changes not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesired operation.

Avis: FCC (Federal Communications Commission) 15.105(a) et 15.21: ce matériel a été testé et est soumis aux restrictions d'utilisation d'un appareil numérique de classe A conformément à la partie 15 de la réglementation de la FCC. Ces restrictions ont été conçues pour assurer une protection satisfaisante concernant les interférences nuisibles produites par l'appareil dans le cadre d'une utilisation commerciale. Dans son fonctionnement, cet équipement génère une fréquence radio qui, dans le cas d'une installation et/ou d'une utilisation non-conforme avec les instructions de la notice explicative, peuvent causer des interférences néfastes avec d'autres fréquences radio. Toute utilisation de cet équipement dans une zone résidentielle risque fortement de causer des interférences nocives qui, le cas échéant, devront être rectifiées à la charge de l'utilisateur. Il est fortement déconseillé de procéder à des modifications ou changements qui ne sont pas explicitement approuvés par une autorité compétente en charge de sa conformité. Dans le cas d'une modification non-autorisée, l'utilisateur pourra voir son autorisation d'utilisation révoquée. Cet appareil est conforme à la partie 15 des règles de la FCC. Son utilisation est soumise aux deux conditions suivantes: (1) Cet appareil ne doit pas causer d'interférences nuisibles, et (2) Cet appareil doit accepter toute interférence reçue, y compris les interférences pouvant entraîner un fonctionnement indésirable.