

User Manual

WIFI Otoscope

Scan the following QR code to install APP:



IOS

Android

Specification:

- 1.WIFI transmission, compatible with Android and IOS systems
- 2.Camera resolution: 1280*720(720P)
- 3.6LED
- 4.600mA battery
- 5.One button photograph

Operation:

1. press the power button about 2S for a long time to start the product. When the power starts, the red wireless signal indicator lights up after 0.5s, and then automatically extinguishes until the WIFI signal is working normally, then it stays bright (no data transmission) or continues to blink (data transmission), and the product is initialized.
2. turn on the mobile phone WIFI switch, search the hot signal, connect the WIFI hotspot named PLX_Camera_xxxx (XXXX is the 4 bit random code), enter the default password: 12345678, after the connection completes, open the mobile phone software Mo-View, enter the view interface.
3. click the middle part of the screen to open / close the control icon. The function of the icon is to take pictures, video, playback, and settings in turn. Click the camera will automatically save the current picture to the mobile phone; click on the video recording video, click again to end video recording, video is automatically saved to the mobile phone; click playback can view the stored photos and videos; click on the settings you can modify the hotspot name and password resolution etc.. Hot names and passwords need to be restarted to take effect.
- 4., the blue indicator is the battery indicator lamp. The product shared 4 lights to indicate the electricity. 1 lights indicate the remaining electricity is 25%, 2 lights indicate the remaining electricity is 50%, and the analogous value is 75%, 100%.
5. long according to the power button about 2s to red wireless signal lights out, product shutdown. The blue power indicator will be extinguished automatically after about 15s.

Warning:

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

NOTE: This device and its antenna(s) must not be co-located or operation in conjunction with any other antenna or transmitter

RF Exposure Statement

To maintain compliance with FCC's RF Exposure guidelines, This equipment should be installed and operated with minimum distance of 5mm the radiator your body. This device and its antenna(s) must not be co-located or operation in conjunction with any other antenna or transmitter