

In the box



N60 headset



Bone wing



User manual



USB
C-Type
charging cable

Specification

Version	Bluetooth 5.0
Wireless Range	10m
Chip Set	CSR8670
Profile	HSP, HFP, A2DP and AVRCP
Weight	50g
Speaker Type	Stereo Bone Conduction Transducer
Frequency Response	20Hz ~ 20KHz
Speaker Sensitivity	100±3dB
Mic Sensitivity	Omni Directional (-42dB±3dB)
Power Supply	Lithium Polymer Battery 210mAH
Charger	Micro USB C-Type 5V
Hearing Aid Usable Time	12 Hour
Music Usable Time	8 Hour
Waiting Time	100 Hour
Charging Time	3 Hour

Buttons and connectors



- 1 Speaker
- 2 LED
- 3 Volume Button
- 4 MIC
- 5 Power/Call/Music
- 6 USB Port (charge)
- 7 Adjusting Band
- 8 Hearing Aid MIC

ETEREO

You must fully charge the headset before using the device for the first time.

- ▶ Connect the USB charging cable to the USB charge port of the headset.
- ▶ The full charge takes about 3hours.

▶ **Power on**

- Press the power button about 5seconds to turn on the headset.

▶ **Power off**

- The full charge takes about 3hours.

While the headset is on,
press the power button to switch between
Bluetooth and hearing aid function.



Pair to your device

1. Put on the headset
2. With power off, press the power button for 5 seconds to enter pairing mode
3. Turn on the Bluetooth function of your Bluetooth device.
4. Select “ETEREO N60” on the Bluetooth device list.

▶ Answering a call

Long press the power button to answer a call.

▶ Music play

When your Bluetooth device is playing music, you can control volume and choose song.

- Next song : Press the volume button(+) for 3 seconds
- Previous song : Press the volume button(-) for 3 seconds
- Adjusting volume : Press the volume button(+)/(-).

"ETEREO N60" APP

- ▶ Smartphone app with many unique functions, to fully utilize Etereo N60.
 - ▶ Enter the serial number on the of the user manual when join (If Serial No. is not entered, there may be a limit to the service use.)
 - ▶ 'Etereo' app features
 - Self frequency fitting
 - Self hearing test
 - Brainwave soothing sound
 - TTS scheduler
 - Simultaneous interpretation
 - Tinnitus Care
- (Details can be found in the app.)



APP (Android/IOS)



FCC compliance information

THIS DEVICE COMPLIES WITH PART 15 OF THE FCC RULES.

Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the

FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the

instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not

occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be

determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the

following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Modifications not expressly approved by the manufacturer could void the user's authority to operate the equipment under FCC

rules.

Quality Assurance

Be sure to bring the warranty and receipt when exchanging or repairing the product.

The warranty period is one year.

Products with defects at the factory will be exchanged free of charge within 10 days with receipt, and defects caused by reasons attributable to consumers will be charged.

For more information, please contact your place of purchase.