



SIREN 4G HUB QUICK START GUIDE

4G Hub V 1.0
Hub Guide V3.0

QUESTIONS?

Contact Siren Support with any questions.

 (415) 769 -3987

 support@siren.care

What's in the box?



Siren Hub

The Siren Hub collects information from your Siren Socks and sends the data to your doctor. Connect Side A to the Siren Hub.

Side A




Side B



Micro USB Cable

This cable connects the Siren Hub to the power adapter.




Power Adapter

Plug the power adapter into an outlet to power the Siren Hub. Connect Side B to the power adapter.

1

Connect Side A


Connect Side A to the Siren Hub.



2

Connect Side B


Connect Side B to the power adapter.



3

Plug in the Siren Hub

Find an outlet in your home where you spend most of your time and plug in the Siren Hub. The Siren Hub light will be solid green when powering on.



4

Siren Hub Successfully Paired

The Siren Hub light will blink blue when successfully connected to the Siren system.



FREQUENTLY ASKED QUESTIONS

Q1. Does the Siren Hub light turn off at night?

Yes, the Siren Hub light is programmed to turn off at night between the hours of **10pm to 8am** local time.

Don't see the answer to your question?

For more information visit <https://siren.care/support/>

Q2. What do the Siren Hub colors mean?

- Solid Green - Powering on
- Solid Blue - Onboarding required
- Blinking Blue - Waiting for socks to connect
- Blinking Orange - Network connection issue
- Blinking Green - Sock(s) are connected

Tip: If the Siren App is open, your socks will disconnect from the Siren Hub and connect directly to the Siren app.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

"This equipment complies with RF radiation exposure limits set forth for an uncontrolled environment.

This equipment should be installed and operated with a minimum distance of 20 centimeters between the radiator and your body."

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.