


High-definition Visual Otscope With LUCKIMAGE MicroEye W1/W2 WIFI box for Operating instructions




Introduction

- This product is an image wireless transmission device. The camera images such as endoscopes and otoscopes are transmitted wireless to the display of the mobile phone, and photos or videos can be saved on the mobile phone.
- There is a button on the front of the device, long press can switch machine.

Instructions

- 1.Scan the QR code to download app Anyscope , or search Anyscope in google play(Android) and App store (IOS) to download the app.
- 2.Insert the camera into the corresponding interface of the device and long press the power button until it turn on.
- 3.Mobile phone connected to will hot spot:Anyscope-xxxxxx.
- 4.Open Anyscope app,You can see the image from the camera.
- 5.App can rotate,take pictures and videos,playback, zoom in and out.
- 6.Converter interface is as follows:



Technical specifications

Converter Parameters	
Battery Capacity	1000mAh
Operate Voltage	4.2V
Operate Current	1A @ 3.7V
Operate Ambient Temperature	-20°- 60°
Operate Time	1h
Size (Model:W1)	85,5x50x10mm
Wireless transmission distance	10m

Converter Parameters	
Battery Capacity	400mAh
Operate Voltage	4.2V
Operate Current	1A @ 3.7V
Operate Ambient Temperature	-20°- 60°
Operate Time	45min
Size (Model:W2)	75x35x16mm
Wireless transmission distance	10m

Frequently Asked Questions

- 1.APP support the minimal version:
 - a)Android 4.4
 - b)IOS 8.0
- 2.If the device can not power on, or turn off automatically, it may be the battery has run out.

Warranty Information

- Warranty time: One year from the date of purchase(The date of purchase has prevailed)
- If the following occurs, this product warranty expires.
 - 1).Other non-equipment design, technology, manufacturing, quality problems cause malfunction or damage.
 - 2).Operate or place the device under abnormal conditions of use or service.
 - 3).Malfunction or damage due to force majeure.

Warning Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.