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
Real-time fat burning monitor

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
COMPONENTS

IMPORTANT!

1. Device will be locked after 300 tests.
If your device is locked, please contact local supplier or manufacturer.
2. Do not use alcohol to clean the device. It will cause device malfunction.
3. The following may affect the measurement results;
 - Consumption of alcohol in the last 24 hours.
 - Drinking coffee, juice or other beverage, mouth spray, brushing teeth, smoking 30 minutes before measuring.
4. If measurement result shows more than 60PPM continually for a long-term, consult with an expert.
5. Device required to be recalibrated or replace sensor every 300 tests or 12 months (whichever comes first)

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
Components Diagram

Function Displays

Breath sampling error
 It is caused by weak or unstable breath sampling.
 Wait few seconds until device is ready again
 then blow through a mouthpiece about 3 seconds.

Low battery warning
 Charge the device using the USB cable.

Battery charging (Blinking)
 It will take less than 1 hour to charge fully.

Battery charging is finished

Re-calibration required (Device blocked to operate)
 Device has measured more than 300 tests and
 re-calibration is required to maintain an accurate
 and sustainable operation.
 Please contact local supplier or manufacturer.

Sensor error
 It is caused by a sensor defect.
 Please contact local supplier or manufacturer.

Bluetooth® mode
 Device is connected with mobile app, KETOCAN

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

How to use

Device needs warm up time to be ready to blow and will determine automatically whether self-diagnosis mode is required or not.

- ①  Put the mouthpiece on Ketoscan mini
(Built-in mouthpiece is usable when you do not have an extra mouthpiece)
- ②  Press and hold the power button for 2 seconds.
Power button is located on the upper left side.
- ③  Warming up will be started and it takes about 60 seconds.
- ④  Device is ready for use.
If there is no operation, device will be turned off automatically after 60 seconds.
- ⑤  Take a deep breath and blow into the mouthpiece steadily and continuously for 3 seconds until the beep sounds twice. If your breath sample is not strong enough or unstable, "X" (sampling error message) will be displayed.
*Do not block the air outlet hole while the breath sampling is being finished
- ⑥  Breath sampling is finished and it's being analyzed.
- ⑦  If required, Self-diagnosis will be operated automatically and it will be shown "OK" on the display. Then device goes to measurement mode automatically and user will need to blow again to see a measurement result.
If device does not require self-diagnosis, measurement result will be displayed.
- ⑧  Check the measurement result.
Figure (PPM) will be shown for 15~20 seconds then device will be turned off automatically.
*0.0~9.9ppm: Displaying until one decimal point
Above 10ppm: No displaying one decimal point

** Beep sound level of Ketoscan mini is adjustable on the mobile App.

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The body's state according to breath acetone concentration

Breath Acetone Concentration (PPM)	Level	Effects
0.0 ~ 1.0	Level 0	Body fat is not burning at all or very small amount is burning now. It is required more intensive diet, exercise and controlling what you eat.
1.1 ~ 1.9	Level 1	
2.0 ~ 2.9	Level 2	It is in its early stage of body fat burning or low level of ketosis state. Body fat is slowly burning in average 0.07 oz (2 grams) per hour.
3.0 ~ 3.9	Level 3	
4.0 ~ 4.9	Level 4	
5.0 ~ 5.9	Level 5	It is the proper level of body fat burning or ketosis state. Body fat is burning in average 0.17 oz (5 grams) per hour.
6.0 ~ 6.9	Level 6	
7.0 ~ 7.9	Level 7	
8.0 ~ 8.9	Level 8	
9.0 ~ 9.9	Level 9	
10 ~ 40	Level 10	It is the best level of body fat burning or ketosis state. Body fat is burning in average 0.35 oz (10 grams) per hour.
41 ~ 60	Level Caution	It is the excessive level of body fat burning or ketosis state. If this state persists, it can be hard on the body.
61~	Level Danger	It is ketoacidosis state. And it will be a serious risk to health. Consult a health professional.

Product Specifications

Type of sensor	High resolution acetone gas sensor	Size	2.12x4.33x0.72 inches / 54x110x18.5mm
Measurement unit	PPM (Parts Per Million)	Weight	0.16lbs / 74g
Display range	0.0 ~ 99PPM	Storage	41~122°F (5~50°C)
Accuracy	10%	Working temperature	50~104°F (10~40°C)
Analyzing time	In 5 seconds	Calibration Interval	Every 300 measurements (even less than 300 measurements but more than 1 year since purchased, it is recommended recalibration or replacing sensor) *Self-diagnosis test does not counting
Power supply	Built-in rechargeable lithium-polymer battery	Battery life	More than 50 measurements
Power charging	USB port 5V 1000mA or more	Approval	FDA Class I
Mouthpiece	Disposable mouthpiece or built-in mouthpiece		

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FCC compliance information

THIS DEVICE COMPLIES WITH PART 15 OF THE FCC RULES.

Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and

(2) This device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.

Increase the separation between the equipment and receiver.

Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

Consult the dealer or an experienced radio/TV technician for help.

Modifications not expressly approved by the manufacturer may void the user's authority to operate the equipment under FCC rules.

FCC compliance information

THIS DEVICE COMPLIES WITH PART 15 OF THE FCC RULES.

Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and

(2) This device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15

of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

· Reorient or relocate the receiving antenna.

· Increase the separation between the equipment and receiver.

· Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

· Consult the dealer or an experienced radio/TV technician for help.

Modifications not expressly approved by the manufacturer may void the user's authority to operate the equipment under FCC rules.

Limited ONE year warranty

Warranty coverage

Sentech USA INC warrants this product against defects in materials and workmanship under normal use for a period of ONE

(1) YEAR from the date of retail purchase by the original end-user purchaser ("Warranty Period"). Except as provided in the preceding sentence or applicable law, Sentech USA INC provides this product "as is". If a defect arises and a valid claim is received by Sentech USA INC within the Warranty Period, at its option, Sentech USA INC will either (1) repair the product,

(2) exchange the product with a product that is new or which has been manufactured from new or serviceable used parts and is at least functionally equivalent to the original product, or (3) refund the purchase price of the product. When a refund is given, the product for which the refund is provided must be returned to Sentech USA INC and becomes Sentech USA INC' property.

This warranty is limited to the original end-user purchaser and is not transferable to, or enforceable by, any subsequent owner.

Disclaimers

Sentech USA INC makes no warranties, express or implied, as to the ability of this product to determine whether, or the extent to which, a test subject's mental or physical functioning, or judgment, may be impaired by blood alcohol including whether the test subject is intoxicated under any definition of that word. Sentech USA INC expressly disclaims any liability for direct, indirect, incidental, special, or consequential damages of any nature under any legal theory. Any act or failure to act based on a reading from this product shall be at the test subject's or user's own risk. Sentech USA INC assumes no responsibility for consequences to or of test subjects who use this product and later are shown to have been under the influence of alcohol or to have had their judgment or any mental or bodily function impaired by alcohol. Correlation between breath alcohol content and blood alcohol concentration depends on many variables, including environmental factors (such as air quality, wind, humidity, temperature, etc.) and health conditions of the test subject. A low BAC reading does not mean that the test subject's physical or mental performance or judgment (including reaction times) can respond to an emergency. The concentration of alcohol in the blood of a test subject cannot be exactly determined by using a breath alcohol screening product. Sentech USA INC does not warrant that the operation of the product will be error-free. Sentech USA INC is not responsible for any consequences arising from failure to follow instructions or relating to the product's use. This warranty does not apply: (a) to consequences caused by accident, abuse, tampering, misuse, flood, fire, earthquake or other external causes; (b) to consequences caused by operating the product outside the permitted or intended uses described by Sentech USA INC; (c) to consequences caused by service performed by anyone who is not a representative of Sentech USA INC or Sentech USA INC Authorized Service Provider; or (d) to a product or part that has been modified to alter functionality or capability without the written permission of Sentech USA INC.

TO THE EXTENT PERMITTED BY LAW, THIS WARRANTY AND THE REMEDIES SET FORTH ABOVE ARE EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, REMEDIES AND CONDITIONS, WHETHER ORAL OR WRITTEN, STATUTORY, EXPRESS OR IMPLIED. AS PERMITTED BY APPLICABLE LAW, SENTECH USA INC SPECIFICALLY DISCLAIMS ANY AND ALL STATUTORY OR IMPLIED WARRANTIES, INCLUDING, WITHOUT LIMITATION, WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE AND WARRANTIES AGAINST HIDDEN OR LATENT DEFECTS. IF SENTECH USA INC CANNOT LAWFULLY DISCLAIM STATUTORY OR IMPLIED WARRANTIES THEN TO THE EXTENT PERMITTED BY LAW, ALL SUCH WARRANTIES SHALL BE LIMITED IN DURATION TO THE DURATION OF THE EXPRESS WARRANTY AND TO THE REPAIR OR REPLACEMENT SERVICE AS DETERMINED BY SENTECH USA INC IN ITS SOLE DISCRETION.

No Sentech USA INC reseller, agent, or employee is authorized to make any modification, extension, or addition to this warranty. If any term is held to be illegal or unenforceable, the legality or enforceability of the remaining terms shall not be affected or impaired.

EXCEPT AS PROVIDED IN THIS WARRANTY AND TO THE MAXIMUM EXTENT PERMITTED BY LAW, SENTECH USA INC IS NOT RESPONSIBLE FOR DIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES UNDER ANY LEGAL THEORY, INCLUDING BUT NOT LIMITED TO LOSS OF USE; LOSS OF REVENUE OR INCOME; LOSS OF ACTUAL OR ANTICIPATED PROFITS (INCLUDING LOSS OF PROFITS FROM A CONTRACT); LOSS OF THE USE OF MONEY; LOSS OF ANTICIPATED SAVINGS; LOSS OF BUSINESS; LOSS OF OPPORTUNITY; LOSS OF GOODWILL, OR LOSS OF REPUTATION.

Some countries, states and provinces do not allow the exclusion or limitation of incidental or consequential damages or allow limitations on how long an implied warranty or condition may last, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary by country, state or province. This Limited Warranty is governed by and construed under the laws of the country in which the product purchase took place.

For warranty & recalibration service
Please contact;

SENTECH USA INC
20280 S.Vermont Ave. #200
Torrance, CA 90502
Tel: 310-526-8231
www.sentech-usa.com
service@sentech-usa.com

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



Pairing the device with the KETOSCAN mobile APP.

1


App. Download

Download mobile app. "KETOSCAN" available at App Store & Google play store.

*Creating an account will be required to keep the data storage online.


2

Bluetooth® pairing mode

Open the mobile App, first then turn on the device. Ketoscan mini will be available for pairing on the device list of App. Once Bluetooth® pairing is completed, it will be paired automatically in the future.

3

Connecting mouthpiece & Starting measurement

Put the mouthpiece into the built-in mouthpiece slot and follow the instructions on the mobile App.

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**** What if the Bluetooth® paring is not available?**

- 1) Turn on the device and search it on device list of App.
- 2) Or turn off the device and App. forcibly then try again.

***Self-diagnosis mode**

Device is checking a sensor status by itself,
and determines whether self -diagnosis mode is required or not.

Device is required to pair with a mobile App. Ketoscan to memorize & track your data.