

bebird R1

60x168mm



Download the APP

IOS
Search for "bebird" in App Store

Android
Search for "bebird" in Google Play




1 OR scan the QR code to download the app



- 2 Open bebird APP Select bebird R1
- 3 Connect Product WiFi "bebird-XXXX"

Accessory replacement

Pull the ear spoon out
Smart visual ear-clean Rod
Model: bebird R1
210-151073



↓ slightly push it in an angle and clip it tightly

Product parameters

Brand	bebird
Product Name	Smart visual ear-clean Rod
Model	bebird R1
Rod weight	13g
Network standard	IEEE 802.11b/g/n
Antenna	Integral Antenna
Working frequency	2.4Ghz
Image transmission rate	30fps
image sensor	CMOS
Working temperature	-10 ~ 50°C
Battery	130mAh lithium battery
Battery life	about 30 minutes
Charging time	0.5h
Input current	DC 5V 300mA
Lens diameter	4.5mm
Optimum focal length	1.5 ~ 2cm
Pixel	3M
Gravity Sensor	3 axis
FCC ID:2AU63-BEBIRD R1	



Attention

- 1、When cleaning lens, please wipe it carefully with professional alcohol cotton swab.
- 2、Please don't use it while you are moving, and protect from being hit by others.
- 3、Please do not put the product in liquid to avoid damage.
- 4、Product built-in rechargeable lithium battery, so it must be charged once a month If it is not used frequently.
- 5、Avoid product exposure to sunshine directly, especially the accessories, in order to avoid softening.
- 6、During using it, the temperature of the product will rise slightly (up to 35 degree), please rest assured to use.

Warranty card

User information	
User name:	Telephone number:
Purchase date:	
Seller:	Telephone number:

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.