



How to set up



Remove the wrap from the product and pull out the tab to activate the battery.

A 1

Stretchable silicone Infinite Loop with no clasps



Slide the strap across to insert it.

A 2

How to take a first measurement

auradevices.io/app




Update to the latest iOS and watchOS. Turn on Bluetooth on your Apple Watch. Download AURA App for your iPhone and Apple Watch.



A 3

To take a measurement:

- Open AURA App on your Apple Watch
- Touch and hold the electrodes with the base of your thumb (shaded zone)
- Wait until the end of the measurement.



A 4


Track all the data in the App for iPhone and Apple Watch.



A 5




How to replace battery or loop




Use included guitar pick or credit card to remove the cover.

B 1




Insert a new CR1632 lithium 3V coin battery with the positive side facing up.

B 2



Replace loops in accordance with the big and small holes on the connector.

B 3



Replace the cover, making sure that it is set correctly. Press down the opposite corners of the cover simultaneously until it clicks.

B 4

Need help?

For support and troubleshooting information visit help.auradevices.io



@auradevices

AURA Strap 2 Info

Safety and Handling

Important safety information
Handle AURA Strap 2 for Apple Watch with care. It contains sensitive electronic components, including batteries, and can be damaged or cause injury if dropped, burned, punctured, crushed, disassembled, or exposed to excessive heat or liquid. Don't use AURA Strap 2 if damaged.

AURA Strap 2, bioimpedance analysis feature and AURA App are not for medical use.

Compatibility
AURA Strap 2 works only with Apple Watch Series 1 or later. It may be used with a 38/40 mm or 42/44 mm version of Apple Watch, depending on your AURA Strap 2 (check it on the package).

Coin cell element
AURA Strap 2 contains a CR1632 Lithium 3V coin battery for getting power. The battery may need to be replaced during the life of the product. You can change it by yourself with the included guitar pick or credit card, or ask for help from an authorized service technician.

IMPORTANT: Do not allow moisture to get inside the case during battery replacement.

Choking hazard
AURA Strap 2 may present a choking hazard or cause other injuries to small children. Keep it away from small children.

Medical device interference
AURA Strap 2 contains components and radios that emit electromagnetic fields. These electromagnetic fields may interfere with pacemakers, defibrillators, or other medical devices. Maintain a safe distance between your medical equipment and AURA Strap 2. Consult your physician and medical device manufacturer for information specific to your medical device. Stop using AURA Strap 2 if you suspect it is interfering with your pacemaker, defibrillator, or any other medical equipment.

Warning
It is strictly prohibited to use the body composition and hydration level monitoring feature if you're pregnant or using a pacemaker or defibrillator.

ing feature if you're pregnant or using a pacemaker or defibrillator.

Skin irritation
AURA Strap 2 can lead to skin infections or irritation if not properly cleaned. Clean AURA Strap 2 regularly with a soft lint-free cloth. If a skin problem develops, discontinue use. If the problem persists, consult a physician.

Electrostatic shock
When using AURA Strap 2 in areas where the air is dehydrated, it is easy to build up static electricity and possible for your wrist to receive a small electrostatic discharge from AURA Strap 2. To minimize the risk of electrostatic discharge, avoid using AURA Strap 2 in arid environments, or touch a grounded unpainted metal object before putting AURA Strap 2 on.

Important handling information
Clean AURA Strap 2 with a soft, dry, lint-free cloth. Do not use liquids or cleaning products.

Support
For support and troubleshooting information, visit help.auradevices.io.

Apple Watch® and watchOS® are trademarks of Apple Inc.

Regulatory Compliance Information

FCC Compliance
This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses, and

can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications to this product not authorized by AURA Devices could void the electromagnetic compatibility (EMC) and wireless compliance and negate your authority to operate the product.

This product has demonstrated EMC compliance under conditions that included the use of compliant peripheral devices and shielded cables between system components. It is important that you use compliant peripheral devices and shielded cables between system components to reduce the possibility of causing interference to radios, televisions, and other electronic devices.

Responsible party (contact for FCC matters only):
AURA Devices Inc. 1013 Centre Road, Suite 403-B Wilmington, DE 19805 USA

Limited Warranty

The AURA Devices One Limited Warranty is a manufacturer's warranty. It provides a warranty where it is not provided by consumer law and/or if its duration less than one year. Consumers have the right to choose whether to claim service under the AURA One-Year Limited Warranty or their consumer law rights. Import-

ant: The AURA Devices One-Year Limited Warranty terms and conditions shall not apply to consumer law claims. For further information about consumer law, please visit the AURA Devices website (www.auradevices.io/warranty) or contact your local consumer organization.

Please note: All claims made under the AURA Devices One-Year Limited Warranty will be governed by the terms set out in this warranty document.

Your AURA Devices-branded hardware product ("Product") is warranted against defects in materials and workmanship for ONE (1) YEAR from the date of original retail purchase ("Warranty Period") when used in accordance with AURA Devices' user manuals (refer to help.auradevices.io). Under this warranty, you will be able to direct your claims to AURA Devices even in situations where you purchased the AURA Devices Product from a third party. If a defect arises during the Warranty Period, AURA Devices, at its option will (1) repair the Product at no charge using new parts or parts that are equivalent to new in performance and reliability, (2) exchange the Product with a product with equivalent functionality formed from new and/or previously used parts that are equivalent to new in performance and reliability or with your consent, a product that is at least functionally equivalent to the Product it replaces. This warranty excludes normal depletion of consumable parts such as batteries unless failure has occurred due to a defect in materials or workmanship and, damage resulting from abuse, accident, modifications, unauthorized repairs or other causes that are not defects in materials and workmanship.

AURA Devices does not warrant, represent or undertake that it will be able to repair or replace any Product under this warranty without risk to and/or loss of information and/or data stored on the Product. In no event shall AURA Devices be liable for (a) loss or damage, which as of the Product's purchase cannot be regarded as being caused by AURA Devices' breach of these warranty terms; or (b) losses caused by the user's fault, loss of data or loss of profits or benefits.

Any limitations of liability in this warranty document shall not apply to (i) death or personal injury pursuant to any mandatory law on product liability; (ii) fraud or fraudulent misrepresentation; (iii) in-

tentional misconduct or gross negligence; (iv) or a culpable breach of major contractual obligations. A damages claim based on a breach of major contractual obligations or gross negligence will be limited to foreseeable damage typical for the sale contract concerned.

To obtain warranty service, contact AURA Devices by email: support@auradevices.io. Proof of purchase may be required to verify eligibility.

AURA Devices will provide warranty service through one or more of the following options:

- (i) Mail-in service. If AURA Devices determines that your Product is eligible for mail-in service, AURA Devices will ask you to send the Product by specific address. AURA Devices will pay for shipping to and from your location. If instructions regarding the method of packaging and shipping the Product are followed.
- (ii) Do-it-yourself (DIY) parts service. DIY parts service allows you to service your Product. If DIY parts service is available in the circumstances, the following process will apply.

(a) Service where AURA Devices requires the return of the replaced Product or part. AURA Devices will ship a replacement Product or component to you with installation instructions, if applicable, and any requirements for the return of the replaced Product or part.

(b) Service where AURA Devices does not require the return of the replaced Product or part. AURA Devices will ship you free of charge a replacement Product or part accompanied by instructions on installation, if applicable, and any requirements for the disposal of the replaced Product or part.

(c) AURA Devices is not responsible for any labor costs you incur relating to DIY parts service. Should you require further assistance, contact AURA Devices by email: support@auradevices.io.

If any term is held to be illegal or unenforceable, it shall be severed from this warranty and the legality or enforceability of the remaining terms shall not be affected.

Software distributed by AURA Devices with or without the AURA Devices brand

(including, but not limited to system software) is not covered by this Warranty. Please refer to the licensing agreement accompanying the software for details of your rights with respect to its use. AURA Devices does not warrant that the operation of the AURA Devices Product will be uninterrupted or error-free. AURA Devices is not responsible for damage arising from failure to follow instructions relating to the AURA Devices Product's use.

This Warranty does not apply: (a) to consumable parts, such as batteries or protective coatings that are designed to diminish over time, unless failure has occurred due to a defect in materials or workmanship; (b) to cosmetic damage, including but not limited to scratches, dents and broken plastic on ports unless failure has occurred due to a defect in materials or workmanship; (c) to damage caused by use with a third party component or product that does not meet the AURA Devices Product's specifications (AURA Devices Product specifications are available at www.auradevices.com under the technical specifications for each product and also available in stores); (d) to damage caused by accident, abuse, misuse, fire, earthquake or other external cause; (e) to damage caused by operating the AURA Devices Product outside AURA Devices' published guidelines; (f) to damage caused by service (including upgrades and expansions) performed by anyone who is not a representative of AURA Devices; (g) to an AURA Devices Product that has been modified to alter functionality or capability without the written permission of AURA Devices; (h) to defects caused by normal wear and tear or otherwise due to the normal aging of the AURA Devices Product; (i) if any serial number has been removed or defaced from the AURA Devices Product; or (j) if AURA Devices receives information from relevant public authorities that the product has been stolen or other security measures designed to prevent unauthorized access to the AURA Devices Product, and you cannot prove in any way that you are the authorized user of the product (e.g., by presenting proof of purchase).

This warranty is governed by and construed under the laws of the country in which the product purchase took place.

For Australian consumers: The rights described in this warranty are in addition

to the statutory rights to which you may be entitled under the Competition and Consumer Act 2010 and other applicable Australian consumer protection laws and regulations. Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure. Repair of the goods may result in loss of data. Goods presented for repair may be replaced by refurbished goods of the same type rather than being repaired. Refurbished parts may be used to repair the goods.

For Canadian Consumers: Residents of Quebec are governed by that province's consumer protection legislation.

For purchases made by consumers in the United Kingdom and Ireland: If a product is defective consumers may, in addition to any other rights which they may have under consumer law in the UK and Ireland, avail themselves of the rights contained in: for products purchased in Ireland: the Sale of Goods Act, 1893 (in particular Sections 12, 13, 14 and 15), the Sale of Goods and Supply of Services Act, 1980 and the European Communities (Certain Aspects of the Sale of Consumer Goods and Associated Guarantees) Regulations 2003 (S.I. No. 11/2003); for products purchased in the UK: the Sale of Goods Act 1979 (in particular Section 12), the Supply of Goods and Services Act 1982 (in particular Section 2) and the Sale and Supply of Goods to Consumers Regulations 2002.

Products Warranty (English) v1.1.0

Quick start guide

Welcome to AURA Strap 2, an accessory that boosts your Apple Watch with hydration levels and body composition tracking.

Push your limits. AURA Devices team

For support and troubleshooting information visit help.auradevices.io



FCC Statement

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.