SUMMARY OF IMPORTANT SAFETY NOTES	
Before using the URGOnight headband, make sure to read and follow these instructions, see all warnings, and keep this document for future reference.	
1. If you are travelling by air, the URGOnight headband must be transported in your carry-on baggage, do not place them in your checked baggage.	
2. Do not use the URGOnight headband while it is charging.	
3. Do not use the URGOnight headband when your head/hair is wet, or near water.	
4. Do not use the URGOnight headband if it has been damaged or altered.	
5. Only use accessories that are specified by the manufacturer.	
6. Do not place the URGOnight headband on open wounds, sores or rashes, or over swollen, red, infected or inflamed skin irritations. If any skin irritations occur in using the URGOnight headband, you should stop using it immediately.	
7. Follow cleaning recommendations outlined in the cleaning section.	
	1

User Manual

Charging

1. Please charge the headband for at least 1 hour before the first use.
2. Only use the charger and cable provided with the URGOnight headband by the manufacturer.
3. The URGOnight headband can only be charged when it is switched off.
4. The URGOnight headband must not be worn when charging, unplug the headband from the charger before use.

Water, moisture, heat

1. Do not use the URGOnight headband near water or when you are wet, especially when your head/hair is wet.
2. Do not immerse the URGOnight headband in liquid of any kind, or placing it in a dishwasher or washing machine.
3. The URGOnight headband is not water resistant but may be wiped down with a dampened cloth for cleaning.
4. Do not leave the URGOnight headband near any heat sources, such as radiators, heat registers, stoves or other apparatus that produce heat, such as audio amplifiers or televisions.

CONFORMITY

The declaration of conformity and any other regulatory documents specific to other countries are available on our website: urgonight.com

SAFETY

The device is compliant according to these standards: EN IEC 62368-1:2018

RADIO & ELECTROMAGNETIC INTERFERENCE

The device is compliant according to Electromagnetic compatibility (EMC) tests:

- EN 55032 :2015/EN 55035 :2017/ETSI EN 301489-1 V2.1.1/ETSI EN 301489-17 V2.2.1/ETSI EN 300328 V2.1.1
- CFR 47 Part 15 Subpart B - Radio frequency devices

WARNING TO USERS IN THE UNITED STATES

Federal Communication Commission Interference Statement 47 CFR Section 15.105(b)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed

PRODUCT

Product Operation

To get started with your URGOnight headband, please refer to the Starter Guide enclosed in the packaging.

Product Description

URGOnight was designed to encourage and teach you how to take control of your sleep by using neurofeedback in an autonomous way. The URGOnight headband measures your brain activity by detecting electrical impulses with 4 dry sensors (electrodes) that provide 4 channels of Electroencephalogram (EEG) data.

The URGOnight mobile application will guide you to set up your headband, take you through the exercises and help you keep track of your progress. URGOnight is intended to be used during the day and has been designed to be easily transported anywhere you like.

Health Considerations

URGOnight is not a medical device and is not designed to treat any medical condition.

If you are diagnosed with one of the following conditions (non-exhaustive list), please consult your doctor prior to using URGOnight:

1. Depression
2. Organic sleep troubles: apnea, hypopnea, restless syndrome, narcolepsy, hypersomnia, Kleine-Levin syndrome, somnambulism, nocturnal eating troubles
3. Pathologies disturbing sleep: cancer, degenerative, autoimmune, inflammatory, chronic or acute, hormonal disturbances
4. For users that have a medical implant, the signal of your URGOnight training may be inaccurate, please consult your doctor.
5. Epilepsy: If you are epileptic, please check with your neurologist if URGOnight is suitable for your condition. As the brainwaves trained with URGOnight are close to the beta waves, that might cause seizures in some epileptic people.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment and meets the FCC radio frequency (RF) Exposure Guidelines. This equipment has very low levels of RF energy that is deemed to comply without testing of specific absorption rate(SAR).

WARNING TO USERS IN THE CANADA


This device complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

Under Industry Canada regulations, this radio transmitter may only operate using an antenna of a type and maximum (or lesser) gain approved for the transmitter by Industry Canada.

To reduce potential radio interference to other users, the antenna type and its gain should be so chosen that the equivalent isotropically radiated power (e.i.r.p) is not more than that necessary for successful communication.

This equipment complies with IC radiation exposure limits set forth for an uncontrolled environment and meets RSS-102 of the IC radio frequency (RF) Exposure rules. This equipment has very low levels of RF energy that is deemed to comply without testing of specific absorption rate (SAR).

Headband

1. **ON/OFF Power Button:** Gently press once to switch the headband on. To switch off, hold down the button for more than 2 seconds.
2. **LED (color indications):**
 - Flashing white: the URGOnight headband is on and searching for mobile phone pairing
 - Still white: Paired with the mobile application, ready to start
 - Still orange: Charging
 - Still orange + green: Fully charged

Physical Symptoms

In rare cases, some people experience blackouts or seizures due to exposure to flashing lights and patterns created by the display on a mobile phone or other devices with screens. If you have experienced in the past any of these symptoms, or have experienced any nausea, involuntary movements, tingling, numbness or vision issues while using these types of applications or devices, you should consult your doctor before using URGOnight or any similar applications. You should immediately cease all such use of URGOnight should the symptoms recur. Nevertheless, avoid prolonged use of URGOnight to avoid any possible fatigue or discomfort, including any eye, muscle or joint strain.

Choking Hazard

URGOnight and its accessories may present a choking hazard. Keep these accessories away from young children.

Radio Signals

The URGOnight headband uses radio signals to connect to wireless networks. Observe signals and notices that prohibit or restrict the use of wireless devices (for example, healthcare facilities, airlines and blasting area). Although the URGOnight headband is designed, tested and manufactured to comply with regulations governing radio frequency emissions, such emissions from the URGOnight headband may negatively affect the operation of other electronic equipment, causing them to malfunction.

Permitted Age & Profile


URGOnight is available for individuals above the age of majority. If the user is a minor, the user must review this documentation with their parent or guardian to make sure that both the user and the parent or guardian, understand it.

Europe - CE

URGOTECH is authorized to apply the CE Mark on the device, thereby declaring conformity to the essential requirements and other relevant provisions of Directive 2014/30/EU, 2014/35/EU, 2014/53/EU.

ENVIRONMENT

Disposal of Waste Electrical and Electronic Equipment

1. You may consult and comply with all local laws and regulations for the proper disposal of batteries and electronic devices in your local area.
2. If there are no local laws concerning electronic waste disposal, place the headband and its battery in a waste bin for electronic devices.
3. Batteries may not be disposed of in municipal waste streams and require separate collection.
4. Disposal of a battery into fire or a hot oven, or mechanically crushing or cutting of a battery, can result in an explosion; a battery subjected to extremely low air pressure that may result in an explosion or the leakage of flammable liquid or gas.

- Flashing orange: There might be a problem with your battery, please contact customer support

Note: In order to optimize your battery, the URGOnight headband will automatically switch off after 2 minutes. It will switch back on if you move it or if you gently press the on/off button.

3. **USB connector:** It is positioned under the URGOnight headband.

4. **2 Top Electrodes:** One on each side

5. **2 Bottom Electrodes:** One on each side

URGOnight is not to be used by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction.

Choking Hazard

URGOnight and its accessories may present a choking hazard. Keep these accessories away from young children.

Radio Signals

The URGOnight headband uses radio signals to connect to wireless networks. Observe signals and notices that prohibit or restrict the use of wireless devices (for example, healthcare facilities, airlines and blasting area). Although the URGOnight headband is designed, tested and manufactured to comply with regulations governing radio frequency emissions, such emissions from the URGOnight headband may negatively affect the operation of other electronic equipment, causing them to malfunction.

WARRANTIES & SERVICES

URGOTECH warrants to the original purchaser that your URGOnight headband shall be free from defects in material and workmanship for 1 year. Except where prohibited by applicable law, this warranty is nontransferable and is limited to the original purchaser and the country in which the product was purchased. This warranty gives you specific legal rights, and you may also have other rights, including a longer warranty duration that may vary under local laws.

DISCLAIMER & LIMITATION OF LIABILITY

The URGOnight headband is not a medical device, it is not meant for medical use, clinical diagnosis or therapeutic purposes.

URGOnight is designed for regular and consistent training sessions of 20 minutes, minimum 3 times a week to observe benefits on your sleep. It is recommended to leave at least one night's rest in between 2 sessions to allow enough time to your brain to assimilate what has been learned. URGOTECH assumes no responsibility with regards to the performance of the URGOnight headband.

© 2018 URGOnight. All rights reserved. URGOnight and related trademarks, names and logo are the property of URGOTECH and are registered and/or used in the European Union, US. "Android" is a trademark of Google

Technical Specifications	
Size	One Size Fits All, adjustable (front, sides and back)
Head Circumference	520 mm to 620 mm (20.47" to 24.40")
Electrodes	4 dry gold plated and retractable electrodes for high quality signal measurement and optimal comfort
Weight	160g
Materials	Hypoallergenic Silicone, Fabric
Battery Life	Up to 6 hours, equivalent to 2 weeks of recommended use
Battery Charge	3 hours to reach 100% full charge
Connectivity to Application	Bluetooth 4.0 low energy
URGOnight Application Compatibility	Apple iOS 11 and later or Android 5 and later
	5
	6

Product Cleaning

1. Clean only with a non-abrasive, lint-free cloth. If needed, you can lightly dampen the cloth with water, and soap-free cleanser to clean the headband. Always air dry the band well before putting it back on.
2. Do not use detergents or cleaning gels, which might degrade the surface of the electrodes.

SAFETY & WARNINGS

Read and understand all instructions in this document.

Battery

1. The URGOnight headband is equipped with two Lithium Polymer batteries of 350mAh.
2. If you are travelling by air, you must carry the URGOnight headband in your carry-on baggage, do not place them in your checked baggage.

use while walking down the street or driving. Failure to pay full attention during walking or the operation of your vehicle could result in death, serious injury, or property damage. You assume total responsibility and risk for not following this operating instruction.

Operating Conditions

1. Temperature range: -20°C/-4°F to 30 °C/86°F
2. Charging Temperature: 0°C/32°F and above
3. Humidity: 40 - 95%
4. Electrical Specifications for charging:
 - Voltage input min/max: 4.5V - 5.5V DC (direct current)
 - Voltage nominal input: 5V DC (direct current)
 - Power input: 5W
5. Do not use other charging adapters unless it meets the following specifications:
 - Input : 100 - 240V - 50/60 Hz 0.4A
 - Output : 5V 1000 mA
 - Model : MKC-0501000DEXU
 - Certified: IEC 60950-1 or IEC 62368
6. The charger must be easily accessible as its power plug acts as a disconnecter.

GENCE, URGOTECH, ITS LENSORS AND SERVICE PROVIDERS WILL NOT BE LIABLE UNDER ANY CONTRACT, NEGLIGENCE, STRICT LIABILITY OR OTHER LEGAL, OR EQUITABLE THEORY FOR ANY CLAIMS ARISING FROM OR RELATED TO: (a) ANY INFORMATION AND DOCUMENTATION MADE AVAILABLE TO YOU BY URGOTECH, ITS LENSORS AND SERVICE PROVIDERS; (b) USE OF URGOTECH PRODUCTS, SOFTWARE, OR SERVICES; OR (c) THE INABILITY TO USE OR THE NON-PERFORMANCE OF URGOTECH PRODUCTS, SOFTWARE, OR SERVICES, IN ALL CASES WHETHER OR NOT SUCH DAMAGES WERE FORESEEN OR UNFORSEEN, EVEN IF URGOTECH HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES, IN NO EVENT SHALL THE AGGREGATE LIABILITY OF URGOTECH, ITS LENSORS AND SERVICE PROVIDERS EXCEED THE AMOUNT PAID BY YOU FOR THE PRODUCTS, SOFTWARE, AND SERVICES GIVING RISE TO THE CLAIM.

SUPPORT & SERVICES

If you wish to receive further information or if you need any customer service, please send in a request on our website: urgonight.com or send an email to support@urgonight.com

URGOTECH SAS

15 avenue d'Iena,
75116 Paris, France

Symbol	Meaning
— —	Direct Current (DC Power)

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