EN

FROSTMED

User Manual

Important Safety Information

These warnings are provided in the interests of your safety. Ensure that you understand them all before installing or using FROSTMED. Your safety is of paramount importance. If you are unsure about any of the meanings or these warnings contact the Customer Care Department (Section 8 of this User Manual).

This User Manual applies to the following models:

- FROSTMED PRO G1 and FROSTMED PRO G2
- FROSTMED LITE G1 and FROSTMED LITE G2
- FROSTMED LITE CV G2 and FROSTMED LITE CV G2.

Intended Use

FROSTMED is a cooling device produced by FROSTIME LLC. FROSTMED utilizes compressed liquid CO₂ as coolant instead of ice, dry-ice, gel packs or phase change materials and it is available in the 3-above models.

When placed in an insulated container, FROSTMED is used to transport temperature sensitive biological/medical materials for short distance (the so called “Last Mile”). The desired temperature range is set using a FROSTMED app available for iOS and Android devices.

Cooling lasting time of one set of two fully charged CO₂ canisters varies from 10 up to 40 hours according to several external parameters (as specified in Section 6) which are not under FROSTMED control. See Section 6 for some usage examples.

FROSTMED allows to track, log and report over a period of time the temperature inside the insulated container where it is placed.

FROSTMED intended uses include, but are not limited to:

- Cold chain transport of vaccines, insulin, medicines and specimens requiring specific temperatures.
- Cold chain transport of blood and blood products requiring specific temperatures.
- Cold chain transport of organs requiring specific temperatures.
- Emergency storage of the above in the event of power failures.
- Cold storage of the above for areas afflicted by disaster.
- Cold storage of the above for remote areas with little or no infrastructure.

Prior to Installation

Check FROSTMED for transportation damage. Under no circumstances should a damaged FROSTMED be installed. In the event of damage please contact the place of purchase.

Installation

- Care should be taken when moving FROSTMED.
- It is dangerous to alter the specifications or attempt to modify this product in any way.
- All Canisters should be secure before turning FROSTMED on.
- FROSTMED should be placed flat on the bottom of the given cooler and/or ice chest that it is placed in.

Child Safety

Do not allow children to use the Bluetooth app that controls FROSTMED unless under adult supervision!

During Use

- Do not place ICE inside the cooler when using FROSTMED.
- Containers with flammable gases or liquids can leak at low temperatures. Do not store any containers holding flammable materials, such as spray cans, fire extinguisher refill cartridges etc. with FROSTMED.
- Never submerge FROSTMED.

Maintenance and Cleaning

- Before cleaning, always switch off FROSTMED.
- Wipe all surfaces with an antimicrobial solution (isopropyl alcohol) and a cotton swab.

Servicing

- This product should be serviced by an authorized engineer and only genuine spare parts should be used.
- Under no circumstances should you attempt to repair FROSTMED by yourself. Repairs carried out by inexperienced persons may cause injury or serious malfunctioning.

Disposing of the FROSTMED

When disposing your FROSTMED, use an authorized disposal site. Check with your local Council or Environmental Health Office to see if there are facilities in your area to recycle FROSTMED.

Handling of CO₂

CO₂ should be handled in accordance with the guidelines and restrictions outlined in this manual to ensure safe use.

CO₂ Canister Transportation

- CO₂ canisters are recommended to be transported in an upright position.
- Verify that canisters are completely sealed before transportation.
- DO NOT pull/tamper with pressure release valves.
- CO₂ filled canisters CAN only be shipped via ground or via air freight in accordance with Federal Regulations.

CO₂ Temperature

- Canisters are to be stored in areas that avoid extreme heat and cold.

CO₂ Pressure

- Canisters should not be exposed to areas with extreme variances in air pressure.

CO₂ Refilling Process

- Canisters refilling process is explained during the training by authorized Frostime's personnel. Please contact Frostime for your refilling needs. Different options are available.
- Only Food or Beverage Grade CO₂ should be used to refill FROSTMED CO₂ canisters.
- Confirm that canisters are not filled to more than 80% capacity to allow for liquid/gaseous expansion.

CO₂ Canister Storage and Transportation (empty/full)

- Canisters should always be stored in an upright position.
- Storage should be at room temperature avoiding extreme temperatures (both hot and cold).
- Do not store under direct sunlight.
- Assure that the canisters are kept away from moisture and other weather conditions that could adversely affect the container.
- Storage should be in a normal pressure environment.
- Keep canisters away from pointed objects or anything that might pierce the container.
- In order to avoid denting or valve damage, DO NOT let CO₂ canisters tumble or roll loose. Always secure CO₂ (empty or full) while traveling.

In the very rare case of a frostbite

Don't touch the plate when it is extremely cold. If accidentally it happens:

- Gently warm the area in warm water or with wet heat until the skin appears red and warm.
- If no water is nearby, breathe on the area through cupped hands and hold it next to your body.
- Do not use direct heat from heating pads, radiator, or fires.

- Do not rub or massage the skin.
- Loosely apply dry, sterile dressings.
- Put gauze or clean cotton balls between fingers or toes to keep them separated.
- Consult medical professional for further treatment.

Content protection

FROSTMED coolant is food-grade CO₂ with 99.90% purity rating. FROSTMED is NOT STERILE.

Biological, medical and any other transported products MUST BE protected and isolated from the internal container's environment.

Compliance

IC and FCC Compliance

This device complies with Industry Canada license-exempt RSS standard(s). This device complies with Part 15 of the FCC rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FCC/ISED RF exposure

This equipment complies with FCC radiation exposure limits and Canada radiation RF exposure limits set forth in CFR 47 Sections 2.1091 and Industry Canada RSS-102 for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 20 cm between the radiator and the body. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

FCC Caution: Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

Digital class B notice

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.

- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FDA clearance

FROSTMED devices (PRO G1, PRO G2, LITE G1, LITE G2, LITE CV G1 and LITE CV G2 models) are FDA-cleared as Class II device, 510(K) exempts.

Table of Contents

1	INTRODUCTION	8
2	STARTING FROSTMED PRO, LITE, LITE CV	8
2.1	FROSTMED SET-UP	8
2.2	ACTIVATING THE ELECTRONICS.....	10
2.3	CONNECTING A MOBILE DEVICE.....	10
3	OPERATING FROSTMED	11
3.1	SETTING THE TEMPERATURE	11
3.2	READING THE LED INDICATORS	14
3.3	CHARGING LEVEL OF THE RECHARGEABLE BATTERY PACK.....	15
3.4	RECOMMENDATIONS.....	15
4	REMOVING EMPTY CO₂ CANISTERS.....	16
4.1	REMOVING CO ₂ CANISTERS	16
5	ALERT PROGRAMMING AND NOTIFICATION	16
6	TROUBLESHOOT	16
7.	SERVICE AND SPARE PARTS.....	18
8.	CUSTOMER SERVICE DEPARTMENT.....	18

1. Introduction

Congratulations on the purchase of your new FROSTMED portable cooling system!

FROSTMED is an innovative cooling product that uses liquid Carbon Dioxide (CO₂) to maintain a low temperature inside your insulated container. The FROSTMED's cooling ability is based on the principle that when liquid CO₂ (compressed CO₂) is released into the atmosphere it absorbs a large amount of heat from the surrounding items which become very cold.

FROSTMED is the ideal source of cold when you transport vaccines, medicine, specimens, blood products, beverage, foods and any temperature sensitive items.

You can monitor and control the temperature inside the insulated container in real-time using the FROSTMED smartphone app (available for G1 and G2 models) or using the FROSTMED web application (only available for G2 models).

Thank you from everyone on FROSTMED team!

2 Starting FROSTMED PRO, LITE, LITE CV

2.1 FROSTMED Set-Up

FROSTMED PRO, LITE and LITE CV models use CO₂ canisters to provide coolant.

To set-up your FROSTMED follow the below instructions:

- First open your FROSTMED package.
- The package contains the following items:
 - 1 FROSTMED device (PRO or LITE or LITE CV)
 - 2 CO₂ refillable canisters (2,5 lbs. each for the PRO model, 1 lb. each for the LITE and LITE CV models)
 - 1 rechargeable Lithium-Ion battery pack
 - 1 USB cable
 - 1 canister retaining clamp (for PRO version only. It can be used to better secure canisters in enclosure)
 - 1 carry handle (for PRO model only)
 - 1 warranty
 - 1 User Manual.
- The rechargeable battery pack must be charged before first use. Charge the rechargeable battery pack via a Micro-USB input on the top of the battery pack case (see Figure 2). The battery pack connects to any common USB phone charger (iPhone /Android) via a Micro-USB cable.

Figure 1: Battery Pack

- Insert the rechargeable battery pack in the battery as shown in Figure 2 and Figure 3.

Figure 2: Battery location for PRO model

Figure 3: Battery location for LITE and LITE CV models

- If you purchase the PRO model, unscrew the canister retaining clamp and separate the 2 CO₂ canisters. LITE and LITE CV models don't have the canister retaining clamp and come already separated.
- Orient the wing of the first CO₂ canister's valve to be parallel into the hole of the inside the FROSTMED's base (canister in upside down position).
- Slowly lower the canister into the FROSTMED's base. **NOTE:** If you do not align the canister with the channel inside the base, don't worry, simply rotate the canister until it drops into the channel.
- Firmly push the canister down and hold.
- Rotate the canister one quarter turn in the clockwise direction. You will feel the canister make contact with a dead stop.
- Repeat the above process for thesecond canister and your FROSTMED will be loaded.
- If you purchased the PRO model, secure the 2 canisters together with the canister retaining clamp. Slide the canister retaining clamp over the two canisters with the thumbscrew over the top of the battery slot on FROSTMED PRO base. Turn the thumb screw until the plastic strap is firmly clamping both canisters together.

NOTE: Slight bursts of gas may occur when seating the canisters, this is normal. A continuous hiss without the electronics is a sign of a misalignment or a faulty valve. Consult the troubleshooting section if the problem persists.

2.2 Activating the Electronics

The following instructions apply to FROSTMED PRO, LITE and LITE CV models. To activate the electronics, press the ON/OFF button. The valve will automatically open and begin cooling. You will hear a hissing sound as the gas exits FROSTMED, this is normal.

2.3 Connecting a Mobile Device

This cooling system is controlled by the FROSTMED app. The operator must have the app downloaded to the operator's Bluetooth device (i.e.: smartphone, tablet, etc.) to control the temperature and other FROSTMED features. The app can be found on the [*Google Play, and Apple app stores*].

To connect the FROSTMED via Bluetooth, open the FROSTMED app represented by the following icon:

From the main screen, touch the three vertical dots in the top right corner to open the options. From the dropdown list, select "Scan for FROSTMED".

Figure 4: Scan FROSTMED

Your mobile Bluetooth device is now connected wirelessly to your FROSTMED cooling system.

Figure 5: FROSTMED Connected and Ready to be Operated

3. Operating FROSTMED

3.1 Setting the Temperature

From the main screen of your Frostime app, you can:

- 1) set the desired temperature (see Figure 6) or

2) select the biological product to be transported (Figure 7). In this case FROSTMED automatically set the temperature according to the following temperature ranges:

- a. Blood:
 - Min: +1°C
 - Max: +10°C
 - Setpoint: +°4C
- b. Platelets:
 - Min: +18°C
 - Max: +24°C
 - Setpoint: +21°C
- c. Sample:
 - Min: +12°C
 - Max: +20°C
 - Setpoint: +15°C
- d. Vaccines
 - Min: +2°C
 - Max: +8°C
 - Setpoint: +5°C
- e. Insulin
 - Min: +2°C
 - Max: +8°C
 - Setpoint: +5°C
- f. Corneas
 - Min: -8°C
 - Max: +2°C
 - Setpoint: -3°C
- g. Room temperature specimens
 - Min: +20°C
 - Max: +25°C
 - Setpoint: +23°C
- h. Refrigerated specimens
 - Min: +4°C
 - Max: +8°C
 - Setpoint: +6°C
- i. Manual

If none of the above range applies to your desired temperature range, please use option 1).

To get back to option 1), if currently using a content type, select content type “i. Manual” from the list of content types.

Figure 6: Option 1) Manual Temperature Setting


Figure 7: Option 2) Pre-configured Temperature Setting

The system reads the temperature inside the container and compares it to the FROSTMED Bluetooth App set point temperature. If the container temperature is above or equal to the set point temperature, FROSTMED will close the valve and stop releasing CO₂.

When the temperature in the container rises above the desired temperature, the system will start releasing CO₂ until the container temperature falls below the set point temperature.

At any moment you can control the temperature inside the container using your app (see Figure 9, 10, 11 and 12) or using a web dashboard (only available for G2 models).

Figure 8: Manual Option has been selected

Figure 9: Blood option has been selected


Figure 10: Celsius/Fahrenheit Conversion

Figure 11: Blast function

At any moment you can switch from Celsius to Fahrenheit Degree selecting “C/F” from the main menu on the top right of the app (see Figure 10).

The “Blast” function in Figure 10 and 11 allows to temporary open the valve and cool down the interior of the container in case you need an immediate source of cold.

Recommendation: Use the “Blast” function only if strictly needed. This function uses a higher amount of CO₂ quantity and can reduce the cooling time of the system.

3.2 Reading the LED Indicators

The LED indicators for the FROSTMED models: PRO, LITE and LITE CV, show different states of the system:

- A **Red Light** indicates the system is powered ON and operating/cooling.
- A **slowly blinking Blue Light** indicates that the Bluetooth connection is available.
- A **rapid blinking Blue light** indicates the Bluetooth connection is established.
- A **Green Light** indicates the control valve is turned on and CO₂ will be released.

Figure 12: LED Indicators for the FROSTMED PRO.

3.3 Charging Level of the Rechargeable Battery pack

To determine the level of the charge of the rechargeable battery pack, push the button on the battery pack.

When the rechargeable battery is activated a set of four green LEDs will appear. A 100% charge is equivalent to 4 LED lights. The higher the number of LEDs present, the more the battery is charged

3.4 Recommendations

1. The device is NOT sterile and biological/medical materials MUST BE protected and isolated from the internal cooler's environment.
2. Double check the exposure temperature after transportation:
 - Control the logged temperatures using the suitable app or/and
 - Use a temperature monitor strip as temperature control backup (examples: 3M™ MONITORMARK™ TIME TEMPERATURE INDICATOR 9860C 10°C/50°F 48 HOUR, or Berlinger Mini-tag or equivalent).

4. Removing Empty CO₂ Canisters

4.1 Removing CO₂ Canisters

Rotate the canister a quarter turn counterclockwise. **Note: It is normal to hear a small hiss of gas while removing the canister. If hissing continues after removal of the tank see the trouble shooting section.**

CAUTION: Compressed gases have a large amount of stored energy within them. Removal of the canister at partial pressure may result in unanticipated motion of the canister. Be cautious when removing partially filled canisters.

5. Alert Programming and Notification

It is possible to receive Alerts Notification when:

- The programmed temperature exceeds the desired threshold
- CO₂ level is lower than the desired threshold
- Battery level is lower than the desired threshold.

This feature is available for the User having a Web dashboard included in their package.

To receive the Alert notification, access to the FROSTMED Web Dashboard screen: “Edit Alerts” and program the thresholds, the desired phone number and the email address where receive the alerts. For further details refer to the “Help” in the FROSTMED Web Dashboard.

ATTENTION: Alerts are sent only when FROSTMED is connected to the app via Bluetooth.

6. Troubleshoot

Factors that impact FROSTMED performance:

- Frequency of opening & closing container lid.
- Initial temperature of contents in the container.
- Insolation factor of the cooler.
- Amount of contents in cooler (fuller the cooler, less CO₂ used).
- Size of the cooler.
- Ambient temperature.

Figure 13 shows the lasting time duration with different external temperatures for FROSTMED LITE with a setup temperature at +6°C and mock red blood cells as payload (tests conducted by Canadian Blood Service – August 27, 2019):


Figure 13: Lasting time duration with external temperature

If FROSTMED is not working correctly, please carry out the following necessary checks before contacting your local FROSTMED Customer Service Representative.

No power to the FROSTMED

Check that the rechargeable battery pack in FROSTMED is fully or partially charged.

The temperature in the cooler is too high Ensure that the thermostat is turned to the correct temperature setting relevant to ambient/outside the cooler temperature conditions.

Excessive frost, ice build-up

Check the ambient/outside temperature level. If the ambient/outside temperature is higher than the set point temperature of the app then ice may build up. Ensure that the ambient/outside temperature is not below 10°F/-12°C.

Normal Operating Sounds

Anytime the valve opens up, there is a short hiss sound for about 5 seconds. This is considered normal. If the hissing sound lasts more than 10-15 seconds, refer to 'Noise'.

Noise

Was the noise related to the opening of the cooler lid?

- If so, open the container lid for a few minutes to allow the air pressure to equalize.
- Ensure the container lid has not been left open.
- Check that nothing is preventing the container lid from closing.
- Refer to 'Normal Operating Sounds.'

If after performing all the necessary checks, your FROSTMED still does not operate correctly, please contact your local FROSTMED Service Center Representative.

Please note that it will be necessary to provide proof of purchase for any service calls under warranty.

Customers with FROSTMED under warranty should ensure that the above checks have been made.

If the malfunction is not a mechanical or electrical breakdown, a service fee will be assessed in order to repair FROSTMED.

7. Service and Spare Parts

In the event that the FROSTMED device requires service, or if you wish to purchase spare parts, please contact your local FROSTMED Customer Service Agent.

FROSTMED Customer Service Department

1441 Broadway, Suite 5011

New York, NY 10018

Ph: 646.569.5785

* calls may be recorded for quality purposes.

8. Customer Service Department

For general enquiries or further information on FROSTMED products please contact:

General Inquiries: info@frostmed.com

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