

User Manual

Manual
This operation manual allows you to work safely with the otoscope camera. Please keep this manual for your information. Liability is not for any damage caused by not following this manual.

Maintenance
To clean this device, use a damp cloth or mild water based cleaner. Fully charge the battery before you leave this device for long time to prevent damage to the battery.

Storage and transport
Proper transport can harm the device. Always transport this device in the provided box in order to prevent damage.

Information on disposal
This electronic equipment does not belong into domestic waste, but must be disposed in accordance with the applicable statutory provisions.

General Information

Product Detail

Operation procedure

- Search for "SR-WSP" on APP Store or Google Play Store, download and install the app. Or Scan QR code Download APP.
- Select the spectrum or ear-pick to install as needed.
- Long press the power button to turn on the device.

- Enter the Wi-Fi setting menu of the phone to find and connect with the Wi-Fi named SRnet-Otoscope.
- Open the app across the Wi-Fi connected. Observe the condition within the ear through the otoscope and APP's user interface.

Description

This otoscope camera is a portable, handheld video camera with short focal length, it's compatible with iPhone, iPad, Android devices and Tablet. The high resolution camera provides great video quality, there are 6 LED lights on the tip of the camera which makes the light path and increase the visual field of view.

Technical Data

Description	Data
Network Standard	IEEE 802.11 b/g/n
Antenna	Wave Antenna
Operating Frequency	2.4GHz
Data Rate	802.11a/b/g/n up to 150Mbps
Camera Diameter	3.8mm
Viewing angle	70°
Depth of Field	30mm-30mm
Camera Light Source	6LEDs
Power Input	5V/1A
Battery Type	800mAh Lithium battery
Battery life	Shows

➢ If you cannot turn on the otoscope, or it turns off after a few minutes of working, scan the light of the battery. Please make sure the battery is fully charged.

➢ In order to keep the image quality clear enough and maintain the optimal health & personal hygiene, please follow up these steps:

- Please use the provided cleaning cloth to sanitize the tip of the camera before and after each use. This can prevent ear wax to stick on the tip of the camera and cause blurry images.
- Please keep a focal length of about 3 inch to achieve the best quality image taking the scope for far from your food path. May result in blurry images.
- The ear canal is not perfectly straight, gently adjust the position of the otoscope to view your ear drum clearly.

Warning

- You must read the instruction manual carefully before using.
- You must put the protection cover onto the otoscope before using, then sanitize all the areas that contact the skin.
- Keep out of reach of children. Children are not allowed to use this product.
- Please take great care when using ear picks to clean ear wax. DO NOT force them too far into the ear canal.
- Stop using the otoscope if it overheats.
- When using this product on areas of people around you, be careful that nobody bumps into your hand while holding the device in your ear.
- Stop using this product immediately if you experience any discomfort.
- This product can not replace a doctor or medical provider's diagnosis or treatment. Contact your doctor immediately if you experience any discomfort.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the device.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) the device must accept any interference received, including interference that may cause undesired operation.

IMPORTANT NOTE

Before this equipment has been tested and found to comply with the limits for a Class II digital device, pursuant to part 15 of the FCC Rules, these limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio-frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio-communication devices which can be identified by the user, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into a different electrical outlet than that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC Radiation Exposure Statement:
This equipment complies with FCC radiation exposure limits set forth for an authorized environment.

➢ When connecting your equipment to a power source:

- Check the voltage rating before you connect this equipment to an electrical outlet.

Power Safety

- Use only the AC adapter approved for use with this device.
- Plug the AC adapter and this equipment in grounded area, not to stick tip or on the floor when you use it to charge the battery.
- The AC adapter and this equipment may become hot during charge operation.