

Yoobao



Visual Earwax Removal Tool

USER MANUAL

Thank you for your purchasing. For accurate and safe use, please read the manual carefully before use the product and keep it well.

The illustrations used in the manual only for reference, please refer to the actual product.

Product Introduction

This product is a personal care cleaning product, powered by a built-in lithium battery, and the camera is driven by a rotary switch control circuit. The image is transmitted to the mobile app via WIFI to display the situation in the ear canal in real time, which is convenient for individuals to clean the ear canal. Compared with the previous blind digging, the visual ear scoop is more convenient and safer.


Package List:

1. Visual Earwax Removal Tool
2. User Manual
3. Charging Cable



Visual Earwax Removal Tool User Manual Charging Cable

Product Structure



Functions

1. Turn the switch to the right in the off state to turn it on. When the power is turned on, the white indicator light flashes, and the white light of the camera fill light is always on.
2. The WIFI connection is successful, the power indicator white light is always on, put the product into the ear canal and it can be used.
3. When the battery remains 20%, the power indicator flashes red, the APP pops up a "charging reminder" pop-up window, plug in the power cord, and the low-battery reminder disappears; Charge the product

while it is off, the red light of the power indicator is on; when it is fully charged, the red light is off.

4. Turn the "rotary switch" to the left to power off.

Operation

1. Two Ways to download APP(Support Android 5.1 and IOS 9.0 version or above)


- ① Download the APP from AppStore or Google Play;
- ② Scan the QR code to download the APP. (If use WeChat to scan, need to manually jump to the browser to open)



IOS and Android system Compatible for both

2. Connect Product WiFi

- ① Turn on the product and APP
- ② Enter the WIFI setting
- ③ Connect phone with product WIFI "UDU earpick-xxxx".



3. Open App

Open App and click "Start up", then use it.



Tips:

- ① Disinfect the ear spoon with alcohol cotton pads before using the product;
- ② After use, remove the ear spoon in time and clean it;
- ③ Not recommended for inflamed ear canals.

About charging

- ① Please fully charge the product before first time use.
- ② Please use standard qualified DC5V adapter for charging.
- ③ The red indicator light is always on during charging.
- ④ The red indicator light off when it is fully charged.
- ⑤ The red indicator light flashes slowly when it's in low battery power, and please fully recharge before next time use.

※ When charging doesn't start after 10 seconds, please check if the Type-C cable port is fully plugged in or any poor connection issue.

Device Cleaning and Maintenance

- ① When cleaning, make sure the device is turned off avoiding electrical shock.
- ② Please take off the cleaning head from the instrument and rinse by warm water.
- ③ The instrument body is not waterproof, please wipe it clean with a wet tissue or cloth after each use. Do not rinse the body with water.

- ④ Clean the camera lens with a cotton swab, and avoid hard objects collide with the lens.
- ⑤ It is forbidden to wipe the body with corrosive liquid.
- ⑥ Please store the device in a dry indoor environment.
- ⑦ It is strictly forbidden to disassemble the instrument by yourself. If there occurs any failure or malfunction, please stop using and contact online customer service immediately.

Specifications

Product Name: Visual Earwax Removal Tool
 Model Number: DZ2102
 Input: DC 5V=1A
 Power: about 1W
 Battery capacity: 250mAh
 Working time: about 60 mins
 Product size: 144×14.2×14.2mm
 Executive standard: GB 4706.1-2005
 GB 4706.15-2008

Hazardous Substance

Part Name	Hazardous Substance					
	Lead (pb)	Mercury (Hg)	Cadmium (Cd)	Hexavalent chromium (Cr(VI))	Polybrominated biphenyls (PBB)	Polybrominated diphenylethers (PBDE)
Shell	○	○	○	○	○	○
PCBA	✗	○	○	○	○	○
Battery	○	○	○	○	○	○
LED	○	○	○	○	○	○
Cleaning Head	○	○	○	○	○	○
Type-C Charging Cable	○	○	○	○	○	○
Camera Module	○	○	○	○	○	○
Packaging Material	○	○	○	○	○	○

This table is compiled with the regulations of SJ/T 136 4
 O: Indicates the content of the hazardous substance is below the limitation specified in GB/T 26572.
 X: Indicates that the content of the hazardous substance exceeds the limitation specified in GB/T 26572.
 Note: Components that contain toxic, harmful or elemental elements have no substitutes due to the limitation of the global development level.

Warranty Card

Please keep this warranty card properly. When the product need to be repaired during the warranty period, please fill in the details, and send the warranty card together with the product.

Dear users:
 Nice day!
 Thank you for trusting and supporting our product. Warranty card is attached to each product, please keep it properly for product warranty, we won't reissue if it's lost. Thank you!

User's Name: _____ Phone Number: _____

Address&Postal Codes: _____

Shop's Name: _____ Order Number: _____

Cause of Malfunction: _____

Order Number: _____

Maintenance Record:

Date	Description of Malfunction	Maintenance Point	Result

Dongguan YOQBAO Communication Devices Co., LTD
 Address: Building 7, Huayuan Street, Changlong Village,
 Huangjiang Town, Dongguan, Guangdong, China
 Email:support.us@yoobao.cn
 HQ:1389-V1.0

Visual Earwax Removal Tool
 DZ2102

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction