



by pretty smart labs

Owner's Manual

Important Safety Instructions

THESE ARE THE INSTRUCTIONS FOR THE USE AND MAINTENANCE OF THE PRODUCT. ALL SAFETY AND OPERATION INSTRUCTIONS MUST BE READ BEFORE OPERATING THIS PRODUCT:

1. BEFORE OPERATING THIS PRODUCT, CAREFULLY READ ALL PROVIDED SAFETY AND OPERATION INSTRUCTIONS.
2. BE SURE TO KEEP THESE SAFETY AND OPERATION INSTRUCTIONS FOR FUTURE REFERENCE.
3. ALL WARNINGS ON THE APPLIANCE AND IN THE OPERATION INSTRUCTIONS MUST BE FOLLOWED. PAY ATTENTION TO ALL WARNINGS TO ENSURE YOUR SAFETY.
4. FOLLOW ALL PROVIDED USAGE AND OPERATION INSTRUCTIONS FOR THE PRODUCT. DO NOT SKIP ANY IMPORTANT STEPS OR PROCEDURES.
5. DO NOT USE THIS PRODUCT NEAR WATER OR MOISTURE. ENSURE THAT THE PRODUCT IS OPERATED IN A DRY AND SAFE LOCATION. THE SCANNER IS RATED IP65 - WATER RESISTANT.
6. USE ONLY A DRY CLOTH TO CLEAN THE PRODUCT. DO NOT USE WATER OR OTHER LIQUIDS.
7. DO NOT INSTALL THE PRODUCT NEAR HEAT SOURCES SUCH AS RADIATORS, STOVES, ETC. AND ENSURE THAT THERE IS SUFFICIENT VENTILATION TO PREVENT OVERHEATING.
8. PROTECT THE POWER CORD TO PREVENT IT FROM BEING STEPPED ON OR PINCHED, ESPECIALLY AT THE PLUGS, OUTLETS, AND WHERE THEY EXIT THE APPLIANCE.
9. USE ONLY ACCESSORIES SPECIFIED BY THE MANUFACTURER TO AVOID DAMAGE TO THE PRODUCT.
10. DURING ELECTRICAL STORMS OR WHEN THE PRODUCT IS NOT IN USE FOR LONG PERIODS OF TIME, UNPLUG THE PRODUCT TO PROTECT IT.

11. IF THE PRODUCT HAS BEEN DAMAGED IN ANY WAY, SUCH AS WHEN THE POWER CORD OR PLUG IS DAMAGED, LIQUID HAS BEEN SPILLED OR OBJECTS HAVE FALLEN INTO THE APPLIANCE, IT HAS BEEN EXPOSED TO RAIN OR MOISTURE, DOES NOT FUNCTION NORMALLY, OR HAS BEEN DROPPED AND DAMAGED, REFER ALL MAINTENANCE TO QUALIFIED SERVICE PERSONNEL.
12. DO NOT PLACE NAKED FLAME SOURCES, SUCH AS LIT CANDLES, ON THE PRODUCT.
13. PAY ATTENTION TO ENVIRONMENTAL ASPECTS OF BATTERY DISPOSAL.
14. THIS PRODUCT IS DESIGNED FOR USE IN MODERATE CLIMATES. DO NOT USE IT IN EXTREME CONDITIONS.
15. WARNING: TO REDUCE THE RISK OF FIRE OR ELECTRIC SHOCK, DO NOT EXPOSE THIS PRODUCT TO RAIN OR MOISTURE.
16. BATTERIES (BATTERY PACK OR INSTALLED BATTERIES) SHOULD NOT BE EXPOSED TO EXCESSIVE HEAT.

Table of Contents

INTRODUCTION	01
<hr/>	
WHAT'S IN THE BOX?	01
<hr/>	
SOFTWARE REQUIREMENTS	01
<hr/>	
SCANNER TUTORIAL	01
<hr/>	
SCANNER GENERAL DESCRIPCION	03
<hr/>	
PRODUCT SPECIFICATIONS	03
<hr/>	

INTRODUCTION

CONGRATULATIONS ON YOUR PURCHASE. THIS MANUAL CONTAINS INFORMATION ABOUT SCANNERS. WE ENCOURAGE YOU TO TAKE A FEW MINUTES TO READ THIS MANUAL, WHICH DESCRIBES THE PRODUCT AND STEP-BY-STEP INSTRUCTIONS FOR USING THE PRODUCT AND GETTING STARTED WITH IT. IT IS NECESSARY TO READ AND UNDERSTAND ALL INSTRUCTIONS BEFORE USING THE PRODUCT.

IF YOU HAVE ANY QUESTIONS ABOUT THESE PRODUCTS, THEIR INSTALLATION OR OPERATION, PLEASE CONTACT YOUR DISTRIBUTOR OR THEIR CUSTOMER SERVICE, OR SEND A MAIL: INFO@PRETTYSMARTLABS.COM

WHAT'S IN THE BOX?

1X HOOFLINK SCANNER

1X CABLE USB A TO USB C

1X HANDLE SLEEVE


SOFTWARE REQUIREMENTS

REQUIRES IOS 15.0 OR LATER.

REQUIRES ANDROID 10.0 OR LATER.

SCANNER TUTORIAL

CLICK THE SCANNER'S BUTTON AND WAIT UNTIL THE SCREEN TURNS ON. YOU'LL SEE THE SCREEN DISPLAY "STANDBY".




FIRST, YOU'LL HAVE TO LINK YOUR SCANNER WITH THE HOOFLINK APP. TO DO THAT, OPEN YOUR APP, CLICK THE TOP LEFT ICON TO OPEN THE MENU AND SELECT "MY SCANNERS". ONCE THERE, DOUBLE-CLICK THE SCANNER'S BUTTON TO ENTER TO THE LINKING MODE. A POPUP WILL APPEAR, SHOWING THE SCANNER'S ID. VERIFY THAT IT CORRESPONDS TO THE ID ON THE SCANNER'S SCREEN. CLICK "LINK". A NEW SCREEN WILL APPEAR, WHERE YOU CAN ASSIGN A NAME TO YOUR SCANNER AND A LOCATION (E.G. A BARN).

NOW YOU CAN START USING IT! IN ORDER TO SCAN A CHIP, CLICK THE SCANNER'S BUTTON ONCE AND IT WILL ENTER THE "SCANNING" MODE. THEN, BRING THE SCANNER CLOSE TO THE HORSE'S LIP OR NECK. WHEN IT FINDS THE CHIP YOU WILL HEAR A BEEP NOISE, AND THE SCANNER'S SCREEN WILL SHOW "FOUND ID".




WHEN YOU SEE THIS ICON IN THE APP, IT MEANS THAT YOU NEED TO SCAN A CHIP IN ORDER TO PROCEED.




TO GO FROM "SCANNING" OR "FOUND ID" BACK TO "STANDBY", SIMPLY PRESS THE SCANNER'S BUTTON ONCE.



AND TO MANUALLY EXIT "LINKING" MODE, PRESS THE SCANNER'S BUTTON TWICE.



SCANNER GENERAL DESCRIPCION



PRODUCT SPECIFICATIONS

GENERAL SPECIFICATION: DIMENSIONS

- LCD DISPLAY
- RFID READER
- BLE ANTENNA
- USB C INPUT
- SCANNER DIMENSIONS (W X H X D):
113.6 X 255.69 X 28.5 MM / 4.47 X 10.06 X 1.12 IN
- PACKAGING DIMENSIONS (W X H X D):
138 X 269 X 38.5 MM / 5.43 X 10.59 X 1.51 IN

FCC REGULATORY COMPLIANCE

CAUTION: CHANGES OR MODIFICATIONS NOT EXPRESSLY APPROVED BY THE PARTY RESPONSIBLE FOR COMPLIANCE COULD VOID THE USER'S AUTHORITY TO OPERATE THE EQUIPMENT.

NOTE: THIS EQUIPMENT HAS BEEN TESTED AND FOUND TO COMPLY WITH THE LIMITS FOR A CLASS B DIGITAL DEVICE, PURSUANT TO PART 15 OF THE FCC RULES. THESE LIMITS ARE DESIGNED TO PROVIDE REASONABLE PROTECTION AGAINST HARMFUL INTERFERENCE IN A RESIDENTIAL INSTALLATION. THIS EQUIPMENT GENERATES, USES AND CAN RADIATE RADIO FREQUENCY ENERGY AND, IF NOT INSTALLED AND USED IN ACCORDANCE WITH THE INSTRUCTIONS, MAY CAUSE HARMFUL INTERFERENCE TO RADIO COMMUNICATIONS. HOWEVER, THERE IS NO GUARANTEE THAT INTERFERENCE WILL NOT OCCUR IN A PARTICULAR INSTALLATION.

IF THIS EQUIPMENT DOES CAUSE HARMFUL INTERFERENCE TO RADIO OR TELEVISION RECEPTION, WHICH CAN BE DETERMINED BY TURNING THE EQUIPMENT OFF AND ON, THE USER IS ENCOURAGED TO TRY TO CORRECT THE INTERFERENCE BY ONE OR MORE OF THE FOLLOWING MEASURES:

- REORIENT OR RELOCATE THE RECEIVING ANTENNA.
- INCREASE THE SEPARATION BETWEEN THE EQUIPMENT AND RECEIVER.
- CONNECT THE EQUIPMENT INTO AN OUTLET ON A CIRCUIT DIFFERENT FROM THAT TO WHICH THE RECEIVER IS CONNECTED.
- CONSULT THE DEALER OR AN EXPERIENCED RADIO/TV TECHNICIAN FOR HELP.

THIS DEVICE IS DESIGNED TO MEET THE REQUIREMENTS FOR EXPOSURE TO RADIO WAVES ESTABLISHED BY THE FCC. THESE REQUIREMENTS SET A SAR LIMIT OF 1.6 W/KG AVERAGED OVER 1 GRAM OF TISSUE FOR USE ON THE BODY AND 4.0 W/KG AVERAGED OVER 10 GRAMS OF TISSUE FOR USE ON THE LIMBS. THE HIGHEST SAR VALUE REPORTED UNDER THIS STANDARD DURING PRODUCT CERTIFICATION FOR USE WHEN PROPERLY WORN ON THE BODY AND LIMBS WITH NO SEPARATION.

SUPPLIER'S DECLARATION OF CONFORMITY

UNIQUE IDENTIFIER

TRADE NAME: HOOF LINK SCANNER
MODEL NO.: HL-SC-100-23

RESPONSIBLE PARTY – U.S. CONTACT INFORMATION

US COMPANY NAME: TU MUNDO INC, DBA PRETTY SMART LABS
ADDRESS: 3704 TIMSON CT AUSTIN TEXAS 78731 UNITED STATES
INTERNET CONTACT INFORMATION: INFO@PRETTYSMARTLABS.COM

