

User Manual

Smart earpick

Appdownload method



Please use your mobile browser to scan the QR code to download

YW tool

(Android Google)
(Android China)
(IOS app store)

Apple users

Select the **IOS app Store** option after Scanning the code. Open the **App Store** To download and install. You can also directly Search for **YW-tool** Downloads in the **App Store**.

Android users


After scanning the code, select the **Android Google** option to enter the next page. Click: **Download and install** through third-party droser option s.

YW tool


(Android Google)
(Android China)
(IOS app store)

Connection process


1. When receiving the product, First ensure that the battery is sufficient. Power button for 3 seconds to turn it on, Blue light will light up.
2. Open the phone settings WIFI Function, Connect to the YW tool-xxxx hotspot



3. Open the APP interface click the **Start** button to start using the product.



4. When you encounter this prompt, please do not switch WIFI, just click to **use**.



Product parameters

Product name: Smart earpick

weight: 14gram
Image sensors: CMOS
Battery capacity: 130mAh
Duration of use: 60 minutes/time
Input current: DC5V300mA
Lens pixels: Ultra HD lenses
Connect to a hotspot: YW tool-xxxxxx
Network standards: IEEE802.11b/g/n
Optimal distance: 1.2-1.6CM
Operating frequency: 2.4GHz
Operating temperature: -5-40°



Note

1. Wipe the lens with a cotton swab when the image is not clear.
2. When cleaning the lens, please use a professional alcohol cotton cloth and wipe it carefully
3. Please observe the surroundings before use, and do not use it when someone is running to avoid impact.
4. This product is not suitable for children under 5 years of age.
5. Children are not allowed to use alone, and must be guided by an adult.
6. Please keep the product out of the reach of children.
7. Do not use the endoscope in liquids to avoid damage caused by immersion in water.
8. Lens heating is a normal phenomenon.
9. This product is intended for use as a personal care product only and is not intended for use as a medical device.
10. Please place the product in a cool and dry place, away from fire and corrosive materials.

Charging instructions

1. Please turn off the product while charging.
2. When charging, the indicator light is shown in red.
3. Charging is complete and the red light goes out

Troubleshooting

1. The screen is stuck
 - A. The data transmission distance is too far, and it is best not to use the product too far away from the mobile phone.
 - B. If the power is insufficient, it will also cause lag, please charge it in time.
2. WIFI signal cannot be transmitted/intermittent
 - A. The lower the battery, the weaker the WIFI signal, use it after it is fully charged.
 - B. When the WIFI signal transmission distance is exceeded, the product and the mobile phone should be placed together.
 - C. There are too many WIFI signals around, and the signals that affect the product should be used in areas with less WIFI signals as much as possible.
 - D. Make sure the mobile phone is connected to the product WIFI.
3. No image display.
 - A. Confirm whether the mobile phone is connected to the WIFI hotspot of the product.
 - B. After confirming the connection, if there is no display, close the APP and reopen the product APP.

After-sales protection card

Warranty Notice

In accordance with the relevant regulations of the national quality supervision and inspection "three guarantees" service place, the company will provide you with a high-quality service commitment of 7 days for return, 15 days for replacement and one-year warranty in the case of non-human damage.

Man-made damage includes: obvious falls, scratches, paint loss, defects, water or other liquids, charging with a charging interface with a non-standard voltage, resulting in the unit burning out, disassembling and damaging the machine without permission, incomplete accessories, etc.

7 days package: 7 days from the date of receipt of the product, you are in the case of normal use, if there is a product quality problem. You have the option to return the product (you need to return the full set of products)

15-day replacement: From the date of receipt of the product, within the 8th day to 25 days of normal use, if there is a product quality problem, you can choose to replace or repair. After your replacement or repair, you are entitled to the warranty period calculated from the date of purchase of the original product.

One-year warranty: From the date of receipt of the product, within one year under your normal use, if there is a quality problem, we will provide you with free repair service.

User Information: _____

User's name: _____ Contact number: _____

Recipient address: _____

Selling merchants: _____ Merchant phone: _____

FCC Caution:

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception,

which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.