

---

# How to use The wearable device Measure heart rate and falls and send text messages

---

---

# LIST

---

**Wearable Devices  
(Heart Rate Care)  
How to link – Android**

---

---

**Wearable Devices  
(Heart Rate Care)  
Device Description**

---

---

# Wearable Devices (Smart Band) How to link - Android

---





## How to link - Android

When using a wearable device, you must download the application “Heart Rate Care” from the Google app store to use it.



## How to link - Android


The wearable device and mobile phone must be connected via Bluetooth for the heart rate care application to work properly.



## How to link - Android


When linked to a wearable, data from the wearable device can be retrieved directly from Heart Rate Care.

Depending on your settings, you can send health data via text message to your saved phone number.



## How to link - Android (User management)

After connecting the wearable device and mobile phone, Please access Heart Rate Care to add and configure users.




---

# Wearable Devices (Heart Rate Care) Device Description

---

## Device Description

If possible, wear the wearable device on your left wrist with the ECG measurement part facing upward.



ECG(electrocardiogram meter)

PPG(photoplethysmography)

Wrist strap connection part

## Device Description

Wearable devices can be turned on/off by pressing the power button. And please charge it with type C (USB-C).



## Device Description

<b>COMPANY</b>	<b>WEARCOM INC</b>
<b>Name of the device</b>	<b>ECG and fall management text messaging service device</b>
<b>Model No.</b>	<b>WEARCOM2301</b>
<b>Certification Number</b>	<b>GC220490</b>
<b>Manufacturer/Country of Manufacture</b>	<b>WEARCOM INC / KOREA</b>
<b>RF specification</b>	<b>BLE</b>
<b>Electrical standards</b>	<b>Charge 5 Vdc , Battery 3.7 Vdc</b>

---

## FCC Information to User

---

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

### Caution

THE GRANTEE IS NOT RESPONSIBLE FOR ANY CHANGES OR MODIFICATIONS NOT EXPRESSLY APPROVED BY THE PARTY RESPONSIBLE FOR COMPLIANCE.

SUCH MODIFICATIONS COULD VOID THE USER'S AUTHORITY TO OPERATE THE EQUIPMENT.

IMPORTANT NOTE : FCC RF Radiation Exposure Statement

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment and meets the FCC radio frequency (RF) Exposure Guidelines of the Exposure rules. This equipment has very low levels of RF energy that is deemed to comply without testing of specific absorption rate (SAR).

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

