



### 1. Install two AAA batteries for the first time

Press and hold the battery compartment cover, pull it downward, put in two AAA batteries, and put back the battery compartment cover.

## 2. Pairing

- i. On the remote, push **Power**.
- ii. On the headset, push **Power** and **Option** simultaneously.
- iii. Using the headset touchpad, select **Enable Bluetooth pairing**.
- iv. From the list, select **eSight Controller**.

### Button Functions:

**Power Button:** Short press to turn on the remote.

**Menu Button:** Short press to open/close the menu.

**Back Button:** Short press to go back in a menu or exit a view.

**Spot Zoom Button:** Hold to temporarily zoom out / release to return to the previous zoom level.

**Launcher Button:** Push to open a console to access more features, like Gallery, eReader, or Support, or access eSight system settings.

**Up / Down Arrows** – In live camera view, push to Zoom In (Up) or Zoom Out (Down). In menu view, push to browse up or down.

**Left / Right Arrows** – In live camera view, push to decrease contrast (Left) and increase contrast (Right). In menu view, push to browse left or right.

**Center Button** – In live camera view, push to take a photo. In menu view, push to select an option.

#### FCC Warning Statement

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- - Reorient or relocate the receiving antenna.
- - Increase the separation between the equipment and receiver.
- - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- - Consult the dealer or an experienced radio/TV technician for help.