

Performing System Setup Functions

Use the Setup menu to perform the following functions:

- Check the battery power level
- Put the transmitter into low power standby mode when the transmitter will not be used for temporary periods
- Configuring the Guardian

Checking the Battery

To check the battery power level:

1. From the Main menu press **More**.

2. From the More menu press **Setup**
3. From the Setup menu press **Bat.** The battery status will display as depicted in Figure 5-1.




Figure 5-1. *Battery Status Display*

There are 3 states that can be displayed under battery status. They are:

High Level	Battery voltage above 7.8 volts.
Mid Level	Battery voltage between 6.0 volts and 7.8 volts.
Low Level	Battery voltage below 6.0 volts.

Replace the battery when Battery Status indicates “low battery”.

Putting the Transmitter into Standby Mode

The Guardian transmitter can be put into low power standby mode for temporary storage. For long term storage, remove the battery to protect the Guardian unit. Selecting the standby mode saves up to 60% of the power consumption by turning off the SpO₂, LCD display and transmitter. However, you cannot put the transmitter into standby mode when it is connected to ECG or has a valid SpO₂ signal.

To put the unit into standby mode:

1. Select **More** from the Main menu.
2. Select **Setup**.
3. Select **Off**. The TX Standby menu will be displayed as depicted in Figure 5-2.



Yes will be flashing. To confirm standby, select **Yes**. If you change your mind about selecting standby mode, use the **arrow** key to select **No**. After making a selection, you will return to the Setup menu.

If the transmitter is in standby mode, press and hold the Menu button for 1 second to resume normal operation.

Configuring the Guardian

The Guardian allows you to set certain default values that take effect when the unit is powered on. You can configure transmitter frequency, display mode, and measurement defaults for SpO₂, ECG, and Temperature.

Keep in mind that these functions require that you enter the Conf menu system. Any time you select the Conf option from the menu the key pad will lock. Press and hold the arrow key for 3 seconds to unlock the keypad before proceeding.

Changing Transmitter Frequency

To change the transmitter frequency:

1. From the Main menu select **More**.
2. Select **Setup**.
3. From the Setup menu select **Conf**.
4. From the Conf Setup menu select **Channel**. The Channel menu will be displayed as depicted in Figure 5-3.




Figure 5-3. *Channel Menu*

The current channel setting is displayed in the top row of the display. To change the setting use the right arrow key to highlight either the up or down arrows on the display, depending on whether you want to increment or decrement the setting. Once the arrow you want to use is flashing, press **Enter**. And then use the **Enter** key to change the setting. The channel will change one number at a time when you press and release the **Enter** key. It will change by 10 if you press and hold the **Enter** key. When the desired channel is set, select **Save** and press **Enter**. Once you save the setting, the Main menu will be displayed.

Configuring the Display mode

You can configure the display mode of the Guardian to either operate continuously or to cycle down 90 seconds after the last key is pressed.

Do the following to configure display mode:

1. From the Main menu select **More**.
2. Select **Setup**.
3. From the Setup menu select **Conf**.
4. From the Conf Setup menu select **More**.
5. Select **Disp**. The Display menu will appear as depicted in Figure 5-4.

DSPLAY MENU	
CONT	CYCLE*

Figure 5-4. *Display Menu*

6. Highlight the selected option, either **Cont** or **Cycle** and press Enter. If you put the display into cycle mode, press the Menu key to bring the display up after it has powered off.

Setting Measurement Defaults

The Guardian's default settings are used at power-up. They will appear marked with an asterisk (*) as you navigate through the menus. You can change the default settings for SpO2, ECG, and Temperature from the unit's Default Setup menu.

To access this menu:

1. From Main menu select **More**.
2. Select **Setup**.
3. From Setup menu select **Conf**.

4. From Conf Setup select **More**.
5. Select **Default**.
6. The Default setup menu will appear, as depicted in Figure 5-5.

Figure 5-5. *Default Setup Menu*

Setting SpO₂ Defaults

You can use the SpO₂ default function to set the measurement interval, the filter mode, and parameter off/on mode.

DEFLT SETUP
SPO2 ECG TMP

measurement interval, the filter mode,

1. Display the **Default Setup** menu as described in the section “Setting Measurement Defaults”.
2. Select **SpO₂** . The SpO₂ default menu will display. It is illustrated in Figure 5-6.

SP02 DEFLT
INTERVL MORE

Figure 5-6. *SpO₂ Default Menu*

3. To set the default for the SpO₂ Interval, select **Intervl**. The SpO₂ Interval Default menu will display as depicted in Figure 5-7.

O2	INT /MIN	
CONT *	1	2

Figure 5-7. *SpO₂ Interval Default Menu*

Select Cont to have the Guardian
Guardian monitor at one or two n

SpO ₂	DEFLT
MODE	ON/OFF

Select 1 or 2 to have the

For more information on the characteristics of each interval type, refer to Part 2 of this manual.

4. To set the default for the SpO₂ filter mode, select More from the SpO₂ default menu (see Figure 5-6).
The second page of the SpO₂ Default menu will display. It is depicted in Figure 5-8.

Figure 5-8. *SpO₂ Default Menu*

Select Mode to display the SpO₂ Mode Default menu as depicted in Figure 5-9.

SpO ₂	MODE
NORMAL	SLOW

Figure 5-9. *SpO₂ Mode Default Menu*

Select either Normal or Slow as the default mode. For more information on the characteristics of these settings, refer to Part 3 of this guide.

5. To set the On/Off default for the SpO₂ parameter at power-up, select More from the SpO₂ default menu (see Figure 5-6). The SpO₂ Default menu will display. Refer to Figure 5-8 for an example.

Select On/Off to display the O₂ Off? Confirmation menu. It is depicted in Figure 5-10.

Figure 5-10. *SpO₂ Off? Confirmation Menu*

O ₂	OFF?
YES	NO *

ECG Defaults

Follow the steps in the section “

y the Default Setup menu as

To set ECG Defaults, do the following:

1. Select ECG from the Default Setup menu.

2. The ECG Default menu will display:




Figure 5-11. ECG Default Menu

To set the default lead selection, choose **3LD DFLT**. The 3LD DFLT Menu will be displayed as depicted in Figure 5-12.




Figure 5-12. 3LD DFLT Menu


Select the desired default lead. Refer to part 4 of this guide for specific information on the characteristics of the ECG setting options.

The Guardian transmitter can be used as a stand alone SpO2 device and can be set so that ECG information is not sent to the receiver, as long as there is no ECG cable connected to the patient. If an ECG signal is detected, ECG will be turned off again.

To turn ECG off at power-up or if an ECG signal is detected, the ECG Off? Confirmation menu will display as depicted in Figure 5-13. The ECG Off? Confirmation menu will display as depicted in Figure 5-13.

Figure 5-13. ECG Off? Confirmation Menu

Select **Yes** and press **Enter** to set a default. Once you have made a selection, you will see the menu. Press **Enter** to cancel the operation.



Setting Temperature Display Defaults

You can set a default for the mode in which the Guardian displays temperature.

1. Follow the steps in the section “Setting Measurement Defaults” to display the Default Setup menu as depicted in Figure 5-5.
2. Select TMP to display the TMP Default Menu, as depicted in Figure 5-14.

Figure 5-14. *TMP Default Menu*

3. Select the desired measurement mode and press **Enter**. You will return to the Default Setup menu. Press the Menu key until you return to the Main menu.

