DEX™

A WIRELESS WORLD

USER MANUAL


PHONE-DEX

MODEL: PHONE-DEX USC

WIDEX®
HIGH DEFINITION HEARING

Table of contents

Handset layout	1
Information	2
Parts checklist	2
Installation	2
Connection	2
Battery	3
Range	4
The display	5
Telephone settings	6
Using the handset menu	6
Turn on the telephone for the first time	6
Display settings	7
Sound settings	7
Call settings	8
Operation	9
Calling	9
Directory	10
Caller ID	12
Intercom	13
Find handset	13
Add and register additional handset	13
To join a call	14
Intercom (two or more handsets)	14
Troubleshooting	15
IMPORTANT SAFETY INFORMATION	16
Equipment Approval Information	18
Disclaimer and limitation of liability	21
Warranty and other information	21
Support	23

Handset layout

1. Charge light
2. ▲ scroll up, increase volume,  directory
3. **MUTE / REMOVE** entered digits or characters, or delete entry
4.  / **FLASH**
5. Dialing keypad
6. * switch to tone dialing temporarily when using the pulse service
7.  speakerphone
8. **REDIAL** / press and hold **PAUSE** to insert a dialing pause
9. **INT** Intercom
10. # dialing options for call log
11. ▼ scroll down, decrease volume,  call log
12.  / **CLEAR** silence the ringer temporarily or exit menu
13. **MENU / SELECT** show menu or save an entry or setting

Information

Before using this telephone, please read *“IMPORTANT SAFETY INFORMATION”* on pages 16-17.

The PHONE-DEX is made for use in USA and Canada only.

The PHONE-DEX is a cordless phone that streams sound directly to your hearing aids.

It is incredibly easy to use – the only thing you need to do is hold the phone normally against the ear using the hearing aid.

The PHONE-DEX works as a conventional phone as well. So family and friends can use it.

See *“Operation”* on page 9 for more information.

Parts checklist


Your package contains the following items. Save your sales receipt and original packaging in case it is necessary to ship your telephone for warranty service.

- Battery compartment cover
- Battery
- Telephone line cord
- Telephone base
- Power adapter
- Handset
- User’s manual

Installation


Connection

1. Connect the telephone line cord to the telephone wall jack and to the telephone base (☎).
Note! If you subscribe to digital subscriber line (DSL) high-speed Internet service through your telephone line, you may need to install a DSL filter between the telephone line cord and the telephone wall jack. Contact your DSL service provider for more information.

2. Connect the telephone base power adapter to an electrical outlet not controlled by a wall switch and to the telephone base (⚡).

3. Insert the battery connector in accordance with the markings for the batteries (2.4V NiMH 500mAh) inside the handset battery compartment.
4. Replace the battery cover. Align the cover flat against the battery compartment, then slide it towards the center of the handset until it clicks into place.
5. Place the handset into the base unit with the keypad facing outwards to charge. The charge light is on when charging.

Note!

The power adapter is the disconnect device between the product and wall power. The wall socket outlet must be close to the equipment and easily accessible.

Battery

When the batteries are relatively new (less than 6 months) and fully-charged (10 hours) they will last for approx 100 hours standby time or approx 10 hours call time. Full battery capacity will not be reached until the batteries have been charged 4–5 times. It is not possible to overcharge or damage the batteries by charging them for "too long".

Warning!

Risk of explosion if battery is replaced by an incorrect type. Use only original type of batteries. The guarantee does not cover any damage caused by using the wrong batteries.

Note!

When the battery is low or when there is too great a distance to the base unit the telephone will be unable to connect to a line. Do not place the base unit close to other electrical equipment, in direct sunlight or other strong heat sources! This is to minimize the risk of any interference. It is normal for the device to become warm during charging and this is not dangerous. Dispose of used batteries according to the instructions.

Charging

To keep the battery charged, place the handset in the telephone base or telephone charger when not in use. The battery is fully charged after 10 hours of continuous charging.

If the screen is blank, you need to charge the handset without interruption for at least 30 minutes to give the handset enough charge to use the telephone for a short time.

When the batteries are running low a warning tone is heard and **LOW BATTERY** and  flashes on the display.

After you install your telephone or power returns following a power outage, date and time needs to be set. See “Date / Time” on page 7 for more information.

Range

This cordless telephone operates with the maximum power allowed by the Federal Communications Commission (FCC). Even so, this handset and telephone base can communicate over only a certain distance - which can vary with the locations of the telephone base and handset, the weather, and the layout of your home or office.

When the handset is out of range, the handset displays **OUT OF RANGE OR NO PWR AT BASE**. If there is a call while the handset is out of range, it might not ring, or if it does ring, the call might not connect well when you press  / **FLASH**. Move closer to the telephone base, then press  / **FLASH** to answer the call.

If the handset moves out of range during a telephone conversation, there might be interference. To improve reception, move closer to the telephone base. Coverage can be improved by turning the head so that the handset is facing the base unit.

The display

Indicator	Function
	Battery status. Flashes while charging
	Speakerphone mode activated
	Ringer volume Off
	Voicemail indicator. <i>Note! This function requires a subscription from your service provider!</i> See also “Clear voicemail indicators” on page 8
NEW	New Caller ID Note! <i>This function requires a subscription from your service provider!</i>
11:59 AM PM	Time. See also “Date / Time” on page 7
MUTE	Microphone is muted during a call
07/04	Date. See also “Date / Time” on page 7

Telephone settings

Using the handset menu

The settings described in this chapter are carried out in idle mode.

The handset menu contains feature menus with submenus. The selected menu is marked with >.

To enter the handset menu:

1. Press **MENU / SELECT**.
2. Press **▲ / ▼** to scroll the feature menu:

DIRECTORY (📖)	CALL LOG (📞)	RINGERS	SETTINGS	REGISTRATION
<ul style="list-style-type: none"> • REVIEW • ADD CONTACT • SPEED DIAL 	<ul style="list-style-type: none"> • REVIEW • DEL ALL CALLS 	<ul style="list-style-type: none"> • RINGER VOLUME • RINGER TONE 	<ul style="list-style-type: none"> • LCD LANGUAGE • CLR VOICEMAIL • KEY TONE • SET DATE / TIME • HOME AREA CODE • DIAL MODE 	

3. Press **MENU / SELECT** to enter a menu.
 - To return to the previous menu, press  / **CLEAR**.
 - To return to idle mode, press and hold  / **CLEAR**.

Turn on the telephone for the first time

The first time you turn on your telephone, you need to set date and time. See “Date / Time” on page 7 for more information.

Display settings

LCD language

1. Press **MENU / SELECT**.
2. Select **SETTINGS** → **LCD LANGUAGE**.
3. Select language and press **MENU / SELECT** to save.

Date / Time

1. Press **MENU / SELECT**.
2. Select **SETTINGS** → **SET DATE / TIME**.
3. Press the dialing keys (0-9) to enter the month (mm), date (dd) and year (yy). When done, press **MENU / SELECT**.
 - If needed, press **MUTE / REMOVE** to delete a digit. Press and hold **MUTE / REMOVE** to delete the entire entry.
4. Press the dialing keys (0-9) to enter hour and minute. Press ▲ / ▼ to choose **AM** or **PM**. When done, press **MENU / SELECT** to save.

Sound settings

Ringer volume

1. Press **MENU / SELECT**.
2. Select **RINGERS** → **RINGER VOLUME**:
 - Press ▲ to increase the volume (max 6).
 - Press ▼ to decrease the volume or to silence the ringer.
3. When done, press **MENU / SELECT** to save.

Tip!

You can also press MUTE / DELETE or CLEAR to temporarily silence the ringer when the handset is ringing.

Ringer tone

1. Press **MENU / SELECT**.
2. Select **RINGERS** → **RINGER TONE**. Make sure that the ringer volume is not off.
3. Press ▲ / ▼ to listen to the 10 different ringer tones.
4. When done, press **MENU / SELECT** to save.

Key tone

1. Press **MENU / SELECT**.
2. Select **SETTINGS** ➔ **KEY TONE**.
3. Select **ON** or **OFF** and press **MENU / SELECT** to save.

Call settings

Home area code

1. Press **MENU / SELECT**.
2. Select **SETTINGS** ➔ **HOME AREA CODE**.
3. Press the dialing keys (0-9) to enter home area code.
 - If needed, press **MUTE / REMOVE** to delete a digit. Press and hold **MUTE / REMOVE** to delete the entire entry.
4. When done, press **MENU / SELECT** to save.

Dial mode

1. Press **MENU / SELECT**.
2. Select **SETTINGS** ➔ **DIAL MODE**.
3. Select **TONE** or **PULSE** and press **MENU / SELECT** to save.

Clear voicemail indicators

If you subscribe to a voicemail service offered by your telephone service provider, **NEW VOICEMAIL** appears on the handset when you have new voice mail messages. After listened to the message, the indicators on the handset turns off automatically. You can manually turn off the indicator if you for example access your voicemail while away from home.

1. Press **MENU / SELECT**.
2. Select **SETTINGS** ➔ **CLR VOICEMAIL** and **TURN OFF INDICATOR?** is shown.
3. Press **MENU / SELECT** to save or press  / **CLEAR** to return to the previous menu.

Tip!

This only turns off the indicator, it does not delete your voicemail messages.

Operation

Your telephone is also specially adapted to work with wireless Widex hearing aids. This means that the sound from the telephone is transmitted directly into your hearing aid(s). When you place the telephone near one of your ears wearing a hearing aid, a 'dong' sound in the hearing aid will indicate that the function is activated. Hold the phone close to one ear using the hearing aid.

When the telephone is moved out of range, your hearing aid(s) will return to the listening program chosen as program number 1 (usually the Master program).

If the telephone sound is too weak or too loud, you can use the volume keys ▼ / ▲ to adjust the volume. If you need further information, please consult your hearing health care professional.

For users who do not wear hearing aids, the telephone works as otherwise described in this manual.

Calling

Make a call

- Press  / **FLASH** or  and dial the telephone number.
- Dial the telephone number and press  / **FLASH** or . If needed, delete number with **MUTE / REMOVE**.

Receive a call

When the handset rings, press  / **FLASH**,  or any dialing key (0-9, * or #).

End a call

Press  / **CLEAR** or put the handset in the telephone base.

Use the speakerphone

During a call, press  to switch between speakerphone and normal handset use.

Redial

The handset stores the last 10 telephone numbers dialed.

1. Press **REDIAL / PAUSE**.
2. Press ▼ / ▲ or **REDIAL / PAUSE** repeatedly until the desired entry is displayed.
3. Press  / **FLASH** or  to dial.

Call waiting

If you subscribe to a call waiting service from your telephone service provider, you hear an alert tone when there is an incoming call while you are on another call.

1. Press  / **FLASH** to put the current call on hold to take the new call.
2. Press  / **FLASH** to switch back.

Volume control

During a call, press  /  to adjust the listening volume.

Note!

The handset earpiece volume setting and speakerphone volume setting are independent.

Mute

1. During a call press **MUTE / REMOVE** to turn off your microphone. **MUTE** is displayed when activated.
2. Press **MUTE / REMOVE** to turn on your microphone. The screen temporarily displays **MICROPHONE ON**.

Temporary silent mode

When someone calls, you can press  / **CLEAR** or **MUTE / REMOVE** on the handset to turn off the handset ringsignal. **RINGER MUTED** is displayed. The next call rings automatically at the present volume.

Directory

You can save up to 50 entries in the directory. Only one handset can review the directory at a time. Any modifications made on one handset apply to all.

Save contact

1. Press **MENU / SELECT**.
2. Select **DIRECTORY** → **ADD CONTACT**.
3. Press the dialing keys (0-9) to enter the number (up to 30 digits). When done, press **MENU / SELECT**.
 - If needed, press **MUTE / REMOVE** to delete a digit. Press and hold **MUTE / REMOVE** to delete the entire entry.
4. Press the dialing keys (0-9) to enter name (up to 15 characters). Press the equivalent key several times for the wanted letter. When done, press **MENU / SELECT** to save.
 - If needed, press 0 to add a space.

Tip!

You can copy number from the redial list, when entering number press **REDIAL / PAUSE** and locate the number, then press **MENU / SELECT** to copy. Then press **MENU / SELECT** again to edit number and to enter name.

Make a call from the directory

1. Press  / .
2. Press  /  to scroll between contacts or quick search by pressing the corresponding number key for the first letter.
3. When the wanted contact appears, press  / **FLASH** or  to dial.

Edit number

1. Press  / .
2. Scroll to wanted contact and press **MENU / SELECT**.
3. When **EDIT NUMBER** is displayed, use the dialing keys (0-9) to edit. Press **MENU / SELECT** to continue.
4. When **EDIT NAME** is displayed, use the dialing keys (0-9) to edit.
5. Press **MENU / SELECT** to save.

Delete number

1. Scroll to wanted contact and press **MUTE / REMOVE**.
2. When the screen displays **DELETE CONTACT?**, press **MUTE / DELETE**.

Store a speed dial entry

You can add up to 10 entries.

1. Press **MENU / SELECT**.
2. Select **DIRECTORY** → **SPEED DIAL** and scroll to the wanted speed dial location, then press **MENU / SELECT**.
3. Scroll to wanted directory entry, then press **MENU / SELECT** to save.

Speed dial an entry

1. Press and hold the corresponding dialing key (0-9) to display the wanted speed dial entry.
2. Press  / **FLASH** or  to dial.

Caller ID

The telephone supports caller ID service offered by most telephone service providers. Caller ID allows you to see the name, number, date and time of calls. The call log stores the 50 most recent calls. Any modifications made on one handset apply to all.

Contact your service provider for more information about the service.

The handset shows missed calls as **NEW**, after reviewing them the indicator is no longer shown. If you don't want to review the missed calls one by one, press and hold  / **CLEAR**. All entries are then considered as old and kept in the call log.

Review and dial from the call log

1. Press  / 
2. Scroll to browse through the call log. When the dialed number appears:
 - Press  repeatedly to show different dialing options.
 - Press 1 repeatedly if you need to add 1 in front of the telephone number.
3. Press  / **FLASH** or  to dial.

Save number from the call log

1. Scroll to browse through the call log, when the wanted number is displayed press **MENU / SELECT**.
2. When **EDIT NUMBER** is displayed, use the dialing keys (0-9) to edit. Press **MENU / SELECT** to continue.
3. When **EDIT NAME** is displayed, use the dialing keys (0-9) to edit.
4. Press **MENU / SELECT** to save.

Delete number in the call log

To delete one entry, when the wanted entry displays press **MUTE / REMOVE**.
Delete all entries:

1. Press **MENU / SELECT**.
2. Select **CALL LOG → DEL ALL CALLS**.
3. When the screen displays **DELETE ALL CALLS?**, press **MENU / SELECT**.

Intercom

Your handset and base station is pre-registered. Up to five handsets can be registered to the base station. Each new handset is assigned in sequential order. The number is shown on the display. The base unit should be positioned centrally so that the coverage area is roughly equal for all handsets.



If for some reason, the handset is not registered to the base station (**OUT OF RANGE OR NO PWR AT BASE** is displayed even when the handset is near the base station), register your handset according to the following procedure.

Find handset

Press the find handset icon () on the telephone base once. The handset rings and the display shows ****PAGING****. To stop:

- Press the find handset icon on the telephone base again.
- Press  / **FLASH**,  or any dialing key (0-9, ***** or **#**).
- Place the handset in the telephone base or charger.

Add and register additional handset

1. Place the new handset in the telephone base or charger.
2. On the telephone base: press and hold the find handset icon  for 5-8 seconds.
3. On the handset: remove the handset from the telephone base or charger.
Press **MENU / SELECT**  **REGISTRATION** OR press **#** and **REGISTERING... PLEASE WAIT** is displayed.

When the registration is complete, **HANDSET REGISTERED** is displayed.

If the handset registration is successful, you will hear a confirmation tone. The handset will automatically be allocated the next available handset number. If the handset registration is unsuccessful, the **OUT OF RANGE OR NO PWR AT BASE** will again be displayed.

Replace handset

If you want to replace a handset or reassign the designated handset number, you must unregister all handsets and then individually register each handset.

On the telephone base: press and hold the find handset icon  for about 10 seconds, then press  again.

To join a call

1. When a handset is in a call, press  / **FLASH** or  on another handset.
2. Press  / **CLEAR** to exit the call. The call continues until all handsets hangs up.

Intercom (two or more handsets)

1. Press **INT** on your handset when it is not in use. If having more than three handsets, use the dialing keys to enter a destination handset number.
2. To answer the intercom call, press  / **FLASH**, , **INT** or any dialing keys (0-9, ***** or **#**) on the other handset.
3. To end the call, press  / **CLEAR** or place the handset in the telephone base or charger.

Answer incoming call during intercom call

If you receive a call during an intercom call, an alert tone is heard.

- To answer the incoming call: press  / **FLASH** and the intercom call ends automatically.
- To end the intercom call and *not* answering the incoming call: press  / **CLEAR** and the intercom call ends and the telephone continues to ring.

Transferring calls

During an intercom call, you have the following options:

- Add an destination handset to the outside call in a three-way conversation. Press and hold **INT** on the original handset.
- Transfer a call. Press  / **CLEAR** or place the handset in the telephone base/charger. The destination handset automatically connects to the outside call.
- Switch between the outside call (**OUTSIDE CALL**) and the intercom (**INTERCOM**) call by pressing **INT**.
- End the intercom call by pressing  / **CLEAR** or by placing the handset back in the telephone base/charger. The outside call continues with the original handset.

Troubleshooting

Check that the telephone line cord is undamaged and properly plugged in. Disconnect any additional equipment, extension cords and other telephones. If the problem is resolved, the fault is with some other equipment.

Test the equipment on a known working line (e.g. at a neighbor's house). If the equipment works there, the fault is probably with your telephone line. Please report this to your network operator.

No number shown when ringing

- This function requires a Caller ID subscription from your service provider.

Warning tone/Not possible to make calls

- The batteries may be running low (charge the handset).
- The hand unit may be nearly out of range. Move closer to a base unit.

Telephone does not work

- Check that the power adapter/telephone line cord is connected correctly.
- Check that you have selected Tone Dial, see *System Settings*.
- Check the charge status of the handset batteries.
- Connect another functioning telephone to the network wall socket.

The telephone continues to ring

Some of the telephone's ringers do not follow the line signal, which means that the signal may continue for up to 8 seconds after answering a call.

If the telephone still does not work, please contact the place of purchase for service. Don't forget the receipt or a copy of the invoice.

IMPORTANT SAFETY INFORMATION

Your unit is a technically advanced product and should be treated with the greatest care. Negligence may void the warranty.

When using your unit, basic safety precautions should always be followed to reduce the risk of fire, electric shock and injury to persons, including the following:

1. Read and understand all instructions.
2. Follow all warnings and instructions marked on the product.
3. Unplug this product from the wall outlet before cleaning. Do not use liquid or aerosol cleaners. Use a damp cloth for cleaning.
4. Do not use this product near water such as near a bath tub, wash bowl, kitchen sink, laundry tub or swimming pool, or in a wet basement or shower.
5. Do not place this product on an unstable table, shelf, stand or other unstable surfaces.
6. Slots and openings in the back or bottom of the telephone base and handset are provided for ventilation. To protect them from overheating, these openings must not be blocked by placing the product on a soft surface such as a bed, sofa or rug. This product should never be placed near or over a radiator or heat register. This product should not be placed in any area where proper ventilation is not provided.
7. This product should be operated only from the type of power source indicated on the marking label. If you are not sure of the type of power supply in your home or office, consult your dealer or local power company.
8. Do not allow anything to rest on the power cord. Do not install this product where the cord may be walked on.
9. Never push objects of any kind into this product through the slots in the telephone base or handset because they may touch dangerous voltage points or create a short circuit. Never spill liquid of any kind on the product.
10. To reduce the risk of electric shock, do not disassemble this product, but take it to an authorized service facility. Opening or removing parts of the telephone base or handset other than specified access doors may expose you to dangerous voltages or other risks. Incorrect reassembling can cause electric shock when the product is subsequently used.
11. Do not overload wall outlets and extension cords.
12. Unplug this product from the wall outlet and refer servicing to an authorized service facility under the following conditions:
 13. When the power supply cord or plug is damaged or frayed.
 14. If liquid has been spilled onto the product.
 15. If the product has been exposed to rain or water.
 16. If the product does not operate normally by following the operating instructions. Adjust only those controls that are covered by the operation instructions. Improper adjustment of other controls may result in damage and often requires extensive work by an authorized technician to restore the product to normal operation.

Save these instructions

IMPORTANT SAFETY INFORMATION

17. If the product has been dropped and the telephone base and/or handset has been damaged.
18. If the product exhibits a distinct change in performance.
19. Avoid using a telephone (other than cordless) during an electrical storm. There is a remote risk of electric shock from lightning.
20. Do not use the telephone to report a gas leak in the vicinity of the leak. Under certain circumstances, a spark may be created when the power adapter is plugged into the power outlet, or when the handset is replaced in its cradle. This is a common event associated with the closing of any electrical circuit. The user should not plug the telephone into a power outlet, and should not put a charged handset into the cradle, if the telephone is located in an environment containing concentrations of flammable or flame-supporting gases, unless there is adequate ventilation. A spark in such an environment could create a fire or explosion. Such environments might include: medical use of oxygen without adequate ventilation; industrial gases (cleaning solvents; gasoline vapors; etc.); a leak of natural gas; etc.
21. Only put the handset of your telephone next to your ear when it is in normal talk mode.
22. The power adapters are intended to be correctly oriented in a vertical or floor mount position. The prongs are not designed to hold the plug in place if it is plugged into a ceiling, under-the-table or cabinet outlet.

The advice applies to the unit, battery, power adapter and other accessories. If the telephone is not working as it should, please contact the place of purchase for service. Don't forget the receipt or a copy of the invoice.

Precautions for users of implanted cardiac pacemakers

Cardiac pacemakers (applies only to digital cordless telephones):

Wireless Technology Research, LLC (WTR), an independent research entity, led a multidisciplinary evaluation of the interference between portable wireless telephones and implanted cardiac pacemakers.

Supported by the U.S. Food and Drug Administration, WTR recommends to physicians that:

Pacemaker patients

- Should keep wireless telephones at least six inches from the pacemaker.
- Should NOT place wireless telephones directly over the pacemaker, such as in a breast pocket, when it is turned ON.
- Should use the wireless telephone at the ear opposite the pacemaker. WTR's evaluation did not identify any risk to bystanders with pacemakers from other persons using wireless telephones.

Save these instructions

Equipment Approval Information

FCC, ACTA and IC regulations

This equipment has been tested and found to comply with the requirements for a Class B digital device under Part 15 of the Federal Communications Commission (FCC) rules. These requirements are intended to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This device complies with Part 15 of the FCC rules. Operation is subject to the following two conditions:

(1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Privacy of communications may not be ensured when using this telephone.

To ensure safety of users, the FCC has established criteria for the amount of radio frequency energy that can be safely absorbed by a user or bystander according to the intended usage of the product. This product has been tested and found to comply with the FCC criteria. The handset may be safely held against the ear of the user. The telephone base shall be installed and used such that parts of the user's body other than the hands are maintained at a distance of approximately 8 inches (20 cm) or more.

This Class B digital apparatus complies with Canadian ICES-003.

FCC Part 68 and ACTA

This equipment complies with Part 68 of the FCC rules and with technical requirements adopted by the Administrative Council for Terminal Attachments (ACTA). The label on the back or bottom of this equipment contains, among other things, a product identifier in the format US:AAAEQ##TXXXX. This identifier must be provided to your local telephone company upon request.

The plug and jack used to connect this equipment to premises wiring and the telephone network must comply with applicable Part 68 rules and technical requirements adopted by ACTA. A compliant telephone cord and modular plug is provided with this product. It is designed to be connected to a compatible modular jack that is also compliant. An RJ11 jack should normally be used for connecting to a single line and an RJ14 jack for two lines. See Installation Instructions in the user's manual.

The Ringer Equivalence Number (REN) is used to determine how many devices you may connect to your telephone line and still have them ring when you are called. The REN for this product is encoded as the 6th and 7th characters following the US: in the product identifier (e.g., if ## is 03, the REN is 0.3).

In most, but not all areas, the sum of all RENs should be five (5.0) or less. For more information, please contact your local telephone company.

This equipment may not be used with Party Lines. If you have specially wired alarm dialing equipment connected to your telephone line, ensure the connection of this equipment does not disable your alarm equipment. If you have questions about what will disable alarm equipment, consult your telephone company or a qualified installer.

If this equipment is malfunctioning, it must be unplugged from the modular jack until the problem has been corrected. Repairs to this telephone equipment can only be made by the manufacturer or its authorized agents. For repair procedures, follow the instructions outlined under the Warranty.

If this equipment is causing harm to the telephone network, the telephone company may temporarily discontinue your telephone service. The telephone company is required to notify you before interrupting service. If advance notice is not practical, you will be notified as soon as possible. You will be given the opportunity to correct the problem and the telephone company is required to inform you of your right to file a complaint with the FCC. Your telephone company may make changes in its facilities, equipment, operation, or procedures that could affect the proper functioning of this product. The telephone company is required to notify you if such changes are planned.

If this product is equipped with a corded or cordless handset, it is hearing aid compatible.

If this product has memory dialing locations, you may choose to store emergency telephone numbers (e.g., police, fire, medical) in these locations. If you do store or test emergency numbers, please:

- Remain on the line and briefly explain the reason for the call before hanging up.
- Perform such activities in off-peak hours, such as early morning or late evening.

Industry Canada

This Class B digital apparatus complies with Canadian ICES-003.

Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

The abbreviation, "IC:", before the registration number signifies that registration was performed based on a Declaration of Conformity indicating that Industry Canada technical specifications were met.

The Ringer Equivalence Number (REN) assigned to each terminal device provides an indication of the maximum number of terminals allowed to be connected to a telephone interface. The termination on an interface may consist of any combination of devices subject only to the requirement that the sum of the Ringer Equivalence Numbers of all the devices does not exceed 5.

The Ringer Equivalence Number (REN) for this terminal equipment is 0.1.

This equipment meets the applicable Industry Canada Terminal Equipment Technical Specifications.

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FCC ID: EW780-8938-00

IC: 1135B-80893800

Disclaimer and limitation of liability

Doro AB and its affiliates or licensees assume no responsibility for any damage or loss resulting from the use of this guide. Doro AB and its affiliates or licensees assume no responsibility for any loss or claims by third parties that may arise through the use of this software. Doro AB and its affiliates or licensees assume no responsibility for any damage or loss caused by deletion of data as a result of malfunction, dead battery, or repairs. Be sure to make backup copies of all important data on other media to protect against data loss.

TO THE MAXIMUM EXTENT PERMITTED BY APPLICABLE LAW, UNDER NO CIRCUMSTANCES SHALL DORO AB OR ANY OF ITS AFFILIATES OR LICENSEES BE RESPONSIBLE FOR ANY LOSS OF DATA OR INCOME OR ANY SPECIAL, INCIDENTAL, CONSEQUENTIAL OR INDIRECT DAMAGES HOWSOEVER CAUSED. THE CONTENTS OF THIS DOCUMENT ARE PROVIDED "AS IS." EXCEPT AS REQUIRED BY APPLICABLE LAW, NO WARRANTIES OF ANY KIND, EITHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE MADE IN RELATION TO THE ACCURACY, RELIABILITY, OR CONTENTS OF THIS DOCUMENT. DORO AB RESERVES THE RIGHT TO REVISE THIS DOCUMENT OR WITHDRAW IT AT ANY TIME WITHOUT PRIOR NOTICE.

Warranty and other information

Guarantee

This product is guaranteed for a period of 12 months from the date of purchase. In the unlikely event of a fault occurring during this period, please contact the place of purchase. Proof of purchase is required for any service or support needed during the guarantee period. This guarantee will not apply to a fault caused by an accident or a similar incident or damage, liquid ingress, negligence, abnormal or improper usage, non-maintenance or any other circumstances on the user's part. Furthermore, this guarantee will not apply to any fault caused by a thunderstorm or any other voltage fluctuations. As a matter of precaution, we recommend disconnecting the charger during a thunderstorm.

Except as required by applicable law, batteries are consumables and are not included in any guarantee.

This guarantee does not apply if batteries other than DORO original batteries are used.

EXCEPT FOR THE WARRANTY EXPRESSLY SET FORTH HEREIN, WIDEX USA INC. /WIDEX CANADA LTD. DISCLAIMS AND EXCLUDES ALL REPRESENTATIONS, WARRANTIES AND CONDITIONS, WHETHER WRITTEN OR ORAL, EXPRESS OR IMPLIED, STATUTORY OR OTHERWISE WITH RESPECT TO THE UNIT AND ELEMENTS OR ACCESSORIES THEREOF, INCLUDING, BUT NOT LIMITED TO, WITHOUT LIMITATION, IMPLIED WARRANTIES AND CONDITIONS OF MERCHANTABILITY, DURABILITY, QUALITY OR FITNESS FOR PARTICULAR PURPOSE, INCLUDING ANY AND ALL WARRANTIES AND CONDITIONS FOUND IN APPLICABLE LEGISLATION.

About cordless telephones


- **Privacy:** The same features that make a cordless telephone convenient create some limitations. Telephone calls are transmitted between the telephone base and the cordless handset by radio waves, so there is a possibility that the cordless telephone conversations could be intercepted by radio receiving equipment within range of the cordless handset. For this reason, you should not think of cordless telephone conversations as being as private as those on corded telephones.
- **Electrical power:** The telephone base of this cordless telephone must be connected to a working electrical outlet. The electrical outlet should not be controlled by a wall switch. Calls cannot be made from the cordless handset if the telephone base is unplugged, switched off or if the electrical power is interrupted.
- **Potential TV interference:** Some cordless telephones operate at frequencies that may cause interference to televisions and VCRs. To minimize or prevent such interference, do not place the telephone base of the cordless telephone near or on top of a TV or VCR. If interference is experienced, moving the cordless telephone farther away from the TV or VCR often reduces or eliminates the interference.
- **Rechargeable batteries:** Exercise care in handling batteries in order not to create a short circuit with conducting material such as rings, bracelets, and keys. The battery or conductor may overheat and cause harm. Observe proper polarity between the battery and the battery charger.
- **Nickel-metal hydride rechargeable batteries:** Dispose of these batteries in a safe manner. Do not burn or puncture the battery. Like other batteries of this type, if burned or punctured, they could release caustic material which could cause injury. When the battery in this unit has reached the end of its useful life, please dispose of it according to your local environmental laws and regulations. For information about recycling cell phone batteries, go to: <http://www.rbc.org>

Technical specifications

Transmit frequency	Handset: 1921536-1928448 MHz Telephone base: 1921536-1928448 MHz Widex-Link 106 MHz
Channels	5
Nominal effective range	Maximum power allowed by FCC and IC. Actual operating range may vary according to environment conditions at the time of use.
Power requirements	Handset: 2.4V 500mAh AAA x 2 Ni-MH battery Telephone base voltage: Adaptor Input: 100-120V, AC voltage, 60Hz Adaptor Output: 5V DC @ 300mA

Support

For assistance, please contact an Authorized Widex Hearing Care Professional:

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© 2013 Doro AB. All rights reserved.
Widex USA Inc.