**VIEWWORKS**

# VIVIX-S 4386W Service Manual

## (FXRD-4386WB)

© Vieworks. 2021 All rights reserved.

The specifications and related information in this manual may be changed without notice. Refer to Vieworks Download System (VDS) for the latest version of our manuals.

Under Copyright laws, this manual should not be reproduced, in whole or in part, without the written permission of Vieworks. This document is basically produced in English, and can be produced in the language of the relevant country upon request of the customer.

## Contents

<b>1. Instruction .....</b>	<b>6</b>
<b>1.1 Document Guide .....</b>	<b>7</b>
1.1.1 Precautions .....	7
1.1.2 Target .....	7
1.1.3 Symbols .....	7
1.1.4 Notations .....	7
1.1.5 Contact Us .....	8
<b>1.2 Purpose of Use .....</b>	<b>9</b>
1.2.1 Intended Us .....	9
1.2.2 Patient group .....	9
1.2.3 Contra-Indication .....	9
<b>1.3 Clinical Benefits of Product Use .....</b>	<b>10</b>
<b>1.4 Features .....</b>	<b>11</b>
<b>1.5 Product Usage Guide .....</b>	<b>12</b>
1.5.1 Product Usage .....	12
1.5.2 Disclaimer .....	12
1.5.3 Restrictions on the outdoor use (applied to some nations) .....	13
1.5.4 Product Disposal .....	13
<b>2. Product .....</b>	<b>14</b>
<b>2.1 Product Components .....</b>	<b>15</b>
2.1.1 Detector .....	15
2.1.2 Battery Pack .....	16
2.1.3 Battery Charger .....	16
2.1.4 SCU (System Control Unit) – Optional .....	17
2.1.5 Wireless Power Transmitter – Optional .....	18
2.1.6 Accessories .....	19
<b>2.2 FXRD-4386WB .....</b>	<b>22</b>
2.2.1 Specifications .....	22
2.2.2 Functions .....	23
2.2.3 Deco Sheet .....	25
2.2.4 Wireless Communication .....	26
2.2.5 Use Environment .....	26
<b>2.3 Battery Pack (FXRB-04A) .....</b>	<b>27</b>
2.3.1 Specifications .....	27
<b>2.4 Battery Charger (FXRC-04A) .....</b>	<b>28</b>

2.4.2	Specifications .....	28
2.4.3	Drawings .....	28
2.4.4	Functions .....	29
<b>2.5</b>	<b>SCU mini (FXRS-04A) – Optional .....</b>	<b>31</b>
2.5.1	Specifications .....	31
2.5.2	Drawings .....	31
2.5.3	Functions .....	32
<b>2.6</b>	<b>Wireless Power Transmitter (FXRC-05A) – Optional .....</b>	<b>34</b>
2.6.2	Specifications .....	34
2.6.3	Drawings .....	34
2.6.4	Functions .....	35
<b>2.7</b>	<b>Others .....</b>	<b>36</b>
2.7.1	X-ray Generator (Recommended shooting condition).....	36
2.7.2	Recommended Specifications of Workstation (PC) .....	36
<b>3.</b>	<b>System Configuration.....</b>	<b>37</b>
<b>3.1</b>	<b>Detector Connection Method .....</b>	<b>38</b>
3.1.1	Wireless Connection.....	38
3.1.2	Wired Connection.....	39
<b>3.2</b>	<b>System Configuration.....</b>	<b>41</b>
3.2.1	SCU AP Mode.....	41
3.2.2	External AP Mode .....	43
3.2.3	Detector AP Mode .....	44
3.2.4	Tether Interface Mode.....	46
3.2.5	Detector Stand-Alone Mode .....	47
<b>3.3</b>	<b>Changing the Wireless Setting .....</b>	<b>48</b>
3.3.1	Detector AP / Station.....	48
3.3.2	Preset Switching .....	50
3.3.3	Sync Wireless Settings.....	52
<b>3.4</b>	<b>Generator Interface.....</b>	<b>53</b>
3.4.1	AED (Auto Exposure Detection) Interface .....	53
3.4.2	DR Trigger Interface .....	54
3.4.3	SW Trigger Interface.....	55
<b>3.5</b>	<b>Configuring DR Trigger Interface .....</b>	<b>56</b>
3.5.1	Trigger Interface Way.....	56

3.5.2	Packet Trigger .....	56
3.5.3	Line Trigger .....	57
3.5.4	EXT_INF Port Pin Map .....	58
3.5.5	Input / Output Circuit.....	59
<b>3.6</b>	<b>Exposure Mode.....</b>	<b>60</b>
3.6.1	AED (Auto Exposure Detection) Interface .....	60
3.6.2	DR Trigger Interface .....	61
3.6.3	SW Trigger Interface.....	62
<b>4.</b>	<b>Settings .....</b>	<b>63</b>
<b>4.1</b>	<b>Product Installation .....</b>	<b>64</b>
4.1.1	Connecting SCU mini (FXRS-04A).....	67
4.1.2	Booting up the Detector .....	67
4.1.3	Installing Wireless Power Transmitter (FXRC-05A) .....	68
<b>4.2</b>	<b>Device Setting .....</b>	<b>71</b>
4.2.1	Software Installation.....	71
4.2.2	Setting Detector and PSU / SCU .....	71
<b>4.3</b>	<b>Diagnosis of Devices .....</b>	<b>72</b>
4.3.1	Image Diagnosis.....	72
4.3.2	Diagosing Battery Pack .....	73
4.3.3	Diagnosing Communication Status .....	75
4.3.4	Diagnosing Communication Speed .....	79
4.3.5	Self-Diagnosis .....	80
<b>5.</b>	<b>Maintenance.....</b>	<b>83</b>
<b>5.1</b>	<b>Product Initialization .....</b>	<b>84</b>
5.1.1	SCU Initialization.....	84
5.1.2	Detector Initialization.....	85
5.1.3	Wireless Initialization of Detector .....	86
<b>5.2</b>	<b>Detector Power Save Function .....</b>	<b>87</b>
<b>6.</b>	<b>Troubleshooting .....</b>	<b>89</b>
<b>6.1</b>	<b>Troubleshooting.....</b>	<b>90</b>
6.1.1	Troubleshooting Guide.....	90
6.1.2	Fail to Turn the Detector On .....	90
6.1.3	The Power Switch of SCU or Status LED is not worked normally .....	90
6.1.4	Communication Test is failed.....	91
6.1.5	Communication is connected, but no image is acquired .....	91
6.1.6	Getting Abnormal Images .....	92

6.1.7	When Hours of Battery Use Decreases .....	92
6.1.8	When an Error Number is Displayed on the OLED Display .....	92
<b>7.</b>	<b>Information .....</b>	<b>93</b>
<b>7.1</b>	<b>Service Information .....</b>	<b>94</b>
7.1.1	Product Lifetime .....	94
7.1.2	Regular Inspection and Maintenance .....	94
7.1.3	Repair .....	94
7.1.4	Replacement Support .....	94
7.1.5	Consumables.....	94
<b>7.2</b>	<b>Warranty .....</b>	<b>95</b>
<b>7.3</b>	<b>Revision History .....</b>	<b>96</b>

---

# 1. Instruction

---

This section gives basic information of this manual and products.

Document Guide

Product of Use

Clinical Benefits of Product Use

Features

Product Usage Guide

## 1.1 Document Guide

This Operation Manual explains how to use the **VIVIX-S 4386W** detector made by Vieworks, X-ray interface unit and other peripheral equipment. Use this manual to set up the **VIVIX-S 4386W** detector as well as understand its various functions.

### 1.1.1 Precautions

Read this manual thoroughly to avoid any accident occurred by negligence.

### 1.1.2 Target

This manual is intended for the users who set up and operate the **VIVIX-S 4386W** detector.

### 1.1.3 Symbols

This product should be operated under the safety instructions with the warning or caution symbol in this manual. It is important for you to read and understand the contents to operate the products safely.

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### 1.1.4 Notations

#### Bold Types

Words in bold indicate products terms, or the sentences which are needed to transmit clear meaning to the customers.

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### 1.1.5 Contact Us

- This manual is provided in print format upon request by the customer.
- For comments or inquiries regarding this document and relevant products, contact via email below:

Item	Contents
Department	Customer Support Team at Vieworks
E-mail	<a href="mailto:CustomerSupport@vieworks.com">CustomerSupport@vieworks.com</a>

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## 1.2 Purpose of Use

### 1.2.1 Intended Us

**VIVIX-S 4386W** is a digital flat panel detector that is used for screening and diagnosis of disease of injury. This detector is intended for use by a qualified/trained doctor or technician on both adult subjects for taking diagnostic radiographic exposures of the skull, spinal column, chest, abdomen, extremities, and other body parts.

### 1.2.2 Patient group

VIVIX-S 4386W is used to examine in patients with or suspected of muscle and bone injury, respiratory diseases. It is intended for general patients such as adults, children, and infants, but all radiography should be reviewed by a doctor in charge prior to beginning the examination.

The radiologist should use a proper technique considering the patient's size to decrease the radiation dose when he or she acquires diagnostic images.

### 1.2.3 Contra-Indication

These detectors are not intended to be used for mammography applications.

### 1.3 Clinical Benefits of Product Use

Compared to X-ray exposure to a patient during X-ray imaging using this detector, the benefit of diagnosing a patient's disease or injury is much greater. The digital flat panel detector has superior X-ray image quality compared to the conventional analog film type and enables real-time image transmission and image processing, so that the benefit is much greater in terms of ease of diagnosis and usability compared to electrical risks that may occur.

## 1.4 Features

- The new sensor with 140  $\mu\text{m}$  of pixel pitch produces high-resolution (approx. 18.87 Mega pixels) digital images.
- The Gadox is provided as a scintillator.
  - FXRD-4386WB –2 Gadoxes (Gadolinium Oxysulfide, 430.08 x 430.08mm) applied.
- This product supports **IEEE 802.11n/ac** and has a wireless communications function, Inside AP™, that shooting is possible everywhere.
- The detector connection status and remaining battery level are visible on the OLED Display.
- Makes direct wireless communication with the built-in wireless AP function. (Inside AP™)
- Quick application of various functions with the two buttons on the side of the detector is available.
- The product supports the stable and reliable Auto Exposure Detection (AED) function. (Anytime™)
- The product can be wired with the tether interface cable according to the operational environment.

## 1.5 Product Usage Guide

This chapter provides instructions about the use of the product, disposal, and the liability limit of Vieworks.

### 1.5.1 Product Usage

- 1 Only a physician or a legally certified operator should use the products.
- 2 The equipment should be maintained in a safe and operable condition by maintenance personnel.
- 3 Follow the guidelines in this manual when installing and using this product.
- 4 Use only computers and image display monitors recommended by this manual.
- 5 Use only the dedicated cables provided with the products.
- 6 For details about installation and use of the products, consult your sales representative or a distributor.

### 1.5.2 Disclaimer

- 1 In no event shall Vieworks be liable for damage or loss arising from fire, earthquake, any action or accident by a third party, any intentional or negligent action by users.
- 2 In no event shall Vieworks be liable for damage or loss arising from any trial usage, or other usage under abnormal conditions.
- 3 In no event shall Vieworks be liable for personal physical harm or property damage that is sustained when the instructions are not followed.
- 4 In no event shall Vieworks be liable for any damage arising from moving, alteration, inspection, or repair by a person other than authorized service engineers.
- 5 Vieworks shall not be liable for loss of image data for any reason while using this product.
- 6 Roentgenography, image processing, image reading, and image data storage must be performed in accordance with the laws of the country or religion in which the product is being used.
- 7 The user is responsible for maintaining the privacy of image data made from this product.
- 8 It is the responsibility of the attending physicians to provide medical care services. Vieworks will not be liable for faulty diagnoses.
- 9 Specifications, composition, and appearance of this product may change without prior notice.

### 1.5.3 Restrictions on the outdoor use (applied to some nations)

There are restrictions on the outdoor use of the U-NII Low (5150-5250 MHz) and U-NII Mid (5250-5350 MHz) bandwidths of the WLAN module incorporated in the device in the following Member States: Belgium (BE), Bulgaria (BG), Czech Republic (CZ), Denmark (DK), Germany (DE), Estonia (EE), Ireland (IE), Greece (EL), Spain (ES), France (FR), Croatia (HR), Italy (IT), Cyprus (CY), Latvia (LV), Lithuania (LT), Luxembourg (LU), Hungary (HU), Malta (MT), Netherlands (NL), Austria (AT), Poland (PL), Portugal (PT), Romania (RO), Slovenia (SI), Slovakia (SK), Finland (FI), Sweden (SE) and United Kingdom (UK).

	AT	BE	BG	CY	CZ
	DE	DK	EE	EL	ES
	FI	FR	HR	HU	IE
	IT	LT	LU	LV	MT
	NL	PL	PT	RO	SE
	SI	SK	UK		

### 1.5.4 Product Disposal

Disposal of this product in an unlawful manner may have a negative impact on health and on the environment. When disposing of this product, therefore, be absolutely sure to follow the procedure which is in conformity with the laws and regulations applicable in your area.

#### European Union (and EEA\*) only

This symbol indicates that this product is not to be disposed with your household waste, according to the WEEE Directive (2012/19/EU) and your national law.

This product should be handed over to a designated collection point, e.g., on an authorized one-for-one basis when you buy a new similar product or to an authorized collection site for recycling electrical and electronic equipment (EEE). Improper handling of this type of waste could have a negative impact on the environment and human health due to potentially hazardous substances that are generally associated with EEE. At the same time, your cooperation in the correct disposal of this product will contribute to the effective usage of natural resources. For more information on where you can drop off your waste equipment for recycling, please contact your local city office, waste authority, approved WEEE scheme, or your household waste disposal service.

\*EEA: Norway, Iceland, and Liechtenstein

---

## 2. Product

---

This section gives an instruction about the product components and their specifications.

Product Specifications

FXRD-4386WB

Battery Pack (FXRB-04A)

Battery Charger (FXRC-04A)

SCU mini (FXRS-04A) – Optional

Wireless Power Transmitter (FXRC-05A) – Optional


Others

## 2.1 Product Components

### 2.1.1 Detector

Product Name	Model / Description
--------------	---------------------

**VIVIX-S 4386W**  
**Wireless Detector**

- FXRD-4386WB (8kg)
  - 2 Battery Packs included

 There are two tethered interfaces on the side and the rear of the detector, and connection can be made to only one of them. Available interface for connection may vary by clients.

The design of the deco sheet attached to the detector can be different for each customer.

The battery pack is removable from the detector. The battery pack can be replaced at the end of its life.

The two battery packs are basically equipped with FXRD-4386WB.

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 Do not open the tether interface covered with the instrument. If not, a breakdown or electrical shock may occur.

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### 2.1.2 Battery Pack

Product Name	Image	Model / Description
Battery Pack		<ul style="list-style-type: none"> <li>• FXRB-04A                             <ul style="list-style-type: none"> <li>▫ 3,400mAh</li> <li>▫ 185g</li> </ul> </li> </ul>

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 Wireless power transmitter is not a standard item provided to customers. If necessary, you can purchase it as an option.

• Wireless power transmitter should be wired. The available power voltage is DC 17V±10% (DC 15~19V), up to 24W.

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 Make sure to use the wireless power transmitter only specified for **VIVIX-S 4386W**.

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**2.1.6 Accessories**

Item	Image	Description
Hardware manual CD	 <p>A CD-ROM with the VIVIX logo and text: "Version 2.6 20210510", "VIVIX-S Manual", "Operation Manual", "Service Manual", and "https://clouds.viewworks.com/5001/". It also includes manufacturer information for Viewworks Co., Ltd.</p>	<ul style="list-style-type: none"> <li>• Detector Operation Manual</li> <li>• Detector Service Manual</li> </ul>
Software & Software manual CD (VXvue)	 <p>A CD-ROM with the VIVIX logo and text: "VXvue V1.0.2.6", "Document", "Operation Manual", "Service Manual", "DICOM CS", "Quick-Start Guide", "XIPL User Manual", and "https://clouds.viewworks.com/5001/". It also includes manufacturer information for Viewworks Co., Ltd.</p>	<ul style="list-style-type: none"> <li>• Software file</li> <li>• Software manual</li> </ul>
Cables (Optional)	 <p>A black tether interface cable with a magnetic lock connector on one end and a standard connector on the other.</p>	<ul style="list-style-type: none"> <li>• VW tether interface cable (7m, magnetic lock)                             <ul style="list-style-type: none"> <li>▪ Connect to the tether interface on the side</li> </ul> </li> </ul>
	 <p>A black tether interface cable with a Fix/release holder connector on one end and a standard connector on the other.</p>	<ul style="list-style-type: none"> <li>• N tether interface cable (7m, Fix/release holder)                             <ul style="list-style-type: none"> <li>▪ Connect to the tether interface on the rear</li> </ul> </li> </ul>
	 <p>A black generator interface cable with a circular connector on one end and a standard connector on the other.</p>	<ul style="list-style-type: none"> <li>• Generator interface cable (15m)</li> </ul>
	 <p>A white UTP LAN cable with RJ45 connectors on both ends.</p>	<ul style="list-style-type: none"> <li>• UTP LAN cable (15m, Direct) CAT5E or more than CAT6</li> </ul>

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- Use only cables designated for the **VIVIX-S 4386W** detectors.
- The use of accessories and cables other than those approved and sold by Viewworks Co., Ltd. may result in increased release of electromagnetic waves or decreased stability of the equipment.
- Accessory equipment connected to the analog and digital interfaces must be certified according to the respective IEC standards. All combinations of equipment must follow **IEC 60601-1** system requirements.
- Any person who connects additional equipment to the signal input or signal output ports configures a medical system and is therefore responsible for ensuring that the system complies with the requirements of the system standard **IEC 60601-1**.
- Consult your sales distributor or manufacturer if you have any concerns.

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## 2.2 FXRD-4386WB

### 2.2.1 Specifications


Item	Specifications
<b>Model</b>	• FXRD-4386WB
<b>Image Scan Mode</b>	• Full mode • Partial mode
<b>Image Sensor</b>	• a-Si(Amorphous Silicon) TFT with PIN diode
<b>X-ray Scintillator Type</b>	• Gd <sub>2</sub> O <sub>2</sub> S:Tb (Gadolinium oxysulfide)
<b>Pixel Pitch</b>	• 0.14mm (140μm)
<b>Field of View</b>	• 43cm x 86cm • 43cm x 43cm
<b>Active Area (H x V)</b>	• 430.08mm x 860.16mm • 430.08mm x 430.08mm
<b>Active Array</b>	• 3072 x 6144 pixels • 3072 x 3072 pixels
<b>Effective Area</b>	• 428.4mm x 858.2mm • 428.4mm x 425.6mm
<b>Effective Array</b>	• 3060 x 6130 pixels • 3060 x 3040 pixels
<b>Grayscale</b>	• 16 bit
<b>Spatial Resolution</b>	• Min. 3.5 lp/mm
<b>Image Acquisition Time (Wired)</b>	• Up to 10 sec • Up 5 sec (Exposure time is set to 500ms / Exposure time is excluded)
<b>Image Acquisition Time (Wireless)</b>	• Up to 10 sec • Up to 5 sec (IEEE802.11ac, MiMO 3x3, 5GHz, 80MHz) (Exposure time is set to 500ms / Exposure time is excluded)
<b>Recommended Cycle Time</b>	• At least 12 sec (Optimized wired/wireless environment, except for when Exposure time is set to 500ms and when the software is processing the image)
<b>X-ray Synchronization Control</b>	• AED (Auto Exposure Detection) • DR Trigger (External line trigger) • Software Trigger
<b>Power Supply</b>	• Supplies power via tether interface cable from SCU: DC 24V, Max. 1.0A • Supplies power with the 2 battery packs: DC 10~13.2V, Max. 78.54Wh • Supplies power with the wireless power transmitter: DC 12V, Max. 1.25A
<b>Power Consumption</b>	• Max. 24 W
<b>Operating Time (beginning of a battery life)</b>	• 10hrs (acquires image in every 100 sec.)
<b>Dimensions (H x W x D)</b>	• 900mm x 465mm x 17.0mm
<b>Weight</b>	• 8kg (including 2 battery packs)
<b>Image Transmission</b>	• Wired: Gigabit Ethernet (1000BASE-T) via PoE (Power over Ethernet) • Wireless: IEEE 802.11n/ac (2.4GHz/5GHz), with 3 antennas
<b>Data Transmission Rate (Wired)</b>	• Max. 1Gbps
<b>Data Transmission Rate (Wireless)</b>	• Max. 300Mbps (IEEE802.11n, MIMO 2x2, 5GHz, 40MHz) • Max. 1300Mbps (IEEE802.11ac, MIMO 3x3, 5GHz, 80MHz)

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
## 2.2.2 Functions

### Side

Name	Description
<b>A Tether Interface Connector</b>	Tether interface cable connector. <ul style="list-style-type: none"> <li>• Used for wired connection between a detector and SCU.</li> </ul>
<b>B Magnet for fixing the tether interface</b>	Used for fixing a tether interface cable
<b>C OLED Display</b>	<ul style="list-style-type: none"> <li>• Displays battery status</li> <li>• Displays wired / wireless connection status</li> <li>• Displays sleep mode status</li> </ul>
<b>D AP Button</b>	<ul style="list-style-type: none"> <li>• Changes AP settings button in wireless communication (Change detector AP / STATION or change preset in STATION mode)</li> <li>• Changes OLED screen</li> </ul>
<b>E Power Indicator LED</b>	<ul style="list-style-type: none"> <li>• Displays system power status</li> <li>• Displays system boot status</li> </ul>
<b>F Power Button</b>	<ul style="list-style-type: none"> <li>• System power on/off</li> <li>• Changes OLED screen</li> </ul>
<b>G Antenna for Wireless LAN</b>	Antennas for wireless communication (3ea)
<b>H Charge Status LED</b>	Displays the charge status of the battery
<b>I Xray Transmitting Part</b>	A part in where X-ray is transmitted


## Rear

Name	Description
<b>A Tether Interface Connector</b>	Tightens the tether interface cable. <ul style="list-style-type: none"> <li>Used for wired connection between a detector and SCU.</li> </ul>
<b>B Tether Interface Holder</b>	Used for fixing and releasing a tether interface cable
<b>C Handle</b>	Handles to install/uninstall the detector to the equipment
<b>D Lift Structure</b>	Used when the detector is placed on flat surface
<b>E Wireless Power Receiver</b>	Receives wireless power
<b>F Infrared LED</b>	An infrared transmission that controls the wireless power of the detector
<b>G Battery Pack Cover</b>	The cover needs to be opened and closed when replacing the battery pack.
<b>H Label</b>	An area in where the label is attached.
<b>I Xray Transmitting Part</b>	A part in where X-ray is transmitted.

### 2.2.3 Deco Sheet

Indication Info.	Description
<b>A</b> Center of the detector	Indicates the central position of detector.
<b>B</b> Image starting point	Indicates the starting point of an original image.
<b>C~E</b> Product safety and performance logos	Notices regarding product safety and performance

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- The RF distance between body and product(FXRC-05A) is 20cm.

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## 2.3 Battery Pack (FXRB-04A)

### 2.3.1 Specifications

Item	Specifications
<b>Model</b>	• FXRB-04A
<b>Type</b>	• Lithium Ion Polymer
<b>Nominal Voltage</b>	• DC +11.55V
<b>Nominal Capacity</b>	• 3,400mAh
<b>Number of Cells</b>	• 3S1P (3 Series 1 Parallel)
<b>Battery Life</b>	• Approx. 800 times (One time standard: Fully charged and then fully discharged)
<b>Dimensions (H × W × D)</b>	• 189.0mm × 89.0mm × 6.65mm (max.)
<b>Weight</b>	• 185g (max.)

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## 2.4 Battery Charger (FXRC-04A)

The battery charger quick charges the battery pack provided with the detector.

This device receives power through the provided AC-DC adaptor or the cradle (FXRR-01A).

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


2.4.4 Functions

Front / Side

No.	Item	Description
1	<b>Battery Tray</b>	Battery cover
2	<b>Battery Interface</b>	Battery link connector
3	<b>Battery Fixer</b>	Anchoring device that fixes battery to the tray
4	<b>Battery Tray State Indicator</b>	Displays the charging state of the battery equipped in each battery tray <ul style="list-style-type: none"> <li>• Orange: Charging</li> <li>• Green: Charged</li> </ul>
5	<b>Tether Interface Connector</b>	Supplies power to the battery charger connected to the tether interface cable or equipped with the cradle DC +24
6	<b>DC Power Input</b>	<ul style="list-style-type: none"> <li>• Supplies power to the battery charger (Only when in sole occupancy)</li> </ul>

## Rear

No.	Item	Description
1	Label	Label attached part


## 2.5 SCU mini (FXRS-04A) – Optional

SCU mini synchronizes the image and X-ray signal as locating among the detector, workstation, and the X-ray generator. The SCU Mini has a wired or wireless connection with the detector (Wired connection to the workstation). You can also use the DR Trigger with the X-ray generator.

### 2.5.1 Specifications


Item	Specification
<b>Model</b>	• FXRS-04A
<b>Power Supply</b>	• Input: DC +24V, Max. 2A
<b>Cable Connection Port</b>	• Gigabit Ethernet port (3ea) • PoE (Power over Ethernet) port (1ea) • X-Ray Generator Interface port (1ea)
<b>Wireless Communications</b>	• IEEE 802.11n (2.4 GHz / 5 GHz)
<b>Dimension (H × W × D)</b>	• 210.0mm × 170.0mm × 45.0mm
<b>Antenna</b>	• 140mm (2ea, dual band)
<b>Weight</b>	• 1.2kg

### 2.5.2 Drawings

### 2.5.3 Functions


Front

No.	Name	Description
1	<b>Power switch</b>	Turns on/off the power of SCU mini.
2	<b>Antenna</b>	Assists communications between the detector and SCU mini.
3	<b>Status LED</b>	Indicates status of SCU mini operation and connection. <ul style="list-style-type: none"> <li>• Blinking green: Booting</li> <li>• Green: Completed to boot up</li> <li>• Blue: The detector is connected to the software and ready to communicate.</li> </ul>

Rear

No.	Name	Description
1	EXT_INF	X-ray generator interface connector (D-SUB 15pin, Female)
2	LAN port (Port 1, 2, 3)	Gigabit Ethernet port (1000BASE-T) <ul style="list-style-type: none"> <li>• <b>Port 1:</b> Communication between the workstation and SCU mini.</li> <li>• <b>Port 2, 3:</b> Communication between multiple detectors.</li> </ul>
3	PoE status lamp	Indicates the status of <b>PoE</b> port. <ul style="list-style-type: none"> <li>• Green: 1 Gbps</li> <li>• Orange: 100 Mbps</li> </ul>
4	PoE port	<b>PoE</b> (Power over Ethernet) port (1000BASE-T) <ul style="list-style-type: none"> <li>• Communication between the detector and SCU mini.</li> <li>• Supplies power to the detector.</li> </ul>
5	DC power input port	DC +24V <ul style="list-style-type: none"> <li>• Supplies power to SCU mini.</li> </ul>

## 2.6 Wireless Power Transmitter (FXRC-05A) – Optional

This equipment wirelessly supplies up to 15W of power to the detector. This equipment should be connected to electrical power. The available voltage of power is DC 17V±10% (DC 15~19V), Max. 24W.

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

## 2.6.4 Functions

### Front / Side / Rear

No.	Item	Description
1	<b>Wireless Power Transmitting Antenna</b>	Wireless power transmitter
2	<b>State LED</b>	Displays active state <ul style="list-style-type: none"> <li>• OFF : Wireless power transmitter is off</li> <li>• Green : Wireless power is being transmitted.</li> <li>• Orange : Error state</li> </ul>
3	<b>DC Power Input Port</b>	DC +17V Supplies power to the wireless power transmitter
4	<b>Setscrew Hole</b>	Hole that fixes wireless power transmitter (4ea)
5	<b>Label</b>	Label attached part
6	<b>Infrared Optic Acquisition</b>	Infrared receiver to control wireless power from detector.

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## 2.7 Others

### 2.7.1 X-ray Generator (Recommended shooting condition)

Item	Recommended condition
X-ray Energy Range	40kVp ~ 150kVp
Lifetime Dose	100Gy

### 2.7.2 Recommended Specifications of Workstation (PC)

Item	Minimum Specifications	Recommended Specifications
OS	Windows 7 Professional SP1 (64bit) Windows 8.1 Professional (64bit)	Windows 10 Professional (64bit)
LAN Card	<ul style="list-style-type: none"> <li>• Gigabit Ethernet card (Intel® Series) for detector interface               <ul style="list-style-type: none"> <li>▫ Speed: 1Gbps</li> <li>▫ Jumbo Frames: supports 9k</li> <li>▫ Receive Descriptors: 1024</li> <li>▫ 820.11ac (Wireless)</li> </ul> </li> <li>• LAN Card for network interface (Separate)</li> </ul>	<ul style="list-style-type: none"> <li>• 1000BASE-T Gigabit Ethernet card (Intel® I210 Series) for detector interface               <ul style="list-style-type: none"> <li>▫ Speed: 1Gbps or faster</li> <li>▫ Jumbo Frames: supports 9K</li> <li>▫ Receive Descriptors: 2K or more</li> <li>▫ 820.11ac (Wireless)</li> </ul> </li> <li>• LAN Card supporting Gigabit for network interface (Separate)</li> </ul>
CPU	CPU that supports AVX and SSE	Intel® i5 Quad Core @2.5 GHz or more (or compatible with CPU)
Memory	4GB	8GB or more
HDD or SSD	1TB	2TB or more
Monitor	General Monitor : 1280 x 800	General Monitor : 1920 x 1080, 2560 x 1440 Surface pro 4 : 2160 x 1440
CD-ROM		CD or DVD Reader / Writer

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## 3. System Configuration

---

This section gives information about the various connection / configuration ways among the detector, SCU, workstation, X-ray generator.

You can figure out the most suitable way of system configuration through this instruction.

Detector Connection Methods

Diagram

System Configuration

Changing the Wireless Setting

Generator Interface


Configuring DR Trigger Interface

Exposure Mode

## 3.1 Detector Connection Method

You can organize wired or wireless connection between the **VIVIX-S 4386W** detector and SCU. You can also choose other connection ways suitable for the use environment in case of need.

### 3.1.1 Wireless Connection

- **VIVIX-S 4386W** and PC (Workstation) connect wirelessly to transfer video and data to a PC.
- Using the SCU mini (FXRS-04A), you can use the DR Trigger interface as well as AED (Auto Exposure Detection) interface with the generator interface.
- Instead of using the SCU mini (FXRS-04A) as an AP, you can use an external AP. In this case, DR Trigger cannot be used.
- If your PC (Workstation) supports Wi-Fi, you can connect directly to your PC using the detector as an AP without a separate AP (SCU, external AP). In this case, DR Trigger cannot be used.
- In a wireless connection, the detector is powered by a battery basically.


- We recommend you use the wireless connection method for the following cases.
  - When using the unfixed detector.
  - When the wired communication is in trouble by the interruption of cables.
- The wireless connection is more freely to operate devices than the wired connection.
- Use the wireless connection with a laptop computer to enhance mobility.
- Nous vous recommandons d'utiliser la méthode de connexion sans fil dans les cas suivants.
  - Lors de l'utilisation du détecteur non fixé.
  - Lorsque la communication filaire est en difficulté par l'interruption des câbles.
- La connexion sans fil permet de faire fonctionner des appareils plus librement que la connexion filaire.
- Utilisez la connexion sans fil avec un ordinateur portable pour améliorer la mobility.

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### 3.1.2 Wired Connection

- The **VIVIX-S 4386W** detector is wired to the PC (Workstation) via the SCU and transmits image and data to the PC.
- Using the SCU mini (FXRS-04A), you can use the DR Trigger as well as AED (Auto Exposure Detection) with the generator interface.
- When the detector is wired, the battery charges while the detector is powered by the tether interface cable.

- We recommend you the wired connection in the following cases.
  - If the detector needs a constant power supply.

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
## 3.2 System Configuration

This section gives information about the configuration modes of **VIVIX-S 4386W** detectors, SCU, workstation and the generator.

### 3.2.1 SCU AP Mode

In SCU AP mode, the detector and the SCU mini (FXRS-04A) are connected wirelessly, and the SCU mini and PC are connected by wired.

### SCU mini & Detector

- SCU mini and the detector are connected wirelessly. In this case, SCU operates as **AP** (Access Point) and detector as **STATION**.
- Basically, the detector is powered by a battery.

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### SCU mini & PC (Workstation)

- SCU mini and PC (Workstation) are connected to the LAN cable.

### SCU mini & X-ray Generator


- **DR Trigger** interface is available if the SCU mini and the X-ray generator are wired via the generator interface cable.

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### 3.2.2 External AP Mode

External AP mode uses the external AP instead of the SCU mini (FXRS-04A) in SCU AP mode. The detector and external AP are connected wirelessly, and the SCU and PC are wired.

#### External AP & Detector

- The external AP and **VIVIX-S 4386W** detectors are connected wirelessly. In this case, the external AP and the detector operate as **AP** (Access Point) and **STATION** respectively.
- By default, the detector is powered by a battery.

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### External AP & PC (Workstation)


The external AP and PC are connected to a LAN cable.

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### 3.2.3 Detector AP Mode

Detector AP mode is a method of wirelessly communicating with a PC using a detector as an AP without a separate AP (SCU, external AP).

### Detector & PC (Workstation)

- The PC and **VIVIX-S4386W** detectors are connected wirelessly. In this case, the detector operates as an AP (access point) and the PC operates as a station.
- By default, the detector is powered by a battery.

- PC (Workstation) should support Wi-Fi.
- It cannot be wired to the X-ray generator due to non-use of SCU, and **DR Trigger** interface is unavailable.
- You can use this mode to move or carry the product.

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
### 3.2.4 Tether Interface Mode

The tether interface mode is a wired connection between the detector and the SCU mini (FXRS-04A) or SCU Lite (FXRP-02A), and the SCU and PC are also wired.

#### SCU & Detector

- SCU and the **VIVIX-S 4386W** detectors are connected with the tether interface cable.
- Power is supplied from the SCU via the tether interface cable and the battery is charged.

#### SCU & PC (Workstation)

- SCU and PC are wired using LAN cable.

#### SCU & X-ray Generator


- **DR Trigger** interface is available if the X-ray generator and SCU mini (FXRS-04A) are wired via the generator interface cable.

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### 3.2.5 Detector Stand-Alone Mode


You can use the **Stand-Alone** mode without making a wired or wireless connection.

- Only **AED** trigger interface and software trigger are available. (**DR Trigger** mode is not available.)
- The acquired images are stored in the internal storage of the detector in order and the images can be transmitted to the PC after connecting the detector to the PC.
- The transmitted images are removed automatically from the detector.
- Seules l'interface de déclenchement AED et le déclencheur logiciel sont disponibles. (Le mode de déclenchement DR n'est pas disponible.)
- Les images acquises sont stockées dans la mémoire interne du détecteur dans l'ordre et les images peuvent être transmises au PC après avoir connecté le détecteur au PC.
- Les images transmises sont automatiquement supprimées du détecteur.

### 3.3 Changing the Wireless Setting

In a wireless connection, the **VIVIX-S 4386W** detector provides an external AP button to facilitate changing wireless settings. You can select the function of the AP button in **VIVIX Setup** as follows:

Function	Description
<b>Detector AP / Station</b>	Switch the detector from <b>AP</b> to <b>STATION</b> or from <b>STATION</b> to <b>AP</b> .  If you press the AP button for about 3 seconds while the detector is in STATION, it scans the neighboring APs and automatically changes to the one with the highest signal strength among the detected APs.
<b>Preset Switching</b>	When changing an AP, change the detector setting that is set for use when connecting with the AP.

In addition, there is a function to synchronize wireless settings over a wired connection.

#### 3.3.1 Detector AP / Station

The AP button allows you to easily switch the detector from STATION mode (SCU AP mode or external AP mode) to AP mode, or from AP mode to STATION mode.


##### Use the AP Button

- 1 In the wireless connection, when the detector is STATION, press the AP button for about 3 seconds to switch the detector to the AP. Also, when the detector is an AP, pressing the AP button for about 3 seconds will switch the detector to STATION.
- 2 During switching, the OLED display shows '**AP Mode Start**' or '**STA Mode Start**'. After switching, you can check the changed settings.

## Use VIVIX Setup

- 1 In the Detector Configuration window of the VIVIX Setup, change the AP settings.

The screenshot shows the 'Detector Configuration' window for VIVIX Setup. The 'AP Button Type' is set to 'AP / Station'. Under the 'Access Point' section, the 'Enable' option is set to 'Off'. The 'Country' is set to 'South Korea [KR]'. The 'Frequency' is set to '5 GHz'. The 'Band' is set to '40 MHz'. The 'Channel' is set to '36'. The 'SSID' is 'vivix\_ap' and the 'Key' is '1234567890'.

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
### 3.3.2 Preset Switching


#### AP Button

- Using the AP button, you can easily change the external AP when the detector is in STATION.
  - Ex) You can use this function when you have one detector and want to use it in both ROOM # 1 and ROOM # 2. In ROOM # 1, connect wirelessly to AP①, move to ROOM # 2 and press the AP button of the detector for 3 seconds to disconnect AP① and automatically connect to AP②.

- When the detector is in STATION, pressing the AP button for about 3 seconds will scan the surrounding AP and automatically change to the AP with the largest signal strength.
- In order to change AP automatically, the SSID and KEY of the AP to be changed must be registered as a preset in the detector. If you want to connect the detector to an AP in multiple places, you must register the list of AP devices to be connected as a preset list. AP devices not registered in the preset list cannot be automatically connected even if the AP button is pressed.

- Preset list can be registered as follows using VIVIX Setup program.

- The first two letters of the preset name are displayed as the preset identifier on the OLED default screen. Therefore, it is recommended to set the preset name so that it can be distinguished by the first two letters.
  - Ex) Preset names of X-ray ROOM #1 and X-ray ROOM #2 is set to R1 and R2 each.
- After registering the preset list, it is recommended to use the Preset Load function manually to synchronize with the current settings.
- Preset contains the following information as well as the SSID and KEY values of the AP to be connected. Therefore, other settings can be automatically changed according to the room environment where the external AP is located.

Item	Description
<b>Network</b>	IP address, Netmask, Gateway
<b>Wireless Network</b>	SSID and KEY of SCU / external AP Wireless only option
<b>Power Mode</b>	Whether to use Sleep, Shutdown, and Power control
<b>Exposure Mode</b>	DR Trigger, AED
<b>Image Transmission</b>	The time when transmission times out Preview option

### Preset Switching

- 1 While scanning a neighboring AP, '**AP scan**' is displayed on the OLED display.
- 2 If there is an AP registered in the preset list among the neighboring APs, it is connected to the corresponding AP, excluding the last connected AP, and '**Change AP**' is displayed on the OLED display.
- 3 When the AP connection is complete, the Preset identifier of the AP is displayed blinking on the OLED default screen.
- 4 The blinking of the preset identifier lasts 30 seconds. If you press the AP button again for 3 seconds while blinking, AP scanning will not be performed again, and the AP will be connected to the AP with the highest signal strength among the currently scanned AP devices except the currently connected AP.
- 5 If there are no more accessible APs among the existing AP devices, the Preset identifier will stop blinking. If the AP button is pressed while the flashing stops, the peripheral AP scan is performed again.

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### 3.3.3 Sync Wireless Settings

You can easily switch from tether interface mode to SCU AP mode.

- 1 In the tether interface mode, press the AP button for about 3 seconds to change the detector to the same SSID and KEY as the wired SCU.
- 2 During synchronization, the OLED display shows '**Sync WLAN**'. After switching, you can check the changed settings.

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
## 3.4 Generator Interface

VIVIX-S 4386W detector provides **AED** interface, **DR Trigger** interface, and **SW Trigger** interface as a generator interface method to acquire images by detecting X-ray.

Mode	Description
<b>AED</b>	The detector detects X-ray exposure from the generator automatically and then performs image acquisition without any cable connection.
<b>DR Trigger</b>	The detector and generator receive and send their signal to each other for image acquisition. <ul style="list-style-type: none"> <li>• SCU mini (FXRS-04A) and X-ray generator should be connected to a generator interface cable.</li> </ul>
<b>SW Trigger</b>	The detector and generator receive and send their pre-arranged signal via the workstation for image acquisition

### 3.4.1 AED (Auto Exposure Detection) Interface

If VIVIX-S 4386W detectors are used as the **AED** interface, you can acquire images without connecting the generator to the detector with a generator interface cable. The detector detects the X-ray shot from the generator and automatically starts shooting.

- Make sure to follow the operating environmental condition (Temp: 10°C ~ +35°C).
- Do not give impact to the product. If it receives strong jolt, unwanted images may be acquired without the X-ray exposure because of the malfunction of the **AED** sensor.
- You may not acquire images or horizontal artifacts may occur depending on the external environment such as exposure condition, thickness of object or the use of grid.
- When you set X-ray exposure area to the direction of the detector, the center of the detector should be involved in the X-ray exposure area. Otherwise, you may not acquire an image.
- Assurez-vous de respecter les conditions environnementales de fonctionnement (Temp: 10: ~ +35°C).

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### 3.4.2 DR Trigger Interface


To form **VIVIX-S 4386W** detector as a **DR Trigger** interface, the generator should transmit EXP\_REQ signal to the detector through SCU. The detector sends the EXP\_OK signal to the generator once it is ready for exposure, and the generator synchronizes the signal to shoot X-ray.

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### 3.4.3 SW Trigger Interface

To form **VIVIX-S 4386W** detector as a **SW Trigger** interface, the generator should transmit EXP\_REQ signal to the detector through workstation. The detector sends the EXP\_OK signal to the generator through workstation once it is ready for exposure, and the generator synchronizes the signal to shoot X-ray.

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## 3.5 Configuring DR Trigger Interface

The detector and the X-ray generator should be connected by the generator interface cable to send or receive the synchronizing signal for X-ray exposure. Connect the one end of the generator interface cable to the **EXT\_INT** port of SCU and the other end to the X-ray generator.


In DR Trigger interface, EXP\_REQ and EXP\_OK should be all connected between the SCU and the generator, whereas only EXP\_REQ should be connected in Passive Trigger interface.

### 3.5.1 Trigger Interface Way

SCU Mini (FXRS-04A) can be wirelessly or wired connected to multiple detectors. If connecting multiple detectors to SCU Mini (FXRS-04A), the connection modes of generator interface cable vary by different Trigger interfaces.

Trigger Interface	Description
Packet Trigger	Handles exposure signal by exchanging packet.
Line Trigger	Handles exposure signal by exchanging electrical signal.

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

### 3.5.4 EXT\_INF Port Pin Map

EXT\_INF port is D-SUB 15 pin Female connectors of SCU mini (FXRS-04A) and used for the generator interface.

#### EXT\_INF1 port pin map definition (1 ~ 15)

No.	Signal name	I/O	Type	Color	Ch.	Description
1	EXP_REQ+_A	Input	Contact	Red	A	Receives EXP_REQ
2	EXP_REQ-_A	Input	Contact	Black	A	Returns signal from EXP_REQ+_A
3	EXP_REQ_TTL_A	Input	TTL	Orange	A	Receives EXP_REQ
4	EXP_REQ_GND_A	Input	TTL	Gray	A	Returns signal from EXP_REQ_TTL_A
5	EXP_OK_POWER_A/B	Input	-	Yellow	A/B	Power of TTL signal
6	EXP_OK+_A	Output	-	Green	A	Sends EXP_OK
7	EXP_OK-_A	Output	-	Brown	A	Returns signal from EXP_OK+_A
8	EXP_OK+_B	Output	-	Blue	B	Sends EXP_OK
9	EXP_OK-_B	Output	-	Pink	B	Returns signal from EXP_OK+_B
10	Reserved	-	-	-	-	Reserved for test only.
11	EXP_REQ+_B	Input	Contact	White	B	Receives EXP_REQ
12	EXP_REQ-_B	Input	Contact	Purple	B	Returns signal from EXP_REQ+_B
13	EXP_REQ_TTL_B	Input	TTL	White/Red	B	Receives EXP_REQ
14	EXP_REQ_GND_B	Input	TTL	White/Black	B	Returns signal from EXP_REQ_TTL_B
15	Reserved	-	-	-	-	Reserved pin for test only.

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

### 3.5.5 Input / Output Circuit

#### Exposure Request Input Circuit (EXP\_REQ)


#### Exposure Respond Output Circuit (EXP\_OK)

## 3.6 Exposure Mode

This mode is specialized to expose the general 2D Radiography.

### 3.6.1 AED (Auto Exposure Detection) Interface

#### Signal Processing Order

- 1 X-ray generator exposes X-ray.
- 2 Once the detector automatically detects the X-ray, it acquires images and sends the image data.

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
### 3.6.2 DR Trigger Interface

#### Signal Processing Order

- 1 Detector receives **EXP\_REQ** signal from the generator.
- 2 Detector sends **EXP\_OK** signal to the X-ray generator once it is ready for image acquisition.
- 3 X-ray generator exposes X-ray after checking the **EXP\_OK** signal
- 4 Detector sends the image data after image acquisition.

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
### 3.6.3 SW Trigger Interface

#### Signal Processing Order

- 1 Detector receives EXP\_REQ signal from the generator.
- 2 Detector sends EXP\_OK signal to the X-ray generator once it is ready for image acquisition.
- 3 X-ray generator exposes X-ray after checking the EXP\_OK signal
- 4 Detector sends the image data after image acquisition.

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## 4. Settings

---

This section gives information on installation, setting, and diagnosis to use a detector with SCU.

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
- 1 If you use the **DR Trigger** interface, connect the one end of generator interface cable to the **EXT\_INF** port of SCU, and the other to the X-ray generator.

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

- 2 Connect one end of the LAN cable to one of the LAN ports of SCU mini, and the other to the LAN Card Connector of workstation assigned for data transfer. Stand the antenna of SCU mini upright.
- 3 When in a wired state, use the tether interface cable to link the detector to the SCU mini port.

### Connection to the tether interface on the side

### Connection to the tether interface on the rear

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
4 Connect AC-DC adaptor to the power input port of SCU mini to supply power.

### 4.1.1 Connecting SCU mini (FXRS-04A)

- 1 Turn on the switch in front of the SCU mini (ON).
- 2 Check if Status LED is turned in green.

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### 4.1.2 Booting up the Detector

- 1 Press the power button on the detector for about a second. But Power Control in Settings is set to By SCU, the power is automatically turned on when wired.
- 2 Check that the detector's power LED is lit green.
- 3 Check the battery status (remaining, charged) and communication connection status (wired / wireless, detector AP / STATION, wireless connection strength, etc.) on the detector's OLED display.

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### 4.1.3 Installing Wireless Power Transmitter (FXRC-05A)


#### How to install the FXRC-05A in the detector

- Alignment between equipment:  
Wireless power transmission antennas in the wireless power transmitter and in the detector should bring into a line. Also, the infrared LED of the detector should be aligned with the infrared optic acquisition of the wireless power transmitter. The range of misalignment between 2 equipment should be up to 5mm.
- Separation distance between equipment:  
Wireless power transmitter and the detector should be positioned close together but should not be in contact with each other. The separation distance between these two equipment should be between 1mm and 3mm.
- Equipment materials:  
Wireless power transmitter uses magnetic induction and infrared interface. Therefore, it should have the material that magnetic field and infrared rays can pass through. Excluding the fixing hole, it is recommended to form layer of air or use plastic, which infrared ray can well pass through, as a material.

- If there is misalignment between equipment, the wireless power transmitter becomes less efficient, and heat increases. If misalignment between two equipment goes beyond permissible range ( $\leq 5\text{mm}$ ), the successful power transmission is not guaranteed.

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## 4.2 Device Setting

### 4.2.1 Software Installation

- 1 After connecting all devices, prepare the following software to set, calibrate and operate the detector with PSU or SCU.

Software	Description
VIVIX Device Driver (VDD)	Image filter driver for acquiring images from a detector.
VIVIX Setup	A program for setting and managing the detector / SCU.

- 2 Install **VIVIX Device Driver** and **VIVIX Setup** in sequence.

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- 3 Configure environment for the workstation.

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### 4.2.2 Setting Detector and PSU / SCU

- 1 After executing **VIVIX Setup**, access to the detector and SCU to set each device properly.
- 2 Perform detector calibration to acquire images suitable for the installation environment.
- 3 Take radiographic images to check if the shooting is conducted normally.

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## 4.3 Diagnosis of Devices

Execute the **VIVIX Setup** program to check if there is any problem to operate the detector / SCU after installing and setting devices.

### Diagnosis Items


Items	Description
<b>Image</b>	Diagnoses the acquired images.
<b>Battery Pack</b>	Diagnoses the condition of a battery pack.
<b>Wired/Wireless Communication Status</b>	Diagnoses the status of wired/wireless communication.
<b>Wired/Wireless Communication Speed</b>	Diagnoses the speed of wired/wireless communication.
<b>Self-diagnosis</b>	Diagnoses defects of a detector by self-diagnosis.

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### 4.3.1 Image Diagnosis

- 1 Execute **VIVIX Setup** and move to the **Image** dialog.
- 2 Take an image and check if it has any problem.
- 3 Take a dark image and check if it has any problem.
- 4 Check the effective area and whole area of the image.

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### 4.3.2 Diagosing Battery Pack

#### Check from OLED display of the detector

Check the mounting status / charging status of a battery, and battery levels on the OLED display.

Icon	Description
	No battery pack, Battery capacity is less than 5%
	Charging (Battery capacity changes depending on remaining charge.)
	Charging complete or discharged. (Battery capacity changes depending on remaining charge.)

#### Check the charge status of detector with charge status LED

LED Status	Description
Off	No charging power
Lights up orange	Charging
Lights up green	Charging completed

### Check from VIVIX Setup

- 1 Execute **VIVIX Setup** and go to the **Maintenance** tab of the **Detector Maintenance** dialog.
- 2 Check the mounting status, charging status and remains of a battery.

Battery		
Remain	Low	
Gauge	0.0%	
	#1	#2
Equipped	Not Equipped	Not Equipped
Charging State	Not Charging	Not Charging
Voltage	0.0 V	0.0 V

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### 4.3.3 Diagnosing Communication Status

#### Checking on the OLED display of the detector

1 Check the detector communication status on the initial screen of OLED display.

Classification	Item	Icon	Description
Battery	With or without a battery		No battery packs. Battery capacity is less than 5%.
	Whether the battery is charged		Charging (Battery capacity changes depending on remaining charge.)
	Battery level		Charging complete or discharged. (Battery capacity changes depending on remaining charge.)
Communication	AP (Wireless)		Operates in AP Mode
	STATION (Wireless)		Preset identifier (default is 'Sy' if not using preset switching function)
			Level 5 (Link Quality: 66~70, Very Good)
			Level 4 (Link Quality: 56~65, Good)
			Level 3 (Link Quality: 41~55, Normal)
			Level 2 (Link Quality: 31~40, Bad)
			Level 1 (Link Quality: 1~30, Very Bad)
			Level 0 (Link Quality: 0, Unknown)
		Tether Interface (Wired)	
			Operates in tether interface mode (connected below 100Mbps)
System Power	Sleep Mode		Enters sleep mode

2 On the initial display of the OLED display, pressing the power button or AP button for about 1 second will switch the OLED display to the second screen. On the second screen, check the IP address of the detector.

Classification	Item	Display	Description
Communication	IP address		Detector IP address

3 If you press the power button or AP button on the second screen of the OLED display for about 1 second, the OLED display will switch to the third screen. On the third screen, check the SSID information of the detector that is important for wireless connection.

Classification	Item	Display	Description
Communication	SSID		Displays SSID of the detector's current state (AP or STATION)

4 During certain operations, the OLED display shows the operating status as follows:

Classification	Display	Description
Operating Status		Detector booting up
		Deactivating the detector sleep mode
		Acquiring images
		Saving images in the detector
		Sending images
		Switching to detector AP mode (AP button function: Detector AP/STA)
		Switching to detector station mode (AP button function: Detector AP/STA)
		Initializing wireless setting information
		Syncing wireless information with connected SCU
		Scanning for nearby AP (AP button function: Preset switching)
	Changing external AP (AP button function: Preset switching)	

5 If a detector error is detected, an error number is displayed on the OLED initial screen. Pressing the power button or AP button for about 1 second while the error screen is displayed will switch to the screen showing the battery and communication status.

Classification	Display	Description
Error Status		Cause: FPGA error Expected malfunction: Failed image acquisition or acquired abnormal image
		Cause: ROIC error Expected malfunction: Failed image acquisition or acquired abnormal image
		Cause: Gate IC error Expected malfunction: Image acquisition failure or acquired abnormal image

Cause: AED sensor error  
 Expected malfunction: Image acquisition failed in AED mode

Cause: Internal power failure  
 Expected malfunction: Unintentional power off, failed image acquisition


Cause: Battery count error  
 Expected malfunction: Battery capacity indication error

**Check communication status in VIVIX Setup**


**Wired Communication**

- 1 Connect the detector with a cable.
- 2 Execute **VIVIX Setup** and go to the **Maintenance** tab of the **Detector Maintenance** dialog.
- 3 Check the status of wireless communication from the **Network** item.

**Wireless Communication**

- 1 Connect the detector wirelessly.
- 2 Execute **VIVIX Setup** and go to the **Maintenance** tab of the **Detector Maintenance** dialog.
- 3 Check the status of wireless communication from the **Network** item.

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- In case of using the detector under wireless communication, be sure to check the communication status before starting to use the detector. If the status is bad, the speed of acquiring images will be very slow or you may fail to acquire images.
- Be sure to check the surrounding wireless communication to prevent communication interference. If wireless communication module in the detector has any problems, contact the service engineer in Vieworks.
- En cas d'utilisation du détecteur en communication sans fil, assurez-vous de vérifier l'état de la communication avant de commencer à utiliser le détecteur. Si l'état est mauvais, la vitesse d'acquisition des images sera très lente ou vous risquez de ne pas acquérir d'images.
- Assurez-vous de vérifier la communication sans fil environnante pour éviter les interférences de communication. Si le module de communication sans fil du détecteur rencontre des problèmes, contactez le technicien de maintenance de Vieworks.

#### 4.3.4 Diagnosing Communication Speed

- 1 Execute **VIVIX Setup** and go to the **Maintenance** tab of the **Detector Maintenance** dialog.
- 2 Click on the **Download** button in **Image Transmission Time** and check image transmission speed of the detector.
- 3 Click on the **Start** button in **Throughput Measurement** to check the data traffic per setting time.

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
### 4.3.5 Self-Diagnosis

- 1 Execute **VIVIX Setup** and go to the **Self Diagnosis** tab in the **Detector Maintenance** dialog.
- 2 Check the desired items to diagnose from the **Test Case** list.
- 3 Click on the **Execute** button to perform self-diagnosis.
- 4 Check the status and result of diagnosis for each item in the **Result** dialog.

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## Self-diagnosis items of detector and measures

### Voltage

Item	Form	Expected problem	Measures
MAIN	Decision	Problem with input power	Replace the battery pack. Replace the tether interface cable. Replace the wireless power transmitter
		Problem with system power	Contact a service engineer.

### Battery

Item	Form	Expected problem	Measures
Detection	Decision	The battery is not attached.	Check if a battery is attached or not.
		A defective circuit is connected to a battery pack.	Contact a service engineer.

### Wired

Item	Form	Expected problem	Measures
Detection	Decision	Defective ethernet PHY	Contact a service engineer.
Connection	Decision	The cable is connected incorrectly.	Reconnect the tether interface cable.
		Defective cable	Replace the tether interface cable.
		Defective ethernet PHY	Contact a service engineer.

### Wireless

Item	Form	Expected problem	Measures
Detection	Decision	Defective wireless module	Contact a service engineer.
Connection	Decision	Inconsistent environment of the wireless communication.	Check wireless communication environment.
		Defective wireless module	Contact a service engineer.

### Sensor

Item	Form	Expected problem	Measures
Impact Sensor	Decision	Defective shock sensor	Contact a service engineer.
Temperature	Decision	Defective temperature sensor	Contact a service engineer.
AED	Decision	Defective AED sensor	Contact a service engineer.

### Memory

Item	Form	Expected problem	Measures
Detection	Decision	Impossible to save backup images.	Contact a service engineer.
		Impossible to save logs.	Contact a service engineer.
		The calibration data is inapplicable.	Contact a service engineer.

<b>File system</b>	Decision	Impossible to save backup images.	Contact a service engineer.
		Impossible to save logs.	Contact a service engineer.
<b>Status</b>	Information	N/A	N/A

**IC**

Item	Form	Expected problem	Measures
<b>FPGA</b>	Decision	Impossible to take images from a detector.	Contact a service engineer.
<b>Fuel Gauge</b>	Decision	Impossible to check the remaining of a battery pack.	Contact a service engineer.

**Self-diagnosis items of SCU and measures****Wired**

항목	Form	Expected problem	Measures
<b>Detection</b>	Decision	Defective ethernet PHY	Contact a service engineer.

**Wireless**

Item	Form	Expected problem	Measures
<b>Detection</b>	Decision	Defective wireless module	Contact a service engineer.

**Memory**

Item	Form	Expected problem	Measures
<b>Detection</b>	Decision	Impossible to save logs.	Contact a service engineer.
<b>File system</b>	Decision	Impossible to save logs.	Contact a service engineer.
<b>Status</b>	Information	N/A	N/A

**IC**

Item	Form	Expected problem	Measures
<b>Switching IC</b>	Decision	Impossible to connect the detector and PC.	Contact a service engineer.
<b>Current Controller</b>	Decision	Impossible to block overcurrent when using the wired mode.	Contact a service engineer.

---

## 5. Maintenance

---

This section gives information about maintenance of the product.

Product Initialization  
Detector Power Save Function

## 5.1 Product Initialization

If connection status of the system is not stable or setting value is not correct, you can initialize the products.

### 5.1.1 SCU Initialization

- 1 Click on the **Configuration** button of SCU after running the **VIVIX Setup** program.
- 2 Click on the **Factory Reset** button in the **Configuration** dialogue.
- 3 Wait for SCU to be initialized and rebooted automatically.
- 4 Make sure the SCU is initialized.

#### Default value of SCU initialization

Item	Default Value
<b>Network</b>	
<b>IP Address</b>	169.254.2.100
<b>Subnet Mask</b>	255.255.0.0
<b>Gateway</b>	169.254.2.100
<b>AP</b>	
<b>AP On/Off</b>	ON
<b>Frequency</b>	5GHz
<b>Country</b>	KR
<b>Band</b>	40MHz
<b>Channel</b>	+36
<b>SSID</b>	vivix
<b>Key</b>	1234567890
<b>Trigger</b>	
<b>Method</b>	Packet
<b>Polarity</b>	Auto

## 5.1.2 Detector Initialization

- 1 Click the **Maintenance** button after executing the **VIVIX Setup** program.
- 2 Go to **Detector Maintenance** window > **Maintenance** tab and click the **Factory Reset** button.
- 3 Wait for the detector until it is rebooted automatically.
- 4 Click the **Configuration** button and check if detector initialization is complete.

### Default value of detector initialization

Item	Default Value
<b>Network</b>	
IP Address	169.254.1.10
Subnet Mask	255.255.0.0
Gateway	169.254.2.100
<b>Wireless Network</b>	
SSID	vivix
Key	1234567890
AP Button	AP switching
Wireless Only	OFF
<b>Power Mode</b>	
Sleep	OFF
Shut down	OFF
Power Control	By Detector
<b>AP</b>	
Enable	OFF
Frequency	5GHz
Country	KR
Band	80MHz
Channel	36
SSID	vivix_ap
Key	1234567890
<b>Image Acquisition</b>	
Exposure Mode	AED
<b>Image Transmission</b>	
Timeout	20sec
Option	None

### 5.1.3 Wireless Initialization of Detector

- 1 Connect the detector wirelessly without connecting the tether interface cable.
- 2 Press and hold both the power button and AP button over 3 seconds.
- 3 Initialization will be conducted when the green LED is flickering. After that, the detector will reboot automatically.
- 4 When booting is complete, make sure the detector is initialized.

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## 5.2 Detector Power Save Function

You can extend the battery life by using the power saving function of the detector.

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### Types of Power Save Mode

Item	Description
<b>Normal</b>	The detector can be operated and can take images at any time.
<b>Sleep</b>	The detector cannot be operated. User can take an image by disabling the sleep mode.
<b>Shut Down</b>	The detector has been turned off. User can take an image after the detector is rebooted.

### Setting Items of Power Save Function

Item	Description
<b>Sleep</b>	<p>Sets whether you use the sleep mode of the detector or not.</p> <ul style="list-style-type: none"> <li>• If the detector is not used for the set time, it switches to Sleep mode if the Sleep function is set. (Off / After 10min / After 15min / After 20min / After 25min / After 30min)</li> </ul>
<b>Shut Down</b>	<p>Sets whether you use the shut down function of the detector or not.</p> <ul style="list-style-type: none"> <li>• If the detector is not used during the set time, the detector is turned off. (Off / After 30min / After 60min / After 90min / After 120min)</li> </ul>

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### Entry Condition of Power Save Mode

Item	Description
<b>Normal</b>	-
<b>Sleep</b>	The detector turns to sleep mode if not used for the setting time.
<b>Shut Down</b>	The detector is turned off if not used for the setting time under the sleep mode. If the Sleep function is not used, the detector will shut down if it is not used for a set time.

### Checking Indication in Power Save Mode

Item	Description
<b>Normal</b>	The power indicator LED is on and the OLED display shows the initial screen. <ul style="list-style-type: none"> <li>The power indicator LED is on and the OLED display shows a sleep mode icon.</li> </ul>
<b>Sleep</b>	<ul style="list-style-type: none"> <li>You can check the sleep status in <b>VIVIX Setup</b> or <b>VXvue</b>.</li> <li><b>VIVIX SDK</b> notices <b>Sleep</b> state.</li> </ul>
<b>Shut Down</b>	Power indicator LED and OLED display are off.

**Shows Sleep Mode on OLED Display**

Item	Icon	Description
<b>Sleep</b>		Enters sleep mode
<b>Changing Settings</b>		Cancels sleep mode

**Disabling Power Save Function**

Item	Description
<b>Normal</b>	-
<b>Sleep</b>	<ol style="list-style-type: none"> <li>You can turn off the sleep state in VIVIX Setup or VXvue.</li> <li>You can call the function for turning off the sleep state from VIVIX SDK.</li> <li>You can cancel the sleep state by pressing the detector power button or the AP button.</li> </ol>
<b>Shut Down</b>	Reboot the detector by pressing a power button on the detector.

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## 6.1 Troubleshooting

### 6.1.1 Troubleshooting Guide

When you encounter problems while using the equipment, search for the table below for the problem or error messages and try the solutions. If the problem persists, turn off the detector and consult your sales representative or a distributor. Please refer to the details of the following symptoms or error messages.

- Troubleshooting must be performed by service engineer who is authorized by Vieworks. If an unqualified person performs troubleshooting on the system resulting in damaging the detector, software, or hardware, then the Vieworks or its representative is not responsible for the detector repair regardless of remain warranty. For more detailed information, refer to <7.1 Service Information> and <7.2 Warranty>.
- Le dépannage doit être effectué par un ingénieur de service autorisé par Vieworks. Si une personne non qualifiée effectue un dépannage sur le système entraînant un endommagement du détecteur, du logiciel ou du matériel, alors Vieworks ou son représentant n'est pas responsable de la réparation du détecteur, quelle que soit la garantie. Pour plus d'informations, reportez-vous aux sections <7.1 Service Information> et <7.2 Garantie>.

### 6.1.2 Fail to Turn the Detector On

Category	Description
Symptom	<ul style="list-style-type: none"> <li>• Failed to turn the power of the detector.</li> </ul>
Expected Causes	<ul style="list-style-type: none"> <li>• The battery pack is discharged.</li> <li>• The battery pack is attached wrongly.</li> <li>• The battery pack is broken.</li> <li>• The detector is broken.</li> </ul>
Solutions	<ol style="list-style-type: none"> <li>1 Charge a battery pack</li> <li>2 Attach a battery pack again.</li> <li>3 Replace with another battery pack.</li> <li>4 Replace with another detector.</li> </ol>

### 6.1.3 The Power Switch of SCU or Status LED is not worked normally

Category	Description
Symptom	<ul style="list-style-type: none"> <li>• The power switch of SCU is not working.</li> <li>• The status LED of SCU is not responding.</li> </ul>
Expected Causes	<ul style="list-style-type: none"> <li>• Problems with the connection of power cable.</li> <li>• Power cable is broken.</li> <li>• AC-DC adapter is broken.</li> <li>• LAN cable is broken.</li> <li>• SCU is broken.</li> </ul>

<b>Solutions</b>	1 Check the cable connections (AC-DC adapter, power cable, LAN cable) of the SCU.
	2 Turn on and off SCU.
	3 Replace the power cable.
	4 Replace the AC-DC adapter.
	5 Replace the LAN cable.
	6 Replace SCU.

#### 6.1.4 Communication Test is failed

Category	Description
<b>Symptom</b>	<ul style="list-style-type: none"> <li>The communication connection fails or is lost.</li> <li>A transmission error occurs, and abnormal images or data are sent.</li> </ul>
<b>Expected Causes</b>	<ul style="list-style-type: none"> <li>Network connection problem</li> <li>Network setting problem</li> <li>PC environment setting problem</li> <li>Wireless environment problem</li> <li>SCU (or external AP) failure</li> <li>Detector failure</li> </ul>
<b>Solutions</b>	<ol style="list-style-type: none"> <li>1 Check the connection of network cable between workstation and SCU (or external AP).</li> <li>2 Check if the accurate network cable is used or not. (CAT 5E or 6)</li> <li>3 Set the network information of workstation, SCU (or external AP) and detector again.</li> <li>4 Set the workstation environment again such as firewall setting.</li> <li>5 Check surrounding wireless communication environment. If there is an environment that prevents wireless communication, remove the obstacle or change the wireless communication settings to prevent the interference.</li> <li>6 Reboot the detector and SCU (or external AP)..</li> <li>7 Replace SCU (or external AP).</li> <li>8 Replace the detector.</li> </ol>

#### 6.1.5 Communication is connected, but no image is acquired

Category	Description
<b>Symptom</b>	<ul style="list-style-type: none"> <li>The image cannot be sent to PC even though the X-ray is exposed.</li> </ul>
<b>Expected Causes</b>	<ul style="list-style-type: none"> <li>Problems with the generator interface</li> <li>Problems with the generator interface connection</li> <li>Generator interface cable failure</li> <li>SCU failure</li> <li>Detector failure</li> </ul>
<b>Solutions</b>	<ol style="list-style-type: none"> <li>1 Check the generator interface setting of the detector. (AED mode, DR Trigger mode)</li> </ol>

- 2 In DR trigger mode, check the generator interface connection.
- 3 Replace the generator interface cable.
- 4 Replace the SCU.
- 5 Replace the detector.

### 6.1.6 Getting Abnormal Images

Category	Description
<b>Symptom</b>	<ul style="list-style-type: none"> <li>• The image sent to the PC is abnormal.</li> </ul>
<b>Expected Causes</b>	<ul style="list-style-type: none"> <li>• Problems with calibration data</li> <li>• Problems with setting exposure time</li> <li>• Wrong operating environment</li> <li>• Detector failure.</li> </ul>
<b>Solutions</b>	<ol style="list-style-type: none"> <li>1 Perform the calibration again.</li> <li>2 Make sure that the detector is set to be applied calibration data.</li> <li>3 Set the exposure time to be longer than the X-ray exposure time.</li> <li>4 Check the operating environment (temperature, humidity, X-ray energy range).</li> <li>5 Replace the detector.</li> </ol>

### 6.1.7 When Hours of Battery Use Decreases

Category	Description
<b>Symptom</b>	<ul style="list-style-type: none"> <li>• The fully charged battery pack runs down quickly.</li> </ul>
<b>Expected Causes</b>	<ul style="list-style-type: none"> <li>• Reduced performance (capacity) with prolonged use of the battery pack</li> <li>• Using the battery pack in low or high temperature environments</li> </ul>
<b>Solutions</b>	<ol style="list-style-type: none"> <li>1 If the battery pack is to be used for a long time, replace it with a new battery pack. (The battery pack is a consumable.)</li> <li>2 Use the battery pack at room temperature. In cold or hot environments, the battery pack's capacity decreases.</li> </ol>

### 6.1.8 When an Error Number is Displayed on the OLED Display

Category	Description
<b>Symptom</b>	<ul style="list-style-type: none"> <li>• The error number is displayed on the OLED display.</li> </ul>
<b>Expected Causes</b>	<ul style="list-style-type: none"> <li>• ERR_00 ~ ERR_04: Detector internal IC failure</li> <li>• ERR_05: Battery pack failure or insufficient battery pack</li> </ul>
<b>Solutions</b>	<ol style="list-style-type: none"> <li>1 ERR_00 ~ ERR_04: Replace the detector.</li> <li>2 ERR_05: Check the number and status of batteries installed. If normal, replace with a new battery pack.</li> </ol>

---

## 7. Information

---

This section gives overview information for service and warranty of the product.

Service Information

Warranty

Revision History

## 7.1 Service Information

### 7.1.1 Product Lifetime

The estimated product lifetime is up to seven (7) years under the appropriate regular inspection and maintenance.

### 7.1.2 Regular Inspection and Maintenance

In order to ensure the safety of patients, operating personnel and third parties, and to maintain the performance and reliability of the equipment, be sure to perform regular inspection at least once a year. If necessary, clean up the equipment, adjust, or replace consumables.

There may be cases where overhaul is recommended depending on the conditions. Contact your sales representative or distributor for regular inspections or maintenance.

### 7.1.3 Repair

If a problem cannot be solved even after taking the measures indicated in Troubleshooting and contact your sales representative or a distributor for repairs. Please refer to the name label and provide the following information.

- Product name: **FXRD-4386WB** Detector
- Serial number: 9 digit-number on the product label
- Explanation of problem: Describe it as detailed as possible.

### 7.1.4 Replacement Support

Performance parts (parts required to maintain the functioning of the product) of this product will be stocked for seven years after discontinuance of production, to allow for repair.

### 7.1.5 Consumables

The following consumable can deteriorate because of its characteristics and structure. For purchase of consumables, contact your sales representative or distributor.

- Battery pack: **FXRB-04A**

## 7.2 Warranty

Vieworks Co., Ltd. follows the contract with each customer regarding the product warranty period. If any such product proves defective during this warranty period, Vieworks at its option, either will repair the defective product without charge for parts and labor, or will provide a replacement in exchange for the defective product. In order to obtain service under this warranty, Customer must notify Vieworks of the defect before the expiration of the warranty period and make suitable arrangements for the performance of service. Customer shall be responsible for packaging and shipping the defective product to the service center designated by Vieworks with shipping charges prepaid.

Vieworks shall pay for the return of the product to customer if the shipment is to a location within the country in which Vieworks designated service center is located. Customer shall be responsible for paying all shipping charges, duties, taxes, and any other charges for products returned to any other locations.

This warranty shall not apply to any defect, failure, or damage caused by improper or inadequate maintenance and care. Vieworks shall not be obligated to furnish service under this warranty to repair damage resulting from attempts by personnel other than Vieworks or its representatives to install, repair, or service this product, to repair damage resulting from improper use or connection to incompatible equipment or power source; or to service a product that has been modified or integrated with other products when the effect of such modification or integration increases the time or difficulty of servicing the product

**THIS WARRANTY IS GIVEN BY VIEWORKS WITH RESPECT TO THIS PRODUCT IN LIEU OF ANY OTHER WARRANTIES, EXPRESS OR IMPLIED. VIEWORKS AND ITS VENDOR DISCLAIM ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. VIEWORKS RESPONSIBILITY TO REPAIR OR REPLACE DEFECTIVE PRODUCTS IS THE SOLE REMEDY PROVIDED TO THE CUSTOMER FOR BREACH OF THIS WARRANTY. VIEWORKS AND ITS VENDORS WILL NOT BE LIABLE FOR ANY INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES IRRESPECTIVE OF WHETHER VIEWORKS OR THE VENDOR HAS ADVANCE NOTICE OF THE POSSIBILITY OF SUCH DAMAGES.**

There are no warranties which extend beyond the description mentioned in this document.

**European Union**

This product can be used in which EU members, in accordance with Article 10(10) / or this product can be used in at least one EU country, in accordance with Article 10(2).

**FCC Instructions**

FCC Compliance Statement

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**FCC Interference Statement\_Class A**

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

**FCC Caution**

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

SDoC form will be attached to the manual when placing to the market.

**Industry Canada Statement**

This device complies with RSS-247/RSS-210 of the Industry Canada Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

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**Technical requirements for Wi-Fi operation in the 5 GHz band**

Wi-Fi access points that are capable of allowing your device to operate in the 5.15 - 5.25 GHz band are optimized for indoor use only. If your Wi-Fi network is capable of operating in this mode, please restrict your Wi-Fi use to indoors to not violate Canada ISED regulations to protect mobile satellite services.

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## 7.3 Revision History

Version	Date	Description
1.0	2021-08-03	• Initial release

# VIEWWORKS

### Vieworks Co., Ltd.

#### Headquarter

41-3, Burim-ro 170beon-gil, Dongan-gu, Anyang-si,

Gyeonggi-do, 14055 Republic of Korea

Telephone: +82-70-7011-6161

Fax: +82-31-386-8631

Homepage: <http://www.vieworks.com>

#### Hwaseong Site

25-7, Jeongnamsandan 2-gil, Jeongnam-myeon, Hwaseong-si,

Gyeonggi-do, 18514 Republic of Korea

Telephone: +82-70-7011-6161

Homepage: <http://www.vieworks.com>

### European Authorized Representative: Obelis s.a

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1030 Brussels, BELGIUM

Tel: +(32) 2.732.59.54

Fax: +(32) 2.732.60.03

E-mail: [mail@obelis.net](mailto:mail@obelis.net)