



Instructions for Use

Avalon Fetal Monitor

FM20/30, FM40/50, Avalon CL

Release J.3. with Software Revision J.30.xx

Patient Monitoring

PHILIPS

Table of Contents


1 Introduction	9
<hr/>	
Who this Book is For	9
Confirm Fetal Life Before Using the Monitor	10
Introducing the Avalon Family of Fetal Monitors	11
2 What's New	29
<hr/>	
What's New in Release J.3	29
What's New in Release G.0	32
3 Basic Operation	35
<hr/>	
Supported Measurements	35
Avalon FM20 and FM30	36
Avalon FM40 and FM50	39
Avalon CL Transducer System	40
Getting to Know Your Avalon FM20/FM30	40
Getting to Know Your Avalon FM40/FM50	43
Transducers	45
Getting to Know Your Avalon CL	47
Cableless Transducers	51
CL Pods	55
Operating and Navigating	56
Operating Modes	63
Automatic Screen Layouts	64
Settings	64
Preparing to Monitor	67
Switching the Monitor to Standby	98
After Monitoring	98
Disconnecting from Power	99
Power On/Power Off Behavior	99
Monitoring After a Power Failure	99
Troubleshooting	100
4 FM20/30 Battery Option	101
<hr/>	
External Power Supply M8023A	102
Using Batteries	102
Optimizing Battery Performance	105
5 Alarms	109
<hr/>	
Alarm Mode	110
Nurse Call Systems	110

Visual Alarm Indicators	110
Audible Alarm Indicators	111
Acknowledging Alarms	114
Acknowledging Disconnect INOPs	114
Alarm Reminder	114
Pausing or Switching Off Alarms	115
Alarm Limits	117
Reviewing Alarms	118
Latching Alarms	118
Testing Alarms	120
Alarm Behavior at Power On	121
6 Patient Alarms and INOPs	123
<hr/>	
Alarm Messages	123
Technical Alarm Messages (INOPs)	126
7 Admitting and Discharging	137
<hr/>	
Admit/Discharge on the Monitor	137
New Patient Check	139
OB TraceVue/IntelliSpace Perinatal	139
8 Non-Stress Test Timer	141
<hr/>	
Setting NST Autostart/Autostop	141
Viewing the NST Timer	141
Timer Expiry Notification	142
Accessing the NST Setup Pop-up Keys	142
9 Non-Stress Test Report	145
<hr/>	
Setting Up an NST Report	145
NST Report Status Window	146
NST Criteria	148
Non-Reassuring Report	149
Non-Reactive NST Test	149
10 Cross-Channel Verification (CCV)	151
<hr/>	
Misidentification of Heart Rates	151
Cross-Channel Verification Functionality	152
Overview of Cross-Channel Comparisons	152
Coincidence Examples	154
Recommended Actions for INOP Coincidence	156
11 Monitoring FHR and FMP Using Ultrasound	157
<hr/>	
Technical Description	157
Limitations of the Technology	158
Misidentification of MHR as FHR	158
What You Need	158

Cableless Monitoring - Important Considerations	159
Preparing to Monitor	160
Selecting Fetal Heart Sound	161
Changing the Fetal Heart Sound Volume	161
Fetal Movement Profile	163
Troubleshooting	165
Additional Information	166
Testing Ultrasound Transducers	173
12 Monitoring Twin FHRs	175
Important Considerations	175
Monitoring Twins Externally	176
Monitoring Twins Internally	177
Separating FHR Traces	177
Troubleshooting	181
13 Monitoring Triple FHRs	183
Important Considerations	183
Monitoring Triplets	184
Separating FHR Traces	184
"Standard" Separation Order	184
"Classic" Separation Order	185
Switching Trace Separation On and Off	186
When Trace Separation is On	187
When Trace Separation is Off	188
Troubleshooting	188
14 Fetal Heart Rate Alarms	191
Changing Alarm Settings	191
Changing Signal Loss Delay	192
15 Monitoring Uterine Activity Externally	193
What You Need	193
External Toco Monitoring	194
Toco Sensitivity	195
Troubleshooting	195
Testing Toco Transducers	195
16 Monitoring Uterine Activity Internally	197
What You Need	197
Internal (IUP) Monitoring	199
Troubleshooting	199
17 Monitoring FHR Using DECG	201
Misidentification of MHR as FHR	201
What You Need	203

Making Connections	205
Monitoring DECG	206
Suppressing Artifacts	208
Troubleshooting	209
Testing DECG Mode	209
18 Monitoring Noninvasive Blood Pressure	211
<hr/>	
Introducing the Oscillometric Noninvasive Blood Pressure Measurement	211
Preparing to Measure Noninvasive Blood Pressure	212
Starting and Stopping Measurements	214
Enabling Automatic Mode and Setting Repetition Time	215
Enabling Sequence Mode and Setting Up The Sequence	215
Choosing the Alarm Source	216
Assisting Venous Puncture	216
Calibrating NBP	216
Troubleshooting	217
19 Monitoring Maternal Temperature	219
<hr/>	
Measuring Tympanic Temperature	219
20 Monitoring SpO2	227
<hr/>	
Selecting an SpO2 Sensor	227
Applying the Sensor	227
Connecting SpO2 Cables	228
Measuring SpO2	228
SpO2 Signal Quality Indicator (Fast SpO2 only)	229
Assessing a Suspicious SpO2 Reading	230
Understanding SpO2 Alarms	230
Setting Up Tone Modulation	231
Setting the QRS Volume	231
21 Monitoring Maternal Heart / Pulse Rate	233
<hr/>	
Priority for Maternal Heart / Pulse Rate	233
Misidentification of MHR for FHR	234
MHR from MEKG Electrodes	235
Monitoring MEKG Wave	238
Pulse Rate from Toco MP	240
Pulse Rate from SpO2	241
Adjusting the Heart Rate / Pulse Alarm Limits	241
Average Pulse Rate from Noninvasive Blood Pressure	242
Testing MEKG Mode	242

22	Printing the ECG Waveform	243
23	Paper Save Mode for Maternal Measurements	249
24	Recovering Data	251
	Recovering Traces on Paper	251
	Recovering Traces on an OB TraceVue/IntelliSpace Perinatal System	252
	Recording Stored Data	252
25	Care and Cleaning	255
	General Points	255
	Cleaning and Disinfecting	256
	Cleaning and Disinfecting Monitoring Accessories	257
	Cleaning and Disinfecting the Tympanic Temperature Accessories	258
	Sterilizing	258
26	Maintenance	259
	Inspecting the Equipment and Accessories	259
	Inspecting the Cables and Cords	259
	Maintenance Task and Test Schedule	260
	Recorder Maintenance	260
	Cleaning the Print Head	265
	Disposing of the Monitor	266
27	Accessories and Supplies	267
	Information on Latex	267
	Avalon CL Base Station	268
	Transducers	268
	Fetal Accessories	268
	MECG Accessories	270
	Noninvasive Blood Pressure Accessories	270
	SpO2 Accessories	273
	Tympanic Temperature Accessories	277
	Recorder Paper	277
	Batteries	277
28	Specifications and Standards Compliance	279
	Environmental Specifications	279
	Physical Specifications	281
	Performance Specifications	285
	Recorder Specifications	294
	Compatible External Displays: FM40/FM50 Only	296
	Manufacturer's Information	296
	Trademark Acknowledgment	296
	Regulatory and Standards Compliance	297
	Environment	306



Monitoring After a Loss of Power	307
ESU, MRI and Defibrillation	307
Cardiac Pacemakers and Electrical Stimulators	307
Fast Transients/Bursts	307
Symbols on the System	308
29 Default Settings Appendix	311
<hr/>	
Alarm and Measurement Default Settings	311
Recorder Default Settings	314
Index	317
<hr/>	

Introduction

Who this Book is For

This book is for trained healthcare professionals using the Avalon FM20, FM30, FM40 and FM50 Fetal/Maternal monitors, and the cableless monitoring system Avalon CL. It describes how to set up and use the monitor and transducers. Familiarize yourself with all instructions including warnings and cautions before starting to monitor patients. Read and keep the Instructions for Use that come with any accessories, as these contain important information about application and care and cleaning that is not repeated in this book.

You should be:

- Trained in the use of fetal heart rate (FHR) monitors.
- Trained in the interpretation of FHR traces.
- Familiar with using medical devices and with standard fetal monitoring procedures.

For information on how to configure and service the monitor, see the Configuration Guide and the Service Guides, or contact your authorized service provider.

Your monitor may not have all of the features and options described in this guide. The exact appearance of the monitor may differ slightly from that shown in the illustrations.

In this guide:

- A **warning** alerts you to a potential serious outcome, adverse event or safety hazard. Failure to observe a warning may result in death or serious injury to the user or patient.
- A **caution** alerts you to where special care is necessary for the safe and effective use of the product. Failure to observe a caution may result in:
 - minor or moderate personal injury,
 - damage to the product or other property,
 - possibly in a remote risk of more serious injury.
- **Monitor** refers to the entire fetal/maternal monitor. **Display** refers to the physical display unit. **Screen** refers to everything you see on the monitor's display, such as measurements, alarms, patient data, and so forth.

- FM30** • Whenever a monitor's identifier appears to the left of a heading or paragraph, it means that the information applies to that monitor only. Where the information applies to all models, no distinction is made.

Avalon CTS Whenever one of these identifiers appear to the left of a heading or paragraph, it means that the information applies to that cableless monitoring system. Where the information applies to both **Avalon CL** systems, no distinction is made.

Confirm Fetal Life Before Using the Monitor

Fetal monitoring technology available today is *not always* able to differentiate a fetal heart rate (FHR) signal source from a maternal heart rate (MHR) source in *all* situations. Therefore, you should confirm fetal life *by independent means* before starting to use the fetal monitor, for example, by palpation of fetal movement or auscultation of fetal heart sounds using a fetoscope, stethoscope, or Pinard stethoscope. If you cannot hear the fetal heart sounds, and you cannot confirm fetal movement by palpation, confirm fetal life using obstetric ultrasonography. Continue to confirm that the fetus is the signal source for the FHR during monitoring.

Be aware that:

- an MHR trace can exhibit features that are very similar to those of an FHR trace, even including accelerations and decelerations. Do not rely solely on trace pattern features to identify a fetal source.
- Fetal Movement Profile (FMP) annotations on a fetal trace *alone* may not always indicate that the fetus is alive. The body of a deceased fetus can move and cause the monitor to annotate fetal body movements.

Here are some examples where the MHR can be misidentified as the FHR.

- **When using an ultrasound transducer:**
 - It is possible to pick up maternal signal sources, such as the aorta or other large vessels.
 - Misidentification may occur when the MHR is higher than normal (especially when it is over 100 bpm).
- **When using a fetal scalp electrode:**
 - Electrical impulses from the maternal heart can sometimes be transmitted to the fetal monitor through a recently deceased fetus via the spiral scalp electrode cable, appearing to be a fetal signal source.
 - The recorded MHR (and any artifact) can be misinterpreted as an FHR (especially when it is over 100 bpm).
- **When Fetal Movement Profile (FMP) is enabled:**

FMP annotations in the absence of fetal life may be a result of:

 - Movement of the deceased fetus during or following maternal movement.
 - Movement of the deceased fetus during or following manual palpation of fetal movement (especially if the pressure applied is too forceful).
 - Movement of the ultrasound transducer.
 - The ultrasound transducer detecting a maternal movement source, such as the mother coughing.

See also “Monitoring FHR and FMP Using Ultrasound” on page 155 and “Monitoring FHR Using DECG” on page 195.

To reduce the possibility of mistaking the MHR for an FHR, it is recommended that you monitor both maternal and fetal heart rates. The monitor's cross-channel verification (CCV) facility can help by automatically detecting when an MHR coincides with an FHR. For further details, see “Cross-Channel Verification (CCV)” on page 149.

Introducing the Avalon Family of Fetal Monitors

The Avalon family of fetal monitors consists of the Avalon FM20, FM30, FM40 and FM50. While the FM20/FM30 and the FM40/FM50 have different form factors, the method of operation is very similar for all monitors. The Avalon fetal monitors also share transducers, accessories, software, and are compatible with the Avalon CL, and Avalon CTS Cableless Fetal Transducer Systems.

Intended Use

The Philips Avalon FM20 (M2702A), FM30 (M2703A), FM40 (M2704A) and FM50 (M2705A) Fetal/Maternal Monitors are intended for:

- non-invasive monitoring of fetal heart rates and movements.
- non-invasive monitoring of maternal heart rates, maternal pulse rates, uterine activity, maternal noninvasive blood pressure, maternal oxygen saturation, and maternal temperature.
- invasive monitoring of fetal Direct ECG and intrauterine pressure and for displaying and recording of fetal and maternal ECG waves. (FM30 and FM50 only).
- displaying, storing, and recording patient data and parameter values, and for generating alarms from fetal and maternal parameters.
- transmitting patient data and parameter values to a patient information and surveillance system used by trained health care professionals.
- use in antepartum testing areas, in labor and delivery rooms, and during postpartum recovery in the hospital environment. They are not intended for use in intensive care units or operating rooms.
- monitoring in a bath or shower. (Avalon CL **ableless** transducers Toco⁺ MP, Ultrasound, and ECG/IUP only)
- transport situations in healthcare facilities, for healthcare facilities outside hospitals, such as doctors' offices, and for use in private households. (FM20 and FM30 only)

WARNING

The fetal/maternal monitors are not intended for:

- use during defibrillation, electro-surgery, or magnetic resonance imaging (MRI).
 - ECG measurements on patients connected to electrical stimulator or with cardiac pacemakers.
 - use of the invasive measurements IUP and fetal DECG, use of the patient module (M2738A) and use of the Avalon CL system in domestic establishments and those connected directly to the public low-voltage supply network that supplies buildings used for domestic purposes.
-

WARNING

No modification of the fetal monitors, transducers, and the Avalon CL base station is allowed.

CAUTION

US federal law restricts this device to sale by, or on the order of, a physician.

Connecting the Monitor to AC Mains

WARNING

- Always use the supplied power cord with the earthed mains plug to connect to an earthed AC mains socket. Never adapt the mains plug from the fetal monitor to fit an unearthed AC mains socket.
- Check that the line frequency is correctly configured in the **Global Settings** menu.
- **FM20/FM30 only:** The protective earth conductor is required for EMC purposes. It has no protective function against electric shock. Double and/or reinforced insulation protects this device against electric shock.
- Do not use AC mains extension cords or multiple portable socket-outlets.

Indications for Use

Avalon Fetal/Maternal Monitor FM20

Indicated for use by trained health care professionals whenever there is a need for monitoring of the physiological parameters uterine activity, heart rate, oxygen saturation, noninvasive blood pressure, pulse rate, and temperature of pregnant women and the fetal heart rates of single fetuses, twins, and triplets in labor and delivery rooms, in antepartum testing areas, in private households and during transports in healthcare facilities.

Avalon Fetal/Maternal Monitor FM30

Indicated for use by trained health care professionals whenever there is a need for monitoring of the physiological parameters uterine activity, heart rate, ECG, oxygen saturation, noninvasive blood pressure, and pulse rate and temperature of pregnant women and the fetal heart rates of single fetuses, twins, and triplets in labor and delivery rooms, in antepartum testing areas, in private households and during transports in healthcare facilities.

Avalon Fetal/Maternal Monitor FM40

Indicated for use by trained health care professionals whenever there is a need for monitoring of the physiological parameters uterine activity, heart rate, oxygen saturation, noninvasive blood pressure, and pulse rate and temperature of pregnant women and the fetal heart rates of single fetuses, twins, and triplets in labor and delivery rooms and in antepartum testing areas.

Avalon Fetal/Maternal Monitor FM50

Indicated for use by trained health care professionals whenever there is a need for monitoring of the physiological parameters uterine activity, heart rate, ECG, oxygen saturation, noninvasive blood pressure, and pulse rate and temperature of pregnant women and the fetal heart rates of single fetuses, twins, and triplets in labor and delivery rooms and in antepartum testing areas.

Safety Information

In this guide:

- A warning alerts you to a potential serious outcome, adverse event or safety hazard. Failure to observe a warning may result in death or serious injury to the user or patient.
- A caution alerts you to where special care is necessary for the safe and effective use of the product. Failure to observe a caution may result in minor or moderate personal injury or damage to the product or other property, and possibly in a remote risk of more serious injury.

WARNING

The fetal/maternal monitors are not intended for:

- use during defibrillation, electro-surgery, or magnetic resonance imaging (MRI).
 - ECG measurements on patients connected to electrical stimulator or with cardiac pacemakers.
 - use of the invasive measurements IUP and fetal DECG, use of the patient module (M2738A) and use of the Avalon CL system in domestic establishments and those connected directly to the public low-voltage supply network that supplies buildings used for domestic purposes.
-

WARNING

No modification of the fetal monitors, transducers, and the Avalon CL base station is allowed.

CAUTION

US federal law restricts this device to sale by, or on the order of, a physician.

WARNING

- Always use the supplied power cord with the earthed mains plug to connect to an earthed AC mains socket. Never adapt the mains plug from the fetal monitor to fit an unearthed AC mains socket.
 - Check that the line frequency is correctly configured in the Global Settings menu.
 - **FM20/FM30 only:** The protective earth conductor is required for EMC purposes. It has no protective function against electric shock. Double and/or reinforced insulation protects this device against electric shock.
 - Do not use AC mains extension cords or multiple portable socket-outlets.
-

WARNING

- To avoid magnetic interference affecting the mode of the pacemaker, ensure that the Avalon CL base station does not come into close contact with implanted pacemakers.
 - This equipment generates, uses, and radiates radio-frequency energy, and if it is not installed and used in accordance with its accompanying documentation, may cause interference to radio communications. Operation of this equipment in a residential area may cause interference, in which case the users must take whatever measures may be required to correct the interference.
-

WARNING

Do not change the date and time setting, if the fetal monitor is connected to a Philips OB TraceVue/ IntelliSpace Perinatal system via a LAN -setup. The monitor uses the OB TraceVue/IntelliSpace Perinatal system date and time, including daylight saving time changes. As long as the fetal monitor is connected to the OB TraceVue/IntelliSpace Perinatal system via the LAN-setup (locomotive icon displayed on the monitor's screen), the option to change the date and time settings at the fetal monitor are disabled, this is not valid for RS232 connections, or the connection to other systems.

WARNING

Short range radio connections are subject to interruption due to interference from other radio sources in the vicinity, including microwaves, bluetooth devices, WLAN devices (802.11b,g,n) and cordless phones. Depending on the strength and duration of the interference, the interruption may occur for an extended period. A loss of connection, due to moving out-of-range, interference, or for other reasons, is indicated with a **No Host Monitoring** INOP (here the host is the fetal monitor) on the NBP or SpO₂ Pods, or a **cl NBP Disconnect** or **cl SpO₂ Disconnect** INOP at the fetal monitor. Correct channel configuration is important, refer to the Configuration Guide for details.

CAUTION

Check the fetal monitors housing for damage before you start to monitor as part of your safety precautions.

WARNING

When connecting devices for acquiring measurements, always position cables and NBP tubing carefully to avoid entanglement or potential strangulation.

CAUTION

Using recorder paper that is not approved by Philips can result in accelerated paper fading and can damage the thermal line printhead. This type of damage is not covered by warranty.

WARNING

Never immerse the base station in liquid. You must protect it against water sprays or splashes. Place the base station where there is no chance of contact with, or falling into water or other liquids.

CAUTION

Avoid the use of pulsating water jets in the bath or shower while monitoring, as these can be misinterpreted as an incorrect (or totally artificial) heart rate.

Toco Baseline drift: The accuracy specified for baseline drift cannot be guaranteed for underwater usage. When using transducers under warm water the temperature increase causes a significant baseline change due to internal pressure increase. The depth under water at which the Toco transducer is used also has an effect on the Toco baseline, as the water pressure increases with depth. After immersion, allow one to two minutes for the pressure to stabilize, then adjust the Toco baseline (between contractions), and check it frequently.

When using the transducers underwater, the radio transmission range is reduced, and signal loss may occur.

WARNING

- Always use the supplied power cord with the earthed mains plug to connect the external power supply M8023A (option #E25) to an earthed AC mains socket. Never adapt the mains plug from the power supply to fit an unearthed AC mains socket.
 - Do not use AC mains extension cords or multiple portable socket-outlets. If a multiple portable socket-outlet without an approved isolation transformer is used, the interruption of its protective earthing may result in enclosure leakage currents equal to the sum of the individual earth leakage currents, so exceeding allowable limits.
 - Do not connect any devices that are not supported as part of a system.
 - Any non-medical device placed and operated in the patient's vicinity must be powered via an approved isolation transformer that ensures mechanical fixing of the power cords and covering of any unused power outlets.
-

WARNING

Do not touch the charging contacts for the cableless transducers at the Avalon base station while you are touching the patient.)

CAUTION

Condition the battery with a monitor not currently in use. The monitor switches off automatically when there is no battery power left.

WARNING

Use only **Philips batteries** part number M4605A. Use of a different battery may present a risk of fire or explosion.

Do not open batteries, or dispose of them in fire, or cause them to short circuit. They may ignite, explode, leak or heat up, causing personal injury.

If battery leakage should occur, avoid contact with skin. Refer to qualified and authorized service personnel.

Dispose of used batteries promptly and in an environmentally-responsible manner. Do not dispose of the battery in normal waste containers. Consult your hospital administrator to find out about local arrangements. Do not expose batteries to liquids.

Do not crush, drop or puncture batteries - mechanical abuse can lead to internal damage and internal short circuits which may not be visible externally.

If a battery has been dropped or banged against a hard surface, whether damage is visible externally or not:

- discontinue use
- dispose of the battery in accordance with the disposal instructions above.

Keep batteries out of the reach of children.

CAUTION

Do not disassemble, heat above 100°C (212°F) or incinerate the batteries, to avoid the risk of fire and burns. Keep batteries out of the reach of children and in their original package until you are ready to use them.

If battery leakage should occur, use caution in removing the battery. Avoid contact with skin. Refer to qualified and authorized service personnel.

Do not install or use pre-damaged batteries.

WARNING

Alarm systems of the monitor and those of the connected OB system are independent and not synchronized.

WARNING

In INOP only mode, no fetal/maternal patient alarms are enabled or indicated.

WARNING

Do not rely exclusively on the audible alarm system for fetal monitoring. Adjustment of alarm volume to a low level or off during monitoring may result in a dangerous situation. Remember that the most reliable method of fetal monitoring combines close personal surveillance with correct operation of monitoring equipment.

WARNING

Be aware that the monitors in your care area may each have different alarm settings, to suit different scenarios. Always check that the alarm settings are appropriate before you start monitoring.

CAUTION

In order to ensure that the settings are reset to user defaults for a new patient, always discharge the previous patient from the fetal monitor.

WARNING

Performing ultrasound imaging or Doppler flow measurements together with ultrasound fetal monitoring may cause false FHR readings, and the trace recording may deteriorate.

CAUTION

Never use ultrasound transducers connected to more than one fetal monitor on the same patient.

WARNING

- Ensure that the conductive parts of the fetal scalp electrode and the maternal leg plate electrode do not contact other conductive parts, including earth.
 - Indication of the heart-rate may be adversely affected by the operation of cardiac pacemaker pulses or by cardiac arrhythmias.
 - During ambulant FHR monitoring, the chance of losing the signal or detecting the maternal heart rate is higher than during stationary monitoring. The frequency of the patient's walk may be detected, and mistaken for a FHR signal.
 - Check the mother's pulse periodically during monitoring and compare this with the FHR signal. Beware of mistaking a "doubled" maternal heart rate for FHR. In the case of a dead fetus, there is a risk that the maternal heart rate is monitored and misinterpreted as the fetal heart rate. Therefore, the simultaneous monitoring of maternal heart rate (preferably, the maternal ECG) is encouraged.
 - Do not interpret maternal movements as fetal movements.
 - Artifacts: FMP artifacts are generated during fetal heart rate searching by changing the transducer position, therefore Philips fetal monitors enable the FMP only after detecting a valid heart rate signal for several seconds. FMP is not recommended when the mother is likely to move, and you should disable Fetal Movement Profile (FMP) at the fetal monitor (FMP Off) if the mother is walking.
 - Gaps and maternal heart rate detection can occur:
 - if the transducer is not correctly positioned.
 - due to the pulsation of uterine blood vessels.
 - if the fetus moves.
-

CAUTION

Using ultrasound gel not approved by Philips may reduce signal quality and may damage the transducer. This type of damage is not covered by warranty.

WARNING

Periodically compare the mother's pulse with the signal coming from the monitor's loudspeaker to ensure that you are monitoring fetal heart rate. Do not mistake a doubled or elevated MHR for FHR.

WARNING

Do not catheterize if placenta previa is diagnosed, or if uterine bleeding from an undetermined source is present.

WARNING

Never attempt to connect the fetal scalp electrode to anything other than the correct DECG adapter cable.

WARNING

The fetal/maternal monitor is not a diagnostic ECG device. In particular, the display of fetal/maternal ECG is intended only for evaluating signal quality for fetal/maternal heart rate as derived from the ECG waveform.

When in doubt, it can be used to identify sources of compromised signal quality, such as noise or muscle artifacts. It can subsequently be used to verify the result of measures taken to resolve them (e.g. checking ECG cable connections or adapting the fetal **ArtifactSuppress** configuration).

The safety and effectiveness of the displayed fetal/maternal ECG waveform (i.e. P, QRS and T segments) for evaluation of fetal/maternal cardiac status during labor have not been evaluated.

WARNING

Intravenous infusion: Do not use the NBP cuff on a limb with an intravenous infusion or arterial catheter in place. This could cause tissue damage around the catheter when the infusion is slowed or blocked during cuff inflation.

Skin Damage: Do not measure NBP in cases of sickle-cell disease or any condition where skin damage has occurred or is expected.

Unattended measurement: Use clinical judgment to decide whether to perform frequent unattended blood pressure measurements in cases of severe blood clotting disorders because of the risk of hematoma in the limb fitted with the cuff.

WARNING

Do not use the thermometer in the presence of flammable anesthetics, such as a flammable anesthetic mixture with air, oxygen or nitrous oxide.

WARNING

- Never apply the probe to the patient when the probe is not connected to the unit.
 - Always use a single-use probe cover to limit patient cross-contamination.
 - Measurement errors or inaccurate readings may result when:
 - probe covers other than the specified probe covers are used (see “Tympanic Temperature Accessories” on page 267)
 - Insert the probe slowly and carefully to avoid damage to the ear canal and the tympanic membrane.
 - Inspect the probe for damage, holes, tears, or sharp edges to avoid injuring the skin.
 - Always ensure that the used probe cover is removed before attaching a new probe cover.
-

CAUTION

- Do not immerse the probe in fluids, or drop fluids on the probe.
 - Do not use a probe cover that has been dropped or is damaged.
 - Do not autoclave. To prevent damage to the unit, probe and accessories, refer to the cleaning procedures in the “Cleaning and Disinfecting” on page 246 chapter.
-

CAUTION

Do not use OxiCliq disposable sensors in a high humidity environment, or in the presence of fluids, which may contaminate sensor and electrical connections causing unreliable or intermittent measurements. Do not use disposable sensors when there is a known allergic reaction to the adhesive.

WARNING

Proper Sensor Fit: If a sensor is too loose, it might compromise the optical alignment or fall off. If it is too tight, for example because the application site is too large or becomes too large due to edema, excessive pressure may be applied. This can result in venous congestion distal from the application site, leading to interstitial edema, hypoxemia and tissue malnutrition. Skin irritations or lacerations may occur as a result of the sensor being attached to one location for too long. To avoid skin irritations and lacerations, periodically inspect the sensor application site and change the application site regularly.

Venous Pulsation: Do not apply sensor too tightly as this results in venous pulsation which may severely obstruct circulation and lead to inaccurate measurements.

Ambient Temperature: At elevated ambient temperatures be careful with measurement sites that are not well perfused, because this can cause severe burns after prolonged application. All listed sensors operate without risk of exceeding 41°C on the skin if the initial skin temperature does not exceed 35°C.

Extremities to Avoid: Avoid placing the sensor on extremities with an arterial catheter, an NBP cuff or an intravascular venous infusion line.

WARNING

- For fully conscious pediatric or adult patients, who have a normal function of perfusion and sensory perception at the measurement site:
 - To ensure skin quality and correct optical alignment of the sensor, inspect the application site when the measurement results are suspicious or when the patient complains about pressure at the application site, but at least every 24 hours. Correct the sensor alignment if necessary. Move the sensor to another site, if the skin quality changes.
 - For all other patients:
 - Inspect the application site every two to three hours to ensure skin quality and correct optical alignment. Correct the sensor alignment if necessary. If the skin quality changes, move the sensor to another site. Change the application site at least every four hours.
 - Injected dyes such as methylene blue, or intravascular dyshemoglobins such as methemoglobin and carboxyhemoglobin may lead to inaccurate measurements.
 - Inaccurate measurements may result when the application site for the sensor is deeply pigmented or deeply colored, for example, with nail polish, artificial nails, dye or pigmented cream.
 - Interference can be caused by:
 - High levels of ambient light (including IR warmers) or strobe lights or flashing lights (such as fire alarm lamps). (Hint: cover application site with opaque material.)
 - Another SpO₂ sensor in close proximity (e.g. when more than one SpO₂ measurement is performed on the same patient). Always cover both sensors with opaque material to reduce cross-interference.
 - Electromagnetic interference, especially at perfusion indicator values below 1.0 or signal quality indicator below medium.
 - Excessive patient movement and vibration.
-

WARNING

With pulse oximetry, sensor movement, ambient light (especially strobe lights or flashing lights) or electromagnetic interference can give unexpected intermittent readings when the sensor is not attached. Especially bandage-type sensor designs are sensitive to minimal sensor movement that might occur when the sensor is dangling.

CAUTION

If you measure SpO₂ on a limb that has an inflated noninvasive blood pressure cuff, a non-pulsatile SpO₂ INOP can occur. If the fetal monitor is configured to suppress this alarm there may be a delay of up to 60 seconds in indicating a critical status, such as sudden pulse loss or hypoxia.

WARNING

The fetal/maternal monitor is not a diagnostic ECG device. In particular, the display of fetal/maternal ECG is intended only for evaluating signal quality for fetal/maternal heart rate as derived from the ECG waveform.

When in doubt, it can be used to identify sources of compromised signal quality, such as noise or muscle artifacts. It can subsequently be used to verify the result of measures taken to resolve them (e.g., checking ECG cable connections).

The safety and effectiveness of the displayed fetal/maternal ECG waveform (i.e., P, QRS and T segments) for evaluation of fetal/maternal cardiac status during labor have not been evaluated.

WARNING

The fetal/maternal monitor is NOT intended for use during defibrillation, electro-surgery, or MRI. Remove all transducers, sensors, and accessories before performing electro-surgery, defibrillation, or MRI, otherwise harm to the patient or the user can result.

WARNING

- No alarm is possible when Toco MP or CL Toco⁺ MP transducer is the source of the pulse rate.
 - No QRS tone is audible when Toco MP or CL Toco⁺ MP transducer is the source of the pulse rate.
 - The Toco MP or CL Toco⁺ MP transducer signal is significantly less reliable if the patient is up and moving about, or is pushing during the second stage of labor.
-

WARNING

No alarm is possible when noninvasive blood pressure is the source of the pulse rate.

CAUTION

Only use Philips paper. Using paper other than Philips paper may result in the failure to recover traces.

CAUTION

Ensure that you admit each patient by name, including other patient identification information, and discharge the patient when you have finished monitoring, so that you can identify which trace period (entry in the patient list) refers to which patient.

WARNING

- Do not operate the monitor if it is wet. If you spill liquid on the monitor, contact your service personnel or Philips service engineer.
 - Do not perform underwater monitoring (for example, in a bath or shower) using wired transducers.
 - Place the monitor where there is no chance of contact with, or falling into water or other liquid.
 - Do not dry equipment using heating devices such as heaters, ovens (including microwave ovens), hair dryers and heating lamps.
 - Do not put equipment or accessories in autoclave (for sterilization).
-

CAUTION

Solutions: Do not mix disinfecting solutions (such as bleach and ammonia) as hazardous gasses may result.

Skin contact: To reduce the risk of skin irritations, do not allow a cleaning or disinfecting agent to leave residues on any of the equipment surfaces - wipe it off with a cloth dampened with water, after allowing the appropriate time for the agent to work, or before applying to a patient.

Hospital policy: Disinfect the product as determined by your hospital's policy, to avoid long term damage to the product.

Local requirements: Observe local laws governing the use of disinfecting agents.

Touch display: To clean and disinfect the touch-enabled display, disable the touch operation by switching off the monitor during the cleaning procedure, or by selecting and holding the Main Screen key until the padlock symbol appears on it, indicating that touch operation is disabled. Select and hold again to re-enable touch operation.

WARNING

Do not put device and accessories in autoclave (for sterilization).

WARNING

Schedule: Failure on the part of the responsible individual hospital or institution employing the use of this equipment to implement a satisfactory maintenance schedule may cause undue equipment failure and possible health hazards.

In case of problems: If you discover a problem with any of the equipment, contact your service personnel, Philips, or your authorized supplier.

Electric shock hazard: Do not open the monitor housing. Refer all servicing to qualified service personnel.

WARNING

To avoid contaminating or infecting personnel, the environment or other equipment, make sure you disinfect and decontaminate the monitor appropriately before disposing of it in accordance with your country's laws for equipment containing electrical and electronic parts. For disposal of parts and accessories such as thermometers, where not otherwise specified, follow local regulations regarding disposal of hospital waste.

WARNING

Reuse: Disposable accessories and supplies intended for single use, or single use only, and are indicated as such on their packaging. Never reuse disposable accessories and supplies, such as transducers, sensors, electrodes and so forth that are intended for single use, or single patient use only.

Approved accessories: Use only Philips-approved accessories.

Packaging: Do not use a sterilized accessory if its packaging is damaged.

Protection against electric shocks: The transducers and accessories listed in this chapter are not defibrillator proof.

Electro-Surgery, Defibrillation and MRI: The fetal/maternal monitors are not intended for use during defibrillation, electro-surgery, or MRI. Remove all transducers, sensors, and accessories before performing electro-surgery, defibrillation, or MRI, otherwise harm can result.

WARNING

Explosion Hazard: Do not use in the presence of flammable anesthetics, such as a flammable anesthetic mixture with air, oxygen, nitrous oxide, or in oxygen rich environment. Use of the devices in such an environment may present an explosion hazard.

WARNING

- Do not use additional AC mains extension cords or multiple portable socket-outlets. If a multiple portable socket-outlet is used, the resulting system must be compliant with IEC/EN 60601-1 A1: 2012.
 - Do not connect any devices that are not supported as part of a system.
 - Do not use a device in the patient vicinity if it does not comply with IEC/EN 60601-1 A1: 2012. The whole installation, including devices outside of the patient vicinity, must comply with IEC/EN 60601-1 A1: 2012. Any non-medical device, including a PC running an OB TraceVue system, placed and operated in the patient's vicinity must be powered via a separating transformer (compliant with IEC/EN 60601-1 A1: 2012) that ensures mechanical fixing of the power cords and covering of any unused power outlets.
 - Do not use USB devices with own power supplies unless an appropriate separation device is used (either between USB interface and device or between device and power).
-

CAUTION

- **FM20/FM30 only:** Although this is an electrical Class II device, it has a protective earth conductor which is needed for EMC purposes.
 - Always use the supplied power cord with the three-prong plug to connect the monitor to AC mains. Never adapt the three-prong plug from the power supply to fit a two-slot outlet.
-

CAUTION

The use of accessories, transducers and cables other than those specified may result in increased electromagnetic emissions or decreased electromagnetic immunity of the device.

WARNING

Do not use cordless/mobile phones or any other portable RF communication system within the patient vicinity, or within a 1.0 m radius of any part of the fetal monitoring system.

WARNING

For paced patients: The radiated SRR power of the CL SpO₂ and NBP Maternal Cableless Measurement Devices, the CL Transmitter, and other sources of radio-frequency energy, when used in very close proximity of a pacemaker, might be sufficient to interfere with pacemaker performance. Due to shielding effects of the body, internal pacemakers are somewhat less vulnerable than external pacemakers. However, caution should be exercised when monitoring paced patients.

In order to minimize the possibility of interference, avoid positioning and wearing the Cableless Measurement Devices in very close proximity to a pacemaker. Consult the pacemaker manufacturer for information on the RF susceptibility of their products

CAUTION

Fetal parameters, especially ultrasound and ECG, are sensitive measurements involving small signals, and the monitoring equipment contains very sensitive high gain front-end amplifiers. Immunity levels for radiated RF electromagnetic fields and conducted disturbances induced by RF fields are subject to technological limitations. To ensure that external electromagnetic fields do not cause erroneous measurements, it is recommended to avoid the use of electrically radiating equipment in close proximity to these measurements.

CAUTION

The device should not be used adjacent to, or stacked with, other equipment unless otherwise specified.

WARNING

- **Leakage currents:** If several items of equipment used to monitor a patient are interconnected, the resulting leakage current may exceed allowable limits.
 - **ECG electrodes:** NEVER allow ECG electrodes to contact other electrical conductive parts, including earth.
-

WARNING

The fetal/maternal monitors are NOT intended for use during defibrillation, electro-surgery, or MRI. Remove all transducers, sensors, and accessories before performing electro-surgery, defibrillation, or MRI, otherwise harm to the patient or the user can result.







WARNING

The fetal/maternal monitors are not intended for use for ECG measurements on patients connected to external electrical stimulator or with cardiac pacemakers.

Device Classification and Names within this Manual

New cableless measurements for fetal and maternal monitoring supplement the Avalon fetal monitor family. The **Avalon CL** solution consists of the **Avalon CL base station**, the **Avalon CL transducers**, and the **IntelliVue CL Pods**.

The **IntelliVue CL Pods** are only used for maternal measurements within the **Avalon CL solution**. To differentiate between the **CL transducers** and the **CL Pods**, the **CL Pod** are also referred to as **maternal measurement Pods**. The following table provides an overview of all the devices within the solution.

Avalon FM20/FM30 and FM40/ FM50 Wired Transducers	Avalon CL Base Station Avalon CL Transducers	IntelliVue CL Pods
<p>FM20/ FM30</p>  <p>M2702A and M2703A</p>		
<p>FM40/ FM50</p>  <p>M2704A and M2705A</p>	<p>Avalon CL Base Station</p>  <p>866074</p>	
<p>US transducer wired</p>  <p>M2736A</p>	<p>Avalon CL Toco+ MP Transducer cableless</p>  <p>866075</p>	<p>IntelliVue CL NBP Pod cableless</p>  <p>865216</p>

**Avalon FM20/FM30 and FM40/
FM50
Wired Transducers**

**Toco/ Toco
MP
transducer**

wired



M2734A and M2734B

**Toco+
transducer
with ECG/
IUP
capability**

wired



M2735A

**Patient
Module for
ECG/IUP**



M2738A

**Avalon CL Base Station
Avalon CL Transducers**

**Avalon CL
US
Transducer**

cableless



866076

**Avalon CL
ECG/IUP
Transducer**

cableless



866077

IntelliVue CL Pods

**IntelliVue
CL SpO2
Pod**

cableless



865215

What's New

This section lists the most important new features and improvements to the fetal monitors and their user interface introduced with Release J.3.

You may not have all of these features, depending on the fetal monitor configuration purchased.

What's New in Release J.3

Avalon CL Transducer System

The Avalon CL system provides cableless monitoring with the Avalon FM20/FM30 and FM40/FM50 with the same functionality and performance as the wired measurement devices (e.g. twin and triplets monitoring). The Avalon CL has a straight-forward handling and operating concept. The cableless transducers are assigned by simply docking them at the base station, no further configurations are necessary.

The Avalon CL system includes the following features:

- Cableless monitoring of twins and triplets
- Cableless maternal measurement Pods CL SpO₂ and CL NBP
- Maternal pulse from a CL Toco⁺ MP (Smart Pulse) transducer
- A cableless ECG/IUP transducer measuring IUP and fetal/maternal ECG
- Watertight cableless transducers that can be used to monitor in water
- Patient call that pages an ambulating mother with an audible signal emitted by the worn CL transducers.
- Out of range audible signal emitted by the worn CL transducers to inform an ambulating mother that she is approaching the limit of the active signal area of reach.
- A transducer finder LED on all cableless transducers to help identify the assigned transducer

Support For Use of Maternal Cableless Measurement Devices

The IntelliVue measurement pods are patient-worn, battery powered measurement devices for SpO₂ and NBP. The devices provide measurement values on the built-in display and communicate them to the fetal monitor using the wireless short range radio (SRR) interface of the Avalon CL base station.

Maternal Temperature Measurement

To measure maternal temperature, the new optional Tympanic Temperature device (866149) is available for the Avalon fetal monitors. The measurement data is:

- fully documented at the local recorder and the OB system,
- can be displayed as a numeric on the screen,
- and is printed out on the recorder trace strip.

Manually Entered Maternal Temperature Measurements

Manually measured temperatures can be entered manually into the fetal monitor and stored in the database. They are also displayed as a numeric on screen, and are printed out on the recorder trace if required.

New Design for the User Interface

The user interface for the fetal monitors has been redesigned to bring the presented information into the foreground, letting the structural elements such as keys and frames retreat into the background. Additionally special regard was given to making the "look and feel" similar to that of standard software products.

New SmartKeys

- The **Start ECG** SmartKey is renamed to **Record ECG** .
- With the **Call Patient** SmartKey you can now page patients who are ambulating wearing cableless transducers.
- With the **Tele Info** SmartKey you can call up the **Tele Info** window on the fetal monitor display. In the **Tele Info** window you can control and view the status of the cableless transducers from the connected Avalon CL system.
- With the **Enter Temp** SmartKey a pop-up window opens showing a numeric pad for entering externally measured maternal temperature values.

Coincidence INOP Tone

When the cross-channel verification detects that the signal of the maternal heart rate coincides with the fetal heart rate, the **Coincidence** INOP is issued with a tone at the fetal monitor. The **Coincidence** INOP tone has a configurable delay.

Increased Internal Back-up Memory

The internal back-up memory is now able to store data at least 3 hours.

USB Interface

An optional USB interface allows the use of USB printers, bar code readers, and input devices such as a keyboard or mouse.

Flexible Nurse Call Interface

An optional Flexible Nurse Call interface allows now the connection of a nurse call device to the fetal monitors.

DHCP Support

The DHCP support offers an alternative to BOOTP. DHCP (dynamic host configuration protocol) enables the fetal monitors to request an IP address (internet protocol address) from the connected network (OB TraceVue/IntelliSpace Perinatal) automatically.

Data Export Support

You can now export measurement values from the monitor to other devices via the LAN or RS232 interface.

NBP Configurable Measurement Sequence

Up to four measurement cycles can be set up which will run consecutively. For each cycle you can set the number of measurements and the interval between them. By setting the last cycle to run continuously you can have regular measurements continue after the sequence has run.

Alarms Enhancements

- In addition to the standard blue INOPs, some INOPs can now be configured as red or yellow INOPs to provide a severity indication (**ECG Leads Off**, **NBP Cuff Overpress**, **Cuff Not Deflated**). A small number of INOPs are always yellow or red to indicate a severity corresponding to red and yellow alarms.
- For the invasive pressure measurement, the extreme pressure alarms Extreme High and Extreme Low can be made available for your monitor in Configuration Mode and are additional to the standard High and Low limit alarms.
- The **Review Alarms** window now shows when the monitor was switched on (after being switched off for longer than 1 minute) and any changes made to the Standby and paired status.
- In order to improve alarming on asystole under certain conditions, you can set **Asystole Detect** in Configuration Mode to **Enhanced**. In enhanced mode an asystole alarm will be suppressed for up to five seconds if a valid beat-to-beat Pulse is detected from a Pressure.

Alarm Reminder

In Configuration Mode you can set now an **Alarm Reminder**. The **Alarm Reminder** emits an audible reminder of alarm conditions that remain active after the alarm is acknowledged. This reminder may take the form of a repetition of the alarm tone for a limited time, or an unlimited repetition of the alarm tone (this is the same as a new alarm). The interval between silencing the alarm and sounding the reminder tone can be set to one, two, or three minutes.

Auto Free

In Configuration Mode you can now set an **Auto Free** setting which discharges a patient automatically when the fetal monitor is powered off, or is in standby mode for a set time. Only the demographic patient data is deleted, the trace data is not affected.

What's New in Release G.0

Battery Option and Patient Transport Improvements for the Avalon FM20 and FM30

- For the FM20/FM30, you can now switch between battery-powered and mains-powered operation without interrupting monitoring. The monitor is connected to the AC mains power with the external power supply.
- Traces of multiple patients can be recorded during transport and can automatically be uploaded to OB TraceVue/IntelliSpace Perinatal (Rev. G or higher) when reconnected.

- A bed hanger is available as a mounting option, especially for patient transport purposes within healthcare facilities.

Maternal Pulse from Toco MP Transducer

- The new Toco MP transducer can measure the maternal pulse rate, in addition to the other available sources; MECG, SpO₂, and NBP. With sensors built-into the bottom cover of the Toco MP transducer, and using a measurement technology similar to SpO₂, it automatically gives you an additional maternal pulse source for Cross-Channel Verification (CCV).
- The new Toco MP transducer is standard for the whole product family, Avalon FM20 to FM50.

Non Stress Test (NST) Analysis as Clinical Decision Support (CDS) Application

- The optional **NST Report** Trace Interpretation feature allows you to automatically interpret FHR traces and to generate a printed NST report, equivalent to the NST report functionality in OB TraceVue Rev. F and based on the NICHD guidelines from 1997.
- The application can be adapted to domestic or hospital guidelines. A report of the NST analysis can be printed automatically or on demand.
- This software option is available for the whole product family, Avalon FM20 to FM50.

FHR Sound Source

A newly-connected fetal measurement (FHR or DFHR) automatically becomes the FHR sound source, without the need for manual intervention. If you prefer the previous behavior, this can be changed in Configuration mode (**Audio Select** setting).

FHR Numeric Display

The FHR numeric display shows changes in fetal heart rate faster and is updated more frequently.

Improved FHR and DFHR Label Concept

A DFHR numeric (fetal heart rate from DECG measurement) now also has a unique number (**DFHR1**, **DFHR2**, or **DFHR3**) to allow easier identification of the related trace in OB TraceVue/IntelliSpace Perinatal.

New SmartKeys

You can control fetal heart sound volume directly with the configurable SmartKeys (**FHR Vol. Up** and **FHR Vol.Down**) on the monitor's screen.

Now the **Toco Baseline** SmartKey automatically turns into the **Zero IUP** SmartKey when switching to IUP measurement.

Trace Separation On/Off Operation

Now you can switch the **Trace Separation** on or off in every FHR setup menu.

NBP

- To prevent the NBP measurement from being switched on or off accidentally, this operation is now only available in Configuration Mode.

- You can enter the NBP setup menu by touching the NBP numeric (on the screen) before the first measurement is taken.
- Algorithm enhancements (more tolerant of movement artifacts).

Alarms

- There is an individual Alarm Tone Volume setting for Yellow, Red and Cyan alarms.
- The new alarm settings available in Configuration Mode are:
 - **AutoIncrease Vol** and **IncreaseVolDelay**
 - **AlarmsOffAtStart**
 - **ConfirmAlarmsOff**

For a detailed description of the settings see the Configuration Guide.

CCV INOP

There is a new INOP (technical alarm) **Coincidence** after 1 minute of persistent coincidence warning.

New Demographic Fields

- **Date of Birth** and **Gestational Age** can be entered for complete documentation in an **NST Report** (supporting the optional NST Trace Interpretation software feature).
- A **Middle Name** field is now available in the admission form - whether it appears is configurable.
- Two additional ID fields **Lifetime ID** and **Encounter ID** can also be configured to appear, and their names can be customized to fit hospital requirements. The default label of the **Lifetime ID** is **MRN**.

Recorder

- It can now be configured whether the recorder speed can only be adjusted in Configuration mode, or also in Monitoring mode.
- When recorder speed is changed, a new trace header is printed. Date and/or time changes are annotated in real time recording.

Stored Data Recording

Printing progress is shown while printing a trace from the internal back-up memory.

FHR Sound Volume

Volume steps are optimized in the low range to allow finer FHR sound volume adjustments.

Avalon CTS

An antenna symbol is displayed next to MEKG waves if an Avalon CTS is used.

Basic Operation

This chapter gives you an overview of the monitor and its functions. It tells you how to perform tasks that are common to all measurements (such as entering data, switching a measurement on, changing some monitor settings, and setting up the recorder). The alarms section gives an overview of alarms. The remaining sections tell you how to perform individual measurements, and how to care for and maintain the equipment.



Supported Measurements

The following **Fetal** measurements are supported:

Fetal Monitor or Model	Fetal Heart Rate (FHR) via US (including Twins)	Triple FHR via US	Toco	FHR via Direct ECG (DECG)	Intrauterine Pressure (IUP)
FM20	Standard	Optional	Standard	-	-
	Optional	Optional	Standard	-	-
FM30	Standard	Optional	Standard	Standard	Standard
	Optional	Optional	Standard	Standard	Standard
FM40	Standard	Optional	Standard	-	-
	Optional	Optional	Standard	-	-
FM50	Standard	Optional	Standard	Standard	Standard
	Optional	Optional	Standard	Standard	Standard

3 Basic Operation

The following **Maternal** measurements are supported:

Fetal Monitor or Model	Maternal Heart Rate (MHR) via Maternal ECG Electrodes	Maternal ECG (MECG) Wave	Maternal Pulse from Toco	Non-invasive Blood Pressure with Pulse Rate	Pulse Oximetry (Maternal SpO ₂) with Pulse Rate	Maternal Temperature
FM20	Standard	-	Standard	Optional (built-in)	-	Optional (Interface cable and manual entry)
	Standard	-	Standard	Optional (cableless)	Optional (cableless)	-
FM30	Standard	Standard	Standard	Optional (built-in)	Optional (built-in)	Optional (Interface cable and manual entry)
	Standard	Standard	Standard	Optional (cableless)	Optional (cableless)	-
FM40	Standard	-	Standard	Standard (built-in)	Standard (built-in)	Optional (Interface cable and manual entry)
	Standard	-	Standard	Optional (cableless)	Optional (cableless)	-
FM50	Standard	Standard	Standard	Standard (built-in)	Standard (built-in)	Optional (Interface cable and manual entry)
	Standard	Standard	Standard	Optional (cableless)	Optional (cableless)	-

Avalon FM20 and FM30

This section outlines the capabilities of your monitor.



Avalon FM20

The Avalon FM20 fetal/maternal monitor provides a solution for external fetal monitoring applications, and optional noninvasive maternal vital signs.

You can monitor fetal heart rates (FHRs) externally using ultrasound, uterine activity and maternal pulse using an external Toco transducer, and the maternal heart rate (MHR) with maternal ECG electrodes, and optionally, noninvasive blood pressure and maternal oxygen saturation (SpO₂).

Measurements are displayed on a 6.5-inch color display as numerics. The display is a touchscreen, and you operate the monitor using this touchscreen interface. The integrated recorder documents fetal and maternal measurements as well as the user-defined annotations.

You can connect the monitor to an OB TraceVue/IntelliSpace Perinatal system via the RS232 connection, or over a LAN connection (with OB TraceVue Revision E.00.00 and later, or IntelliSpace Perinatal Revision H.0 and later).

Avalon FM30

The Avalon FM30 fetal/maternal monitor offers a solution for both external and internal fetal monitoring applications, and optional noninvasive maternal vital signs.

The Avalon FM30 shares all the features and capabilities of the Avalon FM20. In addition, you can monitor one FHR internally with a direct fetal electrocardiogram (DECG), uterine activity internally using an intra-uterine pressure (IUP) catheter together with a Toco⁺ transducer or patient module.

The Avalon FM30 carries the **IP** label, indicating that it is capable of intrapartum monitoring.



FM20/30 with Battery Option #E25 Only The battery option for the FM20/30 provides support for the in-transport monitoring of all measurements when disconnected from a power supply. Existing data storage is automatically uploaded to OB TraceVue or IntelliSpace Perinatal after reconnecting it to the system. Trace printing during transport is also possible.

Avalon FM40 and FM50

This section outlines the capabilities of your monitor.



Avalon FM40

The Avalon FM40 fetal/maternal monitor provides a solution for external fetal monitoring applications, and noninvasive maternal vital signs.

You can monitor fetal heart rates (FHRs) externally using ultrasound, uterine activity using an external Toco transducer, and the maternal heart rate (MHR) via maternal ECG electrodes, and non-invasive blood pressure and maternal oxygen saturation (SpO₂).

Measurements are displayed on a 6.5-inch color display as numerics. The display is a touchscreen, and you operate the monitor using this touchscreen interface. The integrated recorder documents fetal and maternal measurements as well as the user-defined annotations.

You can connect the monitor to an OB TraceVue/IntelliSpace Perinatal system with the RS232 connection, or over a LAN connection (with OB TraceVue Revision E.00.00 and later, or IntelliSpace Perinatal Revision H.0 and later).

Avalon FM50

The Avalon FM50 fetal/maternal monitor offers a solution for both external and internal fetal monitoring applications, and noninvasive maternal vital signs.

The Avalon FM50 shares all the features and capabilities of the Avalon FM40. In addition, you can monitor one FHR internally with a direct fetal electrocardiogram (DECG), and uterine activity internally using an intra-uterine pressure (IUP) catheter together with a Toco⁺ transducer or patient module.

The Avalon FM50 carries the **IP** label, indicating that it is capable of intrapartum monitoring.



Avalon CL Transducer System




The Avalon CL Fetal Transducer System lets you monitor the patient continuously with cableless transducers during the antepartum period, labor, and delivery. You can monitor the fetal heart rate (FHR) using noninvasive CL Ultrasound transducers, or invasively using the CL ECG/IUP transducer with the direct electrocardiogram (DECG). The uterine activity can be monitored using an external CL Toco⁺ MP transducer. The fetal and maternal parameters are measured and transmitted via radio frequency from the CL transducers to the CL base station, eliminating the need for patient cables. With the Avalon CL Transducer System you can monitor a single fetus, twins and triplets.

The Avalon fetal monitor (FM20-FM50) connected to the CL base station displays and records the parameters. All the CL transducers are watertight. You can continuously monitor patients in a bath or shower using the CL Toco⁺MP and the CL Ultrasound transducers.

The system should only be used by, or under the direct supervision of, a licensed physician or other health care practitioner, who is trained in the use of FHR monitors and in the interpretation of FHR traces.

Getting to Know Your Avalon FM20/FM30

Overview



Right Side



- 1 On/Off Switch
- 2 Power Connector

with
Battery
Option



- 1 On/Standby Switch
- 2 MSL Connector

Left Side



- 1 SpO₂ Socket (optional)
- 2 Noninvasive Blood Pressure Socket (optional)
- 3 Fetal Sensor Sockets

Each of the fetal sensor sockets accepts any fetal transducer, one Avalon CL or one Avalon CTS Cableless Fetal Transducer System base station, or an event marker.

Bottom



There are five optional interfaces available for the Avalon FM20/30 monitor:

- LAN/RS232 system interface
- Dual PS/2 interface
- MIB/RS232 interface
- Flexible Nurse Call
- USB port

You can use two of the five optional interfaces at the same time.



- 1 LAN/RS232 system interface
- 2 Dual PS/2 system interface

Optional Interfaces	Description
	Flexible nurse call interface card
	USB ports


Rear



- 1 Display Release
- 2 Carrying Handle
- 3 Built-in Stand

Getting to Know Your Avalon FM40/FM50


Front



- 1 Touchscreen color display
- 2 Transparent paper guide with tear-off edge
- 3 Paper eject button
- 4 Power LED
- 5 On/Standby switch
- 6 Recorder paper table
- 7 Fetal sensor sockets
- 8 Noninvasive blood pressure socket
- 9 SpO₂ socket

Connect any fetal sensor or patient module at the fetal sensor sockets, including an Avalon CL or an Avalon CTS via interface cable (with red connector).



Rear




- 1 Reserved for future use: protective earth intended for use in system installations.
- 2 Equipotential grounding point
- 3 Power cord connector
- 4 Loudspeaker
- 5 Slot 01 for optional LAN/RS232 system interface (for connection to an obstetrical information and surveillance system)
- 6 Slot 03 reserved for future use
- 7 Video output (VGA)
- 8 Telemetry interface
- 9 Slot 02 for optional interfaces: *Either* dual PS/2 system interface (A) for mouse and keyboard connection) *Or* MIB interface (B) for external touch screen connection.

Two Avalon CL, or one Avalon CTS can be also connected to the Telemetry interface sockets using the interface cable (with black connector).

Additional Optional Interfaces


Optional Interfaces	Description
	Flexible nurse call interface card
	USB ports

Transducers



- 1 Transducer finder LED - lights up on the transducer providing the measurement source
- 2 "MP" for M2734B "Toco MP" transducers (additionally capable of providing the maternal pulse measurement)
- 3 Belt Button


Toco (M2734A) and Toco MP Transducer (M2734B)



- 1 Cable - connects to any of the four fetal sensor sockets on the monitor


The M2736AA US transducer is identical to the M2736A US transducer, including all specifications

Ultrasound Transducer (M2736A)



- 1 Connector - for connecting ECG/IUP adapter cables (M2735A Toco⁺ transducer only)

Toco⁺ Transducer with ECG/IUP capability (M2735A)



- 1 Butterfly belt clip (shown fitted; for use with belts without button holes)
- 2 Close-up of MECCG adapter cable connected to Toco⁺ transducer
- 3 Close-up of active finder LED



- 1 Connector - for connecting ECG/IUP adapter cables (same as for Toco⁺ transducer)
- 2 Cable - connects to any of the four fetal sensor sockets on the monitor

Patient Module for ECG/IUP (M2738A)

Getting to Know Your Avalon CL

Front



- 1 Transducer docking slots
- 2 On/Standby button with power LED
- 3 Device label
- 4 LED for optional cableless pods
- 5 Docking slots for cableless pods

WARNING

- To avoid magnetic interference affecting the mode of the pacemaker, ensure that the Avalon CL base station does not come into close contact with implanted pacemakers.
 - This equipment generates, uses, and radiates radio-frequency energy, and if it is not installed and used in accordance with its accompanying documentation, may cause interference to radio communications. Operation of this equipment in a residential area may cause interference, in which case the users must take whatever measures may be required to correct the interference.
-


Docking Slot for Cableless Transducers

The Avalon CL base station has three docking slots to mechanically hold the CL transducers. The flexible electronic contacts identify the transducer type after docking. The transducers are charged while docked. The base station has a built-in radio interface with an integrated antenna to communicate with the transducers for assignment, configuration, and for transducer firmware updates. The transducers can be controlled with the base station and the fetal monitor.



On/Standby Button

Pressing the On/Standby button switches the Avalon CL base station between the two modes On and Standby.



1 On/Standby button

If you switch the base station to **On**, the LED button lights up green. The base station is ready for use. The CL transducers are charging. When the transducers are picked up from their docking slot, the base station sets up the radio communication to the CL transducer automatically. Wired transducers connected to the fetal monitor are disabled, and the antenna symbol for the CL transducer is displayed on the screen of the fetal monitor.

If you switch the base station to **Standby**, the LED button turns off. The base station is now in Standby mode. The docked CL transducers continue charging. No radio communication occurs, any existing radio communication is stopped.

If the LED of the On/Standby button turns red, it indicates a technical problem has occurred that needs your attention. Check your monitor for a possible related INOP message.

Docking Slots for Cableless Measurement Pods

The Avalon CL base station has one docking slot to mechanically hold an IntelliVue CL pod. The electronic contacts identify the CL pod type after docking. The CL pod is charged while docked. The LED under the docking slot indicates the battery status of the CL pod. The base station has a built-in short range radio interface with an integrated antenna to communicate with the CL pod for assignment and configuration. The CL pod can be controlled with the base station and the fetal monitor.



1 Docking slot for CL pods

WARNING

Short range radio connections are subject to interruption due to interference from other radio sources in the vicinity, including microwaves, bluetooth devices, WLAN devices (802.11b,g,n) and cordless phones. Depending on the strength and duration of the interference, the interruption may occur for an extended period. A loss of connection, due to moving out-of-range, interference, or for other reasons, is indicated with a **No Host Monitoring** INOP (here the host is the fetal monitor) on the NBP or SpO₂ Pods, or a **cl NBP Disconnect** or **cl SpO₂ Disconnect** INOP at the fetal monitor. Correct channel configuration is important, see the Configuration Guide for details.

Rear and Bottom



1 Name plate



1 Manufacturer label

2 Cable reel


Cableless Transducers

The cableless Avalon CL transducers have a built-in radio interface with an on-board, integrated antenna. The fetal monitor connected to the Avalon CL base station can control the transducer using radio communication. The transducer transmits measured values, alarms, and status information to the fetal monitor.

The cableless transducers are assigned by simply docking them at the Avalon CL base station. If they are picked up from their docking slot at the base station, they are automatically activated.

CAUTION

To ensure the correct assignment of CL transducers to the fetal monitor in use, always dock them at the base station connected to that fetal monitor before use. This is particularly important when you add an additional CL transducer that was previously docked at another base station.




- 1 Transducer finder LED - lights up on the transducer providing the measurement source.
- 2 CL Toco⁺ MP transducers (additionally capable of providing the maternal pulse measurement)
- 3 Belt button



CL Ultrasound transducer

CL ECG/IUP transducer



Radio Range of CL Transducers

The CL transducers have a operating range around the base station of at least 70 m in the line of sight. Obstructions as walls, metal doors, elevators and other environment structures can lead to signal loss. The Tele symbol indicator and the Tele Info window on the monitor's display (see “Screen Elements” on page 56 and “Tele Info Window” on page 92) provide information on the status of the signal strength.

When a patient is approaching the end of range, the US sound is replaced by an artificial QRS sound (like DECG), and the DECG and MECG waves are no longer displayed on the monitor.


Connector Cap for the CL Toco+ MP Transducer

The CL Toco⁺ MP transducer is delivered with a connector cap covering the MECG/DECG/IUP connector. The connector cap is designed to cover the connector, not to protect it from water. (The connector itself is water-proof and may be immersed in water). You can reorder lost connector caps in a kit of 10 (989803184841).



Cableless Transducer LED Indication

The cableless transducers have a multi-color LED that indicates the status of the transducer with specific colors. This LED remains visible when the transducer is correctly attached to the transducer belt (Philips standard belt).



LED Status	Meaning
White	The LED lights up to identify the US transducer among other US transducers, and to easily verify the correct transducer assignment (transducer finder). The transducer finder LED is controlled by the fetal monitor. Click on the numerics to identify the corresponding US transducer. The LED also lights up when the mother is paged with the Call Patient SmartKey.
White one short blink	The LED shortly lights up to indicate that the transducer successfully opened a radio communication with the base station and that it is ready to use.
Green	The LED lights up green when the transducer is fully charged and docked at the base station.
Yellow	The LED lights up yellow when the transducer is charging and docked at the base station.
Red	The LED lights up red when the transducer is out of battery and has to be recharged.
Cyan	The LED lights up cyan to indicate a technical problem that needs your attention. Check your fetal monitor for a related INOP.

CL Transducer Battery

Battery replacement is recommended after 500 charge/discharge cycles *2) or if the battery is older than 4 years whatever is reached first. If the battery of a cableless transducer has aged and an exchange of the battery is highly recommended, a prompt message is displayed at the fetal monitor for ca. 60 seconds. The prompt is repeated whenever the **Tele Info** window is opened until the battery is replaced.

Depending on the transducer type the following messages are displayed:

cl US battery has aged. Replacement strongly recommended

cl Toco battery has aged. Replacement strongly recommended

cl ECG/IUP batt has aged. Replacement strongly recommended


If you are getting this message contact your service personnel so they can replace the CL transducer's battery.

The date of manufacture and number of charge/discharge cycles can be inquired on the connected fetal monitor.

For battery cycle and battery life specification of the CL pods (NBP and SpO₂) please refer to the Service Guide of the Cableless Measurements.

Charge Cycle

The CL transducer batteries are designed for frequent recharging. A complete charging cycle is only reached and counted, when all recharging periods equal a 100% charge (900 mAh equal 8 hours continued operation).



Audio Signal CL Transducers

The Avalon CL transducers have two audio signals:

- If a tone sequence is played, the transducer was triggered by the fetal monitor to page the patient wearing it, or to locate a not docked transducer.
- An INOP tone indicates that the transducer has detected a technical problem (INOP). Check your fetal monitor for a related INOP message.

CL Pods

The two CL Pods provide measurement values for SpO₂ and NBP on the built-in display, and communicate them to other system components using a wireless short range radio (SRR) interface. They are controlled with SRR from the Avalon CL base station and the connected fetal monitor.

The maternal measurement pods are easily assigned by docking them at the Avalon CL base station.

CL NBP Pod



CL SpO₂ Pod



The SpO₂ Pod and the NBP Pod have an LCD display and three keys for basic operation e.g. to assign the device to a patient:




- 1 Integrated LCD display
- 2 Hard keys
- 3 Measurement identifier

For further operational details read the IntelliVue Cableless Measurements Instruction for Use.

Battery Status LED for CL Pods

The CL Pods do not have their own battery status LED. On the Avalon CL base station the battery status LED for the CL Pods is located directly under the docking slot.



1 Battery status LED for CL Pods

The battery status LED shows five different states:

Status	Meaning
Green	The docked CL Pod is fully charged.
Yellow	The docked CL Pod is charging.
Yellow blinking	The communication is established with the docked CL Pod.
Cyan	Indicates that the docked CL Pod or the charging slot has a technical problem that needs your attention. Check your fetal monitor for a related INOP message.
Off	The battery status LED is off, when no CL Pod is docked.

Audio Signal

The Avalon CL base station has audio signals to let you know when a transducer's battery has to be recharged, or if the base station cannot set up radio communication with a CL transducer or a CL Pod. Only active CL transducers and CL Pods emit an audio signal.

- If the base station issues a descending tone sequence, the attempt to set up radio communication to a CL transducer or a CL Pod has failed.
- If an interval tone is played, the battery of one or more the CL transducers or CL Pods need to be recharged.

Operating and Navigating

Your monitor has a touchscreen. Everything you need to operate the monitor, except the on and off switch, is contained on its screen. Most screen elements are interactive. Screen elements include measurement numerics, screen keys, information fields, status indicators, alarms fields, and menus.

Operator Position The typical operator's position is in front of the monitor.

3 Basic Operation

FM40/50 If an optional external touch display is connected to the monitor, you can operate the monitor using the external touch display.



- 1 Monitor Information line
- 2 Measurement area
- 3 Key area


Screen Elements

Monitor Information Line














- 1 LAN connection status indicator only. RS232 system connection is not indicated. The train icon indicates if the fetal monitor is connected to OB TraceVue/IntelliSpace Perinatal, via a LAN cable or not.
- 2 Patient identification
- 3 Date and time
- 4 Bed label (when connected to a Philips OB TraceVue/IntelliSpace Perinatal system)
- 5 Fetal heart sound volume adjust/indicator
- 6 Alarm volume adjust/indicator
- 7 INOP and alarm status area - shows active alarm messages


Measurement Area




- 1 Antenna symbol (indicates a cableless measurement from a connected Avalon CL or Avalon CTS system).
- 2 Configurable alarm limits
- 3 NST test
- 4 Audio source symbol
- 5 Measurement numeric
- 6 Alarms off symbol
- 7 Fetal trace recorder - status indicator
- 8 Avalon CL or Avalon CTS system - status indicator
- 9 Battery status indicator
- 10 Status line - shows status and prompt messages
- 11 Signal quality indicator: good, acceptable, poor
- 12 Fetal heart rate measurement label
- 13 Measurement unit (configurable)
- 14 Cross-Channel Verification symbol (see “Cross-Channel Verification (CCV)” on page 149)
- 15 NST timer, if configured (default is Off)

Screen Details

Icon	Description
	The antenna symbol indicates a cableless measurement (Avalon CL or Avalon CTS).
	Indicates a short range radio measurement (IntelliVue CL Pods).
	Signal quality indicator: 1 Good 2 Acceptable 3 Poor
	Fetal trace recorder - status indicator Fetal recorder is on
	Fetal recorder is off (when Paper Save Mode is off)
	Fetal recorder is off (when Paper Save Mode is on)
	There is a user-solvable recorder error (paper out, paper jam, wrong paper scale set)
	Fetal recorder is defective: call service
	Avalon CL or Avalon CTS system - status indicator: An interface cable is connected to the monitor, a base station is connected, powered on, and at least one cableless transducer has been taken from the base station and is active. Any connected wired transducers are disabled.
	An interface cable is connected to the monitor, the base station is connected, powered on, and cableless transducers are ready to use, but no cableless transducers are currently active (all are still docked at the base station).
	An interface cable is connected to the monitor, but the base station is not connected to the interface cable, or it is disconnected from AC mains, or is in Standby mode.

Icon	Description
	The CL transducer is moved away from the base station approaching the limit of the area of reach. The US sound is replaced by an artificial QRS sound (like DECG). The DECG and MEGC waves are no longer displayed.

Key Area





- 1 SmartKeys - these can vary according to your monitor's configuration
- 2 **Main Screen** - which closes all open menus and windows and returns to main screen
- 3 Scroll to display more SmartKeys
- 4 **Silence** - key which acknowledges all active alarms by switching off audible alarm indicators

Keys

The monitor has three different types of keys.














Permanent Keys

















A permanent key is a graphical key that remains permanently on the screen, giving you fast access to functions.

Key	Name	Function
	Silence	Acknowledges all active alarms by switching off audible alarm indicators.
	Main Screen	Closes all open menus and windows and returns to the main screen.

SmartKeys

SmartKeys are configurable graphical keys, located at the bottom of the main screen. They give you fast access to functions. The selection of SmartKeys available on your monitor depends on your monitor configuration, and on the options purchased.

SmartKey	Name	Function
	Main Setup	Enters main setup menu
	Pause Alarms	Pauses alarm indicators. Pause duration depends on monitor configuration. If the pause duration is infinite, this key is labeled Alarms Off Select again to immediately re-enable alarm indicators
	Start Recordng	Turns the trace recorder on
	Record ECG Wave	Starts printing the MECCG, DECG or both waves, when both are available
	Patient Demogr.	Enters the patient identification menu to admit/discharge
	QuickAdmit	Quick admits the patient for monitoring
	Toco Baseline	Resets Toco baseline
	Zero IUP	Zeros IUP measurement
	Start/ Stop	<ul style="list-style-type: none"> Starts/stops manual noninvasive blood pressure measurement Starts auto series Stops current automatic measurement within series
	Start NBP	<ul style="list-style-type: none"> Starts manual noninvasive blood pressure measurement Starts auto series
	Repeat Time	Sets the time interval between two noninvasive blood pressure measurements
	Sound Vol. Up	Increases the fetal heart rate volume
	Monitor Standby	Enters the Standby mode, suspends monitoring. All numerics and waves disappear from the display. All settings and patient data information are retained

SmartKey	Name	Function
	Recorder Start/ Stop	Turns the trace recorder on or off
	Paper Advance	Advances the paper automatically to the next fold
	Stop Recording	Turns the trace recorder off
	Stored Data Rec	Prints trace data from the monitor's memory
	Enter Notes	Enters notes
	Timer	Enters NST timer window
	Set Marker	Marks an event
	Stop All	Stops all noninvasive blood pressure measurements
	Stop NBP	Stops manual noninvasive blood pressure measurement Stops current automatic measurement within series
	NBP Modes	Access NBP mode selection and setup, with direct start/stop function
	Defaults	Loads User Default
	Sound Vol. Down	Decreases the Fetal Heart Rate volume
	NST Report	Initiates an NST trace interpretation and obtains a Non-stress test (NST) report
	Call Patient	Pages the patient. Only enabled if at least one CL transducer is currently active
	Tele Info	Calls up the Tele Info window
	Enter Temp	Allows the manual input of the patient's temperature

Pop-Up Keys

Pop-up keys are context-sensitive graphical keys that appear automatically on the monitor screen when required. E.g. the **Confirm** pop-up key appears when you need to confirm a change.

Using the Touchscreen

Select screen elements by pressing them directly on the monitor's screen.

Disabling Touchscreen Operation

- 1 To temporarily disable the touchscreen operation of the monitor, press and hold the **Main Screen** permanent key for about three seconds. A red padlock will blink on the **Main Screen** permanent key.
- 2 Press and hold the **Main Screen** permanent key again for about three seconds to re-enable the touchscreen operation.



Operating Modes

When you switch on the monitor, it starts up in Monitoring Mode. To change to a different mode:

- 1 Select the **Main Setup** menu.
- 2 Select **Operating Modes** and select a mode.

Your monitor has four operating modes. Some are passcode protected.

Mode	Description	Password Protected
Monitoring Mode	The Monitoring Mode is the normal operating mode to monitor patients. You can change elements such as alarm limits. When you discharge the patient, these elements return to their default values. You cannot select or change grayed out items. These items are for your information only. To change these items, switch to the Configuration Mode.	no
Demo Mode	The Demo Mode is used for demonstration and training purposes. Do not change into Demo Mode during monitoring. When transducers are connected to the monitor and the recorder is on, a demo trace is recorded. But the demo trace is not transmitted to an information and surveillance system such as OB TraceVue/IntelliSpace Perinatal.	yes

Mode	Description	Password Protected
Configuration Mode	The Configuration Mode is for personnel trained in configuration tasks. You can change and store the default values and patient profiles permanently in the Configuration Mode. These tasks are described in the Configuration Guide. During installation, the fetal monitor is configured for use in your environment. This configuration defines the default settings you work with when you switch on the fetal monitor.	yes
Service Mode	The Service Mode is for trained and authorized service personnel only.	yes

A field displayed at the fetal monitor screen indicates if the monitor is in Demonstration Mode, Configuration Mode, or Service Mode. To change to a different mode, select this field.

Automatic Screen Layouts

Your monitor's preconfigured screen layouts define how measurement information is arranged on the screen. The monitor automatically applies the correct screen layout for the measurements you are monitoring. No user action is required.

Connecting or disconnecting transducers, or switching the noninvasive blood pressure measurement on or off, results in an automatic adjustment of the screen layout. When a measurement is off, its numerics are removed from the monitor's screen. The monitor stops acquiring data and generating alarms for this measurement. If you disconnect a transducer while it is performing a measurement, the monitor issues a disconnect INOP (and in the case of SpO₂, replaces the measurement numeric with a question mark).

Settings

This section describes the various settings available on the monitor.

Active Settings

What the monitor displays, and the way it operates, is controlled by its settings. They determine screen content, layout, high and low alarm limits and so forth.

The "active settings" are the current settings the monitor uses, including any adjustments made by the last user. Active settings are not permanent, but are retained after a loss of mains power.

There are also two preconfigured default settings:

- User Default
- Factory Default

User Default

The User Default is a complete configuration stored in the monitor's long-term memory. You can store the active settings, modified to your preference, in the User Default (in Configuration Mode).

In monitoring mode, you can load the User Default settings to return to your preferred settings:

- 1 Select the **Defaults** SmartKey



- 2 Select **Confirm** in the dialog box to load the User Default.

Factory Default

The Factory Default is a complete configuration pre-defined at the factory. You cannot modify it. In Configuration Mode, you can load the Factory Default as the active settings.

CAUTION

This resets all settings to factory defined values, but be aware that some values will differ from those with which the fetal monitor was originally shipped from the factory (recorder speed and paper scale type will need to be corrected, for instance). After loading the Factory Default, check the settings and, if necessary, change them to the settings you normally use.

You can use the Factory Default as the basis for producing your User Default. See the Service Guide for details.

Global Settings

General monitor configuration settings are stored in the Global Settings. These include settings for line frequency, QRS type, and whether the monitor is automatically reset to the User Default after a power interruption of more than one minute. You can change the Global Settings in Configuration Mode.

Changing Measurement Settings

Each measurement has a setup menu in which you can adjust all of its settings. You can enter a setup menu:

- 1 with the measurement numeric - select the measurement numeric on the screen to enter its setup menu. For example, to enter the **Setup FHR1** menu, select the **FHR1** (fetal heart rate 1) numeric.
- 2 with the **Main Setup** SmartKey - if you want to setup a measurement when the measurement is switched off, use the **Main Setup** SmartKey and select **Measurements**. Then select the measurement name from the pop-up list. With this SmartKey you can access any setup menu in the monitor.

This guide always describes the entry method using the measurement's setup menu. You can use the method you prefer.

Switching the Noninvasive Blood Pressure Measurement On and Off

The noninvasive blood pressure measurement is the only measurement that you can manually switch on and off. To do this:

- 1 Enter the noninvasive blood pressure measurement's setup menu.
- 2 Select **NBP** to toggle between on and off. The screen display indicates the active setting.

Changing Monitor Settings

To change monitor settings such as brightness, or touch tone volume:

- 1 Enter the **Main Setup** menu.



- 2 Select the setting you want to change, or select **User Interface** to enter a sub menu where you can change user interface settings.

Adjusting the Screen Brightness

- 1 Enter the **Main Setup** menu.



- 2 Select **User Interface**.
- 3 Select **Brightness**.
- 4 Select the appropriate setting for the screen brightness. 10 is the brightest, 1 is the least bright. **Optimum** is suitable for most situations.

Adjusting Touch Tone Volume

The touch tone is the tone you hear when you select any field on the monitor screen. To adjust the touch tone volume:

- 1 Enter the **Main Setup** menu.



- 2 Select **User Interface**.
- 3 Select **Touch ToneVolume**, then select the appropriate setting for the touch tone volume: 10 is the loudest and 1 is the quietest. Selecting zero switches the touch tone volume off.

Setting the Date and Time

The current date and time is displayed in its own element in the information section of the monitor screen.

- 1 Select the date and time screen element from the monitor's information line to enter the **Date, Time** menu.
- 2 Select, in turn, the **Year, Month, Day, Hour** (in 24 hour format) and **Minute**, as necessary.
- 3 Select **Store Date, Time** to change the date and time.

WARNING

Do not change the date and time setting, if the fetal monitor is connected to a Philips OB TraceVue/ IntelliSpace Perinatal system via a LAN -setup. The monitor uses the OB TraceVue/IntelliSpace Perinatal system date and time, including daylight saving time changes. As long as the fetal monitor is connected to the OB TraceVue/IntelliSpace Perinatal system via the LAN-setup (locomotive icon displayed on the monitor's screen), the option to change the date and time settings at the fetal monitor are disabled, this is not valid for RS232 connections, or the connection to other systems.

When disconnected from AC power, the monitor retains the date and time setting for at least two months. If the monitor is off longer than two month, and the operating system detects that the date and time settings are invalid, the monitor initiates a "cold" start and sets the date to 1 Jan 1997 and the time to 00:00.

Checking Your Monitor Revision

- 1 Select **Main Setup, Revisions** to open the **Monitor Revision** menu.
- 2 From the **Monitor Revision** menu, select the monitor component for which you need revision information.

Preparing to Monitor

Confirm fetal life before you begin fetal monitoring. Familiarize yourself with the basic operation principles before you start to monitor.

CAUTION

Check the fetal monitors housing for damage before you start to monitor as part of your safety precautions.

After you switch on the monitor:

- 1 Check that you have the correct patient cables and transducers plugged in for the measurement you want to monitor.
- 2 If you use an Avalon CL or Avalon CTS system, check if the cableless transducers are ready and charged (apparent by either a lit up green or yellow LED indicator).
- 3 Admit your patient to the monitor (see "Admitting a Patient" on page 135).
- 4 Check that the alarm limits, alarm and fetal heart rate volumes, patient category, and so forth are appropriate for your patient. Change the settings if necessary.
- 5 Refer to the appropriate measurement section for details of how to perform the measurements you require.
- 6 Start recording.

Switching On: FM20/FM30

- FM20/30**
- 1 Connect the monitor to AC mains and switch the monitor on.
 - The green power-on LED lights up.
 - The monitor performs a self-test as it starts up. **Selftest: OK**, the serial number, and revisions for the software and firmware are printed on the fetal trace paper (if recorder **Auto Start** is configured to **On**).
 - The monitor display comes on.
 - There is a start-up tone from the loudspeaker.
- FM20/30 Battery Option**
- If this option has been chosen, the green power-on LED on both the external power supply and the battery LED indicator will light up.

Switching On: FM40/FM50

- FM40/50**
- 1 Connect the monitor to AC mains.
 - The green LED lights up.
 - 2 Press the On/Standby switch.
 - The monitor performs a self-test as it starts up. **Selftest: OK**, the serial number, and revisions for the software and firmware are printed on the fetal trace paper (if recorder **Auto Start** is configured to **On**).
 - The monitor display comes on.
 - There is a start-up tone from the loudspeaker.

Adjusting the Display Angle (FM20/FM30)

- FM20/30** You can tilt the display on the FM20 and FM30 to one of five different positions, or you can fold it completely down. The tilt/fold mechanism works on a one-way ratchet system. You hear a click as each of the five positions is reached. The screen can be folded back down only after tilting the display forwards as far as it will go.

To tilt the display from the folded position:

- 1 Unlock the display by releasing the catch.



3 Basic Operation

- 2 Lift the display forward. You will hear a click as the first position engages. If you want to tilt the display further, lift the display further forward until you reach the desired angle.



- 3 To fold the display, pull the display forwards as far as it will go



- 4 Then push the display all the way back until it **clicks** shut.



If your monitor is wall-mounted, the display should be folded flat.

Fastening Belts and Transducers

You can use more than one belt if, for example, you are monitoring uterine activity and FHR simultaneously. There are two basic ways to fasten belts and transducers:

- Belts with button fixings.
- Velcro belts together with the butterfly belt clip.

What You Need

- Ultrasound transducer
- Toco MP or CL Toco⁺ MP transducer
- Ultrasound gel
- Transducer belt (and optional butterfly belt clip, if applicable)

Using Belts with Button Fixings

- 1 Place the transducer belt across the bed, so that the fixing button will face away from the mother when it is fastened.
- 2 Lie the patient on the bed and arrange the belt around her until it is tight but still comfortable.
- 3 Fasten the belt by pushing the fixing button through the overlapping section of the belt. Ensure that the fixing button and the loose ends of the belt are at the patient's side.



3 Basic Operation

- 4 When you have positioned a transducer satisfactorily, you can attach it to the belt by pushing the belt button on the transducer through one of the holes in the belt.




- 5 Alternatively, attach the butterfly belt clip to the transducer belt button and use this to attach the transducer to the belt. The clip allows you to slide the transducer for easy repositioning.



Using CL Transducers with a Belt Clip


The Avalon CL transducers have their own belt clip. They can be optional ordered in a kit of 10 (989803184851).



1 Avalon CL belt clip

Using Belt with Velcro Fixings

Insert one end of the belt between the belt guides on one side of the butterfly belt clip, and secure with the velcro fixing. Insert the other end of the belt between the belt guides on the other side of the butterfly belt clip, adjust for the correct tension, then secure with the velcro fixing.



1 Velcro fixing
2 Belt guides
3 Velcro fixing

WARNING

When connecting devices for acquiring measurements, always position cables and NBP tubing carefully to avoid entanglement or potential strangulation.

Repositioning Transducers

A patient possibly wears transducers for long periods without interruption. In rare cases, skin irritations may occur if a transducer is attached to one location for a longer period. To ensure there are no adverse effects on the patient's skin, inspect the transducer application site at least every three hours. If the skin quality changes, move the transducer to another site.

The ultrasound transducer is often repositioned to follow the fetal heart as part of the normal monitoring process, but this is not so for the Toco transducer. Therefore, remember to check its application site (between contractions) at least every three hours.

To reduce the risk of skin irritations, do not allow a cleaning or disinfecting agent to remain on the transducer. Follow all instructions that accompany the specific cleaning and disinfecting agents you are using. Remove agent residues with a cloth dampened in water before applying a transducer to a patient.

See "Care and Cleaning" on page 245 for further information, and a list of approved agents.

Connecting a Transducer to the Monitor



- 1 SpO₂ socket
- 2 Noninvasive blood pressure socket
- 3 Fetal sensor sockets



You can plug a fetal transducer, an ECG/IUP patient module, an Avalon CL or Avalon CTS Cableless Fetal Transducer System interface cable (red connector), or an external event marker into any of the four fetal sensor sockets marked by the fetal icon, or "Fetal Sensors" (depending on geography).



For measuring maternal SpO₂, connect the sensor to the socket marked with the SpO₂ icon or "SpO₂" (depending on geography).



For maternal non-invasive blood pressure, connect the cuff to the socket marked with the NBP icon or "NBP" (depending on geography).

For the FM20 and FM30, you can connect an Avalon CL Cableless Fetal Transducer system interface cable (red connector) to one of the fetal sensor sockets at the left side of the monitor.


For the FM40 and FM50, you can connect an Avalon CL or Avalon CTS Cableless Fetal Transducer System interface cable (black connector) to one of the two dedicated black sockets marked "Tele" at the rear of the monitor, as an alternative to using one of the fetal sensor sockets (red connector) at the front.




- 1 Interface cable to Avalon CL and Avalon CTS Cableless Fetal Transducer System.
- 2 Connect the black connector to one of the two black sockets (marked "Tele") on the rear of the monitor.

What You See on the Monitor

When you connect a transducer or sensor, the measurement numeric appears on the screen.



for FM20/FM30




for FM40/FM50

Fetal heart rate measurements are labeled in the order in which you plug in the transducers for those measurements. It does not matter which fetal sensor socket you use, as the monitor allocates a channel automatically. For instance, when monitoring triplets, the first transducer you connect is automatically allocated a channel, and the measurement is labeled FHR1, the second FHR2, and the third FHR3. See also chapters "Monitoring Twin FHRs" on page 171 and "Monitoring Triple FHRs" on page 179.

When you touch a measurement numeric on the screen, the setup menu for that measurement opens. The fetal sensor socket to which the transducer for this measurement is connected is identified by the transducer position indicator in the blue setup menu header.

3 Basic Operation

The blue finder LED on a wired fetal transducer lights up when you touch the measurement on the screen, allowing you to identify the corresponding transducer.



1 Finder LED

The white finder LED on a cableless fetal transducer lights up when you touch the measurement on the screen, allowing you to identify the corresponding transducer.



1 Finder LED cableless transducer

The recorder prints an annotation showing the date, time, paper speed, and monitoring mode. It repeats this every 10 minutes.

Start Recording

Checking/Setting Paper Scale

You can check the paper **Scale Type** (**US** for USA, or **Internat'l** for other geographies) in the **Fetal Recorder** menu. In Monitoring Mode, you can see these settings (grayed out), but you cannot change them. They can be changed in Configuration Mode, see the Configuration Guide.

Paper Guide FM40/FM50

FM40/50 The recorder in the FM40 and FM50 features a transparent paper guide which:

- facilitates correct alignment of the paper, both during loading and while the recorder is running. See “Loading Paper FM40/FM50” on page 84.
- incorporates a tear-off edge, which not only allows you to tear off the trace paper where you like (not necessarily at a fold), but also helps to avoid paper misalignment while doing so (see “Tearing Off the Paper” on page 80).
- is removable (see “Removing the Paper Guide: FM40/FM50” on page 250).

Switching the Recorder On and Off

In addition to the normal recording of real-time traces, you will sometimes see a trace recovery printout from the monitor's internal backup memory at high speed when the recorder is started. For details, see “Recovering Traces on Paper” on page 241.

For an explanation of the various symbols that can appear on the trace recording, see “Recorder Specifications” on page 284.

To switch the recorder on, select in **Main Setup** the menu item **Fetal Recorder**, or press one of the SmartKeys: **Start/ Stop** or **Start Recording**.



The "recorder on" status indicator is displayed in the bottom right-hand corner of the screen when you switch on the recorder.




The paper advances quickly for 2 cm and then returns to the set speed. Whenever the recorder is switched on, a trace header is printed vertically on the trace paper, containing the following:

- **Selftest: OK:** confirmation that the monitor's self-test completed successfully, and that it is ready to use.
- the software revision and firmware revision
- the serial number
- the time
- the date
- patient name and medical record number (if entered)

3 Basic Operation

- the **Recorder Speed**




- 1 Fetal heart rate label
- 2 Uterine activity label

The current monitoring modes (if any transducers are connected to the monitor) are printed. Whenever a transducer's mode is changed, the following are printed:

- the time
- the date
- trace identification symbols
- the **Recorder Speed**

The monitor prints the time, date, **Recorder Speed**, and monitoring modes in the trace header when first switched on, in a periodic time stamp every ten minutes after, and if the monitoring modes change. The time stamp begins with the symbol shown below. The data is reprinted in the header if the

time and date are locally adjusted, or if an OB system is connected that readjusts the time and date automatically.



- 1 Time stamp printed every ten minutes

The trace records maternal parameters also. When measuring noninvasive blood pressure, the annotation is made at the end of the measurement. If the noninvasive blood pressure measurement repetition time is short, the noninvasive blood pressure numeric may not always be printed.

The recording of notes (see “Entering Notes” on page 88) or time/date information may be interrupted by connecting, or unplugging a transducer, or by a change in measurement-related setting (for example, artifact suppression, Toco sensitivity, or alarm settings).

A new patient admission or a change to the paper scale setting stops all annotations, and prompts a new vertical trace header to be printed.

To switch off the recorder:

- 1 Either select **Start/Stop** from the **Fetal Recorder** menu.



- 2 Or press one of the SmartKeys (depending on configuration): fetal recorder **Start/ Stop** or **Stop Recording**.

If your recorder is configured with **Confirmed Stop** on (a Configuration Mode setting), you will need to confirm that you want to stop the recorder, before it will stop.




When the recorder is off, the "recorder off" status indicator is displayed in the bottom right-hand corner of the screen: When the **Paper Save Mode** is set to **Off** the paper icon shows an **x** mark, and when **Paper Save Mode** is **On** the icon shows a paper trace icon.



Recording Elements

The trace paper of the recorder has much information printed on it. Here is an exemplary trace with the most common print outs and their meaning. Each trace header contains the last name and first name, the patient ID, patient date of birth, the current date and time, patient's bed label, and the gestational age of the pregnancy.



- 1 Name, patient ID, date of birth, current date and time
- 2 **Recorder Speed**
- 3 Gestational age and patient bed label
- 4 FMP - Fetal Movement Profile
- 5 Other measurement for the patient such as temperature, arterial oxygen saturation, and pulse.
- 6 Entered notes such as patient repositioned, or MD notified, or others
- 7 Time stamp
- 8 Trace separation
- 9 **Recorder Speed**
- 10 Maternal ECG
- 11 Direct ECG
- 12 Coincidence of heart rate detected
- 13 HR with Cross Bell
- 14 Alarm Limits
- 15 DFHR2 (second fetal heart rate from DECG)
- 16 FHR1 (first fetal heart rate)

Choosing Paper Speed

You can choose a **Recorder Speed** of 1, 2, or 3 centimeters per minute (cm/min). The default setting is 3 cm/min.

The ACOG technical bulletin on FHR monitoring states that "*accurate pattern recognition is difficult if not impossible at 1 cm/min and that 1 cm/min is only recommended for more economic screening. When FHR abnormalities arise, the faster paper speeds will enhance FHR pattern recognition*".

Additionally, because a change in **Recorder Speed** results in a change in the appearance of an FHR trace, you are advised to ensure ALL monitors in your institution are set to the same speed.

To set the **Recorder Speed** (in Configuration Mode) see the Configuration Guide.

Advancing the Paper

You can advance the paper automatically to the next fold by pressing the **Paper Advance** SmartKey at any time except during a stored data recording. This is also possible using the **Fetal Recorder** menu.



Marking an Event


You can record significant events on the trace paper (for example, when pain medication is administered or when the mother changes position). The mother can use the remote event marker to mark events herself. You connect the remote event marker to any free fetal sensor socket.

To mark an event on the trace paper you can:

- 1 Either select the **Set Marker** SmartKey.




- 2 Or press the button on the remote event marker. The remote event marker is connected to the monitor via any fetal transducer socket.



A small arrow is printed on the heart rate scale on the trace paper.



This reflects exactly when the marker button was first pressed; keeping the button pressed has no influence on the annotation.



Tearing Off the Paper

CAUTION

Never pull on the paper to advance it, as this can cause misalignment of the paper. Always tear off the paper along the perforation.


FM40/FM50 The recorder's paper guide incorporates a tear-off edge, allowing you to tear off the trace paper cleanly where you like (not necessarily at a fold). When **not** using the paper guide, always tear off the paper along the perforation.

To tear off the trace paper after monitoring **using the paper guide**:

- 1 If the recorder is running (the "recorder on" status indicator is displayed), turn off the recorder by selecting the fetal recorder **Start/ Stop** SmartKey or the **Stop Recording** SmartKey.



- 2 Tear off the paper as shown in the picture. To ensure a clean tear, always tear in an upwards motion, as indicated by the arrows. You can start tearing from the left or right (right-handed user shown). You may want to use both hands to guarantee that the paper is not misaligned during the tear off.



- 3 If you wish to tear off the paper at a fold, select the **Paper Advance** SmartKey, wait for the paper to stop, then tear it off.



FM20/FM30 The FM20/FM30 does not have a paper guide. The procedure is the same as described for the FM40/FM50 with the exception that you should advance the paper to a perforation.

If you wish to tear off the paper at a perforation, select the **Paper Advance** SmartKey, wait for the paper to stop, then tear it off.



Paper-Out Indication

Each pack of paper has 150 pages. The monitor issues a paper-out warning in the status line at the bottom of the screen, when there are five pages to go. If you switch on the recorder or press the paper advance key when there are fewer than five pages remaining, it may take two pages before the alarm is activated. Load a new pack in time.

If the recorder runs out of paper, an audible paper-out alarm is sounded, if so configured. See “Loading Paper FM40/FM50” on page 84 and “Loading Paper FM20/FM30” on page 82 to learn how to reload paper.

Fetal traces continue to be recorded into the monitor's backup memory, and can be retrieved and printed completely if new paper is loaded within one hour, when the **Bridge Paperout** setting is enabled in Configuration Mode. See “Recovering Traces on Paper” on page 241 for further information.

Loading Paper FM20/FM30

CAUTION

Using recorder paper that is not approved by Philips can result in accelerated paper fading and can damage the thermal line printhead. This type of damage is not covered by warranty.

FM20/FM30 To load a pack of paper:

- 1 If the recorder is on, press the recorder **Start/ Stop** SmartKey or the **Stop Recording** SmartKey to turn it off before loading a new pack of paper.




- 1 Press the paper table release to unlock the paper drawer and then pull the table forward to open it fully.
- 2 Lift out any remaining paper from the tray.




- 3 Prepare to place the new pack of paper in the tray with the bottom side down. The bottom side is indicated by the word **STOP** printed on the final page of the new pack.

- 4 Unfold the top page of the pack and position the uterine activity scale on the right.



- 5 Slide the pack into the tray.



3 Basic Operation

- 6 Push the paper drawer back until it "clicks" closed.



- 7 Press the recorder **Start/ Stop** SmartKey or the **Start Recording** SmartKey to switch on the recorder.



Annotations of trace information are printed on the trace paper (see “Switching the Recorder On and Off” on page 75 for details).

Loading Paper FM40/FM50

CAUTION

Using recorder paper that is not approved by Philips can result in accelerated paper fading and can damage the thermal line printhead. This type of damage is not covered by warranty.

FM40/FM50 To load a pack of paper:

- 1 If the recorder is on, press the **Start/ Stop Rec** SmartKey or the **Stop Recording** SmartKey to turn it off before loading a new pack of paper.



- 2 Press the paper eject button to open the paper drawer.



- 3 Lift out any remaining paper from the tray.
- 4 Press and hold the paper eject button to partially eject the paper, thus making it easier to remove.



- 5 Hinge the transparent paper guide forward. It is held in the closed position by a small protrusion on each side of the holder.



3 Basic Operation


- 6 **A** - Protrusion holds paper guide in closed position.




- 7 Prepare to place the new pack of paper in the tray with the bottom side down. The bottom side is indicated by the word *STOP* printed on the final page of the new pack.
- 8 Unfold the top page of the pack and position the uterine activity scale on the right.
- 9 Slide the pack into the tray.




- 10 Feed the paper evenly through the paper guide. Do not close the paper guide yet.




- 11 Close the paper drawer.




- 12 Now close the paper guide.



- 13 Press the recorder **Start/ Stop** SmartKey or the **Start Recording** SmartKey to switch on the recorder.



Annotations of trace information are printed on the trace paper (see “Switching the Recorder On and Off” on page 75 for details).



Entering Notes

Your monitor has a set of 15 factory pre-configured notes (see below). The maximum length of one single note is 30 characters. It is possible to edit the notes in Configuration mode (see the Configuration Guide).

To enter a note:

- 1 Press the **Enter Notes** SmartKey to open the **Enter Note** menu.



- 2 Scroll if necessary, then select the note you wish to enter. A confirmation dialog box opens:
- 3 Select **Confirm** to enter the note. The note is then shown in the status line of the display, and is annotated on the fetal trace if the fetal recorder is on.

- 4 By default, notes are printed lengthwise in the direction of the trace, in the space between the FHR grid and the uterine activity grid. If you prefer, you can configure the recorder to print across the trace.
- 5 You can change this in Configuration Mode by changing the **Notes Recording** setting in the **Fetal Recorder** menu from **Along** (default) to **Across** (notes print width wise across the trace).

The following are the pre-configured notes from which to choose:

- 1 **Patient Repositioned**
- 2 **Vaginal Examination**
- 3 **MD Notified**
- 4 **Sitting**
- 5 **On Back**
- 6 **Left Lateral**
- 7 **Ambulating**
- 8 **Tocolytic Given**
- 9 **Membranes Ruptured**
- 10 **Amniotomy**
- 11 **Amniotic Fluid Clear**
- 12 **Amniotic Fluid Not Clear**
- 13 **Oxytocin**
- 14 **Urinary Catheter**
- 15 **Micro Blood Analysis**

Up to two notes can be printed directly, and the monitor can temporarily store up to a further two notes, and these are printed after the first two have been recorded. Any further notes are discarded. For example, if you enter six notes in quick succession, the first two notes you entered are recorded right away, the next two are stored in memory and then printed when the first two have been recorded, and the last two are discarded.

If the printing of two notes happens to coincide with the regular recording of the time stamp that takes place once every ten minutes, the time stamp is delayed until the notes have finished printing.

Signal Quality

During monitoring, if the fetal heart rate signal quality fluctuates, and becomes poor, it does not necessarily mean that the transducer needs repositioning. The fluctuation may be caused by fetal movement. Allow time for the signal to stabilize before deciding whether to reposition the transducer (ultrasound), or apply a new electrode (ECG). For the best trace quality, the signal quality indicator should be full, indicating good signal quality, even though it may be possible to make traces at a lower signal quality level.

Cableless Monitoring

Basics of Cableless Systems

The fetal monitors FM20/FM30 and FM40/FM50 are compatible with the Avalon CL and Avalon CTS Cableless Fetal Transducer system. Regard the following points for cableless monitoring:

- You can connect one Avalon CL with a red connector to a FM20/FM30 or a FM40/FM50 (fetal socket).
- You can connect up to two Avalon CLs with a black connector to a FM40/FM50 (Tele socket).
- You can connect one Avalon CTS system to a FM20/FM30 or FM40/FM50 at a time (either fetal or Tele socket).
- You **cannot** connect Avalon CTS and an Avalon CL at the same time to the same fetal monitor.
- Monitoring multiple pregnancy using cableless transducers is supported by the Avalon CL system only.
- Using a mixture of wired and cableless fetal transducers is not supported. You can use *either* wired *or* cableless fetal transducers.
- If you cannot get sufficient signal quality using Avalon CTS transducers, switch to wired transducers.

Avalon CL Transducer System



Avalon CTS Transducer System



Assigning Cableless Transducers

The cableless transducers of the Avalon CL are assigned by simply docking them at the base station connected to the patient's fetal monitor. See the Avalon CTS Instruction for Use for the assignment of the cableless transducers of this system.

Activating Cableless Transducers

The CL transducers of the Avalon CL are activated by picking them up from the base station. If the CL transducer is activated, an antenna symbol is displayed on the fetal monitor screen next to the numeric of the measurement. See the Avalon CTS Instruction for Use for the activation of the cableless transducers of this system.

Unassigning Cableless Transducers

The CL transducers of the Avalon CL are unassigned by manually removing them in the corresponding setup menu from the group of assigned cableless devices, or by unassigning them directly in the **Tele Info** window. See the Avalon CTS Instruction for Use for the unassignment of the cableless transducers of this system.

Deactivating Cableless Transducers

The CL transducers of the Avalon CL are deactivated by redocking them at a base station. To deactivate all cableless transducers at once press the Standby button of the Avalon CL base station. See the Avalon CTS Instruction for Use for the deactivation of the cableless transducers of this system.

Configuration of Cableless Systems

For the Avalon CL and Avalon CTS Cableless Fetal Transducer system and the corresponding channels have to be configured for the radio communication to work and not to interfere with any other telemetry devices.

The configuration of a connected Avalon CL system is done in Configuration Mode or Service Mode of the fetal monitor. Service Mode functions can be used to identify channel assignment conflicts in the hospital environment. The configuration should be carried out by authorized and qualified service personnel, either by the hospital's biomedical department, or by Philips Support.

Setup OBR	
Revisions	A.00.00
Frequency Band	WMTS
Channel	4
Scan Duration	15 min
Start Scan	
Stop Scan	
Start RSSI Trace	
Stop RSSI Trace	


For a detailed description of the configuration see the Avalon CL Service Guide and the Fetal Monitor Configuration Guide. For details regarding the configuration of a connected Avalon CTS system see the Avalon CTS Service Guide.

Tele Info Window

The **Tele Info** window of the fetal monitors allows you to manage the Avalon CL system. Below is an exemplary view of the **Tele Info** window. Your window may differ depending on the status of the connected system.

To open the **Tele Info** window you can:





- 1 Configure a **Tele Info** SmartKey to open the window,
- 2 click on the TELE icon on the main screen, or
- 3 select the **Main Setup** and then **Tele Info**.



- 1 Transducer icon assigned with the parameter label
- 2 Cableless module icon and equipment ID
- 3 Base station icon with docking indication (the white slot indicates a charging transducer)
- 4 Key Remove
- 5 Key Find
- 6 Key Battery Report
- 7 Tele icon
- 8 Remaining battery time
- 9 Out of battery icon
- 10 Indication of signal quality

Telemetry

When the monitor recognizes a connected Avalon CL or Avalon CTS interface cable (red or black connector), it confirms the recognition with the following status indicators displayed in the lower right-hand corner of the screen:

Indicator	Avalon CTS	Avalon CL
	Avalon CTS interface cable is connected to the monitor, but the Avalon CTS base station is not connected to the interface cable, or it is disconnected from AC mains, or is in Standby mode.	A base station is connected to the monitor, but the base station is in Standby mode.
	Avalon CTS interface cable is connected to the monitor, Avalon CTS base station is connected, powered on, and cableless transducers are ready to use, but no cableless transducers are currently active (all are still docked at the base station).	A base station is connected to the monitor but no cableless transducers and no cableless Pods are currently active. All are still docked at the base station, or the base station holds no transducers.
	Avalon CTS interface cable is connected to the monitor, Avalon CTS base station is connected, powered on, and at least one cableless transducer is active and assigned.	A base station is connected to the monitor, it is on, and minimum one cableless transducer or cableless Pod is active and assigned.
	Not applicable for Avalon CTS.	The CL transducer is moved away from the base station approaching the end of range (at least 70 m). The US sound is replaced by an artificial QRS sound (like DECG). The DECG and MECCG waves are no longer displayed.

NOTE

A cableless transducer is considered still active, even if it has an INOP condition (e.g. **cl US Disconnect**), until it is either deactivated by docking it at the base station, or until it is manually removed from the **Tele Info** window.

CL transducers have priority over wired transducers. If an Avalon CL or Avalon CTS base station is connected to the fetal monitor, and both wired and cableless transducers are connected to the monitor, the wired transducers are disabled whenever one cableless transducer is active. If an CL Pod is activated it does not automatically deactivate all wired transducers, just the one with the corresponding measurement.

To switch back to using wired transducers, redock the cableless transducers at the base station, or switch the base station to Standby mode, and continue monitoring with the wired transducers. If you cannot get sufficient ultrasound signal quality using the Avalon CTS base station, try repositioning the transducers, or switch to wired transducers. Using an Avalon CL base station you should get a better signal by reducing the distance between CL transducer and base station.

- When using an Avalon CL or Avalon CTS you should be aware that FMP is not recommended when the mother is likely to move, and you should disable Fetal Movement Profile (FMP) on the fetal monitor (Fetal Movement Off) if the mother is walking. Maternal movements are likely to create artifact in the FMP output. See also “Switching FMP On and Off” on page 161, the sections “Cableless Monitoring - Important Considerations” on page 156 and “Fetal Movement Profile” on page 160.

- With the Avalon CL Transducer system you can now monitor twins and triplets with CL transducers. The Avalon CTS System does not have this option.

About RF Signal Quality

Signal transmission can be disturbed if:

- the patient is out of range of the receiving area.
- there is interference from another, possibly stronger, RF signal (a broadcasting station, for instance).
- the patient is near material that absorbs electromagnetic waves (for example, metal-reinforced concrete, elevator doors) or the base station antenna is in an enclosed metal rack.

Pop-up Keys

Pop-Up Key	Function
Remove	Selecting the Remove key deactivates and unassigns the selected active transducer or cableless Pod and removes it from the list. This key is disabled if no active device is selected.
Find	Selecting the Find key pages the selected active CL transducer or CL Pod. This key is disabled if no active device is selected.
Battery Report	Selecting the Battery Report key prints a combined battery report on the built in recorder. The report contains at least the following information for each CL transducer or CL Pod currently in the list: <ul style="list-style-type: none"> • Cableless device identification (like in the Tele Info window list) • Battery identification (e.g. part no.) • Battery serial number • Battery manufacturing date • Battery capacity: remaining/when fully charged/battery temperature • Number of cycles This key is only available and enabled in Service mode.

Prompts

The fetal monitor issues prompt messages to certain user interactions. For example if the CL SpO₂ Pod is picked up from the base station to activate it, the monitor displays the prompt message **cl SpO₂ Added** and the equipment label of the Pod.

CL Transducer Assignment

The cableless transducers of the Avalon CL are assigned by simply docking them at the base station connected to the patient's fetal monitor.

Since twins and triplets can be monitored with the Avalon CL base station, it can be delivered with up to seven CL transducers: one CL Toco⁺ MP, three CL US, and three CL ECG/IUP transducers. Only up to four CL transducer can be assigned at one time to the Avalon CL base station. Four assigned CL transducers enable monitoring combinations such as 1 CL Toco⁺ MP transducer + 3 CL US transducers for monitoring triplets externally, or 1 CL Toco⁺ MP transducer + 2 CL US transducer + 1 CL ECG/IUP transducer for monitoring twins.

Three CL transducers can be docked and charged at the Avalon CL base station at the time. If a fourth transducer should be assigned to the Avalon CL base station:

- 1 Remove first the three already docked CL transducers from the base station and place them on the patient.
- 2 Then dock the remaining fourth transducer at the base station and wait until the LED of the fourth transducer signals that it has communication with the base station.
- 3 Place the fourth CL transducer on the patient and check on the fetal monitor in the **Tele Info** window if the base station and the monitor have communication with all four transducers, and if their numerics are displayed on the screen of the fetal monitor.

CL Transducer Unassignment

The cableless transducer are not unassigned by simply removing them from the base station. To unassign a cableless transducer, there are two methods:

- 1 Open the setup menu of the measurement e.g. **FHR1** by pressing the numeric on the screen.
- 2 Select **Remove** from the open setup menu.
- 3 A confirmation window opens with the prompt **Please Confirm** and **Device will be removed**.
- 4 Press the pop-up key **Confirm**. A prompt states then **clDev has been removed**.

or

- 1 Open the **Tele Info** window.
- 2 Select the symbol of the cableless transducer you want to unassign, and select **Remove** from the pop-up menu.
- 3 A confirmation window opens with the prompt **Please Confirm** and **Device will be removed**.
- 4 Press the pop-up key **Confirm**. A prompt states then **clDev has been removed**.

NOTE

When a cableless transducer is unassigned, all the measurements from this transducer are no longer monitored; e.g. if you unassign a CL Toco⁺ MP, the Toco measurement, SpO₂, and MECG are possibly no longer monitored.

CL Pods Assignment

The cableless Pods of the Avalon CL are assigned by simply docking them at the base station connected to the patient's fetal monitor.

CL Pods Unassignment

To unassign a cableless Pods, there are four methods:

- 1 Open the setup menu of the measurement e.g. **SpO₂** by pressing the numeric on the screen.
- 2 Select **Remove** from the open setup menu.
- 3 A confirmation window opens with the prompt **Please Confirm** and **Device will be removed**.
- 4 Press the pop-up key **Confirm**. A prompt states then **clDev has been removed**.

or

- 1 Open the **Tele Info** window.
- 2 Select the symbol of the cableless Pod you want to unassign, and select **Remove** from the pop-up menu.

3 Basic Operation

3 A confirmation window opens with the prompt **Please Confirm** and **Device will be removed**.

4 Press the pop-up key **Confirm**. A prompt states then **clDev has been removed**.

or

Unassign the cableless Pods via it's user interface. See the Cableless Measurement Instruction for Use.

or

Dock the cableless Pod at another base station or charger and therefore assigning it to the new device.

Paging Patients

To page a patient currently not near the fetal monitor and base station select the SmartKey **Call Patient** at the fetal monitor. The base station will then select the cableless transducer with the best signal quality to emit the tone sequence to notify the patient.

Underwater Monitoring

Only the cableless transducers of the Avalon CL and Avalon CTS systems can be used to monitor under water. You can use them to monitor patients in a bathtub or shower. This does not apply to the maternal cableless SpO₂ and NBP Pods, do not immerse the CL Pods into water.

Cableless transmission distances are shorter when monitoring under water. A metal bathtub is likely to further reduce the operating range.

WARNING

Never immerse the base station in liquid. You must protect it against water sprays or splashes. Place the base station where there is no chance of contact with, or falling into water or other liquids.

CAUTION

Avoid the use of pulsating water jets in the bath or shower while monitoring, as these can be misinterpreted as an incorrect (or totally artificial) heart rate.

Toco Baseline drift: The accuracy specified for baseline drift cannot be guaranteed for underwater usage. When using transducers under warm water the temperature increase causes a significant baseline change due to internal pressure increase. The depth under water at which the Toco transducer is used also has an effect on the Toco baseline, as the water pressure increases with depth. After immersion, allow one to two minutes for the pressure to stabilize, then adjust the Toco baseline (between contractions), and check it frequently.

When using the transducers underwater, the radio transmission range is reduced, and signal loss may occur.

Water-proof belts like M1562B are recommended for the cableless transducers when monitoring in water.

Switching the Monitor to Standby

To switch the monitor to Standby:

Either

- 1 Select the **Monitor Standby** SmartKey.



Or

- 1 Enter the **Main Setup** menu using the SmartKey.



- 2 Select **Monitor Standby**.
- 3 Pressing any key or selecting any field on the screen will resume monitoring.

Note

If a Avalon CL base station is connected to your FM 20/30 monitor, do not turn off the monitor if you need to recharge the batteries of the CL transducers. They can only be recharged if the FM20/30 monitor is on. If an Avalon CL base station is connected to an FM40/50, the batteries of the CL transducers can be recharged while the monitor is in Standby mode.

After Monitoring

- 1 Discharge the patient.
- 2 Remove the transducer from the patient and, using a soft tissue, remove any gel from it. Then clean the transducer.
- 3 Dock CL transducers to their base station so they can recharge.
- 4 Tear off the paper at the fold. To avoid misalignment of the recorder mechanism, NEVER pull on the paper to advance it, or try to tear it other than at a fold (unless using the paper guide with the FM40/FM50).
- 5 Switch off the monitor.

Disconnecting from Power

FM20/30 To disconnect the monitor from AC power, switch the monitor off using the On/Off switch located on the right side of the device, or unplug the power cord from the AC mains socket.

FM40/FM50 and FM20/FM30 with Battery Option The On/Standby button does not disconnect the monitor from the AC power source. To disconnect, unplug the power cord from the AC mains socket. Note that if the power cord is unplugged from the AC mains socket before the monitor is put into Standby, a beeper is activated. The beeper warns you if the monitor is accidentally disconnected from AC mains.

Power On/Power Off Behavior

The general rules determining the behavior of the fetal monitors when connected to, or disconnected from power are as follows:

- A fetal monitor that was switched on prior to a temporary power loss, switches on again when power is restored.
- A fetal monitor that was switched off prior to a temporary power loss, remains off when power is restored.
- When AC mains power is lost, a battery powered monitor (FM20/30) continues to run without interruption on battery power.

Monitoring After a Power Failure

- A fetal monitor that was switched on prior to a temporary power loss switches on again when power is restored.
- A fetal monitor that was switched off prior to a temporary power loss remains off when power is restored.
- If the fetal monitor is without power for less than one minute, monitoring will resume with all active settings unchanged.
- If the fetal monitor is without power for more than one minute, the behavior depends on your configuration.
 - If **Automat. Default** is set to **Yes**, the default profile will be loaded when power is restored.
 - If **Automat. Default** is set to **No**, all active settings are retained, if power is restored within 48 hours. The **Automat. Default** setting is made in Configuration Mode.


Troubleshooting

Problem	Possible Causes	Solutions
Light or no trace.	Wrong paper.	Use recommended paper.
	Dirty printhead.	Clean printhead. See “Cleaning the Print Head” on page 255.
	FM20/30 only: Paper misaligned due to drawer not being correctly shut.	Shut the drawer fully, pushing evenly with both hands.
End of paper noted when pack not finished.	Bad paper feed or wrong paper.	Check paper feed and use recommended paper.
Check Paper INOP is displayed.	INOP messages always indicate equipment problems.	See “Patient Alarms and INOPs” on page 121.
FetRec Equip Malf INOP is displayed.		
Paper End INOP is displayed.		
Wrong Paper Scale INOP is displayed.		

FM20/30 Battery Option

You can switch between battery-powered and mains-powered (AC) operation without interrupting monitoring.

The monitor is connected to the AC mains power via the external power supply.




NOTE

The battery option is not available for the FM40/FM50.

External Power Supply M8023A

The external power supply M8023A (option #E25) allows you to operate the fetal monitor from an AC (alternating current) power source of 100 V to 240 V ($\pm 10\%$) and 50/60 Hz ($\pm 5\%$). If this option is used, then the M8023A (option #E25) power supply is included for FM20/30.



- 1 AC power cord. Connect to AC mains socket.
- 2 Measurement Link (MSL) cable. Supplies power to the monitor for operation and for battery charging.
- 3 Power-on LED. The green light is on when the external power supply is connected to AC mains.

WARNING

- Always use the supplied power cord with the earthed mains plug to connect the external power supply M8023A (option #E25) to an earthed AC mains socket. Never adapt the mains plug from the power supply to fit an unearthed AC mains socket.
- Do not use AC mains extension cords or multiple portable socket-outlets. If a multiple portable socket-outlet without an approved isolation transformer is used, the interruption of its protective earthing may result in enclosure leakage currents equal to the sum of the individual earth leakage currents, so exceeding allowable limits.
- Do not connect any devices that are not supported as part of a system.
- Any non-medical device placed and operated in the patient's vicinity must be powered via an approved isolation transformer that ensures mechanical fixing of the power cords and covering of any unused power outlets.

Using Batteries

The Lithium Ion batteries used in your monitor store a large amount of energy in a small package. This allows reliable battery-operated monitoring but also requires care in use and handling of the batteries. Follow the instructions in this chapter and see the Service Guide for further details.

Battery Power Indicators

The battery LED and battery status information on the Main Screen, in combination with INOP messages and prompts, help you keep track of the battery power status. The indicators always show the remaining capacity in relation to the battery's actual maximum capacity, which may lessen as the battery ages. You can see the actual capacity in the **Battery Status** window.

Battery LED

The possible battery LED status on the right side of the monitor is described in the table below.



Battery LED Colors	If the monitor is connected to mains power, this means	If the monitor is running on battery power, this means
Green	Battery power is > 90%	
Yellow	Battery charging (battery power < 90%)	
Red, flashing		Less than 10 minutes power remaining
Red, flashes intermittently	Battery or charger malfunction	Battery malfunction

Battery Status on the Main Screen





Battery status information is permanently displayed on all screens. It shows the status of the battery, with the battery power remaining, with an estimate of the monitoring time this represents.

Battery power gauge: This shows the remaining battery power. It is divided into sections, each representing 20% of the total power. If three sections are filled, as in this example, this indicates that 60% battery power remains. If no data is available from the battery, a question mark is shown in the gauge.

Monitoring Time Available: Below the battery power gauge a time is displayed. This is the estimated monitoring time available with the current battery power. Note that this time fluctuates depending on the system load (the display brightness, the recorder configuration, and how many measurements you carry out).

Battery malfunction symbols: Normal battery function is indicated by the battery power gauge, together with the remaining operating time, on the main screen. You are informed of problems or changes in the status of the battery by the battery status/malfunction indicator. This consists of a blank battery gauge containing a "!" symbol. If the symbol is red, this indicates a critical situation. You can check the specific cause of the problem by looking at the symbol(s) displayed in the **Battery Status** window.

Battery Status Indicator		Battery Malfunction Indicator	
	<p>Alternates with the battery gauge on the main screen.</p> <p>Check in the Battery Status window to see which status symbol is displayed to identify the cause.</p>		<p>The red exclamation mark flashes. Critical battery situation or malfunction.</p> <p>Check in the Battery Status window to see which malfunction indicator is displayed, or refer to the INOP, to identify the cause.</p>

Battery Status Window










- To access the **Battery Status** window and its associated pop-up keys, select the battery status information on the screen, or select **Main Setup, Battery**.

Battery Status		
		TimeToFull: 1:30
Capacity remaining	[mAh]	6072
Capacity fullCharge	[mAh]	6686
Voltage	[V]	11.63
Current	[mA]	406
Temperature	[°C]	25.3

- **Capacity, remaining** tells you how much power is left in the battery.
- **Capacity, fullCharge** tells you how much power the battery can hold when fully charged.
- **TimeToEmpty** tells you approximately how long you can continue to use the monitor without an AC connection. The time span fluctuates depending on the system load (the display brightness and how many measurements you carry out), and the remaining capacity of the battery.
- **TimeToFull** is shown in place of **TimeToEmpty** if the monitor is connected to a power supply, and tells you how much time is left until the battery is charged to 90%. If >10 hrs is shown here, the battery may not charge completely when the monitor is in use.

Battery status/malfunction symbols: If a problem is detected with the battery, an INOP may be issued, and the following symbols are displayed in the **Battery Status** window, where they may be accompanied by a status message providing more details.

Symbols indicating critical situations are colored red.

Battery Status Symbols		Battery Malfunction Symbols	
	Battery is empty		(Red) incompatible battery
	Battery not charging as the temperature is above or below the specified range		(Red) battery malfunction
	Battery requires maintenance		(Red) battery has no power left
	Charging stopped to protect the battery		(Red) battery temperature too high
			(Red) battery is missing, insert battery

Checking Battery Charge

To check the charge status of a battery, refer to the battery power gauge on the screen, or open the **Battery Status** window.

When Battery Lifetime is Expired

When the battery is aged, either after 3 years from manufacturing date or after 500 charge/discharge cycles, it is recommended to replace the battery. To remind you of this, a message will appear in the **Battery Status** window. This message will only be displayed at the appropriate time when the date and time on the monitor is correct.

Optimizing Battery Performance

The performance of rechargeable batteries may deteriorate over time. Battery maintenance as recommended here can help to slow down this process.

Display Brightness Setting

- 1 In the **Main Setup** menu, select **User Interface, Brightness**, then **Optimum**. This selects a level of brightness suitable for most monitoring locations that uses less battery power than brighter settings.
- 2 Ensure that the current level of brightness is suitable for your monitoring location.

Charging the Battery

To charge the battery,

- 1 Connect the monitor to the external power supply (M8023A).
- 2 Charge the battery until it is full, the battery LED is green, and the battery power gauge is filled.

During high load of the monitor (usage of multiple measurements) the battery may not charge. To remedy this:

- reduce the load by removing measurements,
- reduce the screen brightness or stop the recorder.

Internal temperature conditions can also cause the battery to not charge. This is necessary to protect the battery from damage and does not indicate a malfunction. Keep the monitor at room temperature and move it away from heat sources or out of direct sunlight. The battery will resume charging when the temperature is within range again. The battery will charge more quickly if the monitor is switched off.

Conditioning the Battery

You must condition the battery when the "battery requires maintenance" symbol shows on the screen. Do not interrupt the charge or discharge cycle during conditioning.

CAUTION

Condition the battery with a monitor not currently in use. The monitor switches off automatically when there is no battery power left.

To condition the battery,

- 1 Charge the battery until it is completely full. Open the **Battery Status** window and check that the **Batt fully charged** message is displayed.
- 2 Disconnect the monitor from mains power, and let the monitor run until there is no battery power left and the monitor switches itself off.
- 3 Reconnect the monitor to mains power and charge the battery until it is full for use or charge to 50% for storage.

Storing the Battery

The battery should not remain inside the monitor if it is not used for a longer period of time. Batteries should be charged to a maximum of 50% for storage. For battery removal or exchange, contact your service personnel.

NOTE

The battery will discharge over time in a monitor that is not connected to AC via the external power supply (M8023A). The reported values for "remaining capacity" and "runtime" will become less accurate when the battery is stored in this way for a longer period of time (that is, several weeks).

Battery Safety Information

WARNING

Use only **Philips batteries** part number M4605A. Use of a different battery may present a risk of fire or explosion.

Do not open batteries, or dispose of them in fire, or cause them to short circuit. They may ignite, explode, leak or heat up, causing personal injury.

If battery leakage should occur, avoid contact with skin. Refer to qualified and authorized service personnel.

Dispose of used batteries promptly and in an environmentally-responsible manner. Do not dispose of the battery in normal waste containers. Consult your hospital administrator to find out about local arrangements. Do not expose batteries to liquids.

Do not crush, drop or puncture batteries - mechanical abuse can lead to internal damage and internal short circuits which may not be visible externally.

If a battery has been dropped or banged against a hard surface, whether damage is visible externally or not:

- discontinue use
- dispose of the battery in accordance with the disposal instructions above.

Keep batteries out of the reach of children.

CAUTION

Do not disassemble, heat above 100°C (212°F) or incinerate the batteries, to avoid the risk of fire and burns. Keep batteries out of the reach of children and in their original package until you are ready to use them.

If battery leakage should occur, use caution in removing the battery. Avoid contact with skin. Refer to qualified and authorized service personnel.

Do not install or use pre-damaged batteries.

Alarms

The alarm information here applies to all measurements. Measurement-specific alarm information is contained in the sections on individual measurements.

The fetal monitor has two different types of alarm: patient alarms and INOPs.

Patient Alarms

Red and yellow alarms are patient alarms. A red alarm indicates high priority, such as a potentially life threatening situation (for example, SpO₂ below the desaturation alarm limit). A yellow alarm indicates a lower priority alarm (for example, a fetal heart rate alarm limit violation).

INOPs

INOPs are technical alarms. They indicate that the monitor cannot measure and therefore not detect critical conditions reliably. If an INOP interrupts monitoring and alarm detection (for example, **MECG Leads Off**), the monitor places a question mark in place of the measurement numeric and sounds an audible tone. INOPs without this tone indicate that there may be a problem with the reliability of the data, but that monitoring is not interrupted.

Most INOPs are light blue, however there are a small number of INOPs which are always yellow or red to indicate a severity corresponding to red and yellow alarms. The following INOPs can also be configured as red or yellow INOPs to provide a severity indication:

- **MECG Leads Off**
- **<SpO₂ Label> No Pulse**
- **Tele Disconnected**
- **Battery Empty / Replace Battery**

Alarm Delays

There is a delay between a physiological event at the measurement site and the corresponding alarm indication at the monitor. This delay has two components:

- **The general measurement delay time is the time between the occurrence of the physiological event and when this event is represented by the displayed numerical values.** This delay depends on the algorithmic processing and, for certain measurements (SpO₂, EEG and BIS), on the configured averaging time. The longer the averaging time is configured, the longer it takes until the numerical values reflect the physiological event.
- **The time between the displayed numerical values exceeding an alarm limit and the alarm indication on the monitor.** This delay is the sum of the alarm delay configured for the specific measurement plus the system alarm delay. The system alarm delay is the processing time the system needs for any alarm on the monitor to be indicated after the measurement has triggered the alarm. See the performance specifications in “Specifications and Standards Compliance” on page 269 for the system alarm delay specification.

Multiple Alarms

If more than one alarm is active, the alarm messages are shown in the alarm status area in succession. An arrow symbol next to the alarm message informs you that more than one message is active.

The monitor sounds an audible indicator for the highest priority alarm. If more than one alarm condition is active in the same measurement, the monitor announces the most severe.

WARNING

Alarm systems of the monitor and those of the connected OB system are independent and not synchronized.

Alarm Mode

You can configure the alarm mode for your fetal monitor. There are two possible modes:

- **All:** alarms and INOPs are enabled, with all audible and visual indicators active.
- **INOP only:** only INOPs are enabled, with audible and visual indication active. This is the default alarm mode.

WARNING

In **INOP only** mode, no fetal/maternal patient alarms are enabled or indicated.

The alarm status area for yellow and red alarms shows the **INOP only** indication in conjunction with the "Alarms Off" symbol. No alarm limits or alarm off icons are displayed. No fetal/maternal patient alarm settings are available in the setup menus.

Nurse Call Systems

If configured to do so, red, yellow and light blue alarms are indicated on a nurse call system which is connected to the optional nurse call relay.

Visual Alarm Indicators

Alarm message: An alarm message appears in the alarm status area on the second line at the top of the screen indicating the source of the alarm. If more than one measurement is in an alarm condition, the message changes every two seconds, and has an arrow at the side. The background color of the alarm message matches the alarm priority: red for red alarms, yellow for yellow alarms, and light blue for INOPs. The asterisk symbols (*) beside the alarm message match the alarm priority: *** for red alarms, ** for yellow alarms. INOPs are displayed without asterisks.

Depending on how your monitor is configured, it may display alarm limit violation messages:

- in text form, for example **** FHR1 Low** or
- in numeric form, for example ****FHR1 94<110**, where the second number shows the currently set alarm limit, and the first number shows the value at which that alarm limit was violated by the widest margin.

Flashing numeric: The numeric of the measurement in alarm flashes.

Bright alarm limits: If the alarm was triggered by an alarm limit violation, the corresponding alarm limit on the monitor screen is shown more brightly.

Audible Alarm Indicators

The audible alarm indicators configured for your fetal monitor depend on which alarm standard applies in your hospital. Audible alarm indicator patterns are repeated until you acknowledge the alarm by switching it off or pausing it, or until the alarm condition ceases (if audible alarm indication is set to non-latching).

WARNING

Do not rely exclusively on the audible alarm system for fetal monitoring. Adjustment of alarm volume to a low level or off during monitoring may result in a dangerous situation. Remember that the most reliable method of fetal monitoring combines close personal surveillance with correct operation of monitoring equipment.

Alarm Tone Configuration

The audible alarm indicators of your monitor are configurable. In the monitor's Configuration Mode, you can:

- increase the alarm volume of unacknowledged alarms at regular intervals
- change the interval between alarm sounds (ISO/IEC Standard alarms only)
- change the base volume of the red and yellow alarm tones and the INOP tones
- change the alarm sound to suit the different alarm standards valid in different countries.

Standard Philips Alarms

- Red alarms: a high pitched sound is repeated once a second.
- Two-star yellow alarms and yellow INOPs: a lower pitched sound is repeated every two seconds.
- One-star yellow alarms (short yellow alarms): the audible indicator is the same as for yellow alarms, but of shorter duration.
- Standard INOPs: an INOP tone is repeated every two seconds.

ISO/IEC Standard 9703-2 Audible Alarms

- Red alarms and red INOPs: a high pitched tone is repeated five times, followed by a configurable pause.
- Two-star yellow alarms and yellow INOPs: a lower pitched tone is repeated three times, followed by a configurable pause.
- One-star yellow alarms (short yellow alarms): the audible indicator is the same as for yellow alarms, but of shorter duration.
- Standard INOPs: a lower pitched tone is repeated twice, followed by a pause.

Changing the Alarm Tone Volume

The alarm volume symbol at the top right of the monitor screen gives you an indication of the current volume. To change the volume:

- 1 Select the volume symbol.



- 2 The volume scale pops up.



Alarm Volume

- 3 Select the required volume from the volume scale.

When the alarm volume is set to zero (**Off**), the alarm volume symbol shows this symbol:



If you switch the alarm volume off, you will not get any audible indication of alarm conditions.

Minimum Volume for Severe Yellow or Red INOPs

Severe yellow or red INOPs require action to ensure the well-being of the patient. Therefore the minimum volume for the INOP tone is set to at least alarm volume 8, irrespective of the current alarm volume setting. The INOP tone will sound even if the monitor alarm volume is set to zero.

The severe INOPs for which this applies are:

INOP	Factory Default	INOP tone
Battery Empty	cyan	severe (fetal monitor with battery option) hard (CL Pods) soft (CL transducers)
ECG Leads Off	cyan	hard
NBP Cuff Overpress	cyan	severe
Cuff Not Deflated	cyan	severe
!! <SpO₂ Label> No Pulse	cyan	hard

Power Loss Tone

FM20/30 When power is lost - no power is available from the AC power source or from the battery - a beeper will sound. The tone can be silenced by pressing the On/Standby button.
with
Battery
Option,
FM40/50

Acknowledging Alarms

To acknowledge all active alarms and INOPs, select the **Silence** key. This switches off the audible alarm indicators.



A check mark beside the alarm message indicates that the alarm has been acknowledged.

If the condition that triggered the alarm is still present after the alarm has been acknowledged, the alarm message stays on the screen with a check mark symbol beside it.

If the alarm condition is no longer present, all alarm indicators stop and the alarm is reset.

Switching off the alarms for the measurement in alarm, or switching off the measurement itself, also stops alarm indication.

Acknowledging Disconnect INOPs

Acknowledging an INOP that results from a disconnected transducer switches off the associated measurement.

Alarm Reminder

If **Alarm Reminder** is configured on for your monitor, you will get an audible reminder of alarm conditions that remain active after you have acknowledged the alarm. This reminder may take the form of a repetition of the alarm tone for a limited time, or an unlimited repetition of the alarm tone (this is the same as a new alarm). **Alarm Reminder** is available only for yellow and red INOPs and not for standard light blue INOPs.

In Configuration Mode, you can set the interval between silencing the alarm and sounding the reminder tone to one, two, or three minutes.

Pausing or Switching Off Alarms

If you want to temporarily prevent alarms from sounding, for example while you are moving a patient, you can pause alarms. Depending on your fetal monitor configuration, alarms are paused for one, two, or three minutes, or infinitely.

To view the alarm pause setting chosen for your unit:

- 1 Select **Main Setup, Alarms, Alarm Settings.**
- 2 Check the **Alarms Off** setting.

This setting can be changed in Configuration Mode.

There are some settings made in Configuration Mode that can affect the availability of the pause alarms functionality.

- The **Pause Alarms / Alarms Off** permanent key can be removed from the screen to avoid unintentional switching off of alarms. In this case you can only pause alarms or switch alarms off permanently in the **Alarms** menu, under **Main Setup**.
- The **Pause Alarms / Alarms Off** permanent key and the corresponding hardkey can be configured to pause or switch off red and yellow alarms, yellow alarms only, or not to function at all. If they are configured not to function, you cannot pause alarms or switch alarms off permanently at all.

When the alarms off priority is set to **Yellow Only**, the **Pause Alarms / Alarms Off** key becomes the **Pause Yellow / Yellow Al. Off**.

To Pause All Alarms

If you have configured alarms to be paused for one, two, or three minutes, the SmartKey is labeled **Pause Alarms**.

Select the **Pause Alarms** SmartKey to pause all alarms.



Or

- 1 Select **Main Setup.**
- 2 Select **Alarms.**
- 3 Select **Pause Alarms.**

To Switch All Alarms Off

You can switch alarms off permanently if your monitor is configured to allow infinite alarms pause, and the SmartKey is labeled **Alarms Off**.

Select the **Alarms Off** SmartKey.



Or

- 1 Select **Main Setup.**
- 2 Select **Alarms.**

- 3 Select **Alarms Off**.

To Switch Individual Measurement Alarms On or Off

This applies to alarm mode **All**.



- 1 Select the measurement numeric to enter its setup menu.
- 2 Select **Alarms** to toggle between **On** and **Off**.

The alarms off symbol is shown beside the measurement numeric.

While Alarms are Paused or Off

When red alarms are paused or off:

The red Alarms Paused lamp on the monitor front panel is lit.


In the alarm field, the monitor displays the message **Al. Paused x:yy** or **Alarms Off**, together with the alarms paused symbol or the alarms off symbol.

When yellow alarms are paused or off:


The yellow Alarms Paused lamp on the monitor front panel is lit.

In the alarm field, the monitor displays the message **YellowPaused x:yy** or **Yellow Alarms Off**, together with the alarms paused symbol or the alarms off symbol.

For red and yellow alarms:



Alarms Paused



Alarms Off

- No alarms are sounded and no alarm messages are shown.
- INOP messages are shown but no INOP tones are sounded.

The following INOPs are the only exceptions:

NBP Cuff Overpress, **Batt Empty** and **Batt Malfunction** (these INOPs are issued even if alarms are paused or off).

These INOPs switch the alarms on, and the INOP tones are sounded, even if alarms are paused or off. You need to remove the INOP condition first before you can switch the alarm tones off again.

Restarting Paused Alarms

To manually switch on alarm indication again after a pause, select the SmartKey **Pause Alarms** (or **Alarms Off**) again.

Alarm indication starts again automatically after the pause period expires. If the monitor is configured to stay paused infinitely, you must select **Alarms Off** again to restart alarm indication.

Extending the Alarm Pause Time

If your monitor has the function extended alarm pause enabled, you can extend the alarm pause time. Use this to prevent alarms being indicated, for example, while you are washing a patient or carrying out a procedure. Only extend the alarm pause time when you are sure that clinical personnel are available to monitor the patient's condition closely.

To extend the alarm pause time to five or 10 minutes,

- 1 Select one of the alarm fields. This calls up the **Alarm Messages** window.
- 2 Select either the pop-up key **PauseAl. 5 min** or the pop-up key **PauseAl. 10 min**. Each time you select one of these pop-up keys, the Alarm Pause Time is reset to five (or 10) minutes.

Alarm Limits

The alarm limits you set determine the conditions that trigger yellow and red limit alarms.

WARNING

Be aware that the monitors in your care area may each have different alarm settings, to suit different scenarios. Always check that the alarm settings are appropriate before you start monitoring.

Viewing Individual Alarm Limits (Alarm Mode "All" Only)



Alarm Limits

Audio source symbol

You can usually see the alarm limits set for each measurement next to the measurement numeric on the main screen.

If your monitor is not configured to show the alarm limits next to the numeric, you can see them in the appropriate measurement setup menu. Select the measurement numeric to enter the menu and check the limits.

Changing Alarm Limits

To change individual measurement alarm limits using the measurement's setup menu:

- 1 In the measurement's setup menu, select the alarm limit you want to change. This calls up a list of available values for the alarm limit.
- 2 Select a value from the list to adjust the alarm limit.

Reviewing Alarms

To review the currently active alarms and INOPs, select any of the alarm status areas on the fetal monitor screen. The **Alarm Messages** window pops up. All alarms and INOPs are erased from the monitor's alarm history when you discharge a patient, or if you enter Demonstration Mode.

Alarm Messages Window

The **Alarm Messages** window shows all the currently active alarms and INOPs in chronological order, beginning at the top with the most recent. INOPs are shown on the left hand side, and alarms are shown on the right hand side. Any active red alarms are shown first, followed by yellow alarms. Acknowledged alarms or INOPs are shown with the check mark symbol.

The **Alarm Messages** window pop-up keys appear when the window is opened. Selecting the **Review Alarms** pop-up key opens the **Review Alarms** window.

Review Alarms Window

Review Alarms		
5 Apr 16:55:18	Alarms On	1
5 Apr 16:45:15	SpO2 NON-PULSAT.	2
5 Apr 16:44:57	Alarms Silenced	3
5 Apr 16:44:46	** FHR1 HIGH (153>150)	4
5 Apr 16:44:39	** SpO2 LOW (95<99)	5

The **Review Alarms** window contains a list of up to 300 of the most recent alarms and INOPs with date and time information. If configured to do so, each alarm is shown with the alarm limit active when the alarm was triggered and the maximum value measured beyond this limit. The **Review Alarms** window also shows any changes made to the Alarms On/Off or Silence status. Only the main alarms On/Off transitions are logged in the alarm history, and On/Off alarm transitions for individual measurements are not logged.

The information in the **Review Alarms** window is deleted when a patient is discharged.

The **Review Alarms** window pop-up keys appear when the window is opened. Selecting the **Active Alarms** pop-up key opens the **Alarm Messages** window.

Latching Alarms

The alarm latching setting for your monitor defines how the alarm indicators behave when you do not acknowledge them. When alarms are set to non-latching, their indicators end when the alarm condition ends. Switching alarm latching on means that visual and/or audible alarm indications are still displayed or announced by the monitor after the alarm condition ends. The indication lasts until you acknowledge the alarm.

Viewing the Alarm Latching Settings

To see the alarm latching setting for your monitor:

- 1 In the monitor's **Main Setup** menu, select **Alarms**.
- 2 Select **Alarm Settings**, and see the **Visual Latching** and **Audible Latching** settings.

This setting can be changed in Configuration Mode. You should be aware of the settings chosen for your unit. There are three possible choices each for visual and audible latching: **Red**, **Red & Yellow**, and **Off**. The audible latching configuration can never be configured to a higher level than that configured for the visual latching. In other words, the audible latching setting is always the same level, or lower, than the visual latching setting. For example, if visual latching is configured to **Red Only**, then audible latching can only be set to **Red** or **Off**. The following table shows the possible combinations for latching settings:

Possible Combinations for Alarm Latching Settings	
Visual Latching Setting	Audible Latching Setting
Red & Yellow	Red & Yellow
Red & Yellow	Red
Red & Yellow	Off
Red	Red
Red	Off
Off	Off

Alarm Latching Behavior

Alarm Condition		Red and Yellow Measurement Alarms		
Acknowledgment	Presence	Non-latching alarms	Visual and audible latching	Visual latching, audible non-latching
Alarm has not been acknowledged	Alarm condition still present	Alarm tone on Alarm message	Alarm tone on Alarm message Flashing numerics	Alarm tone on Alarm message Flashing numerics
	Alarm condition no longer present	All audible and visual alarm indicators automatically stop	Alarm tone on Alarm message Flashing numerics	Alarm message Flashing numerics Audible alarm indicators automatically stop
Alarm has been acknowledged	Alarm condition still present	Alarm tone off. Alarm message. Audible alarm reminder (if configured)	Alarm tone off Alarm message Flashing numerics Audible alarm reminder (if configured)	Alarm tone off Alarm message Flashing numerics Audible alarm reminder (if configured)
	Alarm condition no longer present	Audible and visual alarm indicators automatically stop	Audible and visual alarm indicators automatically stop	Audible and visual alarm indicators automatically stop

All INOPs except the "unplugged" INOPs are non-latching.

Testing Alarms

In general, to test the functioning of visible and audible alarms, do the following:

- 1 Enable the alarm.
- 2 Set the alarm limits.
- 3 Measure or simulate the parameter that is out of range, or signal loss.
- 4 Verify that the visible and audible alarms are working.

As an example, to test the FHR alarms:

- 1 Connect the US transducer to a fetal sensor socket.
- 2 Enable the FHR alarming (see "Turning Alarms On or Off" on page 185).

Set the high alarm limit and delay to 150 bpm and 60 seconds respectively, and the low alarm limit and delay to 110 bpm and 60 seconds respectively (see "Changing Alarm Limits" on page 185).

- 1 Generate a fetal heart rate of approximately 180 bpm (3 beats per second) for more than one minute.
- 2 Verify the functioning of the visible and audible alarm.

Alarm Behavior at Power On

If the monitor is switched off for longer than one minute and then switched on again, or after a loss of power lasting longer than one minute, the monitor can be configured to restore either the alarm settings from the monitor's User Defaults, or the most recently used alarm settings. When a patient is discharged, the User Defaults of the monitor are always restored independent of the setting. After any of these situations, you should check that the alarm settings are appropriate for your patient. If power is lost for less than one minute, the alarm settings prior to the power loss are restored.

If **AlarmsOffAtStart** is enabled in Configuration Mode alarms will be initially paused or off the next time the monitor is switched on (depending on the configuration setting for **Alarms Off**). Even if it is enabled, this setting only takes effect if the power down time is more than one minute, and the Global Setting **Automat Default** is set to **Yes**.

In order for alarms to be paused or switched off initially:

- the monitor must be switched off for more than one minute, and
- the main alarms were switched off or paused before the monitor was switched off.

Patient Alarms and INOPs

This chapter lists alarms and technical alarms (INOPs) for the fetal monitors irrespective of their priority.

Alarm Messages

Fetal alarms are identified by either "FHR" or "DFHR". All other alarms without these identifiers refer to maternal parameters.

Alarm Message	From	Condition	Indication
*** Brady (Pulse) ***Brady xxx < yyy	SpO ₂	The heart rate from the Pulse signal has fallen below the bradycardia limit. xxx denotes the lowest measured value; yyy is the bradycardia limit.	Numeric flashes and alarm limit is highlighted, red alarm message, alarm tone.
*** Desat *** Desat xx < yy	SpO ₂	The SpO ₂ value has fallen below the desaturation alarm limit. xx denotes the lowest measured value, and yy is the desaturation limit.	Numeric flashes, red alarm message, alarm tone.
** DFHR1 High ** DFHR2 High ** DFHR3 High **DFHR1 xxx>yyy **DFHR2 xxx>yyy **DFHR3 xxx>yyy	FHR (DECG)	The fetal heart rate obtained from DECG has risen above the high alarm limit. xxx denotes the highest measured value, and yyy is the high alarm limit.	Numeric flashes and high limit is highlighted, yellow alarm message, alarm tone.
** DFHR1 Low ** DFHR2 Low ** DFHR3 Low **DFHR1 xxx<yyy **DFHR2 xxx<yyy **DFHR3 xxx<yyy	FHR (DECG)	The fetal heart rate obtained from DECG has fallen below the low alarm limit. xxx denotes the lowest measured value, and yyy is the low alarm limit.	Numeric flashes and low limit is highlighted, yellow alarm message, alarm tone.
*** Extreme Brady ***Brady xxx < yyy	MECG	The maternal heart rate obtained from the maternal ECG has fallen below the extreme bradycardia limit. xxx denotes the lowest measured value, and yyy is the extreme bradycardia limit.	Numeric flashes, red alarm message, alarm tone.

6 Patient Alarms and INOPs

Alarm Message	From	Condition	Indication
*** Extreme Tachy ***Tachy xxx > yyy	MECG	The maternal heart rate obtained from the maternal ECG has risen above the extreme tachycardia limit. xxx denotes the highest measured value, and yy is the extreme tachycardia limit.	Numeric flashes, red alarm message, alarm tone.
** FHR1 High ** FHR2 High ** FHR3 High **FHR1 xxx>yyy **FHR2 xxx>yyy **FHR3 xxx>yyy	FHR (ultrasound)	The fetal heart rate obtained from ultrasound has risen above the high alarm limit. xxx denotes the highest measured value, and yyy is the high alarm limit.	Numeric flashes and high limit is highlighted, yellow alarm message, alarm tone.
** FHR1 Low ** FHR2 Low ** FHR3 Low **FHR1 xxx<yyy **FHR2 xxx<yyy **FHR3 xxx<yyy	FHR (ultrasound)	The fetal heart rate obtained from ultrasound has fallen below the low alarm limit. xxx denotes the lowest measured value, and yyy is the low alarm limit.	Numeric flashes and low limit is highlighted, yellow alarm message, alarm tone.
** HR High **HR xxx>yyy	MECG	The maternal heart rate obtained from the maternal ECG has risen above the high alarm limit. xxx denotes the highest measured value, and yyy is the high alarm limit.	Numeric flashes and high limit is highlighted, yellow alarm message, alarm tone.
** HR Low **HR xxx<yyy	MECG	The maternal heart rate obtained from the maternal ECG has fallen below the low alarm limit. xxx denotes the lowest measured value, and yyy is the low alarm limit.	Numeric flashes and low limit is highlighted, yellow alarm message, alarm tone.
** NBPs High ** NBPd High ** NBPm High	Noninvasive blood pressure	The measured noninvasive blood pressure value is above the high alarm limit. s , d , or m after the label s indicates whether the systolic, diastolic or mean pressure has crossed the limit.	Numeric flashes and high limit is highlighted, yellow alarm message, alarm tone.
** NBPs Low ** NBPd Low ** NBPm Low	Noninvasive blood pressure	The measured noninvasive blood pressure value is below the low alarm limit. s , d , or m after the label s indicates whether the systolic, diastolic or mean pressure has crossed the limit.	Numeric flashes and low limit is highlighted, yellow alarm message, alarm tone.
** Pulse High	SpO ₂	The pulse rate has exceeded the high alarm limit.	Numeric flashes and high limit is highlighted, yellow alarm message, alarm tone.
** Pulse Low	SpO ₂	The pulse rate has dropped below the low alarm limit.	Numeric flashes and low limit is highlighted, yellow alarm message, alarm tone.

Alarm Message	From	Condition	Indication
** SpO₂ High	SpO ₂	The arterial oxygen saturation has exceeded the high alarm limit.	Numeric flashes and high limit is highlighted, yellow alarm message, alarm tone.
** SpO₂ Low	SpO ₂	The arterial oxygen saturation has fallen below the low alarm limit.	Numeric flashes and low limit is highlighted, yellow alarm message, alarm tone.
*** Tachy (Pulse) ***Tachy/P xxx>yyy	SpO ₂	The heart rate from the Pulse signal has exceeded the tachycardia limit. xxx denotes the highest measured value, and yyy is the tachycardia limit.	Numeric flashes, alarm limit is highlighted, red alarm message, alarm tone.
** <iTemp Label> High	iTemp	The temperature value has exceeded the high alarm limit.	Numeric flashes and high limit is high-lighted, yellow alarm lamp, alarm tone
** <iTemp Label> Low	iTemp	The temperature value has dropped below the alarm limit.	Numeric flashes and high limit is high-lighted, yellow alarm lamp, alarm tone

Technical Alarm Messages (INOPs)

Monitor INOPs

INOP Message	Indication	What to do
Bus Master Malfunc	Monitor: INOP tone	There is a problem with the monitor's hardware. Contact your service personnel.
Bus Master Incomp	Monitor: INOP tone	The Bus Master board or firmware is incompatible and cannot be used with this monitor software revision. Contact your service personnel.
Check Flex Texts	Monitor: INOP tone	If this INOP appears, check the monitor and patient settings before you resume monitoring. If the settings are unexpected, there may be a problem with the monitor software. Contact your service personnel.
Check Monitor Func	Monitor: INOP tone	A potential internal problem with the monitor has been detected. Contact your service personnel.
Check Keyboard	Monitor: INOP tone	Perform a visual and functional check of the keyboard. Contact your service personnel.
Check Mouse	Monitor: INOP tone	Perform a visual and functional check of the mouse input device. Contact your service personnel.
Check OBR Config	Monitor: INOP tone	OBR (OB Radio) channel configuration is invalid or settings are not present. Check channel configuration and hardware option.
Check Settings	Monitor: INOP tone	If this INOP appears, check the monitor and patient settings before you resume monitoring. If the settings are unexpected, there may be a problem with the monitor software. Contact your service personnel.
Check TI Config	Monitor: INOP tone	If this INOP appears, the monitor can be used normally, except for the Trace Interpretation feature. Contact your service personnel.
Check Touch Input	Monitor: Display only	Perform a visual and functional check of the touch input device. Contact your service personnel.
Coincidence	Monitor: INOP tone	One or more fetal/maternal heart rates/pulse rates persistently coincide with each other (see "Cross-Channel Verification (CCV)" on page 149). The INOP tone has a configurable delay. The default setting is an instant INOP tone.
Internal.Comm.Malf	Monitor: INOP tone	There is a problem with I2C Bus communication in the monitor. Contact your service personnel.
No Central Monit.	Monitor: INOP tone	There is a problem with the communication to the network. Central monitoring is currently not possible (no patient alarms or information). Check the connection. Contact your service personnel.
Paper End	Monitor: INOP tone	The end of the paper pack is detected. Insert a new pack of paper.

INOP Message	Indication	What to do
Settings Malfunc	Monitor: INOP tone	The monitor cannot use the predefined settings for monitoring. Contact your service personnel.
Speaker Malfunc	Monitor: INOP tone	Contact your service personnel to check the speaker and the connection to the speaker.
SRR Incompatible	Monitor: INOP tone	The connected SRR (Short Range Radio) board or FW cannot be used with this base station SW revision. Contact your service personnel.
SRR Malfunction	Monitor: INOP tone	Malfunction in the SRR (Short Range Radio) module(s). If the INOP persists, contact your service personnel.
SRR Invalid Chan	Monitor: INOP tone	The channel configuration of the Short Range Radio is invalid. Check channel and channel mask configuration.
SRR Interference	Monitor: INOP tone	The short range radio connection has interference from another device. Try using another channel.
Tele Malfunction	Monitor: INOP tone	Malfunction of the connected base station. If the INOP persists contact your service personnel.
Tele Unsupported	Monitor: INOP tone	The base station is not supported by the monitor. Check software revisions.
TimeExpired:NST	Monitor: INOP tone (Message and tone only if Setup NST Timer, Notification is set to Alarm .)	The time has expired for the NST timer. Clearing the timer clears the INOP.
Unsupported LAN	Monitor: INOP tone	There is a problem with the communication to the network and central monitoring is currently not possible. Check the connection. If the INOP persists, switch off the monitor and contact your service personnel.
User I/F Malfunc	Monitor: INOP tone	Perform a visual and functional check of all the monitor input devices. Contact your service personnel.

FM20/FM30 Battery INOPs

INOP Message	Indication	What to do
Batt Empty	Monitor: INOP tone Battery: battery LED flashes during this INOP, alarms cannot be paused or switched off	The estimated remaining battery-powered operating time is less than 10 minutes. Connect the monitor to AC immediately. If the condition persists and the monitor is not connected to mains power, this INOP is re-issued two minutes after you acknowledge it.
Batt Incompat	Monitor: INOP tone	The battery cannot be used with this monitor. Replace it with the correct battery.
Batt Low	Monitor: INOP tone	The estimated battery-powered operating time remaining is less than 20 minutes.
Batt Malfunction	Monitor: INOP tone Battery: battery LED flashes during this INOP, alarms cannot be paused or switched off unless the monitor is connected to mains power.	The monitor cannot determine the battery status. If the INOP persists, replace the faulty battery. If the condition persists and the monitor is not connected to mains power, this INOP is reissued two minutes after you acknowledge it.
Batt Missing	Monitor: INOP tone	No battery found in FM20/30. Contact your service personnel.
Charger Malfunc	Monitor: INOP tone Battery: battery LED may flash	There is a problem with the battery charger in the monitor. Connect the monitor to mains power and contact your service personnel.
Check Batt Temp	Monitor: INOP tone	The temperature of the battery is too high. Check that the monitor is not exposed to heat.

Recorder INOPs

INOP Message	Indication	What to do
Check Paper	Monitor: display only Recorder: Print-out	Check that there is no paper jam, that the print drawer is properly shut, that the paper is loaded with the grid facing upwards, and that the correct Philips paper is being used.
FetRec Equip Malf	Monitor: INOP tone	There is a problem with the fetal recorder hardware. Contact your service personnel.
FetRec Chk Config	Monitor: INOP tone	Recorder Speed and/or Scale Type settings are set to "Unknown" and need to be set to the correct values in Configuration Mode before the recorder can be operated.
Printhead Overheat	Monitor: INOP tone	The printhead is too hot. The recorder stops, the recorder Start/ Stop key is disabled, and remains so until the printhead cools down sufficiently. Wait for the printhead to cool down, then press the recorder Start/ Stop key or the Silence key to clear the INOP.

INOP Message	Indication	What to do
Wrong Paper Scale	Monitor: INOP tone Recorder: Print-out	The grid scale of the paper in the monitor does not match the grid scale configured in the monitor. Make sure that you use the correct paper and scale for your institution: pre-printed: 30-240 in US and Canada, 50-210 in other geographies.

Wired Transducer INOPs

FHR (Ultrasound)

INOP Message	Indication	What to do
FHR1 Equip Malf FHR2 Equip Malf FHR3 Equip Malf	Monitor: INOP tone	There is a problem with the FHR hardware. Contact your service personnel.
FHR1 Signal Loss FHR2 Signal Loss FHR3 Signal Loss	Monitor: INOP tone INOP tone if Alarm Mode is set to ALL. No tone if Alarm Mode is set to INOP only.	The input signal quality is not sufficient to process the measurement. Adjust the position of the transducer to obtain a better signal.
FHR1 Unplugged FHR2 Unplugged FHR3 Unplugged	Monitor: INOP tone	Reconnect the FHR transducer to the monitor. Check all connections are sound.

IUP

INOP Message	Indication	What to do
IUP Equip Malf	Monitor: INOP tone	There is a problem with the IUP hardware. Contact your service personnel.
IUP Unplugged	Monitor: INOP tone	Reconnect the IUP transducer to the monitor. Check all connections are sound.

Toco

INOP Message	Indication	What to do
Toco Equip Malf	Monitor: INOP tone	There is a problem with the Toco hardware. Contact your service personnel.
Toco Unplugged	Monitor: INOP tone	Reconnect the Toco transducer to the monitor. Check all connections are sound.
Pulse(Toco) Malf	Monitor: INOP tone Toco MP Maternal Pulse Measurement	There is a problem with the Toco MP transducer hardware. Contact your service personnel.

Cableless Transducer INOPs

INOP Message	Indication	What to do
cl US Disconnect cl Toco Disconnect cl ECG/IUP Disconn	Monitor: INOP tone CL transducer: tone and cyan LED	The cableless transducer has lost the connection to the monitor for more than 5 sec.
cl US Malfunction cl Toco Malf cl ECG/IUP Malf	Monitor: INOP tone Numeric is replaced by a -?- CL transducer: tone and cyan LED	Malfunction of the transducer
cl US Batt Low cl Toco Batt Low cl ECG/IUP BattLow	Monitor: Display only CL transducer: LED flashes red	The remaining monitoring time with this transducer is below 30 minutes. Charge battery.
cl US Batt Empty !!cl US Batt Empty !!!cl US BattEmpty cl Toco Batt Empty !!clToco BattEmpty !!!clTocoBattEmpty cl ECG/IUP Empty !!cl ECG/IUP Empty !!!clECG/IUP Empty	Monitor: Display only CL transducer: LED flashes red	The remaining monitoring time with this transducer is below 15 minutes. Charge battery.
cl US Batt Malf cl Toco Batt Malf clECG/IUP BattMalf	Monitor: INOP tone CL transducer: tone and cyan LED	Malfunction of the battery system (charger circuit or battery) detected (e.g. persistent communication error, over-voltage, over-current, battery incompatible). Contact your service personnel.
cl US Check Temp cl Toco Chk Temp clECG/IUP Chk Temp	Monitor: INOP tone CL transducer: tone and cyan LED	The temperature of the transducer or its battery is critically high. If the transducer is docked on the base station, charging is interrupted. Check that device is not covered and do not expose device to heat sources. Do not place the base station in a sunny window in a room without active air conditioning. The environment temperature to charge a transducer battery should not exceed 40 °C. If the INOP persists, remove device from patient and contact your service personnel.

DECG and MECG INOPs

INOP Message	Indication	What to do
DFHR1 Equip Malf DFHR2 Equip Malf DFHR3 Equip Malf	Monitor: INOP tone	There is a problem with the DECG hardware. Contact your service personnel.
DFHR1 Leads Off DFHR2 Leads Off DFHR3 Leads Off	Monitor: INOP tone Numeric is replaced by a -?-	One or more DECG lead is not attached. Make sure that all required leads are attached, and no electrodes have been displaced. Check all connections are sound, and that the legplate attachment electrode is properly attached. If the INOP persists, try using another adapter cable, or legplate attachment electrode. If the INOP persists, contact your service personnel.
DFHR1 Signal Loss DFHR2 Signal Loss DFHR3 Signal Loss	Monitor: INOP tone Numeric is replaced by a -?- INOP tone if Alarm us set to ALL. no tone if Alarm Mode is set to INOP only.	The input signal quality is not sufficient to process the measurement. Reapply the fetal scalp electrode.
DFHR1 Unplugged DFHR2 Unplugged DFHR3 Unplugged	Monitor: INOP tone Numeric is replaced by a -?-	Reconnect the DECG transducer to the monitor. Check all connections are sound.
MECG Equip Malf	Monitor: INOP tone	There is a problem with the MECG hardware. Contact your service personnel.
MECG Leads Off	Monitor: INOP tone Numeric is replaced by a -?-	One or more MECG lead is not attached. Ensure that all required leads are attached, and no electrodes have been displaced. Check all connections are sound. If the INOP persists, try using another adapter cable. If the INOP still persists, contact your service personnel.
MECG Unplugged	Monitor: INOP tone Numeric is replaced by a -?-	Reconnect the MECG transducer to the monitor. Check all connections are sound.

NBP INOPs

INOP Message	Indication	What to do
Cuff Not Deflated	Monitor: INOP tone (severe yellow/red INOP tone) Numeric is replaced by a -?- During this INOP, alarms cannot be paused or switched off.	Remove the cuff from the patient. Make sure that the tubing is not kinked or twisted. Try restarting the measurement. You can silence the INOP, but the INOP message remains visible until the next measurement is started or the Stop All SmartKey is selected.

6 Patient Alarms and INOPs

INOP Message	Indication	What to do
NBP Cuff Overpress	Monitor: INOP tone (severe yellow/red INOP tone) Numeric is replaced by a -?- During this INOP, alarms cannot be paused or switched off.	The cuff pressure exceeds the overpressure safety limits. Remove the cuff from the patient. Ensure that the tubing is not kinked or twisted and that the correct patient category is selected. Try restarting the measurement. You can silence this INOP, but the INOP message remains visible until the next measurement is started or the Stop All SmartKey is selected.
NBP Equip Malf	Monitor: INOP tone Numeric is replaced by a -?-	Remove the cuff from the patient. The noninvasive blood pressure hardware is faulty. Contact your service personnel. You can silence this INOP, but the INOP message remains visible until the next measurement is started or the Stop All SmartKey is selected.
NBP Interrupted	Monitor: INOP tone Numeric is replaced by a -?-	Check the tubing and cuff for leakages or kinks. Check that you are using the correct cuff size and placement, and that the correct patient category is selected. Try restarting the measurement. If the INOP occurs repeatedly, contact your service personnel. You can silence this INOP, but the INOP message remains visible until the next measurement is started or the Stop All SmartKey is selected. This INOP arises when the measurement needed longer than the maximum time for inflation, deflation or the total measurement.
NBP Measure Failed	Monitor: INOP tone Numeric is replaced by a -?-	Check that you are using the correct cuff size and placement, and that the correct patient category is selected. Try restarting the measurement. You can silence this INOP, but the INOP message remains visible until the next measurement is started or the Stop All SmartKey is selected. Check the condition and suitability of the patient for noninvasive blood pressure monitoring. Use another cuff to continue measuring.

CL NBP Pod INOPs

INOP Message	Indication	What to do
cl NBP Disconnect	Monitor: INOP tone CL NBP Pod: INOP tone and display message	The CL NBP Pod has lost the SRR connection to the monitor.
cl NBP No Cradle	Monitor: INOP tone CL NBP Pod: display message	The CL NBP Pod is not in the cradle. You can silence this INOP, but the INOP message remains visible until the CL NBP Pod is inserted into the cradle and the next measurement is started or the Stop All SmartKey is selected.
cl NBP Remove	Monitor: INOP tone CL NBP Pod: display message	The temperature of the battery in the CL NBP Pod is too high. Remove the Cableless Measurement Device from the patient and contact service personnel.
cl NBP Batt Empty !!cl NBP BattEmpty !!!cl NBP BatEmpty	Monitor: INOP tone CL NBP Pod: display message	The remaining battery time of the CL NBP Pod is below 30 minutes. Charge the battery.
cl NBP Batt Incomp	Monitor: INOP tone CL NBP Pod: display message	The battery in use with the CL NBP Pod is incompatible. Replace it with one approved for use with the CL NBP Pod.
cl NBP Batt Low	Monitor: INOP tone CL NBP POD: display message	The remaining battery time of the CL NBP Pod is below 2 hours.
cl NBP Batt Malf	Monitor: INOP tone CL NBP Pod: display message	There is a malfunction in the CL NBP Pod's battery system. Contact your service personnel.
cl NBP Batt Temp	Monitor: INOP tone CL NBP Pod: display message	The temperature of the battery in the CL NBP Pod is critically high. Check that the Pod is not covered or exposed to a heat source. If the INOP persists, remove the Pod from the patient and contact your service personnel.
cl NBP Check Batt	Monitor: INOP tone CL NBP Pod: display message	The battery in the CL NBP Pod is nearing the end of its useful life. Only 50 charge/discharge cycles remain. Contact your service personnel to replace the battery.
cl NBP Serv Batt	Monitor: INOP tone CL NBP Pod: display message	The battery in the CL NBP Pod has reached the end of its useful life. It can no longer be charged. Contact your service personnel to replace the battery.

SpO₂ INOPs

INOP Message	Indication	What to do
SpO₂ Equip Malf	Monitor: INOP tone	There is a problem with the SpO ₂ hardware. Contact your service personnel.
SpO₂ Erratic	Monitor: INOP tone Numeric is replaced by a -?-	Check the sensor placement. Try another adapter cable and sensor. If the INOP persists, contact your service personnel.
SpO₂ Extd.Update	Monitor: display only Numeric is replaced by a -?-	The update period of displayed values is extended due to a noninvasive blood pressure measurement on the same limb or an excessively noisy signal.
SpO₂Interference	Monitor: INOP tone Numeric is replaced by a -?-	There is too much interference, caused by a high level of ambient light and/or electrical interference. Cover the sensor to minimize ambient light. If the INOP persists, ensure that the sensor cable is not damaged or positioned too close to power cables.
SpO₂ Low Perf	Monitor: display only Numeric is replaced by a -?-	Accuracy may be compromised due to very low perfusion. Stimulate circulation at sensor site. If INOP persists, change the measurement site.
SpO₂ NoisySignal	Monitor: INOP tone Numeric is replaced by a -?-	Excessive patient movement or electrical interference is causing irregular pulse patterns. Try to reduce patient movement or to relieve the cable strain on the sensor.
SpO₂ No Pulse	Monitor: INOP tone Numeric is replaced by a -?-	Check the perfusion at measurement site. If necessary, stimulate circulation or change measurement site. If the INOP is due to noninvasive blood pressure measurement on the same limb, wait until the measurement is finished.
SpO₂ No Sensor	Monitor: INOP tone Numeric is replaced by a -?-	Ensure the SpO ₂ sensor is connected. If the INOP persists, try another adapter cable and sensor. If you silence this INOP, the measurement will be switched off.
SpO₂ Poor Signal	Monitor: display only Numeric is replaced by a -?-	The signal condition of the SpO ₂ measurement is poor and measurement accuracy may be compromised.
SpO₂ Pulse?	Monitor: INOP tone Numeric is replaced by a -?-	The detectable pulsations of the SpO ₂ signal are outside the specified pulse rate range.
SpO₂ Searching	Monitor: display only Numeric unavailable	SpO ₂ is analyzing the patient signal to derive Pulse, and SpO ₂ values. Please wait until the search analysis is complete.
SpO₂ Sensor Malf	Monitor: INOP tone Numeric unavailable	The SpO ₂ sensor or adapter cable is faulty. Try another adapter cable and sensor. If the INOP persists, contact your service personnel.
SpO₂ Sensor Off	Monitor: INOP tone Numeric unavailable	The SpO ₂ sensor is not properly applied to the patient. Apply the sensor following the instructions supplied by the manufacturer.
SpO₂ Unkn.Sensor	Monitor: display only Numeric is replaced by a -?-	The connected sensor or adapter cable is not supported by the SpO ₂ measurement. Use only specified sensors and cables.

INOP Message	Indication	What to do
SpO₂ Upgrade	Monitor: display only Numeric unavailable	The SpO ₂ measurement is currently in UPGRADE mode. Monitoring is not possible in this mode.

CL SpO₂ Pod INOPs

INOP Message	Indication	What to do
cl SpO₂ Batt Empty !!cl SpO₂ BatEmpty !!!cl SpO₂ BtEmpty	Monitor: INOP tone CL SpO ₂ Pod: display message	The remaining battery time of the CL SpO ₂ Pod is below 30 minutes. Charge the battery.
cl SpO₂ Batt Incmp	Monitor: INOP tone CL SpO ₂ Pod: display message	The battery in use with the CL SpO ₂ Pod is incompatible. Replace it with one approved for use with the CL SpO ₂ Pod.
cl SpO₂ Batt Low	Monitor: INOP tone CL SpO ₂ Pod: display message	The remaining battery time of the CL SpO ₂ Pod is below 2 hours.
cl SpO₂ Batt Malf	Monitor: INOP tone CL SpO ₂ Pod: display message	There is a malfunction in the CL SpO ₂ Pod's battery system. Contact your service personnel.
cl SpO₂ Batt Temp	Monitor: INOP tone CL SpO ₂ Pod: display message	The temperature of the battery in the CL SpO ₂ Pod is critically high. Check that the Pod is not covered or exposed to a heat source. If the INOP persists, remove the Pod from the patient and contact your service personnel.
cl SpO₂ Check Batt	Monitor: INOP tone CL SpO ₂ Pod: display message	The battery in the CL SpO ₂ Pod is nearing the end of its useful life. Only 50 charge/discharge cycles remain. Contact your service personnel to replace the battery.
cl SpO₂ Serv Batt	Monitor: INOP tone CL SpO ₂ Pod: display message	The battery in the CL SpO ₂ Pod has reached the end of its useful life. It can no longer be charged. Contact your service personnel to replace the battery.
cl SpO₂ Disconnect	Monitor: INOP tone CL SpO ₂ Pod: INOP tone and display message	The CL SpO ₂ Pod has lost the SRR connection to the monitor.
cl SpO₂ No Cradle	Monitor: INOP tone CL SpO ₂ Pod: display message	The CL SpO ₂ Pod is not in its cradle.
cl SpO₂ Remove	Monitor: INOP tone CL SpO ₂ Pod: display message	The temperature of the battery in the SpO ₂ Pod is too high. Remove the SpO ₂ Pod from the patient and contact service personnel.

Tympanic Temperature INOPs

INOP Message	Indication	What to do
<iTemp Label> Malfunction	Monitor: INOP tone Numeric is replaced by a -?-	The thermometer or the interface board are defective. Check module on another monitor. Connect another module to this monitor and replace the interface board. If the problem persists, contact your service personnel.
<iTemp Label> Incompat.	Monitor: INOP tone Numeric is replaced by a -?-	The iTemp device firmware is not supported by the monitor software.
<iTemp Label> Meas Failed	Monitor: INOP tone Numeric is replaced by a -?-	Ambient temperature is out of range. If ambient temperature is within valid range, check the thermometer. Try to end the stand-by mode by picking up a new probe cover. If the problem persists, replace thermometer and contact your service personnel.
<iTemp Lbl> NO DEVICE	Monitor: INOP tone Numeric unavailable	The thermometer is not connected to the monitor.
<iTemp Label> Overrange	Monitor: INOP tone Numeric is replaced by a -?-	Temperature measurement is out of range. If temperature of measurement is within valid range, check with another thermometer.

Admitting and Discharging

The fetal monitor can store basic patient demographic information used to identify patients.

Admit/Discharge on the Monitor

This section describes how you admit and discharge patients when using the monitor as a stand-alone device (that is, when not used with an obstetrical information and surveillance system such as OB TraceVue/IntelliSpace Perinatal).

Admitting a Patient

The fetal monitor displays physiological data as soon as a patient is connected. This lets you monitor a patient who is not yet admitted. It is however important to admit patients properly so that you can identify your patient on recordings.

Use the **Patient Demographics** window and its associated pop-up keys to admit and discharge patients.

To admit a patient,

- 1 Select the patient name field or select the **Patient Demogr.** SmartKey to open the **Patient Demographics** window.
- 2 Clear any previous patient data by selecting **Dischrge Patient** and then **Confirm**.
If you do not discharge the previous patient, you will not be able to distinguish data from the previous and current patients, for example, on the recording.
- 3 Select **Admit Patient**.
- 4 Enter the patient information: select each field and use the on-screen keyboard.
If a conventional keyboard is connected to the monitor you can use this to enter patient information:
 - **Last Name:** Enter the patient's last name (family name), for example **Doe**.
 - **First Name:** Enter the patient's first name, for example **Jane**.
 - **MRN:** Enter the patient's medical record number (MRN), for example **12345678**. Depending on your configuration this field may be labeled differently e.g. "Record ID", "Lifetime Id", "Account Number", "Serial Number", or "SSN".
 - **Gestational Age:** Enter the gestational age of the pregnancy. Enter one numeric for the week (0-50) and one numeric for the day (0-6).
 - **Date of Birth:** Enter the patient's date of birth.

- 5 Select **Confirm**. The patient status changes to admitted. If the recorder is running, the recorder stops and immediately restarts to annotate the new patient data.

Quick Admitting a Patient

Use **Quick Admit** to quickly admit a patient using only a limited set of demographic data.

- 1 Select the **Quick Admit** SmartKey.
- 2 Enter the required data (ID fields or last name depending on configuration) with the keyboard or a barcode scanner.
- 3 Select **Enter**.
- 4 In the confirmation window, select **Confirm** to stop monitoring for the previous patient or free up the fetal monitor (if confirmation is configured).
- 5 If the monitor is connected to an Information Center and only the ID field is entered, the patient name may be set to - - - at the Information Center, depending on the configuration.

Complete the rest of the demographic details as soon as possible to fully identify the patient on the network, on the monitor, and on printed reports. To complete the details, open the **Patient Demographics** window and complete all required fields.

Editing Patient Information

To edit the patient information after a patient has been admitted, select the patient name field on the main screen of the fetal monitor to open the **Patient Demographics** window, and enter the required changes.

Discharging a Patient

You should always perform a discharge even if your previous patient was not admitted. A discharge:

- clears the information in the **Patient Demographics** window.
- resets all monitor settings to the settings defined in the User Default.
- advances the paper automatically if the recorder is running.
- stops the fetal recorder.

When a patient is discharged from the monitor, all patient demographic data is deleted (trace data is not affected).

Depending on your configuration, monitoring for a patient may end automatically when the fetal monitor has been powered off or in standby mode for a set time.

To discharge a patient,

- 1 Select the patient name field to display the **Patient Demographics** window and associated pop-up keys.
- 2 Select the pop-up key for **Discharge Patient**.
- 3 Select **Confirm** to discharge the patient.

CAUTION

In order to ensure that the settings are reset to user defaults for a new patient, always discharge the previous patient from the fetal monitor.

NOTE

In order to ensure a continuous record, it is recommended to discharge the patient at the monitor before performing a new patient admission in OB TraceVue/IntelliSpace Perinatal.

New Patient Check

The fetal monitor can be configured to ask you in certain situations:

- after a specified power-off period
- after a specified standby period

whether a new patient is now being monitored. The pop-up window is entitled **Is this a New Patient?** The monitor offers a **Yes** key to discharge the previous patient and begin monitoring a new patient, and a **No** key to continue monitoring with the current patient data and settings.

The time periods for the two conditions can be configured independently.

OB TraceVue/IntelliSpace Perinatal

When the monitor is connected to an OB TraceVue/IntelliSpace Perinatal system over a LAN connection, the OB TraceVue/IntelliSpace Perinatal system acts as the "master" over patient demographic data. All patient and location-related data that is visible on the monitor is set, overwritten or updated by the OB TraceVue/IntelliSpace Perinatal system. See the OB TraceVue/IntelliSpace Perinatal Instructions for Use for details.

Transferring a Patient

When you are transferring a patient from one bed to another, both the monitor from which you are transferring her and the monitor to which you are transferring her must be switched on and connected to the OB TraceVue/IntelliSpace Perinatal network during the patient transfer.

OB TraceVue/IntelliSpace Perinatal via RS232

In contrast to a LAN connection, when the monitor is connected to an OB TraceVue/IntelliSpace Perinatal system over an RS232 connection, the OB TraceVue/IntelliSpace Perinatal system has no control over the monitor's patient admission and discharge functions.

Depending on how OB TraceVue/IntelliSpace Perinatal is configured, *either* the **Last Name, First Name** and the bed label, *or* just the bed label alone, are taken from the OB TraceVue/IntelliSpace Perinatal system. See the OB TraceVue/IntelliSpace Perinatal Instructions for Use for details.

Non-Stress Test Timer

The non-stress test (**NST**) timer shows the elapsed time for the non-stress test. The timer counts up to the time you set for the NST.

Setting NST Autostart/Autostop

You can set the recorder so that it starts automatically (**NST Autostart**) when the NST timer is started, and stops automatically (**NST Autostop**) when the NST is complete (when the set run time has elapsed). As default, **NST Autostart** is **On**, and **NST Autostop** is **Off**.

Viewing the NST Timer

You can configure the timer notification symbol, (the NST label), a progress bar and the elapsed time to be displayed in the top left-hand corner of the screen. By default, the NST timer is not displayed on the screen.



Alternatively, you can view the timer in the **Timers** window.


To open the **Timers** window:

Either

- a. Press the **Timer** SmartKey.

Or

- b. Access the NST pop-up keys (see “Accessing the NST Setup Pop-up Keys” on page 140), and press the **Timers** key.



Timer Expiry Notification

When the timer expires, the color changes from blue to green, you hear a single tone, and a message appears in the status line on the main screen.

The volume of the tone can be set in Configuration Mode.


Accessing the NST Setup Pop-up Keys

You control and set up the NST timer (for example, start, stop, or clear the timer, and set the run time) using a selection of pop-up keys that you access via any one of three possible routes:

- Via the **Timer** SmartKey (Route 1).
- Via the **Main Setup** SmartKey (Route 2).
- Via the NST menu entry in the menu that pops-up when you touch the top left-hand corner of the screen (Route 3).

Via the Timer SmartKey (Route 1)

- ◆ Press the **Timer** SmartKey. The **Timers** window opens, and the pop-up keys for controlling/setting up the NST timer appear (see “Pop-up Keys for NST Timer Setup” on page 141).



Via the Main Setup SmartKey (Route 2)

- 1 Enter the **Main Setup** menu using the SmartKey.
- 2 Select **NST** to enter the **Setup NST** menu. At the same time, the pop-up keys for setting up the NST timer appear (see “Pop-up Keys for NST Timer Setup” on page 141).



Via the NST Field (Route 3)

- ◆ Select the NST field displayed in the top left-hand corner of the screen (when so configured). The pop-up keys for controlling/setting up the NST timer appear (see “Pop-up Keys for NST Timer Setup” on page 141).

Pop-up Keys for NST Timer Setup

Pop-Up Keys	Selecting this pop-up key lets you....	Comments
Start	Start the timer.	
Stop	Stop the timer, allowing either restarting after a pause (Start) or clearing (Clear key).	
Setup NST	Enter the Setup NST menu. From here you can set the run time.	This pop-up key is not available with Route 2, as the Setup NST menu is already open.
Timer	Return to the Timers window.	This pop-up key is not available with Route 1, as the Timers window is already open.

Run Time

The run time can be set from 10 to 60 minutes. To set the run time, you first need to enter the **Setup NST** menu:

- 1 To enter the **Setup NST** menu:

Either

Enter the **Main Setup** menu using the SmartKey. Then select **NST**.

Or

Access the NST pop-up keys (see “Accessing the NST Setup Pop-up Keys” on page 140), and press the **Setup NST** key.

- 2 Select **Run Time**.



Non-Stress Test Report

It is generally accepted that a non stress test (**NST**) allows you to assess fetal well-being. The monitor's NST report process uses fetal ultrasound (but not DECG) heart rate traces and the Maternal Toco trace, generating a printed report when criteria are met and it is an indication of the fetal well-being. The American term Non Stress Test (NST) is used for antepartum testing. The interpretation algorithm and rule set are equivalent to those implemented in OB TraceVue Revision G.xx or IntelliSpace Perinatal H.xx and higher, and are based on the 2008 NICHD guidelines.

An NST report is a diagnostic aid, but it does not replace the clinician's judgment. The interpretation and the appropriate clinical response remain with the clinician.

A fetus normally produces characteristic heart rate patterns. Average baseline variability and acceleration of the FHR in response to fetal movement are considered reassuring signs. This test does not take into account any form of external fetal stimulation.

For every active ultrasound fetal heart rate measurement, one NST report can reside in the monitor's memory. The reports are cleared when you discharge a patient and when you start a new NST report.

When the NST Report option is available and the "NST Report" feature is "on", the NST status for all available ultrasound fetal heart rate measurements are displayed on the screen.

The minimum displayed information is:

- NST identification (by FHR number: 1, 2, 3)
- Current NST status (by color: inverse for "not started yet", white for "running", yellow for "stopped", green for "finished")

Setting Up an NST Report

To setup NST Report functionality:

- 1 Enter the **Main Setup** menu and select the **NST Report** or
- 2 Select the **NST Report** SmartKey.
- 3 Press the "Setup" pop-up key.
- 4 Set your configuration options.

Select from:

- **NST Analysis.** Choose from **On** or **Off**.
This switches the report feature on or off. This is linked to the NST timer. Both must be set to **On** for the NST report to function.
- **Report Recording.** Choose from:
 - **Manual** - press the **Record Report** pop up key to trigger a manual request.
 - **After Recorder Stop** - report is recorded as soon as recorder becomes idle.
 - **Immediately** - If a realtime recording is running, the monitor pauses it. The recording is continued after the report has been recorded.

Average short term variability (STV) value is documented in [bpm] and [ms] if STV is configured as part of the NST Report. This parameter is not considered as reassuring criteria.

NST Report Status Window

The NST Report window displays a detailed overview of the current NST status for any available ultrasound fetal heart rate measurement. You can see:




- NST Status - whether it is ready, ongoing, or the time and date at which it was stopped, or at which it finished.
- Elapsed time - the time that has elapsed since the NST began.
- Accelerations - the number of FHR accelerations detected so far.
- Baseline - the average baseline value.
- Variability - the average variability value.
- Short Term Variability - the current short term variability (STV) value.
- Decelerations - the number of FHR decelerations detected so far.
- FHR Availability - current statistical FHR availability value.
- Sinusoidal - the current status of sinusoidal rhythm detection.

For criteria not yet met, a white arrow symbol marks the overall status on the top line, and also appear against every criterion not yet met. A yellow symbol indicates detection of severe or prolonged decelerations.

The pop-up keys let you perform the following actions:

- **FHR1, FHR2, FHR3** - switch to the window showing the current NST status for the fetal heart rate.
- **Record Report** - print the NST Report on paper.
- **Record Trace** - record the trace episode that belongs to the current report. Depending on device usage, the trace recording might be incomplete.
- **Setup** - open the **Setup NST Report** window.

Example NST Report

Field	Field Content
Report Title, with FHR label and date	NST Report for FHR1 on 12 Oct. 2009
Product Information	Product DE53102345 G.01.70, OB A.04.24, Toco DE52401090, FHR1 DE00002345 A.05.26
Patient Information	Rogers, Alice Age: 27 Gestational Age: Week 34, Day 5
Start time, end time, Elapsed time, configured Run time	Time: 11:34 - 12:06 Elapsed time: 32 minutes Run time: 20 minutes
Overall one-line NST result summary	NST Criteria*: not met
Title	Trace Interpretation Summary
Result Accelerations	Accelerations: 2 at: 11:59 12:02
Result: Contractions	Contractions: 3 at: 11:57 12:00 12:04
Result: Baseline and Variability	Baseline: 125 bpm (Range: 118 bpm - 129 bpm) Variability: 23 bpm (Range: 20 bpm - 24 bpm) Short Term Var.: 0.9 bpm (3.8 msec)
Statistics: FHR availability	FHR available: 95%
Result: Decelerations	Decelerations: 1 at: 11:58   severe  prolonged

Field	Field Content
Result: Sinusoidal Rhythm detected	Sinusoidal: No
Result: Decelerations before Reporting Period	<p>Events before Reporting Period: Decelerations: 1</p> <p>at: 11:38 ↓</p> <p>↓ severe</p> <p>↔ prolonged</p> <p>This field is enabled if there were decelerations between the start of NST and the start of the reporting period.</p>
Guideline / Criteria Information	<p>(*) Interpretation Criteria based on guideline "NICHD 2008, v01"</p> <p>User-defined criteria for CTG tracing:</p> <ul style="list-style-type: none"> • valid FHR for 90% of reporting period • baseline heart rate between 120 and 160 bpm • at least 2 accelerations in 10 min • not more than 1 decelerations • moderate baseline variability (5 - 25 bpm) • short term var. greater than 5.0 msec <p>Additional Criteria:</p> <ul style="list-style-type: none"> • no severe or prolonged decelerations • no sinusoidal pattern in reporting period

NST Criteria

The patient is monitored for a user-definable period of time (10 - 60 minutes in steps of 5 minutes). The test is considered reassuring when the following criteria are met:

- The fetal heart rate is valid at least 90% (this is configurable) of the specified time span.
- The FHR features a user-defined minimum number of accelerations.
- The FHR features a user-defined maximum number of tolerated decelerations, and does not include severe or prolonged decelerations, which are never tolerated.
- The average baseline fetal heart rate lies within the user-defined limits for low heart rate and high heart rate over the whole time span.
- The FHR exhibits a moderate variability (user-defined) for the specified time span.

An NST Report is generated when the reassuring criteria are met the first time in the current monitoring phase. When performing NST with twins or triplets, a separate NST Report is generated for each fetus.

In the printed report, the average short term variability (STV) value is shown in bpm and ms. This is not part of the reassurance criteria.

After the reassurance criteria have been met, the clinician can print the NST Report and then turn the fetal monitor off, or may continue fetal monitoring and print the report at any time.

Non-Reassuring Report

If the reassurance criteria are **not** met when the test has run for 90 minutes, or if you stop anytime during the 90 minute period, then the test is stopped, and a report is generated stating the reassurance criteria have not been met.

Non-Reactive NST Test

If a non-reactive test occurs and you then use acoustic stimulation, you must exercise caution in interpreting the resulting traces, as artificial stimulation is not taken into account when calculating test results.

Cross-Channel Verification (CCV)

Misidentification of Heart Rates

FHR detection by the monitor may not always indicate that the fetus is alive. Confirm fetal life before monitoring, and continue to confirm that the fetus is the signal source for the recorded fetal heart rate (see “Confirm Fetal Life Before Using the Monitor” on page 10).

To reduce the possibility of mistaking the maternal HR or pulse for FHR, or FHR1 for FHR2 or FHR3, it is recommended that you monitor both maternal HR/pulse and the heart rates of all fetuses (see “Monitoring FHR and FMP Using Ultrasound” on page 155, “Monitoring Twin FHRs” on page 171, “Monitoring Triple FHRs” on page 179, and “Monitoring Maternal Heart / Pulse Rate” on page 223).

Here are some examples where the MHR can be misidentified as the FHR, or one FHR for another FHR (twins/triplets).

- **When using an ultrasound transducer:**
 - It is possible to pick up maternal signal sources, such as the maternal heart, aorta, or other large vessels. Especially if the recorded MHR, and any other artifact is over 100 bpm.
 - It is possible to pick up the same fetal heart rate simultaneously with multiple transducers

NOTE

When an ultrasound transducer is connected to the monitor, but not applied to the patient, the measurement may generate unexpected intermittent FHR readings.

- **When Fetal Movement Profile (FMP) is enabled:**

The FMP annotations on a fetal trace **alone** may not always indicate that the fetus is alive. For example, FMP annotations in the absence of fetal life may be a result of:

 - Movement of the deceased fetus during or following maternal movement.
 - Movement of the deceased fetus during or following manual palpation of fetal position (especially if the pressure applied is too forceful).
 - Movement of the ultrasound transducer.
- **When using a scalp electrode (DECG):**
 - Electrical impulses from the maternal heart can sometimes be transmitted to the fetal monitor through a recently deceased fetus via the spiral scalp electrode, appearing to be a fetal signal source.

Cross-Channel Verification Functionality

The cross-channel verification functionality (CCV) of the fetal monitors compares all monitored heart rates (maternal and fetal), and indicates automatically whether any two channels are picking up the same signal, or monitoring similar values.

If the fetal monitor detects that any channels have the same or similar values the **Coincidence** INOP is issued with an **INOP tone** that can have a configurable delay. In addition yellow question marks appear next to the numerics on the touch screen that have the same or similar values. On the recording trace there is also a question mark printed out from the point where recorded traces continuously overlap.

Visual Aids for CCV Detection

!! Coincidence

Coincidence INOP appears on the touch screen of the fetal monitor.



Question mark next to the numerics that show the same or similar values appears on the touch screen of the fetal monitor.



Question mark printed out on the recorder trace from the point where two measured values coincide.

Overview of Cross-Channel Comparisons

Measurements from Transducers

Measurement	Transducer
FHR US	From Ultrasound or CL Ultrasound transducer
DFHR (DECG)	From a fetal scalp electrode
Maternal pulse	From Toco MP, or CL Toco ⁺ MP transducer
Maternal pulse	From SpO ₂ or CL SpO ₂ Pod
Maternal heart rate	From MEEG electrodes

Measurement Comparison Done by the Fetal Monitor for Cross-Channel Verification

	FHR1 (US)	FHR2 (US)	FHR3 (US)	DFHR (DECG)
FHR1 (US)		✓	✓	✓
FHR2 (US)	✓		✓	✓
FHR3 (US)	✓	✓		✓

	FHR1 (US)	FHR2 (US)	FHR3 (US)	DFHR (DECG)
DFHR (DECG)	✓	✓	✓	
Maternal Pulse (Toco MP)	✓	✓	✓	✓
Maternal Pulse (SpO ₂)	✓	✓	✓	✓
Maternal HR (MECG)	✓	✓	✓	✓

Coincidence Examples


Coincidence of Maternal Pulse and FHR

When the maternal pulse and FHR are being monitored, and the measured values are very similar or the same, the coincidence question mark is displayed on the monitor's screen above both of the corresponding numerics (in this case maternal pulse and FHR). Often the signal loss or coincidence happens because the fetal or maternal movement displaced the ultrasound transducer, and a repositioning of the transducer is necessary.



- 1 **Coincidence** INOP
- 2 Coincidence question mark by **FHR1**
- 3 Coincidence question mark by pulse from Toco MP

The coincidence question mark is also printed on the trace paper next to the corresponding FHR and maternal pulse.



- 1 Printed coincidence question mark on trace
- 2 FHR1 and pulse traces


Coincidence of Twins/Triplets FHRs

When both **FHR1** and **FHR2** are being monitored, and the measured values are very similar or the same, the coincidence question mark is displayed on the monitor's screen above both of the corresponding numerics (in this case **FHR1** and **FHR2**).



- 1 **Coincidence** INOP
- 2 Coincidence question mark by **FHR1**
- 3 Coincidence question mark by **FHR2**

The coincidence question mark is also printed on the trace paper next to the corresponding **FHR1** and **FHR2**.



- 1 Printed coincidence question mark on trace
- 2 **FHR1** and **FHR2** traces

Recommended Actions for INOP Coincidence

- 1 Confirm fetal life by palpation of fetal movement or auscultation of fetal heart sounds using a fetoscope, stethoscope, or Pinard stethoscope.
- 2 Manual determination of the maternal pulse and comparison with the fetal heart rate sound signals from the loudspeaker.
- 3 Reposition the transducer, or ensure that the fetal scalp electrode is placed correctly, until you receive a clear signal and the monitor is no longer issuing the **Coincidence** INOP.
- 4 In case of difficulties deriving a stable maternal pulse reading using the Toco MP or CL Toco⁺ MP transducer, use SpO₂ or the CL SpO₂ Pod instead. In case of similar problems with the pulse measurement from SpO₂, use MECG instead. Reasons to switch the method for deriving a maternal pulse or heart rate include: motion artifacts, arrhythmia, and individual differences in pulse signal quality on the abdominal skin (via Toco MP).
- 5 If you cannot hear the fetal heart sounds, and you cannot confirm fetal movement by palpation, confirm fetal life using obstetric ultrasonography.

Monitoring FHR and FMP Using Ultrasound

To monitor a single FHR externally, you use an ultrasound transducer attached to a belt around the mother's abdomen. The ultrasound transducer directs a low-energy ultrasound beam towards the fetal heart and detects the reflected signal. Your monitor can also detect fetal movements and print the fetal movement profile (**FMP**) on the trace. Monitoring using ultrasound is recommended from the 25th week of gestation for non-stress testing or routine fetal monitoring.

WARNING

Performing ultrasound imaging or Doppler flow measurements together with ultrasound fetal monitoring may cause false FHR readings, and the trace recording may deteriorate.

Technical Description

Fetal monitors use the ultrasound Doppler method for externally monitoring the fetal heart rate. Using the Doppler method, the transducer (in transmitter mode) sends sound waves into the body which are then reflected by different tissues. These reflections (Doppler echoes) are picked up by the transducer (in listening mode). These Doppler echoes are amplified and sent to the monitor's speaker through which the fetal heart signal can be heard. In parallel the Doppler echoes are processed through an auto correlation algorithm to determine the fetal heart rate (FHR). The FHR is displayed on the monitor's numeric display and on the recorded trace.

Properly representing the fetal heart rate using a device that derives heartbeats from motion is a formidable task and the limitations of the technology will be discussed shortly. Basic fetal cardiac physiology may contribute to difficulties in obtaining a reliable ultrasound signal.

A heart rate pattern of a fetus is capable of extraordinary variation, ranging from a stable pattern with minimal variation while the fetus is "asleep" to robust accelerations of 40-60 bpm above baseline rate over a few seconds, or exaggerated variability when the fetus is active. Decelerations of the rate 60-80 bpm below baseline may develop even more abruptly than the accelerations. The Beat-to-Beat arrhythmias of the FHR may further exaggerate how much the the FHR variability at the bottom of the variable decelerations, or in the presence of fetal breathing movements, which also tend to lower the fetal heart rate. The recognition of these normal variations in fetal heart rate patterns will greatly assist in the separation of genuine fetal information from the artifact.

Limitations of the Technology

All tissues moving towards or away from the transducer generate Doppler echoes. Therefore, the resulting signal that is provided to the monitor's speaker, and for further fetal heart signal processing, can contain components of the beating fetal heart wall or valves, fetal movements, fetal breathing or hiccup, maternal movements such as breathing or position changes, and pulsating maternal arteries.

The fetal heart signal processing uses an auto correlation algorithm to obtain periodic events such as heart beats. If the signal is erratic such as from a fetal arrhythmia, the ultrasound device may have

trouble tracking the abrupt changes, and may misrepresent the true FHR pattern. Signals such as those from moving fetal limbs are usually very strong, thereby masking the fetal heart signal. During prolonged movements where the fetal heart signal is masked, the FHR appears blank on the numeric display and as a gap on the recorded trace. Fetal position changes, maternal position changes, or uterine contractions can move the fetal heart partly or fully out of the ultrasound beam resulting in signal loss, or even picking up Doppler echoes from pulsating maternal arteries. In these cases a maternal heart rate or sometimes even a rate resulting from the mixture of fetal and maternal signals may be displayed on the monitor's numeric display and on the recorded trace.

In contrast to the timely well-defined R-peak of an ECG signal obtained with a fetal scalp electrode, the ultrasound Doppler signal from a fetal heart consists of multiple components from atria (diastole), ventricles (systole), valves, and pulsating arteries. These components vary depending on fetal and transducer position and angle, and are further modulated by factors such as fetal or maternal breathing. These effects may produce what is called "artifact". Optimal transducer positioning therefore is key to minimizing these effects and thereby minimizing artifact.

Misidentification of MHR as FHR

FHR detection by the monitor may not always indicate that the fetus is alive. Confirm fetal life before monitoring, and continue to confirm that the fetus is the signal source for the recorded heart rate (see "Confirm Fetal Life Before Using the Monitor" on page 10 and "Cross-Channel Verification (CCV)" on page 149).

What You Need

- Ultrasound transducer
- Toco MP or CL Toco⁺ MP transducer
- Ultrasound gel
- Transducer belt (and optional butterfly belt clip, if applicable)

Cableless Monitoring - Important Considerations

When using an Avalon CL or Avalon CTS Fetal Transducer system with your monitor, note the following:

Refer to "Telemetry" on page 93 for general rules regarding the use of cableless transducers from an Avalon CL or Avalon CTS Cableless Fetal Transducer system.

- When using a cableless ultrasound transducer from an Avalon CL or Avalon CTS system to measure the fetal heart rate, note that you cannot use any other ultrasound transducer (whether cableless or wired) at the same time.


CAUTION

Never use ultrasound transducers connected to more than one fetal monitor on the same patient.

- When using an Avalon CL or Avalon CTS you should be aware that FMP is not recommended when the mother is likely to move, and you should disable Fetal Movement Profile (FMP) on the fetal monitor (**Fetal Movement Off**) if the mother is walking. See also "Fetal Movement Profile" on page 160.

- With the Avalon CL Transducer System you can now monitor twins and triplets with cableless transducers. The Avalon CTS System does not have this option.

The wireless symbol appears next to the measurement label, indicating that the measurement is being made by a cableless transducer.



WARNING

- Ensure that the conductive parts of the fetal scalp electrode and the maternal leg plate electrode do not contact other conductive parts, including earth.
 - Indication of the heart-rate may be adversely affected by the operation of cardiac pacemaker pulses or by cardiac arrhythmias.
 - During ambulant FHR monitoring, the chance of losing the signal or detecting the maternal heart rate is higher than during stationary monitoring. The frequency of the patient's walk may be detected, and mistaken for a FHR signal.
 - Check the mother's pulse periodically during monitoring and compare this with the FHR signal. Beware of mistaking a "doubled" maternal heart rate for FHR. If a fetus is dead, there is a risk that the maternal heart rate is monitored and misinterpreted as the fetal heart rate. Therefore, the simultaneous monitoring of maternal heart rate (preferably, the maternal ECG) and the fetal heart rate is encouraged.
 - Do not interpret maternal movements as fetal movements.
 - Artifacts: FMP artifacts are generated during fetal heart rate searching by changing the transducer position, therefore Philips fetal monitors enable the FMP only after detecting a valid heart rate signal for several seconds. FMP is not recommended when the mother is likely to move, and you should disable Fetal Movement Profile (FMP) at the fetal monitor (**Fetal Movement Off**) if the mother is walking.
 - Gaps in maternal heart rate detection can occur:
 - if the transducer is not correctly positioned.
 - due to the pulsation of uterine blood vessels.
 - if the fetus moves.
-

Preparing to Monitor

Prepare for ultrasound monitoring using the list below. The standard procedures in use in your facility determine the sequence of actions.

- 1 Determine fetal position.
- 2 Fasten the belt around the patient.
- 3 Switch on the monitor and the recorder.
- 4 Connect the transducer to a free socket. The signal quality indicator for the heart rate initially displays an invalid signal.
- 5 Apply a thin layer of ultrasound gel to the underside of the transducer.

CAUTION

Using ultrasound gel not approved by Philips may reduce signal quality and may damage the transducer. This type of damage is not covered by warranty.

- 1 Place the transducer on the abdomen, if possible over the fetal back or below the level of the umbilicus in a full-term pregnancy of cephalic presentation, or above the level of the umbilicus in a full-term pregnancy of breech presentation. Work the transducer in a circular motion to ensure the gel layer makes good contact.
- 2 When the sensor is connected correctly and you receive a good signal, the signal quality indicator should be filled out. If an inadequate signal is produced, the signal quality indicator will indicate a poor signal, and no numeric will appear on the screen.
- 3 Adjust the audio volume of the monitor's loudspeaker to a clearly audible level, while moving the transducer over the abdomen. When you have a good signal, secure the transducer in position below the belt.

WARNING

Periodically compare the mother's pulse with the signal coming from the monitor's loudspeaker to ensure that you are monitoring fetal heart rate. Do not mistake a doubled or elevated MHR for FHR.

The ultrasound transducer may warm slightly (less than 1°C/1.8°F above ambient temperature) when applied to the patient. When not applied, the transducer can reach a maximum temperature of 44°C/112.2°F at an air temperature of 40°C/104°F.

Selecting Fetal Heart Sound

You can listen to the fetal heart sound from **one** ultrasound transducer at a time. When the fetal heart sound is selected for an FHR channel, you see the audio source symbol next to the FHR numeric label for that channel.





- 1 FHR1
- 2 Audio source symbol

To select the audio source for an FHR channel:


- 1 Enter the **Setup FHR1** menu for the channel you want to hear (FHR1 used as an example).
- 2 Press **Select Audio**. It may take a few seconds for the audio source symbol to appear.



Changing the Fetal Heart Sound Volume

The FHR volume symbol at the top right of the **Fetal Heart Sound Volume** window gives you an indication of the current volume. To change the volume:

- 1 Select the volume symbol. The volume scale pops up.
- 2 Select the required volume from the volume scale.



- 1 Fetal Heart Sound Volume

Fetal Movement Profile

The Fetal Movement Profile (FMP) parameter detects fetal movements via an ultrasound transducer connected to the monitor. Only the fetus monitored on the FHR1 channel is monitored for FMP.

Once you have enabled FMP (see “Switching FMP On and Off” on page 161), it is triggered automatically whenever:


- You connect an ultrasound transducer.
- A patient is discharged.

When FMP is enabled, the ultrasound transducer detects most fetal body movements. Eye movements are not detected and movement of the feet and hands may not be detected. Positioning or repositioning of the transducer is recorded as fetal movement. Maternal movement, excessive fetal breathing, or fetal hiccups may also be recorded as fetal movement (also in case of fetal demise or during the second stage of labor). You can mark these artifacts on the trace paper using either the remote event marker, or the event marker key as described in “Marking an Event” on page 79. FMP should be interpreted with care, or disabled when the patient is ambulating or during 2nd stage of labor. Ignore these movements when you interpret the FMP. When monitoring twins or triplets, only the fetus monitored on the FHR1 channel is monitored for movement, but be aware that movements recorded for FHR1 may also be caused by movement of the second or third fetus.

The fetal movement profile (FMP) appears as "activity blocks" (see **A** below) along the top of the Toco Scale, the length of each block showing the duration of the activity.

FMP Statistics

FMP statistics are printed every ten minutes.



1 FMP enabled

- 2 FMP started here
- A Indication of current fetal movement.
- B The FMP statistics are presented as two percentage figures:
The first figure shows the percentage of detected fetal movements in the previous ten minutes.
- C The second figure shows the percentage of detected fetal movements since the start of recording.

To mark the start of the FMP statistic, FMP is printed on the paper with an arrow.



The FMP detection activates after about half a minute of steady heart rate signals (signal indicator half-full, or full) to minimize transducer positioning artifact. You will notice this deliberate delay:

- When a new patient is admitted. A patient discharge restarts the FMP statistics from zero.
- When you connect an ultrasound transducer.

Switching FMP On and Off

You can switch FMP on and off from any FHR channel. For example, to set it from the FHR1 channel:

- 1 Enter the **Setup FHR1** menu.
- 2 Select **Fetal Movement** to toggle between **On** and **Off**.
- 3 Return to the main screen.

Troubleshooting

Problem	Possible Causes	Solutions
Erratic trace Erratic display	Fetal arrhythmia	Consider monitoring FHR using DECG after the rupture of membranes.
	Obese patient	
	Transducer position not optimal	Reposition transducer until signal quality indicator shows a good signal (at least half-full).
	Belt loose	Tighten belt
	Too much gel	Remove excess
	Very active fetus	-
	Insufficient gel	Use enough gel to ensure the transducer makes good contact with the mother's skin.
Signal quality indicator is continuously poor	Transducer position not optimal	Reposition transducer until signal quality indicator shows a good signal (at least half-full).
	FHR less than 50 bpm (and the FHR is audible)	If membranes are ruptured, using a fetal scalp electrode (FM30 and FM50 only) allows measurement of FHR down to 30 bpm.
Questionable FHR	Recording MHR by mistake	Reposition transducer. Confirm fetal life.
	Recording periodic signals when the transducer is not applied to the patient	Disconnect all NON-USED ultrasound transducers, as continuous, regular mechanical, or electromagnetic influences can result in an artificial trace.
	Recorded FHR appears to be suspiciously higher, or suspiciously lower, than real FHR. In very rare cases, half- or double-counting of the FHR can occur.	If you have reason to question the validity of the recorded FHR, always verify FHR by independent means (by auscultation, for example). Measure maternal pulse by independent means.
FHR not recorded	FHR is less than 50 bpm or over 240 bpm	If membranes are ruptured, using a fetal scalp electrode (FM30 and FM50 only) allows measurement of FHR down to 30 bpm.
		If FHR is outside of the specified range, verify FHR by independent means.
FHR1 Equip Malf or FHR2 Equip Malf or FHR3 Equip Malf INOP displayed.		See “Patient Alarms and INOPs” on page 121.
FHR1 Signal Loss or FHR2 Signal Loss or FHR3 Signal Loss INOP displayed.		
FHR1 Unplugged or FHR2 Unplugged or FHR3 Unplugged INOP displayed.		
If you suspect the transducer is malfunctioning		Test the transducer.

Additional Information

Artifact in Fetal Heart Rate Measurement

How to detect it and reduce its occurrence using the Avalon Fetal Monitor

The ultrasound derived FHR measurement technique in Avalon fetal monitors, like all other ultrasound fetal monitors' FHR measurement techniques, has limitations that can lead to misrepresentation of the fetal heart rate pattern and potential misinterpretation of the fetal condition.

An incorrect interpretation of the trace may lead to either unnecessary interventions, or to failure to detect fetal distress and the need for intervention. Thus, the on-going evaluation of the recorded trace requires regular confirmation that the trace represents the true FHR. Specific situations requiring such confirmation include the following:

- After starting a measurement or changing a transducer
- After maternal position changes, for example during pushing with contractions
- When the tracing shows abrupt changes in baseline rate, variability, or pattern (decelerations to accelerations) especially in the second stage of labor
- When the baseline maternal heart rate is within about 15 bpm of the FHR
- When the user is unable to determine a baseline rate and variability occurs between consecutive contractions

There are several ways to verify the source and/or accuracy of the recorded fetal heart rate pattern. These include:

Verification of the FHR with:

- An obstetric stethoscope,
- Ultrasound imaging, or
- A fetal scalp electrode

Verification of the maternal heart rate:

- Using pulse oximetry - for a maternal heart rate pattern displayed simultaneously with the FHR (Cross-Channel Verification (CCV) feature),
- Using Maternal ECG - for a maternal heart rate pattern displayed simultaneously with the FHR (CCV feature), or
- Manual determination of the maternal pulse

Whenever possible measure the maternal pulse rate to make use of the monitor's Cross-Channel Verification (CCV) feature, especially during the second stage of labor, or when the maternal pulse is elevated over 100 bpm. The Philips Avalon fetal monitor provides a Toco MP or CL Toco⁺ MP transducer for maternal pulse detection and the creation of a maternal heart rate pattern plotted on the same recorder as the FHR pattern. In case of difficulties deriving a stable maternal pulse reading using the Toco MP or CL Toco⁺ MP transducer, use SpO₂ or MECCG instead.

When either of these parameters is utilized, the monitor will automatically and continuously perform a CCV of the maternal heart rate pattern against the FHR pattern displayed on the monitor. If the patterns and rates are similar, the CCV provides an alarm that both rates are probably from the same source (i.e., they both represent the maternal heart rate pattern and the fetus is not being monitored). Repositioning the ultrasound transducer will usually correct this, but it may be necessary to apply a fetal scalp electrode. Advising the mother to temporarily cease pushing during contractions may help to more rapidly resolve any uncertainty in this situation.

Doubling: The auto correlation algorithm can display a doubled fetal or maternal heart rate if the duration of diastole and systole are similar to each other, and if the heart rate is below 120 bpm. Doubling, usually brief, is accompanied by an abrupt switch of the trace to double the baseline value.

Halving: With fetal tachycardia (above 180 bpm) and some interference from breathing or maternal arteries the auto correlation algorithm may only recognize every second beat resulting in a halved rate for a limited time. If the actual FHR is above the maximum limit of the monitor (240 bpm), the algorithm will also half-count. Halving is accompanied by an abrupt switch of the trace to exactly half the prior baseline value. This switch may simulate a FHR deceleration and be referred to by clinicians as a “false deceleration.”

Switching to maternal heart rate (also referred to as "Maternal Insertion"): The fetal heart can move partly or fully out of the ultrasound beam and the auto correlation algorithm may then pick up and display the maternal heart rate. Depending on the signal mix in the ultrasound signal, switching to the maternal heart rate may mimic several conditions with the potential for erroneous interpretation and response as follows:

- **The switch to the maternal heart rate may simulate a FHR deceleration** (i.e., a decrease of the fetal heart rate, and be referred to by clinicians as a “false deceleration”).
- **The maternal heart rate may simulate a normal fetal heart rate pattern** (i.e., it may mask a FHR deceleration or fetal demise).
Especially during pushing with contractions in the second stage of labor, the maternal heart rate may increase to the point where it may equal or exceed the fetal rate. Here the maternal trace may mimic a normal fetal trace while the fetus may be having decelerations or fetal demise has occurred. This change from fetal to maternal heart rate pattern may not be at all obvious unless CCV is used and represents the most dangerous pitfall of all the artifacts because fetal distress may go unrecognized.
- **The maternal heart rate may simulate a FHR acceleration**, which is an increase of the fetal heart rate.
During expulsive efforts, the maternal heart rate normally accelerates and may be at or above the normal FHR range.
- **The FHR may display gradual appearing decelerations.** Generally, the “false decelerations” described above are abrupt. Rarely, combinations of “noisy/erratic signal” associated with changes in maternal and/or fetal rate or movement will produce more gradual appearing “false decelerations” but these are usually short-lived with an **abrupt** return to an obviously stable FHR baseline.

“Noisy/Erratic” signals: With mixed or weak signals the tracing may reveal very brief episodes of erratic recorded traces. These represent the auto correlation algorithm finding brief sequences of apparent and persistent heartbeats amidst a mixed or weak signal. These erratic recorded traces are commonplace, especially in association with fetal or maternal movement. **During prolonged periods of such noisy/erratic signals, the fetus is not being adequately monitored.**

Drop out: With mixed or weak signals there may be no heart rate tracing at all. These episodes reflect that if the algorithm does not find an apparent and persistent heartbeat amidst a mixed or weak signal, it will not print a heart rate on the tracing. Brief episodes of drop out are commonplace, especially in association with fetal or maternal movement. **During prolonged periods of drop out, the fetus is not being adequately monitored.**

Multiple Fetuses

With multiple fetuses, the potential to experience these artifacts is increased. Positioning of the transducer is even more critical. Ultrasound scanning should be used to help with positioning of individual transducers. See also “Monitoring Twin FHRs” on page 171 and “Monitoring Triple FHRs” on page 179.

Obtaining a Good Heart Signal

To successfully position the ultrasound transducer, first determine the fetal position using palpation. Position the transducer over the strongest audible fetal heart sound from the monitor’s speaker and **wait at least six seconds after each transducer adjustment** to verify a good signal quality displayed on the Signal Quality Indicator and a consistent FHR numeric display. Having determined the position that provides a strong fetal signal, fix the transducer on the abdomen with the belt.

If the quality of the signal or the appearance of the heart rate trace from the ultrasound transducer is questionable, the transducer should be repositioned as described above. Alternatively, the use of an ultrasound scanner will greatly facilitate the determination of the optimal site for the ultrasound heart rate transducer. Factors during the second stage of labor that may influence the quality of the FHR tracing obtained with ultrasound include:

- Uterine contractions
- Changing contour of the maternal abdomen
- Maternal body movement - positioning
- Maternal expulsive efforts - pushing
- Maternal tachycardia/accelerations with contractions
- Fetal decelerations, Fetal tachycardia
- Delayed return of the fetal heart rate from a deceleration
- Descent of the fetus in the birth canal
- Rotation of the fetus in the birth canal

In some cases during the second stage of labor, a good and reliable ultrasound FHR signal may not be obtainable, and the use of a fetal scalp electrode must be considered (fetal ECG).

Heart Rate Sound

The heart rate sound emitted by the device is a representation of movement that, in most cases, permits accurate auscultation of the FHR corresponding to the FHR displayed on the monitor and rate pattern depicted on the trace recording. On occasion, the user may hear a FHR that differs from the FHR display and the recorded trace. This may occur in situations where the fetal heart moves partly out of the transducer ultrasound beam. In these cases, the user may hear the FHR emitted from the monitor’s speaker, even though another periodic signal (usually the maternal heart rate) has become stronger. The auto correlation algorithm will display the stronger maternal heart rate, despite the persistence of a weaker fetal signal. These occurrences are usually very brief and, if persistent, can be addressed by repositioning the transducer.

Signal Quality Indicator

Avalon compared with its predecessor, the Series 50

Signal quality indicator on Avalon fetal monitors:

Instead of a traffic light-like design (red – yellow – green) used on the Series 50, the signal quality on the Avalon fetal monitor is indicated by a triangle on the touch screen that is displayed in one of three ways:




Avalon fetal monitor signal quality indicator display:

- 1 Completely filled triangle, indicating good signal quality (good/full).
- 2 Half-filled triangle, indicating limited signal quality. This condition may indicate a weak or ambiguous signal. If this status persists, reposition the transducer (acceptable/medium).
- 3 Empty triangle, indicating insufficient signal quality. No FHR is displayed on the monitor's numeric display or the recorded trace. If this status persists, reposition the transducer (poor/no signal).

Examples of Artifacts

When monitoring the maternal ECG, a beat-to-beat maternal heart rate trace is printed alongside the FHR recorded trace. When monitoring the maternal SpO₂ derived pulse rate, a filtered and averaged heart rate trace is printed.

Following are recorded trace examples of complaints received regarding inaccurate output from the Avalon monitors. Scaling is 3 cm/min and 30 bpm/cm.




Double-Counting	
Contractions	Excessive, coupling, hypertonus
Artifact	Double-Counting
Comment	Reassuring tracing. The excessive uterine activity should prompt discontinuation of any oxytocic agent.
Remediation	The true fetal rate can be confirmed by auscultation or by fetal scalp electrode.

Half-Counting	
Baseline Rate	120
Baseline Variability	Moderate
Accelerations	Present
Decelerations	Not apparent
Contractions	Minimal
Artifact	Half-counting, noise, drop out
Comment	Reassuring tracing. The half-count at 4-5 minutes into the tracing may simulate a fetal deceleration, but the abruptness and the lack of any compensatory changes when the normal rate returns suggests that this is half-counting. Insertion of the maternal heart rate (see below) may produce a similar pattern. Note also very brief episodes of half-counting, maternal insertion, and signal drop out.
Remediation	Auscultation or the application of a direct scalp electrode, if feasible, will reveal the true fetal heart rate.

Maternal-Switching (Maternal Insertion)	

11 Monitoring FHR and FMP Using Ultrasound

Maternal-Switching (Maternal Insertion)	
Baseline Rate	170 - Tachycardia
Baseline Variability	Moderate
Accelerations	Unable to determine
Decelerations	Absent
Contractions	Absent
Artifact	Maternal insertion, noise
Comment	The fetus has an elevated baseline rate of about 170 bpm with minimal to moderate variability. The ability to assess fetal status is limited because about half of the tracing displays the maternal heart rate.
Remediation	The application of a maternal transducer (ECG or pulse oximeter) will likely resolve any possible confusion with the tracing. Repositioning the transducer may produce a more reliable tracing. Consideration must also be given to applying a fetal scalp electrode.

Noisy/Erratic Signal and Dropout	
	
Baseline Rate	140
Baseline Variability	Moderate
Accelerations	Present
Decelerations	Absent
Contractions	Minimal
Artifact	Noisy signal, drop-out
Comment	Reassuring tracing. Note that there is episodic drop out of the signal with discontinuity of the fetal tracing.
Remediation	Either improving the position of the transducer or the application of a fetal scalp electrode will reduce the amount of artifact in the tracing.

Selection of Literature References on Artifacts

Fetal Monitoring, A Multidisciplinary Approach, Seventh edition

Susan M. Tucker, Lisa A. Miller, David A. Miller

Mosby Inc. 2009, ISBN 978-0-323-08352-2

Signal ambiguity resulting in unexpected outcome with external fetal heart rate monitoring

By Duncan R. Neilson Jr, MD; Roger K. Freeman, MD; Shelora Mangan, RNC, MSN, CNS
 American Journal of Obstetrics & Gynecology, June 2008:

Antepartal and Intrapartal Fetal Monitoring, 3rd Edition (2007)

By Michelle L. Murray, PhD, RNC

Springer Publishing Company, ISBN 0-8261-3262-6

Page 2, Table 2: **Limitations of Continuous EFM**

Item 15: "The US may detect maternal aortic wall movement and the MHR will be printed. A failure to recognize the lack of a FHR may delay appropriate management."

Page 38, "Solving Equipment Problems", Table 3: **The Ultrasound Transducer**

JOGC (Journal of Obstetrics and Gynaecology Canada)

Volume 29, Number 9, September 2007

Chapter 2: **Intrapartum Surveillance**Page S35: **"Methods of Electronic Fetal Monitoring"**

"... Among its disadvantages are the need for readjustment with maternal or fetal movements and the following: the transducer may record the maternal pulse, it may be difficult to obtain a clear tracing in obese women or those with polyhydramnios, artifact may be recorded, and there may be doubling or halving of the fetal heart rate when it is outside of the normal range."

Maternal of Fetal Heart Rate? Avoiding Intrapartum Misidentification

by Michelle L. Murray

JOGNN Clinical Issues, April 2003, 33, 93-104; 2004. DOI: 10.1177/0884217503261161

Figure 9 "The recording is of the MHR with occasional doubling."

Maternal Heart Rate Pattern – A Confounding Factor In Intrapartum Fetal Surveillance

Schifrin BS, Harwell R, Hamilton-Rubinstein T, Visser G:

Prenat Neonat Med 2001; 6:75-82.

Fetal Monitoring in Practice, 2nd Edition 1998

By Donald Gibb, S. Arulkumaran

Butterworth-Heinemann, ISBN 0-7506-3432-2

Page 65, "False or erroneous baseline because of double counting of low baseline FHR"

Page 66, "Bradycardia: fetal or maternal"

Role of Maternal Artifact in Fetal Heart Rate Pattern Interpretation


Klapholz, Henry M, MD; Schiffrin, Barry S. MD; Myrick, Richard RS

Obstetrics & Gynecology, September 1974, Volume 44, Issue 3

Testing Ultrasound Transducers

If any of the following tests fail, repeat the test using another transducer. If the second transducer passes the tests, confirming that the first transducer is defective, contact your service personnel.

If the second transducer also fails the tests, contact your service personnel.



To test an ultrasound transducer:

- 1 Switch on the monitor and the recorder.
- 2 Connect the transducer to the fetal monitor.
- 3 Select the fetal heart sound for this channel.
- 4 Increase the loudspeaker volume to an audible level.
- 5 Holding the transducer in one hand, move your other hand repeatedly towards and then away from the surface.
- 6 Check that a noise is heard from the loudspeaker.

You can test all ultrasound transducers, including the cableless ones as described above.

Monitoring Twin FHRs

The FHRs of twins are externally monitored using two ultrasound transducers. The Avalon CL Transducer system provides the option to monitor twins with cableless transducers. The Avalon CTS system does not have this option.

FM30/50 Twin FHRs are monitored throughout labor and delivery. After rupture of the membranes, you can monitor one twin externally using ultrasound, and the other internally using DECG.

Refer to the appropriate preceding chapters for contraindications, and more information about the available measurement methods.

FHR detection by the monitor does not always indicate that the fetuses are alive. Confirm fetal life before monitoring, and continue to confirm that the fetuses are the signal source for the recorded fetal heart rates. See “Confirm Fetal Life Before Using the Monitor” on page 10 and “Cross-Channel Verification (CCV)” on page 149.

Important Considerations

When monitoring:

- Ensure that you are recording two different fetal heart rates. The cross-channel verification feature alerts you if the two heart rates coincide (if both transducers are recording the same FHR). If this happens, check the trace and if necessary, reposition an ultrasound transducer to detect the second FHR correctly.
- Fetal heart rate measurements are labeled in the **order** in which you plug in the transducers for those measurements. It does not matter which fetal sensor socket you use, as the monitor allocates a channel automatically. For instance, the first transducer you connect is automatically allocated a channel, and the measurement is labeled **FHR1**, the second is labeled **FHR2**, and so on.

If you need to disconnect the transducers measuring the FHR temporarily, with the intention to continue monitoring after the temporary break (for example, if the mother needs to go to the bathroom), it is important that you reconnect the transducers in the same order as you originally connected them to make sure the measurement labels remain consistent.


Upon discharge of the patient all connected transducers are reset from left to right. Example: Only one transducer is still connected to the monitor, it was labeled **FHR2** while monitoring the previous discharged patient; it is now reset to **FHR1**.

- The transducer Finder LED lets you identify at a glance which transducer is monitoring which heart rate channel.
- The fetal sensor socket to which a transducer is connected is identified by the transducer position indicator in the blue setup menu header:

FM20/FM30



FM40/FM50




- The trace recorded for FHR1 is thicker (darker) than that recorded for FHR2. This ensures that the two heart rates are easily distinguishable. The thickness of the recorded trace can be changed in Configuration Mode.
- Remember that only one fetal heartbeat can be heard from the loudspeaker at any time.
- The audio source symbol shows you which fetus you are listening to. To hear the other fetal heartbeat, select the fetal heart rate sound for this channel (see “Selecting Fetal Heart Sound” on page 158).



- Monitor maternal pulse, especially during later stages of labor, to avoid mistaking maternal heart rate for FHR.
- Make sure you are recording the best possible signals by referring to the signal quality indicators and repositioning the transducers if necessary.
- For the Avalon CL transducer system see “Cableless Monitoring - Important Considerations” on page 156.

Monitoring Twins Externally

To monitor twin FHRs externally you need two ultrasound transducers. Follow the procedures described in “Monitoring FHR and FMP Using Ultrasound” on page 155. The transducer Finder LED lets you identify at a glance which transducer is monitoring which FHR channel, and lights when you select the FHR numeric field on the screen.




Example of the screen showing ultrasound monitoring of twin FHRs:



- 1 FHR 1
- 2 Toco parameter
- 3 FHR 2

Monitoring Twins Internally

FM30/50 Monitor one twin using the procedures described in “Monitoring FHR and FMP Using Ultrasound” on page 155. Monitor the second twin using the procedures described in “Monitoring FHR Using DECG” on page 195.



Example of a screen showing twin monitoring using a combination of US and DECG (the fetal heart rate monitored via DECG is labeled "DFHR1"/"DFHR2"/"DFHR3" on the screen):



- 1 DFHR 1
- 2 Toco parameter
- 3 FHR 2

Separating FHR Traces

To help you to interpret traces with similar baselines, you can separate the baselines by an offset of 20 bpm by switching trace separation on. For details of the offset, see “Separation Order Type” on page 174.

Switching Trace Separation On and Off

Connect transducers to the monitor to measure FHR. Depending on the measurement method, you need *either* two ultrasound transducers *or*:

FM30/50 One ultrasound and one Toco⁺ transducer, or one CL ECG/IUP transducer (to monitor DECG).

- 1 Enter the **Main Setup** menu by pressing the **Main Setup**.
- 2 Select **Fetal Recorder**.
- 3 Select **Trace Separation** to toggle between **On** and **Off**.
- 4 Exit the **Main Setup** menu.

Separation Order Type

In Configuration Mode, you can choose between two different ways, **Standard** and **Classic**, for dealing with the trace offsets on the recording (the order in which they are separated) when **Trace Separation** is **On**.

- **Standard:** the FHR2 trace is shifted up by 20 bpm (it is recorded 20 bpm higher than it really is). No offset is ever applied to the FHR1 trace - it stays where it is. (In case of a third FHR, this is shifted down by 20 bpm.)
- **Classic:** the FHR1 trace is shifted up by 20 bpm when there is more than one FHR measurement. No offset is ever applied to the FHR2 trace - it stays where it is. (In case of a third FHR, this is shifted down by 20 bpm.)

When Trace Separation is On

When trace separation is turned on, the recorder prints a dotted line labeled with the two FHRs at the top, and **+20** at the bottom.



Examples of the two methods (**Standard**, **Classic**) for determining the trace separation order are provided here.

"Standard" Separation Order

To make differentiating the traces easier, the trace from the ultrasound transducer connected to the FHR2 channel is separated from that of FHR1 by 20 bpm. In other words, the trace for FHR2 is recorded 20 bpm higher than it really is. The trace for FHR1 is never shifted.


The recorder prints a dotted line labeled **+20** across the FHR scale, to identify the trace for FHR2. The FHR trace is labeled **+20** every 5 cm:



The label for FHR2 is annotated with a black filled **+20**:



The following trace shows trace separation switched on.



Only the FHR2 trace is offset. The numerical FHR value displayed on the monitor remains unchanged. Subtract 20 from the recorded trace for FHR2 to obtain the true FHR2 value. For example, if the recorded trace shows 160, then the true FHR is 140.

"Classic" Separation Order

To make differentiating the traces easier, the trace from the ultrasound transducer connected to the FHR1 channel is separated from that of FHR2 by 20 bpm. In other words, the trace for FHR1 is recorded 20 bpm higher than it really is. The trace for FHR2 is never shifted.


The recorder prints a dotted line labeled **+20** across the FHR scale, to identify the trace for FHR1. The FHR trace is labeled **+20** every 5 cm:



The label for FHR1 is annotated with a black filled **+20**:




The following trace shows trace separation switched on.




Only the FHR1 trace is offset. The numerical FHR value displayed on the monitor remains unchanged. Subtract 20 from the recorded trace for FHR1 to obtain the true FHR1 value. For example, if the recorded trace shows 160, then the true FHR is 140.

When Trace Separation is Off

To indicate that trace separation is switched off, a dotted line labeled **+0** prints across the FHR scale.



1 **Standard** trace separation switched off here




2 **Classic** trace separation switched off here

Troubleshooting

Common problems that may occur when monitoring FHR using ultrasound are listed in “Monitoring FHR and FMP Using Ultrasound” on page 155. See also “Monitoring FHR Using DECG” on page 195 for common problems you might encounter when monitoring FHR directly.

The following problem may occur when monitoring twins.

Problem		Possible Cause	Solution
	The question mark is printed repeatedly, and appears on the screen and the INOP Coincidence is issued.	Both transducers are recording the same FHR, or one fetal transducer is recording the MHR.	Reposition an ultrasound transducer. See “Recommended Actions for INOP Coincidence” on page 154.

For more information, see “Additional Information” on page 163.

Monitoring Triple FHRs

If your monitor is equipped with the triplets option, it carries the label.



You can monitor triple FHRs externally using three ultrasound transducers. With the Avalon CL Transducer system you can now monitor triplets with cableless transducers. The Avalon CTS system does not have this option.

Refer to the appropriate preceding chapters for contraindications, and more information about the available measurement methods.

FHR detection by the monitor may not always indicate that the fetuses are alive. Confirm fetal life before monitoring, and continue to confirm that the fetuses are the signal source for the recorded fetal heart rates. See “Confirm Fetal Life Before Using the Monitor” on page 10 and “Cross-Channel Verification (CCV)” on page 149.

Important Considerations


- The procedures and any contraindications that apply for twins monitoring also apply for monitoring triplets. In addition, when monitoring triplets: Be aware that monitoring three FHRs is inherently more difficult than monitoring single or twin FHRs. The nature of the application increases the likelihood that a fetal heart rate is monitored by more than one transducer.
- Ensure that you are recording three different fetal heart rates. Pay particular attention to any coincidence of heart rates detected by the monitor's cross-channel verification feature.
- Fetal heart rate measurements are labeled in the **order** in which you plug in the transducers for those measurements. It does not matter which fetal sensor socket you use, as the monitor allocates a channel automatically. For instance, the first transducer you connect is automatically allocated a channel, and the measurement is labeled **FHR1**, the second is labeled **FHR2**, and the third **FHR3**. If you need to disconnect the transducers measuring the FHR temporarily, with the intention to continue monitoring after the temporary break (for example, if the mother needs to go to the bathroom), it is important that you reconnect the transducers in the same order as you originally connected them to make sure the measurement labels remain consistent. Upon discharge of the patient all connected transducers are reset from left to right. Example: Only one transducer is still connected to the monitor, it was labeled **FHR2** while monitoring the previous discharged patient; it is now reset to **FHR1**.
- The transducer Finder LED lets you identify at a glance which transducer is monitoring which heart rate channel.

- The fetal sensor socket to which a transducer is connected, is identified by the transducer position indicator in the blue setup menu header:

FM20/FM30



FM40/FM50



- The trace recorded for the **FHR3** is thicker (darker) than that recorded for **FHR1**, which is thicker than that for **FHR2**. This ensures that the three heart rates are easily distinguishable. The thickness of the recorded trace can be changed in Configuration Mode.
- Remember that only one fetal heartbeat can be heard from the loudspeaker at any time. The audio source symbol shows you which fetus you are listening to. To hear the other fetal heartbeat, select the fetal heart rate sound for this channel (see “Selecting Fetal Heart Sound” on page 158).



- Monitor maternal pulse to avoid mistaking maternal heart rate for FHR.
- Ensure you are recording the best possible signals by referring to the signal quality indicators and repositioning the transducers if necessary.

For the Avalon CL transducer system see “Cableless Monitoring - Important Considerations” on page 156.

Monitoring Triplets

To monitor triple FHRs you need three ultrasound transducers. Follow the procedures described in “Monitoring FHR and FMP Using Ultrasound” on page 155 and in “Monitoring Twin FHRs” on page 171. The transducer Finder LED lets you identify at a glance which transducer is monitoring which heart rate channel.

Separating FHR Traces

To help you to interpret traces with similar baselines, you can separate the baselines by an offset of 20 bpm by switching trace separation on. For details of the offset, see “Separation Order Type” on page 174.

"Standard" Separation Order

To make differentiating the traces easier, the trace for FHR2 is offset by **+20** bpm, and the trace for FHR3 is offset by **-20** bpm. In other words, the trace for FHR2 is recorded 20 bpm higher than it really is, while the trace for FHR3 is recorded 20 bpm lower than it really is. The trace for FHR1 is never shifted.



The recorder prints a dotted line labeled **+20** across the FHR scale, to identify the trace for FHR2.

-20

The recorder prints a dotted line labeled **-20** across the FHR scale, to identify the trace for FHR3.

+20

The FHR trace is labeled every 5 cm.


-20

+20

The label for FHR2 is annotated with **+20** and the FHR3 label is annotated with **-20**.

-20

The following trace shows triplets with **Trace Separation** on, and using **Standard** separation order.



The traces for FHR2 and FHR3 are offset. The numerical FHR values displayed on the monitor remain unchanged. Subtract 20 from the recorded trace for FHR2 to obtain the true FHR2. For example, if the recorded trace shows 160, then the true FHR is 140. Similarly, add 20 to the recorded trace for FHR3 to obtain the true FHR3.

"Classic" Separation Order

To make differentiating the traces easier, the trace for FHR1 is offset by **+20** bpm, and the trace for FHR3 is offset by **-20** bpm. In other words, the trace for FHR1 is recorded 20 bpm higher than it really is, while the trace for FHR3 is recorded 20 bpm lower than it really is. The trace for FHR2 is never shifted.

+20

The recorder prints a dotted line labeled **+20** across the FHR scale, to identify the trace for FHR1.

-20

The recorder prints a dotted line labeled **-20** across the FHR scale, to identify the trace for FHR3.

+20

The FHR trace is labeled every 5 cm.


-20

+20

The label for FHR1 is annotated with **+20** and the FHR3 label is annotated with **-20**.

-20

The following trace shows triplets with **Trace Separation** on, and using **Classic** separation order.



The traces for FHR1 and FHR3 are offset. The numerical FHR values displayed on the monitor remain unchanged. Subtract 20 from the recorded trace for FHR1 to obtain the true FHR1. For example, if the recorded trace shows 160, then the true FHR is 140. Similarly, add 20 to the recorded trace for FHR3 to obtain the true FHR3.

Switching Trace Separation On and Off

- 1 Connect three ultrasound transducers to the monitor to measure FHR.
- 2 See “Switching Trace Separation On and Off” on page 173 for details of how to switch trace separation on or off.


When Trace Separation is On

When trace separation is turned on, the recorder prints a dotted line labeled with the three FHRs at the top, and ± 20 at the bottom. Examples of the two methods (**Standard**, **Classic**) for determining the trace separation order are provided here.


± 20

When Trace Separation is Off

To indicate that trace separation is switched off, a dotted line labeled **+0** prints across the FHR scale.



- 1 **Standard** trace separation switched off here




2 **Classic** trace separation switched off here

Troubleshooting

Common problems that may occur when monitoring FHR using ultrasound are listed in “Monitoring FHR and FMP Using Ultrasound” on page 155.

The following problem may occur when monitoring triplets.

Problem	Possible Cause	Solution
 <p>The question mark is printed repeatedly, and appears on the screen and the INOP Coincidence is issued.</p>	<p>More than one transducer is recording the same FHR, or a fetal transducer records the same heart rate as the MHR.</p>	<p>Reposition one or more ultrasound transducer, as appropriate. See “Recommended Actions for INOP Coincidence” on page 154</p>

For more information, see “Additional Information” on page 163.

Fetal Heart Rate Alarms

Fetal heart rate (FHR) alarms can give both audible and visual warning of a non-reassuring fetal condition. **Your monitor must be configured to alarm mode All to enable the FHR alarms** (see “Alarms” on page 109).

Changing Alarm Settings

When you do any of the following actions for any FHR measurement channel, this applies for all active FHR measurements, both ultrasound and DECG:

- Turning FHR alarms on or off.
- Changing alarm limits.
- Changing alarm delays.
- Changing signal loss delay.

The monitor retains these settings, even when switched off. The alarm limits are printed on the trace every few pages if alarms are on.

Turning Alarms On or Off

- 1 Connect either an ultrasound or a DECG transducer to a free socket on the monitor.
- 2 Enter the Setup menu for a connected FHR measurement.
- 3 Select **Alarms** to toggle between **On** and **Off**.

Changing Alarm Limits

- 1 Connect either an ultrasound or a DECG transducer to a free socket on the monitor.
- 2 Enter the Setup menu for a connected FHR measurement.
- 3 To change the high alarm limit, select **High Limit** and select the alarm limit from the pop-up list.
- 4 To change the low alarm limit, select **Low Limit** and select the alarm limit from the pop-up list.

Changing Alarm Delays

You can change the alarm delays if the **Alarm Mode** is set to **All**.

- 1 Connect either an ultrasound or a DECG transducer to a free socket on the monitor.
- 2 Enter the Setup menu for a connected FHR measurement.

- 3 To change the high alarm limit delay time, select **High Delay** and select the delay time (in seconds) from the pop-up list.
- 4 To change the low alarm limit delay time in seconds, select **Low Delay** and select the delay time (in seconds) from the pop-up list.

Changing Signal Loss Delay

The signal loss delay is the configurable delay before a **Signal Loss** INOP is issued. You can change the delay:

- 1 Connect either an ultrasound or a DECG transducer to a free socket on the monitor.
- 2 Enter the Setup menu for a connected FHR measurement.
- 3 Select **SignalLoss Delay** and select the signal loss INOP delay time (in seconds) from the pop-up list.

Monitoring Uterine Activity Externally

You can measure uterine activity externally using a Toco transducer. You can also use a Toco⁺, Toco MP, or a CL Toco⁺ MP transducer for the same purpose, although they also have more (ECG/IUP and Pulse) capabilities.

The external Toco transducer measures the frequency, duration and relative strength of contractions, but not their absolute intensity. Amplitude and sensitivity depend on various factors such as the position of the transducer, the belt tension and the size of the patient.

What You Need



Toco, Toco MP, or Toco⁺ transducer



CL Toco⁺ MP transducers (additionally capable of providing the maternal pulse measurement)



CL ECG/IUP transducer




Abdominal Belt (disposable shown)

External Toco Monitoring

Prepare for Toco monitoring using the list below. The standard procedures in use in your facility determine the sequence of actions.

- 1 Fasten the abdominal transducer belt around the patient.
- 2 Connect the Toco transducer to a free socket on the monitor. The Toco baseline is automatically reset. The Toco display shows 20. "**Toco**", indicating external uterine measurement, is printed on the trace at intervals.
- 3 Place the transducer on the patient's fundus to ensure the optimum recording of uterine activity.
- 4 Reset the Toco baseline as necessary (see "Resetting the Toco Baseline" on page 188), but not during a contraction.

The following example trace shows two contractions.



Resetting the Toco Baseline

- 1 Press the **Toco Baseline** SmartKey. This resets the Toco baseline to 20 on the display and trace.
- 2 If the Toco value is negative for more than five seconds, the Toco baseline is automatically reset to 0 units.

Toco Sensitivity

If the Toco sensitivity is too high, and the Toco trace exceeds the paper scale, you can reduce the Toco sensitivity to 50%. The default setting is 100%.

To change the Toco sensitivity:

- 1 Enter the **Setup Toco** menu.
- 2 Select **Gain** to toggle between **100%** and **50%**.

Troubleshooting

External Toco Monitoring		
Problem	Possible Causes	Solutions
Quality of the trace deteriorates or the Toco baseline varies	The belt is incorrectly fastened and is too slack or too tight, or the belt has lost its elasticity.	The belt must be tight enough to ensure good contact between the patient's skin and the entire surface of the transducer without causing discomfort. Ensure you are using the correct belt. Adjust it as necessary.
	Fetal movement	Check if the belt is correctly fastened and adjust it as necessary. Reposition the transducer and reset the Toco baseline if necessary.
	Maternal respiration superimposed on trace	Check if belt is not too loose.
	Maternal movement/change of position	Following maternal movement, reset Toco baseline
Toco sensitivity is too high (above 100 units) Toco trace is exceeding the paper scale	Physical transmission of pressure from the uterus to the sensor is much higher than the average value.	Check if the belt is too tightly fastened. Select 50% Toco sensitivity.
Toco Equip Malfunction is displayed.		See "Patient Alarms and INOPs" on page 121.
Toco Unplugged is displayed.		See "Patient Alarms and INOPs" on page 121.
If you suspect the signal from the transducer.		Test the Transducer (see "Testing Toco Transducers" on page 189 below).

Testing Toco Transducers

If any of the following tests fail, repeat the test using another transducer. If the second transducer passes the tests, confirming that the first transducer is defective, contact your service personnel.

If the second transducer also fails the tests, contact your service personnel.



To test a Toco transducer:

- 1 Switch on the monitor and the recorder.
- 2 Connect the transducer to the fetal monitor.
- 3 Gently apply pressure to the pick-up button.
- 4 Check if after a few seconds the values on the display and paper show this change in pressure.


You can test all Toco transducers, including the cableless ones, as described above.

Monitoring Uterine Activity Internally

FM30/50 You can monitor intrauterine pressure (IUP) using an intrauterine catheter together with a patient module, Toco⁺, or the CL ECG/IUP transducer, after rupture of the membranes and the cervix is sufficiently dilated.

What You Need


Illustration 1 shows the complete connection chain from the IUP catheter to the fetal monitor using the patient module:



- 1 Disposable Koala IUP Catheter (M1333A)
- 2 Reusable Koala IUP Adapter Cable (9898 031 43931)
- 3 Patient Module (M2738A)


16 Monitoring Uterine Activity Internally

Illustration 2 shows the complete connection chain from the IUP catheter to the fetal monitor using the Toco⁺ transducer:



- 1 Disposable Koala IUP Catheter (M1333A)
- 2 Reusable Koala IUP Adapter Cable (9898 031 43931)
- 3 Toco⁺ transducer (M2735A)

Illustration 3 shows the complete connection chain from the IUP catheter to the fetal monitor using the CL ECG/IUP transducer:



- 1 Disposable Koala IUP Catheter (M1333A)
- 2 Reusable Koala IUP Adapter Cable (9898 031 43931)
- 3 CL ECG/IUP transducer (866077)
- 4 Avalon CL base station (866074)

Internal (IUP) Monitoring

Read the instructions that accompany the intrauterine catheter and the adapter cable before you start monitoring. Zero the monitor when instructed.

WARNING

Do not catheterize if placenta previa is diagnosed, or if uterine bleeding from an undetermined source is present.

Prepare for IUP monitoring using the list below. The standard procedures in use in your facility determine the sequence of actions.

- 1 Perform a complete clinical evaluation.
- 2 Catheterize after membrane rupture. Insert the catheter according to its accompanying instructions.
- 3 Connect the catheter to the socket on the patient module.
- 4 Connect the patient module to a free socket on the monitor. The monitor is automatically zeroed. The IUP display shows 0. "IUP", indicating internal measurement, is printed at intervals on the trace.
- 5 Zero the monitor (see "Zeroing the Monitor" on page 193).
- 6 If you suspect the catheter is not responding appropriately, flush as directed in the catheter's Instructions for Use. A pressure spike appears on the trace if you flush after connecting the transducer to the monitor.

Zeroing the Monitor

Zero the monitor by selecting the **Zero IUP** SmartKey, or selecting **Zero IUP** in the **Setup IUP** menu. This resets the display and trace to 0. If you do not zero the monitor properly, the pressure trace may exceed the paper scaling.



Selecting the IUP Unit

You can select between mmHg (default) and kPa for the IUP unit.

- 1 Enter the **Setup IUP** menu.
- 2 Press **Unit** to toggle between mmHg and kPa.

Troubleshooting

Internal (IUP) Monitoring		
Problem	Possible Causes	Solutions
No change in pressure during contraction.	Dry environment or possible extra-ovular placement of sensor tip.	Refer to catheter Instructions for Use.
Only pressure peaks can be seen (baseline not visible).	Zero adjustment is incorrect.	Zero the system.
Trace is a straight line.	Transducer is defective.	Remove and touch the catheter. If the trace does not show up and down movements, use a new transducer.
	Catheter blocked.	Refer to catheter Instructions for Use.
Trace is superimposed with noise.	End of catheter is in the uterine wall.	Refer to catheter Instructions for Use.
IUP Equip Malf INOP is displayed.		See “Patient Alarms and INOPs” on page 121.
IUP Unplugged INOP is displayed.		

Monitoring FHR Using DECG

FM30/50 This chapter describes how to monitor a single fetal heart rate via direct ECG (DECG), using a spiral fetal scalp electrode in the intrapartum period.

Read and adhere to the instructions that accompany the fetal scalp electrode, the DECG adapter cable, and the attachment electrode. Pay attention to all the contraindications, warnings, and for the DECG adapter cable, the cleaning and disinfection procedures.

Before starting to monitor, first define the fetal position, and ensure that it is suitable for DECG monitoring.

Misidentification of MHR as FHR

Confirm fetal life before monitoring, and continue to confirm that the fetus is the signal source for the FHR during monitoring. Here are two examples where the MHR can be misidentified as the FHR when using a fetal scalp electrode:

- Electrical impulses from the maternal heart can sometimes be transmitted to the fetal monitor through a recently deceased fetus via the spiral scalp electrode, appearing to be a fetal signal source.
- The recorded MHR, and any artifact, can be misinterpreted as a FHR especially when it is over 100 bpm.

To reduce the possibility of mistaking the MHR for FHR, monitor both maternal and fetal heart rates (see “Monitoring Maternal Heart / Pulse Rate” on page 223). The monitor's cross-channel verification (CCV) facility can help by automatically detecting when the same heart rate is being recorded by different transducers. See “Confirm Fetal Life Before Using the Monitor” on page 10 and “Cross-Channel Verification (CCV)” on page 149.

If the **Coincidence** INOP is issued at the fetal monitor if you are measuring FHR with DECG:

- 1 Confirm that the scalp electrode is placed correctly.
- 2 Confirm fetal life by palpation of fetal movement or auscultation of fetal heart sounds using a fetoscope, stethoscope, or Pinard stethoscope.
- 3 If you cannot hear the fetal heart sounds, and you cannot confirm fetal movement by palpation, confirm fetal life using obstetric ultrasonography.
- 4 In case of difficulties deriving a stable maternal pulse reading using the Toco MP or CL Toco⁺ MP transducer, use SpO₂ instead. In case of similar problems with the pulse measurement from SpO₂, use MECCG instead. Reasons to switch the method for deriving a maternal pulse or heart rate include: motion artifacts, arrhythmia, and individual differences in pulse signal quality on the abdominal skin (via Toco⁺ MP).


What You Need

You can measure fetal DECG using the equipment combinations shown in the following illustrations.

WARNING


Never attempt to connect the fetal scalp electrode to anything other than the correct DECG adapter cable.

Illustration 1 shows the complete connection chain from the fetal scalp electrode to the fetal monitor using the Toco⁺ transducer.




- 1 Fetal Scalp Electrode, single spiral (9898 031 37631)
- 2 Fetal Scalp Electrode, double spiral, Europe only, not for USA (9898 031 37641)
- 3 DECG Adapter Cable (9898 031 37651) with Pre-gelled Attachment Electrode (9898 031 39771)
- 4 Toco+ transducer (M2735A)

Illustration 2 shows the equivalent chain using the CL Toco⁺ MP or CL ECG/IUP transducer.



- 1 Fetal Scalp Electrode, single spiral (9898 031 37631)
- 2 Fetal Scalp Electrode, double spiral, Europe only, not for USA (9898 031 37641)
- 3 DECG Adapter Cable (9898 031 37651) with Pre-gelled Attachment Electrode (9898 031 39771)
- 4 CL Toco⁺ MP or CL ECG/IUP transducer (866077)
- 5 Avalon CL base station (866074)

Illustration 3 shows the equivalent chain using the patient module.



- 1 Fetal Scalp Electrode, single spiral, (9898 031 37631)
- 2 Fetal Scalp Electrode, double spiral, Europe only, not for USA (9898 031 37641)
- 3 DECG Adapter Cable (9898 031 37651) with Pre-gelled Attachment Electrode (9898 031 39771)
- 4 Patient Module (M2738A)

Making Connections

WARNING

Follow the instructions supplied with each of the monitoring accessories you are using.

Prepare for DECG monitoring using the list below. The standard procedures in use in your facility determine the sequence of actions.

- 1 If changing monitoring mode from US to DECG, first disconnect the US transducer.
- 2 Depending on the equipment you are using, ensure that the Toco⁺ transducer, CL Toco⁺ MP, the CL ECG/IUP transducer, or the patient module is connected to the fetal monitor.
- 3 Attach the fetal scalp electrode to the fetus, following the instructions supplied with the fetal scalp electrode.
- 4 Attach a pre-gelled attachment electrode to the DECG adapter cable, following the instructions supplied with the DECG adapter cable.
- 5 Fix the attachment electrode to the mother's thigh, following the instructions supplied with the attachment electrode.
Depending on the equipment you are using, connect the red connector plug on the DECG adapter cable to the red connector on the Toco⁺ transducer, CL Toco⁺ MP, the CL ECG/IUP transducer, or the patient module.
- 6 Connect the fetal scalp electrode to the DECG adapter cable.

You are now ready to begin monitoring DECG.

WARNING

The fetal/maternal monitor is not a diagnostic ECG device. In particular, the display of fetal/maternal ECG is intended only for evaluating signal quality for fetal/maternal heart rate as derived from the ECG waveform.

When in doubt, it can be used to identify sources of compromised signal quality, such as noise or muscle artifacts. It can subsequently be used to verify the result of measures taken to resolve them (e.g. checking ECG cable connections or adapting the fetal **ArtifactSuppress** configuration).


The safety and effectiveness of the displayed fetal/maternal ECG waveform (i.e. P, QRS and T segments) for evaluation of fetal/maternal cardiac status during labor have not been evaluated.

Monitoring DECG

To simultaneously measure DECG and MEGC, you need the CL ECG/IUP transducer or the patient module for DECG, and a Toco+, CL Toco⁺ MP, or CL ECG/IUP transducer for MEGC (see “Monitoring Maternal Heart / Pulse Rate” on page 223). Alternatively, you can monitor the maternal pulse rate via pulse oximetry (see “Pulse Rate from SpO₂” on page 230). You can also monitor maternal pulse with the Toco MP or CL Toco⁺ MP transducer. In any case where you would use a Toco⁺ or Toco MP transducer, you can also monitor with a CL Toco⁺ MP transducer.

- 1 Switch on the recorder.
- 2 The heart rate monitored via DECG is labeled **DFHR1** / **DFHR2** / **DFHR3** on the screen. If configured, the DECG wave is displayed automatically on the screen, labeled **DECG** and fetal. If MEGC is being monitored, both waves are displayed, with the DECG wave above the MEGC wave. The MEGC wave is labeled **MECG** and maternal.

- 3 Check the artifact suppression setting and change it if necessary (see “Suppressing Artifacts” on page 201).



- 1 Measurement label (FHR1)
- 2 Measurement label (Toco)
- 3 Measurement label (DFHR2)
- 4 Measurement label (FHR3)
- 5 Measurement label (SpO₂)
- 6 Measurement label (HR)
- 7 1mV scale bar
- 8 MECG wave with maternal label
- 9 DECG wave with fetal label
- 10 Measurement label (Temp) maternal temperature
- 11 Measurement label NBP

NOTE

The 1mV scale bar for the DECG and MECG wave is not displayed on the screen if you monitor DECG or MECG with an Avalon CTS system. The Avalon CTS system does not provide a scaled ECG.

WARNING

Periodically compare the mother's pulse with the signal coming from the monitor's loudspeaker to ensure that you are monitoring fetal heart rate. If the MHR coincides with the FHR, do not misinterpret the MHR as the FHR (see also “Confirm Fetal Life Before Using the Monitor” on page 10 and “Cross-Channel Verification (CCV)” on page 149).

Suppressing Artifacts

When the monitor's artifact suppression is on, instantaneous heart rate changes of 28 bpm or more, however caused, are not recorded. Fetal arrhythmia will also be suppressed. **If you suspect fetal arrhythmia, switch artifact suppression off.** When artifact suppression is off, all recorded fetal heartbeats within the specified range are shown. The default setting is **On** (artifacts are suppressed).

To change the setting:

- 1 Enter the **Setup DFHR1** menu.
- 2 Select **ArtifactSuppress** to toggle between artifact suppression **On** (artifacts are suppressed) and **Off** (no artifact suppression, use this setting if you suspect fetal arrhythmia).

When artifact suppression is off, **Artifact Suppression Off** is annotated on the trace recording.

Printing the Waveform

You can print the DECG wave onto the trace paper. Refer to “Printing the ECG Waveform” on page 233.

Troubleshooting

Problem	Possible Cause	Solutions
DFHR1 Equip Malf DFHR2 Equip Malf DFHR3 Equip Malf	Malfunctioning equipment	See “Patient Alarms and INOPs” on page 121.
DFHR1 Leads Off DFHR2 Leads Off DFHR3 Leads Off Numeric is displayed with a -?;- INOP tone See also “Patient Alarms and INOPs” on page 121. Erratic trace Erratic display	Fetal scalp electrode detached at connector Poor or no contact between leg attachment electrode and mother No contact between the DECG adapter cable and the leg attachment electrode No contact between the fetal scalp electrode connector and the DECG adapter cable No ECG signal Poor contact between the reference electrode and the mother	Reconnect the fetal scalp electrode If the wave is configured to be displayed on the monitor, you can observe if the ECG signal is clear, or if it shows interruptions and noise Check all connections Disconnect and reconnect the connector several times Check all connections Disconnect and reconnect the connector several times If problem persists, use a new fetal scalp electrode Check for fetal demise Use a new fetal scalp electrode if necessary
Signal quality indicator continuously shows a poor signal	Fetal arrhythmia	Use a new fetal scalp electrode if necessary
DFHR1 Signal Loss DFHR2 Signal Loss DFHR3 Signal Loss	No signal	See “Patient Alarms and INOPs” on page 121.
DFHR1 Unplugged DFHR2 Unplugged DFHR3 Unplugged	No connection	See “Patient Alarms and INOPs” on page 121.

Testing DECG Mode

See the monitor's Service Guide.

Monitoring Noninvasive Blood Pressure

This fetal monitor uses the oscillometric method for the noninvasive blood pressure measurement. In adult mode the blood pressure measurements determined with this device comply with the American National Standard for Electronic or Automated Sphygmomanometers (ANSI/AAMI SP10-1992) in relation to mean error and standard deviation, when compared to intra-arterial or auscultatory measurements (depending on the configuration) in a representative sample population. For the auscultatory reference, the fifth Korotkoff sound was used to determine the diastolic pressure.

A physician must determine the clinical significance of the measurement information.

Introducing the Oscillometric Noninvasive Blood Pressure Measurement

Oscillometric devices measure the amplitude of pressure changes in the occluding cuff as the cuff deflates from above systolic pressure. The amplitude suddenly increases as the pulse breaks through the occlusion in the artery. As the cuff pressure decreases further, the pulsations increase in amplitude, reach a maximum (which approximates to the mean pressure), and then diminish.

Studies show that, especially in critical cases (arrhythmia, vasoconstriction, hypertension, shock), oscillometric devices are more accurate and consistent than devices using other noninvasive measuring techniques.

WARNING

Intravenous infusion: Do not use the NBP cuff on a limb with an intravenous infusion or arterial catheter in place. This could cause tissue damage around the catheter when the infusion is slowed or blocked during cuff inflation.

Skin Damage: Do not measure NBP in cases of sickle-cell disease or any condition where skin damage has occurred or is expected.

Unattended measurement: Use clinical judgment to decide whether to perform frequent unattended blood pressure measurements in cases of severe blood clotting disorders, because of the risk of hematoma in the limb fitted with the cuff.

Measurement Limitations

Measurements are impossible with heart rate extremes of less than 40 bpm or greater than 300 bpm, or if the patient is on a heart-lung machine.

The measurement may be inaccurate or impossible:

- with excessive and continuous patient movement such as during contractions
- if a regular arterial pressure pulse is hard to detect
- with cardiac arrhythmias

- with rapid blood pressure changes
- with severe shock or hypothermia that reduces blood flow to the peripheries
- with obesity, where a thick layer of fat surrounding a limb dampens the oscillations coming from the artery
- on an edematous extremity.

Measurement Methods

There are three measurement methods:

- **Manual** - measurement on demand. Results are displayed for up to one hour.
- **Auto** - continually repeated measurements (between one and 120 minute adjustable interval). You can make a manual measurement between two measurements in Auto Mode.
- **Sequence**- up to four measurement cycles which will run consecutively, with number of measurements and interval between them configurable for each cycle.

Reference Method

The measurement reference method can be Auscultatory (manual cuff) or Invasive (intra-arterial). For further information, see the Application Note supplied on the monitor documentation DVD.

To check the current setting, select **Main Setup, Measurements, NBP**, and check whether the **Reference** setting is set to **Auscultatory** or **Invasive**. This setting can be changed in Configuration Mode.

Preparing to Measure Noninvasive Blood Pressure

If possible, avoid taking measurements during contractions because the measurement may be unreliable and may cause additional stress for the patient.

- 1 Connect the cuff to the air tubing.
- 2 Plug the air tubing into the red NBP connector. Avoid compression or restriction of pressure tubes. Air must pass unrestricted through the tubing.
- 3 Make sure that you are using a Philips-approved correct sized cuff and that the bladder inside the cover is not folded or twisted.

A wrong cuff size, and a folded or twisted bladder, can cause inaccurate measurements. The width of the cuff should be in the range from 37% to 47% of the limb circumference. The inflatable part of the cuff should be long enough to encircle at least 80% of the limb.

- 4 Apply the cuff to a limb at the same level as the heart. If it is not, you must use the measurement correction formula to correct the measurement.


The marking on the cuff must match the artery location. Do not wrap the cuff too tightly around the limb. It may cause discoloration, and ischemia of the extremities. Inspect the application site regularly to ensure skin quality and inspect the extremity of the cuffed limb for normal color, warmth and sensitivity. If the skin quality changes, or if the extremity circulation is being affected, move the cuff to another site or stop the blood pressure measurements immediately. Check more frequently when making automatic or stat measurements.

Correcting the Measurement if Limb is not at Heart Level

To correct the measurement if the limb is not at heart level, to the displayed value

Add 0.75 mmHg (0.10 kPa) for each centimeter higher or	Deduct 0.75 mmHg (0.10 kPa) for each centimeter lower or
Add 1.9 mmHg (0.25 kPa) for each inch higher.	Deduct 1.9 mmHg (0.25 kPa) for each inch lower.

Understanding the Numerics



Depending on the numeric size, not all elements may be visible. Your monitor may be configured to display only the systolic and diastolic values. If configured to do so, the pulse from NBP will display with the NBP numeric.

Aging Numerics

The measured NBP value, together with the corresponding pulse rate if this is switched on, will be displayed for one hour. After that the values are regarded as invalid and are no longer displayed. During this time, measurement values may be grayed out or disappear from the screen after a set time, if configured to do so. This avoids older numerics being misinterpreted as current data. The time can be set in Configuration Mode. In **Auto** Mode the measurement values may disappear more quickly (to be replaced by new measurement values), if the repeat time is set to less than one hour.

NBP and SpO₂ numerics from intermittent measurements can be configured to be grayed out or to disappear from the screen after a set time. This avoids older numerics being misinterpreted as current data. The time can be set individually for SpO₂ and NBP in Configuration Mode.

Alarm Sources

If you have parallel alarm sources, the sources are displayed instead of the alarm limits.

NBP Measurement Timestamp

Depending on your configuration, the time shown beside the NBP numeric can be:

- the time of the most recent NBP measurement, also known as the "timestamp", or
- the time until the next measurement in an automatic series, displayed with a graphic representation of the remaining time, as shown here.



- The NBP timestamp will normally show the **completion** time of the NBP measurement.

NBP Measurement Start Time






In **Auto** or **Sequence** mode the monitor is configured to synchronize the measurements in a measurement series to an "easy-to-document" time. For example, if you start the first measurement at 08:23, and the **Repeat Time** is set to **10 min**, the monitor automatically performs the next measurement at 8:30, then 8:40 and so on, unless it has been configured to **NotSynchron..**

During Measurements

The cuff pressure is displayed instead of the units and the repeat time. An early systolic value gives you a preliminary indication of the systolic blood pressure during measurement.

Starting and Stopping Measurements

Use the Setup menu or the SmartKeys to start and stop measurements.

Action to be performed	Setup menu	SmartKeys
Start/Stop manual measurement	Start/Stop	
Start Auto series		
Stop current automatic measurement		Start/ Stop
Start manual measurement	-	
Start Auto series		Start
Stop manual measurement	-	
Stop current automatic measurement		Stop
Stop automatic, or manual measurement AND series	Stop All NBP	 Stop All

Enabling Automatic Mode and Setting Repetition Time

- 1 In the **Setup NBP** menu, select **Mode**.
- 2 Toggle between **Auto** and **Manual**, if necessary, to pick the measurement method.
- 3 If making an automatic measurement, select **Repeat Time** or press the **Repeat Time** SmartKey and set the time interval between two measurements.



NOTE

Be aware that a combination of a recorder speed of less than 3 cm/min and a repetition time of less than five minutes can result in not all noninvasive blood pressure measurements being recorded on the fetal trace. For example, if the recorder speed is set to 1 cm/min and the repetition time is set to two minutes, due to the low speed setting, the recorder will only be able to record every other noninvasive blood pressure measurement. This affects only the local fetal trace recording, and all measurements are displayed as normal on the monitor's screen.

Enabling Sequence Mode and Setting Up The Sequence

- 1 In the **Setup NBP** menu, select **Mode** and select **Sequence** from the pop-up menu.
- 2 Select **Setup Sequence** to open the **Setup Sequence** window.
Up to four measurement cycles can be setup which will run consecutively. For each cycle you can set the number of measurements and the interval between them. If you want to run less than four cycles in a sequence, you can set the number of measurements for one or more cycles to **Off**.
- 3 Select each sequence in turn, and select the number of measurements and the time interval between the measurements.
- 4 To have measurements continue after the sequence, set the number of measurements for your last cycle to **Continuous** and this cycle will run indefinitely.

CAUTION

Be aware that, if none of the cycles are set to **Continuous**, NBP monitoring will end after the last measurement of the cycle.

When the NBP measurement mode is set to **Sequence**, the repetition time for **Auto** Mode cannot be changed.

Choosing the Alarm Source

You can monitor for alarm conditions in systolic, diastolic, and mean pressure, either singly or in parallel. Only one alarm is given, with the priority of mean, systolic, diastolic.

Menu option	Pressure value monitored
Sys.	systolic
Dia.	diastolic
Mean	mean
Sys & Dia	systolic and diastolic in parallel
Dia & Mean	diastolic and mean in parallel
Sys & Mean	systolic and mean in parallel
Sys&Dia&Mean	all three pressures in parallel


If mean is not selected as alarm source (**Sys.**, **Dia.**, or **Sys & Dia** selected), but the fetal monitor can only derive a mean value, mean alarms will nevertheless be announced using the most recent mean alarm limits. Check that the mean alarm limits are appropriate for the patient, even when not using mean as the alarm source. When no value can be derived an **NBP Measure Failed** INOP will be displayed.

Assisting Venous Puncture

You can use the cuff to cause sub-diastolic pressure. The cuff deflates automatically after a set time if you do not deflate it.

- 1 In the **Setup NBP** menu select **VeniPuncture**.
- 2 Puncture vein and draw blood sample.
- 3 Reselect **VeniPuncture** to deflate the cuff.

During measurement, the display shows the inflation pressure of the cuff and the remaining time in venous puncture mode.



- 1 Cuff pressure
- 2 Venous puncture measurement mode
- 3 Time left in venous puncture mode

Calibrating NBP

NBP is not user-calibrated. NBP pressure transducers must be verified at least once every two years by a qualified service professional, and calibrated, if necessary. See the Service Guide for details.

Troubleshooting

Problem	Possible Causes	Solutions
Cuff will not inflate	Monitor is in Service or Configuration Mode	
	Technical defect	Call service
	Cuff tubing not connected	Connect cuff tubing
High or low values measured (against clinical expectations)	Contraction occurring	Wait until contraction has finished
	Patient talking before or during measurement	Allow patient to rest quietly, then try again after three to five minutes
	Incorrect cuff size or cuff not at heart level	Check cuff size, level, and position
	Noninvasive blood pressure reference method set incorrectly	Check the reference method configured (auscultation or intra-arterial) and correct if necessary in Configuration Mode
Displays zeros for systolic and diastolic values. Measurement automatically repeats	Severe vasoconstriction at cuff site	Move cuff to another limb, check for shock, or verify blood pressure using another method
	Erratic blood pressure fluctuations due to arrhythmias or rapid-acting drugs or contractions	Try again, if unsuccessful, verify blood pressure using another method Wait until contraction has finished.
	Excessive patient movement or convulsions	Restrain movement or verify blood pressure using another method
NBP Cuff Overpress INOP is displayed		See “Patient Alarms and INOPs” on page 121.
NBP Equip Malf INOP is displayed		
NBP Interrupted INOP is displayed		
NBP Measure Failed		

Monitoring Maternal Temperature

Measuring Tympanic Temperature

The tympanic thermometer measures the patient's temperature in the ear using infrared technology.



The result of this measurement can be automatically adjusted to correspond to a different body reference site. The result is displayed on the screen of the thermometer and transmitted to the monitor.

The thermometer is used with single-use probe covers for infection control during measurement.

WARNING


Do not use in the presence of flammable anesthetics, such as a flammable anesthetic mixture with air, oxygen or nitrous oxide.







Place the thermometer into its base station, when it is not in use. The base station allows flexible mounting of the thermometer at the point of care. The base station is connected to the monitor's MIB/RS232 interface (optional) via a cable. It has storage space for up to 32 probe covers.

Thermometer Display and Controls




The tympanic thermometer has a liquid crystal display. The display shows the patient's temperature in numerics and guides you with symbols through the measurement process.



Functional Keys	Description of Use
	Press the eject key to eject the probe cover. The eject key symbol is shown on the display when a measurement has been taken and transmitted.
	Press the change unit key after a measurement to switch between °C and °F.
	Press the start measurement key when you are ready to take a patient's temperature.
	The pulse timer key can be used to time vital signs you take manually. Press and hold the pulse timer key to enter timer mode. Press the pulse timer key again to start the timer. The thermometer will issue a one beep at 15 seconds, two beeps at 30 seconds, three beeps at 45 seconds, and four beeps at 60 seconds.

Status Screens

The thermometer performs an internal test at every start-up to verify that the system components are functioning properly. It measures the ambient temperature. During start-up and measurement, the thermometer screens communicate the current status.

Images	Description
	Ambient Temperature above specified range
	Ambient Temperature below specified range
System errors	
	<p>System error 12 - there is a problem with the settings. Contact your service personnel to have them check the settings and reset them, if necessary.</p> <p>If the display shows any other system error, then reset the thermometer by picking up a probe cover. If the system error does not clear, contact your service personnel.</p>

Making a Temperature Measurement

WARNING

Inaccurate measurement results can be caused by:

- incorrect application of the thermometer.
- anatomical variations in the ear.
- build up of earwax in the ear.
- excessive patient movement during the measurement.
- absent, defective or soiled probe covers.
- probe covers other than the specified probe covers.

- 1 Ensure that the base station is connected via the appropriate cable to the connector on the monitor.
- 2 Remove the thermometer from the base station.
The thermometer is latched to the base station to avoid an accidental fall when the base station is moved. To pick up the thermometer from the base station move it slightly up and then lift it from the base station to release the latch.
- 3 Press the eject key on the thermometer to discard any probe that may have been left on the thermometer from a previous use.

19 Monitoring Maternal Temperature

- 4 Pick up a new probe from the probe container on the base station.
- 5 Inspect the probe cover to make sure it is fully seated (no space between cover and tip base) and there are no holes, tears, or wrinkles in the plastic film.
- 6 Place the thermometer with the probe in the ear canal, sealing the opening with the probe tip. For consistent results, ensure that the probe shaft is aligned with the ear canal.
- 7 Press and release the start-up key gently.
- 8 Wait until you hear the three beeps.
- 9 Remove the probe from the ear.

The temperature values are displayed both on the thermometer itself and on the connected monitor.



- 10 Check that the correct temperature label for the measurement site is displayed: **iTrect**, **iToral**, **iTcore**, or **iTtymp**
- 11 Press the eject key to eject the probe cover into a suitable waste receptacle.
- 12 Return the thermometer to the base station.

The thermometer switches to stand-by mode after 30 seconds when not used.

Possible INOPs

Images	Description
	Patient temperature above measurement range.
	Patient temperature below measurement range.

WARNING

- Never apply the probe to the patient when the probe is not connected to the unit.
- Always use a single-use probe cover to limit patient cross-contamination.
- Measurement errors or inaccurate readings may result when probe covers other than the specified probe covers are used (see “Tympanic Temperature Accessories” on page 267)
- Insert the probe slowly and carefully to avoid damage to the ear canal and the tympanic membrane.
- Inspect the probe for damage, holes, tears, or sharp edges to avoid injuring the skin.
- Always ensure that the used probe cover is removed before attaching a new probe cover.

CAUTION

- Do not immerse the probe in fluids, or drop fluids on the probe.
- Do not use a probe cover that has been dropped or is damaged.
- Do not autoclave. To prevent damage to the unit, probe and accessories, refer to the cleaning procedures in the “Cleaning and Disinfecting” on page 246 chapter.

Body Reference Sites and Monitor Labels

The tympanic thermometer measures the patient's temperature in the ear. The thermometer can be configured to adjust the result of the measurement to correspond to a different body reference site. The measurement label displayed on the monitor corresponds to the body reference site that is configured. The following body reference sites are available:

Body Reference Site	Label on Monitor
Ear temperature (no adjustment)	iTt ymp
Oral temperature	iToral
Core temperature	iTcore
Rectal temperature	iTrect

The body reference site can be selected in the Biomed mode of the tympanic thermometer. See the Service Guide for more information.

The measured maternal temperature is only transmitted to a connected OB TraceVue/IntelliSpace Perinatal, when the tympanic thermometer is configured to the **iTt ymp** body reference site, and the **iTt ymp** label is displayed on the fetal monitor screen and printed out on the trace paper.

Entering Temperature Manually

A temperature measurement can be entered manually.

- 1 Press the SmartKey **Enter Temp.**



- 2 A numeric pad opens.
- 3 Enter the temperature values.
- 4 Press the **Enter** key.

The aging time (time the entered value is valid) can be configured in the Configuration Mode.

Monitoring SpO₂

FM30/40/50 The pulse oximetry measurement (SpO₂) is intended for use with maternal patients.

Philips pulse oximetry uses a motion-tolerant signal processing algorithm, based on Fourier Artifact Suppression Technology (FAST). It provides two measurements:

- Oxygen saturation of arterial blood (SpO₂) - percentage of oxygenated hemoglobin in relation to the sum of oxyhemoglobin and deoxyhemoglobin (functional arterial oxygen saturation).
- Pulse rate - detected arterial pulsations per minute. This is derived from the SpO₂ value, and is one of four sources of the maternal heart/pulse rate used for cross-channel verification (see “Monitoring Maternal Heart / Pulse Rate” on page 223 and “Cross-Channel Verification (CCV)” on page 149).

Selecting an SpO₂ Sensor

See “Accessories and Supplies” on page 257 for a list of sensors, and the patient population and application sites for which they are appropriate.

Familiarize yourself with the Instructions for Use supplied with your sensor before using it.

CAUTION

Do not use OxiCliq disposable sensors in a high humidity environment, or in the presence of fluids, which may contaminate sensor and electrical connections causing unreliable or intermittent measurements. Do not use disposable sensors when there is a known allergic reaction to the adhesive.

Applying the Sensor

- 1 Follow the SpO₂ sensor's Instructions for Use, adhering to all warnings and cautions.
- 2 Remove colored nail polish from the application site.
- 3 Apply the sensor to the patient. The application site should match the sensor size so that the sensor can neither fall off, nor apply excessive pressure.
- 4 Check that the light emitter and the photodetector are directly opposite each other. All light from the emitter must pass through the patient's tissue.

WARNING

Proper Sensor Fit: If a sensor is too loose, it might compromise the optical alignment or fall off. If it is too tight, for example because the application site is too large or becomes too large due to edema, excessive pressure may be applied. This can result in venous congestion distal from the application site, leading to interstitial edema, hypoxemia and tissue malnutrition. Skin irritations or lacerations may occur as a result of the sensor being attached to one location for too long. To avoid skin irritations and lacerations, periodically inspect the sensor application site and change the application site regularly.

Venous Pulsation: Do not apply sensor too tightly as this results in venous pulsation which may severely obstruct circulation and lead to inaccurate measurements.

Ambient Temperature: At elevated ambient temperatures be careful with measurement sites that are not well perfused, because this can cause severe burns after prolonged application. All listed sensors operate without risk of exceeding 41°C on the skin if the initial skin temperature does not exceed 35°C.

Extremities to Avoid: Avoid placing the sensor on extremities with an arterial catheter, an NBP cuff or an intravascular venous infusion line.

Connecting SpO₂ Cables

Connect the sensor cable to the color-coded socket on the monitor. If you are using a disposable sensor, plug the sensor into the adapter cable and connect this to the monitor. Connect reusable sensors directly to the monitor.

CAUTION

Extension cables: Do not use more than one extension cable (M1941A). Do not use an extension cable with Philips reusable sensors or adapter cables with part numbers ending in -L (indicates "long" cable version).

Electrical Interference: Position the sensor cable and connector away from power cables, to avoid electrical interference.

Measuring SpO₂

During measurement, ensure that the application site:

- has a pulsatile flow, ideally with a signal quality indicator of at least medium.
- has not changed in its thickness (for example, due to edema), causing an improper fit of the sensor.

WARNING

- For fully conscious maternal patients, who have a normal function of perfusion and sensory perception at the measurement site:
 - To ensure skin quality and correct optical alignment of the sensor, inspect the application site when the measurement results are suspicious or when the patient complains about pressure at the application site, but at least every 24 hours. Correct the sensor alignment if necessary. Move the sensor to another site, if the skin quality changes.
- For all other patients:
 - Inspect the application site every two to three hours to ensure skin quality and correct optical alignment. Correct the sensor alignment if necessary. If the skin quality changes, move the sensor to another site. Change the application site at least every four hours.
- Injected dyes such as methylene blue, or intravascular dyshemoglobins such as methemoglobin and carboxyhemoglobin may lead to inaccurate measurements.
- Inaccurate measurements may result when the application site for the sensor is deeply pigmented or deeply colored, for example, with nail polish, artificial nails, dye or pigmented cream.
- Interference can be caused by:
 - High levels of ambient light (including IR warmers) or strobe lights or flashing lights (such as fire alarm lamps). (Hint: cover application site with opaque material).
 - Another SpO₂ sensor in close proximity (e.g. when more than one SpO₂ measurement is performed on the same patient). Always cover both sensors with opaque material to reduce cross-interference.
 - Electromagnetic interference, especially at perfusion indicator values below 1.0 or signal quality indicator below medium.
 - Excessive patient movement and vibration.

SpO₂ Signal Quality Indicator (Fast SpO₂ only)

The SpO₂ numeric is displayed together with a signal quality indicator (if configured and enough space is available) which gives an indication of the reliability of the displayed values.

The level to which the triangle is filled shows the quality of the signal; the indicator below shows a medium signal quality, the signal quality is at a maximum when the triangle is completely filled.



Assessing a Suspicious SpO₂ Reading

Traditionally, pulse rate from SpO₂ was compared with heart rate from ECG to confirm the validity of the SpO₂ reading. With newer algorithms, such as FAST-SpO₂, this is no longer a valid criteria because the correct calculation of SpO₂ is not directly linked to the correct detection of each pulse.

When pulse rate is very low, or strong arrhythmia is present, the SpO₂ pulse rate may differ from the heart rate calculated from ECG but this does not indicate an inaccurate SpO₂ value.

WARNING

With pulse oximetry, sensor movement, ambient light (especially strobe lights or flashing lights) or electromagnetic interference can give unexpected intermittent readings when the sensor is not attached. Especially bandage-type sensor designs are sensitive to minimal sensor movement that might occur when the sensor is dangling.

Understanding SpO₂ Alarms

This refers to SpO₂ specific alarms. See the “Alarms” on page 109 chapter for general alarm information. SpO₂ offers high and low limit alarms, and a high priority desat alarm. You cannot set the low alarm limit below the desat alarm limit.

CAUTION

If you measure SpO₂ on a limb that has an inflated noninvasive blood pressure cuff, a non-pulsatile SpO₂ INOP can occur. If the fetal monitor is configured to suppress this alarm there may be a delay of up to 60 seconds in indicating a critical status, such as sudden pulse loss or hypoxia.

Alarm Delays

There is a delay between a physiological event at the measurement site and the corresponding alarm at the monitor. This delay has two components:

- The general system delay time is the time between the occurrence of the physiological event and when this event is represented by the displayed numerical values. This delay depends on the algorithmic processing and the averaging time.
- The time between the displayed numerical values crossing an alarm limit and the alarm indication on the monitor. This delay is the combination of the configured alarm delay time plus the general system alarm signal delay time.

Adjusting the SpO₂ Alarm Limits

In the **Setup SpO₂** menu:

- Select **High Limit** then choose the upper alarm limit.
- Select **Low Limit** then choose the lower alarm limit.

Adjusting the Desat Limit Alarm

The Desat alarm is a high priority (red) alarm notifying you of potentially life threatening drops in oxygen saturation.

- 1 In the **Setup SpO₂** menu, select **Desat Limit**.
- 2 Adjust the limit.

Adjusting the Pulse Alarm Limits

See “Adjusting the Heart Rate / Pulse Alarm Limits” on page 230.

Setting Up Tone Modulation

If tone modulation is on, the QRS tone pitch lowers when the SpO₂ level drops. Remember, the QRS tone is derived from either heart rate (from MECG) or pulse (from built-in SpO₂) depending on which is currently displayed (see “Priority for Maternal Heart / Pulse Rate” on page 223).

NOTE

Pulse from CL SpO₂ and Toco MP does not provide a QRS tone.

Setting the QRS Volume

In the **Setup SpO₂** menu, select **QRS Volume** and set the appropriate QRS tone volume.

Monitoring Maternal Heart / Pulse Rate

You can monitor the maternal heart/pulse rate using one of four sources:

- Maternal pulse from Toco MP or CL Toco⁺ MP transducer (pulse rate)
- SpO₂ (pulse rate)
- Maternal heart rate (MHR) via MEECG electrodes
- NBP (pulse rate)

Maternal heart / pulse rates derived from Toco MP or CL Toco⁺ MP, SpO₂ and MEECG are **continuous** measurements, and are compared against the FHR for cross-channel verification. Average pulse rate derived from noninvasive blood pressure is an **intermittent** measurement, and is therefore not used for cross-channel verification.

Priority for Maternal Heart / Pulse Rate

Priority	Maternal Heart / Pulse Rate Source	Alarming	Used for CCV	Provides QRS Tone
1	HR from MEECG measurement	Yes	Yes	Yes
2	Pulse from SpO ₂ measurement	Yes	Yes	Yes
3	Pulse from CL SpO ₂ Pod measurement	Yes	Yes	No
4	Pulse from Toco MP measurement cableless or cabled	No	Yes	No
5	Pulse from CL NBP Pod measurement	No	No	No
6	Pulse from NBP measurement	No	No	No

Only one Maternal Heart Rate / Pulse Rate numeric will be displayed and recorded at a time (see priority table). If higher-priority measurements are connected but temporarily not providing valid numerics, lower-priority numerics may be displayed and recorded instead.


Misidentification of MHR for FHR

To reduce the possibility of mistaking the MHR for FHR, it is recommended that you monitor both maternal and fetal heart rates (see “Confirm Fetal Life Before Using the Monitor” on page 10 and “Cross-Channel Verification (CCV)” on page 149).

MHR from MECG Electrodes


You can measure MHR using the equipment combinations shown in the following illustrations.

Illustration 1 shows the complete connection chain from the foam electrodes applied to the patient to the fetal monitor using the patient module.




- 1 Patient Module (M2738A)
- 2 MECEG Adapter Cable (M1363A)
- 3 Pre-gelled Foam Electrodes (40493A/B/C/D/E)

Illustration 2 shows the equivalent chain using the Toco⁺ transducer.



- 1 Toco⁺ Transducer (M2735A)
- 2 MECEG Adapter Cable (M1363A)
- 3 Pre-gelled Foam Electrodes (40493A/B/C/D/E)

Illustration 3 shows the equivalent chain using a CL Toco⁺ MP transducer or a CL ECG/IUP transducer.




- 1 Avalon CL base station (866074)
- 2 CL Toco⁺ MP or CL ECG/IUP transducer (866075 or 866077)
- 3 MECC Adapter Cable (M1363A)
- 4 Pre-gelled Foam Electrodes (40493A/B/C/D/E)

To simultaneously measure DECG and MECC, you can use a Toco+, a CL Toco⁺ MP, or a CL ECG/IUP transducer for MECC. For measuring DECG you need a CL ECG/IUP transducer or a patient module (see also “Monitoring FHR Using DECG” on page 195).

Applying Electrodes

To derive the MHR (when you do not want to view the MECCG waveform), you can place the electrodes just below the outer end of the clavicle near each shoulder.



1 MECCG Electrodes

Making Connections

WARNING

Follow the instructions supplied with each of the monitoring accessories you are using.

Prepare for monitoring MHR using the list below. The standard procedures in use in your facility determine the sequence of actions.

- 1 Depending on the equipment you are using, ensure that *either* the Patient Module *or* the Toco⁺ transducer is connected to the fetal monitor. If you are using cableless monitoring use the CL ECG/IUP or the CL Toco⁺ MP transducer.
- 2 Connect a pre-gelled foam electrode to each of the two leads on the MECCG adapter cable.
- 3 Apply the foam electrodes to the patient, following the instructions supplied with the foam electrodes.
- 4 Depending on the equipment you are using, connect the pink connector plug on the MECCG adapter cable to the pink connector on *either* the Patient Module *or* the Toco⁺ transducer. If you are using cableless monitoring connect the MECCG Adapter cable to the connectors of the CL ECG/IUP or the CL Toco⁺ MP transducer.

You are now ready to monitor MHR.

Monitoring MHR

- 1 Switch on the recorder.
- 2 The maternal heart rate is labeled **HR** on the screen.

Monitoring MEEG Wave

WARNING

The fetal/maternal monitor is not a diagnostic ECG device. In particular, the display of fetal/maternal ECG is intended only for evaluating signal quality for fetal/maternal heart rate as derived from the ECG waveform.

When in doubt, it can be used to identify sources of compromised signal quality, such as noise or muscle artifacts. It can subsequently be used to verify the result of measures taken to resolve them (e.g., checking ECG cable connections or adapting the fetal ArtifactSuppress configuration).

The safety and effectiveness of the displayed fetal/maternal ECG waveform (i.e., P, QRS and T segments) for evaluation of fetal/maternal cardiac status during labor have not been evaluated.

FM30/50 When measuring MEEG with the Avalon FM30 or FM50, the MEEG waveform, along with the heart rate numeric, is displayed on the screen when using a Toco⁺, CL Toco⁺ MP transducer, a patient module, or a CL ECG/IUP transducer.

WARNING

The fetal/maternal monitor is NOT intended for use during defibrillation, electro-surgery, or MRI. Remove all transducers, sensors, and accessories before performing electro-surgery, defibrillation, or MRI, otherwise harm to the patient or the user can result.


Applying Electrodes

To obtain a satisfactory maternal ECG waveform you must use the RA to LL (lead II) position of the standard 5-lead ECG.



- 1 Place the RA electrode (A) directly below the clavicle and near the right shoulder.
- 2 Place the LL electrode (B) on the left lower abdomen.

Viewing the Waveform on the Screen




- 1 Measurement label (FHR1)
- 2 Measurement label (Toco)
- 3 Measurement label (DFHR3)
- 4 Measurement label (FHR1)
- 5 Measurement label (SpO₂)
- 6 Measurement label (HR)
- 7 1mV scale bar
- 8 MECCG wave with maternal label
- 9 DECCG wave with fetal label
- 10 Measurement label (Temp) maternal temperature
- 11 Measurement label NBP

NOTE

The 1mV scale bar for the DECCG and MECCG wave is not displayed on the screen if you monitor DECCG or MECCG with an Avalon CTS system. The Avalon CTS system does not provide a scaled ECG.

For the FM30/50, the **MECCG** wave is displayed automatically on the screen, labeled **MECCG**. If **DECCG** is also being monitored (FM30/50), and the DECCG Wave is configured to **On**, both waves are displayed, with the DECCG Wave above the MECCG wave. The DECCG Wave is labeled **DECCG**.

Troubleshooting

Problem	Possible Causes	Solutions
MECG Leads Off displayed. Numeric is displayed with a -?- for 10 seconds; INOP tone See also “Patient Alarms and INOPs” on page 121	One or more MECG leads is not attached.	Make sure that all required leads are attached If the wave is configured to be displayed on the monitor, you can observe if the ECG signal is clear, or if it shows interruptions and noise
	Bad electrical contact	Check positioning of the electrode, ensuring that none are displaced Check electrodes and replace if necessary
	Electrodes defective	
 prints repeatedly	The ultrasound transducer is measuring maternal pulse.	Reposition the ultrasound transducer. See “Recommended Actions for INOP Coincidence” on page 154
MECG Equip Malf displayed	Equipment malfunctions	See “Patient Alarms and INOPs” on page 121.
MECG Unplugged	Equipment not connected	

Printing the Waveform

You can print the MECG wave onto the trace paper. Refer to “Printing the ECG Waveform” on page 233.

Pulse Rate from Toco MP

The maternal pulse is taken from the Toco MP or CL Toco⁺ MP transducer when SpO₂ or MECG measurements are not used or have signal loss. When the pulse rate is very low, or strong arrhythmia is present, the pulse rate measured by the Toco MP or CL Toco⁺ MP transducer may differ from the heart rate calculated from MECG. If the mother is moving about, or began pressing during the second stage of labor, this can cause longer gaps in the recording of the maternal pulse signal. In this case use the SpO₂ or MECG measurement to derive the maternal heart rate.

WARNING

- No alarm is possible when Toco MP or CL Toco⁺ MP transducer is the source of the pulse rate.
- No QRS tone is audible when Toco MP or CL Toco⁺ MP transducer is the source of the pulse rate.
- The Toco MP or CL Toco⁺ MP transducer signal is significantly less reliable if the patient is up and moving about, or is pushing during the second stage of labor.

NOTE

In rare cases it is possible to pick up a fetal signal source. When a Toco MP or CL Toco⁺ MP transducer is connected to the monitor but not applied to the patient, the measurement may generate unexpected intermittent pulse readings.

Pulse Rate from SpO₂

If you are not monitoring MHR via MECCG electrodes, but you are monitoring SpO₂, the maternal pulse rate is derived from the SpO₂ measurement. The pulse numeric is labeled **Pulse** on the screen.

WARNING

No QRS tone is audible when the CL SpO₂ Pod is the source of the pulse rate.

Adjusting the Heart Rate / Pulse Alarm Limits

To adjust the pulse alarm limits for SpO₂:

- 1 In the **Setup SpO₂** menu, select **Pulse (SpO₂)**. This opens the **Setup Pulse (SpO₂)** menu.
- 2 Ensure **Pulse (SpO₂)** is **On**. Select **Pulse (SpO₂)** to toggle between **On** and **Off**.
- 3 Set the pulse alarm limit:
 - Select **High Limit** then choose the upper alarm limit for tachycardia from the pop-up list.
 - Select **Low Limit** then choose the lower alarm limit for bradycardia from the pop-up list.

To adjust the pulse alarm limits for MECCG:

- 1 In the **Setup MECCG** menu, select **MECCG/Pulse Alarms**. This opens the **Setup Pulse (MECCG)** menu.
- 2 Ensure **Pulse MECCG** is **On**. Select **Pulse MECCG** to toggle between **On** and **Off**.
- 3 Set the pulse alarm limit.
 - Select **High Limit** then choose the upper alarm limit for tachycardia from the pop-up list.
 - Select **Low Limit** then choose the lower alarm limit for bradycardia from the pop-up list.

Average Pulse Rate from Noninvasive Blood Pressure

WARNING

No alarm is possible when noninvasive blood pressure is the source of the pulse rate.

When you are measuring noninvasive blood pressure, the monitor can also calculate the average pulse rate. This occurs in either manual or automatic mode, when neither MEEG, SpO₂ nor pulse from Toco MP or CL Toco⁺ MP transducer are measured. The value is displayed on the screen, and printed on the trace. It is not the actual pulse value, but an average pulse rate, taken during the most recent noninvasive blood pressure measurement. The value is updated after each successive measurement. If you need a continuous measurement, you should monitor using MEEG, SpO₂ or pulse from Toco MP or CL Toco⁺ MP transducer.

Testing MEEG Mode

See the monitor's Service Guide.

Printing the ECG Waveform

FM30/50 You can print the ECG wave onto the trace paper. If you are monitoring both DECG and MEEG, both waves will be printed. The start of the wave recording is annotated above the wave with **MEEG** for Maternal ECG, with **DECG** for Direct fetal ECG, and with **25 mm/sec** below the wave.

WARNING

The fetal/maternal monitor is not a diagnostic ECG device. In particular, the display of fetal/maternal ECG is intended only for evaluating signal quality for fetal/maternal heart rate as derived from the ECG waveform.

When in doubt, it can be used to identify sources of compromised signal quality, such as noise or muscle artifacts. It can subsequently be used to verify the result of measures taken to resolve them (e.g., checking ECG cable connections or adapting the fetal **ArtifactSuppress** configuration).

The safety and effectiveness of the displayed fetal/maternal ECG waveform (i.e., P, QRS and T segments) for evaluation of fetal/maternal cardiac status during labor have not been evaluated.

The ECG waveform is printed along the bottom of the heart rate grid, and the three different possibilities look like this:

DECG waveform on its own



- 1 DECG
- 2 Paper speed

MECG waveform on its own



- 1 **MECG**
- 2 Paper speed

DECG and MECG waveforms




- 1 **DECG**
- 2 Paper speed

When the recorder is on, there are two choices for printing the ECG wave:

- **Separate:** This recording mode gives you a six-second ECG strip on the fetal trace paper in fast printout mode. The real-time fetal trace recording is temporarily interrupted while the ECG strip

prints. A new trace header is printed out to mark where the MEGC print-out starts and another FHR header mark signals when the fetal trace resumes.


The following trace shows the MEGC waveform:



- 1 FHR1 trace interrupted
 - 2 **MEGC** header
 - 3 Paper speed
 - 4 FHR1 trace resume
- **Overlap:** This recording mode gives you a delayed six-second snapshot of the maternal and/or direct fetal ECG for documentation on the fetal strip, but without interrupting the fetal trace. It takes 5 minutes to print this six-second snapshot at a recorder speed of 3 cm/min. It is documented as if it was recorded at 25 mm/s.

22 Printing the ECG Waveform

The following trace shows both the DECG and MEGC waveforms:



- 1 FHR1 and FHR2
- 2 **DECG** header
- 3 **MECG** header
- 4 Paper speed

To make your choice:

- 1 Enter the **Main Setup** menu.



- 2 Select **Fetal Recorder** to enter the **Fetal Recorder** menu.
- 3 Select **ECG Wave** to toggle between **Separate** and **Overlap**.

To print the ECG wave(s):

Either

Select the **Record ECG Wave** SmartKey.



Or

- 1 Enter the **Main Setup** menu by selecting the SmartKey.



- 2 Select **Fetal Recorder** to enter the **Fetal Recorder** menu.
- 3 Select **Record ECG Wave**.

Or

- 1 Select the **ECG Wave**.
- 2 Select **Record ECG Wave** in the ECG wave menu.

Paper Save Mode for Maternal Measurements

Your monitor's recorder features a Paper Save Mode, where maternal vital signs are recorded using less paper than during a normal trace recording.

When Paper Save Mode is enabled, and if the recorder is stopped, it will start automatically to print data from maternal measurements as they occur, and then stop again to save paper. You enable Paper Save Mode in Configuration Mode (default is off).

- A header is printed first before the measurements are recorded. A new header is also printed when there is a date change at midnight.
- Each NBP measurement is recorded. The time when the measurement **ended** is recorded.
- Each Temperature measurement is recorded. The time when the measurement **ended** is recorded.
- Other maternal parameters (SpO₂, maternal heart rate or Pulse) are recorded every five minutes. The rules described in the section “Priority for Maternal Heart / Pulse Rate” on page 223 apply.
- Paper Save Mode recording stops if there are no valid maternal measurements for more than one hour, and a message will notify you that there are no active parameters. Paper Save Mode recording will restart automatically when another valid measurement is made.

Event	Paper Save Mode Reactivation
One of the maternal measurements (see above) is valid again.	yes
The recorder is turned on and off again or a report has been recorded (e.g. NST Report).	yes
The Paper Advance function is used.	yes
The Paper Save Mode setting is set off and on again.	yes
ADT information has changed (e.g. because patient information has been completed or updated).	no
The monitor is restarted (e.g. by switching it off and on again).	yes
The date has changed (e.g. at midnight).	no

Recovering Data

The monitor stores trace data, including annotations, in its internal backup memory. This allows the monitor to recover trace data that would otherwise be lost under certain circumstances. This trace recovery data can be automatically retrieved and printed in the event of the paper running out, or automatically transmitted to an OB TraceVue/IntelliSpace Perinatal system (LAN connection only), allowing continuity of data.

The fetal trace printed from the trace recovery data contains all data from the real-time trace, with the exception of the maternal heart rate, the pulse numeric and the MEKG wave.

Note that the data in the memory is cleared when a software upgrade is performed.

CAUTION

Only use Philips paper. Using paper other than Philips paper may result in the failure to recover traces.

Recovering Traces on Paper

The monitor is able to recover traces by printing them out at a high speed from the monitor's backup memory. If the monitor runs out of paper, or if the paper drawer is open, the exact time when this happens is annotated in the backup memory. If the **Bridge Paperout** setting is set to **On** (default), when new paper is loaded and the recorder is started, a trace recovery printout of the data recovered from the backup memory is automatically printed out at high speed (up to 20 mm/s), starting from the time noted in the backup memory. This ensures that no data is lost. A minimum of one hour of trace recovery data can be printed out from the backup memory. When the trace recovery printout has finished, the recorder automatically switches back to continue recording the current trace at the normal speed.

Note the following:

- If you press the fetal recorder **Start/ Stop** SmartKey during a trace recovery printout, the recording stops and the next recording following a recorder restart will be a normal, real-time trace.



After switching the monitor off and then back on again, or following a power failure, the time of the last **Check Paper** INOP or paper-out detection is lost, and therefore any trace recovery data in the backup memory is no longer available to print. The next recording made following a restart of the recorder is a normal, real-time trace.

- The change back to a real-time recording from a trace recovery printout prompts the recording to restart. A new vertical trace header annotation consisting of the time, date and paper speed is printed, letting you see where the trace recovery printout ends, and where the real-time trace continues.
- There can be a gap of up to 30 seconds between the trace recovery printout and the beginning of the real-time trace.

Recovering Traces on an OB TraceVue/IntelliSpace Perinatal System

The trace recovery data stored in the monitor's backup memory can also be uploaded at high speed to an OB TraceVue/IntelliSpace Perinatal system connected over the LAN interface (OB TraceVue Revision E.00.00 or later, and IntelliSpace Perinatal H.0 or later).

When the OB TraceVue/IntelliSpace Perinatal system reconnects to the fetal monitor and detects that there is trace recovery data in the monitor's backup memory that has not yet been transmitted to the system, this data is transferred at high speed to the system. No user action is required.

The exact length of the recovered trace will vary depending on the amount of trace information, but will cover **at least one hour** of trace data, regardless of how many parameters are being measured.

To recover traces on an OB TraceVue/IntelliSpace Perinatal system, the following applies:

- The trace data in the monitor's internal memory must relate to a specific patient in the OB TraceVue/IntelliSpace Perinatal system. In other words, there were no discharge events made on the monitor that would change the patient context.
- The patient must have an open episode. No data will be uploaded if the patient is not admitted to OB TraceVue/IntelliSpace Perinatal. For this reason, it is not possible to use the monitor to collect patient data offline for later transmission to OB TraceVue/IntelliSpace Perinatal.
- Current online trace data is held back until the fast upload is complete.

Recording Stored Data

If the recorder is not running, you can choose to print trace data from the monitor's memory at any time. You can see a list of all stored traces, showing patient identification and trace period, in the **Stored Data Recording** window, from which you can choose one entry at a time.

CAUTION

Ensure that you admit each patient by name, including other patient identification information, and discharge the patient when you have finished monitoring, so that you can identify which trace period (entry in the patient list) refers to which patient.

Trace storage can be triggered by:

- Discharging a patient
- Powering on the monitor
- Entering Standby
- Entering Service Mode

Traces are not available for periods the monitor was switched off, in Service mode, in Standby, or if the trace period was shorter than one minute. To delete all stored trace periods press the **Erase All** key. The speed of the printout depends on the configured recorder speed and on the amount of trace data available. The fetal trace printed from the trace data contains all data from the real-time trace, with the exception of the maternal heart rate, the pulse numeric and the ECG wave.

Information for scale type, trace separation and recorder speed are not stored in the trace memory, but is applied when the stored recording starts. While the stored recording is printing, all functions are disabled, except that for stopping the recorder.

To start a stored data recording:

Either

- 1 Select the **Stored Data Rec** SmartKey.



Or

- 1 Enter the **Main Setup** menu using the SmartKey.



- 2 Select **Fetal Recorder** to open the **Fetal Recorder** menu.
- 3 Select **Stored Data Rec** to open the **Stored Data Recording** window.
- 4 Select an entry for a patient.
- 5 Select **All** to print all stored trace data for the selected entry, or select one of the choices on the other pop-up keys to print only a specified portion of the entry (for example, **Last 15 min** for the last 15 minutes of trace data).

To delete all stored trace periods:

Either

- 1 Select the **Stored Data Rec** SmartKey.



Or

- 1 Enter the **Main Setup** menu using the SmartKey.



- 2 Select **Fetal Recorder** to open the **Fetal Recorder** menu.
- 3 Select **Stored Data Rec** to open the **Stored Data Recording** window.
- 4 Select the **Erase All** key to delete all stored trace periods listed.
- 5 Select the **Confirm** key.

The current patient's entry is at the top of the list. The oldest entry at the bottom of the list has no start time specified, as part of the data originally stored may have been over-written by the current patient's data. The first part of the data, including the information for the start time, is no longer accessible.

It may be that you only see one entry (the current patient's data) in the **Stored Data Recording** window if that patient was monitored for a period long enough to erase any earlier entries.

If you make a stored data recording for an old entry (that is, not for the current patient), the recorder performs a fast trace printout of the stored data, advances the paper to the next paper fold, then stops.

If you make a stored data recording for the current patient, the recorder performs a fast trace printout of the stored data, and then reverts automatically to recording the real-time trace.

Care and Cleaning

Use only the Philips-approved substances and methods listed in this chapter to clean or disinfect your equipment. Warranty does not cover damage caused by using unapproved substances or methods.

Philips makes no claims regarding the efficacy of the listed chemicals or methods as a means for controlling infection. Consult your hospital's Infection Control Officer or Epidemiologist. For comprehensive details on cleaning agents and their efficacy refer to "Guideline for Disinfection and Sterilization in Healthcare Facilities" issued by the U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control, Atlanta, Georgia, 2008. See also any local policies that apply within your hospital, and country.

General Points

The transducers and patient modules are sensitive instruments. Handle them with care.

Keep your monitor, transducers, patient modules, cables and accessories free of dust and dirt. After cleaning and disinfection, check the equipment carefully. Do not use if you see signs of deterioration or damage. If you need to return any equipment to Philips, **always** decontaminate it first before sending it back in appropriate packaging.

Observe the following general precautions:

- Always follow carefully and retain the instructions that accompany the specific cleaning and disinfecting substances you are using.
- Always dilute according to the manufacturer's instructions or use lowest possible concentration.
- Do not allow liquid to enter the case.
- Do not immerse the monitor in liquid. Protect it against water sprays or splashes.
- Do not pour liquid onto the system.
- Never use abrasive material (such as steel wool or silver polish).
- Never use bleach.

WARNING

- Do not operate the monitor if it is wet. If you spill liquid on the monitor, contact your service personnel or Philips service engineer.
- Do not perform underwater monitoring (for example, in a bath or shower) using wired transducers.
- Place the monitor where there is no chance of contact with, or falling into water or other liquid.
- Do not dry equipment using heating devices such as heaters, ovens (including microwave ovens), hair dryers and heating lamps.
- Do not put equipment or accessories in autoclave (for sterilization).

Cleaning and Disinfecting

Clean and disinfect the Avalon FM20, FM30, FM40 and FM50 fetal monitors and the transducers M2734A, M2734B, M2735A, M2736A, and M2738A (including ECG adapter cables) and the Avalon CL base station and the cableless transducers after each use. Clean equipment before disinfecting. For other accessories, see “Cleaning and Disinfecting Monitoring Accessories” on page 247.

Clean with a lint-free cloth, moistened with warm water (40°C/104°F maximum) and soap, a diluted non-caustic detergent, tenside, or phosphate based cleaning agent (see “Recommended Disinfectants” on page 247). Do not use strong solvents such as acetone or trichloroethylene. After cleaning, disinfect using only the approved disinfecting agents listed.

CAUTION

Solutions: Do not mix disinfecting solutions (such as bleach and ammonia) as hazardous gasses may result.

Skin contact: To reduce the risk of skin irritations, do not allow a cleaning or disinfecting agent to leave residues on any of the equipment surfaces - wipe it off with a cloth dampened with water, after allowing the appropriate time for the agent to work, or before applying to a patient.

Hospital policy: Disinfect the product as determined by your hospital's policy, to avoid long term damage to the product.

Local requirements: Observe local laws governing the use of disinfecting agents.

Touch display: To clean and disinfect the touch-enabled display, disable the touch operation by switching off the monitor during the cleaning procedure, or by selecting and holding the **Main Screen** key until the padlock symbol appears on it, indicating that touch operation is disabled. Select and hold again to re-enable touch operation.

Take extra care when cleaning the screen of the monitor because it is more sensitive to rough cleaning methods than the housing. Do not permit any liquid to enter the monitor case and avoid pouring it on the monitor while cleaning. Do not allow water or cleaning solution to enter the measurement connectors. Wipe around and not over connector sockets, or those of the Toco⁺, CL Toco transducer, ECG and IUP Patient Modules, CL ECG/IUP transducer and adapter cables.

Wash soiled reusable belts with soap and water. Water temperature must not exceed 60°C/140°F.

Recommended Disinfectants

We recommend that you use one of the following disinfectants:

Product Name	Product Type	Ingredients
Isopropanol	liquid	Isopropanol 80%
Bacillol® AF	liquid, spray	100 g concentrate contains: Propan-1-ol 45.0 g; Propan-2-ol 25.0 g; Ethanol 4.7 g.
Bacillol®25	liquid	Ethanol 100 mg/g Propan-2-ol (= 2-Propanol) 90 mg/g; Propan-1-ol (= 1-Propanol) 60 mg/g
Meliseptol®	spray	50% 1-Propanol
Accel TB RTU	liquid	0.5% accelerated hydrogen peroxide
Oxivir® Tb Cleaner Disinfectant	spray	0.5% accelerated hydrogen peroxide
Oxivir® Tb Wipes	wipes	0.5% accelerated hydrogen peroxide
Carpe Diem ^{TM/MC} Tb Ready-to-Use General Virucide, Bactericide, Tuberculocide, Fungicide, Sanitizer	spray	0.5% accelerated hydrogen peroxide
Carpe Diem ^{TM/MC} Tb Wipes	wipes	0.5% accelerated hydrogen peroxide
Super Sani-Cloth Germicidal Disposable Wipes	wipes	isopropanol 55% quaternary ammonium chlorides 0.5%
SANI-CLOTH® PLUS Germicidal Disposable Wipes	wipes	isopropanol 15% quaternary ammonium chlorides 0.25%
SANI-CLOTH® HB Germicidal Germicidal Disposable Wipes	wipes	isopropanol < 0.15% quaternary ammonium chlorides 0.14%

Cleaning and Disinfecting Monitoring Accessories

To clean, disinfect and sterilize reusable transducers, sensors, cables, leads, and so forth, refer to the instructions delivered with the accessory.

Do not allow a cleaning or disinfecting agent to leave residues on any of the equipment surfaces. Wipe residues off, after allowing the appropriate time to for the agent to work, with a cloth.

Cleaning and Disinfecting the Tympanic Temperature Accessories

Probe and Thermometer Body

- 1 Wipe the thermometer body clean with a damp cloth. The water temperature should not exceed 55°C (130°F). Do not soak, rinse, or submerge the thermometer under water.
You may add a mild detergent to the water.
- 2 Clean the probe tip with a lint free swab. If the probe tip is soiled, clean it with a dampened swab.
- 3 After you have removed all foreign matter, clean the thermometer lens at the end of the probe tip with a lint free swab or lens wipe. The thermometer lens must be free from fingerprints and/or smudges for proper operation.
- 4 Thoroughly dry all surfaces before using the equipment.

CAUTION

Do not use cleaners and disinfectants such as Spray-Nine™, PhisoHex™, Hibiclens™, or Vesta-Syde™ as they may result in damage to the thermometer case.

Occasional use of a 10:1 water and hypochlorite mixture or a damp isopropyl alcohol wipe or Cidex™ or ManuKlenz™ or VIROX™ or CaviWipes™ cleansing agents is acceptable, however, prolonged or repeated use of these chemicals may result in damage to the thermometer case and display area.

Use of a cloth or sponge is recommended for cleaning. Never use an abrasive pad or an abrasive cleaner on the thermometer.

The thermometer is non-sterile. Do not use ethylene oxide gas, heat, autoclave, or any other harsh method to sterilize this thermometer.

Sterilizing

Sterilization is not recommended for this monitor, related products, accessories or supplies unless otherwise indicated in the Instructions for Use that accompany the accessories and supplies.

WARNING

Do not put device and accessories in autoclave (for sterilization).

Maintenance

WARNING

Schedule: Failure on the part of the responsible individual hospital or institution employing the use of this equipment to implement a satisfactory maintenance schedule may cause undue equipment failure and possible health hazards.

In case of problems: If you discover a problem with any of the equipment, contact your service personnel, Philips, or your authorized supplier.

Electric shock hazard: Do not open the monitor housing. Refer all servicing to qualified service personnel.

Inspecting the Equipment and Accessories

You should perform a visual inspection **before each use**, and in accordance with your hospital's policy. With the monitor switched off:

- 1 Examine unit exteriors for cleanliness and general physical condition. Make sure that the housings are not cracked or broken, that everything is present, that there are no spilled liquids that may have entered the housing, and that there are no signs of abuse.
- 2 Inspect all accessories (transducers, sensors and cables, and so forth). Do not use a damaged accessory.
- 3 Avoid frequent drops of the transducers.
- 4 Switch the monitor on and make sure the display is bright enough. If the brightness is not adequate, contact your service personnel or your supplier.

Batteries Preventive Maintenance

For the FM20/30 with a battery option refer to “Using Batteries” on page 102.

Inspecting the Cables and Cords

- 1 Examine all system cables, the power plug and cord for damage. Make sure that the prongs of the plug do not move in the casing. If damaged, replace it with an appropriate power cord.
- 2 Inspect the cables, leads and their strain reliefs for general condition. Make sure there are no breaks in the insulation. Make sure that the connectors are properly engaged at each end to prevent rotation or other strain.
- 3 Carry out performance assurance checks as described in the monitor's Service Guide.

Maintenance Task and Test Schedule

The following tasks are for Philips-qualified service professionals. All maintenance tasks and performance tests are documented in detail in the service documentation supplied on the monitor's documentation DVD.

Ensure that these tasks are carried out as indicated by the monitor's maintenance schedule, or as specified by local laws, whichever comes sooner. Contact a Philips-qualified service professionals if your monitor needs a safety or performance test. Clean and disinfect equipment to decontaminate it before testing or maintaining it.

Maintenance and Test Schedule	Frequency
Visual Inspection	Before each use.
Clean and disinfect the equipment	After each use.
Safety checks according to IEC 60601-1, and where applicable, to national standards	At least once every two years, or as specified by local laws. After any repairs where the power supply has been replaced (by an authorized service professional). If the monitor has been dropped, it must be repaired/checked by an authorized service agent.
Performance assurance for all measurements	At least once every two years, or if you suspect the measurement values are incorrect.
Noninvasive blood pressure calibration	At least once every two years, or as specified by local laws.
Clean the thermal printhead	At each paper pack change, or every 500 m of paper run.

Recorder Maintenance

Removing the Paper Guide: FM40/FM50

FM40/50 The paper guide is removable, and you can use the recorder without it. When **not** using the paper guide, ALWAYS tear off the paper along the perforation to avoid possible paper misalignment (see “Tearing Off the Paper” on page 80).

To remove the paper guide:



- 1 Press the paper eject button to open the paper drawer.



- 2 Hinge the transparent paper guide forward.





A Protrusion holds paper guide in closed position.



3 Release the paper guide from one side of the holder.



4 Then remove the paper guide.



Refitting is a reversal of the removal procedure.

Storing Recorder Paper

Recorder paper is not intended for long-term archival storage. Another medium should be considered if this is required.

Dyes contained in thermal papers tend to react with solvents and other chemical compounds that are being used in adhesives. If these compounds come into contact with the thermal print, the print may be destroyed over time. You can take the following precautionary measures to help avoid this effect:

- Store the paper in a cool, dry and dark place.
- Do not store the paper at temperatures over 40°C (104°F).
- Do not store the paper where the relative humidity exceeds 60%.
- Avoid intensive light (UV light), as this may cause the paper to turn gray or the thermal print to fade.
- Avoid storing the thermal paper in combination with the following conditions:
 - Papers that contain organic solvents. This includes papers with tributyl and/or dibutyl phosphates, for example recycled paper.
 - Carbon paper and carbonless copy paper.
 - Products containing polyvinyl chlorides or other vinyl chlorides for example (but not exclusively) document holders, envelopes, letter files, divider sheets.
 - Detergents and solvents, such as alcohol, ketone, ester and others, including cleaning and disinfecting agents.
 - Products containing solvent-based adhesives such as (but not exclusively) laminating film, transparent film or labels sensitive to pressure.

To ensure long lasting legibility and durability of thermal printouts, store your documents separately in an air-conditioned place and use:

- only plasticizer-free envelopes or divider sheets for protection.
- laminating films and systems with water-based adhesives.

Using such protective envelopes cannot prevent the fading effect caused by other, external agents.

Cleaning the Print Head

To clean the recorder's thermal print-head:



FM20/30



FM40/50

- 1 Switch off the monitor.
- 2 Open the paper drawer, and remove the paper if necessary, to gain access to the thermal print head.
- 3 Gently clean the thermal print head with a cotton swab or soft cloth soaked in isopropyl alcohol.

Disposing of the Monitor

WARNING

To avoid contaminating or infecting personnel, the environment or other equipment, make sure you disinfect and decontaminate the monitor appropriately before disposing of it in accordance with your country's laws for equipment containing electrical and electronic parts. For disposal of parts and accessories such as thermometers, where not otherwise specified, follow local regulations regarding disposal of hospital waste.

You can disassemble the monitor and the transducers as described in the Installation and Service Manual. You will find detailed disposal information on the following web page: <http://www.healthcare.philips.com/main/about/Sustainability/Recycling/pm.wpd>



Do not dispose of waste electrical and electronic equipment as unsorted municipal waste. Collect it separately, so that it can be safely and properly reused, treated, recycled, or recovered.

The Recycling Passports located there contain information on the material content of the equipment, including potentially dangerous materials which must be removed before recycling (for example, batteries and parts containing mercury or magnesium).

Accessories and Supplies

All accessories listed for the fetal monitor may not be available in all geographies. You can order parts, accessories and supplies from Philips supplies at www.healthcare.philips.com or consult your local Philips representative for details. All accessories and supplies listed here are reusable, unless indicated otherwise.

WARNING

Reuse: Disposable accessories and supplies intended for single use, or single patient use only, and are indicated as such on their packaging. Never reuse disposable accessories and supplies, such as transducers, sensors, electrodes and so forth that are intended for single use, or single patient use only.

Approved accessories: Use only Philips-approved accessories.

Packaging: Do not use a sterilized accessory if its packaging is damaged.

Protection against electric shocks: The transducers and accessories listed in this chapter are not defibrillator proof.

Electro-Surgery, Defibrillation and MRI: The fetal/maternal monitors are not intended for use during defibrillation, electro-surgery, or MRI. Remove all transducers, sensors, and accessories before performing electro-surgery, defibrillation, or MRI, otherwise harm can result.

Information on Latex

All Philips transducers and accessories are latex-free, unless indicated otherwise in the following tables.

Avalon CL Base Station

CL Base Station	Part Number
Avalon CL Base Station can either be ordered with the option K30 (red connector) or K40 (black connector)	866074

Transducers

Transducer	Part Number
Avalon Toco Transducer	M2734A
Avalon Toco ⁺ Transducer for Toco, DECG, MECG or IUP monitoring	M2735A
Avalon Toco MP Transducer for Toco and Maternal Pulse	M2734B
Avalon Ultrasound Transducer	M2736A
Avalon Ultrasound Transducer USA	M2736AA
ECG/IUP Patient Module (for DECG, MECG or IUP)	M2738A
Avalon CL Toco ⁺ MP Transducer for use with the Avalon CL base station	866075
Avalon CL Ultrasound Transducer for use with the Avalon CL base station	866076
Avalon CL ECG/IUP Transducer for use with the Avalon CL base station	866077
Connector Caps for Avalon CL Cableless Smart Transducers (kit of 10) for use with the Avalon CL Toco ⁺ MP transducer	989803184841
Belt Clip for Avalon CL Cableless Smart Transducers (kit of 10)	989803184851
CL SpO ₂ Pod for use with the Avalon CL base station	865215
CL NBP Pod for use with the Avalon CL base station	865216
Remote Event Marker	989803143411

Fetal Accessories

Accessory	Description	Part Number
Belt (reusable, gray, water resistant)	32 mm wide, 15 m roll	M4601A1
	60 mm wide, 5 belts	M4602A
	60 mm wide, 15 m roll	M4603A
	50 mm wide, 5 belts	M1562B


Accessory	Description	Part Number
Belt (reusable, brown, contains latex)	32 mm wide, 15 m roll	1500-06281
	50 mm wide, 5 belts	M1562A
	60 mm wide, 5 belts	1500-0642
	60 mm wide, 15 m roll	1500-0643
Belt (disposable, yellow, water resistant)	60 mm wide, pack of 100	M2208A
Ultrasound gel	12 Bottles	40483A
	5 liter refill (with dispenser) for 40483A Shelf life: 24 months max.	40483B
Belt buttons (kit of 10)		M1569A
Butterfly belt clip (kit of 6)		989803143401
Avalon CL Connector Caps	Connector Caps for Avalon CL Cableless Smart Transducers (kit of 10)	989803184841
Avalon CL Belt Clip	Belt Clip for Avalon CL Cableless Smart Transducers (kit of 10)	989803184851
Avalon CL Battery Replacement Kit		989803184861
Cable Management Kit		989803148841
Belt Clips for Avalon Smart Transducers (kit of 6)		989803143401
DECG Accessories: New Philips DECG Solution <i>(NOT compatible with QwikConnect Plus Solution accessories)</i>	DECG reusable legplate adapter cable (with flushing port)	989803137651
	DECG leg attachment electrode for DECG legplate adapter cable	989803139771
	DECG fetal scalp electrode: single spiral, worldwide availability	989803137631
	DECG fetal scalp electrode: double spiral, Europe only. Not for USA	989803137641
Disposable Koala IUP catheter		M1333A
Reusable Koala IUP adapter cable		989803143931

DECG Accessories: Component Compatibility

Use the following pictorial guide to check component compatibility for DECG accessories.

CAUTION

Do not mix accessories from the New Philips DECG Solution (marked 1) with those from the QwikConnect Plus Solution (marked 2).



MECG Accessories

Accessory	Part Number
MECG reusable adapter cable	M1363A
Foam ECG electrodes, snap-fit, for MECG Adapter Cable (disposable)	40493D/E

Noninvasive Blood Pressure Accessories

The following accessories are approved for use with the fetal monitor:

Adult Multi-Patient Comfort Cuffs and Disposable Cuffs

Maternal Patient Category	Limb Circumference	Bladder Width	Disposable cuff Part No.	Reusable cuff Part No.	Tubing
Adult (Thigh)	42.0 - 54.0 cm	20.0 cm	M1879A	M1576A	M1598B (1.5 m) or M1599B (3.0 m)
Large Adult	34.0 - 43.0 cm	16.0 cm	M1878A	M1575A	
Adult	27.0 - 35.0 cm	13.0 cm	M1877A	M1574A	
Small Adult	20.5 - 28.0 cm	10.5 cm	M1876A	M1573A	
Cuff kit of 4 adult sizes				M1578A	

Adult EasyCare Reusable Cuffs

Maternal Patient Category (color)	Limb Circumference	Bladder Width	Part No.	Tubing
Adult Thigh (grey)	45.0 - 56.5 cm	21.0 cm	M4559B	M1598B (1.5 m) or M1599B (3.0 m)
Adult Thigh (grey) pack of 5 cuffs	45.0 - 56.5 cm	21.0 cm	M4559B5	
Large Adult X-Long (burgundy)	35.5 - 46.0 cm	17.0 cm	M4558B	
Large Adult X-Long (burgundy) pack of 5 cuffs	35.5 - 46.0 cm	17.0 cm	M4558B5	
Large Adult (burgundy)	35.5 - 46.0 cm	17.0 cm	M4557B	
Large Adult (burgundy) pack of 5 cuffs	35.5 - 46.0 cm	17.0 cm	M4557B5	
Adult X-Long (navy blue)	27.5 - 36.5 cm	13.5 cm	M4556B	
Adult X-Long (navy blue) pack of 5 cuffs	27.5 - 36.5 cm	13.5 cm	M4556B5	
Adult (navy blue)	27.5 - 36.5 cm	13.5 cm	M4555B	
Adult (navy blue) pack of 5 cuffs	27.5 - 36.5 cm	13.5 cm	M4555B5	
Small Adult (royal blue)	20.5 - 28.5 cm	10.6 cm	M4554B	
Small Adult (royal blue) pack of 5 cuffs	20.5 - 28.5 cm	10.6 cm	M4554B5	
Cuff kits containing one small adult, one adult, one large adult and one thigh cuff.			864288	
Cuff kits containing one small adult, one adult, one adult X-long, one large adult, one large adult X-long and one thigh cuff			864291	

Adult Single Patient Soft Cuffs

Maternal Patient Category	Limb Circumference	Bladder Width	Part No.	Tubing
Adult (Thigh)	45.0 - 56.5 cm	20.4 cm	M4579B	M1598B (1.5 m) or M1599B (3.0 m)
Large Adult X-Long	35.5 - 46.0 cm	16.4 cm	M4578B	
Large Adult	35.5 - 46.0 cm	16.4 cm	M4577B	
Adult X-Long	27.5 - 36.5 cm	13.1 cm	M4576B	
Adult	27.5 - 36.5 cm	13.1 cm	M4575B	
Small Adult	20.5 - 28.5 cm	10.4 cm	M4574B	

IntelliVue CL NBP Pod Accessories

Description	Limb Circumference Range	Bladder Width	Contents	Product No.
Mobile CL Reusable Small Adult Cuff	21 - 27 cm	10.5 cm	1 cuff	989803163171
Mobile CL Reusable Adult Cuff	26.0 - 34.5 cm	13.0 cm	1 cuff	989803163191
Mobile CL Reusable Large Adult Cuff	33.5 - 45.0 cm	16.0 cm	1 cuff	989803163211
Mobile CL Single-Patient Small Adult Cuff	21 - 27 cm	10.5 cm	20 cuffs	989803163181
Mobile CL Single-Patient Adult Cuff	26.0 - 34.5 cm	13.0 cm	20 cuffs	989803163201
Mobile CL Single-Patient Large Adult Cuff	33.5 - 45.0 cm	16.0 cm	20 cuffs	989803163221
Mobile CL NBP Cradle Kit	-	-	20 cradles	989803163251
Mobile CL Extension Air Hose, 1.0 m	-	-	1 extension air hose	989803163131
Mobile CL NBP Battery Kit	-	-	1 Battery 1 disassembly tool 1 front housing	989803163261
Telemetry Pouch with window	-	-	50 pouches	989803137831
Telemetry Pouch with window	-	-	4 boxes of 50 pouches	989803140371
White Telemetry Pouch with Snaps; box of 50	-	-	50 pouches	989803101971 (9300-0768-050)
White Telemetry Pouch with Snaps; 4 boxes of 50	-	-	4 boxes of 50 pouches	989803101981 (9300-0768-200)

SpO2 Accessories

Some Nellcor sensors contain natural rubber latex which may cause allergic reactions. See the Instructions for Use supplied with the sensors for more information. M1901B, M1903B and M1904B are not available in USA from Philips. Purchase Nellcor OxiCliq sensors and adapter cables directly from Tyco Healthcare. Some sensors may not be available in all countries.

Do not use more than one extension cable with any sensors or adapter cables. Do not use an extension cable with Philips reusable sensors or adapter cables with part numbers ending in -L (indicates "Long" version).

All listed sensors operate without risk of exceeding 41°C on the skin if ambient temperature is below 37°C.

Make sure that you use only the accessories that are specified for use with this device, otherwise patient injury can result.

Philips SpO2 Accessories

Description	Comments	Product Number
Philips reusable sensors		
Adult sensor (2.0 m cable), for patients over 50 kg. Any finger, except thumb.	No adapter cable required.	M1191B
M1191A with longer cable (3.0 m)		M1191BL
Small adult, pediatric sensor (1.5 m cable) for patients between 15 kg and 50 kg. Any finger except thumb. Use only on adult patients with FM30/40/50.		M1192A
Ear sensor (1.5 m cable) for patients more than 40 kg. Use only on adult patients with FM30/40/50.		M1194A
Adult clip sensor (3 m cable) for patients over 40 kg. Any finger except thumb.		M1196A
Adult clip sensor (2 m cable) for patients over 40 kg. Any finger except thumb.		M1196S
Adult sensor (0.45 m), for patients over 50 kg. Any finger except thumb.		Requires M1943A (1.0 m) or M1943AL (3.0 m) adapter cable.
Small adult, pediatric sensor (0.45 m cable) for patients between 15 kg and 50 kg. Any finger except thumb. Use only on adult patients with FM30/40/50.	M1192T	
Adult clip sensor (0.9 m cable) for patients over 40 kg. Any finger except thumb.	M1196T	

27 Accessories and Supplies

Description	Comments	Product Number
Special Edition (SE) Adult sensor (3 m cable), for patients over 50 kg. Any finger except thumb.	No adapter cable required.	M1191ANL
Special Edition (SE) Small adult, pediatric sensor (1.5 m cable) for patients between 15 kg and 50 kg. Any finger except thumb. Use only on adult patients with FM30/40/50.	SE sensors work with FM30/40/50, as well as with OxiMax-compatible SpO ₂ versions of other Philips monitors.	M1192AN
Special Edition (SE) Ear sensor (1.5 m cable) for patients more than 40 kg.		M1194AN
Philips disposable sensors. Not available in the USA.		
Identical to OxiMax MAX-A	Requires M1943A (1.0 m) or M1943AL (3.0 m) adapter cable	M1904B
Identical to OxiMax MAX-P		M1903B
Identical to OxiMax MAX-N		M1901B
Philips disposable sensors. Available worldwide.		
Adult/Pediatric finger sensor (0.45 m cable) Use only on adult patients with FM30/40/50.	Requires M1943A (1.0 m) or M1943AL (3.0 m) adapter cable	M1131A
Adult/Infant/Neonatal (0.9 m cable) for patients > 40 kg. Any finger except thumb. Use only on adult patients with FM30/40/50.		M1133A
Adult/Infant/Neonatal (0.9 m cable) for patients > 40 kg. Any finger except thumb. Adhesive-free Use only on adult patients with FM30/40/50.		M1134A

Nellcor SpO2 Accessories

Product Number	Description	Comments
NELLCOR disposable sensors (must be ordered from Nellcor)		
OxiMax MAX-A	Adult finger sensor (patient size > 30 kg)	Requires M1943A (1.0 m) or M1943AL (3.0 m) adapter cable.
OxiMax MAX-AL	OxiMax MAX-A with long cable	
OxiMax MAX-P	Pediatric foot/hand sensor (patient size 10-50 kg) Use only on adult patients with FM30/40/50.	
OxiMax MAX-N	Adult finger or neonatal foot/hand sensor (patient size > 40 kg or < 3 kg) Use only on adult patients with FM30/40/50.	
Oxisensor II D-25	Adult sensor (patient size > 30 kg)	Requires M1943A (1.0 m) or M1943AL (3.0 m) adapter cable.
Oxisensor II D-20	Pediatric sensor (patient size 10-50 kg) Use only on adult patients with FM30/40/50.	
Oxisensor II N-25	Neonatal/Adult sensor (patient size < 3 kg or > 40 kg) Use only on adult patients with FM30/40/50.	
OxiCliq A	See OxiMax MAX-A	Requires M1943A (1.0 m) or M1943AL (3.0 m) adapter cable together with OC3 adapter cable.
OxiCliq P	See OxiMax MAX-P Use only on adult patients with FM30/40/50.	
OxiCliq N	See OxiMax MAX-N Use only on adult patients with FM30/40/50	

Masimo SpO2 Accessories

Product Number	Description	Philips Part Number
MASIMO LNOP reusable sensors (No adapter cable required)		
LNOP DC-I	Adult Finger Sensor (> 30 kg)	989803140321
LNOP DC-IP	Pediatric Finger Sensor (10 - 50 kg) Use only on adult patients with FM30/40/50	989803140331
LNOP YI	Multi-Site Sensor (> 1kg) Use only on adult patients with FM30/40/50	n/a
LNOP TC-I	Ear Sensor (> 30 kg) Use only on adult patients with FM30/40/50	989803140341
MASIMO LNCS reusable sensors (No adapter cable required)		
LNCS DC-I	Adult Finger Sensor (> 30 kg)	989803148281
LNCS DC-IP	Pediatric Finger Sensor (10 - 50 kg) Use only on adult patients with FM30/40/50	989803148291
LNCS TC-I	Ear Sensor (> 30 kg) Use only on adult patients with FM30/40/50	989803148301
MASIMO LNOP disposable adhesive sensors		
LNOP Adt	Adult Sensor (> 30 kg)	989803140231

27 Accessories and Supplies

Product Number	Description	Philips Part Number
LNOP Adtx	Adult Sensor (> 30 kg)	n/a
LNOP Pdt	Pediatric Sensor (10 - 50 kg) Use only on adult patients with FM30/40/50	989803140261
LNOP Pdtx	Pediatric Sensor (10 - 50 kg) Use only on adult patients with FM30/40/50	n/a
LNOP Neo-L	Neonatal Sensor (< 3kg) or Adult adhesive Sensor (> 40 kg) Use only on adult patients with FM30/40/50	989803140291
MASIMO LNCS disposable adhesive sensors (No adapter cable required)		
LNCS Adtx	Adult Sensor (> 30 kg)	989803148231
LNCS Pdtx	Pediatric Finger Sensor (10 - 50 kg) Use only on adult patients with FM30/40/50	989803148241
LNCS Neo-L	Neonatal Foot Sensor (<3 kg) or Adult Finger Sensor (>40 kg) Use only on adult patients with FM30/40/50	989803148271

IntelliVue CL SpO₂ Pod Accessories

All listed sensors operate without risk of exceeding 41°C on the skin, if the initial skin temperature does not exceed 35°C.

Ensure that you use only the accessories that are specified for use with this device, otherwise patient injury can result.

Description	Contents	Product No.
Mobile CL 20 single patient SpO ₂ Sensors and Cradles for use on adult patients >10 kg	20 Single-Patient Mobile CL DSpO ₂ -1A Sensors 20 Single-Patient Wristbands 20 Single-Patient Cradles pre-assembled	989803165941
Mobile CL 20 single patient SpO ₂ Sensors for use on adult patients > 10 kg	20 Single-Patient Mobile CL DSpO ₂ -1A Sensors	989803165921
Mobile CL reusable SpO ₂ sensor and Cradles for use on adult patients > 15 kg	1 Reusable Mobile CL RSpO ₂ -1A Sensor 20 Single-Patient Cradles with pre-attached Wristbands	989803165931 ¹
Mobile CL 20 SpO ₂ Cradles (single patient)	20 Single-Patient Cradles with pre-attached Wristbands	989803165951
Mobile CL 50 SpO ₂ Wristbands (single patient)	50 Single-Patient Wristbands	989803165961
Mobile CL SpO ₂ Battery Kit	1 Battery 1 disassembly tool 1 front housing	989803168861

¹ May not be available in all geographies

Extension / Adapter Cables

Description	Comments	Product No.
Extension cable (2 m)	For use with Philips reusable sensors and adapter cables.	M1941A
Adapter cable (1.1 m cable)	Adapter cable for Philips/Nellcor disposable sensors.	M1943A
Adapter cable (3 m cable)		M1943AL
Adapter Cable for OxiCliq sensors	Available from Nellcor.	OC 3
LNOP MP12	LNOP MP Series Patient Cable (3.6 m) Adapter Cable for Masimo LNOP sensors	M1020-61102
LNC MP10	LNCS MP Series Patient Cable (3.0 m) Adapter Cable for Masimo LNCS sensors	989803148221

Tympanic Temperature Accessories

Description	Part No.
Temperature probe	989803180831
Disposable probe cover with CE marking (22 boxes @ 96)	989803179611
Disposable probe cover (22 boxes @ 96)	989803179381

Recorder Paper

Supplied in cases of 40 packs. Each pack has 150 numbered pages. Single use. Use the paper specified here.

Geography	FHR Scale	Grid Color	Scale Units	Highlighted 3 cm Lines?	Product No.
USA/Canada and Asia	30 - 240	Red/Orange	mmHg	Yes	M1910A
Europe/Japan	50 - 210	Green	mmHg and kPa	No	M1911A
Japan	50 - 210	Green	mmHg	Yes	M1913A
Japan	50 - 210	Green*	mmHg	Yes	M1913J

* Bradycardia and tachycardia alarm ranges are shaded.

Batteries

Description	Comment	Product No.
Smart Battery 10.8 V, 6000 mAh, Lithium Ion	For Avalon FM20 or Avalon FM30 with battery option #E25	M4605A
Avalon CL Cableless Smart transducer Battery Replacement Kit	Consists of one Philips Lithium Ion Battery (Part No. 453564107871), a tool to open and close the cableless transducer for battery replacement and two replacement O-ring seals.	989803184861

Specifications and Standards Compliance

The monitors are intended to monitor a mother and her fetus(es), which from an electrical safety point of view, are one person.

Environmental Specifications

The monitor may not meet the given performance specifications if stored and used outside the specified temperature and humidity ranges.

Avalon CL Base station 866074		
Temperature Range	Operating	0°C to 45°C (32°F to 113°F)
	Storage/Transportation	-20°C to 60°C (-4°F to 140°F)
Humidity Range	Operating	<95% relative humidity @ 45°C/113°F
	Storage/Transportation	<90% relative humidity @ 60°C/140°F
Altitude Range	Operating	-500 to 3000 m/-1640 to 9840 ft.
	Storage/Transportation	-500 to 13100 m/-1640 to 43000 ft.

Monitor (M2702A/M2703A); Interface Cable for Avalon CTS (M2731-60001 and M2732-60001)		
Temperature Range	Operating	Without battery option: 0°C to 45°C (32°F to 113°F)
		With battery option/charging: 0°C to 35°C (32°F to 95°F)
		With battery option/fully charged: 0°C to 40°C (32°F to 104°F)
	Storage/Transportation	-20°C to 60°C (-4°F to 140°F)
Humidity Range	Operating	<95% relative humidity @ 40°C/104°F
	Storage/Transportation	<90% relative humidity @ 60°C/140°F
Altitude Range	Operating	-500 to 3000 m/-1640 to 9840 ft.
	Storage/Transportation	-500 to 13100 m/-1640 to 43000 ft.

Transducers (M2734A/M2734B/M2735A/M2736A/M2738A)		
Temperature Range	Operating	0°C to 40°C (32°F to 104°F)
	Storage/Transportation	-20°C to 60°C (-4°F to 140°F)
Humidity Range	Operating	<95% relative humidity @ 40°C/104°F
	Storage/Transportation	<90% relative humidity @ 60°C/140°F

28 Specifications and Standards Compliance

Transducers (M2734A/M2734B/M2735A/M2736A/M2738A)

Altitude Range	Operating	-500 to 3000 m/-1640 to 9840 ft.
	Storage/Transportation	-500 to 13100 m/-1640 to 43000 ft.

Avalon CL Transducers (866075/866076/866077)

Temperature Range	Operating	0°C to 40°C (32°F to 104°F)
	Charging	0°C to 35°C (32°F to 95°F)
	Storage/Transportation	-20°C to 60°C (-4°F to 140°F)
Humidity Range	Operating	<95% relative humidity @ 40°C/104°F
	Storage/Transportation	<90% relative humidity @ 60°C/140°F
Altitude Range	Operating	-500 to 3000 m/-1640 to 9840 ft.
	Storage/Transportation	-500 to 13100 m/-1640 to 43000 ft.

SpO₂ Sensors

Operating Temperature Range	0°C to 37°C (32°F to 98.6°F)
-----------------------------	------------------------------

WARNING

Explosion Hazard: Do not use in the presence of flammable anesthetics, such as a flammable anesthetic mixture with air, oxygen, nitrous oxide, or in oxygen rich environment. Use of the devices in such an environment may present an explosion hazard.

Physical Specifications

Avalon CL Base Station		
Dimension and Weight	Size mm/(in): width x height x depth	349 x 74 x 183 mm (349 mm-5% max. width 350 mm) (13.7 x 2.9 x 7.2 in 13.7 in -5% max. width 13.8 in)
	Weight	0.985 kg /2.2 lbs
Electrical Class	When connected to M2702A/ M2703A/M2704A/M2705A Class II equipment using a red connector at the fetal sockets.	When connected to M2704A/M2705A Class I equipment using a black connector in the telemetry socket.
Electrical Power Source	External (powered by fetal monitor)	
Mode of Operation	Continuous operation	
Interface cable connector type	Connector color	Red for FM20/30 left side or FM40/50 front fetal connector socket
		Black for FM40/50 rear telemetry connector socket
Interface cable length	1,5 m /4.11 ft	
Ingress Protection	IP31	

Monitor Physical Specifications		M2702A/M2703A	M2704A/M2705A
Power	Supply Voltages	100 VAC to 240 VAC \pm 10%	
	Supply Frequency Range	50 Hz/60 Hz	
	Power Consumption (current)	0.7 - 0.4 A (M2702A/ M2703A) 1.3 - 0.7 A (M8023A#E25)	1.3 - 0.7 A
Dimensions and Weight	Size (without options) mm/ (in): width x height x depth	286 x 133 x 335 \pm 1% (11.3 x 5.2 x 13.2 in \pm 1%)	420 x 172 x 370 \pm 5% (16.5 x 6.8 x 14.6 in \pm 5%)
	Weight	< 5.1 kg/11.2 lbs	< 9.0 kg/19.8 lbs
	Degree of Protection Against Electrical Shock	Type CF	
Electrical Class	Class II equipment	Class I equipment	
Electrical Power Source	External (AC) without battery option #E25 Internal (LiIo battery) if with battery option #E25	External (AC)	
Mode of Operation	Continuous operation		
Water Ingress Protection Code	-	IP X1 (provided recorder drawer is shut)	

28 Specifications and Standards Compliance

Monitor Physical Specifications		M2702A/M2703A	M2704A/M2705A
Global Speed		6.25 mm/sec, 12.5 mm/sec, 25 mm/sec, 50 mm/sec	
Startup Time	Time taken from switching on the monitor to seeing the first parameter labels	< 30 seconds	

Transducers (M2734A/M2734B/M2735A/M2736A/M2738A)				
Shock Resistance		Withstands a 1 m drop to concrete surface with possible cosmetic damage only		
Water Ingress Protection Code	M2734A&B/35/36A	IP 68 (immersion up to 1 m water depth for 5 hours)		
	M2738A	IP 67 (immersion up to 0.5 m water depth for 30 minutes)		
Dimensions and Weight	M2734A&B/35/36A	Size (diameter)	83 mm/3.27 in	
		Weight (without cable)	< 220 g/7.8 oz.	
	M2738A	Maximum size mm/ (in): width x height x depth	50 x 28 x 135 (2.0 x 1.1 x 5.3 in)	
		Cable length	2.5 m	
		Weight	< 150 g/5.3 oz.	
Degree of Protection Against Electrical Shock		Type CF		
Transducer Identification		Optical Signal Element (Finder LED). Not M2738A		

Avalon CL Transducers (866075/866076/866077)				
Shock Resistance		Withstands a 1.5 m drop to concrete surface with possible cosmetic damage only.		
Water Ingress Protection Code		IP 68 (immersion up to 1 m water depth for 5 hours)		
Dimensions and Weight	Avalon Cableless Toco ⁺ MP Transducer 866075	Size (diameter/height)	75.5 mm/ 36.6 mm 2.9 / 1.4 in	
		Weight	123 g ±10%/4.3 oz ±10%	
	Avalon Cableless US Transducer 866076	Size (diameter/height)	75.5 mm/ 36.6 mm 2.9 / 1.4 in	
		Weight	131 g ±10%/4.6 oz ±10%	
	Avalon Cableless ECG/IUP Transducer 866077	Size (diameter/height)	75.5 mm/ 36.6 mm 2.9 / 1.4 in	
		Weight	121 g ±10%/4.2 oz ±10%	
Degree of Protection Against Electrical Shock		Type CF		
Transducer Identification		Optical Signal Element (Finder LED).		

Interface Cable for Avalon CTS (M2731-60001 and M2732-60001)		
Shock Resistance	Withstands a 1m drop to concrete surface with possible cosmetic damage only	
Water Ingress Protection Code	IP X1	
Dimensions and Weight	Maximum size mm/(in): width x height x depth	55 x 28 x 50 (2.2 x 1.1 x 2.0 in)
	Cable length	2.5 m
	Weight	< 200 g/7.0 oz.

M8023A (Option #E25) External Power Supply Weight and Dimensions	
Maximum Weight	700 g \pm 10% (1.5 lb \pm 10%)
Size (W x H x D)	207 x 70 x 113 mm / 8.1 x 2.8 x 4.4 in

Interface Specifications		
Network	Standard	100-Base-TX (IEEE 802.3 Clause 25)
	Connector	RJ45 (8 pin)
	Isolation	Basic insulation (reference voltage: 250 V; test voltage: 1500 V)
MIB/RS232	Standard	IEEE 1073-3.2-2000
	Connectors	RJ45 (8 pin)
	Mode	Software-controllable BCC (Rx/D/TxD cross over) or DCC (Rx/D/TxD straight through)
	Power	5 V \pm 5 %, 100 mA (max.)
	Isolation	Basic insulation (reference voltage: 250 V; test voltage: 1500 V)
USB Interface	Standard	USB 2.0 full-speed (embedded host)
	Connectors	USB series "Standard A" receptacle
	Power	Low power port 4.4V min; max. load for all ports together 500 mA
	Isolation	none
RS232 (Standard)	Connectors	RJ45 (8-pin)
	Power	none
	Insulation	Basic insulation (reference voltage: 250 V; test voltage: 1500 V)
RS232 (Independent display interface option)	Connectors	RJ45 (8-pin)
	Power	none
	Isolation	none
Basic Nurse Call Relay	Connectors	Modular Jack 6P6C, active open and closed contact
	Contact	\leq 100 mA, \leq 24 V DC
	Isolation	Basic insulation (reference voltage: 250 V; test voltage: 1500 V)
	Delay	<[Configured Latency +0.5] sec

Interface Specifications		
Flexible Nurse Call Relay ¹	Connectors	20 pin MDR (Mini D-Ribbon), active open and closed contacts
	Contact	<=100 mA, <=24 V DC
	Isolation	Basic insulation (reference voltage: 250 V; test voltage: 1500 V)
	Delay	<[Configured Latency +0.5] sec

1 The power loss indication functionality of the Nurse Call Relay board is not supported with fetal monitors.

Avalon CL Radio Interface Specifications		
Short Range Radio Interface	Type	Internal SRR interface
	Technology	IEEE 802.15.4
	Frequency Band	2.4 GHz ISM (2.400 - 2.483 GHz)
	Modulation Technique	DSSS (O -QPSK)
	Effective radiated power	max. 0 dBm (1 mW)
	Range	ca. 5 m without any physical obstructions as walls and doors
OBR (WMTS)	Frequency Band	608 - 614 MHz
	Effective radiated power	<20 mW (Base Station) <1 mW (CL transducers)
	Range	min. 70 m
OBR (ISM)	Frequency Band	433.05 -434.79 MHz
	Effective radiated power	<10 mW (Base Station) <1 mW (CL transducers)
	Range	min. 70 m
OBR (T108)	Frequency Band	920.6 - 923.4 MHz
	Effective radiated power	<20 mW (Base Station) <10 mW (CL transducers)
	Range	min. 70 m

Performance Specifications

Battery Specifications

Performance Specifications		
Avalon FM20/30 Battery Option #E25	Operating Time (with new, fully charged battery)	Basic monitoring configuration: >2 hours (Display Brightness: 70%, Recorder: "On" at 3cm/min, NBP: Auto Mode at 15min, 2 US Transducers, 1 Toco+ with MECG, 1 Patient Module with DECG)
	Charge Time	When monitor is off: approx. 6 hours When monitor is in use: more than 10 hours (depending on monitor configuration).

Performance Specification		
Avalon CL Transducer Battery	Operating time	With a new and fully charged battery min. 8 hours
	Charge time	From a "low battery" indication to a "fully charged" indication min. 3 hours

Fetal / Maternal Specifications

Performance Specifications		
Ultrasound		
Measurement Method		Ultrasound Pulse Doppler
Measurement Range	US	50 to 240 bpm
Resolution	Display	1 bpm
	Printer	1/4 bpm
Jitter @ 200 bpm		≤ 3 bpm
Display Update Rate		1 / second
US Intensity (M2736A/AA)	Average output power	$P = (7.4 \pm 0.4) \text{ mW}$
	Peak-negative acoustic pressure	$p_- = (40.4 \pm 4.3) \text{ kPa}$
	Output beam intensity (I_{ob}) (= spatial average - temporal average intensity)	$I_{sata} = (2.38 \pm 0.59) \text{ mW/cm}^2$
	Spatial-peak temporal average intensity	$I_{spta} = (15.0 \pm 3.2) \text{ mW/cm}^2$
	Effective radiating area @ -12 dB	$(3.11 \pm 0.74) \text{ cm}^2$
	Thermal index (TI) and mechanical index (MI) are always below 1.0.	

28 Specifications and Standards Compliance

Performance Specifications		
US Intensity CL (866076)	Average output power	$P = (12.4 \pm 0.4) \text{ mW}$
	Peak-negative acoustic pressure	$p_- = (49.1 \pm 5.2) \text{ kPa}$
	Output beam intensity (I_{ob}) (= spatial average - temporal average intensity)	$I_{sata} = (2.77 \pm 0.56) \text{ mW/cm}^2$
	Spatial-peak temporal average intensity	$I_{spta} = (21.1 \pm 5.1) \text{ mW/cm}^2$
	Effective radiating area @ -12 dB	$A_{-12dB} = (4.47 \pm 0.89) \text{ cm}^2$
	Thermal index (TI) and mechanical index (MI) are always below 1.0.	
Signal Quality Indication	Poor Quality	empty
	Acceptable Quality	half-full
	Good Quality	full
Beat-to -Beat (max.) for Ultrasound		28bmp
US Frequency		1 MHz \pm 100 Hz
US Signal range		3.5 μVpp to 350 μVpp @ 200 Hz
US Burst	Repetition Rate	3.0 kHz
	Duration	$\leq 100 \mu\text{s}$
FMP Signal Range @ 33 Hz		200 μVpp to 40 μVpp

Performance Specifications		
Toco		
Measurement Method		Strain Gauge Sensor Element
Sensitivity		1 unit = 2.5 g
Resolution	Display	1 unit
	Printer	1/4 unit
Measurement Range		400 units
Signal Range		0 to 127 units
Maximum Offset Range		-300 units
Baseline Setting		20 units
Update Rate	Display	1 / second
	Printer	~ 4 / seconds
Auto Offset Correction		3 seconds after connecting the transducer, the Toco value is set to 20 units
Auto Zero Adjust		Toco value is set to zero following a negative measurement value for 5 seconds

Performance Specifications	
Maternal Pulse from Toco	
Emitted Light Energy	$\leq 15 \text{ mW}$

Performance Specifications	
Wavelength Range	780 to 1100 nm
Range	40 to 240 bpm
Resolution	1 bpm
Display Update Rate	1 / s
Accuracy	± 2% or 1 bpm, whichever is greater
Update Rate	every 4 seconds

Performance Specifications		
IUP		
Measurement Method	Passive Resistive Strain Gauge Elements	
Measurement Range	-100 to +300 mmHg	
Signal Range	-99 to 127 mmHg or -13.2 to 16.9 kPa	
Resolution	Display	1 mmHg
	Printer	1/4 mmHg
Sensitivity	5 μ V/V/mmHg	
Offset Compensation	+100 to -200 mmHg	
Accuracy (not including sensor accuracy)	± 0.5% per 100 mmHg	
Update Rate	Display	1 / second
	Printer	~4 / seconds
Auto Offset Correction	3 seconds after connecting the transducer, the IUP value is set to 0 mmHg	

Performance Specifications		
ECG		
Performance Specifications Type	DECG	Single Lead ECG (derived from Fetal Scalp Electrode)
	MECG	Single Lead ECG (derived from RA and LA electrodes)
Measurement Range	30 to 240 bpm	
Resolution	Display	1 bpm (display update rate 1/s)
	Recorder	1/4 bpm
Wave Speed	6.25 mm/sec, 12.5 mm/sec, 25 mm/sec, 50 mm/sec	
Accuracy	± 1 bpm or 1%, whichever is greater (non-averaging)	
Beat-to -Beat (max.)	MECG: 28 bpm DECG: 28 bpm (with Artifact Suppression On)	
Differential Input Impedance	> 15M Ω	

Performance Specifications		
Electrode Offset Potential Tolerance		± 400 mV
Filter Bandwidth		0.8 to 80 Hz
INOP Auxiliary Current (Leads Off Detection)		< 100 µA
Input Signal Range	DECG	20 µVpp to 6 mVpp
	MECG	150 µVpp to 6 mVpp
Dielectric Strength		1500 Vrms
Defibrillator Protection		None
ESU Protection		None
Pace pulse detection		None

WARNING

The fetal/maternal monitor is not a diagnostic ECG device. In particular, the display of fetal/maternal ECG is intended only for evaluating signal quality for fetal/maternal heart rate as derived from the ECG waveform.

When in doubt, it can be used to identify sources of compromised signal quality, such as noise or muscle artifacts. It can subsequently be used to verify the result of measures taken to resolve them (e.g., checking ECG cable connections or adapting the fetal **ArtifactSuppress** configuration).

The safety and effectiveness of the displayed fetal/maternal ECG waveform (i.e., P, QRS and T segments) for evaluation of fetal/maternal cardiac status during labor have not been evaluated.

Fetal Heart Rate (Ultrasound/DECG) Alarm Specifications			
FHR Alarm Limits	Range	Bradycardia (low limit)	60 to 200 bpm adjustable in 10 bpm steps Default: 110 bpm
		Tachycardia (high limit)	70 to 210 bpm adjustable in 10 bpm steps Default: 170 bpm
FHR Alarm Delay	Range	Bradycardia (low limit) Delay	10 to 300 seconds in steps of 10s Default: 240s
		Tachycardia (high limit) Delay	10 to 300 seconds in steps of 10s Default: 300s
		Signal Loss Delay	10 to 30 seconds in steps of 10s

MECG Alarm Specifications	Range	Adjustment
MECG Alarm Limits	High Range: 31 to 240 Default: 120 bpm	1 bpm steps (30 to 40 bpm) 5 bpm steps (40 to 240 bpm)
	Low Range: 30 to 235 Default: 50 bpm	
Tachycardia	Difference to high limit: 0 to 50 bpm Default: 20 bpm	5 bpm steps
	Clamping at: 150 to 240 bpm Default: 200 bpm	5 bpm steps
Bradycardia	Difference to low limit: 0 to 50 bpm Default: 20 bpm	5 bpm steps
	Clamping at: 30 to 100 bpm Default: 40 bpm	5 bpm steps

Maternal ECG Supplemental Information as required by IEC 60601-2-27

Heart Rate Averaging Method	The maternal heart rate is computed by averaging the 12 most recent RR intervals. If each of 3 consecutive RR intervals is greater than 1200 ms (i.e. rate less than 50 bpm), then the 4 most recent RR intervals are averaged to compute the HR.
Display Update Rate	2 s
Response Time of Heart Rate Meter to Change in Heart Rate	HR change from 80 to 120 bpm: 10 sec. HR change from 80 to 40 bpm: 14 sec.

Noninvasive Blood Pressure

Complies with IEC 80601-2-30:2010/EN 80601-2-30:2010.

Performance Specifications		
Measurement Ranges	Systolic	30 to 270 mmHg (4 to 36 kPa)
	Diastolic	10 to 245 mmHg (1.5 to 32 kPa)
	Mean	20 to 255 mmHg (2.5 to 34 kPa)
Accuracy		Max. Std. Deviation: 8 mmHg (1.1 kPa) Max. Mean Error: ± 5 mmHg (± 0.7 kPa)
Pulse Rate	Range	40 to 300 bpm
	Accuracy (average over noninvasive blood pressure measurement cycle)	40 to 100 bpm: ± 5 bpm
		101 to 200 bpm: $\pm 5\%$ of reading 201 to 300 bpm: $\pm 10\%$ of reading
Measurement Time		Typical at HR > 60bpm Auto/manual: 30 seconds (adult) Maximum time: 180 seconds (adult)

Performance Specifications	
Cuff Inflation Time	Typical for normal adult cuff: Less than 10 seconds
Initial Cuff Inflation Pressure	165 ±15 mmHg
Auto Mode Repetition Times	1, 2, 2.5, 3, 5, 10, 15, 20, 30, 45, 60 or 120 minutes
Venipuncture Mode Inflation	
Inflation Pressure	20 to 120 mmHg (3 to 16 kPa)
Automatic deflation after	170 seconds

Measurement Validation: In adult mode, the blood pressure measurements determined with this device comply with the American National Standard for Electronic or Automated Sphygmomanometers (ANSI/AAMI SP10:2002/A1:2003) in relation to mean error and standard deviation, when compared to intra-arterial or auscultatory measurements (depending on the configuration) in a representative population. For the auscultatory reference the 5th Korotkoff sound was used to determine the diastolic pressure.

Alarm Specifications	Range	Adjustment
Systolic	Adult: 30 to 270 mmHg (4 to 36 kPa)	10 to 30 mmHg: 2 mmHg (0.5 kPa) > 30 mmHg: 5 mmHg (1 kPa)
Diastolic	Adult: 10 to 245 mmHg (1.5 to 32 kPa)	
Mean	Adult: 20 to 255 mmHg (2.5 to 34 kPa)	

Overpressure Settings	Adjustment
> 300 mmHg (40 kPa) > 2 sec	not user adjustable

SpO₂

Complies with EN/ISO 80601-2-61 2011 (except alarm system; alarm system complies with IEC 60601-1-8:2006).

Measurement Validation: The SpO₂ accuracy has been validated in human studies against arterial blood sample reference measured with a CO-oximeter. Pulse oximeter measurements are statistically distributed, only about two-thirds of the measurements can be expected to fall within the specified accuracy compared to CO-oximeter measurements.

Display Update Period: Typical: 2 seconds, Maximum: 30 seconds. Max. with noninvasive blood pressure INOP suppression on: 60 seconds.

SpO ₂ Performance Specifications		
SpO ₂ The specified accuracy is the root-mean-square (RMS) difference between the measured values and the reference values	Range	0 to 100%
	Accuracy	<p>Philips Reusable Sensors: M1191A/B, M1191AL/BL, M1191ANL, M1192A, M1192AN = 2% (70% to 100%) M1191T, M1192T, M1194A, M1194AN, M1196A, M1196T = 3% (70% to 100%)</p> <p>Philips Disposable Sensors with M1943A(L): M1131A, M1901B, M1903B, M1904B = 3% (70% to 100%) M1133A, M1134A = ±2% (70% to 100%)</p> <p>NellcorPB® Sensors with M1943A(L): MAX-A, MAX-AL, MAX-P, MAX-N, D-25, D-20, N-25, OxiCliq A, P, N = 3% (70% to 100%)</p> <p>Masimo Reusable Sensors® with LNOP MP12 or LNC MP10: LNOP DC-I, LNOP DC-IP, LNOP YI, LNCS DC-I, LNCS DC-IP: 2% (70% to 100%) LNOP TC-I, LNCS TC-I: 3.5% (70% to 100%)</p> <p>Masimo Disposable Sensors® with LNOP MP12 or LNC MP10: LNOP Adt, LNOP Adtx, LNOP Pdt, LNOP Pdtx, LNCS Adtx, LNCS Pdtx: 2% (70% to 100%) LNOP Neo-L, LNCS Neo-L: 3% (70% to 100%)</p>
	Resolution	1%
Pulse	Range	30 to 300 bpm
	Accuracy	±2% or 1 bpm, whichever is greater
	Resolution	1 bpm
Sensors	Wavelength range	500 to 1000 nm. Information about the wavelength range can be especially useful to clinicians (for instance, when photodynamic therapy is performed).
	Emitted Light Energy	≤ 15mW
Pulse Oximeter Calibration Range		70% - 100%

SpO ₂ Alarm Specifications	Range	Adjustment	Delay
SpO ₂	50 to 100%	1% steps	(0, 1, 2, 3,... 30) + 4 seconds
Desat	50 to Low alarm limit	1% steps	
Pulse	30 to 300 bpm	1 bpm steps (30 to 40 bpm) 5 bpm steps (40 to 300 bpm)	max. 14 seconds

SpO ₂ Alarm Specifications	Range	Adjustment	Delay
Tachycardia	Difference to high limit 0 to 50 bpm	5 bpm steps	max. 14 seconds
	Clamping at 150 to 300 bpm	5 bpm steps	
Bradycardia	Difference to low limit 0 to 50 bpm	5 bpm steps	max. 14 seconds
	Clamping at 30 to 100 bpm	5 bpm steps	

Tympanic Temperature

Complies with:

- EN 12470-5 (Clinical thermometers - Part 5:2003: Performance of infra-red thermometers)
- ASTM E1965-98 (Infrared Thermometers for Intermittent Determination of Patient Temperature)

with minor exceptions as noted below.

Performance Specifications

Performance Specifications		
Temperature Resolution	0.1 °C or 0.1 °F	
Response Time	less than 2 seconds	
Temperature Calibrated Accuracy Specifications (out of the Factory)		
Ambient Temperature	Target Temperature	Accuracy
25°C (77°F)	37.7°C to 38.9°C (98.4°F to 102°F)	±0.1°C (±0.2°F)
16°C to 33°C (60.8° to 91.4°F)	33°C to 42°C (91.4°F to 107.6°F)	±0.2°C (±0.4°F)
Temperature Calibrated Accuracy Specifications (after recalibration using Genius 2 Checker/Calibrator)		
Ambient Temperature	Target Temperature	Accuracy
16°C to 33°C (60.8° to 91.4°F)	36°C to 39°C (96.8°F to 102.2°F)	±0.2°C (±0.4°F)
16°C to 33°C (60.8° to 91.4°F)	Less than 36°C or greater than 39°C (Less than 96.8°F or greater than 102.2°F)	±0.3°C (±0.5°F)
<p>ASTM laboratory requirement for IR thermometers in the display range 37°C to 39°C (98°F to 102°F) is ±0.2°C (±0.4°F), whereas for mercury-in-glass and electronic thermometers, the requirement per ASTM standards E667-86 and E1112-86 is ±0.1°C (±0.2°F).</p> <p>Clinical accuracy characteristics and procedures are available from Covidien llc on request. To verify the accuracy, use a certified black body as specified in EN 12470-5:2003, Annex C, or use a Genius 2 Checker/Calibrator - available from Covidien llc under part number 303097.</p> <p>Clinical repeatability: meets section A.5 of EN 12470-5:2003 (E) per Covidien llc technical report. Data is available from Covidien llc on request.</p>		

Displayed Temperature Measurement Range		
Mode	Range °C	Range °F
Ear	33.0 to 42.0	91.4 to 107.6
Oral (ear + 0.6 °C)	33.6 to 42.0	92.5 to 107.6
Core (ear + 1.04 °C)	34.0 to 42.0	93.2 to 107.6
Rectal (ear + 1.16 °C)	34.2 to 42.0	93.6 to 107.6
Caution: ASTM E1965-98 specifies 34.4°C to 42.2°C (94°F to 108°F)		

Ambient Temperature Range		
Mode	Range °C	Range °F
Operating 10 to 95% RH, non-condensing	16 to 33	60.8 to 91.4
Storage up to 95% RH, non-condensing	-25 to 55.0	-13 to 131
<p>Caution: EN12470-5 specifies 16.0°C to 35.0°C (60.8° F to 95° F), 10 to 95% RH, non-condensing ASTM E1965-98 specifies 16.0°C to 40.0°C (60.8° F to 104° F), up to 95% RH, non-condensing</p> <p>Storing the thermometer outside the specified temperature/humidity range might adversely affect measurement accuracy. Check the calibration after storage in uncertain conditions.</p>		

Tympanic Temperature Alarm Specifications	
Range	33.0°C to 42.0° C (91° F to 108° F)
Adjustment	0.5° C steps (33° C to 35° C) 0.1° C steps (35° C to 42° C) 1.0° F steps (91° F to 95° F) 0.2° F steps (95° F to 108° F)





Physical Specifications















Thermometer	
Dimensions	190 mm x 43 mm x 55 mm (±3 mm)
Cable length	60 ±5 cm (spiral cable relaxed) 250 ±15 cm (spiral cable extended)
Weight (including cable)	180 ±10 g




Base Station	
Dimensions	205 mm x 65 mm x 75 mm (±3 mm)
Weight (excluding cable)	400g ±10 mm

Recorder Specifications

Built-in Thermal Array Fetal Trace Recorder			
Mechanism	Thermal Array Recorder		
Paper & Printing	Type	Standard Z-fold paper	
	Standard Speeds (real-time traces)	3 cm/min, 2 cm/min, 1cm/min	
	Fast Print Speed (stored traces)	Max. 20 mm/s Print speed is variable and depends on the print load	
	ECG Wave Print Speed (not real-time)	Emulated 25 mm/s Print speed is variable and depends on the print load	
	Paper Advance	20 mm/s	
	Sensing	Optical Reflex Sensor for black page marks	
Accuracy @ 3 cm/min, 2 cm/min, 1 cm/min	±5 mm/page		
Usable Print Width	128 mm		
Resolution	8 dots/mm (200 dpi)		
Time Delay to see trace on paper	<30s @ 1 cm/min		
Trace Separation Offset for FHR (Ultrasound and DECG)	Twin	Standard	FHR2 +20 bpm
		Classic	FHR1 +20 bpm
	Triplet	Standard	FHR2 +20 bpm FHR3 -20 bpm
		Classic	FHR1 +20 bpm FHR3 -20 bpm

Recorder Symbols	
Symbol	Description
	Parameter is capable of alarming and alarms were enabled at the time of printing the annotation. Low / high limit are printed surrounding the symbol.
	Parameter is capable of alarming but alarms were disabled at the time of printing the annotation. (Note: There is no alarm related annotation at all if a parameter does not have alarming capability.)
	Alarm is on (printed next to measurement label)
	Upper and lower alarm limit (printed next to measurement label)

Recorder Symbols	
Symbol	Description
	FMP detection is on
	Beginning of the date/time annotation
	Warning (INOP)
	Measurement from a cableless transducer (printed next to measurement label)
	Measurement from a cableless maternal measurement Pod (SpO ₂ or NBP)
	Pulse from SpO ₂
	Pulse from Toco MP
	Pulse from NBP
	Trace separation +20 bpm (in label)
	Trace separation -20 bpm (in label)
	Trace separation Off (in trace)
	Trace separation +20 bpm (in trace)
	Trace separation -20 bpm (in trace)
	Trace separation +20 bpm and -20 bpm (in trace)

Recorder Symbols	
Symbol	Description
	Coincidence of heart rates is detected
	Marker
	Special wave, with different speed and scale (for example, fast printout of MECCG wave on FM30)

Compatible External Displays: FM40/FM50 Only

External displays can be connected with a maximum cable run of 10 m.

Compatible Display Specifications	External XGA Display (M8031B)	External SXGA Display (M8033C)
Resolution*	1024 x 768	1280 x 1024 pixel
Refresh frequency	60 Hz or 75 Hz	60 Hz
Useful screen	Depends on size of display	
Pixel size		

* The video output of the Avalon FM40/FM50 has VGA resolution.

Manufacturer's Information

You can write to Philips at this address:

Philips Medizin Systeme Boeblingen GmbH

Hewlett-Packard-Str. 2

71034 Boeblingen

Germany

Visit our website for local contact information at: www.healthcare.philips.com

© Copyright 2013. Koninklijke Philips N.V. All Rights Reserved.

Trademark Acknowledgment

Oxisensor™ II, Oxi-Cliq™, and OxiMax™ are trademarks of Tyco Healthcare Group LP, Nellcor Puritan Bennett Division.

Regulatory and Standards Compliance

The monitor is in conformity with the requirements of the European Medical Devices Directive 93/42/EEC and bears the CE marking:



The monitor is classified into Class IIb to Annex IX rule 10.

Safety and Performance

The monitor complies with the following major international safety and performance standards:

- IEC 60601-1:2005+A1:2012 / EN 60601-1:2006+AC:2010
- IEC 60601-1-6:2010 / EN 60601-1-6:2010
- IEC 60601-1-8:2006 / EN 60601-1-8:2007+AC:2010
- IEC 60601-2-27:2011+AC:2012 / EN 60601-2-27:2006+AC:2006
- IEC 80601-2-30:2009+C1:2010 / EN 80601-2-30:2010
- IEC 60601-2-37:2007 / EN 60601-2-37:2008
- IEC 60601-2-49:2011 / EN 60601-2-49:2001
- ISO 80601-2-61:2011 / EN 80601-2-61:2011
- ANSI/AAMI ES60601-1+C1:2009+A2:2010
- CAN/CSA C22.2#60601-1-08+TC2:2011
- JIS T 1303 2005
- AS/NZS 3200.1.0-1998

The possibility of hazards arising from hardware and software errors was minimized in compliance with ISO 14971:2012, IEC 60601-1:2005+A1:2012 / EN 60601-1:2006+AC:2010.

Alarm sounds are compliant with Standard ISO/IEC 9703-2 and IEC 60601-1-8:2006 / EN 60601-1-8:2007+AC:2010.

Safety Tests Fetal Monitor

All the safety tests and procedures required after an installation or an exchange of system components are described in your monitor's Service Guide. These safety tests are derived from international standards but may not be sufficient to meet local requirements.

WARNING

- Do not use additional AC mains extension cords or multiple portable socket-outlets. If a multiple portable socket-outlet is used, the resulting system must be compliant with IEC/EN 60601-1 A1: 2012.
 - Do not connect any devices that are not supported as part of a system.
 - Do not use a device in the patient vicinity if it does not comply with IEC/EN 60601-1-1. The whole installation, including devices outside of the patient vicinity, must comply with IEC/EN 60601-1 A1: 2012. Any non-medical device, including a PC running an OB TraceVue/IntelliSpace Perinatal system, placed and operated in the patient's vicinity must be powered via a separating transformer (compliant with IEC/EN 60601-1 A1: 2012) that ensures mechanical fixing of the power cords and covering of any unused power outlets.
 - Do not use USB devices with own power supplies unless an appropriate separation device is used (either between USB interface and device or between device and power).
-

During the installation the fetal monitor is configured for your environment. This configuration defines your custom default settings you work with when you switch on your fetal monitor. See the Fetal Monitor's Service Guide and the Configuration Guide for details on how to configure your fetal monitor.

Electromagnetic Compatibility (EMC)

The device and its accessories, listed in the accessories section, comply with the following EMC standards:

- EN 60601-1-2+AC:2010 (IEC 60601-1-2:2007)

Take special precautions regarding electromagnetic compatibility (EMC) when using medical electrical equipment. You must operate your monitoring equipment according to the EMC information provided in this book. Before using the device, assess the electromagnetic compatibility of the device with surrounding equipment.

This ISM device complies with Canadian ICES-003:2012. Cet appareil ISM est conforme à la norme NMB-003 du Canada.

CAUTION

- **FM20/FM30 only:** Although this is an electrical Class II device, it has a protective earth conductor which is needed for EMC purposes.
 - Always use the supplied power cord with the three-prong plug to connect the monitor to AC mains. Never adapt the three-prong plug from the power supply to fit a two-slot outlet.
-

CAUTION

The use of accessories, transducers and cables other than those specified may result in increased electromagnetic emissions or decreased electromagnetic immunity of the device.

WARNING

Do not use cordless/mobile phones or any other portable RF communication system within the patient vicinity, or within a 1.0 m radius of any part of the fetal monitoring system.

WARNING

For paced patients: The radiated SRR power of the CL SpO₂ and CL NBP Maternal Cableless Measurement Devices, and other sources of radio-frequency energy, when used in very close proximity of a pacemaker, might be sufficient to interfere with pacemaker performance. Due to shielding effects of the body, internal pacemakers are somewhat less vulnerable than external pacemakers. However, caution should be exercised when monitoring paced patients.

In order to minimize the possibility of interference, avoid positioning and wearing the Cableless Measurement Devices in very close proximity to a pacemaker. Consult the pacemaker manufacturer for information on the RF susceptibility of their products

EMC Testing

CAUTION

Fetal parameters, especially ultrasound and ECG, are sensitive measurements involving small signals, and the monitoring equipment contains very sensitive high gain front-end amplifiers. Immunity levels for radiated RF electromagnetic fields and conducted disturbances induced by RF fields are subject to technological limitations. To ensure that external electromagnetic fields do not cause erroneous measurements, it is recommended to avoid the use of electrically radiating equipment in close proximity to these measurements.

Reducing Electromagnetic Interference

CAUTION

The device should not be used adjacent to, or stacked with, other equipment unless otherwise specified.

The product and associated accessories can be susceptible to interference from continuous, repetitive, power line bursts, and other RF energy sources, even if the other equipment is compliant with EN 60601-1-2 emission requirements. Examples of other sources of RF interference are other medical electrical devices, cellular products, information technology equipment, and radio/television transmissions.

When electromagnetic interference (EMI) is encountered, for example, if you can hear spurious noises on the fetal monitor's loudspeaker, attempt to locate the source. Assess the following:

- Is the interference due to misplaced or poorly applied transducers? If so, re-apply transducers correctly according to directions in this book or in the Instructions for Use accompanying the accessory.
- Is the interference intermittent or constant?
- Does the interference occur only in certain locations?
- Does the interference occur only when in close proximity to certain medical electrical equipment?

Once the source is located, there are a number of things that can be done to mitigate the problem:

- 1 Eliminating the source. Turn off or move possible sources of EMI to reduce their strength.
- 2 Attenuating the coupling. If the coupling path is through the patient leads, the interference may be reduced by moving and/or rearranging the leads. If the coupling is through the power cord, connecting the system to a different circuit may help.
- 3 Adding external attenuators. If EMI becomes an unusually difficult problem, external devices such as an isolation transformer or a transient suppressor may be of help. Your service provider can be of help in determining the need for external devices.

Where it has been established that electromagnetic interference is affecting physiological parameter measurement values, a physician, or a suitably qualified person authorized by a physician, should determine if it will negatively impact patient diagnosis or treatment.

System Characteristics

The phenomena discussed above are not unique to this system but are characteristic of fetal patient monitoring equipment in use today. This performance is due to very sensitive high gain front end amplifiers required to process the small physiological signals from the patient. Among the various monitoring systems already in clinical use, interference from electromagnetic sources is rarely a problem.

Electromagnetic Emissions and Immunity

The EMC standards state that manufacturers of patient-coupled equipment must specify immunity levels for their systems. See Tables 1 to 4 for this detailed immunity information. See Table 5 for recommended minimum separation distances between portable and mobile communications equipment and the product.

Immunity is defined in the standard as the ability of a system to perform without degradation in the presence of an electromagnetic disturbance.

Caution should be exercised in comparing immunity levels between different devices. The criteria used for degradation are not always specified by the standard and can therefore vary with the manufacturer.

In the table below, the term "device" refers to the Avalon FM20/30/40/50 fetal monitor together with its accessories. The table gives details of the electromagnetic emissions, and how these are classified, for the device, and the electromagnetic environments in which the device is specified to technically function.

Table 1 - Guidance and Manufacturer's Declaration: Electromagnetic Emissions		
Emissions Test	Compliance	Avoiding Electromagnetic Interference
Radiofrequency (RF) emissions	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations and flicker IEC 61000-3-3	complies	
RF emissions CISPR 11 For the Avalon FM20/30 fetal monitor with all accessories except the IUP/ECG patient module M2738A.	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage supply network that supplies buildings used for domestic purposes ¹ .
RF emissions CISPR 11 For the Avalon FM40/FM50 with all accessories. For the Avalon FM20/30 fetal monitor whenever used with the IUP/ECG patient module M2738A. For the Avalon CTS Interface Cable (M2731-60001/M2732-60001) whenever used with the Avalon CTS Cableless Fetal Transducer System. For the Avalon CL Base Station with cableless transducers whenever used with the fetal monitors.	Class A	The device is suitable for use in all establishments other than domestic establishments and those directly connected to the public low-voltage supply network that supplies buildings used for domestic purposes.

¹ Note that the device is not intended for home use.

Electromagnetic Immunity

The monitor is suitable for use in the specified electromagnetic environment. The user must ensure that it is used in the appropriate environment as described below.

Table 2 - Guidance and Manufacturer's Declaration: Electromagnetic Immunity			
Immunity Test	IEC 60601-1-2 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8kV air	± 6 kV contact ± 8kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial and/or hospital environment
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial and/or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (> 95% dip in U_T) for 0.5 cycles	<5% U_T (> 95% dip in U_T) for 0.5 cycles	Mains power quality should be that of a typical commercial and/or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device is powered from an uninterruptible power supply.
	40% U_T (60% dip in U_T) for 5 cycles	40% U_T (60% dip in U_T) for 5 cycles	
	70% U_T (30% dip in U_T) for 25 cycles	70% U_T (30% dip in U_T) for 25 cycles	
	< 5% U_T (> 95% dip in U_T) for 5 sec	< 5% U_T (> 95% dip in U_T) for 5 sec	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial and/or hospital environment
Key: U_T is the AC mains voltage prior to application of the test level.			

Radio Compliance Notice

Avalon CL with WMTS

Operation of this equipment requires the prior coordination with a frequency coordinator designated by the FCC for the Wireless Medical Telemetry Service.

This device complies with Part 15 of the FCC Rules and RSS-210 of Industry Canada. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Any changes or modifications to this equipment not expressly approved by Philips Medical Systems may cause harmful radio frequency interference and void your authority to operate this equipment.

Installation of this telemetry device is permitted in hospitals and health care facilities only. This device shall not be operated in mobile vehicles (including ambulances and other vehicles associated with health care facilities). The installer/user of this device shall ensure that it is at least 80 km from the Dominion Radio Astrophysical Observatory (DRAO) near Penticton, British Columbia. The coordinates of DRAO are: latitude N 49° 19' 15", longitude W 119° 37' 12". For medical telemetry systems not meeting this 80 km separation (e.g. the Okanagan Valley, British Columbia) the installer/user must coordinate with, and obtain the written concurrence of, the Director of DRAO before the equipment can be installed or operated. The Director of DRAO may be contacted at 250-497-2300 (telephone) or 250-497-2355 (fax). (Alternatively, the Manager, Regulatory Standards, Industry Canada, may be contacted.)

Le présent appareil est conforme aux CNR d'Industrie Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes: (1) l'appareil ne doit pas produire de brouillage, et (2) l'utilisateur de l'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement.

L'utilisation de cet appareil de télémétrie est permise seulement dans les hôpitaux et établissements de soins de santé. Cet appareil ne doit pas être mis en marche dans des véhicules (y compris les ambulances et autres véhicules associés aux établissements de santé). La personne qui installe/utilise cet appareil doit s'assurer qu'il se trouve à au moins 80 km de l'Observatoire fédéral de radioastrophysique (OFR) de Penticton en Colombie-Britannique. Les coordonnées de l'OFR sont : latitude N 49° 19' 15», longitude O 119° 37' 12 ". La personne qui installe/utilise un système de télémétrie médicale ne pouvant respecter cette distance de 80 km (p. ex. dans la vallée de l'Okanagan (Colombie-Britannique), doit se concerter avec le directeur de l'OFR et obtenir de sa part une autorisation écrite avant que l'équipement ne puisse être installé ou mis en marche. Le directeur de l'OFR peut être contacté au 250-497-2300 (tél.) ou au 250-497-2355 (télécopieur). (Le Directeur des Normes réglementaires d'Industrie Canada peut également être contacté).

Avalon CL with T108

Japanese Radio Law and Japanese Telecommunications Business Law Compliance.

This device is granted pursuant to the Japanese Radio Law (電波法)
and the Japanese Telecommunications Business Law (電気通信事業法).

This device should not be modified (otherwise the granted designation number will become invalid).

Finding Recommended Separation Distances

In the following table, P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, and d is the recommended separation distance in meters (m).

Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter.

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.



Interference may occur in the vicinity of equipment marked with this symbol:

Table 3 - Guidance and Manufacturer's Declaration: Electromagnetic Immunity		
Conducted RF Immunity Test EN/IEC 61000-4-6		
IEC 60601-1-2 Test Level over 150 kHz to 80 MHz	Compliance Level	Electromagnetic Environment Guidance: Recommended Separation Distance (d) (in Meters, at Frequency Range Tested) for Ultrasound and ECG Measurements
3.0 V _{RMS}	3.0 V _{RMS}	$d = 1, 2\sqrt{P}$
Key: d = Recommended separation distance in meters (m)		
P = maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer		
$V1$ = Tested compliance level (in Volts) for the Conducted RF Immunity test IEC 61000-4-6		
The device meets the compliance level of 3.0 V _{RMS} according to IEC 60601-1-2 over the specified test frequency range. Over the frequency range 150 kHz to 80 MHz, the recommended separation distance in meters (d) is found by the following equation:		
$d = \left(\frac{3,5}{V1}\right)\sqrt{P}$	For a compliance level of 3.0 V _{RMS} :	$d = 1, 2\sqrt{P}$

Table 4 - Guidance and Manufacturer's Declaration: Electromagnetic Immunity		
Radiated RF Immunity Test EN/IEC 61000-4-3		
IEC 60601-1-2 Test Level over 80 MHz to 2.5 GHz	Compliance Level	Electromagnetic Environment Guidance: Recommended Separation Distance (d) (in Meters, at Frequency Range Tested) for Ultrasound and ECG Measurements
3.0 V/m	3.0 V/m	Over 80 MHz to 800 MHz: $d = 1, 2\sqrt{P}$
		Over 800 MHz to 2.5 GHz: $d = 2, 3\sqrt{P}$

Table 4 - Guidance and Manufacturer's Declaration: Electromagnetic Immunity		
Radiated RF Immunity Test EN/IEC 61000-4-3		
IEC 60601-1-2 Test Level over 80 MHz to 2.5 GHz	Compliance Level	Electromagnetic Environment Guidance: Recommended Separation Distance (<i>d</i>) (in Meters, at Frequency Range Tested) for Ultrasound and ECG Measurements
Key: <i>d</i> = Recommended separation distance in meters (m)		
<i>P</i> = maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer		
<i>E1</i> = Tested compliance level (in Volts/meter) for the Radiated RF Immunity test IEC 61000-4-3		
The device meets the compliance level of 3.0 V _{RMS} according to IEC 60601-1-2 over the specified test frequency range.		
Over the frequency range 80 kHz to 800 MHz, the recommended separation distance in meters (<i>d</i>) is found by the following equation:		
$d = \left(\frac{3,5}{E1}\right)\sqrt{P}$	For a compliance level of 3.0 V _{RMS} :	$d = 1, 2\sqrt{P}$
Over the frequency range 800 kHz to 2.5 GHz, the recommended separation distance in meters (<i>d</i>) is found by the following equation:		
$d = \left(\frac{7,0}{E1}\right)\sqrt{P}$	For a compliance level of 3.0 V _{RMS} :	$d = 2, 3\sqrt{P}$

Field strengths from fixed transmitters, such as base stations or radio (cellular, cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, it should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

If you require further information or assistance, please contact Philips Support.

Recommended Separation Distances from Other RF Equipment

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment and the device as recommended below, according to the maximum output power of the communications equipment.

Table 5 - Separation Distance (d) in Meters According to Frequency of Transmitter at IEC 60601-1-2 Test Compliance Level			
Rated Maximum Output Power (P) of Transmitter (in Watts)	150 kHz to 80 MHz $d = \left(\frac{3.5}{V1}\right) \sqrt{P}$	80 MHz to 800 MHz $d = \left(\frac{3.5}{E1}\right) \sqrt{P}$	800 MHz to 2.5 GHz $d = \left(\frac{7.0}{E1}\right) \sqrt{P}$
0.01 W	0.1 m	0.1 m	0.23 m
0.1 W	0.4 m	0.4 m	0.7 m
1 W	1.2 m	1.2 m	2.3 m
10 W	3.8 m	3.8 m	7.3 m
100 W	12.0 m	12.0 m	23.0 m

Radio Frequency Radiation Exposure Information

The radiated output power of the device is far below the FCC radio frequency exposure limits. Nevertheless, the device shall be used in such a manner that the potential for human contact during normal operation is minimized.

Environment

Before operation, make sure that the fetal monitor is free from condensation. This can form when equipment is moved from one building to another, and is exposed to moisture and differences in temperature.

Use the monitor in an environment which is reasonably free from vibration, dust, corrosive or explosive gasses, extremes of temperature, humidity, and so forth. It operates within specifications at ambient temperatures between 0 and +45°C (32°F to 113°F). Ambient temperatures that exceed these limits can affect the accuracy of the system, and can damage the components and circuits.

Ambient temperature ranges for storage are -20°C to +60°C (-4°F to 140°F) for the monitor, and -40°C to +60°C (-40°F to 140°F) for transducers.

The transducers are watertight to a depth of 1.0 m for at least five hours (rated IP 68).

WARNING

- **Leakage currents:** If several items of equipment used to monitor a patient are interconnected, the resulting leakage current may exceed allowable limits.
 - **ECG electrodes:** NEVER allow ECG electrodes to contact other electrical conductive parts, including earth.
-

Monitoring After a Loss of Power

If the monitor is without power for **less** than one minute, monitoring will resume with all active settings unchanged. If the monitor is without power for **more** than one minute, the behavior depends on your configuration. If **Automat. Default** is set to **Yes**, the User Defaults will be loaded when power is restored. If **Automat. Default** is set to **No**, all active settings are retained, if power is restored within 48 hours. The **Automat. Default** setting is made in Configuration mode.

FM20/30 with Battery Option, FM40/50 When power is lost - no power is available from the AC power source or from the battery - a beeper will sound. The tone can be silenced by pressing the On/Standby switch.

ESU, MRI and Defibrillation

WARNING

The fetal/maternal monitors are NOT intended for use during defibrillation, electro-surgery, or MRI. Remove all transducers, sensors, and accessories before performing electro-surgery, defibrillation, or MRI, otherwise harm to the patient or the user can result.

Cardiac Pacemakers and Electrical Stimulators














WARNING













The fetal/maternal monitors are not intended for use for ECG measurements on patients connected to external electrical stimulator or with cardiac pacemakers.











Fast Transients/Bursts

The equipment will return to the previous operating mode within 10 seconds without loss of any stored data.

Symbols on the System

	<p>This caution/warning symbol indicates that you should consult the Instructions for Use (this document), and particularly any warning messages.</p>		<p>Connection direction indicator FM20/FM30 with battery option</p>
	<p>This symbol indicates that you should consult the Instructions for Use (this document).</p>		<p>Mouse connection indicator (optional)</p>
	<p>Power-On/Off Switch - FM20/FM30 without Battery Option</p>		<p>Keyboard connection indicator (optional)</p>
	<p>Power-On/StandBy Switch - FM40/FM50 and FM20/30 with Battery Option</p>		<p>Serial/MIB connector (optional)</p>
	<p>Power-On LED</p>		<p>USB interface (optional)</p>
	<p>Electrical Class II equipment, in which the protection against electric shock relies on double or reinforced insulation (FM20/FM30).</p>	<p>Video </p>	<p>Analog interface indicator for connection to any analog video display (VGA resolution) FM40/FM50</p>
	<p>Fetal Sensor Socket symbol.</p>	<p>IP 67</p>	<p>Ingress Protection code according to IEC 60529. The IUP/ECG patient module (M2738A) is rated IP 67 (protection against dust, access to hazardous parts, and the effects of continuous immersion in water to a depth of 0.5 meter for 30 minutes).</p>

	SpO ₂ Socket symbol.	IP 68	Ingress Protection code according to IEC 60529. All transducers (excluding M2738A) are rated IP 68 (protection against dust, access to hazardous parts, and the effects of continuous immersion in water to a depth of 1.0 meter for five hours).
	Noninvasive Blood Pressure Socket symbol.	IP X1	Ingress Protection code according to IEC 60529. The monitors and interface cable for the Avalon CTS (M2731-60001/ M2732-60001) are rated IP X1 (protection against water dripping vertically only).
	Symbol indicating the monitor has the triplets option.		Type CF equipment, not defibrillation proof.
	Symbol indicating the monitor is capable of intrapartum monitoring.		Indicates location of service number
	Button to open paper drawer/paper eject. (FM40/FM50).		Indicates location of serial number
	Protective earth terminal (FM40/FM50).		Indicates location of catalog number
	Equipotential grounding point (FM40/FM50).		Indicates location of the date of manufacture and/or name and address of manufacturer

<p>Tele</p>	<p>Socket for connecting Avalon CTS interface cable M2732-60001 or Avalon CL interface cable (with black connector, FM40/FM50)</p>		<p>Identifies the year and month of manufacture.</p>
	<p>Symbol indication for non-ionizing radiation.</p>		<p>Symbol indicating separate collection for waste electrical and electronic equipment.</p>
	<p>China ROHS</p>		<p>CSA US mark</p>
	<p>GOST-R mark</p>		<p>R&TTE Compliance Association: "Class 2 Device" Mark</p>
<p>FCC ID</p>	<p>Federal Communications Commission: FCC ID xxxx</p>	<p>CMIIT ID</p>	<p>Chinese Radio marking: CMIIT ID (China Ministry of Industry and Information Technology)</p>
	<p>Japanese Radio marking: Radio mark + "PMS Japan" (in Japanese) + Prod. No. + [R]-symbol + ID "Taiwan Radio Label" (NCC Logo) + ID</p>		<p>Korea radio mark: KC logo, KCC ID number, and Conformity assessment information</p>
<p>ISM</p>	<p>Industrial, Scientific, & Medical radio frequency band (Avalon CL frequency band used e.g. in the EU)</p>	<p>WMTS</p>	<p>Wireless Medical Telemetry Service (Avalon CL frequency band used e.g. in North America)</p>
<p>T108</p>	<p>Association Of Radio Industries And Businesses T108 (Avalon CL frequency band used e.g. in Japan)</p>	<p>CE 0366</p>	<p>CE marking accompanied by the VDE NB number 0366</p>
<p>IC: 8888X-XXXXX8</p>	<p>IC-ID (Industry Canada ID) One IC-ID labeling for each built in radio: OBR, SRR</p>		<p>Taiwan Radio Label (NCC Logo) + ID</p>

Default Settings Appendix

This appendix documents the most important default settings of your fetal monitor and the Avalon CL Base Station with the cableless transducers as they are delivered from the factory. For a comprehensive list and explanation of default settings see the Configuration Guide supplied with your fetal monitor. The monitor's default settings can be permanently changed in Configuration Mode.

Alarm and Measurement Default Settings

Alarm Defaults Settings

Alarm Setting	Choice	Default
Alarm Mode	INOP only, All	INOP only
Alarm Volume	0..10	5
Alarms Off	1 min, 2 min, 3 min, Infinite	3 min
Alarm Text	Standard/Enhanced	Standard
Visual Latching	Red&Yellow/Red/Off	Off
Audible Latching	Red Only/Off	Off
Alarm Sounds	Traditional/ISO	Traditional
Alarm Low	0..10	4

Fetal / Maternal Defaults Settings		
FHR (Ultrasound/DECG)	Alarms On/Off Default	On
	Default Color for FHR Numeric	Orange
Toco	Default color for Toco numeric	Green
IUP	Default IUP Scale Unit	mmHg
	Default color for IUP numeric	Green
Maternal Heart Rate HR Measurement	Default Color for MECCG Numeric	Red

Coincidence Default Settings

Coincidence Default Settings	
Coincidence Tone	immediately

NBP Default Settings

Factory Default Settings	
Mode	Manual
Repeat Time	15 min
Alarms from	Sys.
Low Limit	90 / 50 (60)
High Limit	160 / 90 (110)
VP Pressure	mmHg
Done Tone	Off
Veni Puncture	60 mmHg
Start Time	Synchronized
NBP	On
Alarms	On
Color	Red
Reference	Auscultatory

CL NBP Default Settings

NBP Setting	Factory Defaults
Mode	Auto
Repeat Time	10 min
NBP	On
VP Pressure	60 mmHg
Reference	Auscultatory
Unit	mmHg
Done Tone	Off
Start Time	Synchronized
Aging Time	10 min
Color	Red

SpO₂ Default Settings

SpO ₂ Factory Default Settings	
Desat Limit	80
Low Limit	90
Low Limit	100
High Limit	20 seconds
Desat Delay	10 seconds
Low Alarm Delay	10 seconds
High Alarm Delay	10 seconds
NBP Alarm Suppr.	On
Alarms	On
Color	Cyan

Pulse Default Settings	
Pulse (SpO ₂)	On
High Limit	120 bpm
Low Limit	50 bpm
Bradycardia: Difference to Low Limit	20 bpm
Bradycardia: Clamp	40 bpm
Tachycardia: Difference to High Limit	20 bpm
Tachycardia: Clamp	200 bpm

CL SpO₂ Default Settings

SpO ₂ Settings	Factory Defaults
Repeat Time	15 min
Signal Quality	On
NBP Alarm Suppr.	On
Color	

Tympanic Temperature Default Settings

pTemp Settings	Factory defaults H10/H20/H40	H30 (deviations from H10/H20/H40)
	Adult	Adult
Low Limit	36°C (96.8°F)	
High Limit	39°C (102.2°F)	
Unit	°C	
Alarms	On	
Color	Green	

Manually Entered Values Default Settings

Manual Entered Values	Default Settings
Label	Temp
Unit	°C
Color	white
Interval	1h
Msmt	off

Recorder Default Settings

Setting	Choice	Default
Recorder Speed	1, 2, or 3 cm/min	3 cm/min
Scale Type	US, Internat'l	US
Trace Style FHR1	Thin, Medium, Thick, Extra Thick	Thick
Trace Style FHR2		Medium
Trace Style FHR3		Extra Thick
Trace Style Toco		Thick
Trace Style HR		Thin
Wave Style ECG		Thin
ECG Wave printing choice		Separate, Overlap
Notes Recording	Along, Across	Along
Auto Start	Off, On	Off
Confirmed Stop		Off
Bridge Paperout		On
Paper Save Mode		Off
NST Autostart		On
NST Autostop		Off
Trace Separation		Off
Separation Order		Standard, Classic

Setting	Choice	Default
Intensity	4 (medium)	n/a
Cal. Offset	5	n/a

Index

A

accessories 267, 268, 270, 271, 272, 273, 277
 cl NBP pod accessories 272
 cl SpO2 pod accessories 276
 fetal accessories 268
 fetal recorder accessories 277
 latex information 267
 MECG accessories 270
 noninvasive blood pressure accessories 270, 271, 272
 adult cuffs 272
 multi-patient comfort cuff kits 271
 reusable cuffs 271
 single-hose disposable cuffs 272
 SpO2 accessories 273
 Philips sensors (disposable) 273
 Philips sensors (reusable) 273
 transducer accessories 268

ACOG technical bulletin 80

actions after monitoring 98

active alarms 109

adjusting brightness 66

adjusting display 68

adjusting screen brightness 66

admitting a patient 137, 138
 editing information 138

aging numerics 213

alarms 109, 110, 111, 113, 114, 115, 116, 117, 118, 120, 121, 123, 191, 216, 230, 285, 289, 290, 311
 acknowledging alarms 114
 active alarms 109
 alarm behavior at on/off 121
 alarm latching 119
 alarm limit delay 285
 alarm reminder 114
 alarm standards 111
 audible alarm indicators 111
 yellow alarms 109

artifact suppression 208
 artifact suppression and fetal arrhythmia 208
 explanation of artifact suppression 208
 switching on and off artifact suppression 208

automatic noninvasive blood pressure repeat time 215

B

baseline measurement 177, 184, 194

Battery Option 101
 using batteries 101

belt actions 70, 72

belt fastening 70

belt fixing button 70

C

Cableless Monitoring 91
 Avalon CL 40
 configuration cableless system 92
 important considerations 159
 LED indication 53
 TELE window 93
 Telemetry 94
 Underwater monitoring 97

calibration 216, 260
 calibrating noninvasive blood pressure 260

cautions 9

changing FHR alarm delays 191

changing FHR alarm limits 191

cleaning monitor actions 255, 256, 257, 265
 cleaning method 256
 cleaning monitoring accessories 257
 general cleaning requirements 255
 infection control 255
 printhead cleaning 265

configuration mode 63

configuring alarm tone 111

confirming fetal life 10, 158, 175, 183

connecting power 12

cross-channel verification 10, 175, 183, 234
 CCV and triplets 183
 CCV and twins 175

D

data recovery 251

date from OB TraceVue 66

date setting 66

DECG 205, 209, 285
 DECG INOPs 209
 DECG specifications 285
 DECG testing 209
 DECG troubleshooting 209

default alarm settings 311

default IUP settings 285

default settings 285
 FHR default settings 285

IUP default settings 285

MECG default settings 285

Toco default settings 285
 user defaults 65

defibrillation precautions 238, 267, 307

demo mode 63

disabling touchscreen 63

discharging a patient 138

disinfecting 255, 256, 257
 infection control 255
 recommended substances 257

display 68, 213, 285
 adjusting display 68
 display noninvasive blood pressure 213
 fetal display specifications 285

disposal of electronic waste 266

disposal of monitor 266

dyshemoglobins 228
 intravascular dyshemoglobins 228

E

early systolic blood pressure 213

ECG 205, 235, 285
 ECG specifications 285

electrical safety tests 298

electrical surgery precautions. See ESU 228, 238, 267, 307

electrodes 237
 applying electrodes 237

EMC 298
 electromagnetic compatibility 298
 electromagnetic emissions 301
 electromagnetic interference 300
 how to reduce electromagnetic interference 300
 EMC and compliant accessories 298
 EMC precautions 298
 EMC standards 298

entering notes 61

ESU precautions 238, 267, 307

event marker remote 80

extension cable for SpO2 228

external monitoring 157, 176, 184, 194
 external monitoring - Toco 194
 external monitoring - triplets 184
 external monitoring - twins 176
 external monitoring - ultrasound 157

F

- factory defaults 289, 290, 294
 - noninvasive blood pressure factory defaults 289
 - recorder factory defaults 294
 - SpO2 factory defaults 290
- FAST 227
 - Fourier Artifact Suppression Technology 227
- fetal accessories 268
- fetal arrhythmia 208
- fetal demise 10
- fetal display specifications 285
 - IUP - fetal display specifications 285
 - Toco - fetal display specifications 285
 - US - fetal display specifications 285
- fetal heart rate alerting. See FHR 191
- fetal heart rate. See FHR 157
- fetal movement 163
- fetal movement profile. See FMP 157
- FHR 10, 158, 165, 175, 177, 183, 184, 191, 192, 234, 285
 - alerting 191
 - cross-channel verification 234
 - misidentification 10
 - trace separation 177, 184
- FHR alarms 285
- finder LED 45, 73, 175, 183
- finding monitor revision 67
- flexible nurse call 30, 42, 44
- FMP 163, 164
 - FMP and twins 163
 - FMP statistics 163
- functional arterial oxygen saturation 227

G

- global settings 65

H

- high priority alarms 109

I

- infection control 255
 - infection and sterilizing 255
- INOPs 109, 114
 - CCV INOP 33
 - coincidence tone 30
 - disconnect INOPs 114
 - INOP alarms 109
 - INOPs indicators 109
 - silencing INOPs 114
- interference 298, 300
 - how to reduce interference 300
 - interference from RF equipment 298

- internal monitoring 177, 199, 201
- intrauterine pressure. See IUP 197
- intravascular dyshemoglobins 228
 - intravascular dyshemoglobins effects on (SpO2) 228
- ISO/IEC standard 9703-2 111
- IUP 197, 199, 285
- IUP specifications 285

L

- latching alarms 118, 119
 - latching alarm behavior 119

M

- mains power (AC) 12, 99
 - connecting to mains power 12
 - disconnecting from mains power 99
- maintenance 259, 260
 - maintenance schedule 259, 260
 - visual inspection 259
- maternal ECG 208, 238, 240
 - electrode position 238
 - printing waveform 208, 240
 - viewing waveforms 238
- maternal heart rate 233, 234, 235
 - maternal heart rate and cross-channel verification 234
 - maternal heart rate and priorities 233
 - maternal heart rate from MCEG 235
 - sources of maternal heart rate 233
- maternal puls from Toco 32
- Maternal Temperature 30
- measurement settings 65
 - changing measurement settings 65
 - entering setup menu 65
- measurements 64, 65, 285
 - setting up measurements 65
 - switching on and off measurements 64
- MCEG 208, 235, 237, 238, 240, 242, 243, 270, 285
 - accessories 270
 - alarm limits 285
 - connection illustration 235
 - default settings 285
 - electrodes 237
 - specifications 285
 - waveform printing 208, 240
- MCEG alarm limits 285
- methemoglobin (SpO2) 228
- MHR 158, 233, 235, 237
 - MHR from MCEG 235
 - MHR priorities 233
 - MHR sources 233
- monitor settings 10, 40, 66, 68, 98, 251
 - changing monitor settings 66

- major parts and keys 40
- power on/off behavior 99
- Standby 98
- switching on 68

- monitoring mode 63
- monitoring triplets and FHR 183
- monitoring twins and FHR 175
- MR imaging and the SpO2 transducer 228
- MRI precautions 238, 267, 307

N

- NBP. See noninvasive blood pressure 211
- NIBP. See noninvasive blood pressure 211
- Noninvasive Blood Pressure (NBP) 211, 212, 213, 214, 215, 216, 217, 242, 260, 270, 271, 272, 289
 - calibration interval 260
 - comfort cuff kits 271
 - factory defaults 289
 - how the measurement works 211
 - measurement 212, 214
 - limitations 212
 - methods 212
 - starting 214
 - stopping 214
 - NBP accessories 270
 - NBP adult cuffs 272
 - NBP automatic mode 215
 - NBP numerics 213
 - NBP pediatric cuffs 272
 - NBP performance specifications 289
 - NBP reusable cuffs 271
 - NBP single-hose disposable cuffs 272
 - NBP site inspection 212
 - NBP troubleshooting 217
 - NBP venous puncture 216
 - preparing to measure NBP 212
- non-medical devices 298
 - non-medical devices in patient vicinity 298
- Non-Stress Test timer. See NST timer 143
- NST Report
 - NST analysis 32
 - Trace Interpretation 145, 148
- NST timer 141, 142, 143

O

- OB TraceVue 38, 39, 252
 - connection to OB TraceVue 38, 39
 - LAN connection 252
- operating modes 63
 - configuration mode 63
 - demo mode 63
 - monitoring mode 63
 - service mode 63
- operating temperatures 306

oscillometric measurement method 211

P

Paging Patient 97

paper 80, 81, 82, 83, 85, 264, 277

default paper speed 80

loading paper 83, 85

paper out INOP 82

paper speed 80

changing 80

default 80

defaults 80

setting 80

paper, when to reload 82

recorder paper 277

removing paper 83, 85

storing paper 264

tearing off a trace 81

Paper Advance key 61

Paper-Out indication 82

passcode protection 63

patient alarms 123

Patient Demographics window 137

patient module 205, 235

paused alarms 115, 116

restarting paused alarms 116

pausing alarms 115

performance specifications 285, 289, 290

noninvasive blood pressure

specifications 289

SpO2 specifications 290

physiological alarms 109

pop-up keys 95

power 12, 99

connecting power 12

disconnecting from 99

power connection 12

power failure 99

pulse 242

pulse from noninvasive blood

pressure 242

pulse rate 233

R

recorder 76, 100, 264, 265, 277, 294

factory defaults 294

paper storage 264

printhead cleaning 265

recorder paper 277

recorder specifications 294

recorder speed 33

recorder troubleshooting 100

Recorder Elements 79

recording stored data 61

recovering traces 251

trace recovery printout 251

recycling 266

red alarms 109

regulatory compliance 297

remote event marker 80

restarting paused alarms 116

reviewing alarms 118

reviewing alarms window 118

S

safety 260, 297

maintenance interval 260

safety standards 297

screen layouts 64

selftest 68, 76, 120

sensor 227

disposable SpO2 sensors 227

selecting SpO2 sensors 227

separating FHR traces 177, 184

separating triplet traces 184

separating twin traces 177

service mode 63

settings 64, 65, 66

active settings 64

changing settings 66

global settings 65

monitor settings 66

user default settings 65

signal loss delay 192

signal quality 90, 230

signal quality during monitoring 90

signal quality of SpO2 230

silencing alarms 114

specifications 279, 281, 285, 294

built-in recorder specifications 294

DECG specifications 285

ECG specifications 285

environmental specifications 279

monitor 279

SpO2 sensors 279

transducers 279

IUP specifications 285

measurement specifications 285

MECG specifications 285

performance specifications 285

Toco specifications 285

transducer specifications 281

US specifications 285

SpO2 227, 228, 230, 273, 290

accessories 273

alarm specifications 290

connecting the cables 228

disposable sensors 227

factory defaults 290

FAST technology 227

performance specifications 290

Philips sensors (disposable) 273

Philips sensors (reusable) 273

signal quality 230

SpO2 alarms 230, 290

standards 297, 298

EMC standards 298

safety standards 297

standards compliance 297

Standby 61

sterilizing 255, 258

infection control 255

storage temperatures 306

stored data recording 61

suppressing artifacts 208

suspended alarm 116

suspicious SpO2 reading 230

switching on and off 164

switching on and off alarms 116

switching on monitor 68

switching on recorder 76

symbols on the system 308

T

technical alarm messages. See INOPs 126

temperatures 306

operating temperatures 306

storage temperatures 306

testing 120, 209, 242, 298

DECG testing 209

MECG testing 242

safety testing 298

testing alarms 120

time 66

setting time 66

time from OB TraceVue 66

Toco 193, 194, 195, 205, 237, 285

default Toco settings 285

testing a Toco transducer 195

Toco baseline 194

Toco display specifications 285

Toco monitoring 193

Toco sensitivity 195

Toco specifications 285

Toco MP transducer 193, 233

Toco+ transducer 193

touch tone volume 66

touchscreen operation 56, 63

trace actions 81, 251, 252

separation on/off 33

tearing off traces 81

trace recovery on OB TraceVue 252

trace recovery on paper 251

transducer 45, 70, 73, 173, 175, 183, 193, 195,

- 281
- connecting transducer to monitor 73
- fixing transducer to belt 70
- repositioning transducers 73
- Toco transducer 193
- Toco+ transducer 193
- transducer belt clip 70
- transducer specifications 281
- transducer testing 173, 195
 - Toco 195
 - ultrasound 173
- triplets 183, 184, 185, 188
 - cross-channel verification 183
 - external monitoring 184
 - importance of monitoring MHR 183
 - monitoring triplets 183
 - offsetting baselines 184
- troubleshooting 100, 165, 181, 188, 195, 199, 209, 217, 240
 - DECG 209
 - FHR 165
 - IUP 199
 - MECG 240
 - troubleshooting FHR 165
- twins 163, 175, 176, 177, 178, 179, 181
 - cross-channel verification 175
 - external monitoring 176
 - importance of monitoring MHR 175
 - invasive monitoring 177
 - monitoring FHR 175
 - monitoring twins 175
 - offsetting baselines 177
 - twins and FMP 163

U

- ultrasound 157, 160, 173, 285
 - applying gel 160
 - testing a transducer 173
 - ultrasound display specifications 285
 - ultrasound monitoring 157
 - ultrasound specifications 285
- understanding screens 64
- USB 30
- using the adapter cable 205
- uterine activity 193, 194, 195, 197, 199, 205, 237
 - external monitoring 193
 - internal monitoring 197
 - monitoring uterine activity 194, 197, 199, 205, 237
 - troubleshooting uterine activity 195, 199

V

- velcro belt 72
- venous puncture 216
- visual alarm indicators 110

- voltage setting 12
- volume 66, 113
 - alarm volume 113
 - IUP, zeroing 199
 - touch tone volume 66
- warnings 9

Part Number xxxxxxxxxxxx
Printed in Germany 06/13



PHILIPS