

OMRON®

Wheeze Detector

WheezeScan

(HWZ-1000T)

Instruction Manual

Rx ONLY



All for Healthcare

Table of Contents

Introduction	2
Safety Instructions.....	2
Intended Use	2
Receiving and Inspection	2
Important Safety Information.....	3
1. Know Your Device	7
1.1 Contents.....	7
1.2 Part Names	8
2. Installing Batteries	9
3. Detecting Wheeze	10
3.1 Taking a Measurement.....	10
3.2 Checking Result.....	14
4. Managing Results with Your Smart Device	15
5. Maintenance	16
5.1 Maintenance	16
5.2 Storage.....	16
5.3 Cleaning	16
5.4 Disposal	16
6. Troubleshooting	17
7. Specifications	19
8. FCC Statement	21
9. Limited Warranty	22
10. Guidance and Manufacturer’s Declaration	23
11. How the Device Works	25

Introduction

Thank you for purchasing the OMRON WheezeScan. WheezeScan uses innovative listening technology to identify "wheeze" or "no wheeze" quickly and accurately. WheezeScan is a convenient tool to help parents to manage their children's asthmatic symptoms at home. The cause of these symptoms can be bronchitis, bronchiolitis or bronchial asthma.

Safety Instructions

This instruction manual provides you with important information about WheezeScan. To ensure safe and proper use of this device, READ and UNDERSTAND all of these instructions. **If you do not understand these instructions or have any questions, contact 1-800- 634-4350 before attempting to use this device.**

Intended Use

The OMRON WheezeScan is intended to detect the presence of wheezing in infants, toddlers and children from 4 months to 7 years of age.

The WheezeScan shows "Wheeze" or "No Wheeze" as the result of measurement. The obtained results should be discussed with your physician.

Environment of use: home, clinic

Receiving and Inspection

Remove this device and other items from the packaging and inspect for damage. If this device or any other items are damaged, DO NOT USE and contact 1-800-634-4350.

Important Safety Information

Read the Important Safety Information in this instruction manual before using this device.

Follow this instruction manual thoroughly for your safety. Keep for future reference. For specific information about patient's symptoms, **CONSULT WITH A PHYSICIAN.**

WARNING

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

- NEVER diagnose or treat a patient based only on the result from this device. ONLY a physician is qualified to diagnose and treat patient's respiratory disease.
- ONLY apply this device on a patient whose illness was diagnosed as bronchitis, bronchiolitis or bronchial asthma.
- Consult with a physician before using this device. This device is intended to be used under the direction of a physician or licensed healthcare professional.
- If a patient is having signs and symptoms of chest tightness, shortness of breath, coughing or wheezing, a wheeze may not be detected because:
 - the breathing sound is too loud or too soft;
 - the wheeze is at extremely small amplitude, extremely short duration, extremely high frequency or extremely low frequency;
 - an asthma attack occurs without wheezing.

Always take action as prescribed by a physician, no matter what the result from this device is. This device also may incorrectly detect any noise as a wheeze because various sounds are occurring in patient's body such as:

- voice or cry of a patient;
- nasal obstruction sound;
- sputum sound etc.
- DO NOT use in case of severe asthma attack or respiratory failure. Consult with a physician immediately.
- Take a measurement again if you are concerned about the result. This device can detect wheeze only while taking a measurement. Wheeze can occur after taking a measurement. If you are still concerned after taking a measurement again, consult with a physician.
- DO NOT use this device on a chest that is injured or that is under medical treatment.
- DO NOT disassemble or attempt to repair this device. This may cause an inaccurate result. This product contains small parts that may cause a choking hazard if swallowed by infants, toddlers or children.

Data Transmission

- This product emits radio frequencies (RF) in the 2.4 GHz band. DO NOT use this product in locations where RF is restricted, such as on an aircraft or in hospitals. Turn off the **Bluetooth®** feature in this device and remove batteries when in RF restricted areas. For further information on potential restrictions refer to documentation on the Bluetooth usage by the FCC.

Battery Handling and Usage

- Keep batteries out of the reach of infants, toddlers or children.

CAUTION

Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient, or cause damage to the equipment or other property.

- ONLY apply this device on a patient who is between 4 months old and 7 years old.
- ONLY parents who are advised by a physician, caregivers, school nurses and physicians are allowed to use this device. DO NOT let patient use this device on his/her own.
- An inaccurate result may happen when:
 - clothes touch the device (ALWAYS place the sound sensor on the patient's skin directly);
 - the device is used on crying, speaking or moving patients;
 - the sound sensor is placed on an incorrect place on the body;
 - the sound sensor is placed on the body too strongly or too lightly;
 - the sound sensor is not adhered to the body;
 - the patient's chest has too much muscle or fat.
- Stop using this device and consult with a physician if you or the patient experience skin irritation or discomfort.
- DO NOT use this device for any purpose other than detecting wheeze.
- DO NOT use in a location where there is moisture or a risk of water splashing on this device. This may damage this device.
- DO NOT wash or immerse the device in water.
- DO NOT use chemicals other than alcohol to clean the sound sensor.
- DO NOT drop or subject this device to strong shocks or vibrations.
- DO NOT use this device in places with high or low humidity or high or low temperatures. Refer to "Operating conditions" in section 7.
- This device may not work when the sound sensor is damaged. If it happens, contact customer service.

Important Safety Information

- This device does not have a waterproof structure. Contact customer support immediately if the device is exposed to or submerged in water.
- Contact customer support if the device receives a strong shock such as being dropped.
- During operation, make sure that no mobile device or any other electrical device that emits electromagnetic fields is within 30 cm of this device. This may result in incorrect operation of the device and/or cause an inaccurate indication.
- DO NOT use this device with other medical electrical (ME) equipment simultaneously. This may result in incorrect operation of the device and/or cause an inaccurate indication.
- Ensure that this device has acclimated to room temperature before taking a measurement. Taking a measurement after an extreme temperature change could lead to an inaccurate indication. OMRON recommends waiting for approximately 2 hours for the monitor to warm up or cool down when the device is used in an environment within the temperature specified as operating conditions after it is stored either at the maximum or at the minimum storage temperature. For additional information on operating and storage/transport temperatures, refer to section 7.
- DO NOT keep touching the sound sensor for more than 1 minute.

Unusable Environment

- This device cannot be used in an environment with noise at 60 dB or above, such as:
 - near working vacuum cleaners, hair dryers or washing machines
 - in moving vehicles
 - in noisy public or private areasIn any of these places, an error may occur while taking a measurement.

Data Transmission

- DO NOT replace batteries while readings are being transferred to your smart device. This may result in incorrect operation of this device and failure to transfer readings.

Battery Handling and Usage

- DO NOT insert batteries with their polarities incorrectly aligned.
- ONLY use 2 "AAA" alkaline batteries with this device. DO NOT use other types of batteries. DO NOT use new and used batteries together.
- Remove batteries if this device will not be used for a long period of time.

Important Safety Information

- If battery fluid should get in patient's or your eyes, immediately rinse with plenty of clean water. Consult with a physician immediately.
- If battery fluid should get on patient's or your skin, wash the skin immediately with plenty of clean, lukewarm water. If irritation, injury or pain persists, consult with a physician.
- DO NOT use batteries after their expiration date.
- Periodically check batteries to ensure they are in good working condition.
- ONLY use batteries specified for this device. Use of unsupported batteries may damage and/or may be hazardous to this device.

1. Know Your Device

1.1 Contents

Device



Storage case




2 "AAA" alkaline
batteries



Instruction manual


Quick reference guide

1.2 Part Names



A	Sound sensor	H	Connection indicator
B	Display	I	[Connection] button
C	Error symbol		
D	Battery symbol Low: flashes / Depleted: lights		
E	Power indicator		
F	[Power] button		
G	Battery compartment		

2. Installing Batteries


When the “” symbol flashes, it is recommended that you replace the batteries. To replace batteries, turn the device off and remove batteries. Then replace with 2 new alkaline batteries.

Note

- Replacing batteries does not delete previous results recorded in the device.
- The supplied batteries may have a shorter life span than new batteries.
- Disposal of used batteries should be carried out in accordance with local regulations.


1 Open the battery cover.

Turn the battery cover anticlockwise to open. DO NOT pull the battery cover strongly. The rubber may be torn off.




2 Insert 2 “AAA” alkaline batteries as indicated in the battery compartment.

Confirm the polarities correctly aligned.



3 Close the battery cover.

Turn the battery cover clockwise to close.



3. Detecting Wheeze

3.1 Taking a Measurement

Do not use the device in locations exposed to extreme temperatures, humidity, direct sunlight, dust or corrosive vapors such as bleach.

CAUTION

This device cannot be used in an environment with noise at 60 dB or above, such as:


- near working vacuum cleaners, hair dryers or washing machines;
- in moving vehicles;
- in noisy public or private areas

In any of these places, an error may occur while taking a measurement.

- 1 Remove phlegm and nasal mucus, if any, before taking a measurement.**



- 2 Let the patient sit down or lay face up.**



3.1 Taking a Measurement

3 Remove clothes to expose patient's chest. Being naked from the waist up is recommended.

⚠ WARNING

DO NOT place the sound sensor on clothes and DO NOT let clothes touch the device. This may cause an inaccurate result. ALWAYS place the sound sensor directly on the patient's skin.

4 Check where to place the sound sensor.


The sound sensor should be placed on the **right chest** under the right collarbone.

For appropriate location on the body, see the illustration shown below.

To avoid damage to the sound sensor, wipe off baby oil or any chemical products from the patient's skin before taking a measurement.

⚠ WARNING

DO NOT use this device on a chest that is injured or that is under medical treatment.






3.1 Taking a Measurement

5 Press the button.

The device will beep and the power indicator will light green.

Note


- When the power is turned on, all indicators will light up for 2 seconds. Then, the previous measurement result (refer to section 3.2) will be displayed for 2 seconds.
- If you do not start a measurement within 30 seconds after turning the power on, “Err” symbol will appear. To restart the device, press and hold the  button to turn the power off and press it again to turn the power on.
- Ensure that both the parent and the patient DO NOT talk while the measurement is being taken.



3.1 Taking a Measurement

6 Apply the sound sensor to the patient's skin.

Place and adhere the sound sensor firmly.



7 The device will beep twice and begin measuring.

When the sound sensor is correctly placed on the body, the device gives twice the sound "beep" and the measurement starts automatically.

After the sound sensor is applied to the skin, the sensor will flash white while taking the measurement.



Note

- Ensure that the patient DOES NOT talk, cry or move while the measurement is being taken.
- Ensure that the patient continues to breathe normally.

8 Wait until the measurement is completed.

When the measurement is completed, the device beeps 3 times. The maximum measurement time is 30 seconds.

Note

KEEP the sensor in close contact and DO NOT move the device until the device beeps 3 times. Otherwise the result may be inaccurate.

3.2 Checking Result

After the measurement, the result will be shown as pictured below.

If wheezing is not detected, “No Wheeze” will light white.

If wheezing is detected, “Wheeze” will light orange.





Wheezing is NOT detected



Wheezing is detected

Note

- The device automatically turns off 3 minutes after completing a measurement. You can turn the power off manually by pressing and holding the  button.
- If you want to check the last result, press the  button to turn the power on. You will see the last result displayed.
- If the device can not detect the sound correctly, it gives 2 longer beeps, “Err” symbol appears and the sound sensor flashes or lights orange. Refer to section 6 for more information.

4. Managing Results with Your Smart Device

This device can pair with your smart device and you can manage the measurement results with your smart device using the “OMRON Asthma Diary” app. Review the list of compatible smart devices at:


www.omron-healthcare.com/asthma-diary-compatibility

Follow the steps below to pair the device with your smart device.


Note

If you do not pair the device, measurements you take will not be stored in the WheezeScan's memory.

1 Enable Bluetooth on your smart device.



2 Download and install the “OMRON Asthma Diary” app onto your smart device.



3 Open the app on your smart device and follow the instructions.

For more information about the app, refer to “Help” section of the app.

5. Maintenance

5.1 Maintenance

To protect the device from damage, please follow the directions below:

- Store the device in a clean, safe location.
- Changes or modifications not approved by the manufacturer will void the user warranty.

CAUTION

- DO NOT disassemble or attempt to repair this device. This may cause an inaccurate result. This product contains small parts that may cause a choking hazard if swallowed by infants, toddlers or children.

5.2 Storage

Put the device in the storage case, and store in a clean, safe location.

Do not store the device:

- If the device is wet.
- In locations exposed to extreme temperatures, humidity, direct sunlight, dust or corrosive vapors such as bleach.
- In locations exposed to vibrations or shocks.

5.3 Cleaning

- Use a soft cloth moistened with sanitizing alcohol to clean the sound sensor.
- Use a soft dry cloth or a soft cloth moistened with mild (neutral) detergent to clean the device other than the sound sensor, and then wipe with a dry cloth.
- DO NOT use any abrasive or volatile cleaners.
- DO NOT use gasoline, thinners or similar solvents to clean the device.

CAUTION


- DO NOT wash or immerse the device in water.

5.4 Disposal




Dispose of the device, battery and accessories in accordance with local regulations. Unlawful disposal may cause environmental pollution.

6. Troubleshooting

If any of the below problems occur, check to make sure that no other electrical device is within 30 cm. If the problem persists, refer to the table below.

Problem	Possible cause	Solution
<p>The sound sensor flashes orange and “Err” symbol appears.</p>  <p>The device is a white, cylindrical sound sensor. It has a circular top with a mesh-like surface. Below the top, there is a small rectangular display area showing the text 'WHEEZE', 'NO WHEEZE', and 'Err'. The device is shown from a front-facing perspective.</p>	<p>The measurement did not start within 30 seconds after the power was turned on.</p>	<p>Place the sound sensor on the patient's body within 30 seconds after the power is turned on.</p>
	<p>The patient, parent or anyone near the patient was speaking.</p>	<ul style="list-style-type: none">• DO NOT speak while the measurement.• DO NOT use this device on a crying, speaking or moving patient.
	<p>The device was used on clothes during the measurement.</p>	<p>DO NOT put this device on clothes. ALWAYS place this device on the patient's skin directly.</p>
	<p>The sound sensor was applied on an incorrect location on the body.</p>	<p>Apply the sensor on appropriate location (See sub-section 3.1 for measurement location).</p>
	<p>The sound sensor was not placed and adhered firmly on the skin.</p>	<p>KEEP the sensor in close contact with the patient and DO NOT move the device until you hear 3 beep sounds.</p>

6. Troubleshooting

Problem	Possible cause	Solution
<p>The sound sensor lights orange and “Err” symbol appears.</p> 	<p>There is an external noise around the device.</p>	<ul style="list-style-type: none"> • Take a measurement in a quiet place. • DO NOT talk while taking a measurement. • DO NOT apply this device on a crying patient.
<p>The sound sensor does NOT flash or light orange but “Err” symbol appears.</p>	<p>An internal error occurred.</p>	<p>Contact customer service. Refer to section 9.</p>
 flashes	<p>Batteries are low.</p>	<p>Replacing all batteries with 2 new batteries is recommended. Refer to section 2.</p>
 lights	<p>Batteries are depleted or battery polarities are not properly aligned.</p>	<ul style="list-style-type: none"> • Immediately replace all batteries with 2 new batteries. Refer to section 2.
<p>No power. Nothing appears.</p>		<ul style="list-style-type: none"> • Confirm the batteries are properly inserted with polarities correctly aligned. Refer to section 2.
<p>Results change when taking measurements repeatedly.</p>	<p>The patient’s condition is always changing.</p>	<p>NEVER diagnose or treat yourself based on the result of measurement. Consult with a physician if you are concerned.</p>

If you are not able to solve your problem, contact customer service. Refer to section 9.

7. Specifications

Product name	WheezeScan
Model	XXXXXXX ^[REF] HWZ-1000T-Z
Patient age range	From 4 months to 7 years old
Power source	2 "AAA" alkaline batteries 3 V (2 x 1.5 V)
Battery life	Approx. 450 measurements (using new alkaline batteries)
Durable period (service life)	5 years
Transmission method	Bluetooth ® Low Energy
Wireless communication	Frequency range: 2.4 GHz (2400 - 2483.5 MHz) Modulation: GFSK Effective radiated power: < 20 dBm
IP classification	IP 22
Operating conditions	+41 °F to +104 °F (+5 °C to +40 °C) / 15 to 90 % RH (non-condensing) / 700 to 1060 hPa
Storage / transportation conditions	-4 °F to +140 °F (-20 °C to +60 °C) / 10 to 90 % RH (non-condensing)
Weight	Approx. 2 oz. (55 g) (not including batteries)
Dimension	Approx. 1.3" (w) x 1.1" (d) x 4.7" (l) (33 x 27 x 121 mm)
Applied part Protection against electric shock	Type BF (sound sensor) Internally powered ME equipment
Operation Mode	Continuous operation
Memory	Stores up to 150 readings
Maximum temperature of the applied part	Lower than +118.4 °F (+48 °C)










7. Specifications

- These specifications are subject to change without notice.
- IP classification is degrees of protection provided by enclosures in accordance with IEC 60529. This device is protected against solid foreign objects of 12.5 mm diameter and greater such as a finger, and against oblique falling water drops which may cause issues during a normal operation.
- Operation mode is classification in accordance with IEC 60601-1.

About a wireless communication interference

This product operates in an unlicensed ISM band at 2.4 GHz. In the event this product is used near other wireless devices such as microwave and wireless LAN, which operate on the same frequency band as this product, there is a possibility that interference may occur. If interference occurs, stop the operation of the other devices or relocate this product away from other wireless devices before attempting to use it.

Symbols Description

	Applied part - Type BF Degree of protection against electric shock (leakage current)	IP XX	Ingress protection degree provided by IEC 60529
	Manufacturer		Temperature limitation
	Humidity limitation		Atmospheric pressure limitation
	Necessity for the user to consult this instruction manual		Power on / off
	The direction of correct polarity of the battery	Rx ONLY	The product is a medical device and Federal Law (USA) restricts this device to sale by or on the order of a physician
	To indicate generally elevated, potentially hazardous, levels of non-ionizing radiation, or to indicate equipment or systems e.g. in the medical electrical area that include RF transmitters or that intentionally apply RF electromagnetic energy for diagnosis or treatment.		

8. FCC Statement

FCC CAUTION

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with Part 15 of FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Note:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This transmitter must not be co-located or operated in conjunction with any other antenna or transmitter.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment and meets the FCC radio frequency (RF) Exposure Guidelines. This equipment has very low levels of RF energy that are deemed to comply without testing of specific absorption ratio (SAR).

9. Limited Warranty

Your WheezeScan, excluding batteries, is warranted to be free from defects in materials and workmanship appearing within 5 years from the date of purchase, when used in accordance with the instructions provided with the device.

We will, at our option, replace without charge any device covered by the above warranty. Replacement is our only responsibility and your only remedy under the above warranty.

To obtain warranty service contact Customer Service by calling 1-800-634-4350 for the address of the inspection center and the return shipping and handling fee.

Enclose the original printed receipt. Include a letter, with your name, address, phone number, and description of the specific problem. Pack the product carefully to prevent damage in transit. Because of possible loss in transit, we recommend insuring the product with return receipt requested.

THE FOREGOING IS THE SOLE WARRANTY PROVIDED BY OMRON IN CONNECTION WITH THIS PRODUCT, AND OMRON HEREBY DISCLAIMS ANY OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. IMPLIED WARRANTIES AND OTHER TERMS THAT MAY BE IMPOSED BY LAW, IF ANY, ARE LIMITED IN DURATION TO THE PERIOD OF THE ABOVE EXPRESS WARRANTY.

OMRON SHALL NOT BE LIABLE FOR LOSS OF USE OR ANY OTHER SPECIAL, INCIDENTAL, CONSEQUENTIAL OR INDIRECT COSTS, EXPENSES OR DAMAGES.

This warranty provides you with specific legal rights, and you may have other rights that vary by jurisdiction. Because of special local requirements, some of the above limitations and exclusions may not apply to you.

FOR CUSTOMER SERVICE

Visit our website at:
Call toll free:

OmronHealthcare.com
1-800-634-4350

10. Guidance and Manufacturer's Declaration

OMRON WheezeScan

Information for Accompanying Documents in the Scope of IEC60601-1-2:2014

Important information regarding Electromagnetic Compatibility (EMC)

This device manufactured by OMRON HEALTHCARE Co., Ltd. conforms to IEC60601-1-2:2014 Electromagnetic Compatibility (EMC) standard. Nevertheless, special precautions need to be observed:

- The use of accessories other than those specified or provided by OMRON could result in increased electromagnetic emission or decreased electromagnetic immunity of the device and result in improper operation.
- During operation, the use of the device adjacent to or stacked with other device should be avoided because it could result in improper operation. In case such use is necessary, the device and other device should be observed to verify that they are operating normally.
- During operation, portable RF communications device (including peripherals such as antenna cables and external antennas) should be used no closer than 12 inches (30 cm) to any part of the device. Otherwise, degradation of the performance of the device could result.
- Refer to further guidance below regarding the EMC environment in which the device should be used.

Table 1 - EMISSION Limits and Compliance

Phenomenon	EMISSION Limits	Compliance
Conducted and radiated RF EMISSIONS	CISPR 11	Group1, Class B

Table 2 - IMMUNITY TEST LEVELS

Phenomenon	Basic EMC standard	IMMUNITY TEST LEVELS
Electrostatic discharge	IEC 61000-4-2	±8 kV contact ±2 kV,±4 kV,±8 kV,±15 kV air for enclosure port and patient coupling port
Radiated RF electromagnetic fields	IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz for enclosure port
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See table 3
Conducted disturbances induced by RF fields	IEC 61000-4-6	6 Vrms 150 kHz to 80 MHz 80 % AM at 1 kHz for patient coupling port
Rated power frequency magnetic fields	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz for enclosure port

10. Guidance and Manufacturer's Declaration

Table 3 - Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications device

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 to 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 to 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704 to 787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745	800 to 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1720	1700 to 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1845	2400 to 2570	Bluetooth, WLAN, 802.11 b/g/n , RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
1970						
5240	5100 to 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5500	5800	a/n	Pulse modulation 217 Hz	0.2	0.3	9
5785						

11. How the Device Works

What is asthma?

Asthma is a condition in which the airways narrow and swell and produce extra mucus. This can make breathing difficult and trigger coughing, wheezing and shortness of breath.

Asthma can't be cured, but its symptoms can be controlled. Because asthma often changes over time, it's important that you work with your doctor to track the signs and symptoms and adjust the treatment as needed.

Asthma symptoms

Asthma symptoms vary from person to person. A person may have infrequent asthma attacks, have symptoms only at certain times — when exercising, or have symptoms all the time.

Asthma signs and symptoms include:

- Coughing
- Wheezing
- Runny nose or other cold symptoms
- Awakening due to coughing or difficulty breathing
- Being less active than usual

What is wheeze?

Wheezing is a high-pitched whistling sound made while breathing. It's often associated with difficulty breathing. Wheezing may occur during breathing out (expiration) or breathing in (inhalation).

How does WheezeScan work?

WheezeScan acts much like the stethoscope the doctor uses when listening to the lungs. By placing the sensor on the patient's chest for 30 seconds of normal breathing, the device can catch and then analyze the breathing sounds for the presence of wheezing.

The **Bluetooth**® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. and any use of such marks by OMRON HEALTHCARE Co., Ltd. is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries and regions. App Store is service mark of Apple Inc.

Google Play and the Google Play logo are trademarks of Google LLC.

Other trademarks and trade names are those of their respective owners.



OMRON HEALTHCARE Co., Ltd.

53, Kunotsubo, Terado-cho, Muko, Kyoto, 617-0002 JAPAN

Distributed by:

OMRON HEALTHCARE, INC.

1925 West Field Court, Lake Forest, IL 60045 U.S.A.

OmronHealthcare.com

© 2020 OMRON HEALTHCARE, INC.

Made in China

For questions, call: 1-800-634-4350