

BIOMONITOR IV
Technical Manual

BIOMONITOR IV

Technical Manual

BIOMONITOR IV Injectable Cardiac Monitor



BIOMONITOR IV
X-Ray identification

Radiopaque Identification

A radiopaque identification code is visible on standard X-ray, and identifies the injectable cardiac monitor:



CAUTION

Federal (U.S.A.) law restricts this device to sale by, or on the order of, a physician (or properly licensed practitioner).

Contents

Chapter 1: Device Description	1
Chapter 2: Indications	3
Chapter 3: Contraindications	5
Chapter 4: Warnings and Precautions	7
4.1 MR SAFETY INFORMATION	7
4.1.1 Patient Pre-MR Conditions	7
4.1.2 MR Scanner Limitations	7
4.1.3 Restrictions during the MR Scan	8
4.2 Implanted Pacemakers and Defibrillators	8
4.3 Medical Therapy	8
4.4 Storage and Handling	9
4.5 Home Monitoring	10
4.6 Electromagnetic Interference (EMI)	10
4.6.1 Home and Occupational Environments	10
4.6.2 Cellular Phones	11
4.6.3 Hospital and Medical Environments	12
4.7 Insertable Cardiac Monitor Explant and Disposal	12
Chapter 5: Programmable Parameters	13
5.1 Parameters	13
5.1.1 ProgramConsult	14
5.1.2 Atrial Fibrillation (AF)	16
5.1.3 Tachycardia	21
5.1.4 Bradycardia	22
5.1.5 Sudden Rate Drop (SRD)	23
5.1.6 Pause Duration	24
5.1.7 Symptom	24
5.1.8 Resting Rate Period	25

5.2	Home Monitoring (HM)	26
5.2.1	HM PID	26
5.2.2	Home Monitoring	26
5.2.3	Time of Transmission	26
5.2.4	Periodic Subcutaneous Electrocardiogram (sECG)	27
5.2.5	Daily Presenting Electrocardiogram (sECG)	27
5.2.6	Last Message	27
5.2.7	Episode Recording/Transmission	27
5.3	Patient Data	28
5.3.1	ID	28
5.3.2	First / Last Name	28
5.3.3	Date of Birth	29
5.3.4	Gender	29
5.3.5	Date of Implant	29
5.3.6	Hospital, City	29
5.3.7	Physician	29
5.3.8	NYHA	29
5.3.9	Symptom	30
5.3.10	Etiology	30
5.3.11	Remark	30
5.4	Sensing Settings	31
5.4.1	SensingConsult	32
5.4.2	Sensing High-Pass Filter	33
5.4.3	Input Signal Polarity	34
5.4.4	Input High-Pass Filter	34
Chapter 6:	Diagnostics	35
6.1	Diagnostics Overview	35
6.2	General Statistical Information	35
6.3	Activity	35
6.3.1	Rate Trends	36
6.3.2	Rate Histogram	36
6.3.3	Activity Trend	37
6.4	AF Details	37
6.4.1	AF Trends	37
6.4.2	AF Time of Occurrence	38
6.4.3	AF Duration	38
6.4.4	Ventricular Rate During AF	39
6.4.5	Ectopy count	39

6.5	Sensing	40
6.5.1	R-wave Trend	40
6.5.2	Noise Duration Trend	40
Chapter 7: Other Functions/Features		41
7.1	Home Monitoring	41
7.2	Transmission of Information	41
7.3	Patient Device	41
7.4	Transmitting Data	41
7.5	Types of Report Transmissions	42
7.6	Description of Transmitted Data	42
7.7	Patient Data Memory	43
7.8	Position Indicator	44
Chapter 8: Product Storage and Handling		45
8.1	Sterilization and Storage	45
Chapter 9: Follow-up Procedures		47
9.1	General Considerations	47
9.2	Real-time sECG Transmission	47
9.3	Follow-up Page	47
9.4	Recordings	48
9.5	sECG	49
9.5.1	Atrial Fibrillation	49
9.5.2	Tachycardia	50
9.5.3	Bradycardia	50
9.5.4	Pause	51
9.5.5	Symptom	51
Chapter 10: Elective Replacement Indication (ERI)		53
Chapter 11: Insertion/Removal		55
11.1	Opening the Sterile Container	55
11.2	Insertion	55

11.3	Removal	58
11.4	Reasons to Remove an Injectable Cardiac Monitor	58
Chapter 12: Remote Assistant III		59
12.1	General Information on the Remote Assistant III	59
12.2	Remote Assistant III Functional Testing	59
12.3	Getting to Know the Remote Assistant III	59
12.4	Triggering a Manual Recording	60
12.5	Battery LED Indicator Explained	60
12.6	Signal Transmission LED Explained	61
Chapter 13: Technical Data		63
13.1	Parameters	63
13.1.1	Atrial Fibrillation	63
13.1.2	Tachycardia	63
13.1.3	Bradycardia	63
13.1.4	Sudden Rate Drop (SRD)	64
13.1.5	Pause Duration	64
13.1.6	Symptom	64
13.1.7	Resting Rate Period	64
13.1.8	Home Monitoring	65
13.2	Materials in Contact with Human Tissue	65
13.3	Electrical Data/Battery	65
13.4	Mechanical Data	65
Chapter 14: Order Information		67

This page left intentionally blank

This page left intentionally blank

Chapter 1: Device Description

BIOMONITOR IV is a programmable, subcutaneous injectable cardiac monitor able to record subcutaneous ECGs (sECGs) and other physiological parameters.

The BIOMONITOR IV is designed to automatically record the occurrence of arrhythmias in a patient. Arrhythmia may be classified as atrial fibrillation (AF), bradyarrhythmia, pause, sudden rate drop, or tachycardia. In addition, the BIOMONITOR IV can be activated by the patient to record cardiac rhythm during symptomatic episodes.

Note - The BIOMONITOR IV subcutaneous ECG may differ from a surface ECG due to differences in electrode separation and device placement in the body.

BIOMONITOR IV detects a subcutaneous ECG from a pair of electrodes. These signals are filtered in two different ways. For detection of QRS complexes, the signals are filtered with a passband of 10-40 Hz in order to suppress T-waves, artifacts, and baseline drift at low frequencies, and myopotentials and EMI at high frequencies. The resulting signal is appropriate for QRS detection as other components of the signal have been suppressed. This signal naturally does not have a typical ECG morphology due to the bandpass. For waveform display (real-time streaming sECG with the physician's programmer and snapshots for review by the physician), a different passband is utilized to retain signal features that may have diagnostic value. This passband is 0.5 – 40 Hz, which is designed to retain morphological features of a typical ECG while still rejecting large low frequency artifacts and baseline drift.

The BIOMONITOR IV system consists of three main components:

1. BIOMONITOR IV insertable cardiac monitor - The BIOMONITOR IV is a small, leadless device that is typically inserted under the skin, in the chest. The device uses two electrodes on the body of the device to continuously monitor the patient's subcutaneous ECG. BIOMONITOR IV can store up to 96 minutes (67 min minimum) of subcutaneous ECG (sECG) recordings from both automatically detected arrhythmias and from manually triggered symptom episodes. When a patient experiences symptoms, the sECG recordings can be manually triggered by placing the Remote Assistant III over the BIOMONITOR IV. The injectable cardiac monitor is provided preloaded in an insertion tool. An incision tool is also provided.
2. BIOTRONIK Renamic / Renamic Neo Programmer – The programmer is used to set up the BIOMONITOR IV to detect arrhythmias. It also allows one to view, save, or print the stored information.
3. BIOTRONIK CardioMessenger® Smart is a telemetry patient device that forwards the data from the BIOMONITOR IV to BIOTRONIK's Home Monitoring Service Center.

BIOMONITOR IV may be used with BIOTRONIK Home Monitoring® technology, which is an automatic, wireless, remote monitoring system for management of patients with insertable cardiac monitors. When active, Home Monitoring enables the exchange of information about a patient's cardiac status from the implant to the Home Monitoring Service Center (HMSC) where the physician may log in to view the data. The HMSC can be used to provide the physician with advanced reports from the

Chapter 1

BIOMONITOR IV Technical Manual

implanted device and process them into a graphical and tabular format that is accessible via the internet platform HMSC. This information may help the physician optimize the therapy process, possibly providing earlier notification of clinically relevant events to help guide future therapy.

BIOTRONIK conducted the TRUST study to evaluate the safety and effectiveness of Home Monitoring. With the TRUST study, BIOTRONIK was able to show the following with regards to Home Monitoring:

- BIOTRONIK Home Monitoring information may be used as a replacement for device interrogation during in-office follow-up visits.
- A strategy of care using BIOTRONIK Home Monitoring with office visits when needed has been shown to extend the time between routine, scheduled in-office follow-ups of BIOTRONIK implantable devices in many patients. Home Monitoring data is helpful in determining the need for additional in-office follow-up.
- BIOTRONIK Home Monitoring provides early detection of arrhythmias.
- BIOTRONIK Home Monitoring provides early detection of silent, asymptomatic arrhythmias.
- Automatic early detection of arrhythmias and device system anomalies by BIOTRONIK Home Monitoring allows for earlier intervention than conventional in-office follow-ups.
- BIOTRONIK Home Monitoring allows for improved access to patient device data compared to conventional in-office follow-ups since device data is automatically collected and reported on a daily basis.

The implanted device's Home Monitoring function can be used for the entire operational life of the implanted device (prior to ERI).

NOTE: When ERI mode is reached, this status is transmitted and Home Monitoring® will be discontinued after two weeks.

Chapter 2: Indications

The BIOMONITOR IV is indicated to detect the following cardiac arrhythmias:

- Atrial fibrillation
- Bradycardia
- Sudden rate drop
- Tachycardia
- Pause

The BIOMONITOR IV is indicated for use in:

- Patients with clinical syndromes or situations at increased risk of cardiac arrhythmias
- Patients who experience transient symptoms that may suggest a cardiac arrhythmia

The device has not been tested for and it is not intended for pediatric use.

This page left intentionally blank

Chapter 3: Contraindications

There are no known contraindications for the insertion of the BIOMONITOR IV. However, the patient's particular physical or medical condition may dictate whether or not a subcutaneous, chronically inserted device can be tolerated.

This page left intentionally blank

Chapter 4: Warnings and Precautions

Consult the technical manuals for information about other devices used with the BIOMONITOR IV, including the CardioMessenger Smart or Renamic/Renamic Neo Programmer, and related accessories. Please keep technical manuals for later use.

Certain therapeutic and diagnostic procedures may cause undetected damage to an insertable cardiac monitor (ICM), resulting in malfunction or failure at a later time. Please note the following warnings and precautions:



MR Conditional - The cardiac monitor is labeled and certified MR conditional.

4.1 MR SAFETY INFORMATION

Conditions for an MR scan are provided below. Failure to adhere to provided patient position or scan time limitations may result in excessive tissue warming during an MR scan.

Cardiac data recorded by the implanted device during an MR scan may include artifacts that are due to the MR scan and not the patient's cardiac function. Exercise care when interpreting any such data.

4.1.1 Patient Pre-MR Conditions

The following requirements must always be fulfilled in order to perform an MR scan using BIOTRONIK's BIOMONITOR IV:

- There are no other active or abandoned cardiac implants (e.g., lead extensions, lead adapters or abandoned leads) in the patient's body.
- Other active or passive implants are permitted if they are identified as MR conditional by the manufacturer.

NOTE: An MR scan is permitted only if the product-specific conditions are met for all implants and if no metal implantable device longer than 5 cm is in the vicinity of the BIOTRONIK cardiac monitor within a distance of less than 4 cm.

- The device is located in the patient's chest area.

4.1.2 MR Scanner Limitations

The MR scanner has to meet the following conditions:

- Use of a clinical MR system with a cylindrical bore and a static magnetic field strength of 1.5 tesla or 3.0 tesla
- The slew rate of the MR scanner's gradient fields should not exceed 200 T/m/s per axis.
- Local receive-only coils may be used without restriction.
- Local circularly polarized or quadrature volume transmit-only or transmit-and-receive coils are only permitted for the head (including cervical spine) and extremities.

Chapter 4

BIOMONITOR IV Technical Manual

- Local circularly polarized or quadrature transmit-only or transmit-and-receive coils are NOT permitted for the thorax (the area from the shoulders to lower ribs), pelvis, or abdominal regions.
- Maximum spatial gradient of the static magnetic field specification must be ≤ 30 T/m (3,000 gauss/cm).
- Under worst case conditions, the BIOMONITOR IV is expected to produce a maximum temperature rise of $<4.5^{\circ}$ C after 30 minutes of continuous scanning.
- Image artifact and distortion can result from the presence of the BIOMONITOR IV device within the field of view. Image artifact and distortion resulting from the presence of the device within the field of view must be considered when selecting the field of view and imaging parameters. These factors must also be considered when interpreting the MR images.

4.1.3 Restrictions during the MR Scan

The following conditions must be met during the MR scan:

- The mean specific absorption rate (SAR) for the whole body as displayed by the MR scanner must not exceed 4.0 W/kg.
- The head absorption rate displayed by the MR scanner must not exceed 3.2 W/kg.

4.2 Implanted Pacemakers and Defibrillators

The BIOMONITOR IV is not intended for use in patients with an implanted pacemaker or defibrillator. If the patient has a co-implanted pacemaker or defibrillator, the automatic detection of arrhythmic episodes in the BIOMONITOR IV may be affected by the paced heart rhythm.

4.3 Medical Therapy

Before applying one of the following procedures, a detailed analysis of the advantages and risks should be made. Following the procedures, insertable cardiac monitor function must be checked.

Therapeutic Diathermy Equipment - Use of therapeutic diathermy equipment is to be avoided for insertable cardiac monitor patients due to possible heating effects of the insertable cardiac monitor and at the implant site. If diathermy therapy must be used, it should not be applied in the immediate vicinity of the insertable cardiac monitor.

Transcutaneous Electrical Nerve Stimulation (TENS) - Transcutaneous electrical nerve stimulation may interfere with insertable cardiac monitor function and is therefore not recommended. If necessary, the following measures may reduce the possibility of interference:

- Place the TENS electrodes as close to each other as possible.
- Place the TENS electrodes as far from the insertable cardiac monitor as possible.

Defibrillation - The following precautions are recommended to minimize the inherent risk of insertable cardiac monitor operation being adversely affected by defibrillation:

- The paddles should not be placed directly over the implant.
- The paddles should be placed anterior-posterior or along a line perpendicular to the axis formed by the insertable cardiac monitor.
- The energy setting should not be higher than required to achieve defibrillation.
- After defibrillation, evaluate the BIOMONITOR IV for proper function.

Radiation - Insertable cardiac monitor electronics may be damaged by exposure to radiation during radiotherapy. To minimize this risk when using such therapy, the insertable cardiac monitor should be protected with local radiation shielding.

Lithotripsy - Lithotripsy treatment should be avoided for insertable cardiac monitor patients since electrical and/or mechanical interference with the insertable cardiac monitor is possible. If this procedure must be used, the greatest possible distance from the point of electrical and mechanical strain should be chosen (25 cm minimum) in order to minimize a potential interference with the insertable cardiac monitor.

Ablation/Electrocautery - Ablation/Electrocautery - Position the grounding pad so that the current path does not pass through or near the device. When possible, a bipolar electrocautery system should be used. After ablation or electrocautery, evaluate the BIOMONITOR IV for proper function.

Transurethral resection of the prostate - It is recommended that the cautery ground plate be placed under the buttocks or around the thigh, but not in the thoracic area where the current pathway could pass through or near the cardiac monitor.

Hyperbaric Oxygen Therapy (HBOT) - Hyperbaric oxygen therapy (HBOT) for patients with BIOTRONIK CRM devices is not recommended due to the potential for damage or impaired function of the implant after exposure. The physician should conduct a risk-benefit analysis if HBOT treatment is necessary. The device specified in this manual has been tested to be in compliance with ISO 14708-2:2012, where the device is exposed to 40 cycles of ambient pressure up to 450 kPa (4.5 bar).

Therapeutic Ultrasound - The BIOMONITOR IV should not be exposed to therapeutic levels of ultrasound energy, as the active implantable medical device can inadvertently concentrate the ultrasound field and cause harm.

4.4 Storage and Handling

Failure to adhere to storage and handling recommendations may result in device damage or malfunction.

Storage (temperature) - Recommended storage temperature range is -10° to 45°C (14°-113°F). Exposure to temperatures outside this range may result in insertable cardiac monitor malfunction (see Section 8.1).

Handling - Do not drop. The monitor is preloaded into the insertion tool. If the tool is dropped onto a hard surface, return it to BIOTRONIK (see Section 8.1).

FOR SINGLE USE ONLY - Do not resterilize the insertable cardiac monitor, incision tool or insertion tool; they are intended for one-time use.

Device Packaging - Do not use the device if the packaging is wet, punctured, opened or damaged because the integrity of the sterile packaging may be compromised. Return the device to BIOTRONIK.

Storage - Store the device in a clean area, away from sources of disturbance to avoid damage to the device.

Temperature Stabilization - Allow the device to reach room temperature before programming or implanting the device. Temperature extremes may affect the initial device function.

Use Before Date - Do not implant the device after the USE BEFORE DATE because the device sterility and longevity may be compromised.

Sharp - Packaging includes an incision tool that is sharp and should be handled with care.

4.5 Home Monitoring

All BIOMONITOR IV devices can be used with BIOTRONIK's Home Monitoring® system. The Home Monitoring system enables wireless automatic transmission of information about a patient's cardiac status from the implanted device to the physician remotely.

Programming Overview

BIOTRONIK's Home Monitoring system is designed to notify clinicians in less than 24 hours about changes to the patient's condition or status of the implanted device. Updated data may not be available if:

- The patient's CardioMessenger® is unplugged or damaged and is not able to connect to the Home Monitoring system through an active cellular network.
- The CardioMessenger cannot establish a connection to the implanted device.
- The cellular network is not operational or the patient lives in a geographical area not covered by cellular networks.
- The Home Monitoring Service Center is off-line (upgrades are typically completed in less than 24 hours).

Patient's Ability - Use of the Home Monitoring system requires the patient and/or caregiver to follow the system instructions and cooperate fully when transmitting data.

If the patient cannot understand or follow the instructions because of physical or mental challenges, another adult who can follow the instructions will be necessary for proper transmission.

Use in Cellular Phone Restricted Areas - The CardioMessenger (transmitter/receiver) should not be utilized in areas where cellular phones are restricted or prohibited (i.e., commercial aircraft). Cellular network outages (including poor signal strength) prevent reliable connections.

4.6 Electromagnetic Interference (EMI)

The operation of any insertable cardiac monitor can be affected by certain environmental sources generating signals that resemble cardiac activity. In some cases the disturbance sources can couple sufficient energy to damage the insertable cardiac monitor.

BIOTRONIK insertable cardiac monitors have been designed to significantly reduce susceptibility to disturbance sources. However, due to the variety and complexity of sources creating interference, there is no absolute protection against disturbance sources. Generally, it is assumed that disturbance sources produce only minor effects, if any, in insertable cardiac monitor patients. If the patient presumably will be exposed to one of the following environmental conditions, then the patient should be given the appropriate warnings.

4.6.1 Home and Occupational Environments

The following equipment (and similar devices) may affect normal insertable cardiac monitor operation: electric arc welders, electric melting furnaces, radio/television and radar transmitters, power generating facilities, high voltage transmission lines, electrical ignition systems (also of gasoline powered devices) if protective hoods, shrouds, etc., are removed, electrical tools, anti-theft devices of shopping centers and electrical appliances, if not in proper condition or not correctly grounded and encased.

Patients should exercise reasonable caution in avoidance of devices which generate a strong electric or magnetic field. Some potential EMI sources include:

- **High Voltage Power Transmission Lines** - High voltage power transmission lines may generate enough EMI to interfere with insertable cardiac monitor operation if approached too closely.
- **Home Appliances** - Home appliances normally do not affect insertable cardiac monitor operation if the appliances are in proper condition and correctly grounded and encased.
- **Communication Equipment** - Communication equipment such as microwave transmitters, linear power amplifiers, or high-power amateur transmitters may generate enough EMI to interfere with insertable cardiac monitor operation if approached too closely.
- **Commercial Electrical Equipment** - Commercial electrical equipment such as arc welders, induction furnaces, or resistance welders may generate enough EMI to interfere with insertable cardiac monitor operation if approached too closely.
- **Electrical Appliances** - Electric hand-tools and electric razors (used directly over or on the skin covering the insertable cardiac monitor) have been reported to cause insertable cardiac monitor disturbances.
- **Electronic Article Surveillance (EAS)** - Equipment such as retail theft prevention systems may interact with the insertable cardiac monitor devices. Patients should be advised to walk directly through and not to remain near an EAS system longer than necessary.

4.6.2 Cellular Phones

Recent studies have indicated there may be a potential interaction between cellular phones and insertable cardiac monitor operation. Potential effects may be due to the radio frequency signal when the phone is within close proximity (within 6 inches [15 centimeters]) to the insertable cardiac monitor.

Based on testing to date, effects resulting from an interaction between cellular phones and the insertable cardiac monitors have been temporary. Simply moving the phone away from the inserted device will return it to its previous state of operation.

To minimize such interactions, patients having an inserted cardiac monitor who operate a cellular phone should:

- Maintain a minimum separation of 6 inches (15 centimeters) between a hand-held personal cellular phone and the inserted device. Portable and mobile cellular phones generally transmit at higher power levels compared to hand held models. For phones transmitting above 3 watts, maintain a minimum separation of 12 inches (30 centimeters) between the antenna and the inserted device.
- Patients should hold the phone to the ear opposite the side of the inserted device. Patients should not carry the phone in a breast pocket or on a belt over or within 6 inches (15 centimeters) of the inserted device as some phones emit signals when they are turned ON but not in use (i.e., in the listen or standby mode). Store the phone in a location opposite the side of the cardiac monitor.

Chapter 4

BIOMONITOR IV Technical Manual

4.6.3 Hospital and Medical Environments

Electrosurgical Cautery - Electrosurgical Cautery may inhibit insertable cardiac monitor sensing operation. If use of electrocautery is necessary, the current path (ground plate) should be kept as far away from the insertable cardiac monitor as possible.

Lithotripsy - Lithotripsy may damage the insertable cardiac monitor. If lithotripsy must be used, do not focus the beam near the insertable cardiac monitor.

External Defibrillation - External defibrillation may damage the insertable cardiac monitor. Attempt to minimize current flowing through the insertable cardiac monitor by following the precautions.

High Radiation Sources - High radiation sources such as cobalt 60 or gamma radiation should not be directed at the insertable cardiac monitor. If a patient requires radiation therapy in the vicinity of the insertable cardiac monitor, place lead shielding over the device to prevent radiation damage.

4.7 Insertable Cardiac Monitor Explant and Disposal

Device Incineration - Never incinerate an insertable cardiac monitor. Be sure the insertable cardiac monitor is explanted before a patient who has died is cremated (see Section 11.2).

Explanted Devices - Return all explanted devices to BIOTRONIK.

Chapter 5: Programmable Parameters

For a complete list of programmable parameters and the available settings, see Sections 5 or 13. Refer to the programmer manual for additional information.

5.1 Parameters

The Diagnostics/Home Monitoring page under Parameters, shown in Figure 1, allows set-up of the recording criteria of the device to include atrial fibrillation (AF), tachycardia, bradycardia, sudden rate drop, pause, and symptom events. With the BIOMONITOR IV, there are preconfigured programming sets that can be selected under program sets to automatically manage these recording parameters based on patient indication through the ProgramConsult. This feature, as well as the other programming options for the recording criteria will be discussed in detail below.




Figure 1: Main parameter screen

Chapter 5

5.1.1 ProgramConsult

The ProgramConsult menu is located under Program sets and allows for the selection of a preconfigured programming set based on common patient indications. Selecting one of these options provides programming assistance by bundling the suggested recording criteria parameters to a single choice. The preconfigured programs, as well as their associated parameters are summarized in Table 1.



Figure 2: ProgramConsult menu

Parameter	Syncope	Palpitations	AF Monitoring	Cryptogenic Stroke
Atrial fibrillation (AF)	ON	ON	ON	ON
AF sensitivity	Low	Medium	Medium	Medium
RR variability limit	12	12	12	12
Confirmation time	10	6	6	2
Bigeminy rejection	Aggressive	Aggressive	Standard	Aggressive
Ectopy rejection	ON	ON	ON	ON
AF recording	ON	ON	ON	ON
AF transmission	ON	ON	ON	Detect only
Tachycardia	ON	ON	ON	ON
Tachy limit	160	180	180	180
Tachy counter	16	32	48	48
Tachy recording	ON	ON	ON	ON
Tachy transmission	ON	ON	ON	ON
Bradycardia	ON	ON	ON	ON
Brady zone limit	35	30	30	30
Brady duration	20	30	30	30
Brady recording	ON	ON	ON	ON
Brady transmission	ON	ON	ON	ON
Sudden rate drop (SRD)	ON	OFF	OFF	OFF
SRD rate decrease	50			
SRD sensitivity	Low			
SRD recording	ON			
SRD transmission	ON			
Pause	ON	ON	ON	ON
Pause duration	3	5	5	5
Pause recording	ON	ON	ON	ON
Pause transmission	ON	ON	ON	ON

Table 1: ProgramConsult programs and parameters

5.1.2 Atrial Fibrillation (AF)

The Atrial Fibrillation Menu (Figure 3) allows the user to program AF detection ON or OFF as well as set the detection criteria.




Figure 3: AF menu

This section includes:

- Atrial Fibrillation Parameters
- Atrial Fibrillation (AF): ON or OFF
- AF Sensitivity: Low, Medium, High

Allows the user to select preset criteria for AF determination. These are shown in the following Table 2. The user can also select non-preset options to make an Individual program.

Low is the least sensitive setting. This setting requires greater instability and more intervals for the device to declare AF. Conversely, the High setting requires less instability and fewer events to declare AF.




Figure 4: AF sensitivity

Parameter	Range	Low	Medium (Default)	High
Atrial fibrillation detection	OFF, ON	On	On	On
Detection/termination window size	8/16, 16/24, 24/32	16/24	8/16	8/16
Detection intervals	5...(2)...23	11	5	5
Number of detection windows	1...(1)...4	3	2	1
Termination intervals	1...(2)...7	5	1	1
Number of termination windows	1...(1)...4	2	2	3
RR variability limit	6, 9, 12, 15, 18 %	12	12	12
Confirmation time	0.5, 1...(1)...6,10,20,30 min	6	6	6
Bigeminy rejection	OFF, Standard, Aggressive	Standard	Standard	Standard
Ectopy rejection	OFF, ON	ON	ON	ON

Table 2: Parameter summary for AF sensitivity

RR Variability

This parameter represents the maximum percentage of variation between V_s - V_s cycle lengths to be considered stable by the device. The smaller the value, the greater the likelihood of AF being declared. Intervals greater than the RR Variability value from the mean cycle length will be considered AF intervals and count towards the detection and termination threshold.




Figure 5: RR variability limit

Detection/Termination Window Size

The number of cycle lengths used to determine detection and termination of AF. Figure 6 shows the selectable values for Detection and Termination. For example, a Detection value of 8 means the device monitors groups of 8 cycle lengths to determine the RR variability by comparing each of those 8 events to the variability limit value.

If the number of events that are determined to be unstable exceed the programmed Detection interval value, AF suspicion criterion is met.

A Termination criterion of 16 means the device is monitoring groups of 16 events. If the number of unstable events is greater than the programmed Termination value, the rhythm will continue to be considered unstable (AF).

The Detection/Termination windows are not sliding windows, but consecutive windows.




Figure 6: Detection/termination window

Number of Detection Windows

Number of consecutive windows that are required to be determined unstable for the device to start the confirmation time.




Figure 7: Number of detection windows

Number of Termination Windows

Number of consecutive windows that are required to be determined stable for the device to terminate the episode.

Number of termination windows

1	Close
2	
3	
4	

Figure 8: Number of termination windows

Detection Intervals

The Detection Intervals represents the number of intervals that must be unstable within the programmed detection window for the rhythm to be considered unstable.

Detection intervals

5	7	9	Close
11	13	15	
17	19	21	
23			

Figure 9: Detection intervals

Termination Intervals

The Termination Intervals represents the maximum number of unstable events within the programmed termination window allowed to terminate an AF episode. If more than the programmed number of unstable intervals are present, the device will continue to declare an AF event active.




Figure 10: Termination Intervals

Confirmation time

The time before a recording of the AF event occurs. If the events are detected but do not reach the confirmation time period (suspicion phase), the event will not be counted.




Figure 11: Confirmation time

Bigeminy Rejection

Short-long interval patterns indicative of bigeminy may be detected as AF. The bigeminy rejection parameter is designed to recognize these periodic interval patterns. When enabled to Standard, the bigeminy rhythms are recognized and prevented from triggering AF detections. The Aggressive setting is provided for more comprehensive filtering of complex bigeminy rhythms. The default setting for the bigeminy rejection parameter is Standard.

NOTE: If bigeminy rejection is set to Aggressive, the AF sensitivity may be reduced.

Ectopy Rejection

An ectopic event is declared when a short interval is followed by a long interval. When Ectopy Rejection is set to ON, an event that is determined to be an ectopic event will not be considered unstable for the purpose of AF detection or confirmation.

5.1.3 Tachycardia

BIOMONITOR IV may be programmed to record tachycardia events using a rate limit and counter for criteria. Both the tachy limit and tachy counter criteria must be met for an event to be classified as a tachy episode. An event meeting the criteria would record an sECG and update the counters on the Diagnostics section of the device.

Tachy Limit

This parameter value represents the lower rate limit required to be considered a tachy episode.

Tachy limit		
100	110	120
130	140	150
160	170	180
190	200	

Figure 12: Tachy limit

Tachy Counter

This parameter value represents the count limit for tachycardia classification. This is an up/down counter. Each event slower than the tachy limit decreases the count by 1, while each event faster than the tachy limit increments the counter by 1.




Figure 13: Tachy counter

5.1.4 Bradycardia

Brady Zone Limit

Rates determined to be below the programmed Brady Zone Limit will be classified as a bradycardia event. In addition to the rate limit, the rate must also meet the Brady Duration limit. This prevents single slow events from being classified as a bradycardia episode.




Figure 14: Brady zone limit

Brady Duration

The Brady Duration is the time over which the average heart rate is assessed. When the average heart rate is below the programmed bradycardia zone limit for the device, bradycardia is confirmed.




Figure 15: Brady duration

5.1.5 Sudden Rate Drop (SRD)

SRD Rate Decrease

This parameter value represents the percentage in rate decrease that triggers a Sudden Rate Drop event. The device compares the average rate of the most recent events (rate-drop intervals) and compares it to the average rate of the previous events (baseline intervals).




Figure 16: SRD rate decrease

SRD Sensitivity

This parameter programs preset value setting for baseline intervals and rate-drop intervals.




Figure 17: SRD sensitivity

Baseline Intervals

This parameter value represents the number of averaged intervals to determine a baseline rate for sudden rate drop determination.




Figure 18: Baseline intervals

Chapter 5

BIOMONITOR IV Technical Manual

Rate-drop Intervals

This parameter value represents the number of averaged intervals to determine a change in the heart rate. It uses the most recent events and determines the average rate of those events to determine the rate-drop rate value.




Figure 19: Rate-drop intervals

5.1.6 Pause Duration

Pause Duration

The minimum total duration in seconds between R waves for the device to declare a Pause event.




Figure 20: Pause duration

5.1.7 Symptom

Symptom

This is an ON/OFF feature which allows a patient to record an sECG by placing the Remote Assistant III over the device and pressing the button.




Figure 21: Symptom

5.1.8 Resting Rate Period

Start Resting Period

The Start Resting Period is the time the device starts collecting heart rate information for the Rate trend diagnostic. The default time is 2 A.M. for data recording. This time was chosen to reduce the chance of patient activity interfering with data collection.

Start resting period				
00:00	01:00	02:00	03:00	04:00
05:00	06:00	07:00	08:00	09:00
10:00	11:00	12:00	13:00	14:00
15:00	16:00	17:00	18:00	19:00
20:00	21:00	22:00	23:00	

Figure 22: Resting period duration

Resting Period Duration

The Resting Period Duration parameter is the time duration the data is collected for the Resting Rate trend diagnostic. BIOMONITOR IV collects resting heart rate values in 10-min blocks of time. The lowest average collected over the recording period is used as the statistical point for that given day.

Resting period duration		
1	2	3
4	5	6
7	8	9
10	11	12

Figure 23: Resting period duration

Chapter 5

5.2 Home Monitoring (HM)

The availability of parameters and parameter values is determined by the software used for programming/interrogating the injectable cardiac monitor.

5.2.1 HM PID

The HM PID is the Product Identification Number. This is a unique ID number for each product and is used when registering a patient to the Home Monitoring Service Center (HSMC).

5.2.2 Home Monitoring

Programmable ON or OFF.

5.2.3 Time of Transmission

By default, the BIOMONITOR IV will transmit all data and a daily trend report between 1:00 A.M. and 2:00 A.M. daily. Transmission time is also programmable by the user and is based on a 24-hour clock. It is important to keep in mind that the programmer updates the BIOMONITOR IV time based on the programmer time. If the programmer time is different than the local time, the transmission time may be different than expected.

Time of transmission						
Std.						00:00
00:30	01:00	01:30	02:00	02:30	03:00	03:30
04:00	04:30	05:00	05:30	06:00	06:30	07:00
07:30	08:00	08:30	09:00	09:30	10:00	10:30
11:00	11:30	12:00	12:30	13:00	13:30	14:00
14:30	15:00	15:30	16:00	16:30	17:00	17:30
18:00	18:30	19:00	19:30	20:00	20:30	21:00
21:30	22:00	22:30	23:00	23:30		

Figure 24: Time of transmission

5.2.4 Periodic Subcutaneous Electrocardiogram (sECG)

The BIOMONITOR IV can send sECGs up to 60 s each (40 s minimum) in length periodically based on user preference. This allows the user to assess sECGs routinely, even when no events have occurred. The schedule for the Periodic sECG is configurable through the Home Monitoring Service Center website.

When the option of Selection is made, the user can enter up to five specific dates on which to send a periodic sECG. Following the last programmed periodic sECG, the device will revert to sending sECGs every 30 days until new dates are entered.

5.2.5 Daily Presenting Electrocardiogram (sECG)

The BIOMONITOR IV sends one sECG of up to 6 seconds in length daily, which will be presented on the Status Summary page on the Home Monitoring Service Center website. This allows the user to assess a short sECG for the last received monitoring period. The daily presenting sECG is not configurable by the user.

5.2.6 Last Message

Message Type

This box shows the last message type created by the device.

Message Created On

This parameter shows the date and time the last message was created. The clock time is based on a 24-hour clock.

Send Test Message

When performing the “Send Test Message” function, a note will appear with the following message:

“Please remove programmer head for 10 seconds to allow implant to send test message. Afterwards, please interrogate to update status.”

Once the OK button is pressed, a “programming was successful” message will appear on the bottom- left corner of the screen.




Figure 25: Last message

5.2.7 Episode Recording/Transmission

This section provides an overview of which triggers are currently programmed ON for recordings, and also what triggers are set to transmit. Only if the device has recordings enabled for a particular trigger can the HM transmission option be selectable.

Trigger type	Detection	Recording	Transmission
Atrial fibrillation	Medium	ON	ON
Tachycardia [bpm]	180	ON	ON
Bradycardia [bpm]	40	ON	ON
Sudden rate drop [%]	OFF		
Pause [s]	3	ON	ON
Symptom	ON	ON	Nightly

Figure 26: HM episode trigger

Chapter 5

BIOMONITOR IV Technical Manual

The user can modify the Recording and Transmission options. The only programmable options are ON and OFF with the exception of Atrial Fibrillation transmissions and symptom transmissions. Atrial Fibrillation transmissions also include a Detect Only option. This option will transmit just the sECG for detection and not the one for termination. Symptom transmissions also include an On Demand option. This option, if enabled, will allow for up to two symptom sECG recordings to be transmitted per day immediately after the recording is taken by the patient if they are in range of their CardioMessenger.




Figure 27: Atrial fibrillation transmission options




Figure 28: Symptom transmission options

5.3 Patient Data

This section allows the user to add patient, physician, hospital and other information. This information is stored in the device and can be accessed with any compatible programmer.

The data in this section can be modified at any time.

5.3.1 ID

This section allows the user to input up to a 12-digit alphanumeric code to serve as a patient identifier. This may be a medical records number or a study number if the patient is enrolled in a study.

5.3.2 First / Last Name

These sections allow the user to input the patient's first and last name into the memory of the device. This is a free text box, allowing up to 20 characters for the first name, as well as for the last name.

Enter the patient's name and select the enter key.

5.3.3 Date of Birth

This section allows the user to input the patient's birth date. The birth date is entered as MM/DD/YYYY. When initially accessed, the current day will be displayed. The date can be changed using the following methods:

- Selecting the keypad icon to the left of the OK button will bring up a number keypad allowing the user to manually input the date.
- The day can be selected simply by touching the appropriate day on the screen.
- Pressing the month will bring up a listing of the 12 months, and the user can select the appropriate month.
- Selecting the year will bring up a numeric keypad, allowing the user to enter a year.
- The double arrow will change the year by one value each time it is touched. The left double arrows decrease the value and the right double arrows increase the value.
- The single arrow will change the month by one. The left arrow decreases the value and the right arrow increases the value.

Once the date is entered, select the OK button.

5.3.4 Gender

This section allows the user to select the patient's gender.

5.3.5 Date of Implant

The implantation date is entered by the user.

5.3.6 Hospital, City

The hospital name and city name can be added. As with entering the patient's name, up to 20 characters are available to add hospital and city information.

5.3.7 Physician

The physician's name can be added. As with the patient's name, up to 20 characters are available to add physician information. It is a good idea to add the physician's first name also to help prevent confusion.

5.3.8 NYHA

This refers to the New York Heart Association classification. A value can be entered if it is known.




Figure 29: NYHA

Chapter 5

BIOMONITOR IV Technical Manual

5.3.9 Symptom

This section allows the user to select one or multiple symptoms related to the patient. Selecting a symptom will result in a check mark appearing in the box to the left. Once completed, press the OK button. The selection(s) will appear on the main patient page.

Symptom	
<input type="checkbox"/>	A1 Rhythm unspecified
<input type="checkbox"/>	A2 Uncoded
<input type="checkbox"/>	B1 Syncope
<input type="checkbox"/>	B2 Dizzy spells
<input type="checkbox"/>	B3 Bradycardia
<input type="checkbox"/>	C1 Tachycardia
<input type="checkbox"/>	D1 Prophylactic
<input type="checkbox"/>	D2 Heart failure
<input type="checkbox"/>	D3 Cerebral dysfunction
<input type="checkbox"/>	D4 Angina pectoris
<input type="checkbox"/>	D5 After electrical resuscitation

Figure 30: Symptom

5.3.10 Etiology

This section allows the user to select an etiology related to the patient.

Etiology	
<input type="checkbox"/>	A1 Unspecified
<input type="checkbox"/>	A2 Uncoded
<input type="checkbox"/>	B1 Unknown
<input type="checkbox"/>	B2 Conduction tissue fibrosis
<input type="checkbox"/>	C1 Ischemic
<input type="checkbox"/>	C2 Post infarction
<input type="checkbox"/>	D1 Congenital
<input type="checkbox"/>	E1 Surgical complication
<input type="checkbox"/>	E2 Surgical
<input type="checkbox"/>	E3 Ablation
<input type="checkbox"/>	E4 Drug related
<input type="checkbox"/>	F1 Carotid sinus syndrome
<input type="checkbox"/>	F2 Vasovagal syndrome
<input type="checkbox"/>	F3 Orthostatic hypotonia
<input type="checkbox"/>	G1 Endocarditis
<input type="checkbox"/>	G2 Myocarditis
<input type="checkbox"/>	G3 Valvular heart disease
<input type="checkbox"/>	G4 Cardiomyopathy
<input type="checkbox"/>	G5 Hypertrophic cardiomyopathy
<input type="checkbox"/>	G6 Dilated cardiomyopathy
<input type="checkbox"/>	G7 Heart transplantation
<input type="checkbox"/>	G8 Radiation therapy

Figure 31: Etiology

5.3.11 Remark

This section allows the user to input a remark for the patient up to 42 characters in length.

5.4 Sensing Settings

This section allows the user to change the SensingConsult, sensing filter settings, and stored/real-time signal choices.




Figure 32: Sensing setting

Chapter 5

BIOMONITOR IV Technical Manual

5.4.1 SensingConsult

SensingConsult is a feature that allows the user to select a sensing profile to match the patient condition. Selecting one of these options results in the device automatically adjusting the threshold decay, reduction time, and threshold percentages to optimize sensing for that particular patient presentation. A preview image of the sensing profile is included on the left side of the screen.




Figure 33: SensingConsult options

5.4.2 Sensing High-Pass Filter

This section allows the user to change the filter setting of the signal. The higher the value, the more of the baseline and T-wave signals are removed from the sensing signal. This may be used if the baseline signal wanders or if oversensing from T-waves or P-waves occurs.




Figure 34: High-pass filters

Chapter 5

5.4.3 Input Signal Polarity

It is possible to select normal or inverted polarity for the sECG channel input. This will take effect for stored sECGs and streamed real-time sECGs.




Figure 35: Input signal polarity selection screen

5.4.4 Input High-Pass Filter

The sECG channel input contains a selectable analog high-pass filter with a cutoff value of 0.05 Hz or 0.5 Hz. This will take effect for stored sECGs and streamed real-time sECGs.




Figure 36: High-pass filter selection screen

Chapter 6: Diagnostics

6.1 Diagnostics Overview

BIOMONITOR IV can store a variety of statistical information. The various statistics consist of such features as rate histograms, rate trends, and activity trends, which are described in the following sections.

AF Details

- AF trends
- AF time of occurrence
- AF duration
- Ventricular rate during AF
- Ectopy Count (as displayed on the Renamic/Renamic Neo Programmer)*

Activity

- Rate trends
- Rate histograms
- Activity trend

Sensing

- R-wave trend
- Noise duration trend

6.2 General Statistical Information

The BIOMONITOR IV statistics modes are always in operation and cannot be selected OFF.

The counters within the statistic features are reset each time the injectable cardiac monitor is permanently programmed.

The histogram information is a 240-day duration. Afterwards, the oldest data are overwritten. Ongoing episodes are not counted.

6.3 Activity

The Activity diagnostic provides information related to heart rate, heart rate at rest, variability, rate histograms and activity.

Data is collected for the most recent 240 days. The user can look at information for a specific day by using the left/right arrows on the lower left screen or by simply touching on the screen. The date is listed at the bottom of the graph with the data results at the top of the graph.

* On the Home Monitoring website, this will be displayed as PAC/PVC burden and PVC count.

6.3.1 Rate Trends

Heart rate trends provide information related to heart rate, mean heart rate at rest and heart rate variability. Data is collected for the most recent 240 days. The user can look at information for a specific day by using the left/right arrows on the lower left screen or by simply touching on the screen. The date is listed at the bottom of the graph with the data results at the top of the graph.

Heart rate information is based on the daily average heart rate and is displayed as a single data point for the day.

BIOMONITOR IV collects resting heart rate values in 10-minute blocks of time during the mean heart rate at rest recording time. The lowest average collected over the recording period is used as the statistical point for that given day.

Heart rate variability is calculated using SDANN. Data is collected in five minute windows and calculated as a single daily data point.



Figure 37: Rate trends

6.3.2 Rate Histogram

The Rate histogram, shown in Figure 38, provides the percentage of activity in each rate bin for the BIOMONITOR IV. Rate bins are divided into 10 bpm increments.




Figure 38: Rate histograms

6.3.3 Activity Trend

The Activity trend, shown in Figure 39, displays the daily percentage of activity as detected by the motion sensor of the device.



Figure 39: Activity trend

6.4 AF Details

6.4.1 AF Trends

The AF trends diagnostic provides information related to the number and duration in hours of AF events on a daily basis.




Figure 40: AF trends

6.4.2 AF Time of Occurrence

The time of occurrence, shown in Figure 41, summarizes the times of day that atrial tachyarrhythmia episodes began and is broken into three-hour time blocks. Knowing the time of day when atrial tachyarrhythmias begin may help determine whether a particular event will precipitate the tachyarrhythmia.

The total number of events is listed at the bottom of the graph.




Figure 41: Time of occurrence

6.4.3 AF Duration

AF duration shows the length of each AF episode in time bins and provides a percentage of the episodes which occur in each time bin versus the total number of episodes. Ongoing episodes are not counted on the graph.




Figure 42: AF duration

6.4.4 Ventricular Rate During AF

The ventricular rate during AF graph provides the mean and the maximum heart rate during AF. Large differences in the mean and maximum rates may indicate an irregular ventricular response during the AF while small differences may imply that ventricular rate is more stable during AF.




Figure 43: Ventricular rate during AF

6.4.5 Ectopy count

Ectopy count shows the number of ectopic events (AES, VES) per day, not including ectopy events which occurred during suspected or confirmed AF.



Figure 44: Ectopy count

6.5 Sensing

6.5.1 R-wave Trend

The R-wave trend provides average daily R-wave measurement values for up to 240 days.



Figure 45: R-wave trend

6.5.2 Noise Duration Trend

The noise duration trend provides the amount of noise sensed daily by the device, expressed as a percentage of time per day by the BIOMONITOR IV. A high percentage of noise events could interfere with the BIOMONITOR IV's ability to detect arrhythmias.




Figure 46: Noise duration trend

Chapter 7: Other Functions/Features

BIOMONITOR IV injectable cardiac monitors offer many additional functions and features to assist the physician in the care of the patient.

7.1 Home Monitoring

Home Monitoring enables the exchange of information about a patient's cardiac status from the cardiac monitor to the physician. Home Monitoring can be used to provide the physician with reports from the BIOMONITOR IV and can process them into graphical and tabular format called a Home Monitoring Report. This information helps the physician optimize the diagnostic process, as it allows the patient to be scheduled for additional clinical appointments between regular follow-up visits if necessary.

7.2 Transmission of Information

The cardiac monitor transmits information with a small transmitter, which has a range of about 6 feet (2 meters). The patient's cardiac monitor data are sent daily to the corresponding patient device and periodic sECGs are sent in configurable intervals when Home Monitoring is programmed ON. The minimal distance between the cardiac monitor and the patient device must be 8 inches (20 cm).

7.3 CardioMessenger

The CardioMessenger is designed for use in or away from the home. Power is supplied by a standard wall plug. The patient device can be placed on the patient's nightstand or within 6 ft of where data transmission is to occur. CardioMessengers are cellular-capable and allow for transmission of data over a cellular network. For additional information about the CardioMessenger, please refer to its manual.

7.4 Transmitting Data

The cardiac monitor's information is digitally formatted by the BIOTRONIK Home Monitoring Service Center and processed into a concise report called a Home Monitoring Report. The Home Monitoring Report, which is adjusted to the individual needs of the patient, contains current and previous cardiac monitor data. The Home Monitoring Report is sent to the attending physician over the Internet. For more information on registering for Home Monitoring, contact your BIOTRONIK sales representative.

The password protected BIOTRONIK Home Monitoring website can be accessed at the following URL: **www.biotronik-homemonitoring.com**

An online help menu is available in order to assist with the use of the Home Monitoring website.

Chapter 7

BIOMONITOR IV Technical Manual

Additionally, the attending physician may register to be informed of the occurrence of an Event Triggered Message through email or SMS (i.e., mobile phone) with a brief text message. If registered for Internet availability, the patient's detailed cardiac monitor data can then be viewed by logging onto the Home Monitoring website.

7.5 Types of Report Transmissions

When the Home Monitoring function is activated, the transmission of a report (Home Monitoring Report) from the implant can be triggered as follows:

- Trend report—the time period (daily) initiates the report.
- Event report—the BIOMONITOR IV detects certain events, which initiate a report.

Trend Report

The time of the report transmission is programmable. For periodic messages, the time can be set anywhere between 00:00 and 23:30 hours. It is recommended to select a time between 0:00 and 4:00.

The length of the time interval (monitoring interval) is preset to “daily.” For each monitoring interval, a data set is generated in the cardiac monitor and the transmission is initiated at the designated time.

Event Report

When certain cardiac and technical events are detected by the cardiac monitor, a report transmission is automatically triggered. This is described as an “event message” as part of the daily transmission.

The following clinical and technical events initiate a Home Monitoring message transmission:

- Event recording
- ERI detected

NOTE: The attending physician can go to the Home Monitoring website to change or modify the events they wish to be informed of.

7.6 Description of Transmitted Data

The Monitoring Interval

The monitoring interval is the time period since the last periodic message was transmitted. In a periodic report, the monitoring interval since the previous periodic report is 24 hours.

The following data are transmitted for the Home Monitoring Report by the Home Monitoring system, when activated. In addition to the medical data, the serial number of the BIOMONITOR IV is also transmitted.

Device Status & Home Monitoring Settings

Contains device and message identifying values that pertain to the cardiac monitor and Home Monitoring:

- Implantation Date
- Device Status
- Remaining capacity for ERI calculation (done by the Service Center)

- Last follow-up
- Device Serial Number
- Message Creation Date/Time
- Device settings

Physiologic data

- Heart rate
- Heart rate variability
- Patient activity
- Temperature
- Daily statistics (including PAC/PVC burden)
- Event recordings

The temperature measurement feature alerts the physician to increased average temperature over 1-5 days compared to a 30-day baseline. This may indicate fever and allow physicians to determine whether further screening or follow-up is needed. Comparison of temperature measurements using this device and a conventional thermometer in swine warmed to increase their temperature by 1°C to mimic fever showed agreement within 0.1°C. Comparison of subcutaneous temperature measurements in humans to a conventional thermometer has not been performed.

Event recordings include up to 6/day for daily transmission and up to 2/day for optional on-demand patient triggered recordings transmitted immediately.

7.7 Patient Data Memory

Individual patient data can be stored in the injectable cardiac monitor's memory. The stored data is automatically displayed upon each interrogation. The patient data memory contains the following data categories:

- Patient ID (Code)
- Patient Name
- Date of Birth
- Gender
- Symptom
- Etiology
- Physician
- Implantation Date
- NYHA Class
- Hospital
- City
- Phone
- Remark

Chapter 7

BIOMONITOR IV Technical Manual

Symptom and etiology are specified using the European PASSPORT code system. The PASSPORT code is an identification system of two character codes that represent specific conditions. A listing of the codes available with definitions is displayed on the screen of the programmer when patient data is selected.

When the patient data screen is entered symptom or etiology may be entered, and can be accessed following interrogation to check code definition.

When the patient data screen is printed, the date of last follow-up is automatically given on the print-out.

7.8 Position Indicator

The position indicator facilitates positioning of the programmer head. The programmer optically and acoustically indicates whether the programmer head is in communication with the injectable cardiac monitor.

CAUTION

EMI – Computerized systems are subject to EMI or “noise”. In the sources of such disturbance, telemetry communication may be interrupted and prevent programming.

Chapter 8: Product Storage and Handling

8.1 Sterilization and Storage

The injectable cardiac monitor is shipped in a cardboard box, equipped with a quality control seal, and product information label. The label contains the model specifications, technical data, serial number, expiration date, and sterilization and storage information of the injectable cardiac monitor. The monitor is preloaded in the insertion tool and an incision tool is also included in the package.

The injectable cardiac monitor and its accessories have been sealed in a container and gas sterilized with ethylene oxide. To assure sterility, the container should be checked for integrity prior to opening. If a breach of sterility is suspected, return the injectable cardiac monitor to BIOTRONIK.

CAUTION

Storage (temperature) – Recommended storage temperature range is -10° to 45°C (14°-113°F). Exposure to temperatures outside this range may result in cardiac monitor malfunction.

Handling – Do not drop. – If an unpackaged injectable cardiac monitor is dropped onto a hard surface, return it to BIOTRONIK.

SHARP – Packaging includes an incision tool that is sharp and should be handled with care.

CAUTION

FOR SINGLE USE ONLY – Do not resterilize the injectable cardiac monitor or accessories packaged with the cardiac monitor, they are intended for one-time use.

Device Packaging – Do not use the device if the packaging is wet, punctured, opened or damaged because the integrity of the sterile packaging may be compromised. Return the device to BIOTRONIK.

Storage – Store the device in a clean area, away from sources of electromagnetic interference (EMI) to avoid damage to the device.

Use Before Date – Do not implant the device after the USE BEFORE DATE because the device may have reduced longevity.

If a replacement injectable cardiac monitor is needed, contact your local BIOTRONIK representative.

This page left intentionally blank

Chapter 9: Follow-up Procedures

9.1 General Considerations

The injectable cardiac monitor follow-up serves to monitor and provide information related to the patient's rhythm.

The follow-up intervals are, therefore, primarily determined by medical judgment.

The following notes are meant to stress certain product features, which are of importance for follow-up visits. For detailed information on follow-up procedures and medical aspects, please refer to the pertinent medical literature.

NOTE: In order to enable full device functionality, including statistics functions and ERI detection, transmit a permanent program after insertion by pressing the [Transmit/Program] button.

9.2 Real-time sECG Transmission

The injectable cardiac monitors provide real-time transmission of the subcutaneous electrogram (sECG) to the programmer. The sECGs may be transmitted to the programmer via the programming head positioned over the inserted monitor. They are then displayed together with surface ECG and markers on the programmer screen and printed on the ECG recorder.

9.3 Follow-up Page

The follow-up page shown in Figure 47 provides information including the last follow-up date, the battery status, number of the diagnostics recordings and Home Monitoring status.




Figure 47: Follow-up page

Chapter 9

The ECG and sECG signal display may be adjusted to make viewing easier by pressing on the icon shown in Figure 47.

9.4 Recordings

The Recordings page provides a list of stored episodes since the last time it was cleared. Information includes the time and date of the event, the duration, the type of event, mean heart rate and a sECG link to the recording.

No.	Date	Duration	Trigger	MR	View
30	04/10/2017 08:46	--:--:--	Brady	48	
29	04/10/2017 08:29	-	AF termination	68	
28	04/09/2017 06:55	-	Pause		
27	04/09/2017 05:17	00:12:24	Brady	37	
26	04/09/2017 03:33	-	SRD	75/42	
25	04/08/2017 18:12	* 1 days	AF detection	66	
24	04/08/2017 16:49	-	AF termination	75	
23	04/08/2017 16:19	-	Symptom	183	
22	04/08/2017 16:17	00:32:08	AF detection	72	
21	04/08/2017 14:47	-	AF termination	70	
20	04/08/2017 06:30	08:17:34	AF detection	65	
19	04/07/2017 17:12	-	Pause		
18	04/07/2017 17:09	-	AF termination	55	
17	04/07/2017 17:03	00:00:44	Brady	39	
16	04/07/2017 17:01	00:08:18	AF detection	48	
15	04/02/2017 22:15	* 00:12:58	Brady	38	
14	04/02/2017 20:11	-	AF termination	72	
13	04/02/2017 06:13	-	Symptom	183	

[Set Diagnostic and Home Monitoring parameters](#)

Print	Help	Restart
-------	------	---------

Figure 48: Recordings page

9.5 sECG

BIOMONITOR IV can store up to 96 minutes (67 minutes minimum) of sECGs. The types of sECG recording include Tachycardia, Bradycardia, SRD, AF, Pause and Symptom events.

If sECG snapshots of all arrhythmia types are available, the minimum number of each type of snapshot in the device is the following:

Arrhythmia Type	Number of Snapshots	Episode recording scheme
AF	3	Oldest, newest, longest
Tachycardia	3	Oldest, newest, longest
Bradycardia	3	Oldest, newest, longest
SRD	3	Two newest, oldest
Pause	3	Two newest, oldest
Symptom	4	Four most recent

Table 3: Minimum number of episode snapshots

Examples of the different recordings are provided in the following sections. Figures 49-53 are for demonstration purposes only and are not clinically derived.

9.5.1 Atrial Fibrillation

Figure 49 shows an example of an atrial fibrillation recording with the sECG and marker channels. The device will record both when the episode meets detection criteria and when it meets termination criteria. The user can scroll through the Holter and print only a section or the entire recording may be printed.



Figure 49: Atrial fibrillation sECG

Chapter 9

BIOMONITOR IV Technical Manual

9.5.2 Tachycardia

Figure 50 shows an example of a tachycardia recording. The black vertical bar indicates when the Tachy criteria were met.




Figure 50: Tachycardia sECG

9.5.3 Bradycardia

Figure 51 shows an example of a bradycardia recording.




Figure 51: Bradycardia sECG

9.5.4 Pause

Figure 52 shows a Pause recording.




Figure 52: Pause sECG

9.5.5 Symptom

Figure 53 shows a recording from a Symptom event.




Figure 53: Symptom sECG

This page left intentionally blank

Chapter 10: Elective Replacement Indication (ERI)

The service time of BIOMONITOR IV may vary based on several factors, including battery properties, storage time, programmed parameters and circuit operating characteristics. Service time is the time from beginning of service (BOS) to the Elective Replacement Indication (ERI). To assist the physician in determining the optimum time for injectable cardiac monitor replacement, an elective replacement indicator is provided that is activated when the battery cell capacity drops to a predetermined level. The following table defines the different service cycles (at standard settings at 37°C). The beginning of the replacement cycle is displayed on the programmer after injectable cardiac monitor interrogation and appears on the printout. Table 4 shows the service cycle definitions.

Abbreviation	Service Cycle	Definition
BOS	Beginning of Service	Normal service cycle; battery in good condition
ERI	Elective Replacement Indication	Identifies the time of elective replacement indication
EOS	End of Service	Identifies the end of the elective replacement indication period.

Table 4: Service cycle definitions

Table 5 shows the expected longevity (in years) from BOS to ERI for the BIOMONITOR IV injectable cardiac monitors. The programmer software for the BIOMONITOR IV injectable cardiac monitors provides a fuel gauge to provide information related to the battery status.

Injectable Cardiac Monitor	Standard (BOS - ERI) in Years
BIOMONITOR IV	5

Conditions: 6 months shelf life, 1 daily automatic sECG HM upload, and 2 patient-triggered symptoms sECG uploads per month.

Table 5: Nominal BIOMONITOR IV longevity

Note: Nominal longevity of BIOMONITOR IV is reduced by 2 months when patient-triggered symptom sECG upload is configured to 'On Demand'.

Note: At a maximum shelf-storage time of 18 months (UBD), nominal longevity is reduced by approximately 6 months.

The remaining minimum service time is provided in Table 6 below.

Monitor Status	ERI to EOS in Months
ERI	2

Table 6: Remaining minimum service time


All service intervals, including the above-cited nominal injectable cardiac monitor longevity, are based on considerations that include the battery discharge behavior and the hybrid circuit properties including current consumption and replacement indicator.

This page left intentionally blank

Chapter 11: Insertion/Removal

11.1 Opening the Sterile Container


The BIOMONITOR IV is preloaded in the insertion tool and is packaged with an incision tool in a single container sterilized with ethylene oxide.



Peel off the sealing paper of the outer container as indicated by the arrow.

11.2 Insertion

The FIT OneStep Tool allows an “injection-like” insertion of the implant using a single tool. It is used for forming the device tunnel and subsequent subcutaneous delivery of the BIOMONITOR IV implant. The BIOMONITOR IV implant is provided preloaded into the blue tunneling end of FIT OneStep tool, which has a rigid clam-shell design. There is a small window over the BIOMONITOR IV implant to allow the physician to see the implant in the FIT OneStep tool. The incision tool and FIT OneStep tool are intended for single use. See Figure 54.



Chapter 11

BIOMONITOR IV Technical Manual

Figure 54: Incision and FIT OneStep tools

BIOMONITOR IV has been developed to be inserted in a close-fitting subcutaneous tunnel, preferably in or around the left side of the chest. Recommended locations are those areas close to the heart where the implant will be exposed to minimal movement from body positional changes or from arm movement. Suitable implant locations are shown below in Figure 55. In position A, a location between the suprasternal notch and the left nipple is shown. Position B shows an implant location of approximately 45° with respect to the midline. The choice of placement location is to be decided by the physician, on the basis of individual patient anatomy and comfort, as well as cosmetic considerations. The insertion process consists of four (4) intuitive steps: Incision, Tunneling, Unlocking and Retraction, see Table 7.








Figure 55: Two recommended positions for the placement of BIOMONITOR IV

	<p>Step 1 Local anesthetic agent is injected at the selected anatomical position, both along the incision line, and along the length of the planned tunnel. After an appropriate delay to allow the local anesthetic agent to take effect, the incision tool is used to make an incision through the skin. The physician is advised to consider the patient's anatomy when making the incision.</p>
	<p>Step 2 The FIT OneStep Tool with the preloaded BIOMONITOR IV implant is then advanced within a sub-dermal plane until the skin reaches the insertion stopping point, to create a tunnel for the BIOMONITOR IV implant.</p>
	<p>Step 3 Once the tunneling part of the tool is fully inserted, the knob at the proximal end of the handle is turned counter-clockwise to the unlocked icon.</p>
	<p>Step 4 Holding the outer white portion of the handle stationary against the incision, retract the blue inner portion by pulling back and away from the white portion. The implant will remain in place within the tunnel.</p> <p>ICM insertion has been associated with a small risk of device migration and loss through the incision. To help promote healing and device integrity, closure of the incision, in addition to skin dressing, should be considered for at risk patients. The protection of the wound from environmental influences finalizes the insertion procedure for BIOMONITOR IV.</p>

11.3 Removal

Removed cardiac monitors or accessories may not be reused. Removed cardiac monitors can be sent either to the local BIOTRONIK representative or the BIOTRONIK home office for expert disposal. If possible, the removed cardiac monitor should be cleaned with a sodium-hyperchlorite solution of at least 1% chlorine and, thereafter, washed with water prior to shipping.

The injectable cardiac monitor should be removed before the cremation of a deceased patient.

CAUTION

Device Incineration – Never incinerate a removed cardiac monitor. Be sure the insertable cardiac monitor is explanted before a patient who has died is cremated.

Removed Devices – Return all removed devices to BIOTRONIK.

11.4 Reasons to Remove an Injectable Cardiac Monitor

An injectable cardiac monitor may be removed emergently or at a physician’s discretion at any time subsequent to an implant procedure. Reasons for removal include, but are not limited to: patient death; loss of sensing; inability to program/interrogate the inserted monitor; infection, ERI/ EOS (normal or premature); system upgrade; physician preference for another injectable cardiac monitor model; and/or other reason(s) which may or may not be known to the injectable cardiac monitor manufacturer. Complications related to other portions of the injectable cardiac monitor system (i.e., patient) may also result in injectable cardiac monitor removal.

Table 8 summarizes some of the more common reasons for injectable cardiac monitor removal.

Source	Cause	Possible Effect
Battery	Premature depletion or other cause(s) resulting in excessive battery current drain.	Inability to program/interrogate; sensing difficulty.
Circuitry	Electrical parameter changes due to shorts, opens, or component parametric drift Electromagnetic Interference (EMI) from large power tools, industrial equipment, electrocautery, defibrillation, radiation therapy, RF ablation therapy, etc.	Reversion to “Elective Replacement” or electrical reset parameters; inability to program/interrogate; other damage to circuit components resulting in permanent or temporary parameter changes.
Patient	Normal medical complication	Infection
	Body rejection phenomena	Fluid accumulation; migration; erosion.
	Physician preference	Upgrade to an implantable cardiac pacemaker or implantable cardioverter defibrillator.

Table 8: Common reasons to remove an injectable cardiac monitor

Chapter 12: Remote Assistant III

12.1 General Information on the Remote Assistant III

The Remote Assistant III is an accessory for the BIOMONITOR IV that allows patients to manually trigger the recording of an sECG.

Place the Remote Assistant over the inserted BIOMONITOR IV and press the trigger key on the Remote Assistant III to send a signal. If the BIOMONITOR IV successfully receives the signal, it will record and store an ECG. The Remote Assistant III indicates a successful or unsuccessful recording using the behavior of the signal transmission LED before it automatically turns off.

The Remote Assistant III is powered by two batteries that are non-replaceable due to the design of the device.

12.2 Remote Assistant III Functional Testing

The Remote Assistant III is supplied for immediate use. However, check the Remote Assistant III for visible damage before use. A manual function test may also be performed:

- Hold the Remote Assistant III and press the trigger key. The signal transmission LED briefly lights up yellow.
- Then observe the behavior of the signal transmission LED.


If the LED lights up continuously green or flashes yellow and turns off after a short time, the Remote Assistant III is undamaged, and you can use it immediately.

If the LED exhibits other or no behavior, do not use the device. Contact BIOTRONIK (1-800-547-0394) for further instructions.

12.3 Getting to Know the Remote Assistant III

Item	Description
A	Trigger Key
B	Signal Transmission LED
C	Battery Indicator LED

Figure 56: Remote Assistant III layout



12.4 Triggering a Manual Recording

1. The patient holds the Remote Assistant III over his or her chest as close as possible to the location where the BIOMONITOR IV was implanted. Make sure that the Remote Assistant III lies with its backside flat on the chest without a gap or finger in between.



2. The trigger key is pressed and the Remote Assistant III is kept over the implanted cardiac monitor for at least three seconds.




The device emits an acoustic signal, and the signal transmission LED lights up continuously yellow for approximately three seconds.

3. If the recording has been successfully triggered in the BIOMONITOR IV, the signal transmission LED lights up continuously green for a maximum of 30 seconds before the device automatically turns off. If the signal transmission LED flashes yellow, no recording was triggered in the BIOMONITOR IV. In this case, re-position the Remote Assistant III and try again.


Note: If the BIOMONITOR IV is programmed to 'On Demand' transmission of Symptom episodes, the patient is required to keep the CardioMessenger in close proximity (6ft) for 3h after triggering a recording with Remote Assistant III.

12.5 Battery LED Indicator Explained

Battery indicator LED: 

LED behavior	Explanation
LED is not lit	The battery in the Remote Assistant III has sufficient capacity.
LED flashes yellow	The battery's capacity is only sufficient for another approx. 20 trigger attempts.
LED is continuously lit yellow	The device is defective or the battery's capacity is depleted.

12.6 Signal Transmission LED Explained

Signal transmission LED: 

LED behavior	Explanation
LED is continuously lit yellow	The trigger key has been pressed.
LED flashes yellow	Recording in the BIOMONITOR IV was not successfully triggered.
LED is continuously lit green	Recording in the BIOMONITOR IV was triggered successfully.

This page left intentionally blank

Chapter 13: Technical Data

13.1 Parameters

13.1.1 Atrial Fibrillation

OFF, **ON**

AF Sensitivity

Low, **Medium**, High

Bigeminy Rejection

OFF, **Standard**, Aggressive

RR variability

6%, 9%, **12%**, 15%, 18%

Confirmation Time

0.5, 1, 2, 3, 4, 5, **6**, 10, 20, 30

Recording

OFF, **ON**

Transmission

OFF, **ON**, Detect Only

13.1.2 Tachycardia

OFF, **ON**

Tachy Limit

100....(10).... **180**, 190, 200

Tachy Count

8, 12, **16**, 20, 24, 32, 48

Recording

OFF, **ON**

Transmission

OFF, **ON**

13.1.3 Bradycardia

OFF, **ON**

Brady zone limit

30....(5)....**40**....(5)....80 bpm

Brady duration

5, **10**....(5)....30 seconds

Chapter 13

BIOMONITOR IV Technical Manual

Recording

OFF, **ON**

Transmission

OFF, **ON**

13.1.4 Sudden Rate Drop (SRD)

OFF, ON

SRD rate decrease (%)

20, 30 40, **50**, 60, 70

SRD Sensitivity

Low, **Medium**, High

Recording

OFF, ON

Transmission

OFF, ON

13.1.5 Pause Duration

OFF, **ON**

Pause duration

2, **3**....(1)....10 seconds

Recording

OFF, **ON**

Transmission

OFF, **ON**

13.1.6 Symptom

OFF, **ON**

Recording

OFF, **ON**

Transmission

OFF, **Nightly**, On Demand

13.1.7 Resting Rate Period

Start resting period (hh:mm)

00:00, 01:00, **02:00**....(01:00)....23:00

Resting period duration (hours)

1, 2, 3, **4**...(1)...12

13.1.8 Home Monitoring

OFF, **ON**

Time of transmission (hh:mm)

Std.: 00:00....(00:30)....23:30

13.2 Materials in Contact with Human Tissue

Device Coating: Silicone

Electrodes Coating: Iridium

13.3 Electrical Data/Battery

NOTE: At 37° C

Parameter	BIOMONITOR IV
Power source	Li-CFx
Battery voltage at BOS	3.0 V

13.4 Mechanical Data

Model	Size	Mass	Volume
BIOMONITOR IV	77.5 x 8.6 x 4.6 mm Rigid portion 47.5 x 8.3 x 4.3 mm	4.0 g	1.9 cc Rigid portion 1.7 cc

This page left intentionally blank

Chapter 14: Order Information

Injectable Cardiac Monitor Type	Order Number
BIOMONITOR IV	471155

FCC Statement:

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with part 95 of the FCC Rules. This device may not interfere with stations operating in the 400.150-406.000 MHz band in the Meteorological Aids, Meteorological Satellite, and Earth Exploration Satellite Services and must accept any interference received, including interference that may cause undesired operation.

This device will be registered with the Federal Communications Commission under the following number: FCC ID: QRI-BM2610P2

Transceiver technical data:

- RF telemetry (Home Monitoring)
 - Operating frequency range: 402 ... 405 MHz
 - Max. transmission power: < 25 μ W (-16 dBm)
- PGH-Telemetry
 - Operating frequency range: 9 ... 90 kHz
 - Max. field strength: < 30 dB μ A/m (at 10 m)

Chapter 14

BIOMONITOR IV Technical Manual

BIOTRONIK, Inc.
6024 Jean Road
Lake Oswego, OR 97035-5369
(800) 547-0394 (24-hour)
(800) 291-0470 (fax)
www.biotronik.com

M4234-A 02/23
© 2023 BIOTRONIK, Inc. All rights reserved.
MN101 3/7/23

