

SD2 Ultrasonic Pocket Doppler User Manual

About this Manual
P/N: 01.54.459245
MPN: 01.54.459245010
Release Date: Mar. 2022
© Copyright EDAN INSTRUMENTS, INC. 2022. All rights reserved.

Statement

This manual will help you understand the operation and maintenance of the product better. It is reminded that the product shall be used strictly complying with this manual. User's operation failing to comply with this manual may result in malfunction or accident for which Edan Instruments, Inc. (hereinafter called EDAN) cannot be held liable.

EDAN owns the copyrights of this manual. Without prior written consent of EDAN, any materials contained in this manual shall not be photocopied, reproduced or translated into other languages. Materials protected by the copyright law, including but not limited to confidential information such as technical information and patent information are contained in this manual, the user shall not disclose such information to any irrelevant third party. The user shall understand that nothing in this manual grants him, expressly or implicitly, any right or license to use any of the intellectual properties of EDAN. EDAN holds the rights to modify, update, and ultimately explain this manual.

Responsibility of the Manufacturer

EDAN only considers itself responsible for any effect on safety, reliability and performance of the equipment if:
Assembly operations, extensions, re-adjustments, modifications or repairs are carried out by persons authorized by EDAN, and
The electrical installation of the relevant room complies with national standards, and
The instrument is used in accordance with the instructions for use.

EDAN will make available on request circuit diagrams, component part lists, descriptions, calibration instructions, or other information that will assist service personnel to repair those parts of the equipment that are designated by EDAN as repairable by service personnel.

Product Information

Product Name: Ultrasonic Pocket Doppler

Model: SD2, SD2 Pro, SD2 Plus, SD2 Lite, SD2 Basic

Terms Used in this Manual

This guide is designed to give key concepts on safety precautions.

WARNING

WARNING label advises against certain actions or situations that could result in personal injury or death.

CAUTION

CAUTION label advises against actions or situations that could damage equipment, produce inaccurate data, or invalidate a procedure.

NOTE

A **NOTE** provides useful information regarding a function or a procedure.

● Safety Precautions

CAUTION

Federal (U.S.) Law restricts this device to sale by or on the order of a physician

NOTE:

This user manual is written to cover the maximum configuration. Therefore, your model may or may not have some of the parameters and functions described, depending on what you have ordered.

This unit is internally powered equipment, and it is an IEC/EN 60601-1 Type BF applied part. Type BF protection means that the connection between the equipment and personnel complies with permitted leakage currents and dielectric strength of IEC/EN 60601-1.

WARNING and **CAUTION** messages must be observed. To avoid the possibility of injury, observe the following precautions during the operation of the device.

WARNING

- It is to be used by health care professionals and patients on the order of a physician.
- Before the SD2 is prescribed for home use, the user (patient) must be instructed/trained in proper use of the equipment.
- Home fetal heart rate detection has not been shown to prevent the onset of preterm labor nor will it prevent the occurrence of preterm birth.
- The Doppler is a tool to aid the user in hospitals, clinics or at home and should not be used in place of normal fetal detection. It is not intended for treatment or use during labor and delivery.
- Placement of the ultrasound transducer on the abdomen is critical to obtaining the fetal heart beat as opposed to maternal heart beat or other abdominal noise. The user should be trained in proper placement techniques either through acceptable Ob/Gyn training and individual state accreditation, or as being prescribed by such a trained clinician and trained in device placement.
- This device is not explosion-proof and cannot be used in the presence of flammable anesthetics.
- Magnetic and electrical fields are capable of interfering with the proper performance of the device. For this reason, make sure that all external devices operated in the vicinity of this device comply with the relevant EMC requirements. X-ray equipment and magnetic resonance imaging (MRI) devices can emit high levels of electromagnetic radiation.
- We recommend that exposure to ultrasound should be kept as low as reasonably achievable. This is considered to be good practice and should be observed at all times.
- Do not use the device with HF surgical equipment and do not use it in an MRI environment.
- The device is not protected against defibrillation.
- SHOCK HAZARD** - Do not attempt to replace batteries or connect or disconnect a power cord with wet hands.
- The device should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.
- Please ensure the light of SpO2 sensor is fully covered by the finger pulp. Please check the application site every 0.5 to 1 hour. If the skin quality changes, please change another site.
- The medical electrical equipment needs to be installed and put into service according to the EMC information provided in this user manual.
- Portable and mobile RF communications equipment can affect medical electrical equipment; refer to section Recommended Separation Distances.
- Do not service or maintain the device or any accessory which is in use with a patient.

CAUTION

- Refer servicing to qualified personnel.
- Keep the device in a clean environment and avoid vibration during storage.
- Do not sterilize the Doppler.
- Only use SpO2 sensor approved by the manufacturer. Using other SpO2 sensors may compromise device functionality.
- Electromagnetic Interference - Ensure that the environment in which the device is operated is not subject to any source of strong electromagnetic emissions, such as radio transmitters, mobile telephones, etc.
- The device and accessories are to be disposed of according to local regulations after their useful lives. Alternatively, they can be returned to the dealer or the manufacturer for recycling or proper disposal. Batteries are hazardous waste. Refer to IEC61429 for standard disposal when necessary. Do NOT dispose them together with household garbage.
- The loss of pulse signal may occur when the patient has poor peripheral perfusion, and the screen will display "----". When there's measurement beyond range, invalid measurement or no measurement value, it will display "----".

● Introduction

Intended Use/Indications for Use

The SD2 detects fetal heart rate (FHR), maternal oxygen saturation (SpO2) and pulse rate (PR) for pregnant women. The product is only intended for use in hospitals, clinics or at home by health care professionals with relevant expertise and pregnant women who have been trained (by studying instructional videos or product manuals) or have received professional guidance.

Features

You may not have all of these features, depending on the model purchased.

Function	Model	SD2	SD2 Pro	SD2 Plus	SD2 Lite	SD2 Basic
Curve Display		√	√	√	○	○
FHR		√	√	√	√	√
PR		√	√	√	√	√
SpO2		○	○	○	○	○
Bluetooth +APP		○	√	√	√	○
Variable Frequency		○	√	√	√	○
AA Alkaline Batteries		○	√	√	√	○
Charging Function		○	○	○	○	○

NOTE: √ = Standard ○ = Optional × = Not Available


Interface:

There are three display modes. Double click the power touch key to switch the display modes. The displayed results of SpO2 and PR are updated every second.



Numeric Display Interface **FHR Curve Display Interface** **SpO2 Wave Display Interface**

Appearance (Above pictures are just for reference)



Item	Description
1	Fetal heart icon Indicates fetal heart beat and flickers to the fetal heart beat.
2	Perfusion bar Indicates pulse beat
2	FHR numeric 142 Displays fetal heart rate within the range from 30 bpm to 240 bpm.
	PR numeric 75 Displays pulse rate within the range from 30 bpm to 240 bpm
	SpO2 numeric 99 Displays SpO2 within the range from 70% to 100%

	Volume numeric	-5-	Volume numeric is displayed in the center of the screen, and ranges from level 0 to 7.
3	Battery indicator		Displayed on the bottom left corner of the screen and there are five degrees. When the battery indicator is empty and keeps flickering, it indicates that the battery level is extremely low and battery needs to be replaced or charged.
4	Fetal heart signal quality indicator		Displayed on the bottom right corner of the screen and there are three levels, representing poor, acceptable and good signal quality.
5	Volume increase touch key		Touch the key to decrease volume.
6	Volume decrease touch key		Touch the key to increase volume.
7	Ultrasound frequency	2MHz, 3MHz	Indicates the current ultrasound frequency. Touch the volume increase and decrease touch keys simultaneously to switch the frequency. 2 MHz: used for more than 12-week gestation 3 MHz: used for more than 9-week gestation
8	Power touch key		Touch this key for a little while to switch it on or off;

● Basic Operation

NOTE:

To ensure that the Doppler works properly, please read this chapter and Chapter **Safety Precautions** before operation; follow the steps when connecting all the components.

Opening the Package and Checking

Open the package; take out the Doppler and accessories carefully. Keep the package for possible future transportation or storage. Check the components according to the packing list.

- Check for any mechanical damage.
- Check all the cables and accessories.

Installing the Battery

- Remove the battery compartment cover.
- Insert the battery into the compartment carefully. Ensure its anode and cathode terminals are aligned with the anode and cathode marks on the compartment.
- Install the compartment cover.

Recommended Battery Type:

AA Alkaline battery (AA, LR6, 1.5V); AA Rechargeable NI-MH battery (AA, NI-MH, 1.2V)

Charging the NI-MH Batteries

When battery level is low, take the batteries out from the main unit, and charge them with a NI-MH battery charger. The specifications of the provided battery charger are as below. You can also purchase a battery charger that meets the following specifications:
Output: 1.2V-500mA*4AA 500mA*2AAA Battery Capacity: 500mAh-3000mAh

NOTE:

- Rechargeable NI-MH battery is configured, please fully charge the battery after each use of the Doppler to ensure sufficient power during subsequent use.
- Please charge the battery after each transportation or storage.

WARNING

- Turn off the Doppler before removing or replacing the battery.
- Replace alkaline batteries with those of identical specifications provided by the manufacturer or purchased locally. See **Product Specifications** for details about battery specifications.
- If the batteries have been inserted incorrectly, the Doppler will not function or it will be damaged.
- Keep battery away from children. If swallowed, consult a physician at once.
- Do not disassemble or short-circuit batteries.
- Do not recharge alkaline batteries.
- Do not dispose of batteries in fire or water.
- Do not mix with used or other battery type (such as alkaline with carbon zinc).
- Do not solder the batteries directly.
- Do not over-discharge batteries.
- Do not allow metal objects to contact the battery terminals.
- To install or remove batteries, follow the equipment manufacturer's instructions.
- Remove the battery and store it at a cool and dry environment if the Doppler is not used for a long time. Do not keep batteries at temperature of 45°C or above, or at humidity of 75% or above. For rechargeable NI-MH batteries, please charge them at least once every 3 months.
- If rechargeable batteries are used, charge them fully before initial use by using the method introduced in this manual.
- The rechargeable NI-MH batteries should be charged by using the qualified dedicated adapter that conforms to requirements of IEC 60601-1.
- If the time that the Doppler using rechargeable NI-MH batteries becomes much shorter than usual, the battery life is at an end. Replace the batteries with new ones.

Switching On

Touch the power touch key for about 1 second when the Doppler is off, and the Doppler will display the switching on interface before switching to display the test interface.

Switching Off

Touch the power touch key for about 1 second when the Doppler is on, and the Doppler will be switched off. If the Doppler is not in operation or no signal is received for 1 minute, the Doppler will be switched off automatically.

● FHR Detection

SD2 is designed to detect the fetal heartbeat from the 10th week of gestation.

Before applying the Doppler to inspect FHR, you should always check whether the Doppler is in good condition and whether there is evident damage that might affect patient's safety and the device's function. If evident damage is found, stop using it at once and replace it with a good one.

Procedures to Detect FHR:

- Have the patient lie on her back.
- Apply appropriate amount of coupling gel to the ultrasonic transducer head of the Doppler and switch on the Doppler.
- Palpate the patient's abdomen gently to confirm the fetus's position.
- Place the Doppler on the patient's abdomen, and move it around the fetus's position or tilt it until a clear and rhythmic heart sound is heard and FHR numeric is stably displayed.


The coupling gel should not exceed this limit.

This area can be immersed in coupling gel.

- Note:**
- Do not mistake the maternal heart rate for fetal heart rate. The fetal pulse should be different from the maternal pulse, which can be measured at the wrist or neck.
 - Do not wear gloves to touch the keys. If there's water and coupling gel on the fingers, please clean them first or the touching effect will be influenced.

How to Find the Best FH Signal:

- The easiest way: take the position the doctor last detected for FHR as a reference and move the Doppler around the position slowly until the best FH signal is found.
- The fetal heart position may change as the fetus moves inside the uterus. You can confirm the fetal position first according to the position of the uterus fundus in different gestational weeks.
 - At the end of the 12-week gestation, the uterus fundus is about 2-3 fingers' breadth (about 3-4.5 cm) above the symphysis pubis.
 - At the end of the 16-week gestation, the uterus fundus is in the middle between the navel and the symphysis pubis.
 - At the end of the 20-week gestation, the uterus fundus is about 1 finger's breadth (about 1.5 cm) below the navel.
 - At the end of the 24-week gestation, the uterus fundus is about 1 finger's breadth (about 1.5 cm) above the navel.
 - At the end of the 28-week gestation, the uterus fundus is about 3 fingers' breadth (about 4.5 cm) above the navel.
 - At the end of the 32-week gestation, the uterus fundus is in the middle between the navel and the xiphisternum.
 - At the end of the 36-week gestation, the uterus fundus is about 2 fingers' breadth (about 3 cm) below the xiphisternum.



The clearest and loudest fetal heart sound is generally obtained when the Doppler is placed on the fetus's back. Fetal movement is usually the movement of fetal limbs. So, if frequent fetal movement occurs at the right side of the abdomen, the fetus's back is probably at the left side and vice versa. You can find the fetus's back according to fetal movement's position.

If the fetus is in cephalic delivery position, the fetal heart is either on the right side or on the left side below the navel; if the fetus is in breech delivery position, the fetal heart is either on the right side or on the left side above the navel.

Steps to Find Fetal Heart:

Have the patient lie on back and relax >> confirm fetal position by hand >> apply coupling gel to the Doppler >> place the Doppler on patient's abdomen and start looking for the fetal heart >> the fetal heart is found when the Doppler gives out a continuing thumping sound "boom-boom-boom".

CAUTION

- The Doppler's degree of protection against harmful ingress of water is IP22. Do not immerse it in water.
- The Doppler is delicate and sensitive. Please handle it with care and try to avoid dropping on to the ground or any hard surfaces. Any damage caused by dropping is not covered by the warranty.
- Keep the coupling gel away from children. If swallowed, consult a physician at once.

Note:

- The best quality of fetal heart signal is obtained only when the Doppler is placed in the best detection position.
- Do not place the Doppler near positions where placental sound or umbilical blood flow sound is loud.
- If the fetus is in the cephalic position and the mother is supine, the clearest heart sound will normally be found on the midline below the navel. During detection, the pregnant woman's prolonged lying in the supine position should be avoided to reduce the possibility of supine hypotension. Putting a pillow or cushion under the patient's head or feet can be of help.
- It is not possible to obtain accurate FHR unless a clear fetal heart signal is detected. If the calculated FHR is not in accordance with the beat of the fetal heart sound, the fetal heart sound auscultation result shall prevail.
- Whether the fetal heart sound can be obtained or not is related to the skill and qualification of the operator.
- If the fetal heart sound cannot be found for a long time, it is necessary make use of other equipment.
- When applied to the patient, the Doppler may warm slightly (less than 2°C (35.6°F) above ambient temperature). When NOT applied, the Doppler may slightly (less than 5°C (41°F) above ambient temperature).

● SOV Prompt

When detecting fetal heart rate, there are possibilities that maternal HR signal is mistaken for FHR signal. The SOV function can reduce these possibilities. Detect FHR and PR at the same time. When FHR and PR signals overlap, SD2 will issue a sound. At this time, examine the patient or reposition the probe until

fetal heart signal is found.

● SpO2 Detection

SpO2 Plethysmogram measurement is employed to determine the oxygen saturation of hemoglobin in the arterial blood. The SpO2 numeric shows the percentage of hemoglobin molecules which have combined with oxygen molecules to form oxyhemoglobin. The SpO2/PLETH parameter can also provide a pulse rate signal and a plethysmogram wave.

NOTE:

- The device is calibrated to display functional oxygen saturation.
- The monitor does not have specific SpO2 calibration baselines.

3 SpO2 waveform is not proportional to the pulse volume.

4 A Functional tester or simulator cannot be used to assess the SpO2 accuracy. However, it can be used to test the accuracy of a particular calibration curve duplicated by the device, and it turns out the calibration curve meets the accuracy.

Inaccurate measurements can be caused but not limited by:

- Incorrect device application
- high levels of ambient light sources, such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight
- failure to cover the sensor with fingerin high levels of ambient light conditions
- dysfunctional hemoglobins
- low peripheral perfusion
- excessive or violent patient movement
- venous pulsations
- intravascular dyes, such as indocyanine green or methylene blue
- externally applied coloring agents (nail polish, dye, pigmented cream)
- defibrillation
- placement of the device on an extremity with a blood pressure cuff, arterial catheter, or intravascular line
- electromagnetic interference
- Low perfusion

Loss of pulse signal can occur for the following reasons:

- a blood pressure cuff is inflated on the same extremity as the one with the device attached
- there is arterial occlusion proximal to the device
- low peripheral perfusion

● Assessing the Validity of a SpO2 Reading

You can check the quality of the pleth wave and the stability of the SpO2 values to assess whether the sensor functions properly and whether the SpO2 readings are valid. Always use these two indications simultaneously to assess the validity of a SpO2 reading.

NOTE:

- The SpO2 accuracy has been validated in human studies against arterial blood sample reference measured with a CO-oximeter. Pulse oximeter measurements are statistically distributed, only about two-thirds of the measurements can be expected to fall within the specified accuracy compared to CO-oximeter measurements. The volunteer population in the studies composed of local healthy men and women from age 19 to 37, with variations of skin pigmentations. The SpO2 accuracy is as follows: ±2% for 90%-100% and ±4% for 70%-90%.
- The pulse rate accuracy is obtained by comparison to the pulse rate generated with an arterial oxygen simulator (also an electronic pulse simulator).

Generally, the quality of the SpO2 pleth wave reflects the quality of the light signals obtained by the sensor. A wave of poor quality manifests a decline of the signal validity. On the other hand, the stability of the SpO2 values also reflects the signal quality. Different from varying SpO2 readings caused by physiological factors, unstable SpO2 readings are resulted from the sensor's receiving signals with interference. The problems mentioned above may be caused by patient movement, wrong sensor placement or sensor malfunction. To obtain valid SpO2 readings, try to limit patient movement, check the placement of the sensor, measure another site or replace the sensor.

● SpO2 Detection

- Hold the device with one hand and press the thumb on the SpO2 sensor continuously.
- Movement is not recommended during measurement.
- When the signals are stable, read corresponding data from screen.

NOTE: Avoid external light sources such as radiated rays or ultrared rays.

After Detection

- Switch off the Doppler.
- Wipe the remaining gel off the patient and the probe with a clean soft cloth or tissue.

● Mobile Application Software (APP)

SD2 can connect to mobile phones with its Bluetooth function (optional). The SD2 APP has both Android and iOS versions.

iOS APP operating environment:	Android APP operating environment:
A) hardware environment Processor: dual-core Apple A6 RAM: ≥2GB	A) hardware environment CPU: frequency≥1.0GHz RAM: ≥4GB
B) software environment: iOS 13.5.1 and above operating system	B) software environment: Android 8.0 and above operating system
C) network environment: support Bluetooth	C) network environment: support Bluetooth

Scan either of the following QR codes to download the SD2APP, and install and run it as prompted.

Note:

- Your mobile phone may prohibit the installation of "applications from unknown sources". Enter **Settings** to allow the installation first.
- For normal functioning of the APP, please give the APP function-related permissions.
- For how to use the APP, read the instructions in the **About** sub-interface under the **Settings** interface of the APP.

WARNING

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:
1) this device may not cause harmful interference, and
2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
-Reorient or relocate the receiving antenna.
-Increase the separation between the equipment and receiver.
-Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
-Consult the dealer or an experienced radio/TV technician for help.

Radiation Exposure Statement

This device complies with RF radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

NOTE:

- This equipment (SD2) has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.
- Any changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

● Maintenance and Cleaning

Maintenance

The overall check of the Doppler, including safety check and function check, should be performed by qualified personnel every 12 months, and each time after service. And safety check must include current leakage test and insulation test. Besides the above requirements, comply with local regulations on maintenance and measurement.

The accuracy of FHR is determined by the Doppler and cannot be adjusted by user. If you have doubt concerning the accuracy of FHR, verify it with other methods such as using a stethoscope, or contact local distributor or the manufacturer for help.

The Doppler is fragile and must be handled with care. Wipe the remaining gel off the Doppler after each use. These measures can help prolong the Doppler's life. Replace the accessories such as the battery according to use. If any of the accessories are damaged, refer to chapter **Ordering Information** for details and order new ones. Please check the label for the date of manufacture, the service life is 5 years (The service life is limited to the Doppler, not including the replaceable accessories. The only replaceable accessory of SD2 is battery. The frequency of usage is 8 hours/day).

Cleaning

Before cleaning, switch off the Doppler. Clean the exterior surface of the Doppler with a soft, clean cloth dampened with ethanol (75%) or mild neutral detergent thoroughly until no visible contaminants remain. After cleaning, wipe off the cleaning solution with a fresh cloth or paper towel dampened with tap water until no visible cleaning agent remains. Dry the monitor in a ventilated and cool place.

CAUTION

- Do not use strong solvent, such as acetone.
- Never use an abrasive such as steel wool or metal polish.
- Do not remain any solution on the surface after cleaning.

Disinfection

Before disinfection, switch off the Doppler. Wipe the exterior surface of the Doppler with a soft, clean cloth dampened with ethanol (75%) or mild neutral detergent. Wipe off the disinfectant solution with a dry cloth after disinfection if necessary. Dry the Doppler for at least 30 minutes in a ventilated and cool place.

Size:	Length*Width* Height: (48±1) mmx (39±1) mmx (145±3) mm
Weight:	<200g(including batteries)
LCD:	Size: 1.3 inch, 23.4 mm × 23.4 mm Resolution: 240×240(Pixels)
Coupling Gel:	pH: 5.5-8.0 Acoustic Impedance: 1.5x10 ⁶ Pa.s/m - 1.7x10 ⁶ Pa.s/m (35°C/95°F)

Environment	
Working:	Temperature: +5°C ~ +40°C (+41°F ~ +104°F) Humidity: 15% RH ~ 95% RH (non-condensing) Atmospheric Pressure: 70 kPa ~ 106 kPa
Transport and Storage:	Temperature: -25°C ~ +70°C (-13°F ~ +158°F) Humidity: 15% RH ~ 95% RH (non-condensing) Atmospheric Pressure: 70 kPa ~ 106 kPa

Note: The time required for the Doppler to warm from the minimum storage temperature between uses until it is ready for intended use is at least 2 hours; the time required for the Doppler to cool from the maximum storage temperature between uses until it is ready for intended use is at least 2 hours

Performance Specifications

FHR Performance):	(Essential) FHR Measuring Range: 30 bpm ~ 240 bpm Accuracy: ±2 bpm Note: FHR measurement result may not be accurate if the equipment is measuring beyond its measuring range.
FHR Resolution:	1bpm
Audio Output:	Output Power: 2w Background noise: <45dBA
Overall Sensitivity:	>90dB
Auto Power-off:	Power off when the Doppler receives no signal or operation for 1 minute.
Bluetooth:	Transmission Range (Without Obstacles) >: 5m (Indoor range depends on the building's structure and material.)
Intelligent Denoising	Noise generated in static state < 100mV
Shock Resistance	Withstands a 1.3 m drop at least 6 times to concrete surface with possible cosmetic damage only
Ultrasound:	Frequency: (3.0±10%) MHz; (2.0±10%) MHz
	p < 1 MPa
	I _{ob} < 10 mW/cm ²
	I _{spta} < 100 mW/cm ²
	I _{sata} < 10 mW/cm ²
	I _{sppa} < 3.190 W/cm ²
	I _{spta} < 94 mW/cm ²
	Effective Radiating Area: 490mm ² ± 15%
	Ultrasound output power: <49mW
	Working Mode: pulse wave
SpO ₂ :	Pulse wave repetition rate: 5kHz (±5%)
	Pulse wave duration: 75µs (±5%)
	Measurement range: 70%~100%
	Resolution: 1%
SpO ₂ Sensor:	Accuracy: ±2%
	Red light: 660nm
	Infrared light: 905 nm
	Wave length: about 660nm and 900nm
PR:	Emitted light energy: <15mW
	Measurement range: 30~240bpm
	Resolution: 1bpm
	Accuracy: ±2bpm or ±3%, whichever is greater

Battery Specifications	
Specification:	Two AA alkaline batteries Two AA rechargeable NI-MH batteries
Working Duration:	≥24h

Bluetooth Specifications	
Modulation:	GFSK π/4-DQPSK 8DPSK
Frequency:	2400-2483.5MHz
Tolerance Frequency:	± 2ppm
RF output power:	≤ 20dBm (EIRP)
Occupied Channel Bandwidth:	≤ 2MHz
Transmitter Unwanted Emissions:	≤ -30dBm

Low Output Summary Table (For systems whose global maximum value does not exceed 1.0) System: SD2 Ultrasonic Pocket Doppler					
Model (MHz)	I _{spta,3} (mW/cm ²)	TI Type	TI Value	MI	I _{sppa,3} (W/cm ²)
SD2 CD3.0	3.25	TIS	0.032	0.01	0.009
		TIB	0.11		
SD2 CD2.0	4.21	TIS	0.028	0.013	0.012
		TIB	0.086		

Ordering Information			
CAUTION Only the parts supplied by the manufacturers should be used with the Doppler.			
Parts	Part Number	Parts	Part Number
SD2 Doppler	02.06.263378	SD2 Lite Doppler	02.06.263382
SD2 Pro Doppler	02.06.263379	Rechargeable NI-MH battery	21.21.064180
SD2 Plus Doppler	02.06.263380	AA Alkaline Battery	01.21.064086
SD2 Basic Doppler	02.06.263381	Zipper bag	01.56.466428

Ultrasound Intensity and Safety

Ultrasound in Medicine
The use of diagnostic ultrasound has proved to be a valuable tool in medical practice. Given its known benefits for non-invasive investigations and medical diagnosis, including investigation of the human fetus, the question of clinical safety with regards to ultrasound intensity arises.

There is no easy answer to the question of safety surrounding the use of diagnostic ultrasound equipment. Application of the ALARA (As Low As Reasonably Achievable) principle serves as a rule-of-thumb that will help you to get reasonable results with the lowest possible track record of over 25 years of use and no confirmed biological effects on patients or instrument operators, the benefits of the prudent use of diagnostic ultrasound clearly outweigh any risks.

Ultrasound Safety and the ALARA Principle
Ultrasound waves dissipate energy in the form of heat and can therefore cause tissue warming. Although this effect is extremely low with Doppler, it is important to know how to control and limit patient exposure. Major governing bodies in ultrasound have issued statements to the effect that there are no known adverse effects from the use of diagnostic ultrasound, however, exposure levels should always be limited to As Low As Reasonably Achievable (the ALARA principle).

Explanation of MI/TI
MI (Mechanical Index)
Cavitations will be generated when ultrasound wave passes through and contacts tissues, resulting in instantaneous local overheating. This phenomenon is determined by acoustic pressure, spectrum, focus, transmission mode, and factors such as states and properties of the tissue and boundary. This mechanical bioeffect is a threshold phenomenon that occurs when a certain level of ultrasound output is exceeded. The threshold is related to the type of tissue. Although no confirmed adverse mechanical effects on patients or mammals caused by exposure at intensities typical of present diagnostic ultrasound instruments have ever been reported, the threshold for cavitation is still undetermined. Generally speaking, the higher the acoustic pressure, the greater the potential for mechanical bioeffects; the lower the acoustic frequency, the greater the potential for mechanical bioeffects. The AIUM and NEMA formulate mechanical index (MI) in order to indicate the potential for mechanical effects. The MI is defined as the ratio of the peak-rarefactional acoustic pressure (should be calculated by tissue acoustic attenuation coefficient 0.3 dB/cm/MHz) to the acoustic frequency.

$$MI = \frac{P_r}{f \times CMI}$$

$$CMI = 1 \text{ (MPa / MHz)}$$

TI (Thermal Index)
Heating of tissues is caused by absorption of ultrasound when the ultrasound energy is applied. The temperature rise is determined by the acoustic intensity, exposed area and thermophysical properties of the tissue.

In order to indicate the potential for temperature rise caused by thermal effects, the AIUM and NEMA formulate thermal index (TI). It is defined as the ratio of the total acoustic power to the acoustic power required to raise the tissue temperature by 1°C (1.8°F).

According to different thermophysical properties of the tissue, TI is divided into three kinds: TIS, TIB and TIC.

TIS (Soft Tissue Thermal Index): It provides an estimate of potential temperature rise in soft or similar tissues.

TIB (Bone Thermal Index): It provides an estimate of potential temperature rise when the ultrasound beam passes through soft tissue and a focal region is in the immediate vicinity of bone.

TIC (Cranial Bone Thermal Index): It provides an estimate of potential temperature rise in the cranial bones or superficial bones.

Measurement Uncertainties
The uncertainties in the measurements were predominantly systematic in origin; the random uncertainties were negligible in comparison. The overall systematic uncertainties were determined as follows:

1. Hydrophone Sensitivity: ± 12 percent for intensity, ± 6 percent for pressure. Based on the hydrophone calibration report by ONDA. The uncertainty was determined within ±1 dB in frequency range 1-15 MHz.
2. Digitizer: ±0.3 percent for intensity, ± 0.15 percent for pressure.
3. Based on the stated accuracy of the 8-bit resolution of the Agilent DSO6012 Digital Oscilloscope and the signal-to-noise ratio of the measurement.
4. Temperature: ±2.4 percent for intensity uncertainty, ±1.2 percent for pressure uncertainty. Based on the temperature variation of the water bath of ±1°C (1.8°F).
5. Spatial Averaging: ± 3.5 percent for intensity, ± 1.75 percent for pressure.
6. Non-linear Distortion: N/A.

No effects of nonlinear propagation were observed.

Since all the above error sources are independent, they may be added on an RMS basis, giving a total uncertainty of ± 12.73 percent for all intensity values reported, ± 6.37 percent for all the pressure values, ± 12.6 percent for the Mechanical Index, uncertainty of ±12.73% percent for power, ±0.15 percent for center frequency, ±6.87% for the MI.

Prudent Use Statement
Although no confirmed bioeffects on patients caused by exposure from present diagnostic ultrasound equipment have ever been reported, the potential exists that such bioeffects may be identified in the future. Therefore, the ultrasound should be used prudently. High levels of acoustic output and long exposure time should be avoided while acquiring necessary clinical information.

Acoustic Output Reporting Table for Track 1 Acoustic output reporting table for IEC60601-2-37 (IEC60601-2-37, Edition 2.1, 2015-0, table 201.103)					
Transducer Model: SD2, Operating Mode: PW mode					
Index label	MI		TIS		TIC
	At Surface	Below Surface	At surface	Below Surface	
Maximum index value	0.013	0.028	0.086		N/A
Index component value		N/A	0.028	NA	0.086

Acoustic Parameters	$\frac{P_r \times CMI}{f}$ (MPa)	zMI	0.019			
	P (mW)			8.62		N/A
	P _{1x1} (mW)			N/A		
	z _s (cm)				5.12	
	z _b (cm)					4.97
	zMI (cm)	0.50				
	zPIL _a (cm)	0.50				
	f _{awf} (MHz)	2.00	3.00			3.00
	p _{rr} (Hz)	4999.00				
	s _{rr} (Hz)	N/A				
npps	N/A					

Other information	$\frac{P_r \times CMI}{f}$ (MPa)	0.0090					
	I _{spta} (W/cm ²)	3.25					
	I _{spta} at zPIL or zSII (mW/cm ²)	3.72					
	p _r at zPIL (MPa)	0.021					
	Operating control conditions	Fixed	N/A	Fixed	N/A	Fixed	N/A
		Fixed	N/A	Fixed	N/A	Fixed	N/A
		2.00	N/A	3.00	N/A	3.00	N/A

	$\frac{P_r \times CMI}{f}$ at zPIL _a (W/cm ²)	0.012				
	I _{spta} at zPIL _a or zSII _a (mW/cm ²)	4.21				
	I _{spta} at zPIL or zSII (mW/cm ²)	4.67				
	p _r at zPIL (MPa)	0.020				
Operating control conditions	Fixed	N/A	Fixed	N/A	Fixed	N/A
	Fixed	N/A	Fixed	N/A	Fixed	N/A
	2.00	N/A	2.00	N/A	2.00	N/A

Transducer Model: SD2, Operating Mode: PW mode						
Index label	MI		TIS		TIB	
	At Surface	Below Surface	At surface	Below Surface	TIC	
Maximum index value	0.011	0.032	0.11		N/A	
Index component value		N/A	0.032	NA	0.11	

Acoustic Parameters	$\frac{P_r \times CMI}{f}$ (MPa)	zMI	0.019			
	P (mW)			10.96		N/A
	P _{1x1} (mW)			N/A		
	z _s (cm)				4.65	
	z _b (cm)					4.55
	zMI (cm)	0.50				
	zPIL _a (cm)	0.50				
	f _{awf} (MHz)	3.00	3.00			3.00
	p _{rr} (Hz)	4999.00				
	s _{rr} (Hz)	N/A				
npps	N/A					

Standard Parameter Equal Contrast List			
IEC60601-2-37 Standard Parameters			
Parameter	Note	Parameter	Note
P_r	Attenuated Peak-rarefactional Acoustic Pressure	f_{awf}	Center Frequency, Acoustic Working Frequency
P_r	Peak-rarefactional Acoustic Pressure	X	-12dB Output Beam Dimensions
P	Output Power	Y	
z_s	Depth for Soft Tissue Thermal Index	t_d	Pulse Duration
$P_r(Z_s)$	Attenuated Output Power	p_{rr}	Pulse Repetition Frequency (Pulse Repetition Rate)
$I_{ta}(Z_s)$	Attenuated Temporal-average Intensity	d_{eq}	Equivalent Beam Diameter
z_{sp}	Break-point Depth	$I_{p,0}$ at max MI	Attenuated Pulse-average Intensity at the point of Maximum MI
z_b	Depth for Bone Thermal Index	A_{opt}	-12dB Output Beam Area
$I_{p,0}$	Attenuated Pulse-intensity Integral	MI	Mechanical Index
I_p	Pulse-intensity Integral	TIS	Soft Tissue Thermal Index
$d_{eq}(Z_s)$	Equivalent Beam Diameter at the point of z_{sp}	TIB	Bone Thermal Index
		TIC	Cranial-bone Thermal Index

EMC Information

Electromagnetic Emissions

Guidance and manufacturer's declaration – electromagnetic emission
The SD2 Ultrasonic Pocket Doppler is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The SD2 Ultrasonic Pocket Doppler uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The SD2 Ultrasonic Pocket Doppler is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC/EN61000-3-2	Not applicable	
Voltage fluctuations /flicker emissions IEC/EN61000-3-3	Not applicable	

Electromagnetic Immunity

Guidance and manufacturer's declaration –electromagnetic immunity
The SD2 Ultrasonic Pocket Doppler is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV ±15 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient/Burst IEC/EN61000-4-4	±2kV for power supplies ±1kV for input/output lines	Not applicable	Not applicable
Surge IEC/EN61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not applicable	Not applicable
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC/EN61000-4-11	<5%UT(>95% dip in UT) for 0.5cycle 40%UT(60%dip in UT) for 5 cycles 70%UT(30%dip in UT) for 25 cycles <5%UT(>95% dip in UT) for 5s	Not applicable	Not applicable
Power frequency (50Hz/60Hz) magnetic field IEC61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Electromagnetic Immunity

Guidance and manufacturer's declaration –electromagnetic immunity
The SD2 Ultrasonic Pocket Doppler is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC61000-4-6	3 V _{rms} 150 kHz ~ 80 MHz 6V _{rms} in ISM bands between 0.15 MHz and 80 MHz	Not Applicable	Portable and mobile RF communications equipment should be used no closer to any part of the SD2 Ultrasonic Pocket Doppler, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Radiated RF IEC61000-4-3	10V/m 80 MHz ~ 2.7 GHz	10 V/m 80 MHz to 2.7 GHz	$d = 0.35\sqrt{P}$ 80 MHz to 800 MHz $d = 0.7\sqrt{P}$ 800 MHz to 2.7 GHz $d = 6\sqrt{P/E}$ at RF wireless communications equipment bands (Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the SD2 Ultrasonic Pocket Doppler, including cables specified by the manufacturer). Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol:
Immunity to Proximity fields from RF wireless communications equipment IEC 61000-4-3	See Table - Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment for test level	See Table-Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment for test level	

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the SD2 Ultrasonic Pocket Doppler is used

exceeds the applicable RF compliance level above, the SD2 Ultrasonic Pocket Doppler should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the SD2 Ultrasonic Pocket Doppler.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

^c The ISM (industrial, scientific and medical) bands between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. The amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz and 50.0 MHz to 54.0 MHz.

Table-Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment

Test Frequency (MHz)	Brand ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum Power(W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380-390	TETRA 400	Pulse modulation ^{b)} 18Hz	1.8	0.3	27
450	430-470	GMRS 460, FRS 460	FM ^{c)} ±5 kHz deviation 1kHz sine	2	0.3	28
710	704-787	LTE Brand 13, 17	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9
745						
780						
810						
870	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	0.3	28
1720						
1845	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	2	0.3	28
1970						
2450	2400-2570	Bluetooth; WLAN 802.11 b/g/n; RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0.3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9
5500						
5785						

Note: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1m. The 1 m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.
b) The carrier shall be modulated using a 50% duty cycle square wave signal.
c) As an alternative FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case

Recommended Separation Distances

Recommended separation distances between portable and mobile RF communications equipment and the SD2 Ultrasonic Pocket Doppler				
The SD2 Ultrasonic Pocket Doppler is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the SD2 Ultrasonic Pocket Doppler can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the SD2 Ultrasonic Pocket Doppler as recommended below, according to the maximum output power of the communications equipment.				
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)			
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 0.35\sqrt{P}$	800 MHz to 2.7 GHz $d = 0.7\sqrt{P}$	
0.01	/	0.035	0.07	
0.1	/	0.11	0.22	
1	/	0.35	0.7	
10	/	1.11	2.21	
100	/	3.5	7	













For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Troubleshooting

Problem	Possible Cause	Solution
Fail to power on, or shut down shortly after switching on	Battery level is very low. Battery is not installed properly. Fail to switch on the Doppler as instructed. The Doppler has malfunctions.	Replace the batteries or charge the rechargeable NI-MH batteries Re-install the battery. Touch the On/Off touch key for a while to power on the Doppler. Contact service personnel.
Loudspeaker does not work.	Sound volume has been turned down to the lowest level. If the Doppler is configured with Bluetooth, fetal heart sound can be played by mobile phone. The Doppler has malfunctions.	Adjust sound volume to appropriate level. Set to play fetal heart sound by mobile phone or the Doppler on the APP. Contact service personnel.
F		

● **Definition of Symbols**

No.	Symbol	Definition	No.	Symbol	Definition
1		Disposal method	9		General symbol for recovery/recyclable
2		Operating instructions	10		Refer to User Manual (Background: Blue; Symbol: White)
3		Caution	11		MR Unsafe—Keep away from magnetic resonance imaging (MRI) equipment
4		Type BF applied part	12		Non-ionizing electromagnetic radiation
5	P/N	Part Number	13	IP22	IP22 Protected against solid foreign objects of 12.5 mm ϕ and greater. Protection against vertically falling water drops when ENCLOSURE tilted up to 15°
6		Serial Number (Start with H on battery compartment cover)	14	Rx Only	Caution: Federal (U.S.) law restricts this device to sale by or on the order of a physician.
7		Date Manufacture of	15		No SpO2 Alarms
8		Manufacturer	16		