PressureWire™ Aeris™

Guidewire Mounted Sensor

INSTRUCTIONS FOR USE

en: Instructions for use

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hr: Upute za upotrebu

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sl: Navodila za uporabo

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

en: Instructions for use

Content

PressureWire™ Aeris™ packages include:

Sealed non-sterile outer box
Sealed sterile pouch
Packaging coil
PressureWire Aeris
Torque device
Transmitter
Quick Guide
Instructions for use

Store at room temperature in a dry and dark place.

- | | |
|---|---|
| <ol style="list-style-type: none"> 1. Guidewire, refer to label for length 2. Flexible length, 31 cm 3. Shaft 4. Connector 5. Radiopaque tip, 3 cm 6. Sensor element 7. Hydrophilic coating 8. PTFE coating | <ol style="list-style-type: none"> 9. Hydrophobic Polyimide insulators between contact members 10. Torque device 11. PressureWire Aeris Transmitter 12. Function indicators 13. Switch 14. Close Cap 15. Open Cap 16. Cap |
|---|---|

Description

PressureWire is a 0.014" guidewire with an integrated sensor element at the tip to enable measurements of physiological parameters. The guidewire is uniquely paired with a specific transmitter for PressureWire Aeris. PressureWire is available in different lengths. Please refer to the label for information about PressureWire length and thermo compatibility.

Intended use

PressureWire is designed to fit inside a percutaneous catheter for the purpose of directing the catheter through a vessel. The signal output from the sensor is used for calculations and presentation of any physiological parameters, functions or indices based on pressure or temperature, e.g. Fractional Flow Reserve (FFR).

Indications for use

PressureWire is indicated to direct a catheter through a blood vessel and to measure physiological parameters in the heart and in the coronary and peripheral blood vessels.

Contraindications

PressureWire is contraindicated for use in the cerebral vasculature.

Warnings

- No modification of this device is allowed.
- PressureWire is supplied sterile. Discard PressureWire if the pouch is opened or damaged, compromising the sterile barrier. PressureWire is designed for single use only and shall not be reused or resterilized. Use of a non-sterile or resterilized PressureWire may result in, but is not limited to:
 - Local and/or systemic infection
 - Mechanical damage
 - Inaccurate readings
- Observe all PressureWire movements. Whenever PressureWire is moved or torqued, the tip movement should be examined under fluoroscopy. Never push, withdraw or torque PressureWire if it meets resistance or without observing corresponding movement of the tip; otherwise vessel/ventricle trauma may occur.
- Torquing or excessive manipulation of PressureWire in a sharp bend, against resistance or repeated attempts to cross a total vessel occlusion may:
 - Damage and/or fracture PressureWire
 - Cause dissection or perforation of blood vessels
 - Cause vessel spasm
- When introducing PressureWire flush the catheter and administer anticoagulation as for a standard catheterization procedure or clotting may occur.
- Do not use PressureWire in the ventricles if the patient has a prosthetic mechanical or biological valve. It may result in damage to both the prosthesis and PressureWire, which may cause injury or death.
- Use of PressureWire in conjunction with interventional devices with a short monorail may result in folded or fractured PressureWire.
- High frequency surgical device must not be used on a patient at the same time as PressureWire.

Precautions

- PressureWire is a delicate instrument and should be handled carefully.
- Make sure that the transmitter is kept dry, to ensure accurate pressure and/or temperature readings. Inaccurate reading may necessitate device replacement.
- Do not use PressureWire in conjunction with atherectomy catheters. It may damage PressureWire.
- Do not withdraw or manipulate PressureWire in a sharp-edged object. It may result in abrasion of PressureWire coating.
- The accuracy of the diagnostic information is affected by, but not limited to:
 - Failure to achieve maximum coronary and myocardial hyperaemia.
 - Interventional devices, such as balloon catheters, which are positioned as they affect the blood flow.
- PressureWire readings may be affected by defibrillation. Rezero PressureWire after defibrillation use.
- Do not measure pressure when the sensor element of PressureWire is in sharp bends or in contact with atrial or ventricular walls. It might result in pressure artefacts.
- Avoid using PressureWire together with another wire, for so called jailed wire technique, due to difficulty in wire withdrawal causing wire entrapment.

Adverse events

Potential complications which may be encountered during all catheterization procedures include but are not limited to: vessel dissection or occlusion, perforation, embolus, spasm, local and/or systemic infection, pneumothorax, congestive heart failure, myocardial infarction, hypotension, chest pain, renal insufficiency, serious arrhythmias or death.

Direction for use

Use PressureWire in conjunction with 6F (2 mm diameter) guiding catheter. Refer to instructions supplied with any interventional devices to be used in conjunction with PressureWire for their intended uses, contraindications and potential

complications. Note that use of PressureWire in conjunction with uncoated diagnostic catheters may reduce the performance of PressureWire.

CAUTION: If PressureWire gets stuck or is damaged during the procedure, immediately disconnect the transmitter and then replace the damaged PressureWire, including the transmitter, with a new one. Damage may be but is not limited to kinks, bends, sheared polymer coating, no pressure signal or inaccurate pressure signal. Incorrect use may result in vessel/ventricle damage, arrhythmia induction, electric shock, inaccurate pressure signals or inaccurate torque response.

CAUTION: When PressureWire is disconnected from the transmitter during procedure make sure the male connector does not come into contact with conductive surfaces, thus avoiding unintentional connection with other equipment or electric shock.

Preparations

1. Open PressureWire packaging using aseptic technique.
2. Make sure that the packaging coil is flat on the table.
3. Fill the packaging coil with saline solution through the flush port. The saline shall cover the pressure sensor element, approximately 20 cc.

Zeroing PressureWire

4. Activate the receiving instrument according to manufacturer's instruction.
5. Zero PressureWire by turning on the Aeris Transmitter.

CAUTION: Make sure to keep PressureWire stable at all time and avoid any excessive movement during zeroing. The transmitter's function indicator shall display a stable green light when it is successfully zeroed. Please refer to table 5 for further information.

6. Carefully remove PressureWire from the packaging coil.
7. The shapeable guidewire tip may be carefully shaped using standard tip shaping practices.

CAUTION: Do not use a shaping tool with a sharp edge.

Diagnostic procedure

8. Introduce PressureWire into the patient's blood vessels according to standard clinical practice.

CAUTION: Avoid turning off the transmitter during the procedure since PressureWire is automatically zeroed when the transmitter is turned on and PressureWire requires zeroing outside of the body.

9. Use the torque device to slowly steer PressureWire. Position PressureWire sensor element just outside the guiding catheter opening.
10. To ensure accurate pressure reading: Place the aortic pressure transducer at the same height as the patient's heart. Pull back the insertion tool and tighten the hemostatic valve. Make sure that there is no remaining contrast fluid in the catheter.
11. Verify that the pressures from the guiding catheter and PressureWire are equal.
12. If the pressures differ: Equalize the pressure signals according to manufacturer's instructions.

NOTE: If pressure difference before equalization is greater than 20 mmHg or limit according to equipment specification whichever is lowest it is recommended to remove the wire, rezero and repeat from step 8.

13. **Perform measurements.** Use the torque device to slowly steer PressureWire to the desired position and perform pressure measurement according to standard clinical procedures.

NOTE: If it is difficult to reach the desired position, PressureWire may be disconnected for better handling.

14. If applicable, considering the results from the measurements, proceed with interventional procedures.

Interventional procedure

15. Disconnect the transmitter.
16. Remove the torque device by loosening the screw and gently withdraw it.
17. Carefully advance the interventional device over PressureWire guidewire and proceed according to manufacturer's instructions.
18. Put the torque device back on a suitable position on the guidewire.
19. Carefully wipe and dry the guidewire.
20. Take care to gently and fully insert the guidewire into the transmitter and secure the cap.
21. Treat the lesion according to manufacturer's instruction.

How to pause procedure

Perform the following procedure if PressureWire is removed from the patient's blood vessels during any diagnostic or interventional procedure, in order to be reintroduced:

1. Carefully withdraw PressureWire. Make sure that the distal flexible part of the guidewire is kept wet.
2. Carefully wipe and dry the guidewire before it is reconnected to the transmitter.
3. Repeat the procedures from step 8 and forward.

How to verify pressure measurement

1. Carefully withdraw PressureWire until the sensor element is positioned just outside the catheter opening. Verify that the pressures registered by the guiding catheter and PressureWire are equal.

How to finish procedure

1. Carefully withdraw PressureWire. Turn off the Transmitter.

Handling and disposal

1. Handle and dispose PressureWire as biohazardous material and in accordance with medical and applicable local, state and federal law and regulations.
2. The battery can be removed by opening the end cap using e.g. forceps in the slits under the label.

CAUTION: The transmitter is powered by a standard Li-ion battery. Li-ion batteries can contain a small concentration of 1,2-Dimethoxyethane. Do not open the battery or disassemble it. Do not expose to fire or high temperatures (>60°C). Dispose battery according to standard local practice for Li-ion batteries.

Technical specifications

Table 1. Pressure measurement

Operating pressure:	-30 to +300 mmHg
Accuracy ¹⁾ :	≤ 2 mmHg or ± 2%
Temperature sensitivity:	0.3 % of pressure change per °C
Zero drift:	< 7 mmHg/h
Operating time:	3 hours

Table 2. Temperature measurement

Temperature range:	34 – 42°C
--------------------	-----------

Table 3. Radio signal specification

Frequency range:	2.4000-2.4835 GHz (ISM-band)
Type:	Frequency Hopping Spread Spectrum
Radiated power:	1 mW peak, 70 µW average (EIRP)
Range:	Up to 5 m. NOTE: Radio range is reduced by objects and walls, keep transmitter and receiver in line of sight wherever possible.
Delay time:	<28 ms

-
- 1) The contribution to the error in average pressure from any external disturbance source will be ≤ 2 mmHg i.e. when a disturbance source is applied in accordance with IEC 60601-1-2, the average pressure may change up to 2 mmHg. The effect of electromagnetic noise on pressure signal is less than 3 mmHg peak to peak.

Table 4. Environmental conditions

Operating, ambient temperature	15°C to 35°C (59°F to 95°F)) max ± 5°C temperature fluctuation
Operating, relative air humidity	30-75%
Ambient pressure	425-850 mmHg
Transport temperature	- 25°C to +60°C (-13°F to +140°F)
Transport relative humidity	10-95%

Table 5. Transmitter function indicators light

Function indicators	Explanation
Fast blinking green light	Zeroing (initializing)
Green steady light	Transmitter on, normal operation
Alternating green/yellow blinking light	Stability check failed (pulsatile pressure detected). Rezero PressureWire at stable zero pressure.
Yellow steady light	PressureWire disconnected from transmitter. Reinsert PressureWire firmly into transmitter.
Yellow blinking light.	Error. Possible causes can be that the battery level is too low for operation or other internal error. Exchange PressureWire.

Table 6. Platforms compatible with PressureWire Aeris

PressureWire™ Receiver
Ilumien™
QUANTIEN™
For other compatible platforms, check with your sales representative.

Compliance with Regulatory Requirements

Guidance and manufacturer's declaration – electromagnetic emissions

PressureWire system is intended for use in the electromagnetic environment specified below. The customer or the user of PressureWire system should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The PressureWire system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The PressureWire system is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies building for domestic purposes.
Harmonic emissions IEC 61000-3-2	N/A No mains connection	
Voltage fluctuations /flicker emissions IEC 61000-3-3	N/A No mains connection	

Guidance and manufacturer's declaration – electromagnetic immunity

The PressureWire system is intended for use in the electromagnetic environment specified below. The customer or the user of the PressureWire system should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/ output lines	N/A 2) 3)	N/A
Surge IEC 61000-4-5	±1 kV line(s) to line (s) ±2 kV line(s) to earth	N/A 2)	N/A
Voltage dips, short interruptions and voltage variations to power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	N/A 2)	N/A
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
N/A	N/A	N/A	Portable and mobile RF communication equipment should be used no closer to any part of the PressureWire system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A 2)	Recommended separation distance N/A

-
- 2) No mains power input
3) No signal or interconnecting cables

Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	<p>Recommended separation distance $d = 1,2\sqrt{P}$ 80 MHz to 800 MHz $d = 2,3\sqrt{P}$ 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey⁴⁾, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
------------------------------	----------------------------	-------	---

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

NOTE: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Recommended separation distances between portable and mobile RF communications equipment and PressureWire system

The PressureWire system is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the PressureWire system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the PressureWire system as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance to frequency of transmitter (m)		
	150 kHz to 80 MHz* $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0,01	N/A	0,12	0,23
0,1	N/A	0,38	0,74
1	N/A	1,2	2,3
10	N/A	3,8	7,4

- 4) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the PressureWire system is used exceeds the applicable RF compliance level above, the PressureWire system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the PressureWire system.

Rated maximum output power of transmitter (W)	Separation distance to frequency of transmitter (m)		
	150 kHz to 80 MHz* $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
100	N/A	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance (*d*) in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

*Not applicable since the transmitter is not connected to mains.



NOTE: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.



NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

WARRANTY DISCLAIMER

Although PressureWire, hereafter referred to as "product", has been manufactured under carefully controlled conditions, St Jude Medical, hereafter called SJM, has no control over the conditions under which the product is used. SJM, therefore disclaims all warranties, both expressed and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. SJM shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind SJM to any representation or warranty with respect to the product. The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

SYMBOLS WITH EXPLANATIONS

	Guidewire Mounted Sensor.
	Agile Tip.
	Instructions for Use.
	Guidewire mounted pressure sensor, disconnectable transmitter, torque device.
	Caution (Attention, consult accompanying documents).
	Follow instructions for use.
	Use-by date.
	Manufacturer.
	Manufacturing facility.
	European Representative.
	Date of manufacture.

REF	Catalogue number, C12058.
MODEL	Catalogue number, C12058.
	Diameter.
	Length.
LOT	Batch code.
QTY	Quantity.
	Keep away from sunlight.
	Keep dry.
	Temperature limit.
	Humidity limitation.
	Defibrillation-proof Type CF equipment.
	Equipment includes RF transmitter.
	Electrostatic sensitive device.
	Do not use if package is damaged.
	For single use only. Do not reuse.
	Do not resterilize.
STERILE EO	Sterilized using Ethylene Oxide.
	Electronic waste bin. Disposal according to EU's Battery Directive 2006/66/EC.
	Electronic waste bin. Disposal according to EU's WEEE Directive 2002/96/EC.
CE	PressureWire Aeris complies with the Medical Device Directive 93/42/EEC (NB 0086), the Radio and Telecommunication Terminal Equipment Directive 1999/5/EC and with the Restriction of Hazardous Substances Directive 2011/65/EU.
R_x Only	USA only: Federal law restricts this device to sale by or on the order of a Physician.

FCC ID:
IC:8466A-2

FCC identifier for the transmitter.
IC identifier for transmitter.

This device complies with part 15 of the FCC Rules and RSS-210 of the IC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

cETLus valid for Canada and US.
Conforms to ANSI/AAMI Std. ES60601-1
Certified to CAN/CSA Std. C22.2 No. 60601-1

C-tick compliance with AS/NZS 4268.

Japanese RF identifier for the transmitter.


KCC identifier for the transmitter.

Canada 310

This device complies with RSS-310 of Industry Canada. Operation is subject to the condition that this device does not cause harmful interference.

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www.sjm.com

Australian sponsor:

St. Jude Medical Australia Pty Limited

17 Orion Rd.

Lane Cove NSW 2066 AUSTRALIA

+61 2 9936 1200

ST. JUDE MEDICAL[®]

PressureWire Aeries is designed, developed and manufactured by St. Jude Medical. Unless otherwise noted, [™] indicates that the name is a trademark or licensed to, St. Jude Medical, Inc. or one of its subsidiaries. PressureWire, ST. JUDE MEDICAL, the nine-square symbol are trademarks and service marks of St. Jude Medical, Inc. and its related companies.

©2014 St. Jude Medical. All rights reserved.

R_x
Only

CE
0086