**inControl** H·E·B™

# **Premium** **Blood Pressure Monitor**

with Bluetooth® Connectivity

Instruction Booklet for model# BP3NQ1-4BHEB

## Important Product and Safety Information

---

**Follow Instructions for Use. This document provides important product operation and safety information regarding this Blood Pressure Monitor. Please read this document thoroughly before using the device and keep for future reference.**

This blood pressure monitor is an automatic digital blood pressure measuring device for use by adults on the upper arm at home or in the doctor's office. It enables a very fast and reliable measurement of the systolic and diastolic blood pressure, as well as the pulse rate, by way of the oscillometric method. This device detects the appearance of irregular heartbeats during measurement and provides a warning signal when the irregular heartbeat is detected.

### **This blood pressure monitor IS intended to be used:**

- For self-measurement/monitoring of blood pressure and pulse in adults; whereas the person being measured may be the user/operator of the device.
- Within a home healthcare environment.
- With a cuff located upon the user's upper arm; ½ inch above the elbow over the artery as indicated on the cuff.

### **This blood pressure monitor IS NOT intended for use with:**

- Pregnancy
- Preeclampsia
- Children under 12
- Neonatal patients

**WARNING – Potentially hazardous situation that if not avoided may result in serious injury or death.**

- Self-measuring means monitoring, not diagnosis or treatment. Unusual values must always be discussed with your doctor. Under no circumstances should you independently alter the dosages of any drugs prescribed by your physician.
- Consult your physician before using this device if any of the following or similar conditions are present: arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, arterial sclerosis, poor perfusion, diabetes, age, pregnancy, preeclampsia, renal diseases. Motion during measurement, including trembling or shivering may affect the measurement.
- Though not for use with children under 12, ensure that any children around this device are supervised; some parts are small enough to be swallowed and any

present tubes or cables may provide a risk of strangulation.

- Ensure the cuff tubing is not kinked during use as harmful injury may occur due to the effect of blood flow interference caused by high pressure in the cuff not releasing.

### **DO NOT:**

- ✗ Use this device if you think it is damaged or anything appears unusual (e.g. sporadic operation, open packaging upon purchase).
- ✗ Use the displayed pulse for checking the frequency of heart pacemakers as this device is not suitable for this action.
- ✗ Open/modify this device; inaccuracy and/or harmful injury may result.
- ✗ Conduct many frequent measurements as this may result in harmful injury due to blood flow interference.
- ✗ Place the Cuff over a wound as this may cause further injury.
- ✗ Place and pressurize the Cuff over/near any present intravascular access or therapy, or arteriovenous shunt, as this may cause blood flow interference and result in harmful injury.
- ✗ Place and pressurize the Cuff over a limb near the side of a mastectomy as this may cause harmful injury.
- ✗ Use the Blood Pressure Monitor on a limb simultaneously with other medical equipment on the same limb due to possible interference with such medical equipment.
- ✗ Maintain pressure in the Cuff applied to the limb for a prolonged amount of time. Ensure that circulation in the limb is not impaired by checking circulation if prolonged/repeated exposure to pressure occurs.
- ✗ Service the device when being used or when power is supplied. When the device is not in use and power is removed, cuff assembly, batteries and ac adapter (if applicable) may be replaced by the user with Microlife supplied replacements. No other parts/components are accessible.

 **CAUTION – Potentially hazardous situation that if not avoided may result in minor/moderate injury, property damage, and/or damage to the device**

- This device contains sensitive electronics components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g., mobile telephones, microwave ovens). These can lead to temporary impairment of the measuring accuracy. Move the device to another location if interference is determined.
- When not using the Blood Pressure Monitor for extended lengths of time, remove

the batteries to avoid potential battery leakage and damage to the monitor.

- When replacing the batteries, ensure all the batteries are replaced simultaneously to avoid battery damage and potential damage to the monitor. Microlife USA does not recommend using rechargeable batteries.
- To avoid inaccurate measurements and to lessen any discomfort from Cuff pressure, ensure the Cuff is placed correctly on the limb and fits correctly when snug (not tight), as indicated by markings with the Cuff.
- Consult your physician in cases of frequent irregular heartbeat detections.
- This Risk Indicator feature is provided in order to help you understand your potential blood pressure risk. However, this feature is neither a diagnosis nor a substitute for a medical examination. It is important to consult with your physician to determine your risk.

### **DO NOT:**

- × Drop this device or expose it to strong vibrations; sensitive components may be affected resulting in inaccuracies and/or operational issues.
- × Use the Blood Pressure Monitor outside of its specified operation temperature and humidity rating, or if stored outside of its specified storage temperature and humidity rating. Avoid storage in direct sunlight.
- × Use this device in a moving vehicle; inaccurate measurements may result.
- × Use third party accessories. Only use Microlife authorized accessories, such as cuffs or AC adapters, as those not approved for use with the device may provide inaccurate measurements, injury, and/or damage the device.

### **STANDARDS**

In addition to the standards stated in the Instruction Manual:

- This medical device is compliant with medical device and non-invasive blood pressure monitor standards IEC 60601-1, IEC 60601-1-2, IEC 60601-1-11, and AAMI/ANSI/IEC 80601-2-30, and
- Electromagnetic standards IEC 60601—1-2 along with FCC Part 15, and
- Clinical Testing per standard ISO 81060-2:2013 was conducted on blood pressure device using the same measurement technology.

Please note: According to international standards, your monitor should be checked for accuracy every 2 years.

## TYPE

IP20 Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water. Keep Dry. 

Type BF Applied Part

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

## Expected Life

Monitor 5 Years

Cuff 2 Years

## FCC

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Changes or modifications to the product are not approved by Microlife USA and could void the user's authority to operate the equipment under FCC jurisdiction.

This equipment has been **tested and found** to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment **generates**, uses and can radiate radio frequency energy and, if not **installed and** used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment on and off, the user is encouraged to try to correct the interference by **one or more of the following measures: 1) Reorient or relocate the receiving antenna, 2) Increase the separation between the equipment and receiver, 3) Connect the equipment into an outlet on a circuit different from that to which the receiver is connected, 4) Consult the dealer or an experienced radio/TV technician for help.**

This device complies with FCC radiation exposure limits set forth for an uncontrolled environment. End users must follow the specific operating instructions for satisfying RF exposure compliance.

**Trademark Usage:**

- Apple, the Apple logo, iPad, and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.
- Android and Google Play are both trademarks of Google Inc.
- The Bluetooth® Smart word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks in this Blood Pressure Monitor is under license. Other trademarks and trade names are those of their respective owners.

# Automatic Blood Pressure Monitor - Instruction Manual

## Table of contents

### 1. Introduction

---

- 1.1. Your automatic blood pressure monitor
- 1.2. Gentle+™ Measurement System

### 2. How is blood pressure measured?

---

- 2.1. What is blood pressure, and why does it fluctuate?
- 2.2. Measurement guidelines
- 2.3. What can I do to change my blood pressure?

### 3. Components of your blood pressure monitor

---

### 4. How do I get started?

---

- 4.1. Inserting the batteries
- 4.2. Setting the date and time
- 4.3. Tube connection
- 4.4. Select the user

### 5. How do I take a measurement?

---

- 5.1. Preparing to take a measurement
- 5.2. Common errors
- 5.3. Using the cuff
- 5.4. Taking a measurement
- 5.5. Irregular heartbeat detection
- 5.6. Hypertension risk indicator
- 5.7. Memory: reviewing readings
- 5.8. Stopping a measurement
- 5.9. Battery change indicator
- 5.10. Using the AC adapter

### 6. Bluetooth® functions & application (app)

---

- 6.1. Activating your Bluetooth®
- 6.2. Pairing your device

**7. Error messages/troubleshooting**

---

**8. Care and maintenance**

---

**9. Limited warranty**

---

**10. Technical specifications**

---

**11. How to contact us**

---

## 1. Introduction

---

### 1.1. Your automatic blood pressure monitor

Thank you for purchasing a fully automatic blood pressure monitor. Your monitor is designed to provide fast and reliable digital readings of your pulse, and systolic / diastolic blood pressure using the oscillometric method on your upper arm. It offers clinically proven accuracy and has been designed to be user friendly.

Before using your blood pressure monitor, please read this instruction manual carefully to ensure correct use. If you have additional questions regarding blood pressure measurements please contact your doctor.

### 1.2. Gentle+™ measurement system

Gentle+™, our most advanced system to date, combines clinical accuracy and premium features including:

- A quieter, more comfortable measurement experience that reduces the amount of time your arm is under pressure. Less pressure = more comfort!
- An ergonomic cuff design that is easier to use and fits a wider range of arm sizes (for arms 8.7" -16.5").

## 2. How is blood pressure measured?

---

### 2.1. What is blood pressure, and why does it fluctuate?

Your level of blood pressure is determined in the circulatory center of the brain and adjusts to a variety of situations through feedback from the nervous system. To adjust blood pressure, the strength and speed of the heart (Pulse), as well as the width of circulatory blood vessels is altered. Blood vessel width is controlled by fine muscles in the blood vessel walls.

Your level of arterial blood pressure changes periodically during heart activity: During the “blood ejection” (Systole) the value is highest (systolic blood pressure value). At the end of the heart’s “rest period” (Diastole) pressure is lowest (diastolic blood pressure value).

Blood pressure values must lie within certain normal ranges in order to prevent particular diseases.

### 2.2. Measurement guidelines

Blood pressure is very high if your systolic blood pressure is over 160 mmHg and/or your diastolic pressure is above 100 mmHg, **while at rest**. In this case, please consult your physician immediately. Long-term values at this level endanger your health due to continual damage to the blood vessels in your body.

If your systolic blood pressure values are between 140 mmHg and 159 mmHg and/or the diastolic blood pressure values between 90 mmHg and 99 mmHg, consult your physician. Regular self-checks are necessary.

If you have blood pressure values that are too low, (i.e., systolic values under 105 mmHg and/or diastolic values under 60 mmHg), consult your physician.

Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. You can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, keep a record of values along with time of day and date. Show these values to your physician. **Never use the results of your measurements to independently alter the medication prescribed by your physician.**

## Which values are normal?

The following standards for assessing high blood pressure in adults have been established by the American Heart Association (AHA) and American College of Cardiology (ACC) in 2017.

<b>Category</b>	<b>Systolic</b> (mm Hg)	<b>Diastolic</b> (mm Hg)
Normal	<120	and <80
Elevated	120-129	and <80
<b>Hypertension</b>		
Stage 1 Hypertension	130-139	or 80-89
Stage 2 Hypertension	140-179	or 90-119
Hypertensive Crisis	≥180	or ≥120

## Additional information

- This chart reflects a 2017 update to blood pressure standards.
- If your values are mostly normal under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called “labile hypertension.” Consult your doctor.

## 2.3. What can I do to change my blood pressure?

- a) Consult your doctor.
- b) Increased blood pressure values (various forms of hypertension) are associated with considerable health risks over time. Arterial blood vessels in your body are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can result from arteriosclerosis. Furthermore, the heart will become structurally damaged with increased blood pressure values.
- c) There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organ malfunctions. Please consult your doctor for information about the possible origins of your

own increased blood pressure values.

d) There are measures which you can take to reduce and even prevent high blood pressure. These measures must be permanent lifestyle changes.

### **A) Eating habits**

- Strive for a normal weight corresponding to your age. See your doctor for your ideal weight.
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

### **B) Previous illnesses**

- Consistently follow all medical instructions for treating illness such as:
- Diabetes (Diabetes mellitus or sugar diabetes)
- Fat metabolism disorder
- Gout

### **C) Habits**

- Give up smoking completely.
- Drink only moderate amounts of alcohol.
- Restrict your caffeine consumption (coffee, tea, chocolate).



### **D) Your physical condition**

- **After a medical examination, and with your doctor's approval and direction; exercise.**
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your exercise routine. Speak with your doctor regarding the type and extent of exercise that is appropriate for you.

### 3. Components of your blood pressure monitor

---

#### a) Monitor

#### b) Wide range cuff

The cuff is for arm circumference 8.7-16.5" (included) or 22-42 cm.

Additional cuffs are available for purchase by calling 1-800-568-4147.

#### **Please Note:**

Arm circumference should be measured with a measuring tape in the middle of the relaxed upper arm. Do not force cuff connection into the opening. Make sure the cuff connection is not pushed into the AC adapter port. If the cuff is too small, call 1-800-568-4147 for additional information. You may prefer to use one of our wrist type monitors.

## 4. How do I get started?

---

### 4.1. Inserting the batteries

Insert batteries immediately after unpacking the unit. The battery compartment is located on the bottom of the device.

- a) Remove cover.
- b) Insert the batteries (4 x size AA 1.5 V), observing the indicated polarity.
- c) If a battery warning appears in the display, the batteries are discharged and must be replaced.

Low Battery Indicator

### Attention!

- After the low battery indicator appears, the device won't function until the batteries have been replaced.
- Please use AA Long-Life or Alkaline 1.5 V batteries and replace them all at the same time.
- We do not recommend using rechargeable batteries.
- If the blood pressure monitor is not used for long periods, remove the batteries from the device.


## 4.2. Setting the date and time

1. After the new batteries are activated, the year number flashes in the display. You can advance the year by pressing the "M" button. To confirm and then set the month, press the time button.
2. You can now set the month using the "M" button. To confirm and then set the day, press the time button.
3. Please follow the instructions above to set the day, hour and minutes.
4. Once you have set the last minute and pressed the time button, the date and time are set and the time is displayed.
5. If you want to change the date and time, press and hold the time button down for approx. 3 seconds or until the year starts flashing. Now you can enter the new values as described above.

 Memory Button

 Time/Date Button

### 4.3. Tube connection

Insert the black cuff tube connector into the side of the monitor indicated by the drawing of a cuff.

### 4.4. Select the user

This blood pressure monitor is designed to store 99 measurements for each of two users. Before taking a measurement, be certain that the correct user has been selected.

- Press and release the TIME button to toggle between users ("1," or "2").
- Release the button. The user icon in the top left-hand corner of the LCD screen flashes.
- Press the power button to make your selection.

## 5. How do I take a measurement?

---

**Please note: You should always be seated before and during measurement.**

### 5.1. Preparing to take a measurement

- Avoid eating and smoking as well as all forms of exertion directly before measurement. These factors influence the measurement result. Find time to relax by sitting in an armchair in a quiet atmosphere for about 5 minutes before taking a measurement.
- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (normally left).
- Always compare measurements taken at the same time of day, since blood pressure changes during the course of the day, as much as 20-40 mmHg.

### 5.2. Common errors

**Note: Comparable blood pressure measurements always require the same conditions! Conditions should always be quiet.**


- All efforts by the user to support the arm can increase blood pressure. Make sure you are in a comfortable, relaxed position and do not flex any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower or higher than the heart, an erroneously high or low blood pressure will be measured! Each 15cm (6 inch) difference in height between your heart and the cuff results in a measurement error of 10 mmHg!
- Cuffs that are too narrow or too short result in false measurement values. Selecting the correct cuff is extremely important. Cuff size is dependent upon the circumference of the arm (measured in the center). The permissible range is printed on the cuff. If this is not suitable for your use, please call 1-800-568-4147.
- A loose cuff or a sideways protruding air pocket causes false measurement values.

- With repeated measurements, blood accumulates in the arm, which can lead to false results.
- Consecutive blood pressure measurements should be repeated after a 1 minute pause in order to allow the accumulated blood to flow away.




### 5.3. Using the cuff

- Pass the end of the cuff through the flat metal ring so that a loop is formed. The hook and loop material must be facing outward. (Ignore this step if the cuff has already been prepared).
- Slide the cuff over the left upper arm so that the tube is closer to your lower arm.
- Position the cuff on the arm as illustrated. Make certain that the lower edge of the cuff, with the red edge, lies approximately  $\frac{1}{2}$ " (1 to 2 cm) above the elbow and that the tube is closer to the inner side of the arm.

Brachial Artery

**Important:** The red strip on the edge of the cuff (Artery Mark) must lie over the artery which runs down the inner side of the arm.

**TIP:** Align red artery mark to pinky finger.

- Tighten the cuff by pulling the end and close the cuff by affixing the hook and loop material.
- There should be little free space between the arm and the cuff. You should be able to fit 2 fingers between your arm and the cuff. Clothing must not restrict the arm. Any piece of clothing which does must be removed. Cuffs that don't fit properly result in false measurement values. Measure your arm circumference if you are not sure of proper fit.

- f) Lay your arm on a table (palm upward) so the cuff is at the same height as your heart. Make sure the tube is not kinked.
- g) Remain seated quietly for at least 5 minutes before you begin the measurement.

**Comment:**

If it is not possible to fit the cuff to your left arm, it can also be placed on your right arm. However, all measurements should be made using the same arm.

Cuff on  
right arm

## 5.4. Taking a measurement

After the cuff has been appropriately positioned on the arm and connected to the blood pressure monitor, the measurement can begin. Remain still during your measurement, do not flex muscles.

- a) Remain seated quietly for 5 minutes before you begin the measurement.
- b) Press the POWER button to begin. After pressing the POWER button the pump begins to inflate the cuff. The measurement is performed while cuff inflates. The pumping speed can vary, which is completely normal. The increasing cuff pressure will show on the display.
- c) When the device has detected your pulse, the heart symbol in the display begins to blink.
- d) After reaching an individual pressure, the pump stops and the pressure will release. The cuff pressure (large digits) is displayed during the measurement.
- e) When the measurement has been concluded, the air will automatically release from the cuff. The measured systolic and diastolic blood pressure values, and the pulse are displayed.
- f) The measurement results are displayed until you switch the device off. If no button is pressed for 1 minute, the device switches off automatically.

Pumping  
Pressure

Measuring

Systolic

Diastolic


Pulse

Measurement  
complete

## 5.5. Irregular heartbeat detection

This symbol  indicates that certain pulse irregularities were detected during the measurement. In this case, the result may deviate from your normal basal blood pressure—repeat the measurement. In most cases, this is no cause for concern. However, if the symbol appears on a regular basis (e.g. several times a week with measurements taken daily), we advise you to tell your doctor. Please show your doctor the following explanation:

### **Information for the doctor on frequent appearance of the Irregular Heartbeat Symbol**

This instrument is an oscillometric blood pressure monitor device that also analyzes pulse frequency during measurement. The instrument is clinically tested.


If pulse irregularities occur during measurement, the irregular heartbeat symbol is displayed after the measurement. If the symbol appears more frequently (e.g. several times per week on measurements performed daily) or if it suddenly appears more often than usual, we recommend the patient to seek medical advice.

The instrument does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

## 5.6. Hypertension risk indicator


The bars on the left-hand edge of the display show you the range within which the indicated blood pressure value lies. Depending on the height of the bar, the readout value is either within the normal (green), borderline (yellow) or danger (orange, red) range.

The classification is based on standards established by the National Institutes of Health JNC7, 2003..

With each reading, a corresponding black triangle is displayed to the right of the hypertension indicator.

- If the bar is in the green zone, your measurement is "Normal," according to NIH standards.
- If the bar is in the yellow zone, your measurement is "Pre-Hypertension."
- If the bar is in the orange zone, it is "Stage 1 Hypertension."
- If the bar is in the red zone, it is "Stage 2 Hypertension."

## 5.7. Memory: reviewing readings

At the end of a measurement, this monitor automatically stores each result with date and time. This unit stores 99 memories for each of 2 users.

### Viewing the stored values

With the unit off, press the "M" button. The display first shows "A", then shows an average of all measurements stored in the unit. Please note: Measurements for each user are averaged and stored separately. Be certain that you are viewing the measurements for the correct user.


Pressing the "M" button again displays the previous value. To view a particular stored memory, press and hold the "M" button to scroll to that stored reading.

### Memory full

When the memory has stored 99 results for the selected user, a new measured value is stored by overwriting the oldest value.

### Clear all values

If you are sure that you want to permanently remove all stored values, hold down the "M" button (the instrument must have been switched off beforehand) until "CL" appears and then release the button. If you do not want to clear the values, press the START/STOP button. To permanently clear the memory, press the "M" button while "CL" is flashing.

**Individual values cannot be cleared.**

## 5.8. Stopping a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g the patient feels unwell), the power button can be pressed at any time. The device then immediately lowers the cuff pressure automatically.

## 5.9. Battery change indicator

### Batteries almost discharged

When the batteries are approximately 75% used the battery symbol will flash as soon as the instrument is switched on. Although the instrument will continue to measure reliably, you should obtain replacement batteries.

### Batteries discharged— replacements required

When the batteries are discharged, the battery symbol will appear unblinking as soon as the instrument is switched on. You cannot take any additional measurements and must replace the batteries.

1. Open the battery compartment at the back of the instrument by pushing inward at the two arrows and pulling it out.
2. Replace all of the batteries at the same time – ensure they are correctly connected, as shown on the symbols in the compartment.
3. The memory retains all values although date and time must be reset - the year number therefore flashes automatically after the batteries are replaced.
4. To set date and time, follow the procedure described in Section 4.2.

### Note:

Use four new, longlife 1.5V AA batteries. Do not use batteries beyond their expiration date. If the monitor is not going to be used for a prolonged period the batteries should be removed.

**Rechargeable batteries are not recommended for use with this device.**

## 5.10. Using the AC adapter

You may also operate this monitor using the included AC Adapter (output 6V DC/600 mA with DIN plug).

Use only the included AC Adapter to avoid damaging the unit.

1. Ensure that the AC adapter and cable are not damaged.
2. Plug the adapter cable into the AC Adapter port on the blood pressure monitor.
3. Plug the adapter into your electrical outlet.

When the AC adapter is connected, no battery current is consumed.

Note: No power is taken from the batteries while the AC Adapter is connected to the monitor. If electrical power is interrupted, (e.g., by accidental removal of the AC adapter from the outlet) the monitor must be reset by removing the plug from the socket and reinserting the AC adapter connection. Please consult us if you have any questions regarding the AC adapter by calling 1-800-568-4147.

## 6. Bluetooth® functions and application (App)

### Please download the Microlife Connected Health App

(referred to as App in this document) from Apple's App Store® or Google Play® before pairing your devices.:

### Questions?


- To access the tutorial, go to the App home screen (Microlife Connected Health) and press the “Settings & Help” button.
- Press “Tutorial” and select your topic.
- If you have any issues with your setup, please contact customer service toll free: 1-800-568-4147 for additional support

### 6.1 Activating your Bluetooth®

The bluetooth® icon indicator on your blood pressure monitor, located in the left portion of the screen, is designed to provide information about the connection between your device and the monitor. Please note the following:

- The bluetooth® icon on your blood pressure monitor will flash initially when ready to connect with a device, and continue flashing for 2 minutes until a connection occurs.
- To manually activate the Bluetooth® on your monitor, press and hold the power button for approximately 8 seconds or until the Bluetooth® icon starts blinking.
- If a connection does not occur within 2 minutes, the Bluetooth® function on the monitor will automatically turn off.
- Note: At least one blood pressure reading must be taken and recorded on the monitor before the data can be transferred.

Bluetooth icon

## 6.2. Pairing your device:

- Download the App to your phone.
- Create User Profile on App.
- Confirm phone's Bluetooth® is turned on.
- Open the App on your phone.
- Select "Dashboard" from the App's home screen.
- Confirm at least one blood pressure reading has been taken and recorded on the monitor. This needs to occur before the data can be transferred.
- Press and hold the POWER/START button on the monitor for approx. 8 seconds or until the Bluetooth® icon starts blinking.
- On the App press the "sync" button located on the "Dashboard" screen.
- Data in monitor should successfully transfer over to the App on your phone.

Sync

## Setting up your blood pressure monitor:

Confirm the correct user is selected on your blood pressure monitor:


- Confirm the user icon on the bottom right hand side of the monitor screen is set to user 1 or 2.
- To change users on your monitor, see section 4.4.


## 7. Error messages/troubleshooting

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed.

(Example Error no. 1)

Error No.	Possible cause(s) / Solutions	
ERR 1	<p>The systolic pressure was determined but afterwards the cuff pressure fell below 20 mmHg.</p> <p>This situation can occur, for example, if the tube becomes unplugged after the systolic blood pressure has been measured. Additional possible cause: No pulse has been detected.</p>	
ERR 2		<p>Unnatural pressure impulses influence the measurement result. Reason: The arm was moved during the measurement (Artifact).</p>
ERR 3		<p>Inflation of the cuff takes too long. The cuff is not correctly seated or the hose connection is not sealed.</p>
ERR 5	<p>The difference between systolic and diastolic is excessive. Measure again carefully following proper cuff procedures and ensure measurement under quiet conditions. Contact your doctor if this message occurs frequently.</p>	
HI	<p>The cuff pressure is too high. Relax for 5 minutes and repeat the measurement.*</p>	
LO	<p>The pulse is too low (less than 40). Repeat the measurement.*</p>	

**\*If this or any other problem occurs repeatedly, please consult your doctor.**

## Other possible errors and their solutions

If problems occur when using the device, the following points should be checked:

### Malfunction

The display remains blank when the instrument is switched ON although the batteries are in place.

### Remedy

- Check battery installation.
- If the display is unusual, remove the batteries and then exchange them for new ones
- Check polarity

The pressure does not rise although the pump is running.

- Check the connection of the cuff tube and connect properly.

The device frequently fails to measure the blood pressure values, or the values measured are too low or too high.

- Fit the cuff correctly on the arm. (see 5.3)
- Before starting the measurement make sure that the cuff is not too tight and that clothing is not exerting pressure on the arm. Take off articles of clothing if necessary.
- Measure blood pressure again in complete peace and quiet.

Every measurement results in different values, although the device functions normally and normal values displayed.

- Please read the following information and the points listed under "Common sources of error".(5.2)
- Repeat the measurement.
- Blood Pressure changes constantly. The observed readings may accurately reflect your pressure.

The blood pressure values differ from those measured by my doctor.

- Record daily values, consult your doctor.
- Pressure readings in your doctor's office may be higher due to anxiety.
- Your blood pressure fluctuates throughout the day so record the time along with your reading.

After the instrument has inflated the cuff the pressure falls very slowly or not at all. (No reasonable measurement possible).

- Check cuff connections.
- Ensure the unit has not been tampered with.

## 8. Care and maintenance

---

- a) Do not store the unit under direct sunlight, at a high temperature, or in high humidity or dust. Performance may be degraded.
- b) The cuff contains a sensitive air-tight bubble. Handle this cuff carefully and avoid all types of stress through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use gas, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. **Do not wash the cuff in a dishwasher, clothes washer, or submerge in water.**
- d) Handle the tube carefully. Do not pull on it. Do not allow the tubing to kink and keep it away from sharp edges.
- e) **Never open the monitor.** This voids the manufacturer's warranty.

## **9. Limited warranty**

---

Your Blood Pressure Monitor is warrantied for 5 years by Microlife USA Inc, against manufacturer defects for the original purchaser only, from date of purchase.

The 5 year warranty applies to the monitor only. The following accessories are warrantied for 1 year: cuff and adapter. Batteries are not covered by this warranty. There are no user serviceable parts inside.

The warranty does not apply to consequential and incidental damages, or damage caused by batteries, improper handling, and accidents. Professional use, not following the operating instructions, and alterations made to the monitor or accessory by third parties, are also not included in this warranty. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

Microlife USA Inc will investigate your concern. A monitor or accessory as defined by this warranty, determined to be out of specification, will be replaced and shipped to you at no cost. A monitor or accessory as defined by this warranty, determined to be within specification, will be returned to you with a report of findings, at no cost.

Please use the below customer service contact information to reach Microlife USA Inc. regarding any warranty concerns. We ask that you please contact us before sending any product back in order to better identify, and more quickly process, your concern.

## 10. Technical specifications:

---

<b>Weight:</b>	445 g (with batteries)
<b>Size:</b>	150 (W) x 115 (L) x 45 (H) mm
<b>Storage temperature:</b>	-20 to +55°C (-4°F to +131°F)
<b>Humidity:</b>	15 to 90% relative humidity maximum
<b>Operation temperature:</b>	10 to 40°C ( 50°F - 104°F)
<b>Display:</b>	LCD-Display (Liquid Crystal Display)
<b>Measuring method:</b>	oscillometric
<b>Pressure sensor:</b>	capacitive
<b>Measuring range:</b>	
<b>SYS:</b>	60 to 255 mmHg
<b>DIA:</b>	40 to 200 mmHg
<b>Pulse:</b>	40 to 199 per minute
<b>Cuff pressure display range:</b>	0–299 mmHg
<b>Memory:</b>	Automatically stores the last 99 measurements for 2 users (total 198)
<b>Measuring resolution:</b>	1 mmHg
<b>Accuracy:</b>	pressure within $\pm 3$ mmHg or 2% of reading > 200 mmHg pulse $\pm 5$ % of the reading
<b>Power source:</b>	a) 4 AA batteries, 1.5 V b) AC adapter 6 V DC 600 mA (voltage 4.5 V DC to 6 V DC)
<b>Accessories:</b>	Wide-range cuff 22-42 cm (8.7" - 16.5")

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