

ViewPoint Professional Series

VPx Sensor User Manual

Document No. DV3412

September 24, 2021

Rev. Draft



Contents

1.0	Figures.....	3
2.0	Introduction.....	4
2.1	FCC NOTICE.....	5
2.2	Industry Canada.....	6
3.0	Features.....	7
3.1	Buttons.....	7
3.2	Audio Visual Alarm.....	7
3.3	Input Types.....	8
4.0	Menus.....	9
4.1	Home.....	9
4.2	Min Max.....	9
4.3	Clear Min Max.....	10
4.4	Diagnostics Menu.....	10
5.0	New Battery.....	12
5.1	Battery Replacement.....	12
6.0	Ports.....	14
7.0	Specs.....	15
7.1	RTD.....	15
7.2	Temperature and Humidity.....	15
7.3	Power.....	15
8.0	LCD Icon Legend.....	16
9.0	Operation Compatibility.....	18
9.1	DS-VP-PRO-900-S Radio Network Compatibility.....	18
9.1.1	ViewPoint G5 Compatibility.....	18
9.1.2	CheckPoint G4 Compatibility.....	18
10.0	Accessory and Probe List.....	19
11.0	VPx Sensor Models.....	19

1.0 Figures


3-1 VPx Sensor with LCD	7
3-2 VPx Sensor with LCD Buttons.....	7
Table 3-1 Input Types.....	8
4-1 Home Screen.....	9
Table 4-1 Access Min Max Menu Procedure	9
4-2 Min Max Screen	10
Table 4-2 Clear Min Max Procedure	10
Table 4-3 Access Diagnostics Screen Procedure	10
4-3 900 MHz Diagnostics Screen.....	11
4-4 Wi-Fi Diagnostics Screen.....	11
5-1 VPx Sensor Top and Back View	12
5-2 VPx Sensor Internal View (Back Cover Removed)	13
Table 5-1 Reset Battery Level Display Procedure (LCD).....	13
6-0 VPx Sensor Bottom Side View Showing Ports.....	14
Table 6-1 VPx Sensor Bottom Side View Showing Ports	14
Table 6-2 VPx Sensor RTD Probe Pin Assignments	14
Table 6-3 VPx Sensor 4-20 ma, 5v, 10v and 30v Inputs	15
Table 8-1 LCD Icon Legend.....	17
Table 9-1 Input and Accessory Type Part Numbers.....	19
Table 10-1 VPx Models	19

2.0 Introduction

The VPx sensor integrates with Mesa's ViewPoint 1.1 or above software solution. The VPx 900 MHz sensor operates in a 902 to 928 MHz range.

This device complies with Part 15 of the FCC Rules. Sensor operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Ex.



Warning: This unit is not explosion proof and is not rated for intrinsically safe installations.

2.1 FCC NOTICE

WARNING

This equipment has been tested and found to comply with the limits for a class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- * Reorient or relocate the receiving antenna.

- * Increase the separation between the equipment and receiver.

- * Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

- * Consult the dealer or an experienced radio/TV technician for help.

In order to maintain compliance with FCC regulations, shielded cables must be used with this equipment. Operation with non-approved equipment or unshielded cables is likely to result in interference to radio and TV reception. The user is cautioned that changes and modifications made to the equipment without the approval of manufacturer could void the user's authority to operate this equipment.

RF Exposure Notice: To satisfy RF exposure requirements, this device and its antennas must operate with a separation distance of at least 20 cm from all persons and must not be co-located or operating in conjunction with any other antenna or transmitter.


2.2 Industry Canada

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions: (1) This device may not cause interference. (2) This device must accept any interference, including interference that may cause undesired operation of the device.





Le présent appareil est conforme aux CNR d'Industrie Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes: (1) l'appareil ne doit pas produire de brouillage, et (2) l'utilisateur de l'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement.

3.0 Features

3.1 Buttons




3-1 VPx Sensor with LCD

	Cycle	Used for cycling through screens or menu options
	Select	Select menu option
	Previous / Next	Return to previous screen or Next option
	Mute	Mute local Audio Visual Alarms


3-2 VPx Sensor with LCD Buttons

3.2 Audio Visual Alarm


VPx Sensor has local Audible Visual alarms to notify individuals of sensor alarm states even with no access to the Viewpoint software.



LED for Visual Alarm

To mute local Audible and/or visual alarms, press the  button to silence the local alarm **ONLY**. * To perform corrective action for alarm states, do so in the ViewPoint software.

***Note:** This is only possible when the sensor is in an alarm state. It cannot be used to mute future alarms.

Muting the sensor with the  button on the device will only mute that alarm. If, after muting the local alarm, the sensor comes back into and then goes back out of range again, then the local alarm will be triggered again and will need to be muted again if desired.


3.3 Input Types

Parameter	Application
Temperature	-200 °C to + 250 °C (Type K thermocouple, up to +538 °C)
Humidity	0% to 100% (Temp 5 °C to +40 °C)
CO ₂	0% to 20%
O ₂	0% to 25% (use 4-20 mA sensor and third-party sensor)
Door Status	Open / Closed Door
Horizontal Motion	Platelet Incubator
Dry Contact	Normal v. Alarm State (Normally Closed or Open)
Power	100 to 240 V (Detect Power Outages)
DC Voltage	0 to 5 V, 0 to 10 V, 0 to 20 V, 0 to 30 V
4-20 mA	Various sensors with powered 4-20 mA interface
Differential Pressure	-1.0 to + 1.0 inches H ₂ O /-.25 to +.25 / -.5 to +.5
Particle Count	Utilizes 4-20 mA Sensor Output
Leak Detection	Detect water in critical areas

Table 3-1 Input Types

4.0 Menus

4.1 Home



4-1 Home Screen

4.2 Min Max

To get to the Min-Max screen, start at the Home screen,



	Press	Result
Step 1	 x1	Min Max Screen

Table 4-1 Access Min Max Menu Procedure



4-2 Min Max Screen

4.3 Clear Min Max




	Press	Result
Step 1	 x2	Option Screen
Step 2	 x1	Select Clear Min Max
Step 3	 x1	Confirm Clear Min Max

Table 4-2 Clear Min Max Procedure

4.4 Diagnostics Menu

The VPx sensor has onboard diagnostic capabilities that can assist in a variety of setup or troubleshooting scenarios.

To access the “Diagnostics” screen, start at the Home screen:







	Press	Result
Step 1	 x2	Option Screen
Step 2	 x1	Move  up to “Diagnostics”
Step 3	 x1	Diagnostics Screen

Table 4-3 Access Diagnostics Screen Procedure

		
	Diagnostics	
900 MHz	Radio Type:	NET ID:
RTD, Voltage, Current	CH1 Type:	HOP Table:
RTD, Voltage, Current	CH2 Type:	LP:
Transmit Rate	TX Rate:	RSSI:
Last Transmission	Last TX:	HOP CNT:
Next Transmission	Next TX:	Link State:
Read-only Memory Writes	EEPROM Wr:	Link Loss:
Flash Memory Writes	FLASH Wr:	NO ACK:
Data Packets Logged	Logged Pkts:	Retry:
Free Memory	Free Mem:	
Sensor Serial Number	SN:	FW v.
		900 MHz Network Identification
		900 MHz, Mobile Access, SpectraLink
		Link Partner
		Received Signal Strength Indicator
		Hop Count to reach destination
		Is Sensor Connected to Access Point
		Lost Link count
		No Acknowledgement from Access Point
		Connection Retries
		Firmware Version

4-3 900 MHz Diagnostics Screen

		
	Diagnostics	
Wi-Fi	Radio Type:	MAC:
RTD, Voltage, Current	CH1 Type:	LP:
RTD, Voltage, Current	CH2 Type:	RSSI:
Transmit Rate	TX Rate:	AP Auth ERRs:
Last Transmission	Last TX:	DHCP Mode:
Next Transmission	Next TX:	DHCP RSP:
Read-only Memory Writes	EEPROM Wr:	Link State:
Flash Memory Writes	FLASH Wr:	Link Loss:
Data Packets Logged	Logged Pkts:	NO ACK:
Free Memory	Free Mem:	Retry:
Sensor Serial Number	SN:	FW v.
		Sensor's MAC Address
		Link Partner
		Received Signal Strength Indicator
		No. of Errors Communicating with AP
		Dedicated IP or DHCP Mode
		Reserved IP Address on the Network
		Is Sensor Connected to Access Point
		Lost Link count
		No Acknowledgement from Access Point
		Connection Retries
		Firmware Version

4-4 Wi-Fi Diagnostics Screen

5.0 New Battery


Note: Use only 3 3.6 V Lithium-Ion Batteries (P/N 166112).

5.1 Battery Replacement

To replace the (3) 3.6V lithium batteries in the VPx sensor, press down on the two recessed areas on the rear plate of the sensor and pull back the back cover to open the unit.



5-1 VPx Sensor Top and Back View



5-2 VPx Sensor Internal View (Back Cover Removed)


After replacing the batteries, it is necessary to reset the battery level indicator. To reset the battery level indicator, start from the Home screen:

	Press	Result
Step 1	⌂ x2	Option Screen
Step 2	⬆ x2	Move ➡ up to "New Battery"
Step 3	▶ x1	New Battery Confirmation
Step 4	▶ x1	"Confirm" New Battery

Table 5-1 Reset Battery Level Display Procedure (LCD)

6.0 Ports

The VPx Sensor has 4 I/O ports that can accept 2 channels of RTD temperature probes, 2 channels of general-purpose analog probes: 4-20 ma, 5v, 10v, 20v and 30v, a discrete switch/contact input, and one I2C input for a Humidity/Temperature probe.



6-0 VPx Sensor Bottom Side View Showing Ports

VPx Sensor Ports Defined:

Port	Channel	Description
1	1	4-Wire RTD 4-Wire Cryo RTD DC Voltage - 0 to 5 V, 0 to 10 V, 0 to 20 V, 0 to 30 V DC 4-20 ma - Various sensors
2	Discrete	Door Status - Open / Closed Door Dry Contact - Normally Closed or Open
3	2	4-Wire RTD 4-Wire Cryo RTD DC Voltage - 0 to 5 V, 0 to 10 V, 0 to 20 V, 0 to 30 V DC 4-20 ma - Various sensors
4	1 and/or 2	Humidity/Temperature
5	N/A	External Power - 5V 1A AC/DC Wall Mount Adapter (TWA22 – Power Supply)

Table 6-1 VPx Sensor Bottom Side View Showing Ports

Connections for RTD Probes:

Channel	Pin	Description
1 & 2	1	RTD – black wire
	2	RTD – black wire
	3	RTD – white wire
	4	RTD – white wire

Table 6-2 VPx Sensor RTD Probe Pin Assignments

Connections for 4-20 ma, 5v, 10v and 30v Inputs:

Channel	Pin	Description
1 & 2	1	Not connected
	2	4-20 ma in; 5,10,20,30v in
	3	Ground
	4	Not connected

Table 6-3 VPx Sensor 4-20 ma, 5v, 10v and 30v Inputs

7.0 Specs

All specs below are estimated ranges for the particular inputs. Please refer to the probe specific documentation for the exact ranges and tolerances.

7.1 RTD

The VPx sensor utilizes (2) 1,000 Ohm RTDs to cover -200°C to +140°C

4-Wire RTD – (P/N CM-000188)

Range: -90°C to +140°C

Tolerance: $\pm 0.5^\circ\text{C}$ @ -25°C to +45°C ($\pm 1.0^\circ\text{C}$ @ -90°C to -26°C and +46°C to +140°C)

4-Wire Cryo RTD – (P/N CM-000189)

Range: -200°C to +70°C

Tolerance: $\pm 1.0^\circ\text{C}$ over full operating range

7.2 Temperature and Humidity

VPx monitors temperature and humidity with an external probe.

Temperature and Humidity Probe – (P/N 72112)

















Range: +5°C to 40°C and 10% RH to 90% RH

Accuracy: $\pm 0.5^\circ\text{C}$ and $\pm 3\%$ RH

7.3 Power

VPx sensor makes use of a 5V 1A AC/DC Wall Mount Adapter (TWA22 – Power Supply)

8.0 LCD Icon Legend

Icon	Description
	Alarm, above/below alarm limits
	Alarm, limits (Min-Max screen)
	Battery Full
	Battery 2/3
	Battery 1/3
	Battery Low
	A/C Power connected
	Home
	Settings
	Link broken
	Link connected
	Min Max
	Fast Transmit
	Memory 1/3 full
	Memory 2/3 full
	Memory full















	Sound On
	Sound Muted
	Icon for Min/Max screen
	Motion (not moving)
	Motion (Moving)
	Contact
	No Contact
	Door Open
	Door Closed
	Signal Quality (no signal)
	Signal Quality – Poor
	Signal Quality – Fair
	Signal Quality – Good
	Signal Quality – Best

Table 8-1 LCD Icon Legend

9.0 Operation Compatibility

9.1 DS-VP-PRO-900-S Radio Network Compatibility

The DS-VP-PRO-900-S supports two radio networks: “ViewPoint G5 Compatibility” and “CheckPoint G4 Compatibility”. The sensor can be configured for either mode using the VPx Configuration Utility. To use the utility to set the mode, first click the “Read All” button, navigate to the Wi / Radio tab and select an option in the “Compatibility Mode” field.

Radio Settings

Radio Type: Flexband

Network ID: 0000

Compatibility Mode:

Check Point G4 Compatibility

Check Point G4 Compatibility

ViewPoint G5 Compatibility

3 (SpectraLink Upper Band)

Preferred Link SN (0 if not used):

00000000

Apply Radio + Flexband Config.

Click the “Apply Radio + Flexband Config” button.

9.1.1 ViewPoint G5 Compatibility

Select “ViewPoint G5 Compatibility” mode when using the sensor with a ViewPoint radio network and when the sensor will be communicating with a ViewPoint AP.

For “ViewPoint G5 Compatibility “ (G5/G6 Mode) (902.62 to 927.62 MHz):

Hop Table	Frequency Range MHz
Main (0)	906.12 to 924.12
Low (1)	902.62 to 914.87
High (2)	914.87 to 927.62

Make sure to select the hop table that the radio network is using. The default hop table used is 0 – Main. Use the other hop tables in case there is interference issues or to segregate between radio networks.

9.1.2 CheckPoint G4 Compatibility

Select “CheckPoint G4 Compatibility” mode when using the sensor with a CheckPoint radio network and when the sensor will be communicating with a CheckPoint AP.

For “CheckPoint G4 Compatibility “ (G4 Compatibility Mode) (903.000 to 926.000 MHz):

Hop Table	Frequency Range MHz	Num Transmitting Channels
Standard (0)	906.000 to 924.000	58
Full (1)	903.000 to 926.000	58
Low (2)	903.000 to 913.325	58
High (3)	914.773 to 926.000	58

Make sure to select the hop table that the radio network is using. The default hop table used is 0 – Standard. Use the other hop tables in case there is interference issues or to segregate between radio networks.

10.0 Accessory and Probe List

Input/Accessory Type	Part Number
Standard RTD	CM-000188
Cryo RTD	CM-000189
Temp/Humidity (Snap)	72112 (requires CM-000164)
Analog Cable	CM-000284
Door Switch, Motion, Alarm Contact	CM-000183
3.6V Li Batteries	166113
AC Adapter	TWA26

Table 10-1 Input and Accessory Type Part Numbers

11.0 VPx Sensor Models

VPx Sensor Model	Part Number
ViewPoint Professional 900MHz Sensor	DS-VP-PRO-900-S

Table 10-1 VPx Models