DESIGN FOR EARS™

BE1444_005MAN002

 **Bellman & Symfon®**
DESIGN FOR EARS


Visit flash receiver

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BE1442/BE1444

EN




FR

2

Flash receiver BE1442/BE1444

Flash receiver BE1442

Overview

- 1 Rotatable reflector
- 2 Flash
- 3 Visit LED, orange
- 4 Visit LED, green
- 5 Visit LED, yellow
- 6 Visit LED, red
- 7 Connector for power supply
- 8 Telephone input for analogue telephone
- 9 Mute / test button
- 10 Bed shaker output
- 11 Cable holder
- 12 Table stand/battery compartment
- 13 External socket
- 14 Power LED

Flash receiver BE1444

Overview

- 1 Rotatable reflector
- 2 Flash
- 3 Visit LED, orange
- 4 Visit LED, green
- 5 Visit LED, yellow
- 6 Visit LED, red
- 7 Connector for power supply
- 8 Telephone input for analogue telephone
- 9 Mute / test button
- 10 Bed shaker output
- 11 Cable holder
- 12 Table stand/battery compartment
- 13 Smoke alarm menu button
- 14 Power LED
- 15 Smoke alarm supervision LED

EN

Bellman Visit flash receiver, BE1442/BE1444

General

The flash receiver is part of the Visit smart home system that helps you to keep tabs on your home. It alerts you with bright flashes when the doorbell or telephone rings, when your baby needs you or if a fire should occur.

It is activated via radio signals from the Visit transmitters and features colored LED's to differentiate the alarms. You can also connect a bed shaker (optional accessory) and wake up to vibrations if something should happen while you are asleep. For indoor use.

Getting started

Unpacking, installing and testing the unit

- 1 Connect the power supply to the receiver and the mains power. Pull the battery tab at the bottom. Place the receiver on a level surface or mount it on the wall using the BE9075 wall bracket accessory. The receiver should be placed where it best can be seen.
- 2 To test the radio link you need a Visit transmitter. Press the test button/s on the transmitter.
- 3 The receiver starts to flash and lights up the corresponding Visit LED. If a bed shaker is connected, it will vibrate. A short press on the test button repeats the last indication. If nothing happens, see **Troubleshooting**.
- 4 You can connect a Bellman bed shaker, BE1271 (optional accessory) to the bed shaker output (10) . When the flash receiver is activated, the bed shaker will vibrate.
- 5 You can connect the flash receiver to the landline phone jack via the telephone input (8)  with the BE9105 Telephone Flex (optional accessory) and an adapter plug (optional accessory). When the phone rings, the flash receiver will alert with flashes and the yellow LED (5)  will light up.

Function

Radio key

All Bellman Visit units are set to the same radio key as factory default. If you have a neighbour with a similar system, you can change the radio key so that you do not affect each other's systems. Please refer to Installer's guide in page 10. If you want to change the radio keys. All units in a Bellman Visit System must have the same radio key in order to operate as a group.

Indicators and Signals

It is generally the transmitters in the Visit smart home system that determine how the receivers will indicate an alarm. See the description in the relevant transmitter user manual for further information.

The Visit flash receiver has a function which allows you to easily check which alarm was the last one detected. Give a short press on the test button (9) and the last alarm will be repeated.

EN

System indicators

The LEDs     that indicate which transmitter has activated the Bellman Visit Flash Receiver normally have the following meanings.

- Orange LED  indicates a baby monitor.
- Green LED  indicates a door transmitter.
- Yellow LED  indicates a telephone transmitter.
- Red LED  indicates a fire alarm.

If the green  and yellow  LEDs blink alternately, this indicates that the Visit flash receiver is in radio key selection mode. The flash receiver will then wait for a radio signal from a transmitter in the Visit system which will adjust the receiver to the same radio key as the transmitter is set to.

Flash light

When the Visit flash receiver is activated, the flash (2) will show a sharp white light. The flash head (1) can be rotated to point in the direction required. The reflector can, for example, be pointed towards a wall to prevent dazzle. A silicone slip-on top is also available in a variety of colours (art. no. BE9164-BE9167).

Vibration

The Bellman Visit Flash Receiver can power a BE1271 Bellman Bed shaker (optional accessory) connected to the socket (10) . The bed shaker is placed under the pillow to wake the user up when the Visit flash receiver is activated.

Refer to the user manual for the relevant Bellman Visit transmitter for more information about vibration patterns.

Power LED indications

When the flash receiver is connected to mains power, the power LED lights up in green.

The BE1442/BE1444 flash receiver is also equipped with battery backup and the power LED indicates the following:

Power LED 	Status
Green light	Connected to mains power. Backup batteries installed.
Green blinks	Connected to mains power. No backup batteries detected.
Red light	Running on battery backup.
Red blinks	Backup batteries nearly depleted.

For further instruction about replacing the batteries, please refer to Installer's guide "Battery backup" section in Page 11.

Troubleshooting in brief

If	Try this
Nothing happens.	<ul style="list-style-type: none">▪ Check that the power supply unit is connected correctly.▪ Check that there is current in the wall socket.
The receiver is not activated.	<ul style="list-style-type: none">▪ Check the battery in the transmitters.▪ Check that the receiver is not placed too far away from the transmitters by moving it closer to the transmitters.▪ Check that the receiver is set to the correct radio key. For further information see Function/Radio key.
The receiver signals when no transmitter is activated.	<ul style="list-style-type: none">▪ Change the Radio Key on all units in the system. For further information see Appendix/Radio key.

EN

For further information, please refer to the **Installer's guide**.

Installer's guide

10	Introduction
11	Battery backup
12	Smoke alarm supervision
14	Smoke alarm supervision indication
16	Technical information

Introduction


Settings

No adjustments are required for normal use. The relevant descriptions are provided below, if you wish to change a setting for some reason.

Changing the radio key

In order to use several Visit systems close to one another without interference, different radio keys can be set on the different systems. All Visit units from the factory are set to the same radio key, channel 0. This means that all radio key switches on the transmitters are set to the OFF position.

Here is how you change the radio key:

- Open the transmitter cover and move any radio key switch to the up=on position to change the radio key. See **Changing the radio key** for the relevant transmitter.
- Press and hold the mute/test button on the top of the receiver until the green  and yellow  Visit LED's blink alternately. Release the button.
- Press the test button/s on the transmitter within 30 seconds to send the new radio key.
- All Visit LEDs     on the receiver blink 4 times to show that the radio key has been changed. It then returns to normal mode.

Please note: All Visit units must be set to the same radio key in order to operate as a group.

Signal pattern


A signal pattern decides how the receiver in the Visit system indicates during the activation. Changing the transmitter signal switch changes the signal pattern which the receivers display when the transmitter is activated.

See the transmitter user manual for further information.


Battery backup

Getting started

- Remove the battery tab at the bottom of the unit before use.
- Please charge the unit for 24 hours to ensure full backup time.

Replacing batteries

- Remove the screw at the bottom of the unit to exchange the batteries.

How to get batteries


- You can contact your local dealer or direct contact the Bellman & Symfon AB.

BE1442/BE1444

CAUTION

RISK OF EXPLOSION IF BATTERY IS REPLACED BY AN INCORRECT TYPE.

DISPOSE OF USED BATTERIES ACCORDING TO THE INSTRUCTIONS.

EN

- Use only 4x1.2VNiMH 600AAA/900AAA Lixel rechargeable batteries to ensure full backup time.
- Battery Model NO.:
BE1442: 600AAA: LH060-3AH45 or
900AAA: LH095-3AH45
BE1444: 900AAA: LH095-3AH45


Smoke alarm supervision (BE1444 only)

BE1444 has an unique feature when working with smoke alarm transmitter BE1281, which is called smoke alarm supervision function.


To activate the smoke alarm supervision function on the flash receiver, you need to manually pair the smoke alarm transmitter BE1281 by allowing the flash receiver to store the hardware ID of the smoke alarm.

1

- Press and hold the mute/test button on the top of the receiver until the green and yellow Visit LEDs blink alternately. Release the button.

2

- Press and hold the test button on the smoke alarm for more than 3 seconds until it alarms.
- Release the button.
- The flash receiver will indicate that pairing was successful by switching between the LED 🔔 🔋 ✎ 🔥 lights four times.

Indication of smoke alarm ID sequence in BE1444

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● LED lights up

☀ LED flashes

- When successfully paired, the corresponding LED will light up to indicate the smoke alarm hardware ID sequence in BE1444 memory. (Sequence indication: please refer to the Information above.)
- Attach the corresponding label (provided with each smoke alarm transmitter) on the smoke alarm transmitter unit so you can easily identify it later.

BE1444 can maximum store 8 smoke alarm hardware ID. When manually paring the 9th smoke alarm to BE1444, the four LED orange, green yellow and red will blink at the same time for 30 seconds, to indicate the hardware ID storage failure which means the smoke alarm supervision function is not available.

EN

Smoke alarm supervision indication

In case there is faulty on the BE1281 smoke alarm which is manually paired to the BE1444 as per mentioned in Page 12, the Visit LED (6) 🔦 on the BE1444 will blink in red.

To identify the faulty, please refer to the indication on Smoke alarm supervision LED (15) ✂ and the table below.



Status	Indication	Try this
	Smoke alarm supervision LED (15) ✂	
Low battery	▪ orange, blinking	<ul style="list-style-type: none"> ▪ Replace the battery using the correct batteries stated at the smoke alarm label.
Lost connection	▪ red, blinking	<ul style="list-style-type: none"> ▪ Check if the battery in the smoke alarm transmitter is still OK. ▪ Move the units to their original position and verify that the error message disappears.
Hardware failure	▪ red	<ul style="list-style-type: none"> ▪ Clean the smoke alarm transmitter according to the smoke alarm transmitter manual. ▪ If it is still faulty, replace the smoke alarm transmitter.

If there are multiple BE1281 smoke alarm manually paired, just enter the Smoke alarm menu to find the faulty unit.



Removing one smoke alarm transmitter from the system

EN

- Press the Menu button (13) on the backside of the unit.
- The Smoke alarm supervision LED (15) will light up green to indicate you are entering the menu system.
- Press the Menu button (13) until the LED (3) switches to the corresponding smoke alarm transmitter as per page 13.
- Press and hold the Smoke alarm transmitter Menu button (13) for more than 5 seconds.
- Release the button.
- The flash receiver will indicate the successful removal of the smoke alarm transmitter by switching between the LED lights four times.

Removing all smoke alarm transmitters from the system

- Press the Smoke alarm transmitter Menu button (13) on the backside of the unit.
- The Smoke alarm supervision LED (15) will light up green to indicate you are entering the menu system.
- Press and hold the Smoke alarm transmitter Menu button (13) for more than 5 seconds.
- Release the button.
- The flash receiver will indicate the successful removal of the smoke alarm transmitter by switching between the LED lights four times.

Technical information

Power supply

Mains power: 7.5 V DC / 1500 mA.
This product is also designed for IT power distribution system with phase to phase voltage of up to 230V. The socket-outlet shall be installed near the equipment and shall be easily accessible.

Battery back-up: 4x1.2V Lexel 600mAh/900mAh

- Please dispose used batteries at municipal collecting points or dispense them to local stores free of charge.

Radio function

Power consumption: Active: 1500 mA

Radio frequency: 433.92 MHz depending on region

Number of Radio Keys: 64 Radio Keys as standard. Special software can be used to increase these to 256 Radio Keys in increments of 64 per software purchase. Contact the nearest supplier for further information.

Coverage: The normal coverage between a transmitter and receiver in the Bellman Visit System is up to 80 metres with clear line of sight. Coverage is reduced if walls and large objects screen off the signals. Any thick walls constructed of reinforced concrete will greatly affect coverage. The system may also be affected by radio transmitters such as TV transmitters, computers, mobile phones, etc. This means that a unit may work perfectly in one part of the room but not at all in another.

Additional information

Programming pin (16): For service personal only (Output: 3.3V/max. 150mA)
CAUTION: Do not use, unit will be damaged when used.

Activation

Radio: Bellman Visit system
 Via analogue
 telephone network: 26 - 120 V RMS, 15 - 100 Hz.

EN

Self-Maintenance

Information: Test this device once a month

WARNING: The telephone port of the Flasher Receiver must be connected using 26 AWG or greater telecommunications wire.

Output signals

Output signals: Built-in flash light signal >10 Candela
 Warning! Flashes can cause epileptic attacks.
 Vibrator power: 2.0 – 4.0 VDC (max. 500mA)
 For indoor use only: Temperature range 0°C – 40°C
 Dimensions Ø x H: 70 x 140 mm
 Weight: 320 g
 Colour: White

Optional accessories (sold seperately)

Wall Mount Bracket BE9075
 Bellman Bed-shaker BE1271
 Telephone Flex BE9105
 Adapter plug for the appropriate country

FCC ID: WMSBE144X

FCC compliance statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Warning: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

IC compliance statement

Industry Canada




This device complies with RSS-210 of the Industry Canada Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

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TM and © 2017 Bellman & Symfon AB. All rights reserved. All product designs are patented.
Bellman & Symfon Europe AB
Bellman & Symfon Asia Ltd.