

TM-12, TM-12B Product Information

Model no: TM-12, TM-12B

Version: 0.11

Date: 08 May 2023

TM-12, TM-12B User Manual:

The TM-12, TM-12B is a wireless transmitter that supports pad, time pad, and floor mat. It is designed to monitor the status of these devices. The following features are available:

1. Learn code transmission.

The device allows for two ways to generate learn codes:

- 1.1. Turn ON the power of TM-12, TM-12B to trigger learn code transmission. The status LED stay on for about 1s and 3s for normal and long delay of alarm debounce.
- 1.2. Press the pad about 1 sec and then release it.

2. Tracking transmission.

The device transmits tracking codes about every 41s. The tracking code includes information about the pad/mat and battery condition. The LED will flash once with tracking transmission.

3. Pad/mat status change alarm/reset transmission

The device sends bursts of transmission for alarm/reset purpose:

- 1.1. After the pad is pressed and then released.
- 1.2. After the floor mat is pressed

4. Bed/Wheelchair Product type

Every transmission contains the product type of either bed or wheelchair.

5. Battery Low detection

The device enters battery low stage under the following conditions:

- 1.1. Low Cell battery voltage.
- 1.2. When the set operation days count down to zero.

After the battery low stage lasts for 7 days, the device will be locked, and all function will be disabled, even if batteries are replaced.

TM-12, TM-12B Product Information

Specifications:

Model No.	TM-12, TM-12B
Battery Type	Lithium Manganese Dioxide /CR2477 x 2
Battery Life	Min. 1 year (standby)
Operation Frequency	433.92MHz
Modulation type	ASK
ID code combinations	0001 to 4095
Supported pad types	Bed, Chair, or Floor Mat
Nominal range with CBM	30M (open area)
To work with the type of monitors individually	NGM-21, NGM-41
Operation Voltage	3V
Housing	ABS(UL94-HB Flammable Grade)
Dimensions	68.2x54.6x13mm
Net weight (with batteries)	55g
Operating temp	0°C - 40°C
Approvals	FCC

Notice:

Any changes or modifications not expressly approved by the party responsible for compliance could void your authority to operate the equipment.

Note:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced RF technician for help.

TM-12, TM-12B Product Information

This product contains a coin / button cell battery. If the coin / button cell battery is swallowed, it can cause severe internal burns in just 2 hours and can lead to death.

Keep new and used batteries away from children.

If the battery compartment does not close securely, stop using the product and keep it away from children.

If you think batteries might have been swallowed or placed inside any part of the body, seek immediate medical attention.