

User Manual (V0.11)



Features

- Includes standard lanyard with anti-ligature point
- Option to use watch strap or hospital wristband
- Splash resistant
- Low battery alert
- Changeable battery (Type: CR2450)

Functions

- Patient CALL when button is pressed
- Man-down alert in case of patient fall
- Assurance LED after alarm trigger – steady ON during transmission ~5s
- Magnetic reset option – flash 3x for 1s each after reset
- Patient tracking

Programmable Settings

Function	Options	Default
Magnetic reset	On/Off	ON
Man-down alert	On/Off	ON
Tracking signal	Off/30s/1m/2m/5m	OFF
Device ID	0001-9999	Set at factory
Area code	0-7	0

Receiving devices

- UNCG-31 CarePort
- NGM-21 Pegasus/CordFree Monitor
- CBM-03 (if possible) CordLess Monitor – send floor mat code
- CMD-11 Central Monitor Display
- PAG-11/21 Caregiver Alert Pager
- CMEX-11/21 Cordless Expander
- ZP-01 Zone Processor for Protektor/SmartBeds
- ERX-01 Field Receiver for WEAS Systems
- ILB-11 Indication Light Box

Packaging

- Individually packaged in white cardboard box

RONDISH COMPANY LIMITED

Unit G/H, 4/F, Block 1, Kwai Tak Industrial Centre, 15-33 Kwai Tak St., Kwai Chung, N.T.

Tel.: (852) 2541 6685 Fax: (852) 2541 7411 Email: rcl@rondish.com Website: www.rondish.com

FCC/CE Statement

FCC Interference Statement

This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

Notice:

Any changes or modifications not expressly approved by the party responsible for compliance could void your authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

The antenna(s) used for this transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

This device complies with FCC radiation exposure limits set forth for an uncontrolled environment. In order to avoid the possibility of exceeding the Fcc radio frequency exposure limits, human proximity to the antenna shall not be less than 20cm (8inches) during normal operation.

RF exposure statements

This equipment complies with FCC radiation exposure requirement set forth for an uncontrolled environment. End users must follow the specific operating instructions for satisfying RF exposure compliance.

Operate Frequency Bands : 433.92MHz

Caution

Risk of Explosion if battery is replaced by an incorrect type. Dispose of used batteries according to the instructions.

EU declaration of conformity

Rondish Company Limited hereby declares that the device is in compliance with the essential requirements and other relevant provisions of directives 2014/53/EU, 2014/30/EU and 2014/35/EU.

RONDISH COMPANY LIMITED

Unit G/H, 4/F, Block 1, Kwai Tak Industrial Centre, 15-33 Kwai Tak St., Kwai Chung, N.T.

Tel.: (852) 2541 6685 Fax: (852) 2541 7411 Email: rcl@rondish.com Website: www.rondish.com