SERENE®
USER'S MANUAL
Real Time
Baby Sound Sensor

Model
CA-BX

Operating Guide

Thank you for purchasing this product, please read before use
(Keep this manual for future reference)

English

Model CA-BX is an optional accessory for the **CentralAlert™** Visual Alert System. It is a wireless audio sensor for monitoring real-time baby sounds. Any loud baby sound triggers Model CA-BX to send a wireless signal to all **CentralAlert™** units. Each unit activates its own speaker, flasher, bed-shaker, and an appropriate icon for alerting. A 3-position sensitivity slide switch is provided for optimal setting.

Model CA-BX offers a real-time response. This means all **CentralAlert™** functions stop immediately, almost as soon as the baby stops crying. This minimizes continuous false triggering.

Install batteries:

1. Install two AA size alkaline batteries (user provided) and then turn on its power.
2. Verify that the (red) **TRANSMIT** indicator blinks when the **TEST** button is pressed.
3. Verify Room Noise indicators monitor room noise level.
4. This indicates the unit is normal.

Registering Model CA-BX:

This optional accessory must be learned by a Master-Unit in order to work with the system.

1. Press and hold the "**SOS/Test**" button on the Master-Unit for about 5 seconds to set it to the Learn mode. The digital display

English

- shows "**LEAr.**" (Time **SET** switch must be on **RUN** position)
2. Press and hold the "**TEST**" button on Model CA-BX for about 5 seconds until the Master-Unit beeps. (Speaker must be **ON**)
3. Registration is now complete. Press the "**SNOOZE / RESET**" button on the Master-Unit to end the Learn mode.

Testing:


1. Press the "**TEST**" button on Model CA-BX.
2. Confirm that **CentralAlert™** activates all alerting outputs in response.

Installation:

1. Place Model CA-BX on a stable

English

- surface near the baby so that its microphone port is pointing toward the baby (about 2-5 feet away).
2. Set its sensitivity switch so that normal room noise only triggers 1 or 2 indicators.
3. When room noise level or baby's crying reaches the fourth indicator, **TRANSMIT** indicator blinks – sending a wireless signal to activate all **CentralAlert™** units.

Specifications:

1. Advanced Technology: 2.4 Ghz Smart-Code
2. Operating range: Up to 200 ft, open space – actual distance varies according to house construction.
3. Up to 3 Model CA-BX's are supported
4. Battery type: Two AA size alkaline
5. Battery life: up to 10 months
6. Operating temperature: +10 to +30°C
7. Size: 2.5" x 3.8" x 1" Wt: 4 Oz
8. Spec is subject to change

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Serene Innovations, Inc.

14731 Carmenita Road, Norwalk, CA. 90650
Toll Free: 866-376 9271, Fax: 560-483-2087
Website: www.sereneinnovations.com
Email: Customercare@Sereneinnovations.com

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.