

Operator's Manual

Transmitter

ZM-540PA

ZM-541PA



About This Manual

In order to use this product safely and fully understand all its functions, read this manual before using the product. Keep this manual near the instrument or in the reach of the operator and refer to it whenever the operation is unclear.

The manual is only included on the CD. We recommend printing a copy of the electronic data for reference in case of emergency. If you require printed version of the manual, contact your Nihon Kohden representative.

Copyright Notice

The entire contents of this manual are copyrighted by Nihon Kohden. All rights are reserved. No part of this document may be reproduced, stored, or transmitted in any form or by any means (electronic, mechanical, photocopied, recorded, or otherwise) without the prior written permission of Nihon Kohden.

This product stores personal patient information. Manage the information appropriately. Patient names on the screen shots and recording examples in this manual are fictional and any resemblance to any person living or dead is purely coincidental. The contents of this manual are subject to change without notice. If you have any comments or suggestions on this manual, please contact us at: <https://www.nihonkohden.com/>

Contents

| | |
|-----------------------------------------------------------------|-----|
| Unpacking..... | i |
| GENERAL HANDLING PRECAUTIONS..... | I |
| WARRANTY POLICY..... | II |
| Equipment Authorization Requirement..... | III |
| EMC RELATED CAUTION | IV |
| Conventions Used in this Manual and Instrument..... | VI |
| Explanations of the Symbols in this Manual and Instrument | VII |
| Intended Use | 1 |
| General | 1 |
| Receiving Monitor | 3 |
| Panel Description | 4 |
| Front Panel | 4 |
| Rear Panel..... | 6 |
| Top Panel..... | 7 |
| Bottom Panel | 7 |
| Important Safety Information | 9 |
| General..... | 9 |
| Output Signal | 15 |
| Battery | 15 |
| Transmitter Channel Management | 16 |
| For Patients Using Implantable Pacemaker | 16 |
| ECG Monitoring | 17 |
| SpO ₂ Monitoring..... | 18 |
| NIBP Monitoring..... | 22 |
| Maintenance | 26 |
| Preparation on Transmitter | 27 |
| Batteries..... | 27 |
| Handling Batteries | 27 |
| Battery Lifetime..... | 27 |
| Installing and Replacing Batteries | 28 |
| Situations Requiring Battery Replacement | 29 |
| Battery Level Indication | 30 |
| Turning On the Transmitter | 31 |
| Check Items Before Use..... | 32 |
| Check Items After Power On | 32 |
| Check Items After Use..... | 33 |

| | |
|---------------------------------------------------------------------------------------------------------|----|
| Turning Off the Transmitter | 33 |
| Changing the Transmitter Channel | 34 |
| Changing Parameter and System Setup Settings | 35 |
| Notes on Parameter Settings..... | 35 |
| Changing PARAMETER SETUP Settings | 36 |
| Parameter Setup Setting List..... | 36 |
| Displaying the PARAMETER SETUP Screen..... | 37 |
| Changing Parameter Setup Settings | 38 |
| Changing SYSTEM SETUP Settings..... | 44 |
| System Setup Setting List..... | 44 |
| Displaying the SYSTEM SETUP Screen..... | 45 |
| Changing System Setup Settings..... | 46 |
| Initializing Settings | 50 |
| Attaching Electrodes, SpO ₂ Probe and NIBP Cuff to the Patient | 51 |
| Attachment Example..... | 51 |
| Selecting Electrode Leads | 52 |
| Connecting the Electrode Lead to the Transmitter..... | 53 |
| Electrode Position..... | 54 |
| Attaching Electrodes to the Patient and Connecting the Electrode Leads to Disposable Electrodes..... | 57 |
| Checking ECG on the Transmitter Screen | 58 |
| Attaching the SpO ₂ Probe..... | 59 |
| Selecting the SpO ₂ Probe | 59 |
| Connecting the SpO ₂ Probe to the Transmitter | 62 |
| Attaching the Probe to the Patient..... | 63 |
| Attaching the NIBP Cuff..... | 65 |
| Selecting the NIBP Cuff..... | 65 |
| Connecting the NIBP Cuff to the Transmitter | 68 |
| Attaching the NIBP Cuff to the Patient | 68 |
| Locking the Keys on the Transmitter | 71 |
| Monitoring..... | 72 |
| Screen Descriptions..... | 72 |
| Check Electrodes Screen | 73 |
| Numeric and Waveform Screen..... | 74 |
| Waveform Review Screen | 75 |
| Numeric Review Screen | 76 |
| Display Off | 76 |
| Basic Monitoring Operation | 77 |
| Using the Function Key..... | 77 |
| Suspending Alarms on the Receiving Monitor | 78 |
| Pausing Monitoring | 79 |
| Resuming Monitoring after Pause | 81 |

| | |
|----------------------------------------------------------------------|-----|
| Confirming the Patient | 81 |
| Turning the Display Off | 82 |
| Turning the Display On after It was Turned Off | 83 |
| ECG and Respiration Monitoring | 83 |
| Turning ECG Measurement On/Off | 86 |
| Turning Respiration Measurement On/Off..... | 86 |
| Electrode Detachment | 86 |
| SpO ₂ Monitoring..... | 87 |
| Monitoring SpO ₂ during NIBP Measurement..... | 92 |
| NIBP Monitoring..... | 93 |
| Selecting the Initial Cuff Inflation Pressure..... | 93 |
| Selecting the Measurement Mode and Interval..... | 93 |
| Measuring NIBP | 95 |
| Monitoring SpO ₂ during NIBP Measurement..... | 100 |
| Indication and Message List..... | 101 |
| Indication | 101 |
| Messages | 101 |
| Message Display Priority..... | 105 |
| Troubleshooting | 106 |
| Transmitter..... | 106 |
| ECG/Respiration..... | 107 |
| SpO ₂ | 108 |
| NIBP..... | 108 |
| Maintenance | 110 |
| 1. External Check | 111 |
| 2. Transmitter Channel | 111 |
| 3. Transmitting/Receiving Signal | 111 |
| 4. Display | 112 |
| 5. Key Operation..... | 113 |
| 6. ECG Check..... | 114 |
| 7. Respiration Check | 115 |
| 8. SpO ₂ Check | 115 |
| With the AX-300T SpO ₂ Checker | 115 |
| With the AX-410G Medical Instrument Checker..... | 116 |
| 9. NIBP Check | 117 |
| 10. NIBP Cuff for Attaching Transmitter to Patient Arm | 120 |
| Lifetime and Disposal | 121 |
| Disposing of Used Batteries | 121 |
| Battery Lifetime..... | 121 |
| Disposal | 121 |
| Disposing of Electrodes, SpO ₂ Probes and NIBP Cuffs..... | 121 |
| Disposing of Transmitter | 121 |




| | |
|---------------------------------------------------------------------------------------------------|-----|
| Cleaning, Disinfection and Sterilization | 122 |
| Transmitter and Electrode Leads | 122 |
| Cleaning | 122 |
| Disinfection | 123 |
| SpO ₂ Probe | 123 |
| NIBP Cuff | 123 |
| Periodic Inspection | 124 |
| Repair Parts Availability Policy | 124 |
| Specifications | 125 |
| ZM-540PA | 125 |
| Measured Parameters | 125 |
| Transmitted Data..... | 125 |
| Display | 125 |
| Displayed Data..... | 125 |
| ECG | 125 |
| ECG Display and Heart Rate Count | 126 |
| Respiration Measurement..... | 126 |
| SpO ₂ Measurement (ISO 9919: 2005 compliant) | 126 |
| Noninvasive Blood Pressure, NIBP (IEC 60601-2-30: 1999 compliant) | 128 |
| Transmitter | 129 |
| Power Requirements | 129 |
| Dimensions and Weight..... | 129 |
| Environment..... | 129 |
| Safety Standards | 130 |
| Electromagnetic Compatibility | 130 |
| Electromagnetic Emissions | 131 |
| Electromagnetic Immunity | 132 |
| Recommended Separation Distances between Portable and Mobile RF Communications Equipment | 134 |
| Recovery Time after Defibrillation | 134 |
| System Composition for EMC Test..... | 134 |
| ZM-541PA..... | 135 |
| Measured Parameters | 135 |
| Transmitted Data..... | 135 |
| Display | 135 |
| Displayed Data..... | 135 |
| ECG | 135 |
| Respiration Measurement..... | 136 |
| SpO ₂ Measurement (ISO 9919: 2005 compliant) | 136 |
| Noninvasive Blood Pressure, NIBP (IEC 60601-2-30: 1999 compliant) | 138 |

| | |
|---------------------------------------------------------------------------------------------------|-----|
| Transmitter | 138 |
| Power Requirements | 139 |
| Dimensions and Weight..... | 139 |
| Environment..... | 139 |
| Safety Standards | 139 |
| Electromagnetic Compatibility | 140 |
| Electromagnetic Emissions | 140 |
| Electromagnetic Immunity | 141 |
| Recommended Separation Distances between Portable and Mobile RF Communications Equipment | 143 |
| Recovery Time after Defibrillation | 143 |
| System Composition for EMC Test..... | 143 |
| Replaceable Part..... | 144 |
| Options | 144 |
| Transmitter..... | 144 |
| ECG/RESP | 145 |
| SpO ₂ | 145 |
| NIBP..... | 146 |
| Transmission Frequencies | 147 |

Unpacking

Check that all the items are included in the package. If there are any missing items, contact your Nihon Kohden representative.

The name and quantity are described under the illustration.

| | | |
|-------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------|
|  <p>ZM-540PA or ZM-541PA transmitter (1)</p> |  <p>Manuals CD (1)</p> |  <p>Certification (1)</p> |
|-------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------|

GENERAL HANDLING PRECAUTIONS

This device is intended for use only by qualified medical personnel.

Use only Nihon Kohden approved products with this device. Use of non-approved products or in a non-approved manner may affect the performance specifications of the device. This includes, but is not limited to, batteries, recording paper, pens, extension cables, electrode leads, input boxes and AC power.

Please read these precautions thoroughly before attempting to operate the instrument.

- 1. To safely and effectively use the instrument, its operation must be fully understood.**
- 2. When installing or storing the instrument, take the following precautions.**
 - (1) Avoid moisture or contact with water, extreme atmospheric pressure, excessive humidity and temperatures, poorly ventilated areas, and dust, saline or sulphuric air.
 - (2) Place the instrument on an even, level floor. Avoid vibration and mechanical shock, even during transport.
 - (3) Avoid placing in an area where chemicals are stored or where there is danger of gas leakage.
 - (4) The power line source to be applied to the instrument must correspond in frequency and voltage to product specifications, and have sufficient current capacity.
 - (5) Choose a room where a proper grounding facility is available.
- 3. Before Operation**
 - (1) Check that the instrument is in perfect operating order.
 - (2) Check that the instrument is grounded properly.
 - (3) Check that all cords are connected properly.
 - (4) Pay extra attention when the instrument is in combination with other instruments to avoid misdiagnosis or other problems.
 - (5) All circuitry used for direct patient connection must be doubly checked.
 - (6) Check that battery level is acceptable and battery condition is good when using battery operated models.
- 4. During Operation**
 - (1) Both the instrument and the patient must receive continual, careful attention.
 - (2) Turn power off or remove electrodes and/or transducers when necessary to assure the patient's safety.
 - (3) Avoid direct contact between the instrument housing and the patient.
- 5. To Shutdown After Use**
 - (1) Turn power off with all controls returned to their original positions.
 - (2) Remove the cords gently; do not use force to remove them.
 - (3) Clean the instrument together with all accessories for their next use.
- 6. The instrument must receive expert, professional attention for maintenance and repairs. When the instrument is not functioning properly, it should be clearly marked to avoid operation while it is out of order.**
- 7. The instrument must not be altered or modified in any way.**

8. Maintenance and Inspection

- (1) The instrument and parts must undergo regular maintenance inspection at the interval which is specified after the GENERAL HANDLING PRECAUTIONS section.
- (2) If stored for extended periods without being used, make sure prior to operation that the instrument is in perfect operating condition.
- (3) Technical information such as parts list, descriptions, calibration instructions or other information is available for qualified user technical personnel upon request from your Nihon Kohden representative.

9. When the instrument is used with an electrosurgical instrument, pay careful attention to the application and/or location of electrodes and/or transducers to avoid possible burn to the patient.

10. When the instrument is used with a defibrillator, make sure that the instrument is protected against defibrillator discharge. If not, remove patient cables and/or transducers from the instrument to avoid possible damage.

WARRANTY POLICY

Nihon Kohden Corporation (NKC) shall warrant its products against all defects in materials and workmanship for one year from the date of delivery. However, consumable materials such as recording paper, ink, stylus and battery are excluded from the warranty.

NKC or its authorized agents will repair or replace any products which prove to be defective during the warranty period, provided these products are used as prescribed by the operating instructions given in the operator's and service manuals.

No other party is authorized to make any warranty or assume liability for NKC's products.

NKC will not recognize any other warranty, either implied or in writing. In addition, service, technical modification or any other product change performed by someone other than NKC or its authorized agents without prior consent of NKC may be cause for voiding this warranty.

Defective products or parts must be returned to NKC or its authorized agents, along with an explanation of the failure. Shipping costs must be pre-paid.

This warranty does not apply to products that have been modified, disassembled, reinstalled or repaired without Nihon Kohden approval or which have been subjected to neglect or accident, damage due to accident, fire, lightning, vandalism, water or other casualty, improper installation or application, or on which the original identification marks have been removed.

In the USA and Canada other warranty policies may apply.

CAUTION

United States law restricts this device to sale by or on the order of a physician.

Equipment Authorization Requirement

Compliance with FCC Requirements

This device complies with Part 15 of the FCC (Federal Communications Commission) Rules. Operation is subject to the following two conditions:

- (1) this device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

This device complies with Part 95 Subpart H of the FCC Rules to be used in wireless medical telemetry service.

Operation of this equipment requires the prior coordination with a frequency coordinator designated by FCC for the Wireless Medical Telemetry Service.

CAUTION

To comply with the FCC radio frequency (RF) exposure compliance requirements, no change to the antenna or the device is permitted. Any change to the antenna or the device could result in the device, exceeding the RF exposure requirements and void user's authority to operate this device.

- NOTE**
- Use this device only indoors.
 - This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment and meets the FCC radio frequency (RF) Exposure Guidelines as this equipment has very low levels of RF energy. This transmitter must not be co-located or operated in conjunction with any other antenna or transmitter.
 - The devices require registration and deployment by an authorized frequency coordinator. The ASHE (American Society for Healthcare Engineering) has been designated by the FCC to manage the WMTS frequencies. This device has frequency bands which may not be used in some areas. For details, contact your Nihon Kohden representative. For details on the guidelines, refer to the ASHE home page.

EMC RELATED CAUTION

This equipment and/or system complies with IEC 60601-1-2 International Standard for electromagnetic compatibility for medical electrical equipment and/or system. However, an electromagnetic environment that exceeds the limits or levels stipulated in IEC 60601-1-2, can cause harmful interference to the equipment and/or system or cause the equipment and/or system to fail to perform its intended function or degrade its intended performance. Therefore, during the operation of the equipment and/or system, if there is any undesired deviation from its intended operational performance, you must avoid, identify and resolve the adverse electromagnetic effect before continuing to use the equipment and/or system.

The following describes some common interference sources and remedial actions:

1. Strong electromagnetic interference from a nearby emitter source such as an authorized radio station or cellular phone:
Install the equipment and/or system at another location. Keep the emitter source such as cellular phone away from the equipment and/or system, or turn off the cellular phone.
2. Radio-frequency interference from other equipment through the AC power supply of the equipment and/or system:
Identify the cause of this interference and if possible remove this interference source. If this is not possible, use a different power supply.
3. Effect of direct or indirect electrostatic discharge:
Make sure all users and patients in contact with the equipment and/or system are free from direct or indirect electrostatic energy before using it. A humid room can help lessen this problem.
4. Electromagnetic interference with any radio wave receiver such as radio or television:
If the equipment and/or system interferes with any radio wave receiver, locate the equipment and/or system as far as possible from the radio wave receiver.
5. Interference of lightning:
When lightning occurs near the location where the equipment and/or system is installed, it may induce an excessive voltage in the equipment and/or system. In such a case, disconnect the AC power cord from the equipment and/or system and operate the equipment and/or system by battery power, or use an uninterruptible power supply.
6. Warning: Use with other equipment:
When the equipment and/or system is adjacent to or stacked with other equipment, the equipment and/or system may affect the other equipment. Before use, check that the equipment and/or system operates normally with the other equipment.

Caution - continued

7. Warning: Use of unspecified accessory, transducer and/or cable:

When an unspecified accessory, transducer and/or cable is connected to this equipment and/or system, it may cause increased electromagnetic emission or decreased electromagnetic immunity. The specified configuration of this equipment and/or system complies with the electromagnetic requirements with the specified configuration. Only use this equipment and/or system with the specified configuration.

8. Use of unspecified configuration:

When the equipment and/or system is used with the unspecified system configuration different than the configuration of EMC testing, it may cause increased electromagnetic emission or decreased electromagnetic immunity. Only use this equipment and/or system with the specified configuration.

9. Measurement with excessive sensitivity:

The equipment and/or system is designed to measure bioelectrical signals with a specified sensitivity. If the equipment and/or system is used with excessive sensitivity, artifact may appear by electromagnetic interference and this may cause mis-diagnosis. When unexpected artifact appears, inspect the surrounding electromagnetic conditions and remove this artifact source.

10. Use with radiation therapy equipment:



When the equipment and/or system is used in a radiotherapy room, it may cause failure or malfunction due to electromagnetic radiation or corpuscular radiation.

When you bring the equipment and/or system into a radiotherapy room, constantly observe the operation. Prepare countermeasures in case of failure or malfunction.

If the above suggested remedial actions do not solve the problem, consult your Nihon Kohden representative for additional suggestions.

Conventions Used in this Manual and Instrument













Warnings and Cautions

| Level | Description |
|--------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  WARNING | A warning alerts the user to possible injury or death associated with the use or misuse of the instrument. |
|  CAUTION | A caution alerts the user to possible injury or problems with the instrument associated with its use or misuse such as instrument malfunction, instrument failure, damage to the instrument, or damage to other property. |







Explanations of the Symbols in this Manual and Instrument

The following symbols are used with this transmitter. The description of each symbol is shown in the table below.

On Panel

| Symbol | Description | Symbol | Description |
|-----------------------------------------------------------------------------------|---------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------|
|  | Change screen |  | Defibrillation proof type BF applied part |
|  | Attention, consult operator's manual |  | Defibrillation proof type CF applied part |
|  | Moves cursor, scrolls data |  | Serial number |
|  | Direction for attaching battery cover |  | Date of manufacture |
|  | Direction for inserting battery |  | RF transmitter Non-ionizing radiation |
|  | Direct current |  | CSA mark |

On LCD

| Symbol | Description | Symbol | Description |
|-------------------------------------------------------------------------------------|---------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------|
|  | Full battery |  | Battery very weak Cannot measure NIBP Replace battery |
|  | Battery 1/3 full |  | Alarm suspended |
|  | Battery weak Replace battery |  | QRS/pulse sync mark |

Intended Use

General

The ZM-540PA and ZM-541PA transmitters transmit ECG, respiration and pulse waveforms, SpO₂ and NIBP data from a patient to a Nihon Kohden monitor for continuous monitoring. The front LCD displays ECG (or pulse wave), numeric values of monitoring parameters, NIBP measuring mode and interval, messages and battery condition.¹ They also display the compressed waveform and numeric data of the latest 10 minutes.

This transmitter is designed for use by qualified medical personnel in a medical facility such as a hospital or clinic. It is not designed for use outdoors or in a home environment, and is not designed to be operated by the patient themselves.

¹ Essential performance of this transmitter.

The difference between the ZM-540PA and ZM-541PA is the transmission frequency range.

ZM-540PA: 608.0250 to 613.9750 MHz (channel number 9002 to 9478)

ZM-541PA: 1395.0250 to 1399.9750 MHz (channel number E002 to E398)
1427.0250 to 1431.9750 MHz (channel number E502 to E898)

- NOTE
- The transmitter channel can be changed with a QI-901PK channel writer.
 - Read the operator's manual for the receiving monitor together with this manual before use.

WARNING

Do not diagnose a patient based only on data acquired by the transmitter. Overall judgement must be performed by a physician who understands the features, limitations and characteristics of the transmitter and by reading the biomedical signals acquired by other instruments.

WARNING

Do not use this transmitter for monitoring a patient. Monitor the patient on a receiving monitor. This transmitter displays the waveforms and measured values but does not have an alarm function.

⚠ WARNING

Do not use the same transmitter for more than one patient at the same time. Do not connect different sensors from different patients to the same transmitter.

⚠ CAUTION

- Do not use the same channel for different patients. If the same channel is used for two patients, the two patients' data will be lost due to mutual modulation interference, or another patient's data may appear on the receiving monitor screen.
- Do not use two transmitters with adjacent channels in the same hospital. If transmitters with adjacent channels are used, their radio waves interfere with each other.

⚠ CAUTION

Signal loss and artifact may occur because of the multipath cancellation* when using a transmitter.

* Multipath Cancellation (Standing Wave Interference)

When a radio wave reflects off a surface, there may be some points in the room where the reflected and direct waves are exactly out of phase. At these points in the room, the reflected and direct waves cancel each other out so that the signal strength is decreased.

Locations where signal loss occurs are called "null spots". If the transmitter is moving or nearby people or objects are moving, null spots can occur anytime and anywhere.

- NOTE**
- To prevent interference between channels, assign a channel administrator in the hospital and only he or she should manage channel assignment.
 - Use Nihon Kohden parts and accessories to assure maximum performance from your instrument.
 - For stable signal reception, it is recommended to use a diversity antenna system on the receiving monitor. Otherwise, spike noise from transient fading of electric field strength (for example, people moving) may interfere with the transmitter signal and may be mistaken as an arrhythmia on the receiving monitor.

- If the transmitter is used outside the specified environment, its performance cannot be guaranteed.
- For details on antennas and antenna construction, contact your Nihon Kohden representative. You can also refer to the Telemetry System Installation Guide.


Receiving Monitor

Any Nihon Kohden receiving monitor (central monitor with multiple patient receiver) can receive signals from this transmitter as long as the protocol version and channel setting are the same on the receiving monitor and transmitter.

- NOTE
- For details on the receiving monitor and upgrade information, contact your Nihon Kohden representative.
 - The transmitter does not give any patient alarm, only a “no battery” alarm. Patient alarms must be managed on the receiving monitor.

Panel Description

Front Panel



1 LCD

Displays numeric values, ECG or pulse wave, NIBP measuring mode and interval, messages and battery status. For details, refer to the “Screen Descriptions” section.

2 Infrared receiver

Used for upgrading the transmitter software.

3 Screen key

Toggles the screen in the following order.

After power on: Start up → Check electrodes → Numeric and waveform → Waveform review → Numeric review → Display off → Check electrodes ...

After auto display off: Numeric and waveform → Waveform review → Numeric review → Display off → Check electrodes → Numeric and waveform ...


On a SETUP or CHECK screen, this key cancels changing setting or exits the screen.

4 Function key

Depending on the setting on the transmitter, this key suspends alarms, pauses monitoring on the receiving monitor or transmits “Patient confirmed” message.

On a SETUP screen, this key registers the selected setting and moves the cursor to the next setting item.

On a CHECK screen, this key starts or stops maintenance test.



Refer to the WARNING on the next page.

5 NIBP Start/Stop key

Starts/stops NIBP measurement in selected mode.

6 NIBP Interval key

Selects NIBP measurement mode.

7 Lead/Scroll keys

On the numeric and waveform screen, these keys change the ECG lead.

On the waveform review screen, these keys scroll data.

On a SETUP screen, these keys move the cursor.


8 Battery case

Contains three 1.5 V dry cell batteries (AA TYPE).

⚠ WARNING

Close the battery case cover during operation. If the transmitter is used with the battery case cover open, anyone who touches the opened battery case may receive an electrical shock when defibrillation is performed. Touching the opened battery case may cause electrostatic discharge and intermittently interfere with the waveform or data.


Rear Panel




WARNING

This transmitter is not waterproof. If detergent or liquid spills into the transmitter, stop using it and contact your Nihon Kohden representative. If a wet transmitter is used, the patient or operator may receive an electrical shock or injury.

Top Panel



Bottom Panel



⚠ WARNING

Before defibrillation, all persons must keep clear of the bed and must not touch the patient or any equipment or cord connected to the patient. Failure to follow this warning may cause electrical shock or injury.

⚠ WARNING

When performing defibrillation, discharge as far as possible from electrodes, patches and any gel, cream or medicine on the chest of the patient. If there is a possibility that the defibrillator paddle could touch these materials, remove them from the patient. If the defibrillator paddle directly contacts these materials, the discharged energy may cause skin burn to the patient.

⚠ WARNING

When the transmitter is used with an electrosurgical unit (ESU), firmly attach the entire area of the ESU return plate. Otherwise, the current from the ESU flows into the electrodes of the transmitter, causing electrical burn where the electrodes are attached. For details, refer to the ESU manual.

⚠ CAUTION

Do not shake or swing the transmitter while holding the leads or cables connected to the transmitter. The transmitter may come off and injure someone or damage surrounding instruments.

Important Safety Information

General

⚠ WARNING

Never use the transmitter in the presence of any flammable anesthetic gas or high concentration oxygen atmosphere. Failure to follow this warning may cause explosion or fire.

⚠ WARNING

When performing MRI test, remove all electrodes and probe from the patient which are connected to this transmitter. Failure to follow this warning may cause skin burn on the patient. For details, refer to the MRI manual.

⚠ WARNING

Never use the transmitter in a hyperbaric oxygen chamber. Failure to follow this warning may cause explosion or fire.

⚠ WARNING

Do not take this transmitter into the MRI test room. This transmitter is not designed to be used during MRI tests.

⚠ WARNING

Before defibrillation, all persons must keep clear of the bed and must not touch the patient or any equipment or cord connected to the patient. Failure to follow this warning may cause electrical shock or injury.

⚠ WARNING

When performing defibrillation, discharge as far as possible from electrodes, patches and any gel, cream or medicine on the chest of the patient. If there is a possibility that the defibrillator paddle could touch these materials, remove them from the patient. If the defibrillator paddle directly contacts these materials, the discharged energy may cause skin burn to the patient.

⚠ WARNING

When the transmitter is used with an electrosurgical unit (ESU), firmly attach the entire area of the ESU return plate. Otherwise, the current from the ESU flows into the electrodes of the transmitter, causing electrical burn where the electrodes are attached. For details, refer to the ESU manual.

⚠ WARNING

Do not use the same transmitter for more than one patient at the same time. Do not connect different sensors from different patients to the same transmitter.

⚠ WARNING

This transmitter is not waterproof. If detergent or liquid spills into the transmitter, stop using it and contact your Nihon Kohden representative. If a wet transmitter is used, the patient or operator may receive an electrical shock or injury.

⚠ WARNING

Close the battery case cover during operation. If the transmitter is used with the battery case cover open, anyone who touches the opened battery case may receive an electrical shock when defibrillation is performed. Touching the opened battery case may cause electrostatic discharge and intermittently interfere with the waveform or data.

⚠ WARNING

When the signal is unstable, keep the patient under close observation. When the signal is unstable, the monitoring and alarm are not reliable and the receiving monitor cannot detect a sudden change of the patient's condition. This may cause critical changes in the patient condition to be overlooked. Install an appropriate antenna system to ensure stable signal condition.

⚠ WARNING

While the "ALARMS SUSPENDED" message is displayed on the transmitter, all alarms on the receiving monitor are suspended so keep the patient under close observation.

⚠ WARNING

When the patient returns to the bed, turn on the transmitter and check that the monitoring is resumed on the receiving monitor.

⚠ WARNING

If the transmitter is not turned off and monitoring continues for the selected interval, pause monitoring is canceled and monitoring continues. Check that the monitoring is resumed on the receiving monitor.

⚠ WARNING

Do not allow the conductive part of the connector which is connected to the patient to contact other conductive parts including earth. This causes leakage current and incorrect measurement value and leads to wrong diagnosis.

⚠ WARNING

After admitting a patient on the central monitor and attaching electrodes and sensors on the patient and connecting cables to the transmitter, check that there is no error messages and that the waveforms and numeric data are appropriately displayed on the transmitter and central monitor screen. If there is an error message, or waveform or numeric data is not appropriate, check the electrodes and sensors attachment, patient condition and settings on the transmitter and remove the cause.

⚠ CAUTION

Only use Nihon Kohden specified electrodes, electrode leads, SpO₂ probes, and NIBP cuffs. Otherwise, the maximum performance from the transmitter cannot be guaranteed.

⚠ CAUTION

Do not reuse disposable parts and accessories.

⚠ CAUTION

For handling and precautions on electrodes, electrode leads, SpO₂ probes and NIBP cuffs, refer to the manual.

⚠ CAUTION

The measurement values and displayed waveforms on the transmitter and receiving monitor may be different due to timing delay of the display or difference in detection settings.

⚠ CAUTION

Do not shake or swing the transmitter while holding the leads or cables connected to the transmitter. The transmitter may come off and injure someone or damage surrounding instruments.

⚠ CAUTION

Signal loss and artifact may occur because of the multipath cancellation* when using a transmitter.

* Multipath Cancellation (Standing Wave Interference)

When a radio wave reflects off a surface, there may be some points in the room where the reflected and direct waves are exactly out of phase. At these points in the room, the reflected and direct waves cancel each other out so that the signal strength is decreased. Locations where signal loss occurs are called “null spots”. If the transmitter is moving or nearby people or objects are moving, null spots can occur anytime and anywhere.

⚠ CAUTION

Turn off the power of mobile phones, small wireless devices and other devices which produce strong electromagnetic interference around a patient (except for devices allowed by the hospital administrator). Radio waves from devices such as mobile phones or small wireless devices may be mistaken as pulse waves or respiration waves and the displayed data may be incorrect.

⚠ CAUTION

When monitoring respiration is needed, measure respiration with an instrument. The transmitter calculates SpO₂ of arterial blood based on the principle of pulse oximeter and does not measure respiration.

⚠ CAUTION

- Do not use the same channel for different patients. If the same channel is used for two patients, the two patients' data will be lost due to mutual modulation interference, or another patient's data may appear on the receiving monitor screen.
- Do not use two transmitters with adjacent channels in the same hospital. If transmitters with adjacent channels are used, their radio waves interfere with each other.

⚠ CAUTION

When monitoring SpO₂ only, detection of arrhythmia and asystole is not available and arrhythmia alarms such as ASYSTOLE, VF or VT are not available. If the patient requires ECG monitoring, monitor the ECG.

⚠ CAUTION

When monitoring SpO₂ only (without ECG monitoring), turn on both the upper and lower limit alarms for PR and SpO₂ on the receiving monitor. If the patient's pulse is not detected during asystole or other condition, a "CANNOT DETECT PULSE" or "CHECK PROBE" alarm occurs instead of an SpO₂ limit alarm. Furthermore, if the patient has no pulse, noise from probe movement could be misjudged as a pulse and cause an incorrect PR or SpO₂ value to be displayed.

⚠ CAUTION

Always lock the battery cover while using the transmitter. Otherwise the battery cover may come off and correct measurement might not be performed.

Output Signal

⚠ WARNING

Do not use the output signal from the receiving monitor as the synchronization signal for other equipment such as IABP, MRI, echocardiography or defibrillator. There may be time delay between the monitor and the other equipment caused by waveform transmission delay and spike noise may interfere on the output signal and be mistaken as a trigger.

Battery

⚠ CAUTION

Battery replacement must be performed by the operator. When replacing the batteries of a transmitter that is currently used for a patient, disconnect the electrode leads from the transmitter before replacing batteries or do not touch the patient during replacement.

⚠ CAUTION

Refer to the battery and battery charger manuals for details on handling the batteries.

Transmitter Channel Management

WARNING

The following actions must be taken to properly receive the transmitter signal of the correct patient on the receiving instrument. Otherwise, there may be signal loss or signals may mix causing a serious accident, such as monitoring a different patient.

- Assign a channel administrator in the hospital and only he or she should manage channel assignment.
- The channel administrator must manage the channels in the facility so that there is no signal interference.
- When the transmitter channel is changed, the channel administrator must check that the channel on the receiving monitor is also changed and the signal is properly received.
- The channel administrator must replace the channel number label on the transmitter with the new one after changing the channel.

For Patients Using Implantable Pacemaker

WARNING

The bioelectric impedance measurement sensor of a minute ventilation rate-adaptive implantable pacemaker may be affected by transmitter which is connected to the same patient. If this occurs, the pacemaker may pace at its maximum rate and the transmitter may give incorrect data to the monitor. If this occurs, disconnect the electrode leads from the patient or change the setting on the pacemaker by referring to the pacemaker's manual. For more details, contact your pacemaker representative or Nihon Kohden representative.

ECG Monitoring

WARNING

Turn the pacing pulse detection to ON on the receiving monitor when monitoring a pacemaker patient. Otherwise the pacemaker pulse is not rejected. However, even when the pacing pulse detection is set to ON, the pacemaker pulse might not be rejected. When the pacemaker pulse is not rejected, the pacemaker pulse is detected as QRS and false heart rate may be indicated or critical arrhythmia such as asystole may be overlooked. Keep pacemaker patients under close observation.

For the pacemaker pulse rejection capability of the ZM-540PA and ZM-541PA transmitters refer to the “Specifications - ECG” section.

WARNING

Even when the pacing pulse detection is set to ON on the receiving monitor, the pacemaker pulse can be overlooked or detected as QRS. You cannot confirm the pacemaker operation only from the detected pacemaker pulse.

CAUTION

Only use Nihon Kohden specified electrodes and electrode leads. When other type of electrodes or electrode leads are used, the “CHECK ELECTRODES” message may be displayed and ECG monitoring may stop.

CAUTION

When the “CHECK ELECTRODES” message is displayed on the receiving monitor, ECG is not monitored properly and the ECG alarm does not function. Check the electrode, electrode leads, and if necessary, replace with new ones.

SpO₂ Monitoring

⚠ WARNING

SpO₂ measurement may be incorrect in the following cases.

- When the patient's carboxyhemoglobin or methemoglobin increases abnormally.
- When dye is injected in the blood.
- During CPR.
- When measuring at a site with venous pulse.
- When there is body movement.
- When the pulse wave is small (insufficient peripheral circulation).

⚠ WARNING

When monitoring SpO₂ of a patient who is receiving photodynamic therapy, the light from the finger probe sensor may cause a burn.

Photodynamic therapy uses a photosensitizing agent that has a side effect of photosensitivity.

The SpO₂ probes manufactured by Nihon Kohden have two wavelengths with peaks in the range of 650 and 950 nm. The maximum light intensity is less than 5.5 mW.

⚠ WARNING

When not monitoring SpO₂, disconnect the SpO₂ cable from the transmitter. Otherwise, noise from the probe sensor may interfere and incorrect data is displayed on the screen.

⚠ WARNING

- When using the TL-201T finger probe, do not fasten the probe and cable to the finger by wrapping with tape. This may cause burn, congestion or pressure necrosis from poor blood circulation.
- When using probes other than the TL-201T finger probe, to avoid poor circulation, do not wrap the tape too tight. Check the blood circulation condition by observing the skin color and congestion at the skin peripheral to the probe attachment site. Even for short-term monitoring, there may be burn or pressure necrosis from poor blood circulation, especially on neonates or low birth weight infants whose skin is delicate. Accurate measurement cannot be performed on a site with poor peripheral circulation.

⚠ WARNING

Check the circulation condition by observing the skin color at the measurement site and pulse waveform. Change the measurement site every 8 hours for disposable probes and every 4 hours for reusable probes (every 8 hours for TL-631T series probe). The skin temperature may increase at the attached site by 2 or 3°C (4 or 5°F) and cause a burn or pressure necrosis. When using the probe on the following patients, take extreme care and change the measurement site more frequently according to symptoms and degree.

- Patient with a fever
- Patient with insufficient peripheral circulation
- Neonate or low birth weight infant with delicate skin

 **CAUTION**

The disposable probe is not sterilized. Use the disposable probe only for a single patient. Never reuse the disposable probe for another patient because it causes cross infection.

 **CAUTION**

If the attachment site is dirty with blood or bodily fluids, clean the attachment site before attaching the probe. If there is nail polish on the attachment site, remove the polish. Otherwise, the amount of transmitted light decreases, and measured value may be incorrect or measurement cannot be performed.

 **CAUTION**

If the skin gets irritated or redness appears on the skin from the probe, change the attachment site or stop using the probe. Take extreme care for the patients with delicate skin.

 **CAUTION**

While a patient is on medication which causes vasodilation, the pulse waveform may change and in rare cases the SpO₂ value might not be displayed.

⚠ CAUTION

NIBP and SpO₂ can be measured on the same limb, but the SpO₂ monitoring might not be accurate during NIBP measurement. Be careful when reading the SpO₂ values.*

*** Monitoring SpO₂ during NIBP Measurement**

When the SpO₂ probe is attached to the same limb as the NIBP cuff, the blood flow decreases during NIBP measurement and pulse wave cannot be detected and SpO₂ cannot be monitored properly. When “INHIBIT SpO₂ DURING NIBP” on the PARAMETER SETUP screen is set to ON (factory default setting), SpO₂ monitoring is paused during NIBP measurement to avoid SpO₂ alarm occurrence. However, when monitoring SpO₂ on the same limb as the NIBP, be careful when reading SpO₂ values.

⚠ CAUTION

Normal external light does not affect monitoring but strong light such as a surgical light or sunlight may affect monitoring. If affected, cover the measuring site with a blanket.

⚠ CAUTION

When a message indicates a faulty probe, stop monitoring and replace the probe with a new one.

⚠ CAUTION

Handle the probe cable according to the following cautions. Failure to follow these cautions may cause cable discontinuity or short circuit of the probe cable which may cause incorrect measurement data or inability to perform measurement. Also in rare cases, the probe temperature may increase and cause skin burn on the patient. If the probe cable is damaged, replace the probe with a new one.

- Do not pull or bend the probe cable.
- Do not let caster feet run over the probe cable.

 **CAUTION**

Do not use a probe which is deteriorated by aging. Accurate measurement cannot be performed.

 **CAUTION**

Do not use a damaged or disassembled probe. It causes incorrect measurement and may injure the patient.

 **CAUTION**

When the probe is attached on an appropriate site with sufficient thickness and the error message confirming the probe attachment repeatedly appears, the probe may be deteriorated. Replace it with a new one.

NIBP Monitoring

 **WARNING**

Be careful when measuring NIBP on a patient with known bleeding disorders or coagulation. After NIBP measurement, there may be dot hemorrhage, or circulatory disorder by thrombus where the cuff is attached.

⚠ WARNING

NIBP measurement may be incorrect in the following situations.

- When using an ESU
- Body movement
- Small pulse wave
- Too many arrhythmias
- Shaking from an external source
- Rapid blood pressure change
- During CPR
- Slow pulse
- Low blood pressure
- Small pulse pressure
- Cuff is too tight or too loose
- Cuff does not fit the arm
- Cuff is wrapped over thick clothing
- Cuff is deteriorated
- Arterial sclerosis
- Poor perfusion
- Diabetes
- Age
- Pregnancy
- Pre-eclampsia
- Renal disease
- Shivering
- Trembling

⚠ WARNING

When performing NIBP measurements in STAT mode or 5 minute intervals, periodically remove the cuff from the patient for ventilation. The skin temperature may increase at the cuff attachment site by 2 or 3°C (4 or 5°F). When measuring a patient with a fever or peripheral circulation insufficiency, it may cause a burn.

⚠ WARNING

Do not attach the NIBP cuff on a limb which is being used for intravascular access or therapy, or an arterio-venous (A-V) shunt. It may cause reflux of blood or medicinal solution or block injection of medicinal solution due to poor blood circulation.

⚠ WARNING

Do not attach the NIBP cuff on an arm which is the same side as a mastectomy. It may cause circulatory disorder such as swelling from poor blood circulation.

⚠ WARNING

While measuring NIBP, if the NIBP cuff and other medical equipment are attached to the same limb, the medical equipment might not function temporarily.

⚠ WARNING

While measuring NIBP, check the patient condition and confirm that the NIBP cuff does not affect blood circulation.

⚠ WARNING

When performing long term measurement at intervals less than 5 minutes, perform measurements while observing the state of the patient, blood vessels and limb to ensure adequate circulation. Congestion may occur at the measurement site. When performing periodic measurement for a long time, periodically check the circulation condition.

 **CAUTION**

Do not wrap the cuff on an arm or thigh which is used for injection. NIBP measurement on an arm or thigh which is used for injection may cause reflux of blood and stop injection.

 **CAUTION**

Do not wrap the cuff too tight. It may cause poor blood circulation and congestion. If the cuff is wrapped too loosely, the NIBP value may increase.

 **CAUTION**

Do not attach the cuff to the site where there is injury or inflammation. If the skin gets irritated or redness appears on the skin from the cuff, change the attachment site or stop using the cuff. Take extreme care on the patients with delicate skin.

 **CAUTION**

When using an extension hose, check that the extension hose is not bent or squeezed. Otherwise, the cuff might not inflate or deflate. If the cuff cannot deflate, it may cause congestion on the patient at the cuff attachment site.

 **CAUTION**

When performing NIBP measurement repeatedly, have a rest between measurements to recover adequate circulation.

Maintenance

 **CAUTION**

If detergent or liquid spills into the transmitter, stop cleaning or disinfecting it and contact your Nihon Kohden representative. The transmitter needs to be checked for safety and function before use.

 **CAUTION**

Before maintenance, cleaning or disinfection, turn the transmitter power off and remove the batteries. Failure to follow this instruction may result in electrical shock and transmitter malfunction.

 **CAUTION**

Never disassemble or repair the transmitter. If there is any problem with the transmitter, contact your Nihon Kohden representative.

 **CAUTION**

Dispose of Nihon Kohden products according to your local laws and your facility's guidelines for waste disposal. Otherwise, it may affect the environment. If there is a possibility that the product may have been contaminated with infection, dispose of it as medical waste according to your local laws and your facility's guidelines for medical waste. Otherwise, it may cause infection.

Preparation on Transmitter

Batteries

Handling Batteries

⚠ WARNING

Do not handle the batteries with wet hands.

⚠ WARNING

When the transmitter is not in use, remove the batteries. When the batteries are installed, battery power is consumed even when measurement is not performed. When NiMH batteries are left in the transmitter when it is not used, the batteries may overdischarge and leak liquid which makes the batteries unusable and damages the transmitter.

⚠ CAUTION

Refer to the battery and battery charger manuals for details on handling the batteries.

- NOTE**
- Remove the batteries from the transmitter before disposing of it.
 - Use either new alkaline batteries or fully charged rechargeable NiMH batteries. Using unspecified batteries, previously used batteries or batteries that have been stored for long periods may result in short battery life or reduced performance resulting in unstable measurement.

Battery Lifetime

Use three AA type alkaline dry cell batteries. NiMH batteries can also be used but the lifetime of alkaline batteries is longer. The lifetime of NiMH batteries is about 1/2 of alkaline batteries (when fully charged).

With new Nihon Kohden recommended alkaline batteries, the transmitter can continuously measure ECG, respiration, SpO₂ and NIBP for approximately 1 day. The measurement is performed at room temperature, NIBP is measured in auto mode at 60 minute intervals and SpO₂ is measured on an index finger of a male patient with weight 60 kg. Operation time depends on the thickness of the SpO₂ probe attachment site.

Recommended batteries

NiMH secondary: Panasonic eneloop

Alkaline primary: Nihon Kohden Medipower (equivalent to Panasonic LR6 (G))


Installing and Replacing Batteries

CAUTION


Battery replacement must be performed by the operator. When replacing the batteries of a transmitter that is currently used for a patient, disconnect the electrode leads from the transmitter before replacing batteries or do not touch the patient during replacement.

If electrode leads are attached to the patient and the person replacing batteries touches the patient during battery replacement, excess patient leakage current may flow into the patient.


- NOTE**
- Replace all batteries at the same time.
 - Do not use different types of batteries together.
 - Insert the batteries with the correct polarity (+ and -).
 - The capacity of rechargeable NiMH batteries is reduced if the batteries are recharged before they are fully discharged. For details, refer to the battery operator's manual.




1 Remove the battery case cover.



2 Insert three new or fully charged batteries into the battery case observing the correct polarity.



NOTE: Insert the (-) end of the battery first and press it against the spring. If you try to force the (+) end of the battery in first as shown, it will deform the spring and damage the battery and transmitter.



3 Close the cover.

The transmitter is automatically turned on when the batteries are installed and the cover is closed.

⚠ CAUTION

Always lock the battery cover while using the transmitter. Otherwise the battery cover may come off and correct measurement might not be performed.

Situations Requiring Battery Replacement





Replace the batteries when any of the following occurs:

- The transmitter displays the “BATTERY WEAK” message or  icon.
- The transmitter generates a constant alarm (continuous “peep” sound).
- The transmitter LCD does not display anything when the power is turned on.
- The receiving monitor displays a battery replacement message.

NOTE: The battery replacement message or icon on the transmitter disappears after a short time depending on the type of battery. When the battery replacement indication appears, immediately replace the transmitter battery with a new one.

Battery Level Indication

The following icons on the display indicate the battery level. When the screen is turned off, press the Screen key to check the battery level. When "PROTOCOL" of the SYSTEM SETUP screen is set to 57, the battery level indication is transmitted to the receiving monitor.

| Indication | Battery Level | Message on the Receiving Monitor |
|-----------------------------------------------------------------------------------|------------------------------------------------------------------------|----------------------------------------------------------------------------------------|
|  | Fully charged batteries | There is no message on the monitor. |
|  | Batteries are 1/3 full. | |
|  | Batteries are weak. Replace batteries. | Message requiring battery replacement is displayed. |
|  | Batteries are very weak. Cannot measure NIBP. Replace batteries. | |
| No indication | Dead batteries | No signal can be transmitted to the monitor. There is no indication on the monitor. |

Turning On the Transmitter

When the batteries are installed correctly, the power is turned on. A “peep” sounds for one second, the startup screen appears, then the check electrodes screen appears. (There is no “peep” sound when there is no battery power.)

NOTE: Check that the “COMMUNICATION LOSS” message is not displayed on the central monitor.



After checking that the ECG is stable on the check electrodes screen, press the Screen key to display the numeric and waveform screen.



For details on the screen, refer to the “Screen Descriptions” section.

Check Items Before Use

Before turning on the transmitter power, check the following to confirm that the transmitter can be used in normal and safe condition.

Operating Environment

If the transmitter has been stored in conditions such as low temperatures that are outside the usage environment range, place the transmitter in the usage environment for sufficient time for it to adjust to the room temperature before use. If the transmitter may be used constantly, store it within the usage environment range.

Appearance

- There are no damaged or dirty parts on the outside of the transmitter (LCD, keys, sockets, battery case cover, battery case, etc.).
- The transmitter is completely dry.
- The electrodes, electrode lead, SpO₂ probe and NIBP cuff are not broken.

Batteries

- The battery polarity is correct.
- The battery case spring is firmly attached and the battery is not loose.
- The battery case cover is firmly closed.

Channel Setting

- The transmitter channel matches the receiving monitor channel.
- There is no nearby transmitter with the same channel.

Other

The transmitter information is correctly sent to the receiving monitor when the transmitter is turned on or off.

Check Items After Power On

After turning on the power, check the following.

Power On

- The transmitter generates a one second “peep” sound and the startup screen appears.
- The transmitter displays the check electrodes screen.
- The transmitter is not too hot.
- The transmitter does not display the “BATTERY WEAK” message.
- The transmitter does not interfere with the operation of other medical instruments.

Daily Check

- The “signal loss” message is not displayed on the receiving monitor when the transmitter is inside the receiving range of the monitor.
- The battery replacement message is not displayed on the monitor.
- The keys on the transmitter function properly.
- The LCD brightness is appropriate. To adjust brightness, refer to the “Changing SYSTEM SETUP Settings” section.

Check Items After Use

To use the transmitter in safe and optimum condition for next time, check the following.

Before Turning Power Off

- Temporarily changed settings are changed back to the previous settings.
- There was no malfunction on the transmitter.

Storage

- ECG electrode leads, SpO₂ probe and NIBP cuff are cleaned and disinfected.
- If the transmitter got wet, liquid is wiped off and the transmitter is thoroughly dried.
- There are enough consumables, such as disposable electrodes.
- The transmitter power is turned off by removing batteries from the transmitter.
- Dead batteries are disposed of properly.

Turning Off the Transmitter

To turn off the power, remove the batteries. When the power is turned off, the saved waveform and numeric data are deleted.

Changing the Transmitter Channel

The channel of the transmitter can be changed with an optional QI-901PK channel writer.


⚠ WARNING

The following actions must be taken to properly receive the transmitter signal of the correct patient on the receiving instrument. Otherwise, there may be signal loss or signals may mix causing a serious accident, such as monitoring a different patient.

- Assign a channel administrator in the hospital and only he or she should manage channel assignment.
- The channel administrator must manage the channels in the facility so that there is no signal interference.
- When the transmitter channel is changed, the channel administrator must check that the channel on the receiving monitor is also changed and the signal is properly received.
- The channel administrator must replace the channel number label on the transmitter with the new one after changing the channel.

- NOTE
- The software version of the QI-901PK channel writer must be 02-01 or later to change the channel on the transmitter.
 - The channel writer must be used outside the patient environment.

The channel is displayed in the upper left corner of the screen.



Changing Parameter and System Setup Settings

The initial settings on the PARAMETER SETUP and SYSTEM SETUP screens can only be changed before monitoring. Changing these settings during monitoring interrupts monitoring.

NOTE: Changing Parameter and System Setup settings must be done by qualified personnel.

Notes on Parameter Settings

When monitoring NIBP and SpO₂, the following setting must be set as indicated in the table to properly transmit the monitoring data to the receiving monitor. Otherwise, SpO₂ cannot be monitored properly during NIBP measurement.

Some receiving monitors require the software to be upgraded. For details, contact your Nihon Kohden representative.

| SpO ₂ Probe Attachment Site | INHIBIT SpO ₂ DURING NIBP SETTING |
|----------------------------------------------------------|----------------------------------------------|
| Probe attached to the same limb as the cuff ¹ | ON |
| Probe attached to the limb without cuff | OFF |

¹ When the SpO₂ probe is attached to the same limb as the NIBP cuff and the cuff is inflated, the SpO₂ value becomes unstable and SpO₂ or PR alarm may occur.

Changing PARAMETER SETUP Settings


Parameter Setup Setting List

The factory default settings are underlined.

| Setting Item | Description | Settings |
|--------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------|
| ECG ELECTRODES | Select the electrode lead type. | IEC, <u>AHA</u> |
| LEAD TYPE | Select the type of ECG leads. | <u>AUTO</u> , 6 LEADS |
| ECG MEASUREMENT | Turn ECG monitoring on or off. When electrodes are attached to the patient and ECG leads are connected, ECG monitoring starts even when this setting is set to OFF. If this setting is set to OFF, the same setting on the receiving monitor must also be set to OFF. NOTE: When "PROTOCOL" of the transmitter is set to 57 and the receiving monitor is able to receive protocol 57, ECG measurement on the receiving monitor is automatically set to OFF when this setting is set to OFF on the transmitter. | <u>ON</u> , OFF |
| RESP MEASUREMENT | Turn respiration monitoring on or off. When this setting is set to OFF, the same setting on the receiving monitor is automatically set to OFF. | <u>ON</u> , OFF |
| SpO ₂ RESPONSE | Select the SpO ₂ response mode. | FAST, NORMAL, <u>SLOW</u> |
| INHIBIT SpO ₂ DURING NIBP | Turn SpO ₂ monitoring on or off during NIBP measurement. | <u>ON</u> , OFF |
| SELECTABLE INTERVALS (min) | Select the NIBP measurement modes for the mode selection. | <u>MANUAL</u> , STAT, <u>5</u> , <u>10</u> , <u>15</u> , <u>30</u> , <u>60</u> , 120, 240 |
| INITIAL INTERVAL (min) | Select the initial NIBP measurement mode at power on. | <u>MANUAL</u> , STAT, 5, 10, 15, 30, 60, 120, 240 |
| NIBP MODE AFTER STAT (min) | Select the NIBP measurement mode after completing STAT measurement. | MANUAL, <u>5</u> , 10, 15, 30 |
| START/FINISH SOUND | Turn ON or OFF the sound for NIBP measurement start/finish. | START: ON, <u>OFF</u> FINISH: ON, <u>OFF</u> |
| OLD NIBP DATA/AFTER (min) | Select whether to hide or dim the NIBP data after measurement and how long to wait after measurement to dim or hide it. | DATA: HIDE, <u>DIM</u> AFTER: 5, 10, <u>30</u> |
| INITIAL CUFF PRESSURE (mmHg) | Select the NIBP cuff inflation pressure. | 120, 150, <u>180</u> , 210, 240 |

Displaying the PARAMETER SETUP Screen

Cursor




Cursor

MENU screen

Setting item

Setting




PARAMETER SETUP screen - page 1

- 1 Turn off the transmitter by removing one battery.
- 2 While pressing the Function key, turn on the transmitter by inserting the battery. The MENU screen appears.
- 3 Press the ► key to move the cursor to “PARAMETER SETUP”.
- 4 Press the Function key to enter PARAMETER SETUP. The current settings are highlighted.
- 5 Change settings.
 - To move the cursor and select the setting item, press the ◀ or ▶ key then press the Function key.
 - To select and register the setting, press the ◀ or ▶ key then press the Function key.
 - To cancel changing the setting of the selected item, press the Screen key.
- 6 When changing settings on the PARAMETER SETUP screen is complete, press the Screen key to return to the MENU screen.
- 7 Press the ◀ or ▶ key to move the cursor to “EXIT”.
- 8 Press the Function key. The numeric and waveform screen appears.


Changing Parameter Setup Settings

ECG ELECTRODES

Select the electrode lead type.



- 1 On the PARAMETER SETUP screen, press the ► key to move the cursor to “ECG ELECTRODES”.
- 2 Press the Function key. The cursor moves to the selection item.
- 3 Press the ► key to select “IEC” or “AHA”.
- 4 Press the Function key to register the selected setting. The cursor returns to “ECG ELECTRODES”.



LEAD TYPE

Select the type of ECG leads. In normal use, select “AUTO”. When using DIN type lead with 6 electrodes, select “6 LEADS”.

- 1 On the PARAMETER SETUP screen, press the ► key to move the cursor to “LEAD TYPE”.
- 2 Press the Function key.
- 3 Press the ► key to select “AUTO” or “6 LEADS”.
- 4 Press the Function key to register the selected setting. The cursor returns to “LEAD TYPE”.

ECG MEASUREMENT

Turn ECG monitoring on or off. When electrodes are attached to the patient and ECG leads are connected, ECG monitoring starts even when this setting is set to OFF.

If this setting is set to OFF, the same setting on the receiving monitor must also be set to OFF.

NOTE: When "PROTOCOL" of the transmitter is set to 57 and the receiving monitor is able to receive protocol 57, ECG measurement on the receiving monitor is automatically set to OFF when this setting is set to OFF on the transmitter.

- 1** On the PARAMETER SETUP screen, press the ► key to move the cursor to "ECG MEASUREMENT".
- 2** Press the Function key.
- 3** Press the ► key to select "ON" or "OFF".
- 4** Press the Function key to register the selected setting. The cursor returns to "ECG MEASUREMENT".

RESP MEASUREMENT


Turn respiration monitoring on or off. When this setting is set to OFF, the same setting on the receiving monitor is automatically set to OFF.

- 1** On the PARAMETER SETUP screen, press the ► key to move the cursor to "RESP MEASUREMENT".
- 2** Press the Function key.
- 3** Press the ► key to select "ON" or "OFF".
- 4** Press the Function key to register the selected setting. The cursor returns to "RESP MEASUREMENT".

SpO₂ RESPONSE

Select the SpO₂ response mode from FAST, NORMAL or SLOW. For details on the response time, refer to the “Specifications - SpO₂ Measurement (ISO 9919: 2005 compliant)” section in this manual.

NOTE: When measurement condition is unstable due to strenuous movement of the patient, etc., response may become slower in all modes.



- 1 On the PARAMETER SETUP screen, press the ► key to move the cursor to “SpO₂ RESPONSE”. “SpO₂ RESPONSE” is on the second page of the PARAMETER SETUP screen.
- 2 Press the Function key.
- 3 Press the ► key to select “FAST”, “NORMAL” or “SLOW”.
- 4 Press the Function key to register the selected setting. The cursor returns to “SpO₂ RESPONSE”.

PARAMETER SETUP screen - page 2

INHIBIT SpO₂ DURING NIBP

Turn SpO₂ monitoring on or off during NIBP measurement.

When the SpO₂ probe is attached to the same limb as the NIBP cuff and this setting is set to OFF, the pulse may become unstable and SpO₂ or PR alarm may occur. It is recommended to set this setting to ON so that SpO₂ is not measured during NIBP measurement.


When the SpO₂ probe is attached to the other limb from the NIBP cuff, this setting can be set to OFF.

NOTE: When this “INHIBIT SpO₂ DURING NIBP” is set to OFF, refer to the “Monitoring SpO₂ during NIBP Measurement” section.

- 1 On the PARAMETER SETUP screen, press the ► key to move the cursor to “INHIBIT SpO₂ DURING NIBP”. “INHIBIT SpO₂ DURING NIBP” is on the second page of the PARAMETER SETUP screen.
- 2 Press the Function key.
- 3 Press the ► key to select “ON” or “OFF”.
- 4 Press the Function key to register the selected setting. The cursor returns to “INHIBIT SpO₂ DURING NIBP”.

SELECTABLE INTERVALS (min)

When the NIBP INTERVAL key is pressed, the measurement mode changes according to the modes selected in this item. MANUAL mode is already selected for the mode selection.



PARAMETER SETUP screen - page 3

- 1 On the PARAMETER SETUP screen, press the ► key to move the cursor to “SELECTABLE INTERVALS”. “SELECTABLE INTERVALS” is on the third page of the PARAMETER SETUP screen.
- 2 Press the Function key.
- 3 Press the ► key to select or unselect the mode. The selected modes are highlighted.
- 4 Press the Function key to register the selected setting. The cursor returns to “SELECTABLE INTERVALS”.

INITIAL INTERVAL (min)

Select the initial NIBP measurement mode at power on.

- 1 On the PARAMETER SETUP screen, press the ► key to move the cursor to “INITIAL INTERVAL”. “INITIAL INTERVAL” is on the third page of the PARAMETER SETUP screen.
- 2 Press the Function key.
- 3 Press the ► key to select the mode. Only the mode or interval selected for “SELECTABLE INTERVALS” are available.
- 4 Press the Function key to register the selected setting. The cursor returns to “INITIAL INTERVAL”.


NIBP MODE AFTER STAT (min)

Select the NIBP measurement mode after completing the STAT measurement.

- 1 On the PARAMETER SETUP screen, press the ► key to move the cursor to “NIBP MODE AFTER STAT”. “NIBP MODE AFTER STAT” is on the third page of the PARAMETER SETUP screen.
- 2 Press the Function key.
- 3 Press the ► key to select the mode. Only the mode or interval selected for “SELECTABLE INTERVALS” are available.
- 4 Press the Function key to register the selected setting. The cursor returns to “NIBP MODE AFTER STAT”.

START/FINISH SOUND

Turn on or off the sound for NIBP measurement start and finish.



PARAMETER SETUP screen - page 4

- 1 On the PARAMETER SETUP screen, press the ► key to move the cursor to “START/FINISH SOUND”. “START/FINISH SOUND” is on the fourth page of the PARAMETER SETUP screen.
- 2 Press the Function key. The cursor moves to “START”.
- 3 Press the ► key to turn “ON” or “OFF” the sound for NIBP measurement start.
- 4 Press the Function key to register the setting for “START”. The cursor moves to “FINISH”.
- 5 Press the ► key to turn “ON” or “OFF” the sound for NIBP measurement finish.
- 6 Press the Function key to register the selected setting. The cursor returns to “START/FINISH SOUND”.

OLD NIBP DATA/AFTER (min)

Select whether to dim or hide the NIBP data after measurement and how long to wait after NIBP measurement to dim or hide it.

- 1 On the PARAMETER SETUP screen, press the ► key to move the cursor to “OLD NIBP DATA/AFTER”. “OLD NIBP DATA/AFTER” is on the fourth page of the PARAMETER SETUP screen.
- 2 Press the Function key. The cursor moves to “DATA”.
- 3 Press the ► key to select “HIDE” or “DIM” the NIBP data.
- 4 Press the Function key to register the setting for “DATA”. The cursor moves to “AFTER”.
- 5 Press the ► key to select the interval after NIBP measurement to dim or hide.
- 6 Press the Function key to register the selected setting. The cursor returns to “OLD NIBP DATA/AFTER”.

INITIAL CUFF PRESSURE (mmHg)

Select the NIBP cuff inflation pressure.

- 1** On the PARAMETER SETUP screen, press the ► key to move the cursor to “INITIAL CUFF PRESSURE”. “INITIAL CUFF PRESSURE” is on the fourth page of the PARAMETER SETUP screen.
- 2** Press the Function key.
- 3** Press the ► key to select the inflation pressure.
- 4** Press the Function key to register the selected setting. The cursor returns to “INITIAL CUFF PRESSURE”.


Changing SYSTEM SETUP Settings

System Setup Setting List


The factory default settings are underlined.

| Setting Item | Description | Settings |
|-----------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------|
| PROTOCOL | <p>Select the transmitting protocol.</p> <p>57: New protocol. A central monitor with an ORG-9100A, ORG-9110A or ORG-9700A multiple patient receiver whose software version 02-01 or later can receive this protocol.</p> <p>42: Old protocol. A central monitor with an ORG-9100A, ORG-9110A or ORG-9700A multiple patient receiver can receive this protocol.</p> <p>NOTE: When 57 is set, the receiving monitor must be able to receive protocol 57. Otherwise, signals from the transmitter cannot be received.</p> | <u>57</u> , 42 |
| BRIGHTNESS | Select the screen brightness. | <u>DARK</u> , BRIGHT |
| PRESSURE UNIT | Select the unit for NIBP. | <u>mmHg</u> , kPa |
| FUNCTION KEY | <p>Select the function of the Function key.</p> <p>SUSPEND ALARM & PAUSE: Suspends alarm on the receiving monitor for 2 minutes. Pauses monitoring on the transmitter and receiving monitor.</p> <p>SUSPEND ALARM: Suspends alarm on the receiving monitor for 2 minutes.</p> <p>CONFIRM: Displays the "PATIENT CONFIRMED" message on the transmitter screen and transmits the message to the receiving monitor.</p> <p>OFF: No function.</p> <p>NOTE: "SUSPEND ALARM & PAUSE" and "CONFIRM" can only be set when PROTOCOL is set to 57.</p> | <u>SUSPEND ALARM & PAUSE</u> , SUSPEND ALARM, CONFIRM, OFF |
| AUTO RESUME AFTER PAUSE | Select the interval to resume monitoring after PAUSE. | 10 s, <u>30 s</u> , 1 min, 2 min, 3 min |
| SELECTABLE SCREEN TIME OUT PERIOD (min) | Select the display time-out period. | 5, 10, 15, <u>30</u> , 60, 120, 240 |

Displaying the SYSTEM SETUP Screen



MENU screen



SYSTEM SETUP screen - page 1

- 1 Turn off the transmitter by removing one battery.
- 2 While pressing the Function key, turn on the transmitter by inserting the battery. The MENU screen appears.
- 3 Press the ► key to move the cursor to “SYSTEM SETUP”.
- 4 Press the Function key to enter SYSTEM SETUP. The current settings are highlighted.
- 5 Change settings.
 - To move the cursor and select the setting item, press the ◀ or ▶ key then press the Function key.
 - To select and register the setting, press the ◀ or ▶ key then press the Function key.
 - To cancel changing the setting of the selected item, press the Screen key.

The SYSTEM SETUP screen has two pages. To display the second page, press the ► key when the cursor is at “PRESSURE UNIT”.
- 6 When changing settings on the SYSTEM SETUP screen is complete, press the Screen key to return to the MENU screen.
- 7 Press the ◀ or ▶ key to move the cursor to “EXIT”.
- 8 Press the Function key. The numeric and waveform screen appears.

Changing System Setup Settings

PROTOCOL

Select the transmitting protocol. For differences between protocols, refer to the table below.


57: New protocol. A central monitor with an ORG-9100A, ORG-9110A or ORG-9700A multiple patient receiver whose software version 02-01 or later can receive this protocol.

42: Old protocol. A central monitor with ORG-9100A, ORG-9110A or ORG-9700A multiple patient receiver can receive this protocol.


NOTE: When 57 is set, the receiving monitor must be able to receive protocol 57. Otherwise, signals from the transmitter cannot be received.

Differences Between Protocols

| Function | Protocol 42 | Protocol 57 |
|--------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------|--------------|
| Setting ECG MEASUREMENT to OFF on the transmitter automatically turns off the ECG measurement setting on the receiving monitor | No (ECG measurement must be turned off on the receiving monitor) | Yes |
| Pause monitoring on the receiving monitor from the transmitter | No | Yes |
| Transmit "PATIENT CONFIRMED" message | No | Yes |
| Display battery level of the transmitter on the receiving monitor | No | Yes |
| Transmit SpO ₂ messages | Some messages (refer to the "Indication and Message List" section) | All messages |



- 1 Press the ► key to move the cursor to "PROTOCOL".
- 2 Press the Function key.
- 3 Press the ► key to select "42" or "57".



NOTE: FUNCTION KEY (on the second page of the SYSTEM SETUP screen) can be set to "SUSPEND ALARM & PAUSE" or "CONFIRM" only when PROTOCOL is "57". If PROTOCOL is changed to "42", FUNCTION KEY is automatically changed to "OFF".

- 4 Press the Function key to register the selected setting. The cursor returns to "PROTOCOL".

BRIGHTNESS

Select the screen brightness.

- 1 Press the ► key to move the cursor to "BRIGHTNESS".
- 2 Press the Function key.
- 3 Press the ► key to select "DARK" or "BRIGHT".
- 4 Press the Function key to register the selected setting. The cursor returns to "BRIGHTNESS".

PRESSURE UNIT

Select the unit for NIBP.

- 1 Press the ► key to move the cursor to "PRESSURE UNIT".
- 2 Press the Function key.
- 3 Press the ► key to select "mmHg" or "kPa".
- 4 Press the Function key to register the selected setting. The cursor returns to "PRESSURE UNIT".

FUNCTION KEY

Select the function of the Function key. For details on using these functions, refer to “Basic Monitoring Operation” in the “Monitoring” section.


SUSPEND ALARM & PAUSE: Suspends alarm on the receiving monitor for 2 minutes. Pauses monitoring on the transmitter and receiving monitor.

SUSPEND ALARM: Suspends alarm on the receiving monitor for 2 minutes.

CONFIRM: Displays the “PATIENT CONFIRMED” message on the transmitter screen and transmits the message to the receiving monitor.

OFF: No function.

NOTE: “SUSPEND ALARM & PAUSE” and “CONFIRM” can only be set when PROTOCOL is set to 57.



SYSTEM SETUP screen - page 2

- 1 On the SYSTEM SETUP screen, press the ► key to move the cursor to “FUNCTION KEY”. “FUNCTION KEY” is on the second page of the SYSTEM SETUP screen.
- 2 Press the Function key.
- 3 Press the ► key to select the function.
- 4 Press the Function key to register the selected setting. The cursor returns to “FUNCTION KEY”.

AUTO RESUME AFTER PAUSE


Select the interval to resume monitoring after PAUSE. When either of the following conditions is met, monitoring resumes on the receiving monitor.

- Heart rate is properly monitored for the selected interval.
- SpO₂ is properly monitored for the selected interval.
- NIBP is properly measured and the SYS, DIA or MAP value is displayed.

- 1 Press the ► key to move the cursor to “AUTO RESUME AFTER PAUSE”. “AUTO RESUME AFTER PAUSE” is on the second page of the SYSTEM SETUP screen.
- 2 Press the Function key.
- 3 Press the ► key to select the interval.
- 4 Press the Function key to register the selected setting. The cursor returns to “AUTO RESUME AFTER PAUSE”.

SELECTABLE SCREEN TIME OUT PERIOD (min)

Select the display time-out period. If no key is pressed for the selected time, the display is automatically turned off. The selected time is shown on the Select Screen Time Out Period screen. Refer to the “Turning the Display Off” section for details.




SYSTEM SETUP screen - page 2

- 1** Press the ► key to move the cursor to “SELECTABLE SCREEN TIME OUT PERIOD (min)”. It is on the third page of the SYSTEM SETUP screen.
- 2** Press the Function key.
- 3** Press the ◀ or ▶ key to select time-out period.
- 4** Press the Function key to register the selected setting. The cursor returns to “SELECTABLE SCREEN TIME OUT PERIOD (min)”.

Initializing Settings

Do the following procedure to initialize all settings, except for channel, to the factory default settings.


Cursor



1 Turn off the transmitter by removing a battery.

2 While pressing the Function key, turn on the transmitter by inserting the battery. The MENU screen appears.


3 Press the ► key to move the cursor to “SYSTEM INITIALIZE”.



4 Press the Function key to enter the SYSTEM INITIALIZE screen.

5 Press the Function key. A confirmation message appears.

To return to the MENU screen without initializing, press the Screen key.



6 Press the Function key to initialize settings. To cancel initializing, press the ► key. The screen returns to the MENU screen.

Attaching Electrodes, SpO₂ Probe and NIBP Cuff to the Patient

The transmitter can be placed on the bedside. The required length of the electrode leads and SpO₂ probe cable depends on how the transmitter is to be attached to the patient.


NOTE: Monitoring SpO₂ during NIBP Measurement

When the SpO₂ probe is attached to the same limb as the NIBP cuff, the blood flow decreases during NIBP measurement and pulse wave cannot be detected and SpO₂ cannot be monitored properly. When "INHIBIT SpO₂ DURING NIBP" on the PARAMETER SETUP screen is set to ON (factory default setting), SpO₂ monitoring is paused during NIBP measurement to avoid SpO₂ alarm occurrence. However, when monitoring SpO₂ on the same limb as NIBP, be careful when reading SpO₂ values.

When monitoring SpO₂ is important, attach the probe to the limb to which the NIBP cuff or catheter is not attached.

Attachment Example

When transmitter is placed on a bed



NOTE: When placing the transmitter on a bedside, place it on a stable and flat place. If the transmitter falls off, it may be damaged.


Selecting Electrode Leads

CAUTION

Only use Nihon Kohden specified electrodes and electrode leads. When other type of electrodes or electrode leads are used, the "CHECK ELECTRODES" message may be displayed and ECG monitoring may stop.

The following optional electrode leads can be used with the transmitter. Refer to the "Options" section when ordering the electrode leads.

BR-903PA (3 electrodes, AHA, clip type), BR-933PA (3 electrodes, AHA, clip type),
BR-943PA (3 electrodes, AHA, hook type), BR-936PA (6 electrodes, AHA, clip type),
BR-946PA (6 electrodes, AHA, hook type)




Example: BR-903PA

Connecting the Electrode Lead to the Transmitter

Connect the electrode lead to the ECG/RESP socket on the transmitter.

NOTE: Hold the sides of the transmitter when connecting or disconnecting the electrode lead. Pressing the display while holding the transmitter may damage the transmitter.



⚠ CAUTION

Do not shake or swing the transmitter while holding the leads or cables connected to the transmitter. The transmitter may come off and injure someone or damage surrounding instruments.

⚠ CAUTION

Hold the connector of the electrode lead when connecting/disconnecting the electrode lead. If you disconnect the electrode lead by pulling the lead, it damages the electrode lead.

Electrode Position

Follow the physician's instructions for electrode placement when available.

For ECG monitoring, electrodes are attached only on the chest to allow patient movement and obtain continuous stable ECG. The following leads are examples. When also monitoring respiration, refer to the "Electrode Position for Respiration Monitoring" section.


NOTE: The optimum electrode positions for ECG measurement are not always optimum for respiration measurement. Select positions that are suitable for both ECG and respiration measurement or positions which give priority to either ECG or respiration measurement.

Electrode Positions for ECG Monitoring

Six Electrodes

The 6-electrode method with lead II and lead V5 is effective for monitoring myocardial ischemia.

You can improve monitoring accuracy considerably by adding lead V4 to this combination. Va and Vb can be at any position of the standard 12 leads V1 to V6, but V4 and V5 are most appropriate for myocardial ischemic monitoring.




| Electrode Position | Symbol | | Lead Color | |
|-----------------------------------------------------------------------------------------------|--------|-----|--------------|-------------|
| | AHA | IEC | AHA | IEC |
| Left infraclavicular fossa | LA | L | Black | Yellow |
| Right infraclavicular fossa | RA | R | White | Red |
| Below lowest rib on the left anterior axillary line | LL | F | Red | Green |
| Right anterior axillary line at the same level as LL/F | RL | RF | Green | Black |
| Fifth intercostal space on the left midclavicular line. (V4 position of standard 12 leads) | Va | Ca | Brown-blue | White-brown |
| Left anterior axillary line at the same level as Va. (V5 position of standard 12 leads) | Vb | Cb | Brown-orange | White-black |


Lead Position

Standard limb leads


Lead I



Lead II




Lead III




Monopolar limb leads


aV_R lead



aV_L lead




aV_F lead




Monopolar chest leads

V₁ to V₆ leads





Three Electrodes

- Lead MII, which is similar to standard lead II, used when ECG measurement has priority



| Electrode Position | Symbol | | Lead Color | |
|-----------------------------------------------------|--------|-----|------------|--------|
| | AHA | IEC | AHA | IEC |
| Left infraclavicular fossa | LA | L | Black | Yellow |
| Right infraclavicular fossa | RA | R | White | Red |
| Below lowest rib on the left anterior axillary line | LL | F | Red | Green |

- Lead MI, which is similar to standard lead I
- Change F/LL and L/LA of lead MII
- Lead MIII, which is similar to standard lead III
- Change R/RA and L/LA of lead MII.




If the electrode positions above are not available due to chest surgery, attach the electrodes to the root of the limbs or below the clavicles for stable ECG monitoring.

Electrode Positions for Respiration Monitoring

Place the R/RA and F/LL electrodes so that the lungs are between the electrodes.

Position 1


In this position, respiration measurement is available; however, there is a difference in amplitude between different patients.



| R or RA | F or LL |
|-----------------------------|------------------------------------------------------------|
| Right infraclavicular fossa | Fifth intercostal space on the left midclavicular line, V4 |

Position 2


In this position, the waveform amplitude is usually large and the ECG lead is similar to Lead MII. This position can be generally recommended.



| R or RA | F or LL |
|-----------------------------|----------------------------------------------------------|
| Right infraclavicular fossa | Fifth intercostal space on the left midaxillary line, V6 |

Position 3


In this position, the respiration waveform is optimum, but the ECG lead is unusual.



| R or RA | F or LL |
|-------------------------------------------------|----------------------------------------------------------|
| Right midaxillary at the horizontal level of V6 | Fifth intercostal space on the left midaxillary line, V6 |

Position 4

In this position, the respiration measurement is influenced by the impedance variation of the abdomen, so the cardiac pulse wave included in the respiration wave is reduced. Note that the waveform is inverted in phase compared with the chest movement (the waveform goes down during inspiration). It is difficult to measure the ECG at the same time.



| R or RA | F or LL |
|------------------------------------------------|-----------------------------------------------|
| Lowest rib on the right anterior axillary line | Lowest rib on the left anterior axillary line |

Attaching Electrodes to the Patient and Connecting the Electrode Leads to Disposable Electrodes

Prepare the Patient Skin

Shave off excessive body hair.

To reduce skin impedance, clean the electrode site with cream or with a cotton pad moistened with alcohol. Thoroughly dry the skin with a clean cotton pad.

- NOTE**
- For a patient with frequent body movement, rub the sites with Skinpure skin preparation gel. However, do not use Skinpure skin preparation gel on sensitive skin.
 - Do not place electrodes on a wound or on an inflamed, wrinkled or uneven skin surface.


Attaching Electrodes to the Patient

CAUTION


Do not reuse disposable parts and accessories.

- NOTE**
- To maintain good contact between the electrode and skin, check that the paste of the disposable electrode is not dry.
 - When contact between the disposable electrode and skin becomes poor, replace electrodes with new ones immediately. Otherwise, contact impedance between the skin and the electrode increases and accurate ECG cannot be obtained.

Refer to the electrode operator's manual for details.



- 1** Carefully remove the backing paper from the electrode. Avoid touching the adhesive surface.
- 2** Place the electrode on the previously cleaned skin. Pay attention to the electrode lead color and symbol.
- 3** Clip the electrode lead to the electrode.
- 4** Fasten the electrode lead wire with surgical tape with an extra length of wire between the tape and the electrode. This lessens the movement of electrode leads by body movement and helps stable monitoring.



Checking ECG on the Transmitter Screen

After attaching electrodes and connecting ECG leads, check that the electrodes are properly attached to the patient and the ECG waveform is acquired on the check electrodes screen. For details on the screen, refer to the "Screen Descriptions" section.

Attaching the SpO₂ Probe

Selecting the SpO₂ Probe

Select an appropriate probe for the patient.

⚠ CAUTION




Only use Nihon Kohden specified electrodes, electrode leads and SpO₂ probes. Otherwise, the maximum performance from the transmitter cannot be guaranteed.

⚠ CAUTION

Do not use a damaged or disassembled probe. It causes incorrect measurement and may injure the patient.

Reusable Probes




When using a TL-201T finger probe, choose the appropriate cable length for attachment.




| Probe | Cable Length | Patient | Attachment Site |
|-------------------------------------------------------------------------------------------------------------------------------------------|--------------------|---------------------------------|-----------------|
| Finger probe TL-201T  | 0.6 m | Adult or child 20 kg or more | Finger |
| | 1.6 m | | |
| Finger probe TL-631T1, TL-631T3  Attachment tape | TL-631T1: 0.6 m | Adult or child 20 kg or more | Finger or toe |
| | TL-631T3: 1.6 m | | |
| Multi-site probe TL-220T  Attachment tape | 1.6 m | Adult or infant 3 kg or more | Finger or toe |
| | | Neonate 3 kg or less | Instep and sole |

Disposable Probes

⚠ CAUTION

The disposable probe is not sterilized. Use the disposable probe only for a single patient. Never reuse the disposable probe for another patient because it causes cross infection.


| Probe | Patient | Attachment Site |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------|------------------------|
| <p>SpO₂ probe TL-051S, TL-052S</p>  <p style="text-align: right;">40 mm</p> <p style="text-align: center;">Attachment tape</p> <p>Cable length TL-051S: 0.8 m TL-052S: 1.6 m</p> | <p>Adult 50 kg or more</p> | <p>Finger</p> |
| | <p>Neonate 3 kg or less</p> | <p>Instep and sole</p> |
| <p>SpO₂ probe TL-061S, TL-062S</p>  <p style="text-align: right;">35 mm</p> <p style="text-align: center;">Attachment tape</p> <p>Cable length TL-061S: 0.8 m TL-062S: 1.6 m</p> | <p>Adult or child 15 to 50 kg</p> | <p>Finger</p> |
| | <p>Infant 3 to 15 kg</p> | <p>Toe</p> |
| <p>SpO₂ probe TL-271T, TL-271T3</p>  <p>Cable length TL-271T: 0.8 m TL-271T3: 1.6 m</p> | <p>Adult 30 kg or more</p> | <p>Finger or toe</p> |

| Probe | Patient | Attachment Site |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------|------------------------|
| <p>SpO₂ probe TL-272T, TL-272T3</p>  <p>Cable length TL-272T: 0.8 m TL-272T3: 1.6 m</p> | <p>Child 10 to 50 kg</p> | <p>Finger or toe</p> |
| <p>SpO₂ probe TL-273T, TL-273T3</p>  <p>Cable length TL-273T: 0.8 m TL-273T3: 1.6 m</p> | <p>Neonate 3 kg or less</p> | <p>Instep and sole</p> |
| <p>SpO₂ probe TL-274T, TL-274T3</p>  <p>Cable length TL-274T: 0.8 m TL-274T3: 1.6 m</p> | <p>Infant 3 to 20 kg</p> | <p>Finger or toe</p> |

Connecting the SpO₂ Probe to the Transmitter

Connect the probe to the SpO₂ socket on the transmitter.

NOTE: Hold the connector of the probe cable when connecting or disconnecting the probe. Pulling the cable may damage the probe.



⚠ CAUTION

Do not shake or swing the transmitter while holding the leads or cables connected to the transmitter. The transmitter may come off and injure someone or damage surrounding instruments.

⚠ CAUTION

Do not use a damaged or disassembled probe. It causes incorrect measurement and may injure the patient.

Attaching the Probe to the Patient

Attach the probe to the patient by referring to the probe's manual. Make sure that the light emitter and photo detector of the probe face each other at the attachment site.

WARNING

- When using the TL-201T finger probe, do not fasten the probe and cable to the finger by wrapping with tape. This may cause burn, congestion or pressure necrosis from poor blood circulation.
- When using probes other than the TL-201T finger probe, to avoid poor circulation, do not wrap the tape too tight. Check the blood circulation condition by observing the skin color and congestion at the skin peripheral to the probe attachment site. Even for short-term monitoring, there may be burn or pressure necrosis from poor blood circulation, especially on neonates or low birth weight infants whose skin is delicate. Accurate measurement cannot be performed on a site with poor peripheral circulation.

WARNING

Check the circulation condition by observing the skin color at the measurement site and pulse waveform. Change the measurement site every 8 hours for disposable probes and every 4 hours for reusable probes (every 8 hours for TL-631T series probe). The skin temperature may increase at the attached site by 2 or 3°C (4 or 5°F) and cause a burn or pressure necrosis. When using the probe on the following patients, take extreme care and change the measurement site more frequently according to symptoms and degree.

- Patient with a fever
- Patient with insufficient peripheral circulation
- Neonate or low birth weight infant with delicate skin

⚠ WARNING

When monitoring SpO₂ of a patient who is receiving photodynamic therapy, the light from the finger probe sensor may cause a burn. Photodynamic therapy uses a photosensitizing agent that has a side effect of photosensitivity.

The SpO₂ probes manufactured by Nihon Kohden have two wavelengths with peaks in the range of 650 and 950 nm. The maximum light intensity is less than 5.5 mW.

⚠ CAUTION

If the attachment site is dirty with blood or bodily fluids, clean the attachment site before attaching the probe. If there is nail polish on the attachment site, remove the polish. Otherwise, the amount of transmitted light decreases, and measured value may be incorrect or measurement cannot be performed.

⚠ CAUTION

When the probe is attached on an appropriate site with sufficient thickness and the error message confirming the probe attachment repeatedly appears, the probe may be deteriorated. Replace it with a new one.

⚠ CAUTION

If the skin gets irritated or redness appears on the skin from the probe, change the attachment site or stop using the probe. Take extreme care for the patients with delicate skin.

⚠ CAUTION

Do not use a probe which is deteriorated by aging. Accurate measurement cannot be performed.

Attaching the NIBP Cuff

Selecting the NIBP Cuff


Select the NIBP cuff appropriate for the patient.

NOTE: NIBP cannot be measured on neonates using this transmitter.

Reusable Cuffs

To use these cuffs, an optional YN-990P extension hose (1.5 m) is required.


| Reusable Cuff | | Model | Width (cm) | Air Hose Length (cm) |
|---------------------|----------|---------|------------|----------------------|
| YAWARA CUFF2 | | | | |
| For infant | | YP-710T | 5 | 15 |
| For child | | YP-711T | 7 | |
| For adult | Small | YP-712T | 10 | |
| | Standard | YP-713T | 13 | |
| | Large | YP-714T | 16 | |



Disposable Cuffs


⚠ CAUTION
 Disposable cuffs are not sterilized. If necessary, sterilize the cuff using glutaraldehyde solution.

| Disposable Cuff | | Model | Width (cm) | Air Hose Length (cm) |
|-----------------|--------------|---------|------------|----------------------|
| For infant | | YP-810P | 6 | 17 |
| For child | | YP-811P | 8 | |
| For adult | Small | YP-812P | 10 | |
| | Standard | YP-813P | 14 | |
| | Medium large | YP-814P | 15 | |
| | Large | YP-815P | 17 | |
| Extra large | | YP-816P | 18 | |



Extension Hose

⚠ CAUTION
 When using an extension hose, check that the extension hose is not bent or squeezed. Otherwise, the cuff might not inflate or deflate. If the cuff cannot deflate, it may cause congestion on the patient at the cuff attachment site.




YN-990P extension hose, 150 cm

Reference for selecting a cuff

The AHA (American Heart Association) recommends that the cuff width be 40% of the circumference of the upper arm. Refer to the following graph and select the cuff which suits the patient's arm.


- NOTE**
- If a range of arm circumference appropriate for the cuff is prescribed, use a cuff within that range.
 - To obtain accurate measured values, select a wide cuff which can be attached to the upper arm. Measuring with a very narrow cuff may result in measured values higher than the actual values.



Connecting the NIBP Cuff to the Transmitter

To use the YP-810P series disposable cuffs or YP-710T series reusable cuffs, an optional YN-990P extension hose (1.5 m) is required.

NOTE: Connect the joints properly. If there is an air leak, NIBP cannot be measured properly.



- 1 Connect the NIBP cuff to the extension hose.
- 2 Connect the other end of the extension hose to the NIBP socket on the transmitter. Turn the joint clockwise until it clicks.

To disconnect the cuff from the transmitter, turn the hose joint counterclockwise.

Attaching the NIBP Cuff to the Patient

⚠ WARNING

Be careful when measuring NIBP on a patient with known bleeding disorders or coagulation. After NIBP measurement, there may be dot hemorrhage, or circulatory disorder by thrombus where the cuff is attached.

⚠ CAUTION

Do not wrap the cuff on an arm or thigh which is used for injection. NIBP measurement on an arm or thigh which is used for injection may cause reflux of blood and stop injection.

⚠ CAUTION

Do not attach the cuff to the site where there is injury or inflammation. If the skin gets irritated or redness appears on the skin from the cuff, change the attachment site or stop using the cuff. Take extreme care on the patients with delicate skin.

⚠ CAUTION

Do not wrap the cuff too tight. It may cause poor blood circulation and congestion. If the cuff is wrapped too loosely, the NIBP value may increase.

⚠ CAUTION

Do not reuse disposable parts and accessories.

⚠ CAUTION

NIBP and SpO₂ can be measured on the same limb, but the SpO₂ monitoring might not be accurate during NIBP measurement. Be careful when reading the SpO₂ values.*


* Monitoring SpO₂ during NIBP Measurement

When the SpO₂ probe is attached to the same limb as the NIBP cuff, the blood flow decreases during NIBP measurement and pulse wave cannot be detected and SpO₂ cannot be monitored properly. When “INHIBIT SpO₂ DURING NIBP” on the PARAMETER SETUP screen is set to ON (factory default setting), SpO₂ monitoring is paused during NIBP measurement to avoid SpO₂ alarm occurrence. However, when monitoring SpO₂ on the same limb as the NIBP, be careful when reading SpO₂ values.

- NOTE**
- Measuring NIBP at a site other than the upper arm gives different values from those measured at the upper arm. When making diagnosis based on the NIBP values, measure NIBP on an upper arm.
 - To accurately detect the pulsatile flow of the artery, the cuff should be wrapped around a bare upper arm.
 - Do not use an abnormal cuff. The cuff deteriorates from use and cleaning. Before use, check the cuff and confirm that there is no flaw, crack or hole in it. Be careful not to damage the inflation bag. If the inflation bag has a hole or a flaw, it may burst during use. Dispose of an abnormal cuff and replace it with a new one.
 - Refer to the NIBP cuff manual for details.

Cuff Position


When placing the transmitter on a bed, make sure that the hose is not bent.




Place the cuffed upper arm (brachium) at the same height as the patient's heart. If the cuff is not at the same level as the heart, the weight of the blood affects the blood pressure reading. The pressure difference per unit height is 0.7 mmHg/cm. The blood pressure reading decreases when the arm is higher than the heart and increases when lower.

The best measuring condition is when the patient is lying on his/her back with arms and legs relaxed. If the cuff position cannot be on the same level as the heart, the displayed blood pressure reading must be mathematically adjusted.

Placing the Transmitter on the Bed



- 1 Put the cuff on the upper arm so that the ▼ mark of "ARTERY ▼" aligns with the artery of the patient.




- 2 Wrap the cuff so that "INDEX →" comes within the "← RANGE →".
If "Index →" is not within the "← RANGE →", change the cuff size.

Locking the Keys on the Transmitter

To prevent the patient from pressing the keys on the transmitter during monitoring, you can lock the keys.

To lock the keys:

Press the ◀ and ▶ keys at the same time and hold for more than 3 seconds. The “Key locked” screen appears.



When the screen time-out period is set to “1 min” (factory default), the “Key locked” screen is displayed for 5 seconds, then the display turns off if there is no key operation.

When the screen time-out period is set to a certain number of minutes, the “Key locked” screen is displayed for 5 seconds, then it changes to the numeric and waveform screen. If there is no key operation, the display turns off when the remaining time elapses. Refer to the “Turning the Display Off” section for details.

To unlock the keys:

Press the ◀ and ▶ keys at the same time and hold for more than 3 seconds.

Monitoring

⚠ CAUTION

The measurement values and displayed waveforms on the transmitter and receiving monitor may be different due to timing delay of the display or difference in detection settings.


- NOTE
- The waveforms displayed on the check electrodes screen are just to check that the electrode or probe is attached to the patient properly. Do not use the displayed waveforms for other purposes.
 - Do not let the transmitter directly touch the patient's skin. The transmitter temperature rises and this may cause burn to the patient.

Screen Descriptions

When the transmitter is turned on, the startup screen appears, then the check electrodes screen appears to check the electrode attachment.

The screen changes in the following order when the Screen key is pressed.

Check electrodes → numeric and waveform → waveform review → numeric review → display off → check electrodes . . .





Check Electrodes Screen

You can check whether the electrodes are properly attached to the patient and the ECG waveform is acquired.

When 6 leads are used, the I, II, Va and Vb lead waveforms are displayed.

When 3 leads are used, only the lead II waveform is displayed.




When electrodes are detached, the “CHECK ELECTRODES” message and detached electrode position appear on the screen.


NOTE: When ECG measurement is set to OFF on the PARAMETER SETUP screen, the check electrodes screen does not appear.

Numeric and Waveform Screen

Numeric values and waveforms of the monitoring parameters are displayed. You can change the ECG lead with the ◀ and ▶ keys.



When ECG and respiration measurement is turned off




NOTE: The pulse wave amplitude varies according to the ratio of the pulsation component to the entire transmitted IR signal. When the pulsation component ratio is 1%, the pulse wave amplitude is about 5 mm at x1 sensitivity on the screen.

Waveform Review Screen

ECG full disclosure for up to 10 minutes can be saved and reviewed. When ECG measurement is turned off and SpO₂ is monitored, the pulse waveform is saved.

When ECG lead is changed on the numeric and waveform screen, the ECG full disclosure of the changed lead is saved.

The saved data is deleted when the transmitter is turned off.



To scroll the waveform, press the ◀ or ▶ key. The waveform is scrolled by 30 seconds.

Numeric Review Screen

Numeric data of heart rate (or pulse rate when ECG is turned off), SpO₂ and respiration rate for up to 10 minutes are saved at 1 minute intervals.

NOTE: NIBP measured values are not saved.

The saved data is deleted when the transmitter is turned off.

| | HR (bpm) | SpO ₂ (%) | RR (/min) |
|------------|----------|----------------------|-----------|
| Older data | 80 | 97 | 20 |
| - 10[min] | 80 | 97 | 20 |
| - 9[min] | 80 | 97 | 20 |
| - 8[min] | 80 | 97 | 20 |
| - 7[min] | 80 | 97 | 20 |
| - 6[min] | 80 | 97 | 20 |
| - 5[min] | 80 | 97 | 20 |
| - 4[min] | 80 | 97 | 20 |
| - 3[min] | 80 | 97 | 20 |
| - 2[min] | 80 | 97 | 20 |
| Newer data | 80 | 97 | 20 |
| - 1[min] | 80 | 97 | 20 |

Display Off

The display can be turned off any time. Refer to the “Turning the Display Off” section for details.

Select Screen Time Out Period.


10-second time bar

1 min 30 min CONTINUE

RETURN LEFT RIGHT

Basic Monitoring Operation

Using the Function Key



One of the following functions can be assigned to the Function key on the SYSTEM SETUP screen. Refer to the “Changing SYSTEM SETUP Settings” section.

SUSPEND ALARM: Suspends alarms on the receiving monitor, before they occur, for 2 minutes.

PAUSE: Pauses monitoring on the transmitter and receiving monitor.

CONFIRM: Transmits the signal that the patient is confirmed and displays the “PATIENT CONFIRMED” message on the transmitter.

NOTE: To use the Function key for pause or confirm, “PROTOCOL” on the SYSTEM SETUP screen of the transmitter must be set to 57 and the receiving monitor must be able to receive protocol 57. A central monitor with an ORG-9100A, ORG-9110A or ORG-9700A multiple patient receiver whose software version 02-01 or later can receive this protocol.

Suspending Alarms on the Receiving Monitor


WARNING

While the “ALARMS SUSPENDED” message is displayed on the transmitter, all alarms on the receiving monitor are suspended so keep the patient under close observation.

When the FUNCTION KEY is set to “SUSPEND ALARM” or “SUSPEND ALARM & PAUSE” on the SYSTEM SETUP screen, alarms can be suspended for 2 minutes on the receiving monitor before they occur.

To suspend alarms:


- 1 Press the Function key. The “Suspend alarms” confirmation screen appears.



- 2 Press the ► key to suspend alarms.


To cancel suspending alarms and return to the previous screen, press the Screen key.

When the alarms are suspended, the “ALARMS SUSPENDED” message and alarm suspended icon with the remaining minutes in alarm suspension appear on the transmitter screen.



To cancel suspending alarms during 2 minute alarm suspension:

- 1 Press the Function key while the “ALARMS SUSPENDED” message is displayed. The confirmation screen appears.



- 2 Press the ► key to cancel alarm suspension.
Press the Screen key to not cancel alarm suspension.

Pausing Monitoring

When FUNCTION KEY is set to “SUSPEND ALARM & PAUSE” on the SYSTEM SETUP screen, you can pause monitoring on the receiving monitor from the transmitter when the patient cannot be monitored, such as during X-ray examination.


NOTE: To use the Function key for pause, “PROTOCOL” on the SYSTEM screen of the transmitter must be set to 57 and the receiving monitor must be able to receive protocol 57.

Monitoring on the receiving monitor resumes when one of the following conditions is met. When OFF is selected for “AUTO RESUME AFTER PAUSE” on the SYSTEM SETUP screen, monitoring does not automatically resume.

- Heart rate is properly monitored on the transmitter for the interval selected for “AUTO RESUME AFTER PAUSE”
- SpO₂ is properly monitored on the transmitter for the interval selected for “AUTO RESUME AFTER PAUSE”
- NIBP is properly measured and the SYS, DIA or MAP value is displayed on the transmitter

To pause monitoring:

- 1 Press the Function key. The “Suspend alarms” confirmation screen appears.



- 2 Press the Function key for 3 seconds to display the “Pause monitoring” confirmation screen.



- 3 Press the ► key to pause monitoring.
To cancel pause monitoring, press the Screen key.
- 4 Wait about 5 seconds until the “Turn power off” screen appears.



5 Turn off the transmitter.

If the transmitter is not turned off and monitoring continues for the interval set for “AUTO RESUME AFTER PAUSE” on the SYSTEM SETUP screen, pause monitoring is cancelled and monitoring continues.

Resuming Monitoring after Pause

To resume monitoring after pause, check that the electrodes, electrode leads and probe are attached to the patient then turn on the transmitter.

⚠ WARNING


When the patient returns to the bed, turn on the transmitter and check that the monitoring is resumed on the receiving monitor.

⚠ WARNING

If the transmitter is not turned off and monitoring continues for the selected interval, pause monitoring is canceled and monitoring continues. Check that the monitoring is resumed on the receiving monitor.

Confirming the Patient


When the FUNCTION KEY is set to “CONFIRM” on the SYSTEM SETUP screen, you can transmit signal to the receiving monitor to indicate that the patient is confirmed by a medical staff by pressing the Function key.



Turning the Display Off

The display can be turned off any time. To turn off the display:

- 1 Press the Screen key several times until the following screen appears.



- 2 Select the timing for turning off the display with the ◀ or ▶ key. The selected item is highlighted in blue. A 10 second countdown starts. You can select a different time within the 10 second countdown.

1 min (factory default): Turns the display off 1 minute later.

5, 10, 15, 30, 60, 120 or 240 min: Turns the display off when the selected time elapses. To set the time, refer to “SELECTABLE SCREEN TIME OUT PERIOD (min)” in the “Changing System Setup Settings” section.

CONTINUE: Keep the display turned on.

NOTE: If longer than “1 min” is selected, it reduces the battery lifetime.

- 3 Wait 10 seconds until the countdown ends. When the countdown ends, the numeric and waveform screen appears.

Or, press the Screen key before the countdown ends. The check electrodes screen appears.

When the selected time elapses without any key operation, the display turns off automatically. If a key is pressed, the countdown resets.

- When displaying the check electrodes screen after selecting the time-out period, the screen automatically changes to the numeric and waveform screen 2 minutes later if there is no key operation, then the display turns off when the remaining time elapses.
- When displaying another screen after selecting the time-out period, the screen changes to the numeric and waveform screen 1 minute later if there is no key operation, then the display turns off when the remaining time elapses.
- When “1 min” is selected, the display turns off without changing to the numeric and waveform screen.

When the display is turned off automatically or the power is turned off, the setting returns to “1 min” (factory default).

Turning the Display On after It was Turned Off

Press the Screen key. One of the following screens appears.


- The previous screen: The Screen key is pressed within 5 minutes after the display turned off.
- Numeric and waveform screen: The Screen key is pressed more than 5 minutes after the display turned off.
- “Key locked” screen: The Screen key is pressed after the display turned off and the “Key locked” screen was the last screen before the display turned off.

ECG and Respiration Monitoring

When the electrodes are attached and the ECG leads are connected, the heart rate, ECG waveform, respiration rate and respiration waveform appear on the receiving monitor. Refer to the operator's manual of the receiving monitor for details.

When 6 leads are used on this transmitter, up to 8 leads (I, II, III, aVR, aVL, aVF, Va and Vb) of ECG waveforms can be displayed on the receiving monitor. The heart rate is also measured.

When 3 leads are used, one channel ECG waveform of lead II can be displayed on the receiving monitor. Refer to the operator's manual of the monitor for details.



If the receiving monitor does not detect the pacemaker pulse, change the electrode positions. Attach the electrode lead clip to change the monitoring electrode.

⚠ WARNING

The bioelectric impedance measurement sensor of a minute ventilation rate-adaptive implantable pacemaker may be affected by transmitter which is connected to the same patient. If this occurs, the pacemaker may pace at its maximum rate and the transmitter may give incorrect data to the monitor. If this occurs, disconnect the electrode leads from the patient or change the setting on the pacemaker by referring to the pacemaker's manual. For more details, contact your pacemaker representative or Nihon Kohden representative.

⚠ WARNING

Turn the pacing pulse detection to ON on the receiving monitor when monitoring a pacemaker patient. Otherwise the pacemaker pulse is not rejected. However, even when the pacing pulse detection is set to ON, the pacemaker pulse might not be rejected. When the pacemaker pulse is not rejected, the pacemaker pulse is detected as QRS and false heart rate may be indicated or critical arrhythmia such as asystole may be overlooked. Keep pacemaker patients under close observation.

For the pacemaker pulse rejection capability of the ZM-540PA and ZM-541PA transmitters, refer to the "Specifications - ECG" section.

⚠ WARNING

Even when the pacing pulse detection is set to ON on the receiving monitor, the pacemaker pulse can be overlooked or detected as QRS. You cannot confirm the pacemaker operation only from the detected pacemaker pulse.

⚠ WARNING

When the transmitter is used with an electrosurgical unit (ESU), firmly attach the entire area of the ESU return plate. Otherwise, the current from the ESU flows into the electrodes of the transmitter, causing electrical burn where the electrodes are attached. For details, refer to the ESU manual.

⚠ CAUTION

Turn off the power of mobile phones, small wireless devices and other devices which produce strong electromagnetic interference around a patient (except for devices allowed by the hospital administrator). Radio waves from devices such as mobile phones or small wireless devices may be mistaken as pulse waves or respiration waves and the displayed data may be incorrect.

⚠ CAUTION

When the "NOISE" message is displayed on the receiving monitor, ECG is not monitored and the ECG alarm does not function. Remove the cause by checking the electrodes, electrode leads, patient's body movement, EMG and peripheral instruments grounding. Also make sure that an electric blanket is not used.

- NOTE
- Noise generated from an electrosurgical unit may interfere on an ECG waveform, but will not damage it.
 - If an electric blanket is used and incorrect heart rate is displayed on the monitor, turn off the pacing spike detection on the monitor.
 - Turn the pacing spike detection to ON on the receiving monitor when monitoring a pacemaker patient. Pacing pulse is detected by the transmitter and transmitted to the monitor. If the pacing spike detection is turned OFF, QRS and pacemaker spike might not be distinguished and pacemaker failure might not be recognized.
 - If defibrillation may be necessary, set the filter to Monitor or Maximum on the receiving monitor. If defibrillation is performed while the setting is OFF, the waveform recovery may be slow.
 - ECG cannot be monitored on a neonate using this transmitter.

Turning ECG Measurement On/Off

ECG measurement can be turned on or off on the PARAMETER SETUP screen. When electrodes are attached to the patient and ECG leads are connected, ECG monitoring starts even when ECG is turned off.

When “PROTOCOL” on the SYSTEM screen is set to 57:

ECG measurement on the receiving monitor is automatically set to OFF. Also, ECG measurement on the transmitter cannot be turned on or off from the receiving monitor.

When “PROTOCOL” on the SYSTEM screen is set to 42:

If ECG measurement is turned off on the transmitter, ECG measurement on the receiving monitor must also be turned off.

Turning Respiration Measurement On/Off

Respiration measurement can be turned on or off on the PARAMETER SETUP screen. If respiration measurement is turned off, respiration measurement on the receiving monitor is also turned off.

Electrode Detachment

In the following conditions, the “CHECK ELECTRODES” message is displayed on the transmitter and receiving monitor.

- Electrode is detached from skin.
- Electrode lead is disconnected from the electrode.
- Polarization voltage between the electrode and skin is excessively high.

In these cases, check the cause and if necessary, replace electrodes with new ones.


CAUTION

When the “CHECK ELECTRODES” message is displayed on the receiving monitor, ECG is not monitored properly and the ECG alarm does not function. Check the electrode, electrode leads, and if necessary, replace with new ones.

SpO₂ Monitoring

When monitoring starts, SpO₂ and the pulse waveform are sent to the monitor and SpO₂ and the pulse level bar graph are displayed on the transmitter screen. When ECG is not measured, the pulse waveform and pulse rate are also displayed.

Check that pulse wave is displayed and that the “CHECK PROBE” message is not displayed. If the probe is detached, SpO₂, pulse rate and pulse wave are not displayed on the display.



NOTE: During SpO₂ monitoring, the transmitter constantly checks the light intensity of the probe and adjusts the light intensity to maintain optimum measurement condition. During adjustment, the pulse wave becomes flat for about 0.5 seconds.



Measuring Principle

SpO₂ is measured by attaching a probe to a place where light can easily penetrate, such as a finger. The probe has two light emitter diodes and one photo detector.

Light of two wavelengths, 660 nm and 940 nm, are emitted from the emitter diodes. Oxygenated and deoxygenated blood absorb different wavelengths of light as shown in fig. 1 and 2. The light absorption of the skin, tissue, bone and venous blood is constant and can be eliminated. The pulsatile change represents arterial blood only, as shown in fig. 3, and this is measured to determine SpO₂. The amount of unabsorbed light at each wavelength is detected by the photo detector. SpO₂ is calculated from the ratio of absorbed light of the two wavelengths (fig. 2).

SpO₂ can be measured correctly when there is an arterial pulse which causes change in the transmitted light of both wavelengths. When the patient has asystole and no arterial pulse, there is a chance that some external factor such as movement of the SpO₂ probe may cause change in the transmitted light and produce a pulse-like signal which causes an incorrect SpO₂ value to be displayed. This incorrect SpO₂ value depends on the optical properties of the change in transmitted light.




Fig. 1




Fig. 2




Fig. 3

⚠ WARNING

SpO₂ measurement may be incorrect in the following cases.

- When the patient's carboxyhemoglobin or methemoglobin increases abnormally.
- When dye is injected in the blood.
- During CPR.
- When measuring at a site with venous pulse.
- When there is body movement.
- When the pulse wave is small (insufficient peripheral circulation).

⚠ WARNING

Check the circulation condition by observing the skin color at the measurement site and pulse waveform. Change the measurement site every 8 hours for disposable probes and every 4 hours for reusable probes (every 8 hours for TL-631T series probe). The skin temperature may increase at the attached site by 2 or 3°C (4 or 5°F) and cause a burn or pressure necrosis. When using the probe on the following patients, take extreme care and change the measurement site more frequently according to symptoms and degree.

- Patient with a fever
- Patient with insufficient peripheral circulation
- Neonate or low birth weight infant with delicate skin

⚠ WARNING

When monitoring SpO₂ of a patient who is receiving photodynamic therapy, the light from the finger probe sensor may cause a burn. Photodynamic therapy uses a photosensitizing agent that has a side effect of photosensitivity.

The SpO₂ probes manufactured by Nihon Kohden have two wavelengths with peaks in the range of 650 and 950 nm. The maximum light intensity is less than 5.5 mW.

⚠ WARNING

When not monitoring SpO₂, disconnect the SpO₂ cable from the transmitter. Otherwise, noise from the probe sensor may interfere and incorrect data is displayed on the screen.

⚠ CAUTION

Turn off the power of mobile phones, small wireless devices and other devices which produce strong electromagnetic interference around a patient (except for devices allowed by the hospital administrator). Radio waves from devices such as mobile phones or small wireless devices may be mistaken as pulse waves or respiration waves and the displayed data may be incorrect.

⚠ CAUTION

Handle the probe cable according to the following cautions. Failure to follow these cautions may cause cable discontinuity or short circuit of the probe cable which may cause incorrect measurement data or inability to perform measurement. Also in rare cases, the probe temperature may increase and cause skin burn on the patient. If the probe cable is damaged, replace the probe with a new one.

- Do not pull or bend the probe cable.
- Do not let caster feet run over the probe cable.

⚠ CAUTION

Normal external light does not affect monitoring but strong light such as a surgical light or sunlight may affect monitoring. If affected, cover the measuring site with a blanket.

⚠ CAUTION

When the probe is attached on an appropriate site with sufficient thickness and the error message confirming the probe attachment repeatedly appears, the probe may be deteriorated. Replace it with a new one.

⚠ CAUTION

When a message indicates a faulty probe, stop monitoring and replace the probe with a new one.

⚠ CAUTION

While a patient is on medication which causes vasodilation, the pulse waveform may change and in rare cases the SpO₂ value might not be displayed.

⚠ CAUTION

When monitoring SpO₂ only, detection of arrhythmia and asystole is not available and arrhythmia alarms such as ASYSTOLE, VF or VT are not available. If the patient requires ECG monitoring, monitor the ECG.

⚠ CAUTION

When monitoring SpO₂ only (without ECG monitoring), turn on both the upper and lower limit alarms for PR and SpO₂ on the receiving monitor. If the patient's pulse is not detected during asystole or other condition, a "CANNOT DETECT PULSE" or "CHECK PROBE" alarm occurs instead of an SpO₂ limit alarm. Furthermore, if the patient has no pulse, noise from probe movement could be misjudged as a pulse and cause an incorrect PR or SpO₂ value to be displayed.

- NOTE**
- In order to maintain sufficient blood circulation, keep the measurement site warm by covering it with a blanket or something similar. Warming the site is effective, especially for a patient with a small pulse amplitude.
 - When monitoring a patient who has an IABP and SpO₂ cannot be measured, monitor the patient on a wired monitor. If the monitor has a sensitivity mode, set the mode to "MAX".
 - Unlike ECG monitoring, detection of arrhythmia and asystole is not available when monitoring SpO₂ only. If SpO₂ is monitored without ECG, PR and SpO₂ alarms do not occur during asystole because PR and SpO₂ values are not measured when there is no pulse. Furthermore, if the patient has no pulse, noise from probe movement could be misjudged as a pulse and cause an incorrect PR or SpO₂ value to be displayed.


Monitoring SpO₂ during NIBP Measurement

CAUTION

NIBP and SpO₂ can be measured on the same limb, but the SpO₂ monitoring might not be accurate during NIBP measurement. Be careful when reading the SpO₂ values.*

- * When the SpO₂ probe is attached to the same limb as the NIBP cuff, the blood flow decreases during NIBP measurement and pulse wave cannot be detected and SpO₂ cannot be monitored properly. When "INHIBIT SpO₂ DURING NIBP" on the PARAMETER SETUP screen is set to ON (factory default setting), SpO₂ monitoring is paused during NIBP measurement to avoid SpO₂ alarm occurrence. However, when monitoring SpO₂ on the same limb as NIBP, be careful when reading SpO₂ values.

When monitoring SpO₂ is important, attach the probe to the limb to which the NIBP cuff or catheter is not attached.

When SpO₂ monitoring is paused during NIBP measurement, the SpO₂ value just before the start of NIBP measurement and an  mark are displayed on the transmitter for 30 seconds.

When NIBP measurement is not completed after 30 seconds, "-- --" is displayed for the SpO₂ value. The same data also appears on the monitor screen.

- NOTE**
- When continuous SpO₂ monitoring is necessary, attach the probe to the limb to which the NIBP cuff is not attached and set "INHIBIT SpO₂ DURING NIBP" on the PARAMETER SETUP screen to OFF.
 - When the probe is attached to the same limb as the NIBP cuff, set the sync source to a parameter other than SpO₂ on the receiving monitor.
 - When monitoring SpO₂ during STAT NIBP measurement, attach the probe to the limb to which the NIBP cuff is not attached.

NIBP Monitoring

Selecting the Initial Cuff Inflation Pressure

The initial cuff inflation pressure can be changed on the PARAMETER SETUP screen. The default setting is 180 mmHg. To change the setting, refer to the “Changing PARAMETER SETUP Settings” section.


Selecting the Measurement Mode and Interval

Measurement Modes

There are three measurement modes: manual, auto and STAT. The selected mode or interval is displayed on the screen.

The measurement mode and interval can be changed by pressing the NIBP Interval key.

When the key is pressed, the NIBP mode setting screen appears. The measurement modes selected at “SELECTABLE INTERVALS” on the PARAMETER SETUP screen are displayed (key color: white). Select the measurement mode with the ◀ and ▶ keys or NIBP Interval key and press the Function key.



To select the modes to be displayed on the NIBP mode setting screen, refer to the “Changing PARAMETER SETUP Settings” section.

Manual Measurement

In Manual mode, a single NIBP measurement is performed when the NIBP Start/Stop key is pressed.

STAT (Continuous) Measurement

In STAT mode, measurement is continuously repeated for 15 minutes after the NIBP Start/Stop key is pressed.

When the STAT measurement for 15 minutes is completed, the measurement mode automatically changes to the Manual mode or Auto mode of selected interval depending on the "NIBP MODE AFTER STAT" setting on the PARAMETER SETUP screen. The default setting is Manual mode. Refer to the "Changing Parameter Setup Settings" section.

The STAT measurement completes within 15 minutes. When more than 12 minutes elapse from the start of measurement, there will be no more measurement performed and the measurement mode changes to the mode selected for "NIBP MODE AFTER STAT" on the PARAMETER SETUP screen.

Auto Measurement

In Auto mode, measurement is performed automatically at the preset time intervals.

In Auto mode, a single measurement can be performed by pressing the NIBP Start/Stop key between auto measurements.

Measuring NIBP

WARNING

Be careful when measuring NIBP on a patient with known bleeding disorders or coagulation. After NIBP measurement, there may be dot hemorrhage, or circulatory disorder by thrombus where the cuff is attached.

WARNING

NIBP measurement may be incorrect in the following situations.

- When using an ESU
- Body movement
- Small pulse wave
- Too many arrhythmias
- Shaking from an external source
- Rapid blood pressure change
- During CPR
- Slow pulse
- Low blood pressure
- Small pulse pressure
- Cuff is too tight or too loose
- Cuff does not fit the arm
- Cuff is wrapped over thick clothing
- Cuff is deteriorated
- Arterial sclerosis
- Poor perfusion
- Diabetes
- Age
- Pregnancy
- Pre-eclampsia
- Renal disease
- Shivering
- Trembling

⚠ WARNING

When performing NIBP measurements in STAT mode or 5 minute intervals, periodically remove the cuff from the patient for ventilation. The skin temperature may increase at the cuff attachment site by 2 or 3°C (4 or 5°F). When measuring a patient with a fever or peripheral circulation insufficiency, it may cause a burn.

⚠ WARNING

When performing long term measurement at intervals less than 5 minutes, perform measurements while observing the state of the patient, blood vessels and limb to ensure adequate circulation. Congestion may occur at the measurement site. When performing periodic measurement for a long time, periodically check the circulation condition.

⚠ CAUTION


When performing NIBP measurement repeatedly, have a rest between measurements to recover adequate circulation.

- NOTE
- When measuring patients who are conscious, help the patient to relax. Measurement may not be accurate if the patient's arm is tense or if the patient talks.
 - The data for measurement on a leg tends to be higher than measurement on the arm. When making diagnosis based on the NIBP values, measure NIBP on an upper arm.
 - Do not apply pressure to the cuff or air hose. NIBP may not be measured correctly because of noise or NIBP measurement may stop due to the NIBP safety circuit.
 - If there is an abnormal noise generated during measurement, stop using the transmitter and contact your Nihon Kohden representative.
 - Do not measure NIBP of a patient on whom an IABP is being used. Measurement may be incorrect due to the mixing of the patient's own pulse and IABP pulse.
 - NIBP cannot be measured on a neonate using this transmitter.


For Accurate NIBP Measurement

- Make sure the patient is comfortably seated. Do not perform measurement immediately after the patient is seated. Keep the patient in the resting state for five minutes before measurement.
- Make sure the patient has legs uncrossed, feet down on the floor, and back and back arm supported during measurement.

1 Press the NIBP Interval key to display the NIBP mode setting screen.




2 Select the measurement mode by pressing the NIBP Interval key or and keys.



3 Press the Function key.

4 Press the NIBP Start/Stop key to perform measurement.



The cuff is inflated and the inflation pressure is displayed on the screen.



In manual mode: Measurement is performed once.

In STAT mode: Measurement is performed repeatedly for 15 minutes.


In auto mode: The first measurement is performed when the NIBP Start/Stop key is pressed. The second measurement is performed when the current time in the transmitter reaches the selected time interval.

To stop measurement during measurement, press the NIBP Start/Stop key again.

In STAT mode, after completing the STAT measurement, the measurement mode changes to the mode set for “NIBP MODE AFTER STAT” on the PARAMETER SETUP screen.

In auto mode, to stop measurement in auto mode, change the mode to manual. To cancel one measurement, press the NIBP Start/Stop key during measurement.

After the measurement is complete, the measured data is displayed on the screen and is transmitted to the monitor.



When ECG and SpO₂ are not monitored (ECG measurement is turned off and SpO₂ probe is not connected to the transmitter), the pulse rate at the end of NIBP measurement is displayed.


A buzzer can be set to sound at the start and end of NIBP measurement. Refer to the “Changing PARAMETER SETUP Settings” section.

Auto Mode Measurement

When auto mode measurement is selected, “STANDBY” message is displayed on the screen until the NIBP Start key is pressed for the first time.



A time bar appears to indicate the interval between auto mode measurements.



During auto mode measurement, the measurement mode can be changed. During the interval, press the NIBP Interval key to change the mode. When “MANUAL” is displayed for more than one second, the measurement in auto mode is stopped.

Data Display After NIBP Measurement






When the time set at “OLD NIBP DATA” on the PARAMETER SETUP screen elapses after the last measurement, the NIBP data is dimmed or hidden. Whether to dim or hide the old data can also be selected at “OLD NIBP DATA”. Refer to the “Changing PARAMETER SETUP Settings” section.

Monitoring SpO₂ during NIBP Measurement

When the SpO₂ probe is attached to the same limb as the NIBP cuff, the blood flow decreases during NIBP measurement and pulse wave cannot be detected and SpO₂ cannot be monitored properly. When “INHIBIT SpO₂ DURING NIBP” on the PARAMETER SETUP screen is set to ON (factory default setting), SpO₂ monitoring is paused during NIBP measurement to avoid SpO₂ alarm occurrence. However, when monitoring SpO₂ on the same limb as the NIBP, be careful when reading SpO₂ values.

Indication and Message List

Indication


| Indication | Cause | Countermeasure |
|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------|
|  | Fully charged batteries | — |
|  | Batteries are 1/3 full. | |
|  | Batteries are weak. | Replace the batteries. |
|  | Batteries are very weak. Cannot measure NIBP. | |
|  | Alarms on the receiving monitor were suspended by pressing the Function key on the transmitter. | Alarms resume when the suspend interval elapses. To cancel alarm suspension, press the Function key again. |

Messages

When “PROTOCOL” on the SYSTEM screen is set to 57, all messages are transmitted to the receiving monitor. When “PROTOCOL” is set to 42, the messages marked with * are not transmitted to the receiving monitor.

| Message | Cause | Countermeasure |
|------------------|----------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------|
| AIR LEAK | The cuff and extension hose are not properly connected | Connect them properly. |
| | The cuff hose (or extension hose) is not properly connected to the NIBP socket | |
| | The cuff or extension hose is damaged | Replace with a new one. |
| ALARMS SUSPENDED | Alarms on the receiving monitor is suspended by pressing the Function key on the transmitter | Alarms resume when the 2 minute suspend interval elapses. To cancel alarm suspension, press the Function key again. |
| BATTERY WEAK | Dead batteries | Replace batteries. |

| Message | Cause | Countermeasure |
|------------------------------------------|-------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------|
| CANNOT DETECT PULSE* (displayed in blue) | Poor blood circulation for measuring the SpO ₂ value | Check the patient condition, probe attachment or change the attachment site. |
| | The probe is attached too tightly and is obstructing the blood circulation | Reattach the probe. |
| | The probe is not attached to the patient properly | Attach the probe to the patient properly. |
| | “LIGHT INTERFERENCE”, “CHECK PROBE SITE” or “DETECTING PULSE” message is displayed for more than 30 seconds | Refer to the cause and countermeasure for each message in this Messages table and remove the cause. |
| CANNOT DETECT PULSE (displayed in pink) | The patient’s pulse wave is too small to measure NIBP | Measure by palpation or auscultation. |
| | The cuff is not wrapped on the patient properly | Wrap the cuff on the patient properly. |
| CHECK ELECTRODES | Electrode lead is disconnected from the electrode | Firmly connect the electrode lead to the electrode. |
| | Electrode lead is disconnected from the transmitter | Firmly connect the electrode lead to the transmitter. |
| | Electrode lead discontinuity | Replace the electrode lead with a new one. |
| | Electrode is not firmly attached to the skin | Replace the electrode with a new one. |
| | Polarization voltage is abnormally high | |
| CHECK PROBE | The probe is not attached to the patient properly | Attach the probe to the patient properly. |
| | The probe is not attached at the appropriate site | Attach the probe to an appropriate site indicated in the probe manual. |
| | The probe is disconnected from the transmitter | Connect the probe cable to the transmitter. |
| | The probe is past its expiration date | Replace the probe with a new one. |
| CHECK PROBE SITE* | The probe is not attached at the appropriate site | Attach the probe to an appropriate site indicated in the probe manual. |
| | The probe is deteriorated | Replace the probe with a new one. |
| | The probe is past its expiration date | |
| CUFF OCCLUSION | Transmitter malfunction | Immediately remove the cuff from the patient and contact your Nihon Kohden representative. |

| Message | Cause | Countermeasure |
|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| DETECTING PULSE | Searching for the correct pulse wave for SpO ₂ monitoring | Wait until the pulse wave is detected. |
| | The SpO ₂ value cannot be obtained because the waveform is unstable | Attach the probe to the patient properly. |
| | The probe is not attached to the patient properly | |
| HIGH CUFF PRESS | Enormous pressure was applied by the pressure of the cuff | Remove the cause. |
| INFLATION PRESS LOW | Insufficient cuff inflation pressure | Wait for the remeasurement to be performed with increased cuff inflation pressure. |
| LIGHT INTERFERENCE | The SpO ₂ measurement site is under fluorescent light, surgical light, sunlight, or other strong light | Cover the measurement site with a blanket or cloth. |
|  | Considerable body movement | When the message is displayed frequently, check the patient condition and, if necessary, change the attachment site. |
| | The probe is not attached to the patient properly | |
| | SpO ₂ monitoring is paused for NIBP measurement | Wait for NIBP measurement to finish. |
| MEAS TIME OUT | The NIBP measuring time exceeded the specified time due to arrhythmia, body movement, vibration or, cuff or air hose being squeezed | Remove the cause if the cause is body movement, vibration or squeezing of cuff or hose. |
| NIBP MODULE ERROR | Module malfunction | Contact your Nihon Kohden representative. |
| NO NIBP CHANGE BATTERIES | NIBP cannot be measured due to low battery | Replace batteries with new ones. |
| PATIENT CONFIRMED* | Function key is pressed and the "PATIENT CONFIRMED" signal is transmitted to the receiving monitor (when "PATIENT" is assigned as the function for the Function key on the SYSTEM SETUP screen) | — |
| PROBE FAILURE* | The probe is past its expiration date | Replace the probe with a new one. |
| | Probe is damaged or short-circuited | |
| REMEASURING | NIBP is being remeasured due to arrhythmia, body movement, vibration or, cuff or air hose being squeezed | If the message still appears after remeasurement, remove the cause if the cause is body movement, vibration or squeezing of cuff or hose. |
| SAFETY CIRCUIT ERROR | The NIBP safety circuit error | Immediately remove the cuff from the patient and contact your Nihon Kohden representative. |

| Message | Cause | Countermeasure |
|-------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------|
| SAFETY CIRCUIT RUNNING (When this message is displayed, measurement cannot be performed for 40 seconds.) | NIBP measurement stopped by the safety circuit | Check that the hose is not bent or squeezed. |
| | | Wait 40 seconds, then perform remeasurement. If the message still appears, contact your Nihon Kohden representative. |
| SpO ₂ MODULE ERROR* | Transmitter failure | Contact your Nihon Kohden representative. |
| SYS OUT OF RANGE | The systolic blood pressure exceeded the 40 to 280 mmHg range and the blood pressure cannot be measured. | Measure by palpation or auscultation. |
| WEAK PULSE* (displayed in blue) | Poor peripheral circulation | Check the patient condition and change the probe attachment site. |
| | The probe is attached too tightly and is obstructing the blood circulation | Check the probe attachment condition and if necessary, reattach the probe. |
| WEAK PULSE (displayed in pink) | The patient's pulse wave is too small to measure NIBP | Measure NIBP by palpation or auscultation. |
| | The cuff is wrapped too loosely | Wrap the cuff properly. |
| | The cuff size is not appropriate | Use the appropriate cuff. |
| ZEROING | NIBP zero balance is being adjusted | Do not touch the cuff during zeroing. Wait for the message to disappear. |

Message Display Priority

When more than one message condition occurs on the transmitter, only the message with the highest priority is displayed.

| Priority | Message |
|----------|-----------------------------------------|
| High | PATIENT CONFIRMED |
| ↑ | SAFETY CIRCUIT RUNNING |
| | CUFF OCCLUSION |
| | PROBE FAILURE |
| | CHECK ELECTRODES |
| | NIBP MODULE ERROR |
| | SYS OUT OF RANGE |
| | HIGH CUFF PRESS |
| | AIR LEAK |
| | MEAS TIME OUT |
| | CANNOT DETECT PULSE (NIBP) |
| | SpO ₂ MODULE ERROR |
| | CHECK PROBE |
| | CHECK PROBE SITE |
| | CANNOT DETECT PULSE (SpO ₂) |
| | LIGHT INTERFERENCE |
| | REMEASURING |
| | INFLATION PRESS LOW |
| | WEAK PULSE (NIBP) |
| | ZEROING |
| | NO NIBP CHANGE BATTERIES |
| | DETECTING PULSE |
| | WEAK PULSE (SpO ₂) |
| | ALARMS SUSPENDED |
| | |
| Low | BATTERY WEAK |

Troubleshooting

If a problem occurs, use the following tables to find and fix it. If the problem still remains after troubleshooting according to these tables, contact your Nihon Kohden representative.

Transmitter

| Problem | Cause | Countermeasure |
|---------------------------------------------------------------------------------------|----------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Nothing is displayed on the LCD after turning the power on. | Batteries are not installed correctly. The battery polarity is wrong. | Install the batteries correctly. |
| | Batteries are completely discharged. | Replace the batteries with new ones. |
| LCD is difficult to see (too dark or too light). | LCD brightness is not appropriate. | Change the LCD brightness on the SYSTEM SETUP screen. Refer to the "Changing SYSTEM SETUP Settings" section. |
| Nothing is displayed on the receiving monitor after turning the transmitter power on. | The channel of the transmitter and monitor does not match. | Set the correct channel on the monitor. |
| | The software version of the multiple patient receiver or central monitor is old. | Upgrade the multiple patient receiver or central monitor software to receive signal from the transmitter. The software version must be 01-09 or later. |
| | Protocol on the transmitter and monitor does not match. | Set the same protocol on the transmitter and monitor. |
| Signal receiving condition is poor. | Another transmitter with the same channel is used nearby. | Turn the transmitter power off. If the monitor still receives a signal, there is a high probability that another transmitter of the same channel is used nearby. Follow the instruction of your channel administrator and use another transmitter with a different channel. |
| | Signals of another patient are mixing. | Follow the instructions of your channel administrator and use another transmitter of a different channel. |
| | The transmitter is damaged. | Contact your Nihon Kohden representative. |

ECG/Respiration

| Problem | Cause | Countermeasure |
|----------------------------------------------------------------------------------------------------|-----------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------|
| The heart rate is unstable. | Pacing detection setting on the monitor is not correct. | Turn off the pacing detection setting on the monitor. When monitoring a pacemaker patient, turn on pacing detection. |
| The "CHECK ELECTRODES" message appears on the receiving monitor. | Electrode lead is disconnected from the electrode. | Firmly connect the electrode lead to the electrode. |
| | Electrode lead discontinuity | Replace the electrode lead with a new one. |
| | Electrode is not firmly attached to the skin. | Replace the electrode with a new one. |
| | Polarization voltage is abnormally high. | Use Nihon Kohden specified electrodes. |
| ECG baseline is thick (AC hum). | The gel on the electrode is dried out. | Replace the electrode with a new one. |
| | The gel on the electrode is coming off. | |
| | An electric blanket is used. | Cover the blanket with a shield cover. |
| | The hum filter is set to OFF on the monitor. | Set the filter to ON. |
| The heart rate of a patient who is using an electric blanket is unstable on the receiving monitor. | Pacing pulse detection is turned ON on the receiving monitor. | Turn OFF the pacing pulse detection on the receiving monitor. |
| No heart rate or ECG is displayed. | "ECG MEASUREMENT" on the PARAMETER SETUP screen is set to OFF. | If ECG monitoring is necessary, set "ECG MEASUREMENT" to ON. |
| Respiration waveform measurement is unstable. | The gel on the electrode is dried out. | Replace the electrode with a new one. |
| | The gel on the electrode is coming off. | |
| No respiration rate is displayed. | "RESP MEASUREMENT" on the PARAMETER SETUP screen is set to OFF. | If respiration monitoring is necessary, set "RESP MEASUREMENT" to ON. |

SpO₂

| Problem | Cause | Countermeasure |
|-----------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|
| SpO ₂ data is unstable and not reliable. | The probe size is not appropriate for the patient. | Use the appropriate probe for the patient. |
| | Probe attachment condition is poor. The probe is partly detached from the skin. External light is entering the probe. | Firmly attach the probe according to the procedure in the probe operator's manual. |
| | Measurement site is dirty. Patient is wearing nail polish. | Remove dirt and nail polish. |
| | The probe is attached to the same limb that is used for NIBP measurement. | Attach the probe to the opposite limb. Avoid a site where blood circulation condition changes greatly. |

NIBP

| Problem | Cause | Countermeasure |
|--------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------|
| Cuff inflation pressure is less than 10 mmHg. | The cuff hose is not connected to the NIBP socket properly. | Connect the cuff hose to the socket properly. |
| | The cuff is not wrapped around the arm or is wrapped too loosely. | Wrap the cuff around the upper arm. |
| The cuff does not inflate when the NIBP Start/Stop key is pressed. | The cuff hose is not connected to the NIBP socket. | Connect the cuff hose to the socket firmly. |
| | The cuff hose or extension hose may be folded or squeezed when the cuff pressure display on the screen increases quickly but the actual cuff does not inflate. | Check the cuff hose and air hose. |
| Abnormal measurement results are displayed. | The cuff size is not correct. | Select the cuff which fits the patient's limb circumference. |
| | The cuff is not wrapped around the arm correctly. | Wrap the cuff around the upper arm, not too tightly or too loosely. |
| | NIBP data is not correct because of body movement. | Prevent the patient from moving during measurement. |
| | Vibration on the cuff. | Check that nothing is touching the cuff during measurement. Change the measuring site. |
| The cuff is suddenly deflated during inflation. | The NIBP Start/Stop key is pressed during inflation. | — |

| Problem | Cause | Countermeasure |
|------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------|
| Auto mode measurement does not start even when the time interval has passed. | The NIBP Interval key is pressed and the measurement mode is changed. | Check the measurement mode and interval. |
| The cuff suddenly inflates. | The measurement mode is set to auto mode. | Check the time interval. If necessary, stop measurement. |
| Cannot connect cuff to the air hose. | Unspecified cuff is used. | Use a cuff specified by Nihon Kohden. |
| Cannot measure NIBP. | Vibration on the cuff. | Check that nothing is touching the cuff during measurement. |
| | The cuff hose or extension hose is bent or squeezed. | Remove the cause. |
| | The cuff has worn out. | Use a new cuff. |
| Blood congestion occurs. | Measuring over a long period of time at short intervals. | Increase the measuring interval. |
| | | Do not measure NIBP over a long time. |
| Thrombus occurs. | Measuring on a patient with known bleeding disorders or coagulation. | Do not perform NIBP measurement on such a patient. |
| NIBP data on the screen is --- or dark. | The time set for "OLD NIBP DATA" on the PARAMETER SETUP screen elapsed from the last measurement. | When NIBP is measured again, the data is displayed in normal brightness. |
| Three loud pip sounds indicating NIBP measurement cannot be started. | The cuff is not deflated enough to start another measurement. | Wait 30 seconds and measure again. |

Maintenance

To use the transmitter in safe and optimum condition, perform maintenance check at least once every year.

The transmitter contains parts which gradually deteriorate with use. Original performance might not be delivered if any part of the transmitter is deteriorated. Perform regular maintenance checks to assure continued safe operation.

The following units are necessary for some checking items.

- AX-410G medical instrument checker
- AX-300T SpO₂ checker
- Electric or mercury manometer
- 700 mL dummy cuff
- Receiving monitor

 **CAUTION**

Never disassemble or repair the transmitter. If there is any problem with the transmitter, contact your Nihon Kohden representative.

 **CAUTION**

Before maintenance, cleaning or disinfection, turn the transmitter power off and remove the batteries. Failure to follow this instruction may result in electrical shock and transmitter malfunction.

- NOTE**
- The measurement accuracy of the above units must be managed to perform accurate maintenance check.
 - For details on the operation of the above units, refer to the manuals provided with these units.

A maintenance check sheet is provided at the end of this operator's manual. Make a copy of this check sheet before performing maintenance check.


1. External Check

- There is no damage or dirt on the outside of the transmitter.
- The battery case cover is not damaged and can be closed firmly.
- No keys are damaged.
- NIBP socket is not damaged.
- No electrode leads are damaged.
- There are no blood or chemicals on the transmitter.
- The springs in the battery compartment are not damaged or detached.
- Terminals in the battery case are not corroded.

2. Transmitter Channel

- The channel number label on the transmitter is not torn or removed.
- The channel of the transmitter matches the label.


The transmitter channel is displayed in the upper left corner of the screen. The channel number also appears on the startup screen.



3. Transmitting/Receiving Signal

Use the AX-410G medical instrument checker and receiving monitor.

- 1 Connect the medical instrument checker to the transmitter only with the electrode leads.




- 2 Place the transmitter 2 to 3 m from the receiving monitor.
- 3 Set the channel on the receiving monitor to the channel of the transmitter.
- 4 Turn on the transmitter and medical instrument checker.
- 5 Check that the ECG of the transmitter appears on the receiving monitor.
- 6 Turn off the transmitter.
- 7 Check that the ECG disappears from the receiving monitor.


4. Display

Check that there are no dots missing on the screen.


- 1 Turn off the transmitter.
- 2 While pressing the Function key, turn on the transmitter. The MENU screen appears.
- 3 Press the ► key to move the cursor to “MANUAL CHECK” and press the Function key.



- 4 Press the ◀ or ▶ key to move the cursor to “LCD TEST” and press the Function key.



- 5** Each time the ► key is pressed, the screen changes as below. Check that no dots are missing.



- 6** Press the Screen key to return to the MANUAL CHECK screen.


- 7** Press the Screen key again to return to the MENU screen.

5. Key Operation


- 1** Turn off the transmitter.
- 2** While pressing the Function key, turn on the transmitter. The MENU screen appears.
- 3** Press the ► key to move the cursor to “MANUAL CHECK” and press the Function key.



- Press the ◀ or ▶ key to move the cursor to “KEY CHECK” and press the Function key.



- Press each key one at a time and check that the pressed key is highlighted on the screen.




Example:
When the Function key is pressed, the key name is highlighted.

- After checking, press and hold the Screen key to return to the MANUAL CHECK screen.

6. ECG Check

- Connect the medical instrument checker to the transmitter only with the electrode leads.




- Place the transmitter 1 m from the receiving monitor.

- 3 Turn on the transmitter and medical instrument checker.
- 4 Check that the ECG of the transmitter appears on the receiving monitor.

7. Respiration Check

- 1 Connect the medical instrument checker to the transmitter only with the electrode leads.




- 2 Place the transmitter 1 m from the receiving monitor.
- 3 Turn on the transmitter and medical instrument checker.
- 4 Check that the respiration waveform of the transmitter appears on the receiving monitor.

8. SpO₂ Check

SpO₂ can be checked with the AX-300T SpO₂ checker or AX-410G medical instrument checker.

With the AX-300T SpO₂ Checker

- 1 Connect the SpO₂ checker to the transmitter only with the SpO₂ connection cable.



- 2 Place the transmitter 1 m from the receiving monitor.
- 3 Turn on the transmitter and SpO₂ checker.
- 4 Check that the pulse bar graph appears on the transmitter screen.

5 Check that SpO₂ and pulse rate on the transmitter is within the following range.


| SpO ₂ on the SpO ₂ Checker | | Range |
|--------------------------------------------------|---------------|----------------------------------------|
| SpO ₂ ¹ | 97% | 95 to 99%SpO ₂ (±2 digits) |
| | 80% | 78 to 82%SpO ₂ (±2 digits) |
| | 70% | 66 to 74%SpO ₂ (±4 digits) |
| Pulse rate | 60 beats/min | 57 to 63 beats/min (±3% ±1 beat/min) |
| | 120 beats/min | 115 to 125 beats/min (±3% ±1 beat/min) |

¹ The SpO₂ check by the SpO₂ checker is affected by the checker's tolerance to the SpO₂ measuring accuracy of the transmitter. (The measurement accuracy is described in the "Specifications" section.) For details, refer to the SpO₂ checker manual.

6 Check that the SpO₂ and pulse waveform of the transmitter appear on the receiving monitor.

With the AX-410G Medical Instrument Checker

1 Connect the medical instrument checker to the transmitter only with the SpO₂ connection cable.



2 Place the transmitter 1 m from the receiving monitor.

3 Turn on the transmitter and medical instrument checker.

4 Check that the pulse bar graph appears on the transmitter screen.

5 Check that SpO₂ and pulse rate on the transmitter is within the following range.

| Medical Instrument Checker | | Range |
|-------------------------------|---------------------|-----------------------------------------------------------|
| SpO ₂ ¹ | 97%SpO ₂ | 95 to 99%SpO ₂ (97%SpO ₂ ± 2 digit) |
| | 80%SpO ₂ | 78 to 82%SpO ₂ (80%SpO ₂ ± 2 digit) |
| | 70%SpO ₂ | 66 to 74%SpO ₂ (70%SpO ₂ ± 4 digit) |
| PR | 60 bpm | 57 to 63 bpm (120 bpm ± 3% ±1 bpm) |

¹ The SpO₂ check by the medical instrument checker is affected by the checker's tolerance to the SpO₂ measuring accuracy of the transmitter. (The measurement accuracy is described in the "Specifications" section.) For details, refer to the medical instrument checker manual.


- 6 Check that the SpO₂ and ECG of the transmitter appear on the receiving monitor.

9. NIBP Check


Check that the transmitter displays the correct cuff pressure and that there is no air leak. Also check the pressure sensor. The following procedure uses an electric manometer.

- 1 Connect the electric manometer and dummy cuff to the transmitter.


NOTE: Air leaks from the rubber ball during inflation. Use a clip on the air hose of the rubber ball to stop air leaking.




- 2 Turn on the electric manometer.
- 3 While pressing the Function key, turn on the transmitter. The MENU screen appears.
- 4 Press the ► key to move the cursor to “MANUAL CHECK” and press the Function key.



- 5 Press the ◀ or ▶ key to move the cursor to “NIBP HARDWARE TEST” and press the Function key.



- 6 Press the ◀ or ▶ key to move the cursor to “AIR LEAK TEST” and press the Function key.




- 7 Press the ◀ or ▶ key to move the cursor to “AIRLEAK(AUTO)” and press the Function key.




The transmitter inflates the cuff up to about 250 mmHg and measures air leakage from 60 seconds to 120 seconds after inflation.

- 8 Check the following.
 - The value for “AIR LEAK (AUTO)” is below 10 mmHg.
 - The difference between the pressure value displayed on the manometer and transmitter is within ± 6 mmHg.
- 9 Press the Screen key to return to the NIBP HARDWARE TEST screen.
- 10 Press the ◀ or ▶ key to move the cursor to “PRESSURE TEST” and press the Function key.




- 11 Remove the clip from the air hose of the rubber ball.



- 12 By squeezing the rubber ball, increase pressure up to 300 mmHg on the manometer.

13 Check the following.

- The difference between the pressure value displayed on the manometer and “PRESS. (SAFETY)” value on the transmitter is within ± 15 mmHg.
- The difference between the pressure value displayed on the manometer and “PRESS” value on the transmitter is within ± 6 mmHg.



10. NIBP Cuff for Attaching Transmitter to Patient Arm

The NIBP cuff is a consumable. Check the following and when necessary, replace it with a new one.

Appearance

- There are no dirty parts.
- There are no broken stitches on the cuff.
- The label on the cuff is readable.
- The velcro tape on the cuff is not removed and there are no broken stitches.

Inflation bag

- The inflation bag is not torn or damaged.
- There is no water inside the inflation bag.
- The connector on the inflation bag is not damaged.

Lifetime and Disposal

CAUTION

Dispose of Nihon Kohden products according to your local laws and your facility's guidelines for waste disposal. Otherwise, it may affect the environment. If there is a possibility that the product may have been contaminated with infection, dispose of it as medical waste according to your local laws and your facility's guidelines for medical waste. Otherwise, it may cause infection.

Disposing of Used Batteries

Battery Lifetime

Replace the batteries when the battery replacement indication appears on the transmitter. When using rechargeable batteries, recharge them.

Disposal

NOTE: Remove the batteries before disposing of the transmitter.

Before disposing of the batteries, check with your local solid waste officials for details in your area for proper disposal. It may be illegal to dispose of these batteries in the municipal waste stream.

Disposing of Electrodes, SpO₂ Probes and NIBP Cuffs

Refer to the manual for each item.

Disposing of Transmitter

Remove the batteries from the transmitter and dispose of the transmitter following your local laws for disposal.

Cleaning, Disinfection and Sterilization

Transmitter and Electrode Leads

⚠ CAUTION

If detergent or liquid spills into the transmitter, stop cleaning or disinfecting it and contact your Nihon Kohden representative. The transmitter needs to be checked for safety and function before use.

- NOTE**
- The transmitter cannot be sterilized. Sterilizing the transmitter may damage it.
 - When cleaning or disinfecting the transmitter, disconnect it from the patient and remove the batteries beforehand. Otherwise it may result in malfunction.
 - Do not let the NIBP socket get wet during cleaning and disinfection.
 - Do not use a hair dryer etc. to dry the transmitter when it is wet. Doing so may cause the transmitter to deform or become damaged. It may also cause loss of the waterproofing function.
 - Do not use volatile liquids such as thinner or benzine because these will cause the materials to deform or crack.
 - Avoid using flammable disinfectants such as ethanol in a closed place. Ventilate the room if you use flammable disinfectants.
 - Wipe the transmitter with a dry cloth before use and dry it completely after cleaning and disinfecting.

Before cleaning or disinfecting the transmitter, disconnect it from the patient and remove the batteries from the transmitter. Be careful not to let any liquid get inside the transmitter.

Cleaning

Wipe the transmitter and electrode leads with a soft cloth moistened with disinfecting alcohol or neutral detergent diluted with water. After cleaning, dry them completely.

Use cotton swab moistened with neutral detergent diluted with water to clean inside the battery compartment.

Disinfection

⚠ CAUTION

Do not immerse the electrode lead connector in liquid.

Wipe the outside surface of the transmitter and electrode lead with a non-abrasive cloth moistened with any of the disinfectants listed below. For details on the disinfectants, refer to the instruction provided with the disinfectants. Use the recommended concentration.

| <u>Disinfectant</u> | <u>Concentration (%)</u> |
|----------------------------------------|--------------------------|
| Glutaraldehyde solution | 2.0 |
| Alkyldiaminoethylglycine hydrochloride | 0.5 |
| Benzalkonium chloride | 0.2 |
| Benzethonium chloride solution | 0.2 |
| Chlorhexidine gluconate solution | 0.5 |

SpO₂ Probe

Refer to the probe manual.

NIBP Cuff

Refer to the cuff manual.

Periodic Inspection

If the periodic inspection is not performed, degradation or loss of function may go unnoticed and lead to misdiagnosis.

Service personnel should perform the periodic inspection at least once every year. Make sure that the transmitter operates properly and replace the consumables.

If you found abnormalities as a result of inspection and the transmitter is suspected to be faulty, attach an "Unusable" or "Repair request" label to the transmitter and contact your Nihon Kohden representative.

Repair Parts Availability Policy

Nihon Kohden Corporation (NKC) shall stock repair parts (parts necessary to maintain the performance of the instrument) for a period of 8 years from the date of delivery. In that period NKC or its authorized agents will repair the instrument. This period may be shorter than 8 years if a board or part necessary for the faulty section is not available.

Specifications

ZM-540PA

Measured Parameters

| | |
|---------------|-------------------------------------------------------------------|
| Waveforms: | ECG, impedance method respiration, pulse |
| Numeric data: | Heart rate, respiration rate, SpO ₂ , NIBP, pulse rate |

Transmitted Data

| | |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Waveforms: | ECG, respiration, pulse wave |
| Numeric data: | SpO ₂ , NIBP, pulse rate |
| Status information: | Battery replacement, battery level ¹ , alarm suspended, pause monitoring ¹ , patient confirmed ¹ , ECG lead, pacing detection, electrode detachment, electrode impedance ¹ , ECG off ¹ , respiration method (impedance) ¹ , SpO ₂ status, NIBP status, channel ID, time constant (3.2 s), type of transmitter, transmitter code number ¹ , transmitter serial number ¹ |

¹ Transmitted only when the protocol is "57".

Display

| | |
|---------------|--------------------------|
| Display size: | 2.2 inch TFT color LCD |
| Viewing area: | 44.16 (H) × 33.12 (V) mm |
| Resolution: | 320 (H) × 240 (V) dots |

Displayed Data

| | |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Numeric and waveform screen: | ECG (one waveform from lead I, II, III, Va or Vb), heart rate, pulse rate, respiration rate, SpO ₂ , NIBP (systolic, diastolic, MAP), message, battery level, QRS/pulse sync mark, pulse bar graph, NIBP measurement mode and status information, ECG lead |
| Waveform review screen: | ECG or pulse wave of past 10 minutes |
| Numeric review screen: | Heart rate or pulse rate, respiration rate and SpO ₂ at 1 minute interval for past 10 minutes |
| Check electrodes screen: | ECG for checking electrode attachment |

ECG

ECG Measurement

| | |
|---------------------------------------|------------------------------------|
| Channels: | 4 |
| Input dynamic range: | ±10 mV or more |
| Electrode offset potential tolerance: | ±500 mV or more |
| Input impedance: | 5 MΩ or more |
| Common mode rejection ratio: | 95 dB or more |
| | IEC 60601-2-27 50.102.10 compliant |

| | |
|----------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Pacing pulse detection: | amplitude ± 2 to 700 mV, duration 0.1 to 2 ms IEC 60601-2-27: 2005 compliant Based upon pacemaker pulse rejection capability |
| Defibrillation-proof: | ECG input protected against 400 Ws/DC 5 kV IEC 60601-2-27 17.101 compliant |
| ECG recovery time after defibrillation: | within 10 s |
| Electrode condition: | Displays CHECK ELECTRODES message |
| Tall T-wave rejection capability: | Complies with the heights of T-waves from 0 to 1.6 mV IEC 60601-2-27: 2005 50.102.17 compliant |
| Pacemaker pulse rejection capability, without overshoot: | Complies with the amplitudes of pacemaker pulses ± 2 to ± 700 mV and widths 0.1 to 2 ms (As defined by method B of IEC 60601-2-27: 2005 50.102.13, this corresponds to the pacemaker pulses amplitudes and widths of ± 2.8 mV/2 ms to amplitudes ± 45.6 mV/0.1 ms.) |
| Pacemaker pulse rejection capability, with overshoot: | Overshoot amplitudes and time constants of ± 0.12 mV/100 ms to ± 2 mV/4 ms (As defined by method B of IEC 60601-2-27: 2005 50.102.13, this corresponds to the pacemaker pulses amplitudes and widths of ± 2.8 mV/2 ms to amplitudes ± 45.6 mV/0.1 ms.) |

ECG Display and Heart Rate Count

| | |
|---------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| Frequency characteristic: | filter on: 1 to 18 Hz, filter off: 0.05 to 60 Hz |
| Heart rate detection method: | Average |
| QRS detection: | 70 to 120 ms: amplitude ≥ 0.5 mV, rate 30 to 200 beats/min 40 to 120 ms: amplitude ≥ 0.5 mV, rate 30 to 250 beats/min |
| Heart rate counting range: | 0, 15 to 300 beats/min |
| Heart rate counting accuracy ¹ : | ± 2 beats/min, (0, 15 to 300 beats/min) |

¹ Essential performance of this transmitter

Respiration Measurement

| | |
|----------------------------------------------------|---------------------------------------------|
| Measuring method: | Impedance method |
| Measuring lead: | Between R and F |
| Impedance range: | 2 k Ω or less |
| Respiration rate measuring accuracy ¹ : | ± 2 counts/min (at 0 to 150 counts/min) |
| Respiration rate counting range: | 0 to 150 counts/min |

¹ Essential performance of this transmitter

SpO₂ Measurement (ISO 9919: 2005 compliant)

| | |
|------------------------|----------------------------|
| Measuring range: | 0 to 100%SpO ₂ |
| Declared range: | 70 to 100%SpO ₂ |
| Minimum display range: | 1%SpO ₂ |
| Display update cycle: | Every 3 seconds |

| | |
|------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Measuring accuracy (rms) ¹ : | Accuracy assurance temperature: 18 to 40°C |
| Total accuracy including probe: | 80%SpO ₂ ≤ %SpO ₂ ≤ 100%SpO ₂ : ±2%SpO ₂ 70%SpO ₂ ≤ %SpO ₂ < 80%SpO ₂ : ±3%SpO ₂ under 70%SpO ₂ : not specified |
| Accuracy of the transmitter: | 80%SpO ₂ ≤ %SpO ₂ ≤ 100%SpO ₂ : ±1%SpO ₂ 50%SpO ₂ ≤ %SpO ₂ < 80%SpO ₂ : ±2%SpO ₂ under 50%SpO ₂ : not specified |
| Pulse rate measuring range: | 30 to 300 bpm |
| Pulse rate display range: | 30 to 300 bpm |
| Pulse rate accuracy (rms) ¹ : | ±3% ±1 bpm |


¹ Essential performance of this transmitter


- NOTE**
- The SpO₂ measuring accuracy was tested using the TL-201T, TL-260T, TL-271T and TL-631T SpO₂ probes. The testing was performed during induced hypoxia on healthy volunteers (Ethnicity: 10 Caucasians, 2 Africans, 1 Asian and 3 Indians), (Skin: 8 light, 4 medium, 4 dark), (Age: 21 to 34), (5 women and 11 men) under the condition of no motion. Arterial blood was sampled and measured by a CO-oximeter. The difference between SpO₂ measured by the SpO₂ probe and functional SaO₂ measured by a CO-oximeter was calculated using the root-mean-square (rms) method. This measurement accuracy figure represents 2/3 of all test measurements.
 - A pulse oximeter tester that generates simulated signals can be used to check the difference from the design specification, but it cannot be used as a replacement for human signals for testings accuracy.

Response time:


Selectable from “Slow”, “Normal” and “Fast”.

The following graphs show response time for 0.6%/s SpO₂ change.





The following graph shows response time for 10 bpm/s pulse rate change.



Noninvasive Blood Pressure, NIBP (IEC 60601-2-30: 1999 compliant)

| | |
|----------------------------------------------|--------------------------------------------------------------------------------------------------------|
| Measuring method: | Oscillometric |
| Measurement mode: | Manual, STAT (≤ 15 min), Periodic |
| Intended patient type: | Adult, child |
| Measuring range: | 0 to 300 mmHg |
| Pressure display range: | 0 to 300 mmHg |
| Measuring accuracy ¹ : | ± 3 mmHg ($0 \text{ mmHg} \leq \text{NIBP} \leq 300 \text{ mmHg}$) AAMI SP-10: 2002 compliant |
| Cuff inflation time: | ≤ 20 s (700 cc), 0 to 200 mmHg ≤ 15 s (70 cc), 0 to 200 mmHg |
| Pressure retention: | ≤ 5 mmHg (250 cc at 250 mmHg inflation for 10 seconds) |
| Air leakage: | ≤ 3 mmHg/min (700 cc at 300 mmHg inflation) |
| Maximum temperature at cuff attachment site: | 43°C (109°F) |
| Power discontinuity: | Deflate immediately after power down |

| | | |
|---------|------------------------------------------------------|-----------------|
| Safety: | Maximum pressurization value cuff inflation limiter: | 300 to 330 mmHg |
| | Cuff inflation time limiter: | ≤ 180 s |
| | Interval time limiter: | ≤ 30 s |

¹ Essential performance of this transmitter

Transmitter

| | |
|-------------------------------|--------------------------------------------------------------------|
| FCC regulation: | FCC part 95 Subpart H Wireless Medical Telemetry Service (WMTS) |
| Field strength limits: | < 200 mV/m (at 3 m) |
| Undesired emissions: | |
| FCC part 95 95.2379 (a), (b): | ≤ 960 MHz: 200 μV/m (at 3 m) ≥ 960 MHz: 500 μV/m (at 3 m) |
| Antenna: | Internal |
| Transmission channel: | Indicated on the transmitter |
| Transmission frequency range: | 608.0250 to 613.9750 MHz |
| Channel spacing: | 50 kHz or 37.5 kHz (12.5 kHz when interleaved) |
| Modulation: | FSK (frequency shift keying) |
| Type of emission: | F1D |
| Occupied bandwidth: | < 20 kHz |
| Effective radiated power: | 1.0 mW |

Power Requirements

| | |
|------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Rated voltage: | 3.6 V |
| Operating voltage: | 3.2 to 4.8 V |
| Battery type: | Three AA (R6) type NiMH secondary batteries Three AA (R6) type alkaline dry cell primary batteries |
| Battery lifetime (with alkaline batteries, at room temperature): | approximately 1 day (measuring ECG, respiration, SpO ₂ of approximately 60 kg weight adult male patient at the index finger, NIBP at 60 minute intervals) |

Dimensions and Weight

| | |
|-------------|---------------------------------------------------------------------------------------------------------------------------|
| Dimensions: | 114 W × 125 H × 63 D (mm) |
| Weight: | about 340 g (excluding batteries and other accessories) about 410 g (including batteries, excluding other accessories) |

Environment

Operating environment

| | |
|-----------------------|---------------------------|
| Temperature: | 5 to 40°C (41 to 104°F) |
| Humidity: | 30 to 85% (noncondensing) |
| Atmospheric pressure: | 700 to 1060 hPa |

Storage and transport environment

| | |
|-----------------------|-----------------------------|
| Temperature: | -20 to +65°C (-4 to +149°F) |
| Humidity: | 10 to 95% |
| Atmospheric pressure: | 700 to 1060 hPa |

Safety Standards

| | |
|------------------|---------------------------------|
| Safety standard: | CAN/CSA-C22.2 No. 601-1 M90 |
| | CAN/CSA-C22.2 No. 601-1. 1S1-94 |
| | CAN/CSA-C22.2 No. 601-1. 1B-90 |
| | CAN/CSA-C22.2 No. 60601-2-49-04 |
| | CAN/CSA-C22.2 No. 60601-2-27-06 |
| | CAN/CSA-C22.2 No. 60601-2-30-02 |
| | IEC 60601-1:1988 |
| | IEC 60601-1 Amendment 1: 1991 |
| | IEC 60601-1 Amendment 2: 1995 |
| | IEC 60601-2-27: 2005 |
| | IEC 60601-2-30: 1999 |
| | IEC 60601-2-49: 2001 |
| | ISO 9919: 2005 |
| | UL 60601-1: 2003 |

Type of protection against electrical shock:

INTERNALLY POWERED EQUIPMENT

Degree of protection against electrical shock:

ECG and impedance method respiration:

DEFIBRILLATION-PROOF TYPE CF APPLIED PART

SpO₂ and NIBP:

DEFIBRILLATION-PROOF TYPE BF APPLIED PART

Degree of protection against harmful ingress of water:

IPX0 (non-protected)

Degree of safety of application in the presence of a FLAMMABLE ANAESTHETIC MIXTURE WITH AIR, OR WITH OXYGEN OR NITROUS OXIDE:

Equipment not suitable for use in the presence of FLAMMABLE ANAESTHETIC MIXTURE WITH AIR, OR WITH OXYGEN OR NITROUS OXIDE

Mode of operation:

CONTINUOUS OPERATION

Electromagnetic Compatibility

IEC 60601-1-2: 2001

IEC 60601-1-2 Amendment 1: 2004

Electromagnetic Emissions

This ZM-540PA is intended for use in the electromagnetic environment specified below.

The customer or the user of the ZM-540PA should assure that it is used in such an environment.

| Emissions test | Compliance | Electromagnetic environment - guidance |
|-------------------------------------------------------------|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF emissions CISPR 11 | Group 1 | The ZM-540PA uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emissions CISPR 11 | Class B | The ZM-540PA is suitable for use in all establishments, including domestic establishments. |
| Harmonic emissions IEC 61000-3-2 | Not applicable | |
| Voltage fluctuations/ flicker emissions IEC 61000-3-3 | Not applicable | |


Electromagnetic Immunity

This ZM-540PA is intended for use in the electromagnetic environment specified below. The customer or the user of the ZM-540PA should assure that it is used in such an environment.

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
|--------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±6 kV contact ±8 kV air | ±6 kV contact ±8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast transient/ burst IEC 61000-4-4 | ±2 kV for power supply lines ±1 kV for input/output lines | Not applicable | — |
| Surge IEC 61000-4-5 | ±1 kV differential mode ±2 kV common mode | Not applicable | — |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | <5% U_T (>95% dip in U_T) for 0.5 cycles 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 s | Not applicable | — |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 3 A/m | 3 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |
| NOTE: U_T is the AC mains voltage prior to application of the test level. | | | |

Avoiding Electromagnetic Interference (Impedance Respiration)

Impedance respiration measurement is very sensitive and affected by electromagnetic interference. Technological limitations do not allow immunity levels higher than 1 V/m for radiated RF electromagnetic fields. Electromagnetic fields with field strengths above 1 V/m may cause measurement error. Do not use electrically radiating equipment near the impedance respiration measurements.

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p> | <p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p> | <p>3 Vrms</p> <p>3 V/m 80 MHz to 2.5 GHz</p> <p>(1 V/m 80 MHz to 2.5 GHz for respiration)</p> | <p>Portable and mobile RF communications equipment should be used no closer to any part of the ZM-540PA, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance $d = 1.2\sqrt{P}$</p> <p>$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>($d = 3.5\sqrt{P}$ 80 MHz to 800 MHz for respiration $d = 7.0\sqrt{P}$ 800 MHz to 2.5 GHz for respiration)</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey¹, should be less than the compliance level in each frequency range².</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p> |
| <p>NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.</p> | | | |
| <p>NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> | | | |
| <p>¹ Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ZM-540PA is used exceeds the applicable RF compliance level above, the ZM-540PA should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the ZM-540PA.</p> <p>² Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 1 V/m for respiration and 3 V/m for all other functions.</p> | | | |

Recommended Separation Distances between Portable and Mobile RF Communications Equipment

The ZM-540PA is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ZM-540PA can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ZM-540PA as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter (W) | Separation distance according to frequency of transmitter (m) | | |
|-----------------------------------------------|---------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| | 150 kHz to 80 MHz $d = 1.2\sqrt{P}$ | 80 MHz to 800 MHz $d = 1.2\sqrt{P}$ (For respiration: $d = 3.5\sqrt{P}$) | 800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$ (For respiration: $d = 7.0\sqrt{P}$) |
| 0.01 | 0.12 | 0.12 (0.35 ¹) | 0.23 (0.7 ¹) |
| 0.1 | 0.38 | 0.38 (1.1 ¹) | 0.73 (2.2 ¹) |
| 1 | 1.2 | 1.2 (3.5 ¹) | 2.3 (7.0 ¹) |
| 10 | 3.8 | 3.8 (11 ¹) | 7.3 (22 ¹) |
| 100 | 12 | 12 (35 ¹) | 23 (70 ¹) |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

(¹ For respiration)

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Recovery Time after Defibrillation

The transmitter returns to the normal operating mode within 10 seconds after defibrillation. The stored settings are not affected.

System Composition for EMC Test

The ZM-540PA transmitter is tested to comply with IEC 60601-1-2:2001 and Amendment 1:2004 with the following composition. If any part which is not specified by Nihon Kohden is used, the EMC specifications might not comply.

| Units | Cable Length |
|----------------------------|--------------|
| ZM-540PA transmitter | — |
| YP-503P NIBP cuff | 0.15 m |
| BR-906P ECG electrode lead | 0.8 m |
| TL-201T finger probe | 1.6 m |

ZM-541PA

Measured Parameters

| | |
|---------------|-------------------------------------------------------------------|
| Waveforms: | ECG, impedance method respiration, pulse |
| Numeric data: | Heart rate, respiration rate, SpO ₂ , NIBP, pulse rate |

Transmitted Data

| | |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Waveforms: | ECG, respiration, pulse wave |
| Numeric data: | SpO ₂ , NIBP, pulse rate |
| Status information: | Battery replacement, battery level ¹ , alarm suspended, pause monitoring ¹ , patient confirmed ¹ , ECG lead, pacing detection, electrode detachment, electrode impedance ¹ , ECG off ¹ , respiration method (impedance) ¹ , SpO ₂ status, NIBP status, channel ID, time constant (3.2 s), type of transmitter, transmitter code number ¹ , transmitter serial number ¹ |

¹ Transmitted only when the protocol is "57".

Display

| | |
|---------------|--------------------------|
| Display size: | 2.2 inch TFT color LCD |
| Viewing area: | 44.16 (H) × 33.12 (V) mm |
| Resolution: | 320 (H) × 240 (V) dots |

Displayed Data

| | |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Numeric and waveform screen: | ECG (one waveform from lead I, II, III, Va or Vb), heart rate, pulse rate, respiration rate, SpO ₂ , NIBP (systolic, diastolic, MAP), message, battery level, QRS/pulse sync mark, pulse bar graph, NIBP measurement mode and status information, ECG lead |
| Waveform review screen: | ECG or pulse wave of past 10 minutes |
| Numeric review screen: | Heart rate or pulse rate, respiration rate and SpO ₂ at 1 minute interval for past 10 minutes |
| CHECK ELECTRODE screen: | ECG for checking electrode attachment |

ECG

ECG Measurement

| | |
|---------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Channels: | 4 |
| Input dynamic range: | ±10 mV or more |
| Electrode offset potential tolerance: | ±500 mV or more |
| Input impedance: | 5 MΩ or more |
| Common mode rejection ratio: | 95 dB or more |
| Pacing pulse detection: | IEC 60601-2-27 50.102.10 compliant amplitude ±2 to 700 mV, duration 0.1 to 2 ms IEC 60601-2-27: 2005 compliant Based upon pacemaker pulse rejection capability |
| Defibrillation-proof: | ECG input protected against 400 Ws/DC 5 kV IEC 60601-2-27 17.101 compliant |

| | |
|----------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ECG recovery time after defibrillation: | within 10 s |
| Electrode condition: | Displays CHECK ELECTRODES message |
| Tall T-wave rejection capability: | Complies with the heights of T-waves from 0 to 1.6 mV IEC 60601-2-27: 2005 50.102.17 compliant |
| Pacemaker pulse rejection capability, without overshoot: | Complies with the amplitudes of pacemaker pulses ± 2 to ± 700 mV and widths 0.1 to 2 ms (As defined by method B of IEC 60601-2-27: 2005 50.102.13, this corresponds to the pacemaker pulses amplitudes and widths of ± 2.8 mV/2 ms to amplitudes ± 45.6 mV/0.1 ms.) |
| Pacemaker pulse rejection capability, with overshoot: | Overshoot amplitudes and time constants of ± 0.12 mV/100 ms to ± 2 mV/4 ms (As defined by method B of IEC 60601-2-27: 2005 50.102.13, this corresponds to the pacemaker pulses amplitudes and widths of ± 2.8 mV/2 ms to amplitudes ± 45.6 mV/0.1 ms.) |

ECG Display and Heart Rate Count

| | |
|---------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| Frequency characteristic: | filter on: 1 to 18 Hz, filter off: 0.05 to 60 Hz |
| Heart rate detection method: | Average |
| QRS detection: | 70 to 120 ms: amplitude ≥ 0.5 mV, rate 30 to 200 beats/min 40 to 120 ms: amplitude ≥ 0.5 mV, rate 30 to 250 beats/min |
| Heart rate counting range: | 0, 15 to 300 beats/min |
| Heart rate counting accuracy ¹ : | ± 2 beats/min, (0, 15 to 300 beats/min) |

¹ Essential performance of this transmitter

Respiration Measurement

| | |
|----------------------------------------------------|---------------------------------------------|
| Measuring method: | Impedance method |
| Measuring lead: | Between R and F |
| Impedance range: | 2 k Ω or less |
| Respiration rate measuring accuracy ¹ : | ± 2 counts/min (at 0 to 150 counts/min) |
| Respiration rate counting range: | 0 to 150 counts/min |

¹ Essential performance of this transmitter

SpO₂ Measurement (ISO 9919: 2005 compliant)

| | |
|-----------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Measuring range: | 0 to 100%SpO ₂ |
| Declared range: | 70 to 100%SpO ₂ |
| Minimum display range: | 1%SpO ₂ |
| Display update cycle: | Every 3 seconds |
| Measuring accuracy (rms) ¹ : | Accuracy assurance temperature: 18 to 40°C |
| Total accuracy including probe: | 80%SpO ₂ \leq %SpO ₂ \leq 100%SpO ₂ : $\pm 2\%$ SpO ₂ 70%SpO ₂ \leq %SpO ₂ $<$ 80%SpO ₂ : $\pm 3\%$ SpO ₂ under 70%SpO ₂ : not specified |
| Accuracy of the transmitter: | 80%SpO ₂ \leq %SpO ₂ \leq 100%SpO ₂ : $\pm 1\%$ SpO ₂ 50%SpO ₂ \leq %SpO ₂ $<$ 80%SpO ₂ : $\pm 2\%$ SpO ₂ under 50%SpO ₂ : not specified |

Pulse rate measuring range: 30 to 300 bpm
 Pulse rate display range: 30 to 300 bpm
 Pulse rate accuracy (rms)¹: $\pm 3\% \pm 1$ bpm



¹ Essential performance of this transmitter

- NOTE**
- The SpO₂ measuring accuracy was tested using the TL-201T, TL-260T, TL-271T and TL-631T SpO₂ probes. The testing was performed during induced hypoxia on healthy volunteers (Ethnicity: 10 Caucasians, 2 Africans, 1 Asian and 3 Indians), (Skin: 8 light, 4 medium, 4 dark), (Age: 21 to 34), (5 women and 11 men) under the condition of no motion. Arterial blood was sampled and measured by a CO-oximeter. The difference between SpO₂ measured by the SpO₂ probe and functional SaO₂ measured by a CO-oximeter was calculated using the root-mean-square (rms) method according to ISO 9919: 2005. This measurement accuracy figure represents 2/3 of all test measurements.
 - A pulse oximeter tester that generates simulated signals can be used to check the difference from the design specification, but it cannot be used as a replacement for human signals for testings accuracy.


Response time:

Selectable from “Slow”, “Normal” and “Fast”.

The following graphs show response time for 0.6%/s SpO₂ change.



The following graph shows the response time example when pulse rate changes 10 bpm/s.



Noninvasive Blood Pressure, NIBP (IEC 60601-2-30: 1999 compliant)

| | |
|-----------------------------------|--------------------------------------------------------------------------------------------------------|
| Measuring method: | Oscillometric |
| Measurement mode: | Manual, STAT (≤ 15 min), Periodic |
| Intended patient type: | Adult, child |
| Measuring range: | 0 to 300 mmHg |
| Pressure display range: | 0 to 300 mmHg |
| Measuring accuracy ¹ : | ± 3 mmHg ($0 \text{ mmHg} \leq \text{NIBP} \leq 300 \text{ mmHg}$) AAMI SP-10: 2002 compliant |
| Cuff inflation time: | ≤ 20 s (700 cc), 0 to 200 mmHg ≤ 15 s (70 cc), 0 to 200 mmHg |
| Pressure retention: | ≤ 5 mmHg (250 cc at 250 mmHg inflation for 10 seconds) |
| Air leakage: | ≤ 3 mmHg/min (700 cc at 300 mmHg inflation) |
| Power discontinuity: | Deflate immediately after power down |
| Safety: | Maximum pressurization value cuff inflation limiter: 300 to 330 mmHg |
| | Cuff inflation time limiter: ≤ 180 s |
| | Interval time limiter: ≤ 30 s |

¹ Essential performance of this transmitter

Transmitter

| | |
|-------------------------------|----------------------------------------------------------------------------------------------------|
| FCC regulation: | FCC part 95 Subpart H Wireless Medical Telemetry Service (WMTS) |
| Field strength limits: | < 740 mV/m (at 3 m) |
| Undesired emissions: | |
| FCC part 95 95.2379 (a), (b): | below 960 MHz: < 200 $\mu\text{V/m}$ (at 3 m) above 960 MHz: < 500 $\mu\text{V/m}$ (at 3 m) |
| Antenna: | Internal |
| Transmission channel: | Indicated on the transmitter |
| Transmission frequency range: | 1395.0250 to 1399.9750 MHz 1427.0250 to 1431.9750 MHz |
| Channel spacing: | 50 kHz or 37.5 kHz (12.5 kHz when interleaved) |

| | |
|---------------------------|--------------------------------------|
| Modulation: | FSK (frequency shift keying) |
| Type of emission: | F1D |
| Occupied bandwidth: | < 20 kHz |
| Effective radiated power: | 5.0 mW |
| | Can be changed to 1.0 mW if required |

Power Requirements

| | |
|------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Rated voltage: | 3.6 V |
| Operating voltage: | 3.2 to 4.8 V |
| Battery type: | Three AA (R6) type NiMH secondary batteries Three AA (R6) type alkaline dry cell primary batteries |
| Battery lifetime (with alkaline batteries, at room temperature): | approximately 1 day (measuring ECG, respiration, SpO ₂ of approximately 60 kg weight adult male patient at the index finger, NIBP at 60 minute intervals) |

Dimensions and Weight

| | |
|-------------|---------------------------------------------------------------------------------------------------------------------------|
| Dimensions: | 114 W × 125 H × 63 D (mm) |
| Weight: | about 340 g (excluding batteries and other accessories) about 410 g (including batteries, excluding other accessories) |

Environment

Operating environment

| | |
|-----------------------|---------------------------|
| Temperature: | 5 to 40°C (41 to 104°F) |
| Humidity: | 30 to 85% (noncondensing) |
| Atmospheric pressure: | 700 to 1060 hPa |

Storage and transport environment

| | |
|-----------------------|-----------------------------|
| Temperature: | -20 to +65°C (-4 to +149°F) |
| Humidity: | 10 to 95% |
| Atmospheric pressure: | 700 to 1060 hPa |

Safety Standards

| | |
|------------------|---------------------------------|
| Safety standard: | CAN/CSA-C22.2 No. 601-1 M90 |
| | CAN/CSA-C22.2 No. 601-1. 1S1-94 |
| | CAN/CSA-C22.2 No. 601-1. 1B-90 |
| | CAN/CSA-C22.2 No. 60601-2-49-04 |
| | CAN/CSA-C22.2 No. 60601-2-27-06 |
| | CAN/CSA-C22.2 No. 60601-2-30-02 |
| | IEC 60601-1:1988 |
| | IEC 60601-1 Amendment 1: 1991 |
| | IEC 60601-1 Amendment 2: 1995 |
| | IEC 60601-2-27: 2005 |
| | IEC 60601-2-30: 1999 |
| | IEC 60601-2-49: 2001 |

ISO 9919: 2005
 UL 60601-1: 2003

Type of protection against electrical shock:

INTERNALLY POWERED EQUIPMENT

Degree of protection against electrical shock:

ECG and impedance method respiration:

DEFIBRILLATION-PROOF TYPE CF APPLIED PART

SpO₂ and NIBP:

DEFIBRILLATION-PROOF TYPE BF APPLIED PART

Degree of protection against harmful ingress of water:

IPX0 (non-protected)

Degree of safety of application in the presence of a FLAMMABLE ANAESTHETIC MIXTURE WITH AIR, OR WITH OXYGEN OR NITROUS OXIDE:

Equipment not suitable for use in the presence of FLAMMABLE ANAESTHETIC MIXTURE WITH AIR, OR WITH OXYGEN OR NITROUS OXIDE
 CONTINUOUS OPERATION

Mode of operation:

Electromagnetic Compatibility

IEC 60601-1-2: 2001

IEC 60601-1-2 Amendment 1: 2004

Electromagnetic Emissions

This ZM-541PA is intended for use in the electromagnetic environment specified below.

The customer or the user of the ZM-541PA should assure that it is used in such an environment.

| Emissions test | Compliance | Electromagnetic environment - guidance |
|-------------------------------------------------------------|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF emissions CISPR 11 | Group 1 | The ZM-541PA uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emissions CISPR 11 | Class B | |
| Harmonic emissions IEC 61000-3-2 | Not applicable | |
| Voltage fluctuations/ flicker emissions IEC 61000-3-3 | Not applicable | |


Electromagnetic Immunity

This ZM-541PA is intended for use in the electromagnetic environment specified below. The customer or the user of the ZM-541PA should assure that it is used in such an environment.

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
|--------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±6 kV contact ±8 kV air | ±6 kV contact ±8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast transient/ burst IEC 61000-4-4 | ±2 kV for power supply lines ±1 kV for input/output lines | Not applicable | — |
| Surge IEC 61000-4-5 | ±1 kV differential mode ±2 kV common mode | Not applicable | — |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | <5% U_T (>95% dip in U_T) for 0.5 cycles 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 s | Not applicable | — |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 3 A/m | 3 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |
| NOTE: U_T is the AC mains voltage prior to application of the test level. | | | |

Avoiding Electromagnetic Interference (Impedance Respiration)

Impedance respiration measurement is very sensitive and affected by electromagnetic interference. Technological limitations do not allow immunity levels higher than 1 V/m for radiated RF electromagnetic fields. Electromagnetic fields with field strengths above 1 V/m may cause measurement error. Do not use electrically radiating equipment near the impedance respiration measurements.

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
|----------------------------------------------------------------------------|-----------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p> | <p>3 V_{rms} 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p> | <p>3 V_{rms}</p> <p>3 V/m 80 MHz to 2.5 GHz</p> <p>(1 V/m 80 MHz to 2.5 GHz for respiration)</p> | <p>Portable and mobile RF communications equipment should be used no closer to any part of the ZM-541PA, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p>$d = 1.2\sqrt{P}$</p> <p>$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>($d = 3.5\sqrt{P}$ 80 MHz to 800 MHz for respiration $d = 7.0\sqrt{P}$ 800 MHz to 2.5 GHz for respiration)</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey¹, should be less than the compliance level in each frequency range².</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

¹ Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ZM-541PA is used exceeds the applicable RF compliance level above, the ZM-541PA should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the ZM-541PA.

² Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 1 V/m for respiration and 3 V/m for all other functions.

Recommended Separation Distances between Portable and Mobile RF Communications Equipment

The ZM-541PA is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ZM-541PA can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ZM-541PA as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter (W) | Separation distance according to frequency of transmitter (m) | | |
|-----------------------------------------------|---------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| | 150 kHz to 80 MHz $d = 1.2\sqrt{P}$ | 80 MHz to 800 MHz $d = 1.2\sqrt{P}$ (For respiration: $d = 3.5\sqrt{P}$) | 800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$ (For respiration: $d = 7.0\sqrt{P}$) |
| 0.01 | 0.12 | 0.12 (0.35 ¹) | 0.23 (0.7 ¹) |
| 0.1 | 0.38 | 0.38 (1.1 ¹) | 0.73 (2.2 ¹) |
| 1 | 1.2 | 1.2 (3.5 ¹) | 2.3 (7.0 ¹) |
| 10 | 3.8 | 3.8 (11 ¹) | 7.3 (22 ¹) |
| 100 | 12 | 12 (35 ¹) | 23 (70 ¹) |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

(¹ For respiration)

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Recovery Time after Defibrillation


The transmitter returns to the normal operating mode within 10 seconds after defibrillation. The stored settings are not affected.

System Composition for EMC Test

The ZM-541PA transmitter is tested to comply with IEC 60601-1-2:2001 and Amendment 1:2004 with the following composition. If any part which is not specified by Nihon Kohden is used, the EMC specifications might not comply.

| Units | Cable Length |
|----------------------------|--------------|
| ZM-541PA transmitter | — |
| YP-503P NIBP cuff | 0.15 m |
| BR-906P ECG electrode lead | 0.8 m |
| TL-201T finger probe | 1.6 m |

Replaceable Part



| Name | Code No. | Qty |
|--------------------|-------------|-----|
| Battery case cover | 6143-903101 | 1 |

Options

⚠ CAUTION

Only use Nihon Kohden specified electrodes, electrode leads and SpO₂ probes. Otherwise, the maximum performance from the transmitter cannot be guaranteed.

Transmitter

- QI-901PK channel writer

ECG/RESP

| Name | Type | Length (m) | Model | Qty | Supply Code |
|----------------|------------------------------|------------|----------|-----|-------------|
| Electrode lead | 3 electrodes, AHA, clip type | 0.8 | BR-903PA | 1 | K911A |
| | 3 electrodes, AHA, clip type | 0.8 | BR-933PA | 1 | K903B |
| | 3 electrodes, AHA, hook type | 0.8 | BR-943PA | 1 | K904B |
| | 6 electrodes, AHA, clip type | 0.8 | BR-936PA | 1 | K903D |
| | 6 electrodes, AHA, hook type | 0.8 | BR-946PA | 1 | K904D |

SpO₂

| Name | Lead Length (m) | Model | Qty | Supply Code | |
|--------------------------------------------------------|-----------------|----------|--------|--------------------------|-------|
| Finger probe (reusable) | 0.6 | TL-201T | 1 | P225H | |
| | 1.6 | | | P225F | |
| Multi-site probe (reusable) | 1.6 | TL-220T | | P225G | |
| Finger probe (reusable) | 0.6 | TL-631T1 | | P311A | |
| | 1.6 | TL-631T3 | | P311C | |
| SpO ₂ probe (for adult, disposable) | 0.8 | TL-271T | 24 | P203A | |
| | 1.6 | TL-271T3 | | P203E | |
| SpO ₂ probe (for child, disposable) | 0.8 | TL-272T | | P203B | |
| | 1.6 | TL-272T3 | | P203F | |
| SpO ₂ probe (for neonate/adult, disposable) | 0.8 | TL-273T | | P203C | |
| | 1.6 | TL-273T3 | | P203G | |
| SpO ₂ probe (for infant, disposable) | 0.8 | TL-274T | | P203D | |
| | 1.6 | TL-274T3 | | P203H | |
| SpO ₂ probe (for neonate/adult, disposable) | 0.8 | TL-051S | | 5 | P228A |
| | 1.6 | TL-052S | | | P228B |
| SpO ₂ probe (for child/infant, disposable) | 0.8 | TL-061S | P229A | | |
| | 1.6 | TL-062S | P229B | | |
| TL-05X TL-06X foam tape | — | — | 4 × 25 | P260, P260D ¹ | |
| BLUPRO attachment tape | — | — | 3 × 30 | P263, P263A ¹ | |
| Probe fastener | — | YS-093P2 | 30 | P267 | |

¹ There is more than one supply code. The supply code to be used depends on the country where the probe is used. Contact your Nihon Kohden representative for further details.

NIBP

| Name | | Width (cm) | Length (m) | Qty | Model | Supply Code |
|---------------------------------------|--------------|------------|------------|-----|---------|-------------|
| Extension hose | | — | 1.5 | 1 | YN-990P | S903 |
| Reusable cuff | | | | | | |
| NIBP cuff, infant, YAWARA CUFF2 | | 5 | 0.15 | 1 | YP-710T | S951A |
| NIBP cuff, child, YAWARA CUFF2 | | 7 | 0.15 | 1 | YP-711T | S951B |
| NIBP cuff, adult, small, YAWARA CUFF2 | Small | 10 | 0.15 | 1 | YP-712T | S951C |
| NIBP cuff, adult, small, YAWARA CUFF2 | Standard | 13 | 0.15 | 1 | YP-713T | S951D |
| NIBP cuff, adult, small, YAWARA CUFF2 | Large | 16 | 0.15 | 1 | YP-714T | S951E |
| Disposable cuff | | | | | | |
| Disposable cuff, infant | | 6 | 0.17 | 20 | YP-810P | S945C |
| Disposable cuff, child | | 8 | 0.17 | 20 | YP-811P | S945D |
| Disposable cuff, adult | Small | 10 | 0.17 | 20 | YP-812P | S946E |
| | Standard | 14 | 0.17 | 20 | YP-813P | S946F |
| | Medium large | 15 | 0.17 | 20 | YP-814P | S946G |
| | Large | 17 | 0.17 | 20 | YP-815P | S946H |
| | Extra large | 18 | 0.17 | 20 | YP-816P | S946I |

Transmission Frequencies

Channel: 9002 to 9478

| Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) |
|---------------|-----------------|---------------|-----------------|---------------|-----------------|
| 9002 | 608.0250 | 9040 | 608.5000 | 9078 | 608.9750 |
| 9003 | 608.0375 | 9041 | 608.5125 | 9079 | 608.9875 |
| 9004 | 608.0500 | 9042 | 608.5250 | 9080 | 609.0000 |
| 9005 | 608.0625 | 9043 | 608.5375 | 9081 | 609.0125 |
| 9006 | 608.0750 | 9044 | 608.5500 | 9082 | 609.0250 |
| 9007 | 608.0875 | 9045 | 608.5625 | 9083 | 609.0375 |
| 9008 | 608.1000 | 9046 | 608.5750 | 9084 | 609.0500 |
| 9009 | 608.1125 | 9047 | 608.5875 | 9085 | 609.0625 |
| 9010 | 608.1250 | 9048 | 608.6000 | 9086 | 609.0750 |
| 9011 | 608.1375 | 9049 | 608.6125 | 9087 | 609.0875 |
| 9012 | 608.1500 | 9050 | 608.6250 | 9088 | 609.1000 |
| 9013 | 608.1625 | 9051 | 608.6375 | 9089 | 609.1125 |
| 9014 | 608.1750 | 9052 | 608.6500 | 9090 | 609.1250 |
| 9015 | 608.1875 | 9053 | 608.6625 | 9091 | 609.1375 |
| 9016 | 608.2000 | 9054 | 608.6750 | 9092 | 609.1500 |
| 9017 | 608.2125 | 9055 | 608.6875 | 9093 | 609.1625 |
| 9018 | 608.2250 | 9056 | 608.7000 | 9094 | 609.1750 |
| 9019 | 608.2375 | 9057 | 608.7125 | 9095 | 609.1875 |
| 9020 | 608.2500 | 9058 | 608.7250 | 9096 | 609.2000 |
| 9021 | 608.2625 | 9059 | 608.7375 | 9097 | 609.2125 |
| 9022 | 608.2750 | 9060 | 608.7500 | 9098 | 609.2250 |
| 9023 | 608.2875 | 9061 | 608.7625 | 9099 | 609.2375 |
| 9024 | 608.3000 | 9062 | 608.7750 | 9100 | 609.2500 |
| 9025 | 608.3125 | 9063 | 608.7875 | 9101 | 609.2625 |
| 9026 | 608.3250 | 9064 | 608.8000 | 9102 | 609.2750 |
| 9027 | 608.3375 | 9065 | 608.8125 | 9103 | 609.2875 |
| 9028 | 608.3500 | 9066 | 608.8250 | 9104 | 609.3000 |
| 9029 | 608.3625 | 9067 | 608.8375 | 9105 | 609.3125 |
| 9030 | 608.3750 | 9068 | 608.8500 | 9106 | 609.3250 |
| 9031 | 608.3875 | 9069 | 608.8625 | 9107 | 609.3375 |
| 9032 | 608.4000 | 9070 | 608.8750 | 9108 | 609.3500 |
| 9033 | 608.4125 | 9071 | 608.8875 | 9109 | 609.3625 |
| 9034 | 608.4250 | 9072 | 608.9000 | 9110 | 609.3750 |
| 9035 | 608.4375 | 9073 | 608.9125 | 9111 | 609.3875 |
| 9036 | 608.4500 | 9074 | 608.9250 | 9112 | 609.4000 |
| 9037 | 608.4625 | 9075 | 608.9375 | 9113 | 609.4125 |
| 9038 | 608.4750 | 9076 | 608.9500 | 9114 | 609.4250 |
| 9039 | 608.4875 | 9077 | 608.9625 | 9115 | 609.4375 |

| Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) |
|---------------|-----------------|---------------|-----------------|---------------|-----------------|
| 9116 | 609.4500 | 9159 | 609.9875 | 9202 | 610.5250 |
| 9117 | 609.4625 | 9160 | 610.0000 | 9203 | 610.5375 |
| 9118 | 609.4750 | 9161 | 610.0125 | 9204 | 610.5500 |
| 9119 | 609.4875 | 9162 | 610.0250 | 9205 | 610.5625 |
| 9120 | 609.5000 | 9163 | 610.0375 | 9206 | 610.5750 |
| 9121 | 609.5125 | 9164 | 610.0500 | 9207 | 610.5875 |
| 9122 | 609.5250 | 9165 | 610.0625 | 9208 | 610.6000 |
| 9123 | 609.5375 | 9166 | 610.0750 | 9209 | 610.6125 |
| 9124 | 609.5500 | 9167 | 610.0875 | 9210 | 610.6250 |
| 9125 | 609.5625 | 9168 | 610.1000 | 9211 | 610.6375 |
| 9126 | 609.5750 | 9169 | 610.1125 | 9212 | 610.6500 |
| 9127 | 609.5875 | 9170 | 610.1250 | 9213 | 610.6625 |
| 9128 | 609.6000 | 9171 | 610.1375 | 9214 | 610.6750 |
| 9129 | 609.6125 | 9172 | 610.1500 | 9215 | 610.6875 |
| 9130 | 609.6250 | 9173 | 610.1625 | 9216 | 610.7000 |
| 9131 | 609.6375 | 9174 | 610.1750 | 9217 | 610.7125 |
| 9132 | 609.6500 | 9175 | 610.1875 | 9218 | 610.7250 |
| 9133 | 609.6625 | 9176 | 610.2000 | 9219 | 610.7375 |
| 9134 | 609.6750 | 9177 | 610.2125 | 9220 | 610.7500 |
| 9135 | 609.6875 | 9178 | 610.2250 | 9221 | 610.7625 |
| 9136 | 609.7000 | 9179 | 610.2375 | 9222 | 610.7750 |
| 9137 | 609.7125 | 9180 | 610.2500 | 9223 | 610.7875 |
| 9138 | 609.7250 | 9181 | 610.2625 | 9224 | 610.8000 |
| 9139 | 609.7375 | 9182 | 610.2750 | 9225 | 610.8125 |
| 9140 | 609.7500 | 9183 | 610.2875 | 9226 | 610.8250 |
| 9141 | 609.7625 | 9184 | 610.3000 | 9227 | 610.8375 |
| 9142 | 609.7750 | 9185 | 610.3125 | 9228 | 610.8500 |
| 9143 | 609.7875 | 9186 | 610.3250 | 9229 | 610.8625 |
| 9144 | 609.8000 | 9187 | 610.3375 | 9230 | 610.8750 |
| 9145 | 609.8125 | 9188 | 610.3500 | 9231 | 610.8875 |
| 9146 | 609.8250 | 9189 | 610.3625 | 9232 | 610.9000 |
| 9147 | 609.8375 | 9190 | 610.3750 | 9233 | 610.9125 |
| 9148 | 609.8500 | 9191 | 610.3875 | 9234 | 610.9250 |
| 9149 | 609.8625 | 9192 | 610.4000 | 9235 | 610.9375 |
| 9150 | 609.8750 | 9193 | 610.4125 | 9236 | 610.9500 |
| 9151 | 609.8875 | 9194 | 610.4250 | 9237 | 610.9625 |
| 9152 | 609.9000 | 9195 | 610.4375 | 9238 | 610.9750 |
| 9153 | 609.9125 | 9196 | 610.4500 | 9239 | 610.9875 |
| 9154 | 609.9250 | 9197 | 610.4625 | 9240 | 611.0000 |
| 9155 | 609.9375 | 9198 | 610.4750 | 9241 | 611.0125 |
| 9156 | 609.9500 | 9199 | 610.4875 | 9242 | 611.0250 |
| 9157 | 609.9625 | 9200 | 610.5000 | 9243 | 611.0375 |
| 9158 | 609.9750 | 9201 | 610.5125 | 9244 | 611.0500 |

| Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) |
|---------------|-----------------|---------------|-----------------|---------------|-----------------|
| 9245 | 611.0625 | 9288 | 611.6000 | 9331 | 612.1375 |
| 9246 | 611.0750 | 9289 | 611.6125 | 9332 | 612.1500 |
| 9247 | 611.0875 | 9290 | 611.6250 | 9333 | 612.1625 |
| 9248 | 611.1000 | 9291 | 611.6375 | 9334 | 612.1750 |
| 9249 | 611.1125 | 9292 | 611.6500 | 9335 | 612.1875 |
| 9250 | 611.1250 | 9293 | 611.6625 | 9336 | 612.2000 |
| 9251 | 611.1375 | 9294 | 611.6750 | 9337 | 612.2125 |
| 9252 | 611.1500 | 9295 | 611.6875 | 9338 | 612.2250 |
| 9253 | 611.1625 | 9296 | 611.7000 | 9339 | 612.2375 |
| 9254 | 611.1750 | 9297 | 611.7125 | 9340 | 612.2500 |
| 9255 | 611.1875 | 9298 | 611.7250 | 9341 | 612.2625 |
| 9256 | 611.2000 | 9299 | 611.7375 | 9342 | 612.2750 |
| 9257 | 611.2125 | 9300 | 611.7500 | 9343 | 612.2875 |
| 9258 | 611.2250 | 9301 | 611.7625 | 9344 | 612.3000 |
| 9259 | 611.2375 | 9302 | 611.7750 | 9345 | 612.3125 |
| 9260 | 611.2500 | 9303 | 611.7875 | 9346 | 612.3250 |
| 9261 | 611.2625 | 9304 | 611.8000 | 9347 | 612.3375 |
| 9262 | 611.2750 | 9305 | 611.8125 | 9348 | 612.3500 |
| 9263 | 611.2875 | 9306 | 611.8250 | 9349 | 612.3625 |
| 9264 | 611.3000 | 9307 | 611.8375 | 9350 | 612.3750 |
| 9265 | 611.3125 | 9308 | 611.8500 | 9351 | 612.3875 |
| 9266 | 611.3250 | 9309 | 611.8625 | 9352 | 612.4000 |
| 9267 | 611.3375 | 9310 | 611.8750 | 9353 | 612.4125 |
| 9268 | 611.3500 | 9311 | 611.8875 | 9354 | 612.4250 |
| 9269 | 611.3625 | 9312 | 611.9000 | 9355 | 612.4375 |
| 9270 | 611.3750 | 9313 | 611.9125 | 9356 | 612.4500 |
| 9271 | 611.3875 | 9314 | 611.9250 | 9357 | 612.4625 |
| 9272 | 611.4000 | 9315 | 611.9375 | 9358 | 612.4750 |
| 9273 | 611.4125 | 9316 | 611.9500 | 9359 | 612.4875 |
| 9274 | 611.4250 | 9317 | 611.9625 | 9360 | 612.5000 |
| 9275 | 611.4375 | 9318 | 611.9750 | 9361 | 612.5125 |
| 9276 | 611.4500 | 9319 | 611.9875 | 9362 | 612.5250 |
| 9277 | 611.4625 | 9320 | 612.0000 | 9363 | 612.5375 |
| 9278 | 611.4750 | 9321 | 612.0125 | 9364 | 612.5500 |
| 9279 | 611.4875 | 9322 | 612.0250 | 9365 | 612.5625 |
| 9280 | 611.5000 | 9323 | 612.0375 | 9366 | 612.5750 |
| 9281 | 611.5125 | 9324 | 612.0500 | 9367 | 612.5875 |
| 9282 | 611.5250 | 9325 | 612.0625 | 9368 | 612.6000 |
| 9283 | 611.5375 | 9326 | 612.0750 | 9369 | 612.6125 |
| 9284 | 611.5500 | 9327 | 612.0875 | 9370 | 612.6250 |
| 9285 | 611.5625 | 9328 | 612.1000 | 9371 | 612.6375 |
| 9286 | 611.5750 | 9329 | 612.1125 | 9372 | 612.6500 |
| 9287 | 611.5875 | 9330 | 612.1250 | 9373 | 612.6625 |

| Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) |
|---------------|-----------------|---------------|-----------------|---------------|-----------------|
| 9374 | 612.6750 | 9417 | 613.2125 | 9460 | 613.7500 |
| 9375 | 612.6875 | 9418 | 613.2250 | 9461 | 613.7625 |
| 9376 | 612.7000 | 9419 | 613.2375 | 9462 | 613.7750 |
| 9377 | 612.7125 | 9420 | 613.2500 | 9463 | 613.7875 |
| 9378 | 612.7250 | 9421 | 613.2625 | 9464 | 613.8000 |
| 9379 | 612.7375 | 9422 | 613.2750 | 9465 | 613.8125 |
| 9380 | 612.7500 | 9423 | 613.2875 | 9466 | 613.8250 |
| 9381 | 612.7625 | 9424 | 613.3000 | 9467 | 613.8375 |
| 9382 | 612.7750 | 9425 | 613.3125 | 9468 | 613.8500 |
| 9383 | 612.7875 | 9426 | 613.3250 | 9469 | 613.8625 |
| 9384 | 612.8000 | 9427 | 613.3375 | 9470 | 613.8750 |
| 9385 | 612.8125 | 9428 | 613.3500 | 9471 | 613.8875 |
| 9386 | 612.8250 | 9429 | 613.3625 | 9472 | 613.9000 |
| 9387 | 612.8375 | 9430 | 613.3750 | 9473 | 613.9125 |
| 9388 | 612.8500 | 9431 | 613.3875 | 9474 | 613.9250 |
| 9389 | 612.8625 | 9432 | 613.4000 | 9475 | 613.9375 |
| 9390 | 612.8750 | 9433 | 613.4125 | 9476 | 613.9500 |
| 9391 | 612.8875 | 9434 | 613.4250 | 9477 | 613.9625 |
| 9392 | 612.9000 | 9435 | 613.4375 | 9478 | 613.9750 |
| 9393 | 612.9125 | 9436 | 613.4500 | | |
| 9394 | 612.9250 | 9437 | 613.4625 | | |
| 9395 | 612.9375 | 9438 | 613.4750 | | |
| 9396 | 612.9500 | 9439 | 613.4875 | | |
| 9397 | 612.9625 | 9440 | 613.5000 | | |
| 9398 | 612.9750 | 9441 | 613.5125 | | |
| 9399 | 612.9875 | 9442 | 613.5250 | | |
| 9400 | 613.0000 | 9443 | 613.5375 | | |
| 9401 | 613.0125 | 9444 | 613.5500 | | |
| 9402 | 613.0250 | 9445 | 613.5625 | | |
| 9403 | 613.0375 | 9446 | 613.5750 | | |
| 9404 | 613.0500 | 9447 | 613.5875 | | |
| 9405 | 613.0625 | 9448 | 613.6000 | | |
| 9406 | 613.0750 | 9449 | 613.6125 | | |
| 9407 | 613.0875 | 9450 | 613.6250 | | |
| 9408 | 613.1000 | 9451 | 613.6375 | | |
| 9409 | 613.1125 | 9452 | 613.6500 | | |
| 9410 | 613.1250 | 9453 | 613.6625 | | |
| 9411 | 613.1375 | 9454 | 613.6750 | | |
| 9412 | 613.1500 | 9455 | 613.6875 | | |
| 9413 | 613.1625 | 9456 | 613.7000 | | |
| 9414 | 613.1750 | 9457 | 613.7125 | | |
| 9415 | 613.1875 | 9458 | 613.7250 | | |
| 9416 | 613.2000 | 9459 | 613.7375 | | |

Channel: E002 to E398

| Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) |
|---------------|-----------------|---------------|-----------------|---------------|-----------------|
| E002 | 1395.0250 | E040 | 1395.5000 | E078 | 1395.9750 |
| E003 | 1395.0375 | E041 | 1395.5125 | E079 | 1395.9875 |
| E004 | 1395.0500 | E042 | 1395.5250 | E080 | 1396.0000 |
| E005 | 1395.0625 | E043 | 1395.5375 | E081 | 1396.0125 |
| E006 | 1395.0750 | E044 | 1395.5500 | E082 | 1396.0250 |
| E007 | 1395.0875 | E045 | 1395.5625 | E083 | 1396.0375 |
| E008 | 1395.1000 | E046 | 1395.5750 | E084 | 1396.0500 |
| E009 | 1395.1125 | E047 | 1395.5875 | E085 | 1396.0625 |
| E010 | 1395.1250 | E048 | 1395.6000 | E086 | 1396.0750 |
| E011 | 1395.1375 | E049 | 1395.6125 | E087 | 1396.0875 |
| E012 | 1395.1500 | E050 | 1395.6250 | E088 | 1396.1000 |
| E013 | 1395.1625 | E051 | 1395.6375 | E089 | 1396.1125 |
| E014 | 1395.1750 | E052 | 1395.6500 | E090 | 1396.1250 |
| E015 | 1395.1875 | E053 | 1395.6625 | E091 | 1396.1375 |
| E016 | 1395.2000 | E054 | 1395.6750 | E092 | 1396.1500 |
| E017 | 1395.2125 | E055 | 1395.6875 | E093 | 1396.1625 |
| E018 | 1395.2250 | E056 | 1395.7000 | E094 | 1396.1750 |
| E019 | 1395.2375 | E057 | 1395.7125 | E095 | 1396.1875 |
| E020 | 1395.2500 | E058 | 1395.7250 | E096 | 1396.2000 |
| E021 | 1395.2625 | E059 | 1395.7375 | E097 | 1396.2125 |
| E022 | 1395.2750 | E060 | 1395.7500 | E098 | 1396.2250 |
| E023 | 1395.2875 | E061 | 1395.7625 | E099 | 1396.2375 |
| E024 | 1395.3000 | E062 | 1395.7750 | E100 | 1396.2500 |
| E025 | 1395.3125 | E063 | 1395.7875 | E101 | 1396.2625 |
| E026 | 1395.3250 | E064 | 1395.8000 | E102 | 1396.2750 |
| E027 | 1395.3375 | E065 | 1395.8125 | E103 | 1396.2875 |
| E028 | 1395.3500 | E066 | 1395.8250 | E104 | 1396.3000 |
| E029 | 1395.3625 | E067 | 1395.8375 | E105 | 1396.3125 |
| E030 | 1395.3750 | E068 | 1395.8500 | E106 | 1396.3250 |
| E031 | 1395.3875 | E069 | 1395.8625 | E107 | 1396.3375 |
| E032 | 1395.4000 | E070 | 1395.8750 | E108 | 1396.3500 |
| E033 | 1395.4125 | E071 | 1395.8875 | E109 | 1396.3625 |
| E034 | 1395.4250 | E072 | 1395.9000 | E110 | 1396.3750 |
| E035 | 1395.4375 | E073 | 1395.9125 | E111 | 1396.3875 |
| E036 | 1395.4500 | E074 | 1395.9250 | E112 | 1396.4000 |
| E037 | 1395.4625 | E075 | 1395.9375 | E113 | 1396.4125 |
| E038 | 1395.4750 | E076 | 1395.9500 | E114 | 1396.4250 |
| E039 | 1395.4875 | E077 | 1395.9625 | E115 | 1396.4375 |

| Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) |
|---------------|-----------------|---------------|-----------------|---------------|-----------------|
| E116 | 1396.4500 | E155 | 1396.9375 | E194 | 1397.4250 |
| E117 | 1396.4625 | E156 | 1396.9500 | E195 | 1397.4375 |
| E118 | 1396.4750 | E157 | 1396.9625 | E196 | 1397.4500 |
| E119 | 1396.4875 | E158 | 1396.9750 | E197 | 1397.4625 |
| E120 | 1396.5000 | E159 | 1396.9875 | E198 | 1397.4750 |
| E121 | 1396.5125 | E160 | 1397.0000 | E199 | 1397.4875 |
| E122 | 1396.5250 | E161 | 1397.0125 | E200 | 1397.5000 |
| E123 | 1396.5375 | E162 | 1397.0250 | E201 | 1397.5125 |
| E124 | 1396.5500 | E163 | 1397.0375 | E202 | 1397.5250 |
| E125 | 1396.5625 | E164 | 1397.0500 | E203 | 1397.5375 |
| E126 | 1396.5750 | E165 | 1397.0625 | E204 | 1397.5500 |
| E127 | 1396.5875 | E166 | 1397.0750 | E205 | 1397.5625 |
| E128 | 1396.6000 | E167 | 1397.0875 | E206 | 1397.5750 |
| E129 | 1396.6125 | E168 | 1397.1000 | E207 | 1397.5875 |
| E130 | 1396.6250 | E169 | 1397.1125 | E208 | 1397.6000 |
| E131 | 1396.6375 | E170 | 1397.1250 | E209 | 1397.6125 |
| E132 | 1396.6500 | E171 | 1397.1375 | E210 | 1397.6250 |
| E133 | 1396.6625 | E172 | 1397.1500 | E211 | 1397.6375 |
| E134 | 1396.6750 | E173 | 1397.1625 | E212 | 1397.6500 |
| E135 | 1396.6875 | E174 | 1397.1750 | E213 | 1397.6625 |
| E136 | 1396.7000 | E175 | 1397.1875 | E214 | 1397.6750 |
| E137 | 1396.7125 | E176 | 1397.2000 | E215 | 1397.6875 |
| E138 | 1396.7250 | E177 | 1397.2125 | E216 | 1397.7000 |
| E139 | 1396.7375 | E178 | 1397.2250 | E217 | 1397.7125 |
| E140 | 1396.7500 | E179 | 1397.2375 | E218 | 1397.7250 |
| E141 | 1396.7625 | E180 | 1397.2500 | E219 | 1397.7375 |
| E142 | 1396.7750 | E181 | 1397.2625 | E220 | 1397.7500 |
| E143 | 1396.7875 | E182 | 1397.2750 | E221 | 1397.7625 |
| E144 | 1396.8000 | E183 | 1397.2875 | E222 | 1397.7750 |
| E145 | 1396.8125 | E184 | 1397.3000 | E223 | 1397.7875 |
| E146 | 1396.8250 | E185 | 1397.3125 | E224 | 1397.8000 |
| E147 | 1396.8375 | E186 | 1397.3250 | E225 | 1397.8125 |
| E148 | 1396.8500 | E187 | 1397.3375 | E226 | 1397.8250 |
| E149 | 1396.8625 | E188 | 1397.3500 | E227 | 1397.8375 |
| E150 | 1396.8750 | E189 | 1397.3625 | E228 | 1397.8500 |
| E151 | 1396.8875 | E190 | 1397.3750 | E229 | 1397.8625 |
| E152 | 1396.9000 | E191 | 1397.3875 | E230 | 1397.8750 |
| E153 | 1396.9125 | E192 | 1397.4000 | E231 | 1397.8875 |
| E154 | 1396.9250 | E193 | 1397.4125 | E232 | 1397.9000 |

| Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) |
|---------------|-----------------|---------------|-----------------|---------------|-----------------|
| E233 | 1397.9125 | E272 | 1398.4000 | E311 | 1398.8875 |
| E234 | 1397.9250 | E273 | 1398.4125 | E312 | 1398.9000 |
| E235 | 1397.9375 | E274 | 1398.4250 | E313 | 1398.9125 |
| E236 | 1397.9500 | E275 | 1398.4375 | E314 | 1398.9250 |
| E237 | 1397.9625 | E276 | 1398.4500 | E315 | 1398.9375 |
| E238 | 1397.9750 | E277 | 1398.4625 | E316 | 1398.9500 |
| E239 | 1397.9875 | E278 | 1398.4750 | E317 | 1398.9625 |
| E240 | 1398.0000 | E279 | 1398.4875 | E318 | 1398.9750 |
| E241 | 1398.0125 | E280 | 1398.5000 | E319 | 1398.9875 |
| E242 | 1398.0250 | E281 | 1398.5125 | E320 | 1399.0000 |
| E243 | 1398.0375 | E282 | 1398.5250 | E321 | 1399.0125 |
| E244 | 1398.0500 | E283 | 1398.5375 | E322 | 1399.0250 |
| E245 | 1398.0625 | E284 | 1398.5500 | E323 | 1399.0375 |
| E246 | 1398.0750 | E285 | 1398.5625 | E324 | 1399.0500 |
| E247 | 1398.0875 | E286 | 1398.5750 | E325 | 1399.0625 |
| E248 | 1398.1000 | E287 | 1398.5875 | E326 | 1399.0750 |
| E249 | 1398.1125 | E288 | 1398.6000 | E327 | 1399.0875 |
| E250 | 1398.1250 | E289 | 1398.6125 | E328 | 1399.1000 |
| E251 | 1398.1375 | E290 | 1398.6250 | E329 | 1399.1125 |
| E252 | 1398.1500 | E291 | 1398.6375 | E330 | 1399.1250 |
| E253 | 1398.1625 | E292 | 1398.6500 | E331 | 1399.1375 |
| E254 | 1398.1750 | E293 | 1398.6625 | E332 | 1399.1500 |
| E255 | 1398.1875 | E294 | 1398.6750 | E333 | 1399.1625 |
| E256 | 1398.2000 | E295 | 1398.6875 | E334 | 1399.1750 |
| E257 | 1398.2125 | E296 | 1398.7000 | E335 | 1399.1875 |
| E258 | 1398.2250 | E297 | 1398.7125 | E336 | 1399.2000 |
| E259 | 1398.2375 | E298 | 1398.7250 | E337 | 1399.2125 |
| E260 | 1398.2500 | E299 | 1398.7375 | E338 | 1399.2250 |
| E261 | 1398.2625 | E300 | 1398.7500 | E339 | 1399.2375 |
| E262 | 1398.2750 | E301 | 1398.7625 | E340 | 1399.2500 |
| E263 | 1398.2875 | E302 | 1398.7750 | E341 | 1399.2625 |
| E264 | 1398.3000 | E303 | 1398.7875 | E342 | 1399.2750 |
| E265 | 1398.3125 | E304 | 1398.8000 | E343 | 1399.2875 |
| E266 | 1398.3250 | E305 | 1398.8125 | E344 | 1399.3000 |
| E267 | 1398.3375 | E306 | 1398.8250 | E345 | 1399.3125 |
| E268 | 1398.3500 | E307 | 1398.8375 | E346 | 1399.3250 |
| E269 | 1398.3625 | E308 | 1398.8500 | E347 | 1399.3375 |
| E270 | 1398.3750 | E309 | 1398.8625 | E348 | 1399.3500 |
| E271 | 1398.3875 | E310 | 1398.8750 | E349 | 1399.3625 |

| Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) |
|---------------|-----------------|---------------|-----------------|
| E350 | 1399.3750 | E389 | 1399.8625 |
| E351 | 1399.3875 | E390 | 1399.8750 |
| E352 | 1399.4000 | E391 | 1399.8875 |
| E353 | 1399.4125 | E392 | 1399.9000 |
| E354 | 1399.4250 | E393 | 1399.9125 |
| E355 | 1399.4375 | E394 | 1399.9250 |
| E356 | 1399.4500 | E395 | 1399.9375 |
| E357 | 1399.4625 | E396 | 1399.9500 |
| E358 | 1399.4750 | E397 | 1399.9625 |
| E359 | 1399.4875 | E398 | 1399.9750 |
| E360 | 1399.5000 | | |
| E361 | 1399.5125 | | |
| E362 | 1399.5250 | | |
| E363 | 1399.5375 | | |
| E364 | 1399.5500 | | |
| E365 | 1399.5625 | | |
| E366 | 1399.5750 | | |
| E367 | 1399.5875 | | |
| E368 | 1399.6000 | | |
| E369 | 1399.6125 | | |
| E370 | 1399.6250 | | |
| E371 | 1399.6375 | | |
| E372 | 1399.6500 | | |
| E373 | 1399.6625 | | |
| E374 | 1399.6750 | | |
| E375 | 1399.6875 | | |
| E376 | 1399.7000 | | |
| E377 | 1399.7125 | | |
| E378 | 1399.7250 | | |
| E379 | 1399.7375 | | |
| E380 | 1399.7500 | | |
| E381 | 1399.7625 | | |
| E382 | 1399.7750 | | |
| E383 | 1399.7875 | | |
| E384 | 1399.8000 | | |
| E385 | 1399.8125 | | |
| E386 | 1399.8250 | | |
| E387 | 1399.8375 | | |
| E388 | 1399.8500 | | |

Channel: E502 to E898

| Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) |
|---------------|-----------------|---------------|-----------------|---------------|-----------------|
| E502 | 1427.0250 | E544 | 1427.5500 | E586 | 1428.0750 |
| E503 | 1427.0375 | E545 | 1427.5625 | E587 | 1428.0875 |
| E504 | 1427.0500 | E546 | 1427.5750 | E588 | 1428.1000 |
| E505 | 1427.0625 | E547 | 1427.5875 | E589 | 1428.1125 |
| E506 | 1427.0750 | E548 | 1427.6000 | E590 | 1428.1250 |
| E507 | 1427.0875 | E549 | 1427.6125 | E591 | 1428.1375 |
| E508 | 1427.1000 | E550 | 1427.6250 | E592 | 1428.1500 |
| E509 | 1427.1125 | E551 | 1427.6375 | E593 | 1428.1625 |
| E510 | 1427.1250 | E552 | 1427.6500 | E594 | 1428.1750 |
| E511 | 1427.1375 | E553 | 1427.6625 | E595 | 1428.1875 |
| E512 | 1427.1500 | E554 | 1427.6750 | E596 | 1428.2000 |
| E513 | 1427.1625 | E555 | 1427.6875 | E597 | 1428.2125 |
| E514 | 1427.1750 | E556 | 1427.7000 | E598 | 1428.2250 |
| E515 | 1427.1875 | E557 | 1427.7125 | E599 | 1428.2375 |
| E516 | 1427.2000 | E558 | 1427.7250 | E600 | 1428.2500 |
| E517 | 1427.2125 | E559 | 1427.7375 | E601 | 1428.2625 |
| E518 | 1427.2250 | E560 | 1427.7500 | E602 | 1428.2750 |
| E519 | 1427.2375 | E561 | 1427.7625 | E603 | 1428.2875 |
| E520 | 1427.2500 | E562 | 1427.7750 | E604 | 1428.3000 |
| E521 | 1427.2625 | E563 | 1427.7875 | E605 | 1428.3125 |
| E522 | 1427.2750 | E564 | 1427.8000 | E606 | 1428.3250 |
| E523 | 1427.2875 | E565 | 1427.8125 | E607 | 1428.3375 |
| E524 | 1427.3000 | E566 | 1427.8250 | E608 | 1428.3500 |
| E525 | 1427.3125 | E567 | 1427.8375 | E609 | 1428.3625 |
| E526 | 1427.3250 | E568 | 1427.8500 | E610 | 1428.3750 |
| E527 | 1427.3375 | E569 | 1427.8625 | E611 | 1428.3875 |
| E528 | 1427.3500 | E570 | 1427.8750 | E612 | 1428.4000 |
| E529 | 1427.3625 | E571 | 1427.8875 | E613 | 1428.4125 |
| E530 | 1427.3750 | E572 | 1427.9000 | E614 | 1428.4250 |
| E531 | 1427.3875 | E573 | 1427.9125 | E615 | 1428.4375 |
| E532 | 1427.4000 | E574 | 1427.9250 | E616 | 1428.4500 |
| E533 | 1427.4125 | E575 | 1427.9375 | E617 | 1428.4625 |
| E534 | 1427.4250 | E576 | 1427.9500 | E618 | 1428.4750 |
| E535 | 1427.4375 | E577 | 1427.9625 | E619 | 1428.4875 |
| E536 | 1427.4500 | E578 | 1427.9750 | E620 | 1428.5000 |
| E537 | 1427.4625 | E579 | 1427.9875 | E621 | 1428.5125 |
| E538 | 1427.4750 | E580 | 1428.0000 | E622 | 1428.5250 |
| E539 | 1427.4875 | E581 | 1428.0125 | E623 | 1428.5375 |
| E540 | 1427.5000 | E582 | 1428.0250 | E624 | 1428.5500 |
| E541 | 1427.5125 | E583 | 1428.0375 | E625 | 1428.5625 |
| E542 | 1427.5250 | E584 | 1428.0500 | E626 | 1428.5750 |
| E543 | 1427.5375 | E585 | 1428.0625 | E627 | 1428.5875 |

| Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) |
|---------------|-----------------|---------------|-----------------|---------------|-----------------|
| E628 | 1428.6000 | E671 | 1429.1375 | E714 | 1429.6750 |
| E629 | 1428.6125 | E672 | 1429.1500 | E715 | 1429.6875 |
| E630 | 1428.6250 | E673 | 1429.1625 | E716 | 1429.7000 |
| E631 | 1428.6375 | E674 | 1429.1750 | E717 | 1429.7125 |
| E632 | 1428.6500 | E675 | 1429.1875 | E718 | 1429.7250 |
| E633 | 1428.6625 | E676 | 1429.2000 | E719 | 1429.7375 |
| E634 | 1428.6750 | E677 | 1429.2125 | E720 | 1429.7500 |
| E635 | 1428.6875 | E678 | 1429.2250 | E721 | 1429.7625 |
| E636 | 1428.7000 | E679 | 1429.2375 | E722 | 1429.7750 |
| E637 | 1428.7125 | E680 | 1429.2500 | E723 | 1429.7875 |
| E638 | 1428.7250 | E681 | 1429.2625 | E724 | 1429.8000 |
| E639 | 1428.7375 | E682 | 1429.2750 | E725 | 1429.8125 |
| E640 | 1428.7500 | E683 | 1429.2875 | E726 | 1429.8250 |
| E641 | 1428.7625 | E684 | 1429.3000 | E727 | 1429.8375 |
| E642 | 1428.7750 | E685 | 1429.3125 | E728 | 1429.8500 |
| E643 | 1428.7875 | E686 | 1429.3250 | E729 | 1429.8625 |
| E644 | 1428.8000 | E687 | 1429.3375 | E730 | 1429.8750 |
| E645 | 1428.8125 | E688 | 1429.3500 | E731 | 1429.8875 |
| E646 | 1428.8250 | E689 | 1429.3625 | E732 | 1429.9000 |
| E647 | 1428.8375 | E690 | 1429.3750 | E733 | 1429.9125 |
| E648 | 1428.8500 | E691 | 1429.3875 | E734 | 1429.9250 |
| E649 | 1428.8625 | E692 | 1429.4000 | E735 | 1429.9375 |
| E650 | 1428.8750 | E693 | 1429.4125 | E736 | 1429.9500 |
| E651 | 1428.8875 | E694 | 1429.4250 | E737 | 1429.9625 |
| E652 | 1428.9000 | E695 | 1429.4375 | E738 | 1429.9750 |
| E653 | 1428.9125 | E696 | 1429.4500 | E739 | 1429.9875 |
| E654 | 1428.9250 | E697 | 1429.4625 | E740 | 1430.0000 |
| E655 | 1428.9375 | E698 | 1429.4750 | E741 | 1430.0125 |
| E656 | 1428.9500 | E699 | 1429.4875 | E742 | 1430.0250 |
| E657 | 1428.9625 | E700 | 1429.5000 | E743 | 1430.0375 |
| E658 | 1428.9750 | E701 | 1429.5125 | E744 | 1430.0500 |
| E659 | 1428.9875 | E702 | 1429.5250 | E745 | 1430.0625 |
| E660 | 1429.0000 | E703 | 1429.5375 | E746 | 1430.0750 |
| E661 | 1429.0125 | E704 | 1429.5500 | E747 | 1430.0875 |
| E662 | 1429.0250 | E705 | 1429.5625 | E748 | 1430.1000 |
| E663 | 1429.0375 | E706 | 1429.5750 | E749 | 1430.1125 |
| E664 | 1429.0500 | E707 | 1429.5875 | E750 | 1430.1250 |
| E665 | 1429.0625 | E708 | 1429.6000 | E751 | 1430.1375 |
| E666 | 1429.0750 | E709 | 1429.6125 | E752 | 1430.1500 |
| E667 | 1429.0875 | E710 | 1429.6250 | E753 | 1430.1625 |
| E668 | 1429.1000 | E711 | 1429.6375 | E754 | 1430.1750 |
| E669 | 1429.1125 | E712 | 1429.6500 | E755 | 1430.1875 |
| E670 | 1429.1250 | E713 | 1429.6625 | E756 | 1430.2000 |

| Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) | Channel (HEX) | Frequency (MHz) |
|---------------|-----------------|---------------|-----------------|---------------|-----------------|
| E757 | 1430.2125 | E800 | 1430.7500 | E843 | 1431.2875 |
| E758 | 1430.2250 | E801 | 1430.7625 | E844 | 1431.3000 |
| E759 | 1430.2375 | E802 | 1430.7750 | E845 | 1431.3125 |
| E760 | 1430.2500 | E803 | 1430.7875 | E846 | 1431.3250 |
| E761 | 1430.2625 | E804 | 1430.8000 | E847 | 1431.3375 |
| E762 | 1430.2750 | E805 | 1430.8125 | E848 | 1431.3500 |
| E763 | 1430.2875 | E806 | 1430.8250 | E849 | 1431.3625 |
| E764 | 1430.3000 | E807 | 1430.8375 | E850 | 1431.3750 |
| E765 | 1430.3125 | E808 | 1430.8500 | E851 | 1431.3875 |
| E766 | 1430.3250 | E809 | 1430.8625 | E852 | 1431.4000 |
| E767 | 1430.3375 | E810 | 1430.8750 | E853 | 1431.4125 |
| E768 | 1430.3500 | E811 | 1430.8875 | E854 | 1431.4250 |
| E769 | 1430.3625 | E812 | 1430.9000 | E855 | 1431.4375 |
| E770 | 1430.3750 | E813 | 1430.9125 | E856 | 1431.4500 |
| E771 | 1430.3875 | E814 | 1430.9250 | E857 | 1431.4625 |
| E772 | 1430.4000 | E815 | 1430.9375 | E858 | 1431.4750 |
| E773 | 1430.4125 | E816 | 1430.9500 | E859 | 1431.4875 |
| E774 | 1430.4250 | E817 | 1430.9625 | E860 | 1431.5000 |
| E775 | 1430.4375 | E818 | 1430.9750 | E861 | 1431.5125 |
| E776 | 1430.4500 | E819 | 1430.9875 | E862 | 1431.5250 |
| E777 | 1430.4625 | E820 | 1431.0000 | E863 | 1431.5375 |
| E778 | 1430.4750 | E821 | 1431.0125 | E864 | 1431.5500 |
| E779 | 1430.4875 | E822 | 1431.0250 | E865 | 1431.5625 |
| E780 | 1430.5000 | E823 | 1431.0375 | E866 | 1431.5750 |
| E781 | 1430.5125 | E824 | 1431.0500 | E867 | 1431.5875 |
| E782 | 1430.5250 | E825 | 1431.0625 | E868 | 1431.6000 |
| E783 | 1430.5375 | E826 | 1431.0750 | E869 | 1431.6125 |
| E784 | 1430.5500 | E827 | 1431.0875 | E870 | 1431.6250 |
| E785 | 1430.5625 | E828 | 1431.1000 | E871 | 1431.6375 |
| E786 | 1430.5750 | E829 | 1431.1125 | E872 | 1431.6500 |
| E787 | 1430.5875 | E830 | 1431.1250 | E873 | 1431.6625 |
| E788 | 1430.6000 | E831 | 1431.1375 | E874 | 1431.6750 |
| E789 | 1430.6125 | E832 | 1431.1500 | E875 | 1431.6875 |
| E790 | 1430.6250 | E833 | 1431.1625 | E876 | 1431.7000 |
| E791 | 1430.6375 | E834 | 1431.1750 | E877 | 1431.7125 |
| E792 | 1430.6500 | E835 | 1431.1875 | E878 | 1431.7250 |
| E793 | 1430.6625 | E836 | 1431.2000 | E879 | 1431.7375 |
| E794 | 1430.6750 | E837 | 1431.2125 | E880 | 1431.7500 |
| E795 | 1430.6875 | E838 | 1431.2250 | E881 | 1431.7625 |
| E796 | 1430.7000 | E839 | 1431.2375 | E882 | 1431.7750 |
| E797 | 1430.7125 | E840 | 1431.2500 | E883 | 1431.7875 |
| E798 | 1430.7250 | E841 | 1431.2625 | E884 | 1431.8000 |
| E799 | 1430.7375 | E842 | 1431.2750 | E885 | 1431.8125 |

| Channel (HEX) | Frequency (MHz) |
|--------------------------|----------------------------|
| E886 | 1431.8250 |
| E887 | 1431.8375 |
| E888 | 1431.8500 |
| E889 | 1431.8625 |
| E890 | 1431.8750 |
| E891 | 1431.8875 |
| E892 | 1431.9000 |
| E893 | 1431.9125 |
| E894 | 1431.9250 |
| E895 | 1431.9375 |
| E896 | 1431.9500 |
| E897 | 1431.9625 |
| E898 | 1431.9750 |

Maintenance Check Sheet

Hospital/Organization: _____

Service Personnel: _____

Instrument Name: Transmitter

Instrument Model: ZM-540PA, ZM-541PA

Instrument Serial Number: _____

Hardware Revision Number: _____

Software Revision Number: _____

| | | |
|--------------------------------------------------------|------|------|
| 1. External Check | Pass | Fail |
| 2. Transmitter Channel | Pass | Fail |
| 3. Transmitting/Receiving Signal | Pass | Fail |
| 4. Display | Pass | Fail |
| 5. Key Operation | Pass | Fail |
| 6. ECG Check | Pass | Fail |
| 7. Respiration Check | Pass | Fail |
| 8. SpO ₂ Check | Pass | Fail |
| 9. NIBP Check | Pass | Fail |
| 10. NIBP Cuff for Attaching Transmitter to Patient Arm | Pass | Fail |

Overall Judgement

- OK
- Can be used but needs maintenance
- Maintenance required. Cannot be used.



NIHON KOHDEN CORPORATION
 1-31-4 Nishiochiai, Shinjuku-ku Tokyo 161-8560, Japan
 Phone +81 3-5996-8041

NIHON KOHDEN ITALIA S.r.l.
 Via Fratelli Bronzetti 28, 24124 Bergamo, Italy
 Phone +39 035-219543
 Fax +39 035-232546

North and South America

NIHON KOHDEN AMERICA, INC.
 15353 Barranca Parkway, Irvine, CA 92618, U.S.A.
 Toll-free +1-800-325-0283
 Phone +1 949-580-1555
 Fax +1 949-580-1550

UK Responsible Person
NIHON KOHDEN UK LTD.
 Unit 3, Heyworth Business Park,
 Old Portsmouth Road, Peasmarsh,
 Guildford, Surrey, GU3 1AF, UK
 Phone +44 14-8333-1328

NIHON KOHDEN MEXICO S.A. DE C.V.
 Insurgentes Sur 730, Piso 9 Oriente, Col. Del Valle
 C.P. 03100, Delegacion Benito Juarez, Ciudad de Mexico
 Phone +52 55-8851-5550
 Fax +52 55-8851-5580

Asia

**SHANGHAI KOHDEN
 MEDICAL ELECTRONIC INSTRUMENT CORP.**
 No. 567 Huancheng Bei Road
 Shanghai Comprehensive Industrial Development Zone
 Fengxian District, Shanghai 201401, China
 Phone +86 21-5743-6998

NIHON KOHDEN DO BRASIL LTDA.
 Rua Diadema, 89, 1º andar, conjuntos 11 a 17, bairro Mauá
 no Município de São Caetano do Sul, Estado de São Paulo
 CEP 09580-670, Brasil
 Phone +55 11-3044-1700
 Fax +55 11-3044-0463

NIHON KOHDEN SINGAPORE PTE LTD
 1 Maritime Square, #10-34 HarbourFront Centre
 Singapore 099253
 Phone +65 6376-2210
 Fax +65 6376-2264

Europe

EC REP European Representative

NIHON KOHDEN EUROPE GmbH
 Raiffeisenstrasse 10, 61191 Rosbach, Germany
 Phone +49 6003-827-0
 Fax +49 6003-827-599

NIHON KOHDEN INDIA PVT. LTD.
 308, Tower A, Spazedge, Sector 47, Sohna Road
 Gurgaon-122 002 Haryana, India
 Toll-free +91 1800-103-8182
 Phone +91 124-493-1000
 Fax +91 124-493-1029

NIHON KOHDEN DEUTSCHLAND GmbH
 Raiffeisenstrasse 10, 61191 Rosbach, Germany
 Phone +49 6003-827-0
 Fax +49 6003-827-599

NIHON KOHDEN MIDDLE EAST FZE
 JAFZA One Tower A, 19th floor, Office No. 1912
 P.O. Box 261516, Jebel Ali Free Zone, Dubai, U.A.E.
 Phone +971 4-884-0080
 Fax +971 4-880-0122

NIHON KOHDEN FRANCE SARL
 Centre d’Affaires La Boursidière, Bât C – RDC,
 92357 Le Plessis Robinson, France
 Phone +33 1-49-08-05-50
 Fax +33 1-49-08-93-32

NIHON KOHDEN KOREA, INC.
 3F, Cheongok Bldg., 88, Dongmak-ro, Mapo-gu,
 Seoul, 04075, Republic of Korea
 Phone +82 2-3273-2310
 Fax +82 2-3273-2352

NIHON KOHDEN IBERICA S.L.
 Calle Toronga 23, Oficina 1 28043 Madrid, Spain
 Phone +34 917-161-080
 Fax +34 913-004-676

Contact information is accurate as of Sep 2022. Visit <https://www.nihonkohden.com/> for the latest information.

The model and serial number of your device are identified on the rear or bottom of the unit.

Write the model and serial number in the spaces provided below. Whenever you call your representative concerning this device, mention these two pieces of information for quick and accurate service.

Model _____

Serial Number _____

Your Representative

Note for users in the territory of the EEA and Switzerland:

Any serious incident that has occurred in relation to the device should be reported to the European Representative designated by the manufacturer and the Competent Authority of the Member State of the EEA and Switzerland in which the user and/or patient is established.

ZM-540PA, ZM-541PA_0614-902754K



NIHON KOHDEN CORPORATION
1-31-4 Nishiochiai, Shinjuku-ku, Tokyo 161-8560, Japan
Phone +81 3-5996-8041
<https://www.nihonkohden.com/>