


OrthoPulse® Gen 2.0



 Biolux Technology GmbH
Neubaugasse 31
3462 Absdorf, Austria




Patented: orthopulse.com/patents

Charging Case: IP 32

OrthoPulse device: IP 37 FCC ID: 2A6CA-OPI IC: 28421-OPI



Charger case status







-  Charging
-  Error
-  Charging completed

Caution/Attention:

Plug cable to standard USB receptacle 1A 5V

Status Light Guide

In hand

-  Ready for treatment
-  Low battery, not ready
-  Treatment in progress
-  Treatment completed
-  Treatment paused
-  Error, email

support@orthopulse.com

In charging case (plugged in)

-  Charging, ready
-  Charging, not ready
-  Fully charged, ready
-  Restarting
-  Ready for Bluetooth*
-  Syncing



Made in Canada

380-642 Rev. A2