

Features of TAIYINENG:

TAIYINENG is used with tail Ageing yoga Mask to provide energy and proper operation mode for the yoga mask. There are two modes of operation: manual operation and Bluetooth intelligent operation. The operation mode is divided into 5 kinds, and the intensity can be adjusted in 10 grades. I can use the yoga mask I use to deeply activate the skin, tighten it deeply, promote the absorption of active ingredients, make the skin glowing and achieve amazing beauty effect.

Short press to increase strength, long press for three seconds to change the mold test

Indicator light

Short press to reduce strength, long press for three seconds to change the mold test



Open the key

Product specifications:

Name	TAIYINENG
rated voltage	3V DC
Rated current	4mA
working	Low frequency pulse boot
battery	button battery
导入强度	十档设置
Structure design	磁吸+面膜
country	China

Matters needing attention:

- 1, please clean the face before use, and then fully soak the mask essence into the mask cloth, and stick to the face with the pattern outwardly.
2. Taiyi can be combined with the corresponding contact on the decorative side of the mask cloth.
3. During use, please do not use metal materials to contact taiyinan.
- 4, each use time should not exceed 25 minutes, not more than 2 times a day.
5. Please place it out of reach of children and avoid high temperature and dust environment.
6. During use, please stop immediately if you feel uncomfortable.

Use method: choose one of the two

一, manual use:

1. Apply the mask to the surface of your skin, with the patterned side facing out. 2. Turn on the snooze switch. 3. Put it on the corresponding contact on the patterned side of the mask cloth, so that it can be closely combined with the mask. 4. Press + button to adjust to the most comfortable intensity. If the adjustment strength is too high, press - button to adjust back. 5. To change the mold test, long press + or - to change another mode. 6. Please turn off the switch on the laptop when it will stop automatically after 20 minutes.

二. Bluetooth smart (connected to mobile phones) :

1. Apply the mask to the surface of your skin, with the patterned side facing out. 2. Turn on the snooze switch. 3. Put it on the corresponding contact on the patterned side of the mask cloth, so that it can be closely combined with the mask. 4. Turn on the bluetooth of your phone, scan the QR code with wechat, and enter the Mini program. 5. Adjust the intensity to the most comfortable gear. 6. The default mode is integrated mode. If you want to change the mold test, tap the button that needs the mold test. 7. The system will stop automatically after 20 minutes. Please turn off the switch on the laptop after stopping.

Bluetooth connection: Turn on Bluetooth and scan the QR code to enter.

Abnormal connection:

1. Open Settings -- Privacy -- Bluetooth -- turn on the wechat permission
2. Open Settings -- Privacy -- Location permission -- wechat -- Bluetooth permission
3. Enable Positioning -- Hold down positioning -- WLAN and Bluetooth scan -- Enable Bluetooth scan

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by the party responsible for compliance

could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception,

which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help. The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction