




## Owner's Manual

### BodySense Smart Scale Lite




For FAQs and more information, please visit:  
[www.eufylife.com](http://www.eufylife.com)

## What's In the Box



## Inserting the Batteries


- 1 Turn the scale over and place it on a flat surface.
- 2 Open the battery door and insert three AAA batteries. Make sure the positive and negative ends are facing the correct polarity directions marked in the battery compartment.  
 - The digits "00.00" will show on the LED display.
- 3 Close the battery door and wait until the digits "00" are shown on the LED display.



**Lo** will appear on the LED display when the batteries are low. Replace all three batteries at the same time. Do not mix new and used batteries.


## Selecting the Measuring Unit

Press the **Unit** button on the bottom of the scale to cycle through and select the measuring unit: Pound / Stone / Kilogram. The measuring unit can be selected from within the EufyLife app as well.



## Installing the EufyLife App and Pairing

- 1 Download the EufyLife app from the App Store (iOS devices) or Google Play (Android devices).
- 2 Enable Bluetooth on your smartphone or tablet.
- 3 Open the EufyLife app and follow the on-screen instructions to pair the scale with your device.  
 - Once paired, the app and scale will automatically connect the next time the app is opened.



	Status Icon	
	Flashing	Pairing
	Solid	Connected



- This scale is compatible with devices running iOS 8.0 and above or Android 4.3 and above.  
- Do not pair the scale from your smartphone's Bluetooth menu. Use the EufyLife app to pair the scale with your smartphone.

## Measuring with the EufyLife App

Before taking your first measurement, you will need to create and set up a user ID in the EufyLife app. It is possible to create up to 16 user IDs.

- 1 Place the scale on a hard and flat surface.
- 2 Open the EufyLife app and select Smart Scale Lite on the device homepage.
- 3 Step onto the scale and remain still.
  - Your weight and Body Mass Index (BMI) will be analyzed automatically and shown in the app.



To get accurate measurements, it is recommended to step onto the scale barefoot.





## Measuring without the EufyLife App

This scale can also be used as a conventional weight-reading scale. You can simply step on the scale to measure your current weight without opening the EufyLife app.

## Maintenance

Regularly clean the scale with a dry or slightly dampened soft cloth.



- Do not wash the scale with water or immerse it in water.  
- Never use detergents to clean the scale as this may cause it to malfunction or discolor.

## Troubleshooting

Error Message

Error Display	Description	Solution
OL	Weight overload. The device will power off.	Stop using this scale for measurement.
Lo	Low battery. The device will power off.	Replace all three batteries at the same time.

When Measuring

Issue	Cause	Solution
Abnormal measurement results: <ul style="list-style-type: none"> <li>• Too high</li> <li>• Too low</li> <li>• Large difference detected between two recent measurements</li> </ul>	<ul style="list-style-type: none"> <li>• Incorrect posture.</li> <li>• Scale is placed on soft carpeting or uneven surface.</li> </ul>	<ul style="list-style-type: none"> <li>• Step onto the scale barefoot and stand still.</li> <li>• Place the scale on a hard and flat surface.</li> </ul>
LED display does not turn on after stepping onto the scale.	<ul style="list-style-type: none"> <li>• Batteries not installed.</li> <li>• Worn batteries.</li> </ul>	<ul style="list-style-type: none"> <li>• Install three new batteries.</li> <li>• Replace all three batteries at the same time.</li> </ul>

## Specifications

Product Model	T9142
Power Source	4.5 V (3 × AAA batteries)
LED Display	3.1 × 1.3 in / 80 × 32 mm
Measurement Unit	Pound / Stone / Kilogram
Unit Increments	0.1 lbs / 0.1 kg
Measurement Range	11 lbs - 410 lbs / 0 st 11 lbs - 29 st 4 lbs / 5 kg - 186 kg
Operating Environment	Temperature: 10 - 40 °C / 50 - 104 °F Humidity: 20 - 90% RH
Storage Environment	Temperature: -10 - 50 °C / 14 - 122 °F Humidity: ≤ 90% RH
Operating Frequency Range	2402 - 2480 MHz
Radio-Frequency Power	5 dBm
Dimension	12 × 12 × 1 in / 300 × 300 × 25 mm
Net Weight	3.5 lbs / 1.6 kg

## Customer Service



15-Month Limited Warranty



Lifetime Technical Support



support@eufylife.com



+1 (800) 988 7973 (US) Mon-Fri 9:00 - 17:00 (PT)  
+44 (0) 1604 936200 (UK) Mon-Fri 6:00 - 11:00 (GMT)  
+49 (0) 69 9579 7960 (DE) Mon-Fri 6:00 - 11:00  
+81 03 4455 7823 (JP) 月-金 9:00-17:00  
+ 86 400 0550 036 (中国) Mon-Fri 9:00 - 17:30



@EufyOfficial



@EufyOfficial



@EufyOfficial



@EufyOfficial

## Important Safety Instructions

- Read all instructions before using the device. Keep the manual for future reference.
- Only use the device as described in this user manual. Any misuse may cause electric shocks, burns, start a fire or other unexpected hazards.
- Children should not play with the appliance. Cleaning and user maintenance should not be done by children without supervision.
- Be aware that this device is designed for individuals to measure their weight. Any information provided by this device is in no way meant to treat, cure, or prevent diseases or illnesses from happening.
- The device is equipped with a data transmission function that emits electromagnetic energy. Portable Radio Frequency (RF) communications equipment (including peripherals such as antenna cables and external antennas) should not be used within 30 cm (12 inches) of the device, otherwise, performance may be affected.
- Do not modify this device.
- Do not expose this device to extreme temperatures, direct sunlight, moisture, or corrosive environments.
- Do not mix and use different types of batteries together e.g. lithium-ion (Li-ion) and nickel-metal hydride (NiMH) batteries.
- Do not mix new and used batteries.
- To prevent slipping, do not step on the device after bathing or showering.
- This device should not be used when it is adjacent to or stacked on top of another device because it could cause the device to operate improperly. If such use is necessary, this device and the other device should be observed to verify that they are operating normally.
- This device is powered by three AAA batteries. Do not attempt to recharge non-rechargeable batteries as this can cause the batteries to leak or explode.
- This device can be powered by rechargeable batteries.
- Used batteries should be placed in a sealed plastic bag and disposed of safely according to local environmental regulations.
- If the device will not be used for an extended period of time, remove batteries before storing the device in a cool and dry place.

## Usage Tips:

To ensure measurement accuracy, follow the instructions below when you start your measurement:

- Place the scale on a flat, hard surface. Soft surfaces such as carpeting will affect the performance of the scale.
- Step on the platform with bare feet. Stand still until the measurement is complete.
- For the most consistent measurements, weigh yourself at the same time of day and location.

## Notice



This symbol means the product must not be discarded as household waste, and should be delivered to an appropriate collection facility for recycling. Proper disposal and recycling helps protect natural resources, human health and the environment. For more information on disposal and recycling of this product, contact your local municipality, disposal service, or the shop where you bought this product.



This symbol means "Refer to instruction manual or booklet".

## FCC Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Warning: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures: (1) Reorient or relocate the receiving antenna. (2) Increase the separation between the equipment and receiver. (3) Connect the equipment into an outlet on a circuit different from that to which the receiver is connected. (4) Consult the dealer or an experienced radio/TV technician for help.

## RF Warning Statement

The device has been evaluated to meet general RF exposure requirements. The device can be used in portable exposure conditions without restriction.

SM-A560-V01