



KRAKEN

MEDICATION MANAGEMENT SYSTEM

User Manual

iMPAK©

1350 Campus Parkway • Suite 200

Neptune, NJ 07753

Phone 732.751.3475 • Fax 732.481.8525

TABLE OF CONTENT

1. PRODUCT DESCRIPTION	3
2. INDICATIONS FOR USE	3
3. INTENDED USER	3
4. GENERAL SAFETY PRECAUTIONS	3
5. STORAGE AND HANDLING	3
6. WARNINGS	3
7. MAINTENANCE	4
8. PREPARE THE KRAKEN FOR USE	4
9. REPLACE THE BATTERY IN YOUR DAYPAK DEVICE	4
10. UNDERSTANDING YOUR KRAKEN APPLICATION.....	5
11. TROUBLESHOOTING.....	9
12. iMPAK HEALTH CONTACT INFORMATION	10

1. PRODUCT DESCRIPTION

The KRAKEN Medication Management System is a medication adherence device aimed at increasing the compliance of the individual user. The system consists of:

- a. Weekly Sorting Station
- b. Seven (7) daily medication “wells”
- c. DayPak device
- d. Mobile application

2. INDICATION FOR USE

The Kraken System is to be used to track the usage of daily medication to increase the users compliance and enable the user to track and report medication compliance to share with their physicians and family members.

3. INTENDED USER

This device is intended to be used by an individual taking multiple medications at various times on a regular basis.

4. GENERAL SAFETY PRECAUTIONS

- Users should always follow any current medical treatment and/or instructions from their medical provider.
- The device should not be used in conjunction with water, or in a wet environment.
- Do not sterilize this unit with an autoclave or glass sterilizer.
- Disperse any static electricity from your body before using the unit.
- Do not expose the unit to strong shocks or vibrations.
- Do not disassemble, repair, or modify the unit.
- Do not insert battery with polarity reversed.
- Do not use batteries of a type other than that specified for use with the device.
- Do not use for any purpose other than tracking and reporting medication compliance.
- Keep away from pets and children.

5. STORAGE AND HANDLING

Do not store the unit in:

- Locations exposed to direct sunlight.
- Locations subject to high temperatures and high humidity.
- Wet or damp locations where water may get on the unit.
- Dusty locations.
- Near fires or open flames.
- Locations exposed to strong vibration.
- Locations exposed to strong electromagnetic fields.

6. WARNINGS

- This device is not designed or intended for diagnosis or as treatment.
- Patients should always consult their physician if they notice changes in their health.
- This device should never be used as a basis for starting or modifying medical treatment.
- Do not attempt self-diagnosis or self-treatment based on the recorded results.

- Keep out of reach of infants, small children, or anyone incapable of using the device properly.


7. MAINTENANCE

No maintenance of this system is required, except:


- The battery should be replaced when necessary. See instructions below.
- Units should be wiped with a damp cloth when necessary.

8. PREPARE THE KRAKEN SYSTEM FOR USE

- Download KRAKEN application to an iPhone 4S/5, iPad Mini and/or iPad (4th generation). Device must be Bluetooth 4.0 enabled. Note: the iPhone/iPad must be purchased separately.
- Create an *iMpak Connect* account at www.impakconnect.com (if necessary) Information needed can be found in the “about” section of your application.




- Pairing the device-to pair your device simply type in your serial number and password. This information can be found on the white label inside the Daypak Unit



- Fill the sorting Station with 1 week (7 days) of medications

- Load a new daily “well” into the Daypak every day and take all medication daily as prescribed
- Every 1-3 days, log onto the KRAKEN app and pair device to upload medication compliance information.

Send information (based on enrolled program, if applicable) to corresponding medical provider and/or caregiver. This should be done via button in application for “Send data to PHR”. See “Understanding Your App.”




- Receive calls/texts regarding medication compliance.

9. REPLACE THE BATTERY IN YOUR KRAKEN DAYPAK

- Remove “wells” and all medication from Daypak device
- Open Daypak unit, turn battery door ¼ turn to open and reveal battery
- Carefully pry battery out to remove
- Replace with new Lithium Ion CR2032 coin cell battery making sure the positive (+) side is facing up.

10. UNDERSTANDING YOUR KRAKEN APPLICATION

HOME –



COMPLIANCE –

Use this button to change date range



Taken Meds taken at correct time
Missed Missed dosage
Improper Meds taken within an hour of correct time
No meds no medications


SORT -

Click individual day's box to assign medication and time


Click time icon to set entire week

Click here to make changes to medication list






TREND –



SETTING REMINDERS–

To set reminders simply click the “Set or Clear Reminders” button on the “pairing page” (see pairing device) Choose your desired timeslot for the reminder you wish to set and click it. This will bring up a secondary dialog box which will allow you to select the frequency, duration and the time you wish the reminder to occur.

OPTIONS



11. TROUBLESHOOTING

Problem	Solution
I do not see any data showing up on my smart device.	Make sure your Bluetooth is turned on and your battery is still good
I no longer hear a beep or see a light flash when I open a door.	Change your battery (cr2032)
The application won't let me send my data to my physician.	Data can only be sent over 4g or wifi. If you are in an area where signal service is low you will not be able to send your data. You will

	need to try again once you return to an area where the signal service is strong.
I have not received any alerts or calls regarding my program	Please make sure to sync your device with a smart device. (see pairing above)

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and*
- (2) This device must accept any interference received, including interference that may cause undesired operation.*

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna,*
- Increase the separation between the equipment and receiver.*
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.*
- Consult the dealer or an experienced radio/TV technician for help.*

12. IMPAK CONTACT INFORMATION

iMPak Health
 1350 Campus Parkway
 Suite 200
 Neptune, NJ 07753
 732-751-7547
 www.iMPakhealth.com

For support please email support@impakhealth.com

Limited Warranty

Your Kraken Medication Management Device is warranted to be free from manufacturing defects for a period of one year under normal use. This warranty excludes any devices used to read the data from the Kraken Medication Management Device such as Low Energy Bluetooth enabled readers, smart phones, other mobile phones, tablets, or data collection hubs. This warranty extends only to the original retail purchaser. This warranty is complete and exclusive. This warranty does not apply; (a) to consumable parts, such as batteries, unless damage has occurred due to a defect in materials and workmanship; (b) to cosmetic damage, including but not limited to scratches, dents, and broken plastic; (c) to damage caused by use with non iMPAK HEALTH, LLC products; (d) to damage caused by accident, abuse, misuse, liquid contact, fire, earthquake or other external causes; (e) to damage caused by operating the product outside the permitted or intended uses described by iMPAK HEALTH, LLC; (f) to damage caused by service (including upgrades and expansions) performed by anyone who is not a representative of iMPAK HEALTH, LLC; (g) to a product or part that has been modified to alter functionality or capability without the written permission of iMPAK Health, LLC; or (h) if any iMPAK HEALTH, LLC serial number has been removed or erased from the device. The warrantor expressly disclaims liability for incidental, special, or consequential damages of any nature. Any implied warranties arising by the operation of law shall be limited in duration to the term of this warranty. This warranty gives you specific legal rights, and you may have other rights which vary from state to state. As a condition to operation of your warranty, you must register your device within 10 days of purchase by filling out and mailing the enclosed warranty card.

Warranty Card

Please complete the following information to register your device for our limited warranty.

Name: _____ Address: _____

Device Serial Number: _____ Date Purchased/Received: _____

Is this device for your personal use? Yes No If no who will be the primary user? _____

How did you hear about our device?

- | | | |
|------------------------------------|--|--|
| <input type="checkbox"/> Physician | <input type="checkbox"/> Family Member | <input type="checkbox"/> Insurance Company |
| <input type="checkbox"/> Hospital | <input type="checkbox"/> Advertisement | <input type="checkbox"/> Other |

Please Detach and mail to iMPak Health, 1350 Campus Pkwy., Neptune NJ 07753