

Rainbow-Operation manual


Technical Specification:

Model	Rainbow
Bluetooth version	Bluetooth4.0
Battery life	Continues working 24 hours
Connection distance	10M
Battery capacity	150mA Lithium battery
Recharging time	2 hours
Recharging current	100mA
Working current	5mA
Transmit power	-8dbm
Standby time	1 month
Standby current	0.1mA
Deadline for the release of electricity.	3.65V

User instruction:

1 – Charging your device

- a- Using USB cable (included).
- b- Using power cable (not included).
- c- Remove the bottom of the protective plastic cover
- d- The light is red when recharging.
- e- The light is green when charging is done. completed, the indicator



2- Download and launch the Rainbow App

3- Check the Bluetooth connection

- a- Turn On Bluetooth in your smartphone's settings.
- b- Long press button on the device, to boot.
- c- Rainbow is connected to the smartphone when the light is switched on.

4- Ready to start brushing



(Holding steady and horizontally the toothbrush)

- a- Begin brushing in the App interface, Bluetooth automatically connect.
- b- Follow the instructions on your mobile App to calibrate your toothbrush, holding steady and horizontally the grip.

5- Start brushing

- a- Follow the toothbrush movements and the toothbrush location on the App to brush your teeth with Rainbow.

Federal Communications Commission (FCC) Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Warning: Changes or modifications made to this device not expressly approved by **INNOVATIVE PRECISION INSTRUMENTS LIMITED** may void the FCC authorization to operate this device.

Note: The manufacturer is not responsible for any radio or tv interference caused by unauthorized modifications to this equipment. Such modifications could void the user's authority to operate the equipment.