



QUICK START GUIDE

Smart Pillbox

STEP 1 Remove Battery Separator



Use a coin, unscrew counterclockwise the battery cover.



Pull the clear plastic battery separator.

* Please note once battery is enabled, it remains on at all times.

STEP 2 Download App

Search "Memo Health" on

Visit <http://tinylogics.com>



or




STEP 3 Activate Memo Box



Smart Pillbox >


Select "More > Smart Pillbox"



Add New Device

Click "Add New Device" and follow the prompts

STEP 4 Use Memo Box



Indicator light


- No No need to take meds
- Red Double dosing
- Green Time to take meds

Button

1. Press once, to check battery, if no 'beep' = battery dead.
2. Press once, either 'no light' = need not take meds or 'green light' = time to take meds.

* Please note it is not an on/off button

Compartment for separate meds or cartridges



FCC ID: 2AH3P-M1602

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

IC : 21387-M1602

This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

Le présent appareil est conforme aux CNR d'Industrie Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes : (1) l'appareil ne doit pas produire de brouillage, et (2) l'utilisateur de l'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement.

CAN ICES-3(B)/NMB-3(B)